

1/16




Fig. 1

2/16




Fig. 2

3/16




Fig. 3

4/16




Fig. 4

5/16




Fig. 5

6/16

**Screening Checklist for Traditional/Abbreviated Premarket
Notification [510(k)] Submissions**

based on

**Guidance for Industry and FDA Staff
Format for Traditional and Abbreviated 510(k)s**
<http://www.fda.gov/cdrh/ode/guidance/1567.html>

Title	Related Information	Present	Inadequate	N/A
MDUFMA Cover Sheet	Medical Device User Fee Cover Sheet www.fda.gov/oc/mdufma/cover-sheet.html	✓ Sec. 1		
CDRH Premarket Review Submission Cover Sheet	CDRH Premarket Review Submission Cover Sheet www.fda.gov/opacom/morechoices/fdaforms/FDA-3514.pdf	✓ Sec. 2		
510(k) Cover Letter	Appendix A of "Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s" updated November 17, 2005	✓ Sec. 3		
Indications for Use Statement	Device Advice "Content of a 510(k)" Section D www.fda.gov/cdrh/devadvice/314312.html#link_6	✓ Sec. 4		
510(k) Summary or 510(k) Statement	Device Advice "Content of a 510(k)" Section E www.fda.gov/cdrh/devadvice/314312.html#link_7	✓ Sec. 5		
Truthful and Accuracy Statement	Device Advice "Content of a 510(k)" Section G www.fda.gov/cdrh/devadvice/314312.html#link_9	✓ Sec. 6		
Class III Summary and Certification	Class III Summary and Certification Form www.fda.gov/cdrh/manual/stmnciii.html			✓
Financial Certification or Disclosure Statement	FORM FDA 3454, Certification: Financial Interests and Arrangements of Clinical Investigators www.fda.gov/opacom/morechoices/fdaforms/FDA-3454.pdf FORM FDA 3455, Disclosure: Financial Interests and Arrangements of Clinical Investigators www.fda.gov/opacom/morechoices/fdaforms/FDA-3455.pdf Financial Disclosure by Clinical Investigators www.fda.gov/oc/guidance/financialdis.html			✓
Declarations of Conformity and Summary Reports (Abbreviated 510(k)s)	Use of Standards in Substantial Equivalence Determinations www.fda.gov/cdrh/ode/guidance/1131.html FDA Standards program www.fda.gov/cdrh/stdsprog.html Declaration of conformity www.fda.gov/cdrh/devadvice/3145.html#link_9 Required Elements for Declaration of Conformity to Recognized Standard www.fda.gov/cdrh/ode/reqrecstand.html	✓ Sec. 7		
Executive Summary	See section 10 in Chapter II of "Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s" updated November 17, 2005.	✓ Sec. 8		

Fig. 6

7/16

Title	Related Information	Present	Inadequate	N/A
Device Description	See section 11 in Chapter II of "Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s" updated November 17, 2005 www.fda.gov/cdrh/k863.html	✓ Sec. 9		
Substantial Equivalence Discussion	Guidance on the CDRH Premarket Notification Review Program 6/30/86 (K86-3), www.fda.gov/cdrh/k863.html	✓ Sec. 9 Sec. 15		
Proposed Labeling	Device Advice "Content of a 510(k)" Section H www.fda.gov/cdrh/devadvice/314312.html#link_10	✓ Sec. 11		
Sterilization/Shelf Life	Updated 510(k) Sterility Review Guidance (K90-1) www.fda.gov/cdrh/ode/guidance/361.html For reuse of single use devices, see Guidance for Industry and FDA Staff -- Medical Device User Fee and Modernization Act of 2002 Validation Data in Premarket Notification Submissions (510(k)s) for Reprocessed Single-Use Medical Devices www.fda.gov/cdrh/ode/guidance/1216.html	✓ Sec. 12		
Biocompatibility	FDA Blue Book Memo, G95-1, Use of International Standard ISO-10993, "Biological Evaluation of Medical Devices Part 1: Evaluation and Testing" www.fda.gov/cdrh/g951.html	✓ Sec. 13		
Software	Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices www.fda.gov/cdrh/ode/software.html			✓
Electromagnetic Compatibility/Electrical Safety	CDRH Medical Device Electromagnetic Compatibility Program www.fda.gov/cdrh/emc See also IEC 60601-1-2 Medical Electrical Equipment -- Part 1: General Requirements for Safety; Electromagnetic Compatibility -- Requirements and Tests (Second Edition, 2001)			✓
Performance Testing -- Bench	See section 18 in Chapter II of "Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s" updated November 17, 2005	✓ Sec. 14		
Performance Testing -- Animal	See section 19 in Chapter II of "Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s" updated November 17, 2005			✓
Performance Testing -- Clinical	See section 20 in Chapter II of "Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s" updated November 17, 2005 Certification/Disclosure Forms: Financial Interests and Arrangements of Clinical Investigators www.fda.gov/opacom/morechoices/fdaforms/FDA-3454.pdf www.fda.gov/opacom/morechoices/fdaforms/FDA-3455.pdf			✓
Kit Certification	Device Advice http://www.fda.gov/cdrh/devadvice/314c.html			✓

Fig. 7

8/16

TABLE OF CONTENTS

SECTION 1: Medical Device User Fee Cover Sheet	3
SECTION 2: CDRH Premarket Review Submission Cover Sheet	4
SECTION 3: Cover Letter	5
SECTION 4: Indications for Use.....	7
SECTION 5: 510(k) Summary	8
SECTION 6: Truthful and Accurate Statement.....	10
SECTION 7: Class III Summary and Certification	11
SECTION 8: Financial Certification or Disclosure Statement.....	11
SECTION 9: Declaration of Conformity.....	12
SECTION 10: Device Description	13
10.1 Introduction	13
10.2 Implants.....	13
10.3 Abutments	14
10.4 Abutments for Casting	15
10.5 Healing Caps	15
10.6 Surface Treatment	15
10.6.1 Sand Blasting.....	15
10.6.2 Cleaning Process.....	15
SECTION 11: Substantial Equivalence Discussion	16
11.1 Predicate device Information.....	16
11.1.1 MIS - Implant Technologies Ltd. (K040807 (Seven) ; K080162 (Uno))	16
11.1.2 Implant Direct LLC (K030007)	16
11.2 Comparison Table	17
11.3 Discussion of similarities and differences.....	21
SECTION 12: Proposed Labeling	23
12.1 Instructions for Use	23
12.2 Device Labels.....	23
12.3 Risks Warnings and Precautions	23
Caution, consult accompanying documents.....	23
SECTION 13: Sterilization and Shelf Life	26
13.1 Sterilization	26
13.2 Packaging – Microbial Barrier.....	26
13.3 Shelf Life.....	27
SECTION 14: Biocompatibility.....	28
SECTION 15: Software	29
SECTION 16: Electromagnetic Compatibility and Electrical Safety	29
SECTION 17: Performance Testing – Bench.....	30
17.1 Risk Analysis	30
17.2 Static and Fatigue Testing	30
17.3 Validation of Impurity Removal	32
SECTION 18: Performance Testing – Animal.....	34
SECTION 19: Performance Testing – Clinical	34

Fig. 8

9/16

 **WARNING.** If an implant is placed under weight that exceeds its functional capacity, excessive bone loss and/or structural sagging of the implant may occur. Different physiological and anatomical conditions could have an effect on the effectiveness of implant therapy, in particular, the following: inadequate bone quality and quantity, poor or absent oral hygiene monitoring, presence of blood disorders, and non-controlled endocrine problems. The restoration of small diameter implants with angled abutments, in the molar areas, is not recommended. Incorrect manipulation of small components inside of the oral cavity implies a breathing and/or ingestion risk. Forcing the implant into the surgical site, other than into the one decided upon, with the instrument, might cause stripping of the screwdriver's hexagon, stripping of the implant hexagon, cold soldering of the interface, between the insertion instrument and the implant, or stripping of the implant site's walls, a situation that could compromise the effective initial fixation of the implant.

Fig. 9

10/16

TEST REPORT - - GAMMA RADIATION VALIDATION

TESTING SERVICES COMPANY

Report Date: 18/1/09

To - CUSTOMER

From - SERVICE PROVIDER

Microbiology Division

Starting Date - 19/11/08

Laboratory Completion Date - 15/12/08

SAMPLE - Dental Implants

Test ReportGamma Radiation Validation of 22.5 kGy according to
AAMI TIR 33:2005Milouda Report Number - 49259, 49261, 49262, 50850
Order # 1012/08

Authorization of Report

	Name	Title	Signature
Prepared by	Dr. Ronit Ben Avraham	Head of Biology Lab	
Lab technician			
Approved by	Dr. Ronit Ben Avraham	Head of Biology Lab	

Page 2 of 9

Remarks:

1. The laboratory operates in accordance with recognized standards of the International ISO/IEC 17025 in all tests where recognition has been granted.
2. The microbiological tests are recognized and published by the Israeli Board of Health.
3. The results relate to the sample tested only.
4. Laboratory results are to be used in their entirety and no part may be quoted or copied in other documents.
5. Sampling was provided by and is the sole responsibility of the customer.
6. The Target Laboratory Accreditation Authority is not responsible for the test results.

Fig. 10

11/16

We, **COMPANY X** located in **LOCATION**, manufacturers of Dental Implants and accessories, have classified our:

- **Dental Implants & healing Cups & Abutments as a Class IIIb**, in accordance with Rule 8, of Annex II of the Medical Device Directive 93/42EEC "All implantable devices and long term surgically invasive devices are in class IIIb".
- **Dental Implants accessories & surgical kits are Class I** in accordance to rule 5 "All invasive devices with respect to body orifices, other than surgically invasive devices and which are not intended for connection to an active medical device: — are in Class I if they are intended for transient use".
- **Accessories / ball attachments as Class IIa** - For those products that are being sale separately to support the patient and are being remain in the patient mouth for several weeks Can be more than 30 days), such as the healing cup are classified as Class IIa according to rule 5: All invasive devices with respect to body orifices are in Class IIIb if they are intended for long-term use, except if they are used in the oral cavity as far as the pharynx, in an ear canal up to the ear drum or in a nasal cavity and are not liable to be absorbed by the mucous membrane, in which case they are in Class IIa."

Full Name: _____

Position: _____

Date: _____

Signature: _____

Fig. 11

12/16

Application Number: APP-00212	
Client: Global Devices	
Country: Canada	Device: Heart Monitor
Risk Class: A	Application Date: 4/12/2011
Documents:	
1-Executive Summary	Download
2-Company Profile	Download
3-Device EMC	Download
4-Insurance papers	Download
5-Insurance papers-FR	 File not found
6-Financial info	Download

Fig. 12

13/16

Upload to Library:

Select Document: <one of the forms under some module>

Select Language: English

Upload File:

Fig. 13

14/16




Fig. 14

15/16




Fig. 15

16/16




Fig. 16