

IDENTIFICATION OF ATYPICAL ANTIBODIES IN HUMAN BLOOD AND
BLOOD PRODUCTS

TECHNICAL FIELD

5 Described herein are methods for identifying atypical antibodies in blood and blood products.

BACKGROUND

Natural human antibodies have been identified that cause
10 pyrogenic responses during USP pyrogen tests. Some human donors naturally produce these "atypical antibodies," perhaps resulting from exposure to rabbits, rodents, or paracitic insects preying upon such host animals (i.e., fleas). The atypical antibodies are anomalous and uncommon,
15 but are not harmful to humans. Atypical antibodies can cross-react with rabbit white blood cell antigens and cause a pyrogenic response during rabbit pyrogen tests. The pyrogenic response is, however, a "false positive" result because other methods, such as the limulus amoebocyte lysate
20 assays (LAL), showed that suspect plasma samples giving pyrogenic responses in the rabbit assay do not contain endotoxins. In addition, the results of in vitro pyrogen tests (aka monocyte activation tests) indicate the absence of non-endotoxin pyrogens. Accordingly, atypical antibodies
25 in human blood or plasma cause erroneous rabbit pyrogen test results and may result in the disposal of individual or pooled blood or plasma that falsely tests positive as "pyrogenic."

The methods described herein permit the identification of blood or plasma samples containing atypical antibodies that result in false positives in pyrogen assays. This method is advantageous because atypical antibody-containing samples
5 can be eliminated before pooling with other blood or plasma and contaminating the pool. Accordingly, the method reduces manufacturing cost by preventing the unnecessary contamination of blood or plasma pools with samples that appear "pyrogenic" owing to the presence of atypical
10 antibodies. High throughput testing methods described herein permit identification of suspect falsely positive samples. Such samples can be discarded prior to pooling with other samples and prevent tainting the pool with atypical antibodies.

15

SUMMARY

Described herein are methods for identifying atypical antibodies in blood and blood products.

20 One aspect described herein is method for identifying atypical reactive antibodies in a blood products manufacturing process, the method comprising: (a) obtaining a sample of blood or plasma; (b) testing the sample and a control using any one or more of cellular agglutination,
25 fluorescence microscopy, immunoprecipitation, immunodiffusion, immunofluorescence, ELISA, flow cytometry, FACS, or Western blotting; (c) comparing the sample and control testing results; (d) determining whether the sample contains reactive atypical antibodies; and; (e) interdicting

a blood or plasma unit that was the source of the sample if the sample contains reactive atypical antibodies.

In one aspect, there is provided a method for interdicting a blood or plasma unit in a blood products manufacturing process, the method comprising: (a) obtaining a sample of blood or plasma from a blood or plasma unit; (b) testing the sample and a control using any one or more of cellular agglutination, fluorescence microscopy, immunoprecipitation, immunodiffusion, immunofluorescence, ELISA, flow cytometry, FACS, or Western blotting; (c) comparing the sample and control testing results; (d) determining whether the sample contains atypical reactive antibodies that cross reacts with rabbit white blood cell antigens; and (e) interdicting the blood or plasma unit that was the source of the sample if the sample contains said atypical reactive antibodies.

Another aspect described herein is a method for identifying atypical reactive antibodies in a blood products manufacturing process, the method comprising: (a) obtaining a sample of blood or plasma; (b) testing the sample and a control using a cellular agglutination assay; (c) comparing the sample and control testing results; (d) determining whether the sample contains reactive atypical antibodies; and; (e) interdicting a blood or plasma unit that was the source of the sample if the sample contains reactive atypical antibodies.

Another aspect described herein is a method for identifying atypical reactive antibodies in a blood products manufacturing process, the method comprising: (a) obtaining a sample of blood or plasma; (b) testing the sample and a control using flow cytometry; (c) comparing the sample and control testing results; (d) determining whether the sample contains reactive atypical antibodies; and; (e) interdicting a blood or plasma

unit that was the source of the sample if the sample contains reactive atypical antibodies.

Another aspect described herein is a method for identifying atypical reactive antibodies in a blooda plasma products manufacturing process, the method comprising: (a) obtaining a sample of blood or plasma; (b) testing the sample and a control using Western blotting; (c) comparing the sample and control testing results; (d) determining whether the sample

contains reactive atypical antibodies; and; (e) interdicting a blood or plasma unit that was the source of the sample if the sample contains reactive atypical antibodies.

Another aspect described herein is a method for evaluating potential false positive pyrogen test results in a blood products manufacturing process, the method comprising: (a) obtaining a sample of blood or blood product; (b) testing the sample and a control using fluorescence microscopy, immunoprecipitation, immunodiffusion, immunofluorescence, ELISA, flow cytometry, FACS, or Western blotting; (c) testing the sample and a control using the Limulus Amebocyte Lysate (LAL) or similar assay; (d) comparing the sample and control testing results in both assays; (e) evaluating whether the sample gave a false positive pyrogen test result; and (f) interdicting a blood or plasma unit that was the source of the sample if the sample contains reactive atypical antibodies.

In another aspect, there is provided a method for evaluating potential false positive pyrogen test results in a blood products manufacturing process, the method comprising: (a) obtaining a sample of blood or blood product; (b) testing the sample and a control using fluorescence microscopy, immunoprecipitation, immunodiffusion, immunofluorescence, ELISA, flow cytometry, FACS, or Western blotting; (c) testing the sample and a control using the Limulus Amebocyte Lysate (LAL) or similar assay; (d) comparing the sample and control testing results in both assays; and (e) evaluating whether the sample gave a false positive pyrogen test result, wherein the sample containing atypical reactive antibodies that cross react with rabbit white blood cell antigens is indicative of false positive pyrogen test results.

BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1: Microscopy of Rabbit WBC Agglutination Assays. Donor X and control plasma were incubated with rabbit WBCs and then viewed using microscopy. Light (A, B) and phase contrast (C, D) photomicrographs of rabbit white blood cells incubated with Donor X plasma (A, C) or control plasma (B, D). Agglutination was observed with the Donor X plasma (A, C), but not with the plasma control (B, D).

Figure 2: Microscopy of Human WBC Agglutination Assays. Donor X and control plasma were incubated with human WBCs

and then viewed using microscopy. Light (A, B) and phase contrast (C, D) photomicrographs of human white blood cells incubated with Donor X plasma (A, C) or a control plasma (B, D). No agglutination was observed with the Donor X plasma (A, C) or the plasma control (B, D).

Figure 3: Phase Contrast and Fluorescence Microscopy of Donor X Rabbit WBC Agglutination. Donor X and control plasma were incubated with rabbit WBCs, reacted with fluorescein-labeled anti-human IgG, and then viewed using phase contrast or fluorescence microscopy. Panel (A) shows a phase contrast photomicrograph of rabbit WBCs incubated with Donor X plasma. Panel (B) shows a fluorescent photomicrograph of rabbit WBCs incubated with Donor X plasma and then treated with fluorescein-labeled anti-human IgG. Strong fluorescence was observed in the agglutinated cell clusters indicating human IgG from Donor X plasma was responsible for cellular agglutination.

Figure 4: Fluorescence Microscopy of Rabbit WBC Agglutination Assays. Donor X and control plasma were incubated with rabbit WBCs and then reacted with fluorescein-labeled anti-human IgG. Panel (A) shows positive fluorescence of rabbit WBCs. Panel (B) shows the results with control plasma.

Figure 5: Fluorescence Microscopy of Rabbit WBC Agglutination Assays. Donor X and control plasma were incubated with rabbit WBCs, reacted with fluorescein-labeled anti-human IgG, and then viewed using fluorescence

microscopy. Panel (A) shows positive fluorescence of rabbit WBCs incubated with fluorescein-labeled anti-human IgG and a 40× image of the agglutinated cells (bottom panel). These results indicate that Donor X's plasma contains IgGs reactive with rabbit WBC cell-surface antigens. Panel (B) shows the weak fluorescence observed with the control plasma and a phase contrast photomicrograph of the rabbit WBCs showing no agglutination (bottom panel).

10 **Figure 6:** Fluorescence Microscopy of Human WBC Agglutination Assays. Donor X and control plasma were incubated with human WBCs, reacted with fluorescein-labeled anti-human IgG, and then viewed using fluorescence microscopy. Neither (A, Donor X) nor (B, control plasma) have strong fluorescence or
15 agglutination (bottom panels) indicating that Donor X's plasma contains IgGs not reactive to human WBC cell-surface antigens.

Figure 7: Fluorescence Flow Cytometry Histograms. Rabbit
20 WBCs were incubated with Donor X or control plasma, reacted with fluorescein-labeled anti-human IgG, and then analyzed by flow cytometry. Histograms show the relative fluorescence observed. Panels A and B show two separate experiments. Donor X plasma gave the largest signal with a relative
25 fluorescence at least three times that of the control plasma. See Tables 4 and 5 for quantitative results.

Figure 8: Fluorescence Flow Cytometry Histograms. Human WBCs were incubated with Donor X or control plasma, reacted with
30 fluorescein-labeled anti-human IgG, and then analyzed by flow

cytometry. Histograms show the relative fluorescence observed. Panels A and B show two separate experiments. Donor X plasma was not significantly different from the control plasma. See **Tables 6 and 7** for quantitative results.

5

Figure 9: Electrophoresis and Western Blotting Analysis of Samples of Pyrogenic IVIG-C Containing Donor X Plasma Isolates. Lane 1: MW markers; Lane 2: Pyrogenic IGIV-C (Donor X); Lane 3: Non-pyrogenic IGIV-C; Lane 4: Rabbit
10 Serum; Lane 5: Fetal Bovine Serum; Lane 6: Horse Serum. (A) Instant Blue stained gel. (B) Western blot control (no primary antibody) probed with anti-human IgG alkaline phosphatase immunoconjugate (IgG-AP). (C) Western blot using pyrogenic IGIV-C containing Donor X plasma isolates as
15 primary antibody and anti-human IgG-AP. (D) Western blot using non-pyrogenic IGIV-C as primary antibody and anti-human IgG-AP.

Figure 10: Electrophoresis and Western Blotting Analysis of
20 Samples of Pyrogenic IVIG-C Containing Donor X Plasma Isolates. Lane 1: MW markers; Lane 2: Fetal Bovine Serum; Lane 3: Rabbit Serum; Lane 4: Rabbit Serum, diluted 1:5; Lane 5: Rabbit Serum, diluted 1:10; Lane 6: Rabbit Serum, diluted 1:50. (A) Western blot control (no primary
25 antibody) probed with anti-human IgG alkaline phosphatase immunoconjugate (IgG-AP). (B) Western blot using pyrogenic IGIV-C containing Donor X plasma isolates as primary antibody and anti-human IgG-AP. (C) Western blot using non-pyrogenic IGIV-C as primary antibody and anti-human IgG-AP.
30 (D) Instant Blue stained gel.

Figure 11: Donor X's Plasma Cross-reacts with Rat WBCs: Fluorescence Microscopy of Rat WBC Agglutination Assays. Rat WBCs were incubated with Donor X or control plasma, reacted with fluorescein-labeled anti-human IgG, and then viewed using fluorescence microscopy. Panel (A) shows positive fluorescence and agglutination of rat WBCs incubated with fluorescein-labeled anti-human IgG and a 40× image of the agglutinated cells (bottom panel). These results indicate that Donor X's plasma contains IgGs reactive with rat WBC cell-surface antigens. Panel (B) shows the weak fluorescence observed with the control plasma and a 40× image of the rat WBCs showing no agglutination (bottom panel).

15

DETAILED DESCRIPTION

Described herein are methods for identifying atypical antibodies in blood and blood products. An individual, referred to herein as "Donor X," donated plasma that was pooled with other units for the manufacture of biotherapeutic protein products. During processing, the pooled plasma was assayed for pyrogenicity using the USP pyrogenicity rabbit assay. Unexpectedly, the pooled plasma tested positive as pyrogenic in the USP rabbit test. Further assays traced the pyrogenic agent to Donor X's plasma. Limulus amoebocyte lysate (LAL) assays showed that Donor X plasma was not contaminated with bacterial endotoxins. Instead, the assays described herein demonstrated that Donor X's plasma contained atypical antibodies that were responsible for the pyrogenic response.

Specifically, light and fluorescence microscopic examinations showed that Donor X plasma agglutinated rabbit and rat WBCs but not human WBCs. The rabbit WBC cross-reactivity was specific to Donor X because plasma from Donor X's parents, siblings, and children did not react. Fluorescent flow cytometry experiments showed that Donor X plasma contained IgG antibodies reactive with rabbit WBC cell-surface antigens and Western blot experiments confirmed the reactivity of the IgGs with rabbit sera. Collectively, these results suggest that Donor X may have been exposed to rodents or insect vectors of rodents that could have induced humoral immunity cross-reactive with rabbit WBCs. Thus, plasma from some individuals can test positive in USP rabbit pyrogen tests, not because they are contaminated with bacteria, but because they contain atypical antibodies that are cross-reactive with rabbit WBC antigens.

EXAMPLES

Example 1

20 United States Pharmacopia (USP) Pyrogen Assays

The current United States Pharmacopia § 151 outlines the pyrogen assay. The test involves measuring the rise in temperature of rabbits following the intravenous injection of a test solution. This assay is designed to determine whether products can be tolerated by the test rabbit in a dose not to exceed 10 mL per kg injected intravenously within a period of not more than 10 minutes. Initially, three rabbits are injected. If any rabbit shows an individual temperature rise $\geq 0.5^{\circ}\text{C}$, the test is continued using five additional rabbits. If three or more of the

eight rabbits show individual rises in temperature of $\geq 0.5^{\circ}\text{C}$ and/or the sum of the eight individual temperature rises exceeds 3.3°C the material under examination is considered pyrogenic.

5

A sample of Donor X plasma or a pooled sample without any Donor X plasma was diluted 1:100 into 10 mL of sodium chloride solution (0.9% NaCl) and injected into the ear veins of three healthy mature rabbits. The rabbits' temperatures were measured rectally within 10 minutes of injection. Temperature data are shown in **Table 1**. The Donor X-free plasma sample did not induce a temperature increase in any of the rabbits. In contrast, when Donor X plasma was tested, temperature increases of between 1.1 - 1.2°C were measured. Since the total temperature increase for the 3 rabbits was 3.4°C , Donor X plasma was considered pyrogenic and there was no need to extend testing to another 5 rabbits.

20 The immunoglobulins in a Pool containing 0% or 10% Donor X plasma were captured using a protein A column and tested for pyrogenicity. The Donor X-free plasma sample did not induce an increase in temperature but samples containing Donor X plasma were highly pyrogenic. These results indicated that
25 the pyrogenic response in rabbits could be related to immunoglobulins in Donor X plasma.

Table 1: Donor X plasma causes pyrogenic responses in rabbits				
Sample	USP Pyrogen Max Temp Increase (°C)			Total ΔT
	Rabbit 1	Rabbit 2	Rabbit 3	
Plasma Donor X (1:100)	1.1	1.1	1.2	3.4
Pooled Plasma w/o Donor X (1:100)	0.0	0.0	0.0	0.0
Pooled Mab Eluate w/ 10% Donor X	2.0	2.4	1.6	6.0
Pooled Mab Eluate w/o Donor X	0.0	0.1	0.0	0.1

Example 2

A series of experiments were performed using Donor X plasma to better understand the nature of its pyrogenicity.

5

White Blood Cell Agglutination and Microscopy Experiments

Agglutination experiments were performed to evaluate interactions between Donor X plasma and rabbit or human white blood cells (WBC). WBCs were harvested from rabbit and human whole blood by density gradient centrifugation using Histopaque® (Sigma-Aldrich) and suspended in normal buffered saline supplemented with BSA. The rabbit and human WBCs were then incubated with Donor X and control plasma in a 96-well microplate. Following incubation and washing, fluorescent-labeled anti-human IgG was added, and microplates were incubated, washed, and examined microscopically. Each well was examined for agglutination

10

15

using visible light and phase contrast microscopy, and then viewed using fluorescence microscopy (results discussed in the subsequent section).

5 Significant agglutination was observed in test wells containing Donor X plasma and rabbit WBCs (**Figures 1A and 1C**) but not in test wells containing Donor X plasma and human WBCs (**Figures 2A and 2C**). No agglutination was observed in wells containing control plasma and rabbit WBCs
10 (**Figures 1B and 1D**) or control plasma and human WBCs (**Figures 2B and 2D**). These results were reproduced in numerous assays and indicate that Donor X IgG binds to rabbit WBCs , but not to human WBCs. See **Figure 2**.

15 During several of the agglutination experiments, cytotoxicity was observed in samples containing rabbit WBCs and Donor X plasma but not in wells containing control plasma and rabbit WBCs, nor any wells containing human WBCs with Donor X or control plasma. The observation that Donor
20 X plasma is toxic to rabbit WBCs suggested specific binding of Donor X immunoglobulins to these cells. See **Figures 1-3**.

Fluorescent Microscopy Experiments

Fluorescent microscopy experiments were performed in
25 parallel with the agglutination and light microscopy studies described above and the results are presented in **Table 2**. Rabbit WBCs were strongly fluorescent in samples incubated with Donor X plasma (**Figures 4A and 5A**), compared to a relatively weak degree of fluorescence for samples incubated
30 with control plasma (**Figure 4B**) or with only fluorescent-

labeled anti-human IgG. See **Figures 4 and 5**. These findings indicated that the fluorescence observed for rabbit WBCs was specific to IgG present in Donor X. A weak degree of fluorescence was observed for human WBCs incubated with both Donor X and control plasma (**Figure 6**) or with only the fluorescent-labeled anti-human IgG (no Donor X or control plasma added). **Figure 6B**. These findings indicated that the fluorescence observed with human WBCs represented non-specific binding, independent of the presence of Donor X or control IgG.

Table 2: Fluorescent Microscopy Results

Plasma Source	Host WBC	Fluorescent-labeled anti-human IgG	Fluorescence Score (0-4)
Donor X	Rabbit	Yes	3-4
Control	Rabbit	Yes	0-2
None	Rabbit	Yes	0-1
Donor X	Human	Yes	0-1
Control	Human	Yes	0-1
None	Human	yes	0-1

Example 3**Flow Cytometry Experiments**

In order to quantitate antibody binding and the fluorescence observed by microscopy, flow cytometry studies were performed. In these experiments, rabbit and human WBCs were incubated with Donor X and control plasma and washed before adding and incubating with fluorescent-labeled anti-human IgG. The cell samples were washed, resuspended in normal

buffered saline to a concentration ranging from approximately 3×10^6 to 5×10^6 cells/mL and analyzed by flow cytometry.

5 **Figure 7** contains histogram overlays showing relative fluorescent intensity of two different rabbit WBC samples incubated with Donor X or control plasma. Two additional samples were included as assay controls, unstained rabbit WBCs (cell control) and rabbit WBCs treated with
10 fluorescent-labeled anti-human-IgG only. The median fluorescence of the rabbit WBC sample incubated with Donor X plasma was 3264 in experiment 1 and 922 in experiment 2 and, significantly higher than that observed after incubating rabbit WBCs with control plasma (median fluorescence of 499
15 and 175 for experiments 1 and 2, respectively). See **Tables 3 and 4**. Thus, even though there was some overlap between the rabbit WBCs incubated with Donor X and control plasma, there was a distinct difference in fluorescent intensity. These results correlated well with the microscopy study
20 results, discussed above, where WBCs reacted with Donor X plasma and produced significantly stronger fluorescence than those incubated with control plasma.

Table 3: Flow cytometry histogram data showing relative fluorescence of rabbit WBCs incubated with Donor X or control plasma					
Experiment 1	Relative Fluorescence (N = 20,0000)				
Cell Type	Mean	Median	Mode	Mode Count	Std Dev
Donor X Plasma	4,507.4	3,263.8	10,000.1	1,632	3,688.4
Control Plasma	1,508.9	498.9	9,003.1	509	2,404.3
Anti-IgG only	13.4	5.7	3.7	400	176.3
Cell Control	4.3	3.0	3.4	630	49.7

Table 4: Flow cytometry histogram data showing relative fluorescence of rabbit WBCs incubated with Donor X or control plasma		
Experiment 2	Relative Fluorescence (N = 20,0000)	
Cell Type	Mean	Median
Donor X Plasma	1752	922
Control Plasma	501	175
Anti-IgG only	13	1
Cell Control	1	1

The flow cytometry experiments were repeated with human WBCs and the results are shown in Figure 8 and Tables 5 and 6. Human WBCs incubated with Donor X plasma and control plasma had similar histograms, indicating no significant difference in fluorescent intensity. These results correlated with

microscopy study results, discussed previously, in which human WBCs treated with Donor X and control plasma produced a similar degree of fluorescence. The human WBC sample stained with only fluorochrome-labeled anti-human-IgG also showed significant overlap with Donor X and control plasma samples, and thus indicated a significant degree of non-specific binding by the secondary antibody.

Table 5: Flow cytometry histogram data showing relative fluorescence of human WBCs incubated with Donor X or control plasma					
Experiment 3	Relative Fluorescence (N = 20,0000)				
Cell Type	Mean	Median	Mode	Mode Count	Std Dev
Donor X Plasma	1004.9	47.7	16.1	304	1739.0
Control Plasma	820.9	37.0	11.3	294	1472.1
Anti-IgG only	381.8	17.4	9.4	332	785.5
Cell Control	5.6	2.8	2.9	695	2663.1

Table 6: Flow cytometry histogram data showing relative fluorescence of human WBCs incubated with Donor X or control plasma

Experiment 4	Relative Fluorescence (N = 20,0000)	
Cell Type	Mean	Median
Donor X Plasma	174	7
Control Plasma	274	16
Anti-IgG only	76	4
Cell Control	1	1

In summary, the flow cytometry analyses showed significant binding of Donor X immunoglobulins (i.e., IgGs) to rabbit WBCs as compared to control plasma and minimal binding to human WBCs.

Example 4

Supplemental Rabbit Pyrogen Testing

In order to evaluate a possible genetic association for Donor X immunoglobulins and their effect on rabbit pyrogenicity, USP pyrogen assays were performed on sera donated by relatives of Donor X. A sample of Donor X serum was also tested as a control. Because previous studies demonstrated that Donor X plasma produced a significant pyrogen response at dilutions of 1:100, all test samples were diluted 1:100 in sterile normal saline (0.9% NaCl, USP, for injection) prior to rabbit pyrogen testing. An aliquot of each sample was also used in a LAL assay to examine endotoxin contamination as a source of pyrogen response. Pyrogen and LAL results are presented in **Table 7**.

Table 7: USP Pyrogen Test Results for Donor X and Immediate Family Members

	USP Rabbit Pyrogen Results (max. temperature increase in °C)			
Sera Samples (1:100 dilution)	Rabbit 1	Rabbit 2	Rabbit 3	Total ΔT
Donor X	1.0	0.3	0.6	1.9
Father	0.0	0.1	0.0	0.1
Mother	0.0	0.2	0.0	0.2
Sister 1	0.0	0.0	0.0	0.0
Brother	0.0	0.0	0.0	0.0
Sister 2	0.0	0.1	0.0	0.1
Daughter 1	0.0	0.1	0.1	0.2
Daughter 2	0.0	0.2	0.0	0.2
Son	0.0	0.1	0.0	0.1

Serum from Donor X produced a significant temperature increase in two of the three test rabbits, with a total temperature increase of 1.9 °C. This response was consistent with previous testing with Donor X plasma. Serum from Donor X's relatives, including parents, siblings, and children produced no significant temperature increases. Limulus Amoebocyte Lysate assay (LAL) results for all samples were negative, indicating that exogenous endotoxin did not contribute to the rabbit pyrogenic responses.

Example 5**Red Blood Cell Agglutination Studies**

Donor X plasma was tested with rabbit red blood cells in a series of agglutination experiments to determine if Donor X plasma contains immunoglobulins specific for antigens on rabbit RBCs. Incompatibility between Donor X plasma immunoglobulins and rabbit RBCs could potentially cause hemolysis and pyrogenicity. For these studies, Donor X and control plasma were titrated against a suspension of rabbit RBCs. The suspension was observed at three time points: (1) immediately; (2) after a 30-minute incubation at 37 °C; and (3) after anti-human-globulin serum was added.

Both Donor X and control plasma produced strong agglutination of rabbit RBCs at all time points, and an equivalent titer was observed for Donor X and the positive control. Hemolysis was observed at low dilutions of both Donor X and the positive control.

The presence of anti-A, and/or anti-B immunoglobulins present in Donor X and control plasma could potentially cross-react with rabbit RBC antigens with similar epitopes to human A and B antigens. Accordingly, Donor X and control plasma were pre-absorbed with human A and/or B RBCs to remove cross-reacting anti-A and anti-B antibodies. The pre-absorbed plasma was then tested against rabbit RBCs as described above. Both Donor X and control plasma produced strong agglutination, similar to the initial results. No difference in reactivity was observed between Donor X and the control plasma. These results showed the presence of

antibodies in Donor X plasma with broad cross-reactivity to antigens/epitopes on rabbit RBCs. In addition, these results suggest that an RBC-mediated process is not responsible for the pyrogenic response in rabbits.

5

Donor X plasma was also tested for antibodies to human RBC antigens using a RBC antibody-identification panel. Negative results were obtained with all panel cells, confirming that Donor X plasma contains no clinically
10 significant alloantibodies.

Antigen phenotyping was also performed on Donor X RBCs, included typing for RBC antigens belonging to the Rh, Kell, Duffy, Kidd, Lewis, MNS, P, and Lutheran blood group
15 systems. Donor X RBCs were of a common RBC phenotype, and there were no unusual results.

Example 6

Downstream plasma products containing Donor X plasma were
20 assayed to identify the factor responsible for generating the pyrogenic-response. Human Immunoglobulin Globulin, Intravenous containing 10% caprylate/chromatography purified (e.g., IGIV-C 10%, i.e., Gamunex®, Grifols Therapeutics Inc., formerly Talecris Biotherapeutics, Inc.) that was
25 produced from plasma pools containing Donor X plasma was assayed using Western blotting.

Western Blotting

Samples of "pyrogenic" IGIV-C produced from Donor X-
30 containing plasma pools, non-pyrogenic IGIV-C produced from

Donor X-free plasma pools, rabbit serum, fetal bovine serum, and horse serum were run on four 4-20% reducing SDS-PAGE gels. One gel was stained with Instant Blue (**Figure 9A**) while the other three were transferred to PVDF membranes (**Figures 9B-D**). One membrane was reacted with only anti-Human IgG conjugated to alkaline phosphatase (**Figure 9B**). The remaining membranes were reacted with pyrogenic IGIV-C or with non-pyrogenic IGIV-C, and then with the anti-Human IgG alkaline phosphatase conjugate (**Figures 9C-D**). Lane 1: MW markers; Lane 2: Pyrogenic IGIV-C (Donor X); Lane 3: Non-pyrogenic IGIV-C; Lane 4: Rabbit Serum; Lane 5: Fetal Bovine Serum; Lane 6: Horse Serum.

The Instant Blue-stained gel showed that comparable amounts of rabbit serum, fetal bovine serum and horse serum were loaded on the gel. Lanes 4-6 in **Figure 9A**. The membrane stained with only anti-Human IgG showed that the secondary antibody (anti-human IgG) was specific and reacted only with human IgG (pyrogenic and non-pyrogenic IGIV-C). Lanes 2 and 3 in **Figure 9B**. The membrane reacted with pyrogenic IGIV-C and the anti-Human IgG alkaline phosphatase conjugate showed that pyrogenic IGIV-C reacted strongly to rabbit serum and weakly with fetal bovine and horse serum. Lanes 4-6 in **Figure 9C**. The membrane reacted with non-pyrogenic IGIV-C and the anti-Human IgG alkaline phosphatase conjugate showed that non-pyrogenic IGIV-C reacted weakly with the three test serums, including rabbit serum. Lanes 4-6 in **Figure 9D**. Collectively, these results indicate that pyrogenic IGIV-C reacts strongly with rabbit sera, while non-pyrogenic IGIV-C does not.

Samples of fetal bovine serum and various concentrations of rabbit serum were run on four 4-20% SDS-PAGE gels. Three gels were transferred to PVDF membranes (**Figures 10A-C**), and one gel was stained with Instant Blue (**Figure 10D**). One membrane was reacted with only anti-Human IgG conjugated to alkaline phosphatase (**Figure 10A**). The remaining membranes were reacted with pyrogenic IGIV-C or with non-pyrogenic IGIV-C, and then with the anti-Human IgG alkaline phosphatase immunoconjugate (**Figures 10B-C**). Lane 1: MW markers; Lane 2: Fetal Bovine Serum; Lane 3: Rabbit Serum; Lane 4: Rabbit Serum, diluted 1:5; Lane 5: Rabbit Serum, diluted 1:10; Lane 6: Rabbit Serum, diluted 1:50.

The Western blot was negative when the membrane was probed with only anti-Human IgG. **Figure 10A**. However, when pyrogenic IGIV-C was used as the primary Ab, it specifically reacted to undiluted rabbit serum and rabbit serum diluted 1:5 and 1:10. See Lanes 3-5 in **Figure 10B**. The 1:50 dilution of rabbit serum was not reactive (Lane 6). Only undiluted rabbit serum was detected when non-pyrogenic IGIV-C was used as the primary antibody. Lane 4 in **Figure 10C**. The Instant Blue-stained gel showed the relative amounts of fetal bovine serum and rabbit serum that were loaded onto the gels. Lanes 2-3 in **Figure 10D**. Overall, these results indicate that pyrogenic IGIV-C contains ~10-fold more antibodies against rabbit serum than non-pyrogenic IGIV-C.

Example 7**Rat WBC Fluorescence Microscopy and Agglutination**

Rat WBCs were isolated from whole blood by density gradient centrifugation using Histopaque®. Rat WBCs were reacted
5 with Donor X plasma as described for Example 2 above. Donor X plasma produced distinct agglutination and fluorescence with Rat WBCs. See Figure 11.

Example 8**10 Summary of Results**

The results of the microscopy (**Figures 1-6**) and flow cytometry experiments (**Figures 7-8**) indicated that Donor X IgG readily bound to rabbit WBCs. Antibody-binding was demonstrated by both agglutination and fluorescent
15 intensity. Although control plasma produced some fluorescence with rabbit WBCs, the degree of fluorescence was significantly less than that observed with Donor X. Control plasma did not cause agglutination of rabbit WBCs in any of the experiments, while Donor X plasma consistently
20 produced agglutination. Donor X plasma (IgG) did not produce agglutination with human WBCs in any of the experiments, indicating that Donor X IgG did not bind to human WBCs.

25 The binding of Donor X IgG to rabbit WBCs is a likely trigger for rabbit WBC activation and release of endogenous (leukocytic) pyrogens, which caused the observed fever response.

Rabbit pyrogen assays conducted on immediate family of Donor X (i.e., parents, siblings, and children) were uniformly negative. **Table 7.** These results indicated that the unique properties of Donor X IgG, with regard to rabbit temperature response, were not dominant allele-based, but rather antibody-specific to Donor X.

Experiments with rabbit RBCs demonstrated that both Donor X and control plasma contained antibodies with broad cross-reactivity to an antigen on rabbit RBCs. Donor X and control plasma produced very similar reactions with rabbit RBCs, suggesting that an RBC-mediated process is not responsible for the pyrogenic response in rabbits.

Western Blotting experiments showed that that pyrogenic IGIV-C (containing Donor X plasma isolates) reacts strongly with rabbit sera, while non-pyrogenic IGIV-C does not. This indicates that the presence of atypical IgGs from Donor X in the pyrogenic IGIV-C, were responsible for eliciting a pyrogenic response in the USP pyrogen test. In addition, the Western Blotting experiments showed that pyrogenic IGIV-C contains ~10-fold more antibodies against rabbit serum than non-pyrogenic IGIV-C. This indicated a probable stochastic effect from the Donor X atypical antibodies that caused pyrogenic responses.

Experiments with rat white blood cells show that Donor X plasma is capable of cross reacting with and causing agglutination of rat WBCs. These results suggested that direct exposure to rats or indirect exposure by insect

vectors of rodents (e.g., fleas), might have lead to "atypical" IgG-immunoglobulin immunity with cross-reactivity to both rat and rabbit cells.

5 **Example 9**

High-Throughput Assays

High-throughput ELISA, fluorescence, or Western blot experiments are performed by incubating test samples in 96-well, 192-well, or 384-well plates or membranes, washing, 10 blocking, and probing the samples using enzyme- or fluorophore-immunoconjugates and then analyzing the results via fluorometry, luminometry, densitometry, colorimetry, or UV/visible absorbance, among other detection methods. Such high-throughput assays permit in-line analysis of blood or 15 plasma samples or products before, during, and after processing and can eliminate reactive samples, such as those containing atypical immunoglobulins, which may produce false-positive pyrogen results in assays.

CLAIMS

1. A method for interdicting a blood or plasma unit in a blood products manufacturing process, the method comprising:

- (a) obtaining a sample of blood or plasma from a blood or plasma unit;
- (b) testing the sample and a control using any one or more of cellular agglutination, fluorescence microscopy, immunoprecipitation, immunodiffusion, immunofluorescence, ELISA, flow cytometry, FACS, or Western blotting;
- (c) comparing the sample and control testing results;
- (d) determining whether the sample contains atypical reactive antibodies that cross react with rabbit white blood cell antigens; and
- (e) interdicting the blood or plasma unit that was the source of the sample if the sample contains said atypical reactive antibodies.

2. The method of claim 1, wherein the testing step uses a cellular agglutination assay.

3. The method of claim 1, wherein the testing step uses flow cytometry.

4. The method of claim 1, wherein the testing step uses Western blotting.

5. A method for evaluating potential false positive pyrogen test results in a blood products manufacturing process, the method comprising:

- (a) obtaining a sample of blood or blood product;
- (b) testing the sample and a control using fluorescence microscopy, immunoprecipitation, immunodiffusion, immunofluorescence, ELISA, flow cytometry, FACS, or Western blotting;
- (c) testing the sample and a control using the Limulus Amebocyte Lysate (LAL) or similar assay;
- (d) comparing the sample and control testing results in both assays; and
- (e) evaluating whether the sample gave a false positive pyrogen test result, wherein the sample containing atypical reactive antibodies that cross react with rabbit white blood cell antigens is indicative of false positive pyrogen test results.