




Figure 1A

Continue B

e-Dispenser starts dispensation process once payment is successful.

For each product that has been dispensed, a snapshot of the product image will be taken by the system for Pharmacist / Doctor to validate that the dispensed medication is correct

The dispensed medication will then be dropped into the collection box.

Status of dispensation process will be shown on screen until all the products have been dispensed and dropped into the collection box

At this stage, the e-Dispenser system will take a final snapshot of all the products that have been dropped into the collection box

It is mandatory at this stage that the Pharmacist / Doctor has to perform a final validation of the dispensed products based on the final image taken.

Once this final validation has been done by the Pharmacist / Doctor, the collection box will open and instruction will be shown on the screen to ask the customer // patient to collect the medication(s).

After medication (s) has been collected by the customer (patient), Pharmacist / Doctor will request them to show the taken medication(s) in front of the camera to ensure that all medications have been collected.

The collection box door will close by itself automatically after 30 seconds. Any uncollected goods will be dropped into the uncollected bin before the next session starts.

Figure 1B