

Fig. 1B

NME7 (~42KDa)







Fig. 1A

NME1 (~17KDa)



Western Blot Analysis of Cancer and Stem Cells, Probing for Presence of NME1, NME6 and NME7 in the Cell Lysates





Reducing western blot

- 1- lyse cells with RIPA buffer
- 2- run reducing SDS-PAGE (20ul per lane)
- 3- transfer protein to PVDF membrane
- 4- probe membrane with anti NM23-H7 (B-9, Santa Cruz Biotechnology)
- 5- use goat anti mouse-HRP as secondary
- 6- detect protein by chemiluminescence

Reducing western blot

- 1- run reducing SDS-PAGE (20ul of CM per lane)
- 2- transfer protein to PVDF membrane
- 3- probe membrane with anti NM23-H7 (B-9, Santa Cruz Biotechnology)
- 4- use goat anti mouse-HRP as secondary
- 5- detect protein by chemiluminescence

Reducing western blot

- 1- lyse cells with RIPA buffer
- 2- run reducing SDS-PAGE (20ul per lane)
- 3- transfer protein to PVDF membrane
- 4- probe membrane with anti NM23-H7 (B-9, Santa Cruz Biotechnology)
- 5- use goat anti mouse-HRP as secondary
- 6- detect protein by chemiluminescence

Recombinant NME7 novel variant containing NDPK domains A and B, “NME7-AB”, expresses well with high yield in *E. coli* and as the soluble protein

Fig. 4A

FPLC purification of NME7-AB




Fig. 4B

SDS-PAGE of NME7-AB





Fig. 4C
NME7-AB is purified



ELISA shows NME7 Dimerizes MUC1*

MUC1* extra cellular domain peptide immobilized on plate was bound by NME7 to saturation; a second MUC1* peptide with a C-terminal His-tag or Biotin tag was added and visualized by HRP labeled antibody to either His-tag or HRP labeled streptavidin

Fig. 5A










Fig. 5B





T47D + anti-NM23-H7 rabbit polyclonal - t=48h







Fig. 7K

T47D + anti-NM23-H7 rabbit polyclonal - t=96h



Stem Cell Lysates + Corresponding Conditioned Media Probed for Presence of NME-7




Fig. 9

Fig. 10B





Fig. 10A



Cancer Stem Cell Marker Expression - T47D




Fig. 11




Fig. 12

Cancer Stem Cell Marker Expression - DU145




Fig. 13

Breast cancer cells are progressed to a more metastatic state by culturing in the presence of 2i inhibitors or in NME7-AB




Fig. 14A

Breast cancer cells are progressed to a more metastatic state by culturing in the presence of a bacterial NME that has high sequence homology to human NME1 or NME7-AB




Fig. 14B

[illegible]

Fig. 16

NME7 specific peptides for generating antibodies to inhibit NME7 for the treatment or prevention of cancers.

The following peptide sequences are identified as being immunogenic peptides giving rise to antibodies that target human NME7 but not human NME1. The sequences were chosen for their lack of sequence homology to human NME1.

1. LALIKPDA
2. MMMLSRKEALDFHVDHQS
3. ALDFHVDHQS
4. EILRDDAICEWKRL
5. FNELIQFITTGP
6. RDDAICEW
7. SGVARTDASESIRALFGTDGIRNAA
8. ELFFPSSGG
9. KFTNCTCCIVKPHAVSEGLLGKILMA
10. LMAIRDAGFEISAMQMFNMDRVNVEEFYEVYKGVVT
11. EFYEVYKGVVTEYHD
12. EIQQNNATKTFREFCGPADPEIARHLRPGTLRAIFGKTKIQNA
13. YSGPCVAM
14. FREFCGP
15. VHCTDLPEDGLLEVQYFFKILDN
16. IQNAVHCTD
17. TDLPEDGLLEVQYFFKILDN
18. PEDGLLEVQYFFK
19. EIINKAGFTITK
20. MLSRKEALDFHVDHQS
21. NELIQFITT
22. EILRDDAICEWKRL
23. SGVARTDASESIRALFGTDGI
24. SGVARTDASES

Fig. 16 (Continued)

25. ALFGTDGI

26. NCTCCIVKPHAVSE

27. LGKILMAIRDA

28. EISAMQMFNMDRVNVE

29. EVYKGVVT

30. EYHDMVTE

31. EFCGPADPEIARHLR

G T N

Fig. 17

NME7 specific peptides for generating antibodies to inhibit NME7 for the treatment or prevention of cancers.

The following are preferred as they are likely areas that are important for structural integrity or for binding to the MUC1* peptide. Bivalent antibodies where in each variable region would bind to each one of a pair are preferred.

35. ICEWKRL

36. LGKILMAIRDA

37. HAVSEGLLGK

38. VTEMYSGP

39. NATKTFREF

40. AIRDAGFEI

41. AICEWKRLLG PAN

42. DHQSRPFF

43. AICEWKRLLG PAN

44. VDHQSRPF

45. PDSFAS

46. KAGEIIEIINKAGFTITK

20/47

Fig. 18

The following peptide sequences are from human NME1 and were selected for their high homology to human NME7 as well as for their homology to other bacterial NME proteins that are able to mimic its function.

47. MANCERTFIAIKPDGVQRGLVGEIHKRFE

48. VDLKDRPF

49. HGSDSVESAEKEIGLWF

Especially preferred for their high homology to human NME7-A or -B and also to HSP 593 are:

50. ERTFIAIKPDGVQRGLVGEIHKRFE

51. VDLKDRPFFAGLVKYMHS GPVVAMVWEGLN

52. NIIHGSDSVESAEKEIGLWFHPEELV

53. KPDGVQRGLVGEI

21/47

Fig. 19

NME7-AB specific peptides preferred for generating antibodies for the treatment or prevention of cancer.

NME7A peptide 1

MLSRKEALDFHVDHQSC

NME7A peptide 2

SGVARTDASESC

NME7B peptide 1

DAGFEISAMQMFNMDRVNVEC

NME7B peptide 2

EVYKGVVTEYHDMVTEC

NME7B peptide 3

AIFGKTKIQNAVHCTDLPEDGLLEVQYFFC

The Test of Whether or not a Cancer Cell is Metastatic or a 'Cancer Stem Cell' is if 200 or less are able to Form a Tumor in a test Animal; Regular Cancer Cells Require 6M or more.

As Few as 50 NME-p-Induced Metastatic Cancer Cells are Sufficient to Form Tumors in test Animals. Animals that were also injected Daily with Human NME-p Developed Metastatic Cancer rather than just Localized Tumors



Figure 20

MOUSE #	# CELLS IMPLANTED	+ NME7	TUMOR	MLTS	SIZE
1	50	Y	Y	Y	M
2	50	N	Y	Y	S
3	50	Y	Y	Y	M
4	50	N	Y	N	M
5	50	Y	Y	Y	M
6	50	N	N	N	No visible tumor
7	100	Y	Y	Y	S
8	100	N	N	N	N
9	100	Y	Y	Y	S
10	100	N	N	N	No visible tumor
11	100	Y	Y	Y	S
12	1,000	Y	Y	Y	S/M
13	1,000	N	Y	N	S/M
14	1,000	Y	N	N	No visible tumor
15	1,000	N	Y	N	L
16	1,000	Y	Y	N	No visible tumor
17	1,000	N	N	N	No visible tumor
18	10,000	Y	Y	Y	S
19	10,000	N	Y	N	S
20	10,000	Y	Y	Y	L
21	10,000	N	Y	N	L
22	10,000	Y	Y	N	L

Figure 21




Fig. 22A





Fig. 22B

Tumor Analysis: Off-Injection-Site 'Bumps' are Human Breast Tumors



KEY

Number is Mouse #

' = primary tumor

r = remote tumor

RED = injected daily with MVE7-AB

KEY MOUSE # AND NUMBER OF CANCER CELLS ENGRAFTED


50 cells: MICE # 1-6

100 cells: MICE # 7-11

1,000 cells: MICE # 12-17

10,000 cells: MICE # 18-22


Fig. 23



KEY MOUSE # AND NUMBER OF CANCER CELLS ENGRAFTED

- 50 cells: MICE # 1-6
- 100 cells: MICE # 7-11
- 1,000 cells: MICE # 12-17
- 10,000 cells: MICE # 18-22

Fig. 24



Anti-NME7 antibodies generated with peptides A1, A2, B1, B2 & B3 bind to NME7 (left) but not to NME1 (right); C20 is an anti-NME1 antibody




Fig. 26

ELISA assay tests NME7 antibodies for ability to block binding of NME7-AB to MUC1* peptide on surface, but not human NME1 dimers




Fig. 27

Anti-NME7 antibodies generated with peptides A1, A2, B1, B2 & B3 as well as the immunizing peptides themselves inhibit cancer cell growth





Fig. 28

Antibodies generated from NME7-specific peptides inhibit cancer cell growth; peptides themselves also inhibit cancer cell growth



T47D cancer cells. antibody combo 52 (A1), 55 (A2), 56 (B1), 57 (B1), 58 (B2), 59 (B2) at 5ug/mL each 30 ug total . All five peptides were combined and added at 1uM each or 10uM each

Fig. 29

Treating cancer cells with anti-NME7 antibodies inhibits transition to “floater” cells, which PCR shows have greatly increased expression of metastatic markers such as CXCR4; xenograft experiments show that the floater cells form a tumor at extremely low copy number – 50 – and thus fulfill the requirement for being classified cancer stem cells or metastatic cancer cells.

Antibodies	Floater observation
Control IgG	100%
53,55,57 (A1,A2,B1)	70%
53,57 (A1,B1)	50%
61 (B3)	5%

JR observations

The number of “Floater” cells, which are the ones that have higher expression of metastatic markers and that form tumors in animals at extremely low copy number is typically 20% of the amount of plated cells by Day 7. Here, we define 100% as the number of floater cells that results when a control antibody is added. Other percentages of floater cells is relative to the control in which a control IgG antibody was added.


Antibodies	Floater observation
Control IgG	100%
53,55,57 (A1,A2,B1)	65%
53,57 (A1,B1)	40%
61 (B3)	5%

VH observations

Fig. 30

Treating cancer cells with anti-NME7 antibodies inhibits transition to metastatic cancer cells; CXCR4 is metastasis marker for cancers


Fig. 31A



T47D breast cancer cells. 2i are MEK and GSK3-beta inhibitors. NF=non-floaters; FL=floaters. 61= rabbit #61 immunized with B3 peptide derived from NME7. Comb-2 = combination of antibodies from rabbits immunized with A1 and A2. Comb-3 = combination of antibodies from rabbits immunized with peptides A1, A2 and B1.

Treating cancer cells with NME7-AB peptides inhibits transition to metastatic cancer cells; CXCR4 and SOX2 are metastasis marker for cancers

Fig. 31B




72 hour

144 hour

Treating cancer cells with NME7-AB peptides inhibits transition to metastatic cancer cells; CXCR4 and SOX2 are metastasis marker for cancers

Fig. 31C



Cell Line T47D	Medium	Total RNA		EEF1A1	
		ng/ μ L	Yield (μ g)	Threshold	Cycle
Control	RPMI + 10% FBS	395.8	47.49		15.2
2i attached	MM + 2i	387.8	11.64		14.5
2i Floaters	MM + 2i	234.6	7.04		14.9
NME7 attached	MM+NME7 4nM	334.8	10.04		14.7
NME7 Floaters	MM+NME7 4nM	259.3	7.78		16.0
2i+ Antibody B3 rabbit 61 attached	MM + 2i	2.7	0.08		25.9
2i+ Antibody B3 rabbit 61 Floaters	MM + 2i	3.6	0.11		25.0
2i + Antibody Combination 3 attached	MM + 2i	44.7	1.34		17.0
2i+ Antibody Combination 3 Floaters	MM + 2i	39.0	1.17		15.9
2i+ Antibody Combination 2 attached	MM + 2i	46.0	1.38		15.6
2i+ Antibody Combination 2 Floaters	MM + 2i	77.8	2.33		16.3
NME7+ Antibody B3 rabbit 61 attached	MM+NME7 4nM	65.7	1.97		17.0
NME7+ Antibody B3 rabbit 61 Floaters	MM+NME7 4nM	15.8	0.47		19.9
NME7+ Antibody Combination 3 attached	MM+NME7 4nM	32.1	0.96		17.0
NME7+ Antibody Combination 3 Floaters	MM+NME7 4nM	109.3	3.28		16.1
NME7+ Antibody Combination 2 attached	MM+NME7 4nM	134.5	4.03		16.1
NME7+ Antibody Combination 2 Floaters	MM+NME7 4nM	139.5	4.19		18.6

Figure 32

NME7-X1 expression in stem and cancer cells




Figure 33

NME7 isoform expression in stem cells and cancer cells




Figure 34




Figure 35




Figure 36

Fig. 37A

SK-OV3 NME7AB 144 hours -
Floaters




Fig. 37B

OV-90 NME7AB 144 hours - Floaters





Fig. 37C



Various cancer cell lines that express MUC1* and NME7-AB-like species

Fig. 38A




Fig. 38B





Fig. 38C



CHL-1 – Melanoma
A2058 – Melanoma
CAPAN-2 – Pancreatic (+)
PANC-1 – Pancreatic (-)


All samples normalized to
40µg/mL
SDIX – Goat anti-Rabbit
VU4H5 and B-9 – Goat anti
Mouse

BT-474 HER2 positive breast cancer cells express almost no MUC1 or MUC1* until they become Metastatic, which is resistant to Herceptin and other chemotherapy drugs. Blocking MUC1* with siRNA or anti-MUC1* Fab reversed the metastatic transition

Fig. 38D



Co-Immunoprecipitation of MUC1* and NME7 in serum grown T47D breast cancer cells




Co-Immunoprecipitation of MUC1* and NME7 from human embryonic stem cells & iPS cells

Fig. 40A

iPS7 HES3

IgG Ab5 IgG Ab5

Fig. 40C



NME7 (B9)

Fig. 40B

MUC1*
(anti-PSMGFR)

