

THE KIT FOR SCREENING COLORECTAL CANCER AND ADVANCED ADENOMA AND ITS APPLICATION

CROSS-REFERENCE TO RELATED APPLICATIONS

[0001] This application claims priority to Chinese Patent Application Serial Nos. 201810502359.7, filed May 23, 2018, and 201810502387.9, filed May 23, 2018, each of which is herein incorporated by reference in its entirety for all purposes.

FIELD OF THE INVENTION

[0002] The present invention relates to compositions and methods for screening colorectal cancer and advanced adenoma, and other applications.

DESCRIPTION OF TEXT FILE SUBMITTED ELECTRONICALLY

[0003] The contents of the text file submitted electronically herewith are incorporated herein by reference in their entirety: A computer readable format copy of the Sequence Listing (filename: NEWH-017_01WO_SeqList_ST25.txt, date recorded: May 22, 2019, file size: 22 kilobytes).

BACKGROUND OF THE INVENTION

[0004] Colorectal cancer (CRC) is the fourth common cancer in the world, with a mortality only less than lung cancer, liver cancer and stomach cancer. The annual death caused by CRC is nearly 700,000. CRC is a "modernized" disease with a higher incidence in developed countries compared to that in developing countries. In the United States, colorectal cancer remains the second leading cause of death (Clinical Interventions in Aging 2016; 11:967–976). With the improvement of the living standards of people, the incidence and mortality of CRC have been increasing in China since 2000 (CA CANCER J CLIN 2016; 66:115–132). The 5-year survival rates of the patients with early-stage, localized disease (stage I and II) approach 90%, while the survival rate of those with late-stage CRC is only 13.1%. The cost of treatment for patients with late-stage CRC often is enormous and can only alleviate the symptoms of the disease (Clinical Interventions in Aging 2016; 11:967–976).

[0005] The development of CRC is a slow process which is generally asymptomatic and difficult to be detected at the early stage until the tumor grows to a few centimeters in size, which may block the

passage of feces and lead to cramping, pain or visible bleeding. The development of CRC has gone through multistep process involving a series of histological, morphological, and genetic changes that accumulate over time: namely, from health, hyperplasia, small polyps, large polyps, and adenocarcinoma to cancer. Polyps are abnormal cells that grow or accumulate locally within the intestinal mucosa. The dividing cells in the polyps may accumulate enough genetic changes to penetrate the intestinal wall and eventually evolve into CRC. However, only a small number of polyps have evolved into CRC after more than a decade of development. Two main types of the malignant potential polyps are adenomas and sessile serrated polyps (SSPs), each of which developed into CRC with different risks. The risk of adenoma developing into CRC is related to its size. Generally, Adenomas are larger in size and have a greater potential for developing into CRC. Advanced adenomas (AA) refer to ≥ 1 cm in size or with $\geq 25\%$ villous component or high-grade dysplasia of any size. Although only about 10% of the most AA become cancerous, 60%-70% of CRC develop from adenomas, the remaining 25%-35% of CRCs develop from SSPs (Clinical Interventions in Aging 2016; 11:967-976). Therefore, early discovery of CRC and AA, and removal of the lesions, can effectively block the progress of CRC to save the lives of patients, significantly improve the patient's 5-year survival rate, and reduce expensive treatment costs in the late-stage of CRC, which greatly reduces the economic burden of the family and society.

[0006] At present, there are several tests for detection of CRC, mainly including colonoscopy, sigmoidoscopy, CT colonography, fecal occult blood test (FOBT) and fecal immunochemical test (FIT).

[0007] The sensitivity of colonoscopy for detecting CRC is $>95\%$. Its screening interval is every 10 years. The advantage of colonoscopy is high sensitivity, which can inspect the entire colon and remove the lesion at the same time. However, the disadvantage of which is invasive examination and bowel preparation will bring discomfort, and the patient needs to be calmed. There is a risk of bowel perforation and bleeding during colonoscopy. These limitations contribute to low compliance with colonoscopy screening.

[0008] The sensitivity of sigmoidoscopy to detect distal colon is greater than 95%. The screening interval of CRC with sigmoidoscopy is every 5 years in combination with FOBT. The advantage of screening CRC with sigmoidoscopy is high sensitivity, no need for systemic sedation, and the lesion

can be removed simultaneously during the examination. The disadvantage of which is semi-invasive examination, easily causing discomfort during examination and the inspection cost is high.

[0009] CT colonography uses radiation to visualize the colon, which sensitivity is $>90\%$ and performed every 5 years. The advantage of which is so high sensitivity that the entire colon can be observed and no sedation is needed. The disadvantage is that the assay is a semi-invasive examination so that patients will easily feel uncomfortable during the screening process. In addition, the lesions cannot be removed at the same time and radiation safety needs to be considered.

[0010] Overview, the tests above to detect CRC based on imaging has high sensitivity, but they are expensive and the bowel preparation is easy to cause discomfort and other side effects. As a result, the patient compliance is low. Additionally, these assays require professional equipment, and doctors with professional skills and rich experience, which may not be available. As a result, the overall screening/detecting rate is low. In addition, some patients are not suitable for these assays. For example, patients with diabetics have a lower success rate of bowel preparation and higher risk of side effects (J Gastrointest Liver Dis 2010; 19: 369-372, World J Gastrointest Endosc 2013; 5: 39-46).

[0011] FOBT and FIT detect hemoglobin in feces of patients by an enzyme reaction and an immunochemical method respectively, with a sensitivity of 33%-75% and 60%-85% respectively for CRC detection, and the tests are performed every 1 year. Although FOBT and FIT are easy to popularize, noninvasive and low cost, the detection rate of precancerous lesions is low (Clinical Interventions in Aging 2016; 11 967-976).

[0012] During the development of polyps into CRC, mutation and methylation changes in some genes such as APC, KRAS, p53, BRAF, NDRG4, BMP3, etc. are accumulated (Clinical Interventions in Aging 2016; 11:967-976). Therefore, detecting of these mutations or methylation changes helps to detect of CRC and precancerous lesions.

[0013] Zou et al. (Clinical Chemistry 2012; 58: 2375-383) used methylation qPCR to detect methylation levels of BMP3, NDRG4, VIM, and TFPI2 genes in tissue samples. In total 37 cases of CRC tissue samples, 25 adenoma tissue samples and 29 healthy human tissue samples were tested. When the specificity was 95%, the sensitivity of BMP3, NDRG4, VIM and TFPI2 genes for CRC detecting was 84%, 92%, 86%, and 92% respectively, and the sensitivity for adenoma detecting was 68%, 76%, 76%, and 88% respectively. It was shown that the detection of genes methylation in colon cancer tissues has high sensitivity and specificity. However, the tissue sampling method is difficult to

be widely used, because the sampling process leads to certain damage to the patient's body. Therefore, it is not suitable for screening CRC and precancerous lesions in the general population.

[0014] Multitarget stool DNA (mt-sDNA) testing includes methylation and mutation detection of tumor exfoliated cells and hemoglobin detection in stool samples, which is screened every 3 years and has the advantages of high sensitivity, non-invasive and easy to popularize (Clinical Interventions In Aging 2016;11:967–976). As a screening method, mt-sDNA can detect CRC and AA early which greatly improves the survival rate of patients. Imperiale et al. (N Engl J Med 2014; 370:1287-97) established a system based on mt-sDNA for methylation detection of BMP3 and NDRG4 genes, point mutation detection of KRAS gene and fecal hemoglobin detection, and then assessed the risk of CRC and AA according to a logistic regression formula. The sensitivity of CRC and AA detection was 92.3% and 42.4% respectively, and the specificity was 86.6%.

[0015] Mt-sDNA is applied to screen sporadic CRC and AA with a advantage of noninvasive compared to colonoscopy and more sensitive compared to FOBT and FIT, but the sensitivity of AA detection is still far lower than that of CRC (Clinical Interventions in Aging 2016; 11:967-976).

[0016] Currently, the products based on mt-sDNA for detecting CRC or AA, such as Cologuard® are mainly developed for European and American populations. No product for CRC and AA detection in the Asian population is available. Particularly, according to “Summary of safety and effectiveness data (SSED)” of Cologuard® issued by the U.S. Food and Drug Administration (www.accessdata.fda.gov/cdrh_docs/pdf13/P130017b.pdf), the sensitivity of AA detection using Cologuard® in the white population and the African-American population is 42.3% and 42.4%, respectively, but the sensitivity of AA detection using the same product in the Asian population is only 30.8%. Therefore, there remains a need to develop an effective system for CRC and/or AA detection in the Asian population to cope with the current increase in the incidence and mortality of colorectal cancer in Asian countries.

[0017] Although there were many studies on the methods for detecting methylation of BMP3 and NDRG4 genes in stool samples from patients with CRC and AA, there is no detailed and comprehensive researches about the hypermethylated CpG sites in BMP3 and NDRG4 genes of the Asian population (ONCOLOGY LETTERS 2014;8:1751-1756; ONCOLOGY LETTERS 2015;9:1383-1387). Moreover, due to limited sample sizes, previous studies of BMP3 and NDRG4 genes methylations in Asian patients were not helpful to identify methylation sites that are most

relevant to CRC and AA. Therefore, there remains a need to determine the exact location of the hypermethylated CpG sites of BMP3 and NDRG4 genes in the Asian population, and to design and optimize kits based on these methylated CpG sites, making the detection of AA more sensitive.

SUMMARY OF THE INVENTION

[0018] The present disclosure provides DNA sequences comprising of hypermethylated CpG sites in the promoter region of BMP3 and NDRG4 genes.

[0019] The present disclosure also provides preferred primers and probes for detecting methylation of BMP3 or NDRG4 genes, and combinations thereof for detecting methylation of both BMP3 and NDRG4 genes.

[0020] The present disclosure further provides a kit for detecting CRC and AA of the Asian population. The DNA sequences comprising of hypermethylated CpG sites in the promoter region of BMP3 and NDRG4 genes can be used as markers for CRC and/or AA detection in Asian population.

[0021] Compared to other primers and probes, the pairs of preferred primers and probes for detecting methylation levels of BMP3 and/or NDRG4 genes have surprisingly higher sensitivity and specificity to detect tumor tissue, such as CRC and AA, and especially AA. In addition, the combinations of these preferred primers and probes of the present disclosure also achieve surprisingly higher sensitivity and specificity to detect tumor tissue, such as CRC and AA, and especially AA.

[0022] A kit for detecting CRC and AA of Asian population based on above preferred primer and probe combinations is also provided.

[0023] In some embodiments, the kit comprises: (1) a preferred combination of pairs of primer and probe and corresponding qPCR reagents; (2) primers and probes for detecting seven mutations () in a coding region of KRAS gene and corresponding qPCR reagents; (3) reagents for detecting hemoglobin in feces.

[0024] In some embodiments, results obtained from an assay using the kit are corrected and analyzed according to a logistic regression formula. In some embodiments, the formula is used to calculate a value for determining the presence or absence of CRC and/or AA. In some embodiments, the formula is $P = e^K / (1 + e^K)$, where P is a comprehensive index, and $K = a * \Delta Ct1 + b * \Delta Ct2 + c * \Delta Ct3 + d * FIT + X$, wherein e is a natural constant, and a, b, c, d, X are clinical constant. In some embodiments, when the P value is equal or more than a predetermined threshold, the result indicates a positive detection of

CRC and/or AA in the patient. In some embodiments, when the P value is less than the threshold, the result indicates a negative detection of CRC and/or AA in the patient, and the patient is determined to be healthy.

[0025] The present disclosure provides a kit for detecting the presence or the absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof. A patient in need thereof is a patient suspected to have CRC and/or AA, such as a patient having at least one sign of developing CRC and/or AA, or a patient having a risk of developing CRC and/or AA, or a subject having a routine medical checkup but otherwise having no sign or risk.

[0026] In some embodiments, the kit comprises a) a first pair of primers and a first probe for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in a biological sample obtained from the patient. In some embodiments, each of the first pair of primers and first probe comprises a contiguous sequence of at least 16 nucleotides that is identical to, complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 1,

[0027] In some embodiments, the kit comprises b) a second pair of primers and a second probe for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in a biological sample obtained from the patient. In some embodiments, each of the second pair of primers and second probe comprises a contiguous sequence of at least 16 nucleotides that is identical to, complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 2,

[0028] In some embodiments, the first pair of primers and the first probe are selected from the group consisting of:

[0029] i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5;

[0030] ii) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11; and

[0031] iii) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17;

[0032] In some embodiments, wherein the second first pair of primers and the second probe are selected from the group consisting of:

[0033] iv) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8;

[0034] v) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14; and

[0035] vi) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20;

[0036] In some embodiments, the kit comprises:

i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

[0037] In some embodiments, the kit comprises:

i) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

[0038] In some embodiments, the kit comprises:

i) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

[0039] In some embodiments, both the first probe and the second probe comprise a fluorescent donor and an acceptor fluorophore.

[0040] In some embodiments, the first probe and the second probe are TAQMAN® probes.

[0041] In some embodiments, the kit further comprises:

(1) means for detecting the presence or absence of at least one mutation in the KRAS gene in the patient; and

(2) means for detecting the presence or absence of hemoglobin in a biological sample obtained from the patient.

[0042] In some embodiments, the means for detecting the presence or absence of at least one mutation in the KRAS gene in the patient comprises at least one pair of primers capable of amplifying the Exon 12 and/or Exon 13 region of the KRAS gene in a polymerase chain reaction (PCR).

[0043] In some embodiments, the means for detecting the presence or absence of hemoglobin in the biological sample comprises an anti-hemoglobin antibody.

[0044] In some embodiments, the primers are capable of amplifying a KRAS gene region comprising at least one KRAS mutation selected from the group consisting G12D, G12V, G12C, G13D, G12A, G12R, G12S, and G13C.

[0045] In some embodiments, the antibody is a colloidal gold-conjugated antibody.

In some embodiments, the kit further comprises means for amplifying an internal quality control gene. An internal control can detect (1) inhibition contamination from sample or extraction method, (2) detect instrument malfunction, (3) chemistry failures (e.g., expired or degraded kit or components, or false combination of reagents), and (4) human error. In some embodiments, the internal control gene is a positive control, such as a gene in a positive control sample that has been determined to have methylation. In some embodiments, the internal control gene is a negative control, such as a gene in a negative control sample that has been determined not to have methylation.

[0046] In some embodiments, the kit further comprises instructions for use and/or interpretation of a test result obtained by using the kit.

[0047] In some embodiments, the kit further comprises means to detect a complex formed by the antibody and the hemoglobin in the biological sample.

[0048] In some embodiments, the biological sample obtained from the patient is a fecal sample.

[0049] In some embodiments, the kit further comprises a bisulfite reagent, and a container suitable for mixing the bisulfite reagent and the biological sample of the patient, or polynucleotides obtained from the biological sample.

[0050] In some embodiments, instead of using bisulfite, the kit further comprises a methylation sensitive restriction enzyme reagent.

[0051] In some embodiments, the kit further comprises: (1) a positive standard and a negative standard for detecting BMP3 methylation in the biological sample, and (2) a positive standard and a negative standard for detecting NDRG4 methylation in the biological sample.

[0052] In some embodiments, the positive standard for detecting BMP3 methylation comprises a polynucleotide sequence of:

GTTAGTTTGGTTCGGGTGTTTTTAAAAATAAAGCGAGGAGGGAAGGTATAGATAGATTTT
GAAAATATTCGGGTATATACGTCGCGATTTATAGTTTTTTTTAGCGTTGGAGTGGAGA
CGGCGTTCGTAGCGTTTTGCGCGGGTGAGGTTTCGCGTAGTTGTTGGGGAAGAGTTTATT
TGTTAGGTTGCGTTGGGTTAGCGTAGTAAGTGGGGTTGGTCGTTATTTTCGTTGTATTCGG
TCGCGTTTTCGGGTTTCGTGCGTTTTTCGTTTTAG (SEQ ID NO: 67);

[0053] In some embodiments, the negative standard for detecting BMP3 methylation comprises a polynucleotide sequence of:

GTTAGTTTGGTTGGGTGTTTTTAAAAATAAAGTGAGGAGGGAAGGTATAGATAGATTTT
GAAAATATTTGGGTATATATGTTGTGATTTATAGTTTTTTTTTAGTGTTGGAGTGGAGA
TGGTGTGTTGTAGTGTTTTGTGTGGGTGAGGTTTGTGTAGTTGTTGGGGAAGAGTTTATTT
GTTAGGTTGTGTTGGGTTAGTGTAGTAAGTGGGGTTGGTTGTTATTTTGTGTTATTTGGT
TGTGTTTTGGGTTTTGTGTGTTTTTGTGTTTTAG (SEQ ID NO: 68);

[0054] In some embodiments, the positive standard for detecting NDRG4 methylation comprises a polynucleotide sequence of:

TGAGAAGTCGGCGGGGGCGCGGATCGATCGGGGTGTTTTTTAGGTTTCGCGTCGCGGTT
TTCGTTTCGTTTTTCGTTTCGTTTATCGGGTATTTTAGTCGCGTAGAAGGCGGAAGTTACG
CGCGAGGGATCGCGGTTTCGTTTCGGGATTAGTTTTAGGTTTCGGTATCGTTTCGCGGGTCG
AGCGTTTATATTCGTTAAATTTACGCGGGTACGTTTTCGCGGCGTATCGTTTTTAGTT
(SEQ ID NO.: 69).

[0055] In some embodiments, the negative standard for detecting NDRG4 methylation comprises a polynucleotide sequence of:

TGAGAAGTTGGTGGGGGTGTGGATTGATTGGGGTGTTTTTTAGGTTTTGTGTTGTGGTTT
TTGTTTGTTTTTTTGTGTTTATTGGGTATTTTAGTTGTGTAGAAGGTGGAAGTTATGTG
TGAGGGATTGTGGTTTGTGTTGGGATTAGTTTTAGGTTTGGTATTGTTTTGTGGGTTGAGT

GTATATATTTGTTAAATTTATGTGGGTATGTTTTGTGGTGTATTGTTTTAGTT (SEQ ID NO.: 70).

[0056] Also provided is a method for detecting the presence or absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof.

[0057] In some embodiments, the method comprises a) obtaining genomic DNA from a biological sample of the patient.

[0058] In some embodiments, the method further comprises b) treating the genomic DNA of a), or a fragment thereof, with one or more reagents to convert cytosine bases that are unmethylated thereof to uracil or another base that is detectably dissimilar to cytosine in terms of hybridization properties.

[0059] In some embodiments, the method further comprises c) contacting the treated genomic DNA, or the treated fragment thereof, with a first pair of primers for detecting the presence or absence of methylation sites of a gene encoding bone morphogenetic protein 3 (BMP3) in the patient. In some embodiments, the method further comprises contacting the treated genomic DNA, or a fragment thereof, with a second pair of primers for detecting the presence or absence of methylation sites of a gene encoding NDRG family member 4 protein (NDRG4) in the patient.

[0060] In some embodiments, the first pair of primers comprise a contiguous sequence of at least 9 nucleotides that is identical to, complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 1. In some embodiments, the second pair of primers comprise a contiguous sequence of at least 9 nucleotides that is complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 2.

[0061] In some embodiments, the treated genomic DNA or the fragment thereof is either amplified to produce at least one amplificate by the first pair of primers or the second pair of primers, or is not amplified.

[0062] In some embodiments, the method further comprises d) determining the presence or absence of CRC or AA in the patient, based on a presence or absence of said amplificate, the methylation state or level of at least one CpG dinucleotide of the BMP3 gene and the NDRG4 gene in the patient.

[0063] In some embodiments, a quantitative PCR is used to amplify the methylated BMP3 gene in the sample. In some embodiments, a quantitative PCR is used to amplify the methylated NDRG4 gene in the sample.

[0064] In some embodiments, the method also comprises using primers for amplifying a reference gene (a.k.a., normalizer, housekeeping gene, or endogenous control). In some embodiments, a quantitative PCR is used to amplify the reference gene in the sample.

[0065] In some embodiments, the first pair of primers and the first probe are selected from the group consisting of:

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5;
- ii) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11; and
- iii) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17.

[0066] In some embodiments, the second first pair of primers and the second probe are selected from the group consisting of:

- iv) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8;
- v) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14; and
- vi) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20.

[0067] In some embodiments, the method comprises using

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and
- ii) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

[0068] In some embodiments, the method comprises using

- i) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

[0069] In some embodiments, the method comprises using

i) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

[0070] In some embodiments, both the first probe and the second probe comprise a fluorescent donor and an acceptor fluorophore. In some embodiments, the first probe and the second probe are TAQMAN® probes.

[0071] In some embodiments, the method further comprises a step of detecting the presence or absence of at least one mutation in the KRAS gene in a biological sample obtained from the patient.

[0072] In some embodiments, the method further comprises a step of detecting the presence or absence of hemoglobin in a biological sample obtained from the patient. In some embodiments, the step of detecting the presence or absence of hemoglobin in the biological sample comprises using an anti-hemoglobin antibody. In some embodiments, the antibody is a colloidal gold-conjugated antibody.

[0073] In some embodiments, the step of detecting the presence or absence of at least one mutation in the KRAS gene in the patient comprises using at least one pair of primers capable of amplifying the Exon 12 and/or Exon 13 region of the KRAS gene in a polymerase chain reaction (PCR). In some embodiments, the primers are capable of amplifying a KRAS gene region comprising at least one KRAS mutation selected from the group consisting G12D, G12V, G12C, G13D, G12A, G12R, G12S, and G13C.

[0074] In some embodiments, the mutant KRAS gene is amplified by one or more pairs of primers selected from the group consisting of:

(1) a forward primer G12D-F comprising SEQ ID NO.: 35, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

- (2) a forward primer G13D-F comprising SEQ ID NO.: 36, and a reverse primer Kras-R comprising SEQ ID NO.: 42;
- (3) a forward primer G12V-F comprising SEQ ID NO.: 37, and a reverse primer Kras-R comprising SEQ ID NO.: 42;
- (4) a forward primer G12C-F comprising SEQ ID NO.: 38, and a reverse primer Kras-R comprising SEQ ID NO.: 42;
- (5) a forward primer G12S-F comprising SEQ ID NO.: 39, and a reverse primer Kras-R comprising SEQ ID NO.: 42;
- (6) a forward primer G12A-F comprising SEQ ID NO.: 40, and a reverse primer Kras-R comprising SEQ ID NO.: 42; and
- (7) a forward primer G12R-F comprising SEQ ID NO.: 41, and a reverse primer Kras-R comprising SEQ ID NO.: 42.

[0075] In some embodiments, the KRAS probe for the qPCR comprises SEQ ID NO.: 46.

[0076] In some embodiments, the amplification of BMP3 gene is performed in a quantitative PCR (qPCR), and the method further comprises amplifying a first reference gene (i.e., a first reference gene) to determining the Ct value of the BMP3 amplification as $\Delta Ct1$.

[0077] In some embodiments, the amplification of NDRG4 gene is performed in a quantitative PCR (qPCR), and the method further comprises amplifying a second reference gene (i.e., a second reference gene) to determining the Ct value of the NDRG4 amplification as $\Delta Ct2$.

[0078] In some embodiments, the amplification of mutant KRAS gene is performed in a quantitative PCR (qPCR), and the method further comprises amplifying a third reference gene (i.e., a third reference gene) to determining the Ct value of the mutant KRAS amplification as $\Delta Ct3$.

[0079] In some embodiments, the first and the second reference genes are the same. In some embodiments, the same reference gene is a B2M gene.

[0080] In some embodiments, the third reference gene is an ACTB gene. In some embodiments, qPCR primers for amplifying ACTB gene comprise SEQ ID NOs.: 43 and 44, and the probe comprise SEQ ID NO.: 46.

[0081] In some embodiments, the method comprises using (1) a positive standard and a negative standard for detecting BMP3 methylation in the sample, and (2) a positive standard and a negative standard for detecting NDRG4 methylation in the sample.

[0082] In some embodiments, the positive standard for detecting BMP3 methylation comprises a polynucleotide sequence of:

GTTAGTTTGGTTCGGGTGTTTTTAAAAATAAAGCGAGGAGGGAAGGTATAGATAGATTTT
GAAAATATTCGGGTATATACGTCGCGATTTATAGTTTTTTTTAGCGTTGGAGTGGAGA
CGGCGTTCGTAGCGTTTTGCGCGGGTGAGGTTTCGCGTAGTTGTTGGGGAAGAGTTTATT
TGTTAGGTTGCGTTGGGTAGCGTAGTAAGTGGGGTTGGTCGTTATTTTCGTTGTATTCGG
TCGCGTTTCGGGTTCGTGCGTTTTTCGTTTTAG (SEQ ID NO: 67);

[0083] In some embodiments, the negative standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTTGGGTGTTTTTAAAAATAAAGTGAGGAGGGAAGGTATAGATAGATTTT
GAAAATATTTGGGTATATATGTTGTGATTTATAGTTTTTTTTTAGTGTTGGAGTGGAGA
TGGTGTGTTGTAGTGTTTTGTGTGGGTGAGGTTTGTGTAGTTGTTGGGGAAGAGTTTATTT
GTTAGGTTGTGTTGGGTAGTGTAGTAAGTGGGGTTGGTTGTTATTTTGTGTTATTTGGT
TGTGTTTTGGGTTTTGTGTGTTTTTGTGTTTTAG (SEQ ID NO: 68);

[0084] In some embodiments, the positive standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTCGGCGGGGGCGCGGATCGATCGGGGTGTTTTTTAGGTTTCGCGTCGCGGTT
TTCGTTTCGTTTTTTCGTTTCGTTTATCGGGTATTTTAGTCGCGTAGAAGGCGGAAGTTACG
CGCGAGGGATCGCGGTTTCGTTTCGGGATTAGTTTTAGGTTTCGGTATCGTTTCGCGGGTCG
AGCGTTTATATTCGTTAAATTTACGCGGGTACGTTTTCGCGGCGTATCGTTTTTAGTT
(SEQ ID NO.: 69).

[0085] In some embodiments, the negative standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

[0086] TGAGAAGTTGGTGGGGGTGTGGATTGATTGGGGTGTTTTTTAGGTTTGTGTTGT
GGTTTTTGTGTTGTTTTTTGTTTGTGTTTATTGGGTATTTTAGTTGTGTAGAAGGTGGAAGTT
ATGTGTGAGGGATTGTGGTTTGTGTTGGGATTAGTTTTAGGTTTGGTATTGTTTTGTGGGT
TGAGTGTTTATATTTGTTAAATTTATGTGGGTATGTTTTTGTGGTGTATTGTTTTTAGTT
(SEQ ID NO.: 70).

In some embodiments, the method comprises amplifying a quality control standard.

[0087] In some embodiments, a method for detecting the presence or absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof, comprising using a kit of the present disclosure.

[0088] The present disclosure further provide a method for detecting the presence or absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof, comprising:

- a) obtaining an untreated genomic DNA from a fecal sample of the patient;
- b) treating the genomic DNA of a), or a fragment thereof, with one or more reagents to convert cytosine bases that are unmethylated thereof to uracil or another base that is detectably dissimilar to cytosine in terms of hybridization properties;
- c) performing a quantitative PCR (qPCR) using the treated genomic DNA of b) as a template, and determining the Ct value of BMP3 gene in the patient as $\Delta Ct1$;
- d) performing a qPCR using the treated genomic DNA of b) as a template, and determining the Ct value of NDRG4 gene in the patient as $\Delta Ct2$;
- e) performing a qPCR using the untreated genomic DNA as a template, and determining the Ct value of a mutant KRAS gene in the patient as $\Delta Ct3$;
- f) performing a fecal immunochemical test of hemoglobin protein in the fecal sample and determining a score as FIT;
- g) determining the value of K, wherein $K = a*\Delta Ct1 + b*\Delta Ct2 + c*\Delta Ct3 + d*FIT + X$, wherein a, b, c, d, X are clinical constants; and
- h) determining the value of a comprehensive index P, wherein $P = e^K / (1 + e^K)$, wherein e is the natural constant.

[0089] The clinical constants a, b, c, d, and X can be determined by analyzing clinical data distribution among a patient population.

[0090] In some embodiments, when P is equal or more than a predetermined threshold value, the patient is determined to have CRC and/or AA, and when P is less than the predetermined threshold value, the patient is determined to be health.

[0091] In some embodiments, the predetermined threshold value is calculated from clinical data distribution, such as clinical data obtained from patients that have been determined to have CRC and/or AA, and patients that have been determined not to have CRC and/or AA.

[0092] In some embodiments, the qPCR for amplifying BMP3 gene comprises a first pair of primers and a first probe, wherein the first pair of primers and the first probe are selected from the group consisting of:

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5;
- ii) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11; and
- iii) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17.

[0093] In some embodiments, the qPCR for amplifying NDRG4 gene comprises a second pair of primers and a second probe, wherein the second pair of primers and the second probe are selected from the group consisting of:

- iv) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8;
- v) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14; and
- vi) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20.

[0094] In some embodiments, the method comprises using

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the sample, and
- ii) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the sample.

[0095] In some embodiments, the method comprises using

- i) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the sample, and

ii) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the sample.

[0096] In some embodiments, the method comprises using

i) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the sample, and

ii) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the sample.

[0097] In some embodiments, both the first probe and the second probe comprise a fluorescent donor and an acceptor fluorophore. In some embodiments, the first probe and the second probe are TAQMAN® probes.

[0098] In some embodiments, the mutant KRAS gene comprising at least one KRAS mutation selected from the group consisting G12D, G12V, G12C, G13D, G12A, G12R, G12S, and G13C.

[0099] In some embodiments, the fecal immunochemical test comprises a colloidal gold-conjugated antibody.

[00100] In some embodiments, step c) and step d) of the method comprises using B2M gene as a reference gene.

[00101] In some embodiments, the method comprises using

(1) a positive standard and a negative standard for detecting BMP3 methylation in the sample, and

(2) a positive standard and a negative standard for detecting NDRG4 methylation in the sample.

[00102] In some embodiments, the positive standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTCGGGTGTTTTTAAAAATAAAGCGAGGAGGGAAGGTATAGATAGATTTT
GAAAATATTCGGGTTATATACGTCGCGATTTATAGTTTTTTTTAGCGTTGGAGTGGAGA
CGGCGTTCGTAGCGTTTTGCGCGGGTGAGGTTGCGGTAGTTGTTGGGGAAGAGTTTATT
TGTTAGGTTGCGTTGGGTTAGCGTAGTAAGTGGGGTTGGTCGTTATTTTCGTTGTATTCGG
TCGCGTTTCGGGTTTCGTGCGTTTTTCGTTTTAG (SEQ ID NO: 67).

[00103] In some embodiments, the negative standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTTGGGTGTTTTTAAAAATAAAGTGAGGAGGGAAGGTATAGATAGATTTT
GAAAATATTTGGGTATATATGTTGTGATTTATAGTTTTTTTTTAGTGTTGGAGTGGAGA
TGGTGTTTGTAGTGTTTTGTGTGGGTGAGGTTTGTGTAGTTGTTGGGGAAGAGTTTATTT
GTTAGGTTGTGTTGGGTAGTGTAGTAAGTGGGGTTGGTTGTTATTTTGTGTATTTGGT
TGTGTTTTGGGTTTTGTGTGTTTTTGTTTTAG (SEQ ID NO: 68).

[00104] In some embodiments, the positive standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTCGGCGGGGGCGCGGATCGATCGGGGTGTTTTTTAGGTTTCGCGTCGCGGTT
TTCGTTTCGTTTTTCGTTTCGTTTATCGGGTATTTTAGTCGCGTAGAAGGCGGAAGTTACG
CGCGAGGGATCGCGGTTTCGTTTCGGGATTAGTTTTAGGTTTCGGTATCGTTTCGCGGGTCG
AGCGTTTATATTCGTTAAATTTACGCGGGTACGTTTTTCGCGGCGTATCGTTTTTAGTT
(SEQ ID NO.: 69).

[00105] In some embodiments, the negative standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTTGGTGGGGGTGTGGATTGATTGGGGTGTTTTTTAGGTTTTGTGTTGTGGTTT
TTGTTTGTTTTTTTGTTTGTTTATTGGGTATTTTAGTTGTGTAGAAGGTGGAAGTTATGTG
TGAGGGATTGTGGTTTGTGTTGGGATTAGTTTTAGGTTTGGTATTGTTTTGTGGGTTGAGT
GTTTATATTTGTTAAATTTATGTGGGTATGTTTTGTGGTGTATTGTTTTTAGTT (SEQ ID
NO.: 70).

In some embodiments, the method comprises amplifying a quality control standard in the step c) and the step d).

[00106] Also provided is a method for diagnosing and treating a colorectal cancer (CRC) and/or advanced adenoma (AA) in a patient in need thereof, comprising determining the presence or absence of CRC and/or AA in the patient by using a kit of the present disclosure, and treating the patient depends on the presence or absence of CRC and/or AA in the patient.

[00107] Also provided is a method for diagnosing and treating a colorectal cancer (CRC) and/or advanced adenoma (AA) in a patient in need thereof, comprising determining the presence or absence

of CRC and/or AA in the patient by using a method described herein, and treating the patient depends on the presence or absence of CRC and/or AA in the patient.

BRIEF DESCRIPTION OF THE DRAWINGS

[00108] Figure 1A to Figure 1D depict the results of CpG island prediction and relative position of amplicons of two BMP3 and NDRG4 genes. “Y”, “R” is degenerate bases. Figure 1A - The result of CpG islands prediction of promoter region of BMP3 genes; Figure 1B - The relative position of amplicons of BMP3 gene; Figure 1C - The result of CpG islands prediction of promoter region of NDRG4 genes; Figure 1D - The relative position of amplicons of NDRG4 gene.

[00109] Figure 2A depicts the difference in methylation CpG sites of BMP in the white and Asian populations. Figure 2B depicts the difference in methylation CpG sites of NDRG4 gene in the white and Asian populations.

[00110] Figure 3A depicts the analytical sensitivity amplification curve of BMP3 with the primers and probes in the preferred group 1. Figure 3B depicts the analytical sensitivity amplification curve of NDRG4 with primers and probes in the preferred group 1. Figure 3C depicts the analytical sensitivity amplification curve of BMP3 with the primers and probes in the preferred group 2. Figure 3D depicts the analytical sensitivity amplification curve of NDRG4 with primers and probes in the preferred group 2. Figure 3E depicts the analytical sensitivity amplification curve of BMP3 with the primers and probes in the preferred group 3. Figure 3F depicts the analytical sensitivity amplification curve of NDRG4 with primers and probes in the preferred group 3. Figure 3G depicts the analytical sensitivity amplification curve of BMP3 with the primers and probes in the comparative group 1. Figure 3H depicts the analytical sensitivity amplification curve of NDRG4 with primers and probes in the comparative group 1. Figure 3I depicts the analytical sensitivity amplification curve of BMP3 with the primers and probes in the comparative group 2. Figure 3J depicts the analytical sensitivity amplification curve of NDRG4 with primers and probes in the comparative group 2. Figure 3K depicts the analytical sensitivity amplification curve of BMP3 with the primers and probes in the comparative group 3. Figure 3L depicts the analytical sensitivity amplification curve of NDRG4 with primers and probes in the comparative group 3.

[00111] Figure 4A depicts the analytical specificity amplification curve of BMP3 with the primers and probes in the preferred group 1. Figure 4B depicts the analytical specificity amplification curve of

NDRG4 with primers and probes in the preferred group 1. Figure 4C depicts the analytical specificity amplification curve of BMP3 with the primers and probes in the preferred group 2. Figure 4D depicts the analytical specificity amplification curve of NDRG4 with primers and probes in the preferred group 2. Figure 4E depicts the analytical specificity amplification curve of BMP3 with the primers and probes in the preferred group 3. Figure 4F depicts the analytical specificity amplification curve of NDRG4 with primers and probes in the preferred group 3. Figure 4G depicts the analytical specificity amplification curve of BMP3 with the primers and probes in the comparative group 1. Figure 4H depicts the analytical specificity amplification curve of NDRG4 with primers and probes in the comparative group 1. Figure 4I depicts the analytical specificity amplification curve of BMP3 with the primers and probes in the comparative group 2. Figure 4J depicts the analytical specificity amplification curve of NDRG4 with primers and probes in the comparative group 2. Figure 4K depicts the analytical specificity amplification curve of BMP3 with the primers and probes in the comparative group 3. Figure 4L depicts the analytical specificity amplification curve of NDRG4 with primers and probes in the comparative group 3.

[00112] Figure 5A to Figure 5C depicts the amplification curves using the primers and probes in the preferred group 1, preferred group 2, and preferred group 3 respectively for detecting BMP3 methylations in clinical sample. Figure 5D to Figure 5F depict the amplification curves using the primers and probes in the comparative group 1, comparative group 2, and comparative group 3 respectively for detecting BMP3 methylations in the same assay.

[00113] Figure 6A to Figure 6C depicts the amplification curves using the primers and probes in the preferred group 1, preferred group 2, and preferred group 3 respectively for detecting NDRG4 methylations in clinical sample. Figure 6D to Figure 6F depict the amplification curves using the primers and probes in the comparative group 1, comparative group 2, and comparative group 3 respectively for detecting NDRG4 methylations in the same assay.

DETAIL DESCRIPTION

Definitions

[00114] References to “one embodiment”, “an embodiment”, “one example”, and “an example” indicate that the embodiment(s) or example(s) so described may include a particular feature, structure,

characteristic, property, element, or limitation, but that not every embodiment or example necessarily includes that particular feature, structure, characteristic, property, element or limitation. Furthermore, repeated use of the phrase “in one embodiment” does not necessarily refer to the same embodiment, though it may.

[00115] “Nucleic acid” or “oligonucleotide” or “polynucleotide”, as used herein means at least two nucleotides covalently linked together. The depiction of a single strand also defines the sequence of the complementary strand. Thus, a nucleic acid also encompasses the complementary strand of a depicted single strand. Many variants of a nucleic acid may be used for the same purpose as a given nucleic acid. Thus, a nucleic acid also encompasses substantially identical nucleic acids and complements thereof. A single strand provides a probe that may hybridize to a target sequence under stringent hybridization conditions. Thus, a nucleic acid also encompasses a probe that hybridizes under stringent hybridization conditions. Nucleic acids may be single stranded or double stranded, or may contain portions of both double stranded and single stranded sequences. The nucleic acid may be DNA, both genomic and cDNA, RNA, or a hybrid, where the nucleic acid may contain combinations of deoxyribo- and ribo-nucleotides, and combinations of bases including uracil, adenine, thymine, cytosine, guanine, inosine, xanthine hypoxanthine, isocytosine and isoguanine. Nucleic acids may be obtained by chemical synthesis methods or by recombinant methods.

[00116] As used herein the phrase “subject in need thereof” refers to an animal or human subject who is known to have cancer, at risk of having cancer (e.g., a genetically predisposed subject, a subject with medical and/or family history of cancer, a subject who has been exposed to carcinogens, occupational hazard, environmental hazard) and/or a subject who exhibits suspicious clinical signs of cancer (e.g., blood in the stool or melena, unexplained pain, sweating, unexplained fever, unexplained loss of weight up to anorexia, changes in bowel habits (constipation and/or diarrhea), tenesmus (sense of incomplete defecation, for rectal cancer specifically), anemia and/or general weakness). Additionally or alternatively, the subject in need thereof can be a healthy human subject undergoing a routine well-being check up.

[00117] As used herein the term “about” refers to $\pm 10\%$.

[00118] The phrase “consisting essentially of” means that the composition or method may include additional ingredients and/or steps, but only if the additional ingredients and/or steps do not materially alter the basic and novel characteristics of the claimed composition or method.

[00119] As used herein, the singular form “a”, “an” and “the” include plural references unless the context clearly dictates otherwise. For example, the term “a compound” or “at least one compound” may include a plurality of compounds, including mixtures thereof.

[00120] “Stringent hybridization conditions” as used herein mean conditions under which a first nucleic acid sequence (e.g., probe) will hybridize to a second nucleic acid sequence (e.g., target), such as in a complex mixture of nucleic acids. Stringent conditions are sequence-dependent and will be different in different circumstances. Stringent conditions may be selected to be about 5-10° C. lower than the thermal melting point (T_m) for the specific sequence at a defined ionic strength pH. The T_m may be the temperature (under defined ionic strength, pH, and nucleic concentration) at which 50% of the probes complementary to the target hybridize to the target sequence at equilibrium (as the target sequences are present in excess, at T_m , 50% of the probes are occupied at equilibrium). Stringent conditions may be those in which the salt concentration is less than about 1.0 M sodium ion, such as about 0.01-1.0 M sodium ion concentration (or other salts) at pH 7.0 to 8.3 and the temperature is at least about 30° C. for short probes (e.g., about 10-50 nucleotides) and at least about 60° C. for long probes (e.g., greater than about 50 nucleotides). Stringent conditions may also be achieved with the addition of destabilizing agents such as formamide. For selective or specific hybridization, a positive signal may be at least 2 to 10 times background hybridization. Exemplary stringent hybridization conditions include the following: 50% formamide, 5×SSC, and 1% SDS, incubating at 42° C., or, 5×SSC, 1% SDS, incubating at 65° C., with wash in 0.2×SSC, and 0.1% SDS at 65° C.

[00121] “Substantially complementary” as used herein means that a first sequence is at least 60%, 65%, 70%, 75%, 80%, 85%, 90%, 95%, 97%, 98% or 99% identical to the complement of a second sequence over a region of 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100 or more nucleotides, or that the two sequences hybridize under stringent hybridization conditions.

[00122] “Substantially identical” as used herein means that a first and a second sequence are at least 60%, 65%, 70%, 75%, 80%, 85%, 90%, 95%, 97%, 98% or 99% identical over a region of 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100 or more nucleotides or amino acids, or with respect to nucleic acids, if the first sequence is substantially complementary to the complement of the second sequence.

[00123] As used herein the term “diagnosing” refers to classifying pathology, or a symptom, determining a severity of the pathology (e.g., grade or stage), monitoring pathology progression, forecasting an outcome of pathology and/or prospects of recovery.

[00124] The phrase “consisting essentially of” means that the composition or method may include additional ingredients and/or steps, but only if the additional ingredients and/or steps do not materially alter the basic and novel characteristics of the claimed composition or method.

1. DNA sequences comprising hypermethylated CpG sites in the promoter region of BMP3 and NDRG4 genes in Chinese population of CRC and AA

[00125] This invention provides DNA sequences comprising of detailed hypermethylated CpG sites in the promoter region of BMP3 and NDRG4 genes in Asian population (e.g., the Chinese population), which is can be used as marker for CRC and AA detection.

[00126] In some embodiments, provided is the natural sequence of BMP3 gene as follows (5' to 3'), showing potentially methylated sites marked by the superscript “m”:

GCCAGTTTGGC^mCGGGTGTTCCTCCAAAAATAAAG^mCGAGGAGGGAAGGTACAGACAGAT
CTTGAAAACACC^mCGGGCCACACA^mCGC^mCG^mCGACCTACAGCTCTTTCTCAG^mCGTTGG
AGTGGAGA^mCGG^mCGCCCGCAG^mCGCCCTG^mCG^mCGGGTGAGGTC^mCG^mCGCAGCTGCTG
GGGAAGAGCCACCTGTCAGGCTG^mCGCTGGGTCAG^mCGCAGCAAGTGGGGCTGGC^mCG
CTATCT^mCGCTGCACCCGGC^mCG^mCGTCC^mCGGGCTC^mCGTG^mCGCCCT^mCGCCCCAG(SEQ
ID NO.: 65).

[00127] In some embodiments, provided is the natural sequence of NDRG4 gene as follows (5' to 3'), showing potentially methylated sites marked by the superscript “m”:

TGAGAAGT^mCGG^mCGGGGG^mCG^mCGGAT^mCGAC^mCGGGGTGTCCCCCAGGCTC^mCG^mCGT
^mCG^mCGGTCCC^mCGCT^mCGCCCTCC^mCGCC^mCGCCCAC^mCGGGCACCCCAGC^mCG^mCGCAG
AAGG^mCGGAAGCCA^mCG^mCG^mCGAGGGAC^mCG^mCGGTC^mCGTC^mCGGGACTAGCCCCAGG
CC^mCGGCAC^mCGCCC^mCG^mCGGGC^mCGAG^mCGCCCACACC^mCGCCAAACCCA^mCG^mCGGG
CA^mCGCCCC^mCG^mCGG^mCGCAC^mCGCCCCCAGCC(SEQ ID NO.: 66).

[00128] After treating the natural genomic DNA, or a fragment thereof, with one or more reagents to convert cytosine bases that are unmethylated thereof to uracil (e.g., by bisulfite) or another base that

is detectably dissimilar to cytosine in terms of hybridization properties, the converted sequence of BMP3 gene is as follow (5' to 3'), showing potentially methylated sites marked by the superscript "m":
 GTTAGTTTGGT^mCGGGTGTTTTTAAAAATAAAG^mCGAGGAGGGAAGGTATAGATAGATTTT
 GAAAATATT^mCGGGTTATATA^mCGT^mCG^mCGATTTATAGTTTTTTTTTTAG^mCGTTGGAGTGGA
 GA^mCGG^mCGTT^mCGTAG^mCGTTTTG^mCG^mCGGGTGAGGTT^mCG^mCGTAGTTGTTGGGGAAG
 AGTTTATTTGTTAGGTTG^mCGTTGGGTTAG^mCGTAGTAAGTGGGGTTGGT^mCGTTATTT^mCG
 TTGTATT^mCGGT^mCG^mCGTTT^mCGGGTTT^mCGTG^mCGTTTT^mCGTTTTAG (SEQ ID NO.:1).

[00129] After treating the natural genomic DNA, or a fragment thereof, with one or more reagents to convert cytosine bases that are unmethylated thereof to uracil (e.g., by bisulfite) or another base that is detectably dissimilar to cytosine in terms of hybridization properties, the converted sequence of NDRG4 gene is as follow (5' to 3'), showing potentially methylated sites marked by the superscript "m":

TGAGAAGT^mCGG^mCGGGGG^mCG^mCGGAT^mCGAT^mCGGGGTGTTTTTTAGGTTT^mCG^mCGT^mC
 G^mCGGTTTT^mCGTT^mCGTTTTTT^mCGTT^mCGTTTAT^mCGGGTATTTTAGT^mCG^mCGTAGAAGG
^mCGGAAGTTA^mCG^mCG^mCGAGGGAT^mCG^mCGGT^mCGTT^mCGGGATTAGTTTTAGGTT^mCGG
 TAT^mCGTTT^mCG^mCGGGT^mCGAG^mCGTTTATATT^mCGTTAAATTTA^mCG^mCGGGTA^mCGTTTT^m
 CG^mCGG^mCGTAT^mCGTTTTTTAGTT (SEQ ID NO.:2).

[00130] The DNA sequences comprising of detailed hypermethylated CpG sites in the promoter region of BMP3 and NDRG4 genes in Asian population of the present disclosure is particularly useful for detecting CRC and/or AA in Asian population. For example, primers and probes can be designed to target one or more specific methylation sites in BMP3 and/or NDRG4 genes, as tools to determine the BMP3 and/or NDRG4 methylation state and level, therefore to determine the tumor condition in a patient in need thereof.

2. The three pairs of preferred primers and probes for detecting methylation of BMP3 and NDRG4 gene respectively and corresponding reagents.

[00131] This invention provides three pairs of preferred primers and probes for detecting methylation levels of BMP3 and NDRG4 genes, respectively. These primers and probes are designed to target the high-frequency methylated CpG sites in the Asian population (e.g., the Chinese population).

[00132] These specific pairs of preferred primers and probes have surprisingly higher sensitivity and specificity in detecting CRC and AA, especially for AA detection in the Asian population, when compared to those in existing commercial products such as Cologuard®.

The sequences of primers and probes are as follows:

Gene	Group	Primer/Probe	Sequence ID	Sequences (5' to 3')
BMP3	Preferred 1	Forward Primer	SEQ ID NO.:3	TTTGAAAATATTCGGGTTATATACGTCGC
		Reverse Primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
		Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTCG
	Preferred 2	Forward Primer	SEQ ID NO.:9	AATATTCGGGTTATATACGTCGCGA
		Reverse Primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCCAA
		Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTCGTAGCGT
	Preferred 3	Forward Primer	SEQ ID NO.:15	AATATTCGGGTTATATACGTCGCGATT
		Reverse Primer	SEQ ID NO.:16	ACTTACTACGCTAACCCAAACG
		Probe	SEQ ID NO.:17	TAGCGTTGGAGTGGAGACGGCGTTCGTA
NDRG4	Preferred 1	Forward Primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
		Reverse Primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAAATACCCGAT
		Probe	SEQ ID NO.:8	CGTCGCGGTTTTTCGTTCTGTTTTTCGTTCTG
	Preferred 2	Forward Primer	SEQ ID NO.:12	GCGGGTGAGAAAGTCGGC
		Reverse Primer	SEQ ID NO.:13	GTAACCTCCGCTTCTACGC
		Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTTCGTT
	Preferred 3	Forward Primer	SEQ ID NO.:18	CGGTTTTCGTTCTGTTTTTTCG
		Reverse Primer	SEQ ID NO.:19	AACCTAAACTAATCCCGAACGAACC
		Probe	SEQ ID NO.:20	TCGTTATCGGGTATTTTAGTCGCGTAG

[00133] The oligonucleotides of the present disclosure, advantageously permit the extremely specific amplification of the hypermethylated CpG sites in the promoter region of BMP3 or NDRG4 in a biological sample obtained from an Asian patient.

[00134] In some embodiments, provided are oligonucleotides that are partially or completely complementary to a sequence of SEQ ID NOs: 3 to 20.

[00135] In some embodiments, provided are oligonucleotides having one or more modifications compared to a probe sequence, such as SEQ ID NOs: 5, 11, 17, 8, 14, and 20. In some embodiments, the modification can happen at the 5' end and/or 3' end of one of the nucleotide sequences recited in SEQ ID NO.: 5, 11, 17, 8, 14, and 20.

[00136] Examples of modified base moieties which can be used to modify nucleotides at any position on its structure include, but are not limited to: 5-fluorouracil, 5-bromouracil, 5-chlorouracil, 5-

iodouracil, hypoxanthine, xanthine, acetylcytosine, 5-(carboxyhydroxymethyl) uracil, 5-carboxymethylaminomethyl-2-thiouridine, 5-carboxymethylaminomethyluracil, dihydrouracil, beta-D-galactosylqueosine, inosine, N-6-sopentenyladenine, 1-methylguanine, 1-methylinosine, 2,2-dimethylguanine, 2-methyladenine, 2-methylguanine, 3-methylcytosine, 5-methylcytosine, N6-adenine, 7-methylguanine, 5-methylaminomethyluracil, methoxyarninomethyl-2-thiouracil, beta-D-mannosylqueosine, 5'-methoxycarboxymethyluracil, 5-methoxyuracil, 2-methylthio-N6-isopentenyladenine, uracil-5-oxyacetic acid, pseudouracil, queosine, 2-thiocytosine, 5-methyl-2-thiouracil, 2-thiouracil, 4-thiouracil, 5-methyluracil, uracil-5-oxyacetic acid methylester, uracil-S-oxyacetic acid, 5-methyl-2-thiouracil, 3-(3-amino-3-N-2-carboxypropyl) uracil, and 2,6-diaminopurine amongst others.

[00137] Examples of modified sugar moieties which may be used to modify nucleotides at any position on its structure include, but are not limited to: arabinose, 2-fluoroarabinose, xylose, and hexose, or a modified component of the phosphate backbone, such as phosphorothioate, a phosphorodithioate, a phosphoramidothioate, a phosphoramidate, a phosphordiamidate, a methylphosphonate, an alkyl phosphotriester, or a formacetal or analog thereof.

[00138] In some embodiments, an oligonucleotide in a sequence of SEQ ID NOs: 5, 11, 17, 8, 14, and 20 is replaced by a unnatural nucleotide, such as an artificial nucleic acid. Artificial nucleic acids include, but are not limited to, peptide nucleic acid (PNA), Morpholino, locked nucleic acid (LNA), glycol nucleic acid (GNA) and threose-nucleic acid (TNA). Each of these is distinguished from naturally occurring DNA or RNA by changes to the backbone of the molecule.

[00139] In some embodiments, a probe of the present disclosure comprises a label at the 5' and the probe.

[00140] In some embodiments, the label at the 5' of a probe comprises a fluorescent dye, such as a fluorophore. As used herein, fluorophore is a fluorescent chemical compound that can re-emit light upon light excitation. Fluorophores typically contain several combined aromatic groups, or planar or cyclic molecules with several π bonds. Non-protein organic fluorophores include, but are not limited to, xanthene derivative (e.g., fluorescein, rhodamine, Oregon green, eosin, and Texas red); cyanine derivatives (e.g., cyanine, indocarbocyanine, oxacarbocyanine, thiocarbocyanine, and merocyanine), squaraine derivatives and ring-substituted squaraines (e.g., Seta, SeTau, and Square dyes), naphthalene derivatives (e.g., dansyl and prodan derivatives), coumarin derivatives; oxadiazole derivatives (e.g.,

pyridyloxazole, nitrobenzoxadiazole and benzoxadiazole); anthracene derivatives (e.g., anthraquinones, including DRAQ5, DRAQ7 and CyTRAK Orange); pyrene derivatives (cascade blue, etc.), oxazine derivatives (e.g., Nile red, Nile blue, cresyl violet, oxazine 170, etc.; acridine derivatives (e.g., proflavin, acridine orange, acridine yellow, etc.); arylmethine derivatives (e.g., auramine, crystal violet, malachite green); tetrapyrrole derivatives (e.g., porphyrin, phthalocyanine, bilirubin). Particular examples include, but are not limited to, VIC, PET, Texas Red, Cy3, Cy5, FAM(6-Carboxyfluorescein), HEX (6-carboxy-2',4,4',5',7,7'-hexachlorofluorescein), ROX (5(6)-Carboxy-X-rhodamine), JOE (6-carboxy-4',5'-dichloro-21,71-dimethoxyfluorescein), TET (5'-tetrachlorofluorescein phosphoramidite), NED (fluorescein benzoxanthene), TAMRA (6-carboxy-N,N,N,N-tetramethylrhodamine), FITC (fluorescein isothiocyanate). Examples of particular fluorophores that can be used in the probes disclosed herein are known to those of skill in the art and include those provided in U.S. Pat. No. 5,866,366 to Nazarenko et al., such as 4-acetamido-4'-isothiocyanatostilbene-2,2'-disulfonic acid; acridine and derivatives such as acridine and acridine isothiocyanate, 5-(2'-aminoethyl)aminonaphthalene-1-sulfonic acid (EDANS), 4-amino-N-[3-vinylsulfonyl]phenyl]naphthalimide-3,5 disulfonate (Lucifer Yellow VS), N-(4-anilino-1-naphthyl)maleimide, anthranilamide; Brilliant Yellow; coumarin and derivatives such as coumarin, 7-amino-4-methylcoumarin (AMC, Coumarin 120), 7-amino-4-trifluoromethylcoumarin (Coumarin 151); cyanosine; 4',6-diaminidino-2-phenylindole (DAPI); 5',5''-dibromopyrogallol-sulfonephthalein (Bromopyrogallol Red); 7-diethylamino-3-(4'-isothiocyanatophenyl)-4-methylcoumarin; diethylenetriamine pentaacetate; 4,4'-diisothiocyanatodihydro-stilbene-2,2'-disulfonic acid; 4,4'-diisothiocyanatostilbene-2,2'-disulfonic acid; 5-[dimethylamino]naphthalene-1-sulfonyl chloride (DNS, dansyl chloride); 4-dimethylaminophenylazophenyl-4'-isothiocyanate (DABITC); eosin and derivatives such as eosin and eosin isothiocyanate; erythrosin and derivatives such as erythrosin B and erythrosin isothiocyanate; ethidium; fluorescein and derivatives such as 5-carboxyfluorescein (FAM), 5-(4,6-dichlorotriazin-2-yl)aminofluorescein (DTAF), 2',7'-dimethoxy-4',5'-dichloro-6-carboxyfluorescein (JOE), fluorescein, fluorescein isothiocyanate (FITC), QFITC (XRITC), 6-carboxy-fluorescein (HEX), and TET (Tetramethyl fluorescein); fluorescamine; IR144; IR1446; Malachite Green isothiocyanate; 4-methylumbelliferone; ortho cresolphthalein; nitrotyrosine; pararosanilin; Phenol Red; B-phycoerythrin; o-phthalaldehyde; pyrene and derivatives such as pyrene, pyrene butyrate and succinimidyl 1-pyrene butyrate; Reactive Red 4 (CIBACRON™ Brilliant

Red 3B-A); rhodamine and derivatives such as 6-carboxy-X-rhodamine (ROX), 6-carboxyrhodamine (R6G), lissamine rhodamine B sulfonyl chloride, rhodamine (Rhod), rhodamine B, rhodamine 123, rhodamine X isothiocyanate, N,N,N',N'-tetramethyl-6-carboxyrhodamine (TAMRA), tetramethyl rhodamine, and tetramethyl rhodamine isothiocyanate (TRITC); sulforhodamine B; sulforhodamine 101 and sulfonyl chloride derivative of sulforhodamine 101 (Texas Red); riboflavin; rosolic acid and terbium chelate derivatives; LightCycler Red 640; Cy5.5; and Cy56-carboxyfluorescein; boron dipyrromethene difluoride (BODIPY); acridine; stilbene; 6-carboxy-X-rhodamine (ROX); Cy3; Cy3.5, Cy5, Cy5.5, VIC® (Applied Biosystems); LC Red 640; LC Red 705; OregonGreen™, CALRed™, Red640; and Yakima yellow; LightCycler®Cyan500; LightCycler®, Red610; Alexa 647; Alexa 555; 5-(2-aminoethyl)amino-1-naphthalene sulfonic acid (EDANS); tetramethyl rhodamine (TMR); tetramethyl rhodamine isocyanate (TMRITC), fluorescein isocyanate (FITC), χ -rhodamine, derivatives thereof, or any combination thereof, amongst others. More fluorescent dyes are described in US Patent NOs: 5866366, 6818431, 6056859, 9140688, 9581587, 6165765, 6485909, 8158358, 7625723, 7560236, 7867701, 9150912, 7960543, 6555383, 6881570, 8198026, 5625081, 8445291, 9194801, 8835110, 7893227, 9243289, 7427674, 9512493, US Patent Application Publication NOs: 20170152552, 20030170672, 20160281151, 20130084558, 20060281100, 20140234833, 20150072340, 20050089910, 20090081677, 201400240220180171393, 20060188886, 20010018185, 20110151446, and WO/2000/017330A1, WO/2008/030071A1, WO/2013/049631A1, WO/2016/179090A1, WO/2016/123895A1, WO/2003/079022A1, each of which is herein incorporated by reference in its entirety.

[00141] In some embodiments, a probe of the present disclosure comprises a fluorescent donor and an acceptor fluorophore. As used herein, an acceptor fluorophore (e.g., a “fluorescent quencher”), is a fluorophore which absorbs energy from a donor fluorophore, for example in the range of about 400 to 900 nm. Acceptor fluorophores generally absorb light at a wavelength which is usually at least 10 nm higher (such as at least 20 nm higher) than the maximum absorbance wavelength of the donor fluorophore. Acceptor fluorophores have an excitation spectrum which overlaps with the emission of the donor fluorophore, such that energy emitted by the donor can excite the quencher. Any acceptor fluorophores known in the art can be utilized. In a particular example, an acceptor fluorophore is a dark quencher, such as Dabcyl, QSY7 (Molecular Probes), QSY9 (Molecular Probes), QSY21 (Molecular Probes), QSY33 (Molecular Probes), BLACK HOLE QUENCHERS™ (Glen Research,

e.g., BHQ-1, BHQ-2, BHQ-3), ECLIPSE™ Dark Quencher (Epoch Biosciences), DDQ-I, DDQ-II, Dabcyl, Eclipse, or IOWA BLACK™ (Integrated DNA Technologies, e.g., Iowa Black FQ, Iowa Black RQ). More fluorescent quenchers are described in US Patent NOs: 9957546, US9274008, US Patent Publication NOs: 20140295422, 20090042205, 20160281182, 20180142284, 20140147929, and WO/2009/009615A1, WO/2016/160572A1, WO/2016/178953A1, WO/2018/229663A1, WO/2010/051544A2, WO/2013/152220A2, each of which is herein incorporated by reference in its entirety. A quencher can reduce or quench the emission of a donor fluorophore. In such an example, instead of detecting an increase in emission signal from the acceptor fluorophore when in sufficient proximity to the donor fluorophore (or detecting a decrease in emission signal from the acceptor fluorophore when a significant distance from the donor fluorophore), an increase in the emission signal from the donor fluorophore can be detected when the quencher is a significant distance from the donor fluorophore (or a decrease in emission signal from the donor fluorophore when in sufficient proximity to the quencher acceptor fluorophore).

[00142] In some embodiments, primers and probes of the present disclosure are based on fluorescence resonance energy transfer (FRET). Examples of oligonucleotides using FRET that can be used to detect amplicons include linear oligoprobes, such as HybProbes, 5' nuclease oligoprobes, such as TAQMAN® probes, hairpin oligoprobes, such as molecular beacons, scorpion primers and UniPrimers, minor groove binding probes, and self-fluorescing amplicons, such as sunrise primers.

[00143] In some embodiments, primers and/or probes of the present disclosure are labeled by other functional entities, such as biotin, haptens, antigens, chemical groups, radioactive substances, enzymatic markers, etc. The detection of a marked amplification product may be accomplished, for example, using fluorescence methods, chemoluminescence methods, densitometry methods, photometry methods, precipitation reactions, enzymatic reactions including enzymatic reinforcement reactions, SPR ("surface plasmon resonance") methods, ellipsometry methods, measurement of the index of refraction, measurement of reflectance, and similar methods.

[00144] In some embodiments, primers and probes described herein can be used in a quantitative PCR to determine the methylation state and level in BMP3 and/or NDRG4 gene in a patient. In some embodiments, an additional reaction can be included to amplify one or more reference gene. In some embodiments, the reference gene is a gene in a patient whose activity would not be affected by the presence or absence of CRC and AA, and would not be affected by the methylation state and level of

BMP3 and NDRG4. In some embodiments, the reference genes include, but are not limited to β -globin (HBB), telomerase (TERT), glyceraldehyde-3-phosphate dehydrogenase (GAPDH), albumin (ALB), β -actin (ACTB), Beta 2 microglobulin (B2M), and T cell receptor γ (TRG).

[00145] In some embodiments, B2M gene is used as a reference gene in the quantitative PCR for detecting methylation state and level of BMP3/NDRG4.

[00146] In some embodiments, one or more other controls can be introduced, including but not limited to, no-template control (for detecting reagent or equipment contamination and to confirm positive results); no-amplification control (for detecting background fluorescence generated by degraded labelled probes), and positive control (for detecting inhibitors or malfunction, and to confirm that reagents and equipment are working).

[00147] In some embodiments, a qPCR is used for determining if there is amplification of methylated BMP3 gene or methylated NDRG4 gene in the sample. The detected signal from the probe of BMP3 or NDRG4 is quantitated either by reference to a standard curve or by comparing the Ct values to that of a reference gene. Analysis of housekeeping genes is often used to normalize the results. The cycle threshold (Ct) is defined as the number of cycles required for the fluorescent signal to cross a predetermined threshold (e.g., exceeding the background level, such as exceeding the level of amplification in a negative control sample). In some embodiments, the threshold is automatically determined by the software of the qPCR instrument or other suitable methods. In some embodiments, the threshold is set just above (e.g., about 0.01%, 0.1%, 1%, 5%, or 10% higher) the terminal fluorescent value in a negative control sample.

[00148] In some embodiments, when the Ct value associated with BMP3 or NDRG4 amplification in a test sample is no more than ($\leq$) about 35, 34, 33, 32, 31, 30, or less, the sample is determined as containing methylated BMP3 or NDRG4, and the patient has CRC and/or AA (positive result), otherwise the sample is determined as not containing the methylated BMP3 or NDRG4, and the patient does not have CRC or AA (negative result). For the reference gene amplification, when the Ct value associated with a control gene amplification in the sample is not more than ($\leq$) about 34, 33, 32, 31, 30, 29 or less, the reference gene amplification is determined to be positive, otherwise the reference gene amplification is determined to be negative. When the reference gene amplification is determined to be negative, the test result is invalidated.

[00149] In some embodiments, the difference between the Ct value associated with BMP3 is compared to the Ct value associated with the reference gene amplification ($\Delta Ct = Ct_{\text{gene of interest}} - Ct_{\text{reference gene}}$) and referred as $\Delta Ct1$. In some embodiments, when $\Delta Ct1$ is not more than a predetermined critical value ($\leq$ the critical value), then the sample is determined to have BMP3 methylation (positive result), and the patient is determined to have CRC or AA. In some embodiments, when $\Delta Ct1$ is more than a predetermined critical value ($>$ the critical value), then the sample is determined not to have BMP3 methylation (negative result), and the patient is determined to be healthy. In some embodiments, the critical value is the corresponding ΔCt value for a sample comprising 5ng/ μ L nucleotide sequence having a methylation rate of 1%, such as about 8, 9, or 10.

[00150] In some embodiments, the difference between the Ct value associated with NDRG4 is compared to the Ct value associated with the reference gene amplification ($\Delta Ct = Ct_{\text{gene of interest}} - Ct_{\text{reference gene}}$) and referred as $\Delta Ct2$. In some embodiments, when $\Delta Ct2$ is not more than a predetermined critical value ($\leq$ the critical value), then the sample is determined to have NDRG4 methylation (positive result), and the patient is determined to have CRC or AA. In some embodiments, when $\Delta Ct2$ is more than a predetermined critical value ($>$ the critical value), then the sample is determined not to have NDRG4 methylation (negative result), and the patient is determined to be healthy. In some embodiments, the critical value is the corresponding ΔCt value for a sample comprising 5ng/ μ L nucleotide sequence having a methylation rate of 1%, such as about 8, 9, or 10.

[00151] Preferred primers and probes of the present disclosure have surprisingly high sensitivity and specificity in detecting CRC and AA in the Asian population, especially for AA detection, when compared to those in existing commercial products such as Cologuard®.

[00152] As used herein, the term “sensitivity” refers to the rate when patients actually having CRC and/or AA in a given population are correctly detected. In some embodiments, the sensitivity is at least about 65%, at least about 70%, at least about 75%, at least about 80%, at least about 85%, at least about 86%, at least about 87%, at least about 88%, at least about 89%, at least about 90%, at least about 91%, at least about 92%, at least about 93%, at least about 94%, at least about 95%, at least about 96%, at least about 97%, at least about 98%, at least about 99%, or more, or 100%. In some embodiments, the size of the population is at least about 50, 100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 2000, 3000, 4000, 5000, 6000, 7000, 8000, 9000, 10,000, or more.

[00153] As used herein, the term “specificity” refers to the rate when patients actually not having CRC or AA a given population are correctly diagnosed as not having the condition. In some embodiments, the specificity is at least about 85%, at least about 86%, at least about 87%, at least about 88%, at least about 89%, at least about 90%, at least about 91%, at least about 92%, at least about 93%, at least about 94%, at least about 95%, at least about 96%, at least about 97%, at least about 98%, at least about 99%, or more, or 100%. In some embodiments, the size of the population is at least about 50, 100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 2000, 3000, 4000, 5000, 6000, 7000, 8000, 9000, 10,000, or more.

[00154] As demonstrated in the Examples, the preferred primers and probes provide surprisingly high sensitivity and specificity in detecting methylation of BMP3 and NDRG4, therefore lead to surprisingly high sensitivity and specificity in detecting CRC and/or AA. For example, the sensitivity of CRC detection with the three pairs of preferred primers and probes of BMP3 gene is at least 85%; the sensitivity of CRC detection with the three pairs of preferred primers and probes of NDRG4 gene is at least 90%; the sensitivity of AA detection with the three pairs of preferred primers and probes of BMP3 gene is at least 66%; the sensitivity of AA detection with the three pairs of preferred primers and probes of NDRG4 gene is at least 73%. Also, the specificity of CRC and AA detection with the three pairs of preferred primers and probes of BMP3 gene is about 97%-100% (e.g., at least about 97.8%); the specificity of CRC and AA detection with the three pairs of preferred primers and probes of NDRG4 gene is also about 97%-100% (e.g., at least about 97.8%).

3. Three sets of preferred primers-probe combinations for detecting CRC and/or AA in patients.

[00155] The present disclosure also provides three sets of preferred primer and probe combinations for detecting methylation levels of BMP3 and NDRG4 genes, in order to determine the presence or absence of CRC or AA in a patient. These specific combinations have surprisingly higher sensitivity and specificity in detecting CRC and AA in the Asian population, especially for AA detection, when compared to those in existing commercial products such as Cologuard®

The sequences of three sets of preferred primers and probes are as follows:

Combination No.	Gene	Group	Primers/ Probes	Sequence ID	Sequence (5' to 3')
4	BMP3	Preferred 1	Forward Primer	SEQ ID NO.:3	TTTGAAAATATTCCGGGTATATACGTCGC
			Reverse Primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTTCG
	NDRG4	Preferred 1	Forward Primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
			Reverse Primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGTTTTTCGTTTCGTTTTTCGTTTCGT
5	BMP3	Preferred 1	Forward Primer	SEQ ID NO.:3	TTTGAAAATATTCCGGGTATATACGTCGC
			Reverse Primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTTCG
	NDRG4	Preferred 2	Forward Primer	SEQ ID NO.:12	GCGGGTGAGAAAGTCGGC
			Reverse Primer	SEQ ID NO.:13	GTAACCTCCGCCTTCTACGC
			Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTTCGTT
7	BMP3	Preferred 2	Forward Primer	SEQ ID NO.:9	AATATTCCGGGTATATACGTCGCGA
			Reverse Primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCCAA
			Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTCTAGCGT
	NDRG4	Preferred 1	Forward Primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
			Reverse Primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGTTTTTCGTTTCGTTTTTCGTTTCGT

[00156] As demonstrated in the Examples, the preferred combinations of BMP3 and NDRG4 primers-probe sets provide surprisingly high sensitivity and specificity in detecting methylation of BMP3 and NDRG4, therefore lead to surprisingly high sensitivity and specificity in detecting CRC and/or AA. In some embodiments, this sensitivity and specificity are obtained when the same assay also includes KRAS gene analysis and hemoglobin test, as explained in Example 5. For example, the overall sensitivity of CRC detection with the three preferred combinations of primers and probes of BMP3 and NDRG4 genes in combination of the KRAS gene analysis and the hemoglobin test is at least 95%; the sensitivity of AA detection with the three preferred combinations of primers and probes of BMP3 and NDRG4 gene in combination of the KRAS gene analysis and the hemoglobin test is at least 93%; the sensitivity of CRC + AA detection with the three preferred combinations of primers and probes of BMP3 and NDRG4 genes in combination of the KRAS gene analysis and the hemoglobin test is at least 97%; the specificity of CRC+ AA detection with the three preferred combinations of primers and probes of BMP3 and NDRG4 genes in combination of the KRAS gene analysis and the hemoglobin test is at least 97%.

4. The kit for detecting CRC and AA

[00157] The present disclosure provides kits for BMP3/NDRG4 methylation detection, and/or for CRC/AA detection in a patient in need therefore. In some embodiments, the kits are particularly suitable for Asian patients, such as Chinese patients.

[00158] In some embodiments, this kit comprise: (1) at least one of the three sets of preferred primer and probe combinations for detecting CRC and AA (Combination Nos. 4, 5, and 7 in Table 20), and corresponding qPCR reagents; (2) means for detecting seven mutations in the coding region of KRAS gene (i.e., G12D、G13D、G12V、G12C、G12S、G12A and G13R), such as suitable primers and probes, and corresponding qPCR reagents; (3) means for detecting hemoglobin in a stool sample, such as reagents based on FIT technology (e.g., an anti-hemoglobin antibody, and reagents for detecting the complex formed by the antibody and hemoglobin in the stool sample).

[00159] In some embodiments, the kit comprises at least one set of the preferred primer-probe combination below:

(1) The three sets of preferred primer and probe combinations is as follows:

NO.	Gene	Group	Primer/Probe	Sequence ID	Sequence (5' to 3')
4	BMP3	Preferred 1	Forward Primer	SEQ ID NO.:3	TTTGAAAATATTCGGGTATATACGTCGC
			Reverse Primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTTCG
	NDRG4	Preferred 1	Forward Primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
			Reverse Primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGGTTTTCGTTTCGTTTTTCGTTTCGT
5	BMP3	Preferred 1	Forward Primer	SEQ ID NO.:3	TTTGAAAATATTCGGGTATATACGTCGC
			Reverse Primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTTCG
	NDRG4	Preferred 2	Forward Primer	SEQ ID NO.:12	GCGGGTGAGAAGTCGGC
			Reverse Primer	SEQ ID NO.:13	GTAACCTCCGCCTTCTACGC
			Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTCGTT
7	BMP3	Preferred 2	Forward Primer	SEQ ID NO.:9	AATATTCGGGTATATACGTCGCGA
			Reverse Primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCCAA
			Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTTCGTAGCGT
	NDRG4	Preferred 1	Forward Primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
			Reverse Primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGGTTTTCGTTTCGTTTTTCGTTTCGT

(2) Triplex quantitative PCR reagents for detecting methylation level of BMP3 and NDRG4 genes.

[00160] In some embodiments, the kits comprise reagents to perform a multiplex PCR for detecting methylation of BMP3 and NDRG4 at the same time. In some embodiments, the multiplex PCR is a quantitative PCR.

[00161] In some embodiments, the reagents include Taq DNA polymerase. In some embodiments, the final concentration of Taq DNA polymerase is about 2U/reaction. In some embodiments, the reagents include MgCl₂. In some embodiments, the final concentration of MgCl₂ in the reaction is 2mM. In some embodiments, the reagents include dNTPs. In some embodiments, the dNTP has a final concentration of 0.2mM. In some embodiments, the reagents include about 0.5mM to 0.75mM primers of amplifying BMP3, NDRG4 and reference genes. In some embodiments, the reagents include about 0.1mM to 0.25mM probes hybridized to DNA sequences of BMP3, NDRG4 and reference genes. In some embodiments, the reagents further include PCR buffer, such as concentrated PCR buffer (e.g., 5x or 10x), which can be diluted to the final concentration of 1x. In some embodiments, B2M is the reference gene and amplified for quality control in quantitative PCR.

(3) Primers and probes for detecting seven mutations in the coding region of KRAS gene.

[00162] In some embodiments, the kits comprise means for detecting mutations in the KRAS gene. In some embodiments, the kits comprise primers and probes that are designed to amplify and detect seven mutant hotspots in Exon 12 and Exon 13 in the open reading region of KRAS gene, which are G12D, G13D, G12V, G12C, G12S, G12A, and G13R. In some embodiments, the sequences of primers and probes are as follows:

Name	Primer/Probe	Sequence ID	Sequences(5' to 3')
G12D-F	Forward primer	SEQ ID NO.:35	AACTTGTGGTAGTTGGAGGTGA
G13D-F	Forward primer	SEQ ID NO.:36	AACTTGTGGTAGTTGGAGCTGGGGA
G12V-F	Forward primer	SEQ ID NO.:37	AACTTGTGGTAGTTGGAGTTGT
G12C-F	Forward primer	SEQ ID NO.:38	AAACTTGTGGTAGTTGGGGCTT
G12S-F	Forward primer	SEQ ID NO.:39	AAACTTGTGGTAGTTGGTGCTA
G12A-F	Forward primer	SEQ ID NO.:40	AACTTGTGGTAGTTGGAGCAGC
G12R-F	Forward primer	SEQ ID NO.:41	AAACTTGTGGTAGTTGGAGCTC
Kras-R	Reverse primer	SEQ ID NO.:42	GAATGGTCCTGCACCAGTAATATG
ACTB-F	Forward primer	SEQ ID NO.:43	AGGGCTTCTTGTCTTTCCTT
ACTB-R	Reverse primer	SEQ ID NO.:44	CGTGCTCGATGGGGTACTTC
KRAS-P	Probe	SEQ ID NO.:45	AGGCAAGAGTGCCTTGACGATACAGC

ACTB-P	Probe	SEQ ID NO.:46	CGTGATGGTGGGCATGGGTCAGAAGGA
--------	-------	---------------	-----------------------------

(4) The multiplex quantitative PCR reagents for detecting codon mutations of KRAS gene.

[00163] Also provided is a multiplex quantitative PCR system for detecting all seven mutations of KRAS gene. In some embodiments, a multiplex quantitative PCR reaction comprise Taq DNA polymerase with final concentration of 2.5U/reaction, MgCl₂ with final concentration of 1mM, dNTPs with final concentration of 0.1mM, 0.3-0.9μM primers of amplifying KRAS and ACTB genes, 0.05-0.1μM probes hybridized to DNA sequences of KRAS and ACTB genes, and 1×PCR buffer. ACTB is a reference gene and amplified for quality control in quantitative PCR.

(5) The kit of detecting hemoglobin in feces

[00164] In some embodiments, the kits comprise reagents for detecting hemoglobin. In some embodiments, hemoglobin is tested qualitatively by enzyme-linked immune sorbent assay (ELISA).

5. The logistic regression model

[00165] In some embodiments, after having the results of the BMP3/NDRG4 methylation test, the KRAS mutation test, and the hemoglobin test, all results are compiled and subjected to a logistic regression model, in order to determine the presence or the absence of CRC and/or AA in the patients.

[00166] In some embodiments, the method comprises calculating the value of a comprehensive cancer index value P, wherein $P = e^K / (1 + e^K)$, wherein e is the natural constant.

[00167] In some embodiments, K is defined to be $K = a * \Delta Ct1 + b * \Delta Ct2 + c * \Delta Ct3 + d * FIT + X$, wherein a, b, c, d, X are clinical constant. $\Delta Ct1$, $\Delta Ct2$ and $\Delta Ct3$ are Ct value of BMP3, NDRG4 and KRAS subtract to that of reference genes.

[00168] In some embodiments, when P value is equal or more than predetermined threshold, the test result is positive, otherwise negative. Positive result indicates the person is possible to have CRC or AA, otherwise is healthy.

6. Treatment methods

[00169] Methods of the present disclosure in some embodiments comprise treating the patients in need after the patients are classified to having colorectal cancer and/or adenoma. In some embodiments, the treating include, but are not limited to, surgery, chemotherapy, radiation therapy, immunotherapy, palliative care, exercise.

[00170] As used herein the phrase “treatment regimen” refers to a treatment plan that specifies the type of treatment, dosage, schedule and/or duration of a treatment provided to a subject in need thereof (e.g., a subject diagnosed with a pathology). The selected treatment regimen can be an aggressive one which is expected to result in the best clinical outcome (e.g., complete cure of the pathology) or a more moderate one which may relieve symptoms of the pathology yet results in incomplete cure of the pathology. It will be appreciated that in certain cases the treatment regimen may be associated with some discomfort to the subject or adverse side effects (e.g., damage to healthy cells or tissue). The type of treatment can include a surgical intervention (e.g., removal of lesion, diseased cells, tissue, or organ), a cell replacement therapy, an administration of a therapeutic drug (e.g., receptor agonists, antagonists, hormones, chemotherapy agents) in a local or a systemic mode, an exposure to radiation therapy using an external source (e.g., external beam) and/or an internal source (e.g., brachytherapy) and/or any combination thereof. The dosage, schedule and duration of treatment can vary, depending on the severity of pathology and the selected type of treatment, and those of skills in the art are capable of adjusting the type of treatment with the dosage, schedule and duration of treatment.

[00171] In some embodiments, the treatments include, but is not limited to, fluorouracil, capecitabine, oxaliplatin, irinotecan, UFT, FOLFOX, FOLFOXIRI, and FOLFIRI, antiangiogenic drugs such as bevacizumab, and epidermal growth factor receptor inhibitors (e.g., cetuximab and panitumumab).

7. Advantages of the present invention

Without wishing to be bound by any particular theory, the present invention at least has the following advantages:

- (1) The detailed methylated CpG sites of promoter regions of BMP3 and NDRG4 genes in the Chinese population are provided and can be as a biomarker for detecting CRC and/or AA.
- (2) The kit is more suitable for CRC and AA detection in Chinese population than other similar products based on mt-sDNA, such as Cologuard®, because the present invention targets Chinese population specific methylation CpG sites, while the previous products are directed to the white people population.
- (3) The sensitivity and specificity of CRC detection are higher compared to other similar products based on mt-sDNA.
- (4) The sensitivity and specificity of AA detection are significantly improved compared to other

similar products based on mt-sDNA.

- (5) The method is noninvasive and easy for home sampling, therefore leads to good patient compliance, so that it can be widely used as a CRC and AA screening method. The method reduces the incidence and mortality due to CRC in the Asian population.

EXAMPLES

Example 1

Discovery that methylated CpG sites in the promoter regions of BMP3 and NDRG4 genes respectively in Chinese CRC and AA population.

(1) Sample collection

[00172] In total 191 colon FFPE tissue samples were collected from patients with CRC and AA confirmed by colonoscopy, including 50 colorectal cancer tissues and 49 paired adjacent normal tissues, 46 adenomas cancer tissues and 46 paired adjunct normal tissues.

(2) DNA extraction

[00173] Genomic DNAs were extracted from FFPE samples with the TaKaRa MiniBEST FFPE DNA Extraction Kit (catalogue number: 9782). The detailed operation steps are described as follows:

- i. Scrap 30mg paraffin section tissue with a sterile scalpel and remove excess paraffin.
- ii. Put the paraffin section tissue in a 1.5 mL centrifuge tube, and add 500 μ L of buffer DP, Mix and incubate in water at 80°C for 1 minute, and then vortex for 10 seconds. Add 180 μ L of Buffer GL and vortex.
- iii. The mixture was centrifuged at 12,000 rpm for 1 minute at room temperature, and then the solution was formed into two layers (upper oil phase, lower aqueous phase). Add 20 μ L of Proteinase K (20 mg/mL) and 10 μ L of RNase(10 mg/mL) to the lower aqueous phase, and mix thoroughly with pipette up and down gently. Be careful not to disturb the layers. And then water bath at 56°C for 1 hours.
- iv. The solution of the previous step was incubated at 90°C for 30 minutes and cooled to room temperature. And then add 200 μ L of Buffer GB and 200 μ L of 100% ethanol to the solution, vortex for 10 seconds. Centrifuge at 12,000 rpm for 1 minute at room template, and the solution

was formed into two layers (upper oil phase, lower aqueous phase)

- v. Put a spin column in a collection tube, and add lower aqueous phase of the previous step into the spin column. Be careful not to disturb the layers, and centrifuge at 12,000 rpm for 2 minutes at room temperature, and then discard the waste.
- vi. Add 500 μ L of Buffer WA to the spin column, and centrifuge at 12,000 rpm for 1 minute at room temperature, and then discard the waste.
- vii. Add 500 μ L of Buffer WB to the spin column, centrifuge at 12,000 rpm for 1 minute at room temperature, and then discard the waste. And repeat once.
- viii. Put the spin column in the collection tube, centrifuge at 12,000 rpm for 2 minutes at room temperature.
- ix. Put the spin column in a new 1.5 mL centrifuge tube, add 50-100 μ L of sterilized water or elution buffer to the center of the spin column membrane, and then place it at room temperature for 5 minutes.
- x. Centrifuge at 12,000 rpm for 2 minutes at room temperature, and DNA is eluted.
- xi. The eluted DNA is quantified with Nanodrop 2000 fluorometer, and stored at -20°C for use.

(3) The prediction of CpG islands in promoter regions of BMP3 and NDRG4 genes and primer designing for amplicon sequencing.

- i. The prediction of CpG islands of promoter regions of BMP3 and NDRG4 genes.

[00174] The promoter sequences of BMP3 and NDRG4 genes were downloaded, including DNA sequences approximately 1000-1500 bp upstream from transcription start site (TSS) and 5'UTR region. The CpG islands of the sequences were predicted with MethPrimer software (www.urogene.org/methprimer/). As shown in Figure 1A to Figure 1D, two larger ones of the three CpG islands of BMP3 gene are located at approximately 400 bp upstream from TSS and whole 5'UTR region(chr4: 81951752-81952760), and only one CpG island of NDRG4 gene is located at approximately 500 bp upstream from TSS and partial 5'UTR region(chr16:58497061-58497938). The build of human reference genome is GRCh37/hg19.

- ii. Primer designing to amplify sequence of CpG islands predicted of BMP3 and NDRG4 genes.

[00175] Based on the length of the sequences to be amplified and read length of sequencing, four and five pairs of primers were designed for BMP3 and NDRG4 genes, respectively. There is as much overlap as possible between adjacent amplicons so that the CpG islands of the two genes can be sequenced thoroughly. The primers of the two genes are listed in Table 1, and the relative position of amplicons of two genes is shown in Figure 1A to Figure 1D.

Table 1 the primers for amplicon sequencing of BMP3 and NDRG4 genes

Gene	Amplicon	Primer	Sequence ID	Sequence (5' to 3')
BMP3	Amplicon 1	Forward Primer 1	SEQ ID NO.:47	AGTTTGGTGTAAGTTAAGAG
BMP3		Reverse Primer 1	SEQ ID NO.:48	CTAACTCTATTTTAAACRCCA
BMP3	Amplicon 2	Forward Primer2	SEQ ID NO.:49	GTTTTAATTTTGGAAAAGGTAA
BMP3		Reverse Primer 2	SEQ ID NO.:50	ACCTAACAAATAAACTCTTCC
BMP3	Amplicon 3	Forward Primer 3	SEQ ID NO.:51	GAAGGTATAGATAGATTTTGAA
BMP3		Reverse Primer 3	SEQ ID NO.:52	CACCTAACACAACCTTACRAAACT
BMP3	Amplicon 4	Forward Primer 4	SEQ ID NO.:53	GTATTTAGTTATGGTTGGGGYGAGTA
BMP3		Reverse Primer 4	SEQ ID NO.:54	CTCACCTACTACTACCGCCCR
NDRG4	Amplicon 1	Forward Primer 1	SEQ ID NO.:55	AGGTTTTTGAGTTTTTGGTTTTTTT
NDRG4		Reverse Primer 1	SEQ ID NO.:56	CCCTCCAAACCCCTATAAC
NDRG4	Amplicon 2	Forward Primer2	SEQ ID NO.:57	GGATGGGGATGTTTTTGTAG
NDRG4		Reverse Primer 2	SEQ ID NO.:58	RGRGAAACCTAAAAAACACC
NDRG4	Amplicon 3	Forward Primer 3	SEQ ID NO.:59	GYGGAGYGGGTGAGAAGT
NDRG4		Reverse Primer 3	SEQ ID NO.:60	CRAACAACCAAAAACCCCTC
NDRG4	Amplicon 4	Forward Primer 4	SEQ ID NO.:61	GTTYGTTYGGGATTAGTTTGTAGG
NDRG4		Reverse Primer 4	SEQ ID NO.:62	CRCAAAACRAAAAACRAAAC
NDRG4	Amplicon 5	Forward Primer 5	SEQ ID NO.:63	GYGGYGTTTYGTTTTTG
NDRG4		Reverse Primer 5	SEQ ID NO.:64	CRACRACTAAAAATCCCCAA

[00176] DNA samples were subjected to bisulfite treatment as described below.

(4) Bisulfite treatment

- The extracted DNA is put at room temperature to thaw, and the DNA concentration is diluted to 20 ng/ μ L. Add 40 μ L of the diluted DNA to a 1.5 mL centrifuge tube, and then add 4 μ L of 3M NaOH solution, incubate at 42°C for 20 minutes.
- Add 400 μ L of a conversion solution, mix, and incubate at 50°C of 16 hours in the dark.
- Add 550 μ L of binding solution, mix, and transfer the solution to a DNA purification column. Centrifuge at 13,000 rpm for 90 seconds and discard the waste. Centrifuge again for 3 minutes and discard the waste.
- Add 600 μ L of 90% ethanol to the DNA purification column, centrifuge at 13,000 rpm for 90 seconds, and discard the waste. Centrifuge again at 13,000 rpm for 15 seconds.
- Add 300 μ L of desulfurization solution (0.3 M NaOH of 90% ethanol solution) to the DNA

purification, put it at room temperature for 30 minutes. Centrifuge at 13,000 rpm for 90 seconds and discard the waste.

vi. Add 600 μL of 90% ethanol, centrifuge at 13,000 rpm for 90 seconds, and discard the waste.

Repeat this step once, and centrifuge again at 13,000 rpm for 3 minutes.

vii. Put the DNA purification column into a new 1.5 mL centrifuge tube, and add 40 μL of eluate.

Incubate the tube at 50°C for 30 minutes, and centrifuge at 13,000 rpm for 90 seconds. Store the converted DNA solution at -20°C for use.

[00177]

(5) Library preparation and amplicon sequencing

(a) Multiplex PCR amplification

i. Prepare PCR master mixture as follows:

Table 2

No.	reagents	Final concentration
1	2 $\times$ QIAGEN Multiplex PCR Master Mix(QIAGEN, Catalogue NO.: 206143)	1 $\times$
2	10 $\times$ primer mix(2 μM each)	0.2 μM
3	RNase-free Water	Up to 40 μL

ii. Vortex the mixture gently and pipette 40 μL of which to each PCR tubes, and then add 10 μL of bisulfate treated DNA.

iii. Vortex gently and PCR amplification was preformed as follows: one cycle of denatured at 95°C for 15 minutes, 35 cycles of denatured at 94°C for 30 seconds, annealed at 55°C for 90 seconds, and elongated at 72°C for 90 seconds, one cycle of elongated at 72°C for 10 minutes, and finally held at 4°C forever.

iv. 5 μL of PCR products were blended with 6X loading buffer (TakaRa, catalogue No.:9156) and loaded to 1% (w/v) agarose gel with a control of DL2000 DNA marker (TakaRa, catalogue No.:3427Q) , and electrophoresis at 120V for 40 minutes.

v. If there was a non-specific amplification, it was necessary to recovered remaining 45 μL of PCR products after electrophoresis according to the kit instructions (QIAGEN, catalogue No.: 28704).

[00178]

(b) Purification and quantification of PCR products

PCR product was purified with MinElute PCR Purification Kit according to the kit instructions (QIAGEN, Catalogue No.: 28004). Purified PCR product was quantified with Qubit™ dsDNA BR Assay Kit (Catalogue No.: Q32850).

[00179]

(c) Adapter ligation

Adapters were ligated with purified PCR products with NEBNext Quick Ligation Module (NEB, Catalogue Number: E6056L), and the reaction mixture were prepared as follows:

Table 3. Reaction Mixture for adapter ligation

No.	Reagents	Volume (μL)/reaction
1	NEBNext Quick Ligation Reaction Buffer(5×)	10
2	Adapter	5
3	Quick T4 DNA Ligase	5
4	RNase-free Water	10

Add 30 μL of mixture to each PCR tubes and then add 20 μL of purified PCR products respectively.

The solution were mixed gently, incubated at 20°C for 15 minutes and held at 4°C for 10 minutes.

[00180]

(d) Purification of ligation products

- i. Put Agencourt AMPure XP beads (Beckman, Catalogue number: A63882) to room temperature for use.
- ii. Centrifuge the PCR tubes at 280g for 1 minute at 20°C, and transfer the 50 μL of ligation products to a new 96-well PCR plate.
- iii. Vortex the AMPure XP beads for 30 seconds to be dispersed evenly, and add 56 μL of beads to each well of the PCR plate. Pipette the mixture gently for 10 times to mix, and keep still at room temperature for 5 minutes.
- iv. Place the 96-well PCR plate on a 96-well magnetic plate, keep it still for 2 minutes until the supernatant was clear. Keep the 96 well PCR plate on the 96-well magnetic plate and remove the supernatant, and then add 200 μL of freshly prepared 80% ethanol to each wells, and incubate at room temperature for 30 seconds, and remove the supernatant.
- v. Keep the 96-well PCR plate on the 96-well magnetic plate and add 200 μL of freshly prepared

80% ethanol to each wells. And then incubate at room temperature for 30 seconds. Discard the supernatant and remove the residual ethanol with 10 μ L pipette.

- vi. Keep the 96-well PCR plate on the magnetic plate and let the beads dry naturally for 10 minutes.
- vii. Remove the 96-well PCR plate from the magnetic plate and add 27.5 μ L of 10mM Tris (pH 8.5) to each wells of the 96-well plate. Pipette up and down gently to mix beads and Tris for 10 times until the beads were dispersed thoroughly. Keep the plate still at room temperature for 2 minutes.
- viii. Place the 96-well PCR plate on the magnetic plate and still for 2 minutes until the supernatant is clear. Pipette 25 μ L of the supernatant to a new PCR tube and store it at -20°C for use.

[00181]

(e) PCR amplification and product purification

PCR master mixture was prepared as follows and mixed gently. Add 40 μ L of mixture to each PCR tubes and then add 5 μ L of purified ligation products.

Table 4. PCR reaction mixture

NO.	Reagents	Final concentration
1	NEBNext Ultra Q5 Master Mix (NEB, catalogue number: M0544)	1×
2	Forward Primer(10 μ M)	1 μ M
3	Reverse Primer(10 μ M)	1 μ M
4	Nuclease-free water	Up to 40 10 μ L

Vortex the PCR tubes gently and centrifuge briefly, and carries PCR reaction as the following condition: one cycle of denaturing at 98°C for 30 seconds, 8 to 15 cycles of denaturing at 98°C for 10 seconds and elongating at 65°C for 75 seconds, one cycle of denaturing at 65°C for 5 minutes, and finally hold at 4°C. PCR products were purified according to (d) Purification of ligation products.

[00182]

(f) Detection and sequencing of library

The concentration and size distribution of purified PCR products were analyzed with Qubit™ dsDNA BR Assay kit (catalogue number: Q32850) and the Agilent 2100 Bioanalyzer Instrument, respectively. Libraries were pooling with equal molar concentration and sequenced on Illumina HiSeq2500 with read length of PE125.

[00183]

(g) Sequencing data analysis

Raw data was demultiplexed according sample indexes and sequencing reads were mapped to reference gene sequence of BMP3 and NDRG4 with SHRiMP V2.04 software. Methylated CpG sites were identified based on mapping results. And finally, It founded that 26 and 39 CpG sites of BMP3 and NDRG4 genes were hypermethylated at CRC and AA tissues compared to adjacent normal tissues($p < 0.05$), which indicated that these methylated CpG sites can be a DNA biomarker for early diagnosis of CRC and AA.

Table 5 Methylated frequency of CpG sites of BMP3 gene in CRC tissues and adjacent normal tissues.

CpG sites	Chromosome coordinate (hg19)	Location to transcript start site (NM_001201.3)	The number of methylated samples divided by the total number of samples- CRC	Methylation percentage- CRC	The number of methylated samples divided by the total number of samples- Normal	Methylation percentage- Normal	P-value
#1	81952078	-41	34/50	0.68	2/49	0.04	<0.01
#2	81952099	-20	38/50	0.76	2/49	0.04	<0.01
#3	81952135	17	33/50	0.66	1/49	0.02	<0.01
#4	81952146	28	34/50	0.68	2/49	0.04	<0.01
#5	81952149	31	34/50	0.68	1/49	0.02	<0.01
#6	81952151	33	35/50	0.7	2/49	0.04	<0.01
#7	81952172	54	33/50	0.66	0/49	0.00	<0.01
#8	81952186	68	33/50	0.66	1/49	0.02	<0.01
#9	81952189	71	36/50	0.72	1/49	0.02	<0.01
#10	81952193	75	37/50	0.74	0/49	0.00	<0.01
#11	81952198	80	34/50	0.68	1/49	0.02	<0.01
#12	81952205	87	37/50	0.74	2/49	0.04	<0.01
#13	81952207	89	37/50	0.74	2/49	0.04	<0.01
#14	81952218	100	34/50	0.68	1/49	0.02	<0.01
#15	81952220	102	37/50	0.74	2/49	0.04	<0.01
#16	81952255	137	35/50	0.7	0/49	0.00	<0.01
#17	81952266	148	33/50	0.66	1/49	0.02	<0.01
#18	81952285	167	38/50	0.76	1/49	0.02	<0.01
#19	81952293	175	36/50	0.72	1/49	0.02	<0.01
#20	81952302	184	35/50	0.7	1/49	0.02	<0.01
#21	81952306	188	33/50	0.66	1/49	0.02	<0.01
#22	81952308	190	36/50	0.72	0/49	0.00	<0.01
#23	81952313	195	35/50	0.7	0/49	0.00	<0.01
#24	81952320	202	35/50	0.7	2/49	0.04	<0.01
#25	81952324	206	38/50	0.76	2/49	0.04	<0.01
#26	81952330	212	33/50	0.66	2/49	0.04	<0.01

Table 6 Methylated frequency of CpG sites of BMP3 gene in AA tissues and adjacent normal tissues.

CpG sites	Chromosome coordinate (hg19)	Location to transcript start site (NM_001201.3)	The number of methylated samples divided by the total number of samples-AA	Methylation percentage-AA	The number of methylated samples divided by the total number of samples-AA	Methylation percentage-Normal	P value
#1	81952078	-41	30/50	0.60	2/49	0.04	<0.01
#2	81952099	-20	29/50	0.58	1/49	0.02	<0.01
#3	81952135	17	32/50	0.64	2/49	0.04	<0.01
#4	81952146	28	29/50	0.58	2/49	0.04	<0.01
#5	81952149	31	30/50	0.60	2/49	0.04	<0.01
#6	81952151	33	32/50	0.64	0/49	0.00	<0.01
#7	81952172	54	32/50	0.64	2/49	0.04	<0.01
#8	81952186	68	32/50	0.64	2/49	0.04	<0.01
#9	81952189	71	28/50	0.56	0/49	0.00	<0.01
#10	81952193	75	29/50	0.58	1/49	0.02	<0.01
#11	81952198	80	29/50	0.58	0/49	0.00	<0.01
#12	81952205	87	31/50	0.62	1/49	0.02	<0.01
#13	81952207	89	29/50	0.58	0/49	0.00	<0.01
#14	81952218	100	32/50	0.64	1/49	0.02	<0.01
#15	81952220	102	28/50	0.56	1/49	0.02	<0.01
#16	81952255	137	29/50	0.58	0/49	0.00	<0.01
#17	81952266	148	31/50	0.62	0/49	0.00	<0.01
#18	81952285	167	32/50	0.64	1/49	0.02	<0.01
#19	81952293	175	33/50	0.66	0/49	0.00	<0.01
#20	81952302	184	32/50	0.64	1/49	0.02	<0.01
#21	81952306	188	30/50	0.60	0/49	0.00	<0.01
#22	81952308	190	31/50	0.62	1/49	0.02	<0.01
#23	81952313	195	29/50	0.58	1/49	0.02	<0.01
#24	81952320	202	32/50	0.64	1/49	0.02	<0.01
#25	81952324	206	31/50	0.62	2/49	0.04	<0.01
#26	81952330	212	32/50	0.64	1/49	0.02	<0.01

Table 7 Methylated frequency of CpG sites of NDRG4 gene in CRC tissues and adjacent normal tissues.

CpG sites	Chromosome coordinate (hg19)	Location to transcript start site (NM_020465.3)	The number of methylated samples divided by the total number of samples-CRC	Methylation percentage-CRC	The number of methylated samples divided by the total number of samples-normal tissue	Methylation percentage-Normal	P value
#1	chr16:58497349	-200	38/50	0.75	1/49	0.02	<0.01
#2	chr16:58497352	-197	35/50	0.70	2/49	0.04	<0.01
#3	chr16:58497358	-191	35/50	0.70	2/49	0.04	<0.01

#4	chr16:58497360	-189	39/50	0.78	1/49	0.03	<0.01
#5	chr16:58497365	-184	39/50	0.78	1/49	0.03	<0.01
#6	chr16:58497369	-180	39/50	0.77	1/49	0.02	<0.01
#7	chr16:58497388	-161	37/50	0.74	1/49	0.02	<0.01
#8	chr16:58497390	-159	37/50	0.73	1/49	0.03	<0.01
#9	chr16:58497393	-156	40/50	0.79	1/49	0.03	<0.01
#10	chr16:58497395	-154	37/50	0.74	1/49	0.02	<0.01
#11	chr16:58497402	-147	38/50	0.76	1/49	0.03	<0.01
#12	chr16:58497406	-143	39/50	0.78	1/49	0.02	<0.01
#13	chr16:58497414	-135	37/50	0.74	1/49	0.03	<0.01
#14	chr16:58497418	-131	36/50	0.72	0/49	0.00	<0.01
#15	chr16:58497425	-124	36/50	0.72	1/49	0.02	<0.01
#16	chr16:58497438	-111	36/50	0.72	2/49	0.04	<0.01
#17	chr16:58497440	-109	36/50	0.71	1/49	0.02	<0.01
#18	chr16:58497449	-100	37/50	0.73	2/49	0.04	<0.01
#19	chr16:58497458	-91	39/50	0.77	2/49	0.04	<0.01
#20	chr16:58497460	-89	37/50	0.74	1/49	0.03	<0.01
#21	chr16:58497462	-87	38/50	0.75	1/49	0.03	<0.01
#22	chr16:58497470	-79	35/50	0.70	1/49	0.02	<0.01
#23	chr16:58497472	-77	36/50	0.71	2/49	0.04	<0.01
#24	chr16:58497477	-72	36/50	0.72	1/49	0.03	<0.01
#25	chr16:58497481	-68	37/50	0.74	1/49	0.02	<0.01
#26	chr16:58497499	-50	37/50	0.73	1/49	0.02	<0.01
#27	chr16:58497505	-44	40/50	0.79	2/49	0.04	<0.01
#28	chr16:58497510	-39	40/50	0.79	1/49	0.02	<0.01
#29	chr16:58497512	-37	40/50	0.79	2/49	0.04	<0.01
#30	chr16:58497517	-32	35/50	0.70	1/49	0.02	<0.01
#31	chr16:58497521	-28	39/50	0.77	1/49	0.02	<0.01
#32	chr16:58497531	-18	38/50	0.76	1/49	0.03	<0.01
#33	chr16:58497542	-7	36/50	0.71	1/49	0.02	<0.01
#34	chr16:58497544	-5	37/50	0.74	1/49	0.03	<0.01
#35	chr16:58497550	2	37/50	0.74	0/49	0.00	<0.01
#36	chr16:58497556	8	38/50	0.75	1/49	0.03	<0.01
#37	chr16:58497558	10	39/50	0.78	2/49	0.04	<0.01
#38	chr16:58497561	13	36/50	0.71	2/49	0.04	<0.01
#39	chr16:58497566	18	35/50	0.70	1/49	0.03	<0.01

Table 8 Methylated frequency of CpG sites of NDRG4 gene in AA tissues and adjacent normal tissues.

CpG sites	Chromosome coordinate (hg19)	Location to transcript start site (NM_020465.3)	The number of methylated samples divided by the total number of samples-AA	Methylation percentage-AA	The number of methylated samples divided by the total number of samples-normal tissue	Methylation percentage-Normal	P value
#1	chr16:58497349	-200	28/46	0.60	1/46	0.02	<0.01
#2	chr16:58497352	-197	28/46	0.60	1/46	0.02	<0.01
#3	chr16:58497358	-191	27/46	0.59	1/46	0.03	<0.01
#4	chr16:58497360	-189	29/46	0.62	1/46	0.03	<0.01
#5	chr16:58497365	-184	30/46	0.65	2/46	0.04	<0.01
#6	chr16:58497369	-180	28/46	0.61	2/46	0.04	<0.01
#7	chr16:58497388	-161	29/46	0.64	1/46	0.03	<0.01
#8	chr16:58497390	-159	29/46	0.64	1/46	0.03	<0.01
#9	chr16:58497393	-156	30/46	0.65	1/46	0.02	<0.01
#10	chr16:58497395	-154	28/46	0.60	2/46	0.04	<0.01
#11	chr16:58497402	-147	29/46	0.62	1/46	0.03	<0.01
#12	chr16:58497406	-143	29/46	0.62	1/46	0.03	<0.01
#13	chr16:58497414	-135	28/46	0.61	1/46	0.02	<0.01
#14	chr16:58497418	-131	30/46	0.65	1/46	0.02	<0.01
#15	chr16:58497425	-124	29/46	0.64	2/46	0.04	<0.01
#16	chr16:58497438	-111	29/46	0.64	1/46	0.02	<0.01
#17	chr16:58497440	-109	29/46	0.62	2/46	0.04	<0.01
#18	chr16:58497449	-100	29/46	0.64	1/46	0.03	<0.01
#19	chr16:58497458	-91	29/46	0.62	1/46	0.02	<0.01
#20	chr16:58497460	-89	27/46	0.59	1/46	0.03	<0.01
#21	chr16:58497462	-87	29/46	0.62	1/46	0.02	<0.01
#22	chr16:58497470	-79	29/46	0.64	1/46	0.02	<0.01
#23	chr16:58497472	-77	30/46	0.65	2/46	0.04	<0.01
#24	chr16:58497477	-72	30/46	0.65	2/46	0.04	<0.01
#25	chr16:58497481	-68	28/46	0.60	2/46	0.04	<0.01
#26	chr16:58497499	-50	29/46	0.63	1/46	0.02	<0.01
#27	chr16:58497505	-44	28/46	0.61	1/46	0.02	<0.01
#28	chr16:58497510	-39	27/46	0.59	1/46	0.02	<0.01
#29	chr16:58497512	-37	30/46	0.65	2/46	0.04	<0.01
#30	chr16:58497517	-32	29/46	0.62	2/46	0.04	<0.01
#31	chr16:58497521	-28	29/46	0.64	1/46	0.02	<0.01
#32	chr16:58497531	-18	29/46	0.62	2/46	0.04	<0.01
#33	chr16:58497542	-7	29/46	0.64	1/46	0.02	<0.01
#34	chr16:58497544	-5	29/46	0.63	1/46	0.02	<0.01
#35	chr16:58497550	2	29/46	0.62	1/46	0.02	<0.01
#36	chr16:58497556	8	28/46	0.61	1/46	0.03	<0.01

#37	chr16:58497558	10	30/46	0.65	1/46	0.03	<0.01
#38	chr16:58497561	13	29/46	0.63	1/46	0.02	<0.01
#39	chr16:58497566	18	30/46	0.65	1/46	0.02	<0.01

EXAMPLE 2

Comparison of differential methylation CpG sites of BMP3 and NDRG4 genes related to CRC among different races.

[00184] We analyzed the methylation microarray data of BMP3 and NDRG4 genes in TCGA database (Illumina Human Methylation 450 data) and found that there are five significantly different methylation CpG sites in the promoter region of NDRG4 genes between the white and Asian populations (Figure 2 and Table 9). In order to further verify the differential methylation CpG sites, tissue and blood samples were collected from 106 CRC and AA patients in China. The DNAs were extracted and treated with bisulfate. The promoter regions of BMP3 and NDRG4 genes were amplified and sequenced. We analyzed the sequencing data and found indeed that hypermethylation CpG sites of Asian population are different from that of the white in TCGA database, and the hypermethylation CpG sites are different between CRC and AA tissue samples. According to these different methylation CpG sites, we developed a detection kit specifically for CRC and AA screening in the Asian population (e.g. Chinese population).

Table 9 the difference in methylation CpG sites of BMP3 and NDRG4 gene in the white and Asian populations.

NO.	Gene	Methylation Site	One-way Anova	Coordinate (hg19)	Gene region	Methylation percentage(White>Asian)
1	BMP3	cg19675063	p = 0.001735, (f = 4.408)	81975265	Exon3	Y
2	NDRG4	cg27147718	p = 0.00004290, (f = 6.558)	58496542	promoter	Y
3	NDRG4	cg04190807	p = 0.01761, (f = 3.036)	58497230	promoter	N
4	NDRG4	cg00687686	p = 0.01657, (f = 3.072)	58497236	promoter	N
5	NDRG4	cg04942472	p = 0.01955, (f = 2.973)	58497239	promoter	N
6	NDRG4	cg01466678	p = 0.02327, (f = 2.867)	58497395	promoter	N
7	NDRG4	cg05469759	p = 0.03612, (f = 2.599)	58498456	Intron1	N
8	NDRG4	cg00262031	p = 0.02585, (f = 2.803)	58498574	Intron1	N
9	NDRG4	cg06650115	p = 0.01221, (f = 3.256)	58498585	Intron1	N
10	NDRG4	cg04005075	p = 0.01938, (f = 2.978)	58498636	Intron1	N
11	NDRG4	cg09324514	p = 0.0004509, (f = 5.194)	58498710	Intron1	N
12	NDRG4	cg16812519	p = 0.01862, (f = 3.002)	58498754	Intron1	N
13	NDRG4	cg26824423	p = 0.007544, (f = 3.543)	58498818	Intron1	N
14	NDRG4	cg05333442	p = 0.003665, (f = 3.970)	58533743	Intron3	N
15	NDRG4	cg01343363	p = 0.01099, (f = 3.319)	58533808	Intron3	Y
16	NDRG4	cg27113419	p = 0.006302, (f = 3.649)	58533979	Intron3	Y
17	NDRG4	cg05725404	p = 0.01145, (f = 3.294)	58534157	Intron3	N

EXAMPLE 3**Screening of primers and probes for BMP3 and NDRG4 genes****(1) Designing and selection of primers and probes of BMP3 and NDRG4 genes.**

[00185] qPCR Primers and probes were designed based on the methylated CpG sites of BMP3 and NDRG4 genes. Three pairs of preferred primers and probes were identified. The preferred primers and probes are compared with several other candidate primes and probes of BMP3 and NDRG4 genes with positive controls and negative controls.. The information of primers and probes is shown in Table 10.

Table 10 the information of preferred and remaining primers and probes

Gene	Group	Primer/Probe	Primer/Probe ID	Sequence (5' to 3')
BMP3	Preferred 1	Forward primer	SEQ ID NO.:3	TTTGAAAATATTCGGGTATATACGTCGC
BMP3	Preferred 1	Reverse primer	SEQ ID NO.:4	ATAAACTCTTCCCAACAACCTACGCGAA
BMP3	Preferred 1	Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTTCG
NDRG4	Preferred 1	Forward primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTTAGGTTTC
NDRG4	Preferred 1	Reverse primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
NDRG4	Preferred 1	Probe	SEQ ID NO.:8	CGTCGCGGTTTTCTGTTCTGTTTTCTGTTCTG
BMP3	Preferred 2	Forward primer	SEQ ID NO.:9	AATATTCGGGTATATACGTCGCGA
BMP3	Preferred 2	Reverse primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCAA
BMP3	Preferred 2	Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTCGTAGCGT
NDRG4	Preferred 2	Forward primer	SEQ ID NO.:12	GCGGGTGAGAAGTCGGC
NDRG4	Preferred 2	Reverse primer	SEQ ID NO.:13	GTAACCTCCGCCTTCTACGC
NDRG4	Preferred 2	Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTCTGTT
BMP3	Preferred 3	Forward primer	SEQ ID NO.:15	AATATTCGGGTATATACGTCGCGATT
BMP3	Preferred 3	Reverse primer	SEQ ID NO.:16	ACTTACTACGCTAACCCAACG
BMP3	Preferred 3	Probe	SEQ ID NO.:17	TAGCGTTGGAGTGGAGACGGCGTTCGTA
NDRG4	Preferred 3	Forward primer	SEQ ID NO.:18	CGGTTTTCTGTTCTGTTTTTTCG
NDRG4	Preferred 3	Reverse primer	SEQ ID NO.:19	AACCTAAACTAATCCCGAACGAACC
NDRG4	Preferred 3	Probe	SEQ ID NO.:20	TCGTTTATCGGGTATTTTAGTCGCGTAG
BMP3	Control 1	Forward primer	SEQ ID NO.:21	GAGTGGAGACGGCGTTCGTA
BMP3	Control 1	Reverse primer	SEQ ID NO.:22	CCACTTACTACGCTAACCCAACG
BMP3	Control 1	Probe	SEQ ID NO.:23	CGGGTGAGGTTTCGCGTAGTTGTTGGG
NDRG4	Control 1, Control 2, Control 3	Forward primer	SEQ ID NO.:24	GGGTGTTTTTTAGGTTTCGCGTC
NDRG4	Control 1	Reverse primer	SEQ ID NO.:25	CGTAACTCCGCCTTCTACGC
NDRG4	Control 1, Control 2, Control 3	Probe	SEQ ID NO.:26	ACGCGACTAAAATACCCGATAAACGAACGAAA AAACGAAC
BMP3	Control 2	Forward primer	SEQ ID NO.:27	TTAGGTTGCGTTGGGTTAGCG

Gene	Group	Primer/Probe	Primer/Probe ID	Sequence (5' to 3')
BMP3	Control 2	Reverse primer	SEQ ID NO.:28	ACTCCGAAAACGCAAAAAACCG
BMP3	Control 2	Probe	SEQ ID NO.:29	ATTCGGTCGCGTTTCGGGTTTCGTGC
NDRG4	Control 2	Reverse primer	SEQ ID NO.:30	GACCCGCGAAACGATACCG
BMP3	Control 3	Forward primer	SEQ ID NO.:31	TATTCGGGTTATATACGTCCG
BMP3	Control 3	Reverse primer	SEQ ID NO.:32	CTTACTACGCTAACCCAACG
BMP3	Control 3	Probe	SEQ ID NO.:33	CCCAACAACCTACGCGAACCTCACCCG
NDRG4	Control 3	Reverse primer	SEQ ID NO.:34	TCGCGCGTAACTTCCGCCTT

(2) Comparison preferred primers and probes with controls with positive and negative controls of BMP3 and NDRG4 genes.

[00186] Standard samples with different methylation ratios were formed by spiking positive control DNAs into negative control DNAs of BMP3 and NDRG4 genes respectively in different ratio (Table 11). Analytical sensitivity was compared by amplify the standard sample DNA with preferred and control primers and probes of BMP3 and NDRG4 genes respectively. Analytical specificity was compared by amplify negative control DNA of BMP3 and NDRG4 genes with preferred and control primers and probes of BMP3 and NDRG4 genes respectively, and the quantity of negative control DNAs is 10^4 copies, 10^5 copies, 10^6 copies, 10^7 copies and 10^8 copies per reaction.

Table 11 standard sample of BMP3 and NDRG4 genes

Methylation ratio	Copy number of methylated DNA	Copy number of total DNA
1	10^5	10^5
1/10	10^4	10^5
1/10 ²	10^3	10^5
1/10 ³	10^2	10^5
1/10 ⁴	10	10^5

[00187] The master mixture was prepared as Table 5 and quantitative PCR reaction condition is denaturation at 95°C for 1 minutes firstly, and then 50 cycles of denaturation at 95°C for 20 seconds and elongation at 60°C for 1 minutes.

Table 12 the final concentration of reagents of methylation qPCR

No.	Reagents	Final concentration
1	10×PCR buffer	1×
2	MgCl ₂	2mM
3	dNTPs	0.2mM

4	Taq DNA polymerase	2U
5	Each primer	0.75mM
6	Each probe	0.25 mM
7	DNA template	2 μ L
8	Ultrapure water	Up to 50 μ L

[00188] As shown in Table 13 and Figure 3A to Figure 3L, $1/10^4$ methylated DNA can be detected using the three pairs of preferred primers and probes, but three pairs of controls cannot be detected.

[00189] As shown in Table 13 and Figure 3A to Figure 3L, the preferred three pairs of BMP3 and NDRG4 primers and probes of the present invention are capable of stably detecting methylation levels as low as one in ten thousand, while three pairs of comparative primers and probes are incapable of detecting one in ten thousand methylation levels.

Table 13 analytical sensitivity comparison of preferred and control primers and probes of BMP3 and NDRG4 genes.

Gene	Methylated ratio	Copy number of methylated DNA	Copy number of total DNA	Preferred primers and probes			Control primers and probes		
				Preferred 1	Preferred 2	Preferred 3	Control 1	Control 2	Control 3
BMP3	1	10^5	10^5	Y	Y	Y	Y	Y	Y
	1/10	10^4	10^5	Y	Y	Y	Y	Y	Y
	1/10 ²	10^3	10^5	Y	Y	Y	Y	Y	Y
	1/10 ³	10^2	10^5	Y	Y	Y	Y	N	N
	1/10 ⁴	10	10^5	Y	Y	Y	N	N	N
NDRG4	1	10^5	10^5	Y	Y	Y	Y	Y	Y
	1/10	10^4	10^5	Y	Y	Y	Y	Y	Y
	1/10 ²	10^3	10^5	Y	Y	Y	Y	Y	Y
	1/10 ³	10^2	10^5	Y	Y	Y	N	N	Y
	1/10 ⁴	10	10^5	Y	Y	Y	N	N	N

Y: detected, N: Undetermined

[00190] Analytical sensitivity amplification curves of BMP3 and NDRG4 are shown in Figure 3A to Figure 3L for each of the preferred combination, compared to control groups.

Primers and probes	BMP3	NDRG4
Preferred 1	Figure 3A	Figure 3B
Preferred 2	Figure 3C	Figure 3D
Preferred 3	Figure 3E	Figure 3F
Control 1	Figure 3G	Figure 3H
Control 2	Figure 3I	Figure 3J

Control 3	Figure 3K	Figure 3L
-----------	-----------	-----------

[00191] As shown in Table 14 and Figure 4A to Figure 4L, the three pairs of preferred primers and probes for BMP3 and NDRG4 genes have no amplification signals for different concentrations of unmethylated DNA, while the comparative primers and probes exhibit different degrees of non-specific amplification.

Table 14 Comparison of analytical specificity between preferred and control primers and probes.

BMP3 (Copies per reaction)	NDRG4 (Copies per reaction)	Preferred primers and probes			Control primers and probes		
		Preferred 1	Preferred 2	Preferred 3	Control 1	Control 2	Control 3
10 ⁸	10 ⁸	N	N	N	Y	Y	Y
10 ⁷	10 ⁷	N	N	N	N	Y	Y
10 ⁶	10 ⁶	N	N	N	N	N	N
10 ⁵	10 ⁵	N	N	N	N	N	N
10 ⁴	10 ⁴	N	N	N	N	N	N

Y: no amplification, N: non-specific amplification

Analytical specificity amplification curves of BMP3 and NDRG4

Primers and probes	BMP3	NDRG4
Preferred 1	Figure 4A	Figure 4B
Preferred 2	Figure 4C	Figure 4D
Preferred 3	Figure 4E	Figure 4F
Control 1	Figure 4G	Figure 4H
Control 2	Figure 4I	Figure 4I
Control 3	Figure 4K	Figure 4L

Example 4

Validation of preferred and comparative methylated primers and probes is performed with fecal samples of BMP3 and NDRG4 genes.

[00192] Methylation level of BMP3 and NDRG4 genes in 81 fecal samples was detected with the three pairs of preferred and comparative primers and probes. The three comparative primers and probes are:

Comparative 1: BMP3 forward primer SEQ ID NO.:21, BMP3 reverse primer SEQ ID NO.: 22, and BMP3 probe SEQ ID NO.: 23; NDRG4 forward primer SEQ ID NO.: 24, NDRG4 reverse primer SEQ ID NO.: 25, and NDRG4 probe SEQ ID NO.: 26;

Comparative 2: BMP3 forward primer SEQ ID NO.:27, BMP3 reverse primer SEQ ID NO.: 28, and BMP3 probe SEQ ID NO.: 29; NDRG4 forward primer SEQ ID NO.: 24, NDRG4 reverse primer SEQ ID NO.: 30, and NDRG4 probe SEQ ID NO.: 26;

Comparative 3: BMP3 forward primer SEQ ID NO.:31, BMP3 reverse primer SEQ ID NO.: 32, and BMP3 probe SEQ ID NO.: 33; NDRG4 forward primer SEQ ID NO.: 24, NDRG4 reverse primer SEQ ID NO.: 34, and NDRG4 probe SEQ ID NO.: 26.

[00193] The information of fecal samples is showed in Table 15.

Table 15 the statistics of eighty-one fecal samples

Colonoscopy	Numbers
Normal	46
adenomas	15
Colorectal cancer	20
total	81

[00194] Fecal DNAs were extracted from samples by following the method described below.

(1) Extraction of fecal DNA

- i. Add 40 mL of lysate to 4-6g fecal samples, and vortex thoroughly, and then incubate at 50°C for 16 hours.
- ii. Then centrifuge at 5000rpm for 10 minutes. Pay attention to the weighing balance before centrifugation. After the end of the centrifugation, carefully remove the centrifuge tubes and do not shake it vigorously.
- iii. Pipette 9 mL of supernatant into a new 50 mL centrifuge tube, and then add 1 mL of extracting adjuvant, 60 μ L of magnetic beads and 10 mL of isopropanol. Vortex for 10 seconds, and incubate at 65°C for 20 minutes. Mix up and down once per 5 minutes during incubating.
- iv. After the incubation, put the 50 mL centrifuge tube on a magnetic stand and keep steady for 3 minutes until the supernatant is clear and discard the supernatant.
- v. Remove the 50 mL centrifuge tube from magnetic stand and add 12 mL of washing solution. Vortex until the magnetic beads fall off the tube well, and keep steady for 3 minutes, and put the 50 mL centrifuge tube back to magnetic stand for 3 minutes until the supernatant is clear and discarded.
- vi. Add 15 mL of 80% ethanol solution, vortex until beads fall off the tube well and keep steady for 3 minutes. Put the tube back into magnetic stand and keep steady for 3 minutes until the supernatant is clear and discarded.
- vii. Repeat the previous step once.

- viii. Pipette the bottom residual liquid, keep the tube open, and incubate at 65°C for 5 minutes. Remove the tube from the magnetic stand until the beads is dry and add 1.5 mL of preheated eluent I . Pipette the beads from the well to eluent I with 1000 μ L pipette, and transfer the mixture to a 2mL centrifuge tube and close the lid, and then incubate at 65°C for 5 minutes.
- ix. Centrifuge at 13000 rpm for 3 minutes, pipette 600 μ L of supernatant to a new 1.5 mL of tube, and then add 600 μ L of binding solution and vortex thoroughly.
- x. Transfer 600 μ L of the above mixture to a DNA purification column, centrifuge at 13,000 rpm for 1 minutes, and discard the waste.
- xi. Deal with the remaining mixture as previous step, and centrifuge at 13,000 for 2 minutes.
- xii. Add 600 μ L of 90% ethanol solution to the DNA purification column, centrifuge at 13,000 rpm for 1 minute and discard the waste.
- xiii. Repeat the previous step 2 times.
- xiv. Centrifuge at 13,000 rpm for 3 minutes, and put the DNA purification column into a new 1.5 mL centrifuge tube. Open the DNA purification column and incubate at 65°C for 5 minutes to dry.
- xv. Drop 100 μ L of preheated eluate II to the middle of the DNA purification column, close the lid, and then incubate at 65°C for 5 minutes. Centrifuge at 13,000 rpm for 2 minutes. Obtain the eluted DNA solution and store it at 2-8°C for use. Long-term storage should be kept at -25°C to -15°C.

[00195]

(2) Bisulfite treatment

The detailed operation steps according to example 3.

[00196]

(3) qPCR

QPCR master mixture was prepared as follows:

Table 16

NO.	Reagent	Final concentration
1	10 $\times$ PCR buffer	1 $\times$
2	MgCl ₂	2 mM
3	dNTP s	0.2 mM

4	Taq DNA Polymerase	2 U
5	Each primer	0.75 μ M
6	Each probe	0.25 μ M
7	DNA template	2 μ L
8	Ultrapure water	Up to 50 μ L

qPCR reaction condition is one cycle of denaturing at 95°C for 2 minutes, and 50 cycles of denaturing at 95°C for 20 seconds and elongating at 95°C for 1 minutes. B2M gene is as a reference gene for quality control of qPCR reaction.

(4) Results

[00197] As shown in Table 17, the sensitivity of CRC and AA detection in 81 fecal samples with the three pairs of preferred primers and probes of BMP3 and NDRG4 genes are up to 85.0%-95.0% (CRC by BMP3 methylation), 66.7%-73.3% (AA by BMP3 methylation), and 90.0%-95.0% (CRC by NDRG4 methylation) and 73.3%-86.7% (AA by NDRG4 methylation), respectively. In addition, the specificity of CRC and AA detection using either BMP3 methylation data or NDRG4 methylation data is both about 97.8%-100.0%.

[00198] However, the sensitivity of CRC and AA detection in eighty-one fecal samples with the three pairs of comparative primers and probes of BMP3 and NDRG4 genes are 85.0%-90.0% (CRC by BMP3 methylation), 46.7%-60.0% (AA by BMP3 methylation), and 90.0%-95.0% (CRC by NDRG4 methylation) and 66.7%-73.3% respectively (AA by NDRG4 methylation). Also, the overall specificity of CRC and AA detection using either BMP3 methylation data or NDRG4 methylation data is up about 91.3%-93.5% and 93.5%-95.7%, respectively.

[00199] It can be seen that the preferred primers and probes are superior to the comparative primers and probes in the detection of the clinical samples, especially for AA detection.

Table 17 the result statistics of clinical samples detection

BMP3										
Item	Preferred 1			Preferred 2			Preferred 3			
Samples grouping	CRC	AA	Normal	CRC	AA	Normal	CRC	AA	Normal	
Number of positive samples detected	19/20	11/15	0/46	18/20	11/15	1/46	17/20	10/15	1/46	
Sensitivity	95.0%	73.3%	/	90.0%	73.3%	/	85.0%	66.7%	/	
Specificity	100.0%			97.8%			97.8%			
Total compliance rate	93.8%			91.4%			88.9%			
Amplification curve	Figure 5A			Figure 5B			Figure 5C			
Item	Comparative 1			Comparative 2			Comparative 3			
Samples grouping	CRC	AA	Normal	CRC	AA	Normal	CRC	AA	Normal	
Number of positive samples detected	18/20	9/15	3/46	17/20	8/15	4/46	17/20	7/15	4/46	
Sensitivity	90.0%	60.0%	/	85.0%	53.3%	/	85.0%	46.7%	/	
Specificity	93.5%			91.3%			91.3%			
Total compliance rate	86.4%			82.7%			81.5%			
Amplification curve	Figure 5D			Figure 5E			Figure 5F			

NDRG4										
Item	Preferred 1			Preferred 2			Preferred 3			
Samples grouping	CRC	AA	Normal	CRC	AA	Normal	CRC	AA	Normal	
Number of positive samples detected	19/20	13/15	0/46	19/20	12/15	0/46	18/20	11/15	1/46	
Sensitivity	95.0%	86.7%	/	95.0%	80.0%	/	90.0%	73.3%	/	
Specificity	100.0%			100.0%			97.8%			
Total compliance rate	96.3%			95.1%			91.4%			
Amplification curve	Figure 6A			Figure 6B			Figure 6C			
Item	Comparative 1			Comparative 2			Comparative 3			
Samples grouping	CRC	AA	Normal	CRC	AA	Normal	CRC	AA	Normal	
Number of positive samples detected	18/20	10/15	2/46	18/20	10/15	3/46	19/20	11/15	3/46	
Sensitivity	90.0%	66.7%	/	90.0%	66.7%	/	95.0%	73.3%	/	
Specificity	95.7%			93.5%			93.5%			
Total compliance rate	88.9%			87.7%			90.1%			
Amplification curve	Figure 6D			Figure 6E			Figure 6F			

Table 18 detection results of 81 fecal samples

No.	Gender	Age	Colonoscopy	Ct BMP3						NDRG4					
				Preferred primers and probes			Comparative primers and probes			Preferred primers and probes			Comparative primers and probes		
				1	2	3	1	2	3	1	2	3	1	2	3
1	Man	50	CRC-stage I	33.9	32.1	37.5	32.7	35.8	36.2	U	33.6	35	33.9	38.8	40.1
2	Man	51	CRC-stage II	34.9	32.2	29.1	41.3	37.6	39.9	39.4	38.4	33.8	34	38.1	40.3
3	Female	56	CRC-stage III	35.1	33.1	31.4	34.9	41.2	35.1	30.7	30.4	33.2	34.2	42.6	35
4	Man	72	CRC-stage I	35.7	U	U	33.4	40.1	36	37.2	38.9	U	U	U	36.7
5	Female	75	CRC-stage II	31.9	36.6	32.7	37.9	36.3	34.5	35.4	37.4	34.5	36.4	33	40.5
6	Female	55	CRC-stage II	37.5	35.1	29.5	43	39.8	39.7	34.4	35.9	35.6	35.1	37.4	42
7	Female	50	CRC-stage II	36.2	32.2	38.8	39.7	37.2	37.3	33.6	37.4	38.8	42.3	36.4	37.5
8	Man	62	CRC-stage IV	U	U	U	U	U	U	36.1	34.5	39.9	41	41.8	40.4
9	Female	74	CRC-stage I	35.7	37.7	29.3	38	36	44	29.6	32.3	36.6	42.6	36.8	33
10	Female	73	CRC-stage II	31.7	29.7	33.9	32.6	42	32.1	36.5	35.7	39.7	39.3	33.5	31.8
11	Man	52	CRC-stage II	34	40	35.9	39.5	38.1	42.4	36.9	U	36.5	42.6	37.8	36.6
12	Female	57	CRC-stage II	34.5	29.1	29.3	34.6	42.6	41.3	37.7	31.8	35.6	42.3	39	42.6
13	Man	51	CRC-stage I	37.6	34	33.2	38.8	U	U	31.9	35.4	35.9	39.1	35.1	34.6
14	Man	60	CRC-stage II	36.2	31.7	35.6	40.1	42.2	33.1	29.6	31.7	29.1	33.8	36.4	36.5
15	Female	49	CRC-stage III	31.1	37	29.4	41.9	35.8	34.3	34.4	33.6	35.6	40.3	39.2	39
16	Man	49	CRC-stage IV	34	35.5	U	41.8	35.9	37.5	37.6	32.3	32.3	38.1	U	31.3
17	Man	48	CRC-stage I	38.1	37.1	35.8	39	40.7	38.3	29.4	30.7	U	U	37.1	U
18	Man	54	CRC-stage II	36.1	33.7	35.2	U	U	U	29.4	37	33.5	34.7	34	41.5
19	Female	60	CRC-stage III	29.2	36.4	33.6	38.4	36.5	33.7	32.9	30.1	35.1	33.5	38.5	41.3
20	Female	62	CRC-stage IV	36.4	34.7	38.5	39.5	37.9	38.9	36.4	29.4	33.1	35.8	43.1	42.8
21	Female	50	AA	38.7	U	U	35.3	U	U	34.6	37	U	U	U	40.9
22	Female	56	AA	U	35.6	30	40.9	38.6	35	34.6	U	29.3	42.7	33.4	U
23	Man	58	AA	U	U	29.6	39.3	38	41.1	39.3	35.1	39.5	39.4	U	37.4

No.	Gender	Age	Colonoscopy	CI BMP3						NDRG4					
				Preferred primers and probes			Comparative primers and probes			Preferred primers and probes			Comparative primers and probes		
				1	2	3	1	2	3	1	2	3	1	2	3
24	Man	67	AA	36.3	32.8	34.2	U	U	U	33.9	33	32.7	36	35.5	U
25	Man	51	AA	35.9	33.7	33.9	36.4	38	39	32.7	33.9	U	U	35.8	40.4
26	Female	71	AA	29.3	30	35.4	40.3	37.7	34.4	30.7	34.9	38.3	33.8	41.1	33
27	Man	65	AA	U	30.3	U	U	U	U	34.8	34	36.8	41.9	42.7	36.1
28	Man	62	AA	37.5	29.2	U	U	U	U	U	U	36.8	43.1	43.9	32.8
29	Female	45	AA	31	37.6	37.4	37.2	38.5	37.8	30.8	31.7	30.8	U	U	U
30	Man	67	AA	29.7	34.3	29.8	31.3	U	U	33.9	U	U	U	U	U
31	Female	53	AA	36.3	U	U	U	35.4	34.6	29.5	34.9	37.7	39.6	40	42.9
32	Female	68	AA	38.4	37.3	30.4	U	U	U	U	37.2	U	U	U	35.1
33	Female	50	AA	36.9	39	32.1	36.6	U	U	39.5	29.1	30.7	37.1	37.6	36.6
34	Female	48	AA	U	U	U	U	37.3	35.8	38.3	32.2	32.4	36.2	36	43.9
35	Man	50	AA	35.5	29.1	32.7	39.4	42.7	U	29.8	33.8	32.9	36.4	38.2	42
36	Man	45	Normal	U	U	U	U	U	U	U	U	U	U	U	U
37	Female	50	Normal	U	U	U	U	U	U	U	U	U	U	U	U
38	Man	59	Normal	U	U	U	U	U	U	U	U	U	U	U	U
39	Man	73	Normal	U	U	U	U	U	U	U	U	U	U	U	U
40	Man	47	Normal	U	U	U	U	U	U	U	U	U	U	U	U
41	Man	71	Normal	U	U	U	U	U	U	U	U	U	U	U	U
42	Female	45	Normal	U	U	U	U	U	U	U	U	U	U	U	U
43	Female	62	Normal	U	U	U	U	U	U	U	U	U	U	U	U
44	Man	55	Normal	U	U	U	U	U	U	U	U	U	U	U	U
45	Female	54	Normal	U	U	U	U	U	U	U	U	U	U	U	U
46	Man	63	Normal	U	U	U	39.4	U	U	U	U	U	U	U	U
47	Female	53	Normal	U	U	U	U	U	U	U	U	U	41.2	U	U

No.	Gender	Age	Colonoscopy	CI BMP3						NDRG4					
				Preferred primers and probes			Comparative primers and probes			Preferred primers and probes			Comparative primers and probes		
				1	2	3	1	2	3	1	2	3	1	2	3
48	Man	67	Normal	U	U	U	U	37.2	41.2	U	U	U	U	U	U
49	Man	69	Normal	U	U	U	U	U	U	U	U	U	U	U	U
50	Man	54	Normal	U	U	U	U	U	U	U	U	U	U	U	U
51	Man	50	Normal	U	U	U	U	U	U	U	U	U	U	U	U
52	Female	65	Normal	U	U	U	U	U	U	U	U	U	40.3	41	U
53	Female	56	Normal	U	U	U	U	U	U	U	U	U	U	U	U
54	Man	67	Normal	U	35.6	U	U	41.3	U	U	U	U	U	U	U
55	Man	70	Normal	U	U	U	U	U	U	U	U	U	U	U	U
56	Female	75	Normal	U	U	U	U	U	U	U	U	U	U	U	U
57	Female	50	Normal	U	U	U	U	U	U	U	U	U	U	U	U
58	Female	66	Normal	U	U	U	U	U	U	U	U	U	U	U	U
59	Man	49	Normal	U	U	U	U	U	41.8	U	U	U	U	U	U
60	Man	71	Normal	U	U	U	40.8	U	U	U	U	U	U	U	U
61	Man	70	Normal	U	U	U	U	U	U	U	U	U	40.5	U	U
62	Man	73	Normal	U	U	U	U	U	U	U	U	U	U	U	U
63	Man	48	Normal	U	U	U	U	38.5	U	U	U	U	U	U	U
64	Man	52	Normal	U	U	U	U	U	U	U	U	U	U	U	39.8
65	Man	71	Normal	U	U	U	U	U	U	U	U	U	U	U	U
66	Female	60	Normal	U	U	U	U	U	41.9	U	U	35.9	38.4	U	U
67	Female	65	Normal	U	U	U	U	U	U	U	U	U	U	U	U
68	Female	74	Normal	U	U	U	U	U	U	U	U	U	U	U	U
69	Man	71	Normal	U	U	U	U	U	U	U	U	U	U	U	U
70	Female	73	Normal	U	U	U	U	U	U	U	U	U	U	U	U
71	Man	49	Normal	U	U	U	U	U	U	U	U	U	U	U	U

No.	Gender	Age	Colonoscopy	CI BMP3						NDRG4					
				Preferred primers and probes			Comparative primers and probes			Preferred primers and probes			Comparative primers and probes		
				1	2	3	1	2	3	1	2	3	1	2	3
72	Female	69	Normal	U	U	U	U	U	U	U	U	U	U	U	U
73	Man	59	Normal	U	U	U	42.7	41.9	42.5	U	U	U	U	U	U
74	Man	71	Normal	U	U	37.1	U	U	U	U	U	U	U	U	U
75	Female	54	Normal	U	U	U	U	U	U	U	U	U	U	38.9	35.2
76	Man	71	Normal	U	U	U	U	U	U	U	U	U	U	U	U
77	Female	63	Normal	U	U	U	U	U	U	U	U	U	U	U	U
78	Man	74	Normal	U	U	U	U	U	U	U	U	U	U	U	U
79	Man	62	Normal	U	U	U	U	U	U	U	U	U	U	U	U
80	Female	75	Normal	U	U	U	U	U	U	U	U	U	U	U	U
81	Man	48	Normal	U	U	U	U	U	U	U	U	U	U	U	U

Notes U: undetermined

Example 5**Screening of methylation primer and probe combinations of BMP3 and NDRG4 genes****(1) Fecal DNA extraction**

[00200] Fecal DNA is extracted according to the protocol of example 3.

(2) Methylation detection of BMP3 and NDRG4 genes

[00201] The combinations of primer and probe of BMP3 and NDRG4 genes is shown as follow. qPCR was performed with the nine combinations respectively according to the steps of example 3.

Table 19. Combinations of primer and probe

Combination No.	BMP3	NDRG4
1	Preferred 3	Preferred 1
2	Preferred 3	Preferred 2
3	Preferred 3	Preferred 3
4	Preferred 1	Preferred 1
5	Preferred 1	Preferred 2
6	Preferred 1	Preferred 3
7	Preferred 2	Preferred 1
8	Preferred 2	Preferred 2
9	Preferred 2	Preferred 3

The sequences of primers and probes of nine combinations of BMP3 and NDRG4 genes are as follows:

Table 20. Sequences of primers and probes in the combinations

Combination No.	Gene	Group	Primer/Probe	Sequence ID	Sequence (5' to 3')
1	BMP3	Preferred 3	Forward primer	SEQ ID NO.:15	AATATTCGGGTTATATACGTCGCGATT
			Reverse primer	SEQ ID NO.:16	ACTTACTACGCTAACCCAACG
			Probe	SEQ ID NO.:17	TAGCGTTGGAGTGGAGACGGCGTTCGTA
	NDRG4	Preferred 1	Forward primer	SEQ ID NO.:6	ATCGATCGGGTGTTTTTAGGTTTC
			Reverse primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGGTTTTTCGTTTCGTTTTTCGTTTCGT
2	BMP3	Preferred 3	Forward primer	SEQ ID NO.:15	AATATTCGGGTTATATACGTCGCGATT
			Reverse primer	SEQ ID NO.:16	ACTTACTACGCTAACCCAACG
			Probe	SEQ ID NO.:17	TAGCGTTGGAGTGGAGACGGCGTTCGTA
	NDRG4	Preferred 2	Forward primer	SEQ ID NO.:12	GCGGGTGAGAAGTCGGC
			Reverse primer	SEQ ID NO.:13	GTAACCTCCGCCTTCTACGC
			Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTTCGTT
3	BMP3	Preferred 3	Forward primer	SEQ ID NO.:15	AATATTCGGGTTATATACGTCGCGATT
			Reverse primer	SEQ ID NO.:16	ACTTACTACGCTAACCCAACG

Combination No.	Gene	Group	Primer/Probe	Sequence ID	Sequence (5' to 3')
	NDRG4	Preferred 3	Probe	SEQ ID NO.:17	TAGCGTTGGAGTGGAGACGGCGTTCGTA
			Forward primer	SEQ ID NO.:18	CGGTTTCGTTTCGTTTTTCG
			Reverse primer	SEQ ID NO.:19	AACCTAAATAATCCCGAACGAACC
			Probe	SEQ ID NO.:20	TCGTTTATCGGGTATTTTAGTCGCGTAG
4	BMP3	Preferred 1	Forward primer	SEQ ID NO.:3	TTTAAAAATATTCGGGTATATACGTCGC
			Reverse primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTCG
	NDRG4	Preferred 1	Forward primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
			Reverse primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGGTTTTTCGTTTCGTTTTTCGTTTCGT
5	BMP3	Preferred 1	Forward primer	SEQ ID NO.:3	TTTAAAAATATTCGGGTATATACGTCGC
			Reverse primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTCG
	NDRG4	Preferred 2	Forward primer	SEQ ID NO.:12	GCGGGTGAGAAATCGGC
			Reverse primer	SEQ ID NO.:13	GTAACCTCCGCCTTCTACGC
			Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTTCGTT
6	BMP3	Preferred 1	Forward primer	SEQ ID NO.:3	TTTAAAAATATTCGGGTATATACGTCGC
			Reverse primer	SEQ ID NO.:4	ATAAACTCTTCCCCAACAACTACGCGAA
			Probe	SEQ ID NO.:5	AGCGTTGGAGTGGAGACGGCGTTCG
	NDRG4	Preferred 3	Forward primer	SEQ ID NO.:18	CGGTTTCGTTTCGTTTTTCG
			Reverse primer	SEQ ID NO.:19	AACCTAAATAATCCCGAACGAACC
			Probe	SEQ ID NO.:20	TCGTTTATCGGGTATTTTAGTCGCGTAG
7	BMP3	Preferred 2	Forward primer	SEQ ID NO.:9	AATATTCGGGTATATACGTCGCGA
			Reverse primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCCAA
			Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTCGTAGCGT
	NDRG4	Preferred 1	Forward primer	SEQ ID NO.:6	ATCGATCGGGGTGTTTTTAGGTTTC
			Reverse primer	SEQ ID NO.:7	CCTTCTACGCGACTAAAATACCCGAT
			Probe	SEQ ID NO.:8	CGTCGCGGTTTTTCGTTTCGTTTTTCGTTTCGT
8	BMP3	Preferred 2	Forward primer	SEQ ID NO.:9	AATATTCGGGTATATACGTCGCGA
			Reverse primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCCAA
			Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTCGTAGCGT
	NDRG4	Preferred 2	Forward primer	SEQ ID NO.:12	GCGGGTGAGAAATCGGC
			Reverse primer	SEQ ID NO.:13	GTAACCTCCGCCTTCTACGC
			Probe	SEQ ID NO.:14	TAGGTTTCGCGTCGCGGTTTTTCGTT
9	BMP3	Preferred 2	Forward primer	SEQ ID NO.:9	AATATTCGGGTATATACGTCGCGA
			Reverse primer	SEQ ID NO.:10	GCAACCTAACAAATAAACTCTTCCCCAA
			Probe	SEQ ID NO.:11	TGGAGTGGAGACGGCGTTCGTAGCGT
	NDRG4	Preferred 3	Forward primer	SEQ ID NO.:18	CGGTTTCGTTTCGTTTTTCG
			Reverse primer	SEQ ID NO.:19	AACCTAAATAATCCCGAACGAACC
			Probe	SEQ ID NO.:20	TCGTTTATCGGGTATTTTAGTCGCGTAG

[00202] QPCR master mixture used in the reaction is as follows:

Table 21. QPCR master mixture for BMP3/NDRG4

NO.	Reagents	Final concentration
1	10× PCR buffer	1×
2	MgCl ₂	2 mM
3	dNTP s	0.2 mM
4	Taq DNA Polymerase	2 U
5	Each primer	0.75 μM
6	Each probe	0.25 μM
7	DNA template	2 μL
8	Ultrapure water	Up to 50 μL

[00203] The qPCR reaction condition is one cycle of denatured at 95°C for 2 minutes, and 50 cycles of denatured at 95°C for 20 seconds and elongated at 60°C for 1 minute.

(3) Variants detection of KRAS gene

[00204] Seven mutation hotspots in codon 12 and 13 of KRAS gene were detected. The seven mutant are G12D, G13D, G12V, G12C, G12S, G12A and G13R and sequences of which primers and probes are in Table 22.

Table 22 Primers and probes to detect seven mutation of KRAS gene

Primer Name	Primer/Probe	Sequence ID	Sequence(5' to 3')
G12D-F	Forward primer	SEQ ID NO.:35	AACTTGTGGTAGTTGGAGGTGA
G13D-F	Forward primer	SEQ ID NO.:36	AACTTGTGGTAGTTGGAGCTGGGGA
G12V-F	Forward primer	SEQ ID NO.:37	AACTTGTGGTAGTTGGAGTTGT
G12C-F	Forward primer	SEQ ID NO.:38	AAACTTGTGGTAGTTGGGGCTT
G12S-F	Forward primer	SEQ ID NO.:39	AAACTTGTGGTAGTTGGTGCTA
G12A-F	Forward primer	SEQ ID NO.:40	AACTTGTGGTAGTTGGAGCAGC
G12R-F	Forward primer	SEQ ID NO.:41	AAACTTGTGGTAGTTGGAGCTC
Kras-R	Reverse primer	SEQ ID NO.:42	GAATGGTCCTGCACCAGTAATATG
ACTB-F	Forward primer	SEQ ID NO.:43	AGGGCTTCTTGTCTTTCCTT
ACTB-R	Reverse primer	SEQ ID NO.:44	CGTGCTCGATGGGGTACTTC
KRAS-P	Probe	SEQ ID NO.:45	AGGCAAGAGTGCCTTGACGATACAGC
ACTB-P	Probe	SEQ ID NO.:46	CGTGATGGTGGGCATGGGTCAGAAGGA

QPCR master mixture used in the reaction is as follows:

Table 23. QPCR master mixture for KRAS

No.	Reagents	Final concentration
1	G12A-F	0.72 μ M
2	G12C-F	0.60 μ M
3	G12D-F	0.72 μ M
4	G12R-F	0.48 μ M
5	G12S-F	0.90 μ M
6	G12V-F	0.72 μ M
7	G13D-F	0.48 μ M
8	Kras-R	0.90 μ M
9	Kras-P	0.10 μ M
10	ACTB-F	0.30 μ M
11	ACTB -R	0.30 μ M
12	ACTB -P	0.05 μ M
13	5 $\times$ PCR Buffer,-Mg ²⁺	1 $\times$
14	MgCl ₂	1.0 mM
15	dNTPs	0.1 mM
16	Taq DNA Polymerase	2.5 U
17	DNA template	2 μ L
18	Ultrapure water	Up to 50 μ L

[00205] QPCR reaction condition is one cycle of denatured at 95°C for 5 minutes and 45 cycles of denatured at 95°C for 15 seconds, annealed at 71°C for 60 seconds and then elongated at 55°C for 50 seconds.

[00206] The reaction quality control of qPCR was carried out using ACTB as a reference gene.

(4) Fecal hemoglobin test

[00207] Fecal hemoglobin was detected with fecal immunochemical Test (FIT), and the result is positive or negative.

(5) Generate a score with a formula

[00208] Bring the Ct value of qPCR detection of BMP3, NDRG4 and KRAS genes, and positive and negative result of fecal hemoglobin test into the logistic regression formula as follow:

$$P = e^K / (1 + e^K)$$

[00209] Where: P is a comprehensive index, $K = a \cdot \Delta Ct1 + b \cdot \Delta Ct2 + c \cdot \Delta Ct3 + d \cdot FIT + X$, e is a natural constant, a, b, c, d, X are clinical constant. $\Delta Ct1$, $\Delta Ct2$, and $\Delta Ct3$ are Ct value of target genes subtract to that of reference genes.

[00210] The test result is positive if P value is equal or more than a predetermined threshold, otherwise is negative. Positive result indicates that the subject is possible to have CRC or AA.

(6) Test result

[00211] The eighty-one fecal samples of example 3 were detected and the results of different combinations of primers and probes of BMP3 and NDRG4 genes are shown in Table 24.

Table 24 the detection result of 81 fecal samples

Combination No.	BMP3	NDRG4	Sensitivity (CRC+AA)	Specificity (CRC+AA)	Number of samples tested positive in CRC	Number of samples tested positive in AA	Number of samples tested negative in normal	Sensitivity (CRC)	Sensitivity (AA)
1	Preferred 3	Preferred 1	94.3%	97.8%	19/20	14/15	45/46	95.0%	93.3%
2	Preferred 3	Preferred 2	91.4%	97.8%	18/20	14/15	45/46	90.0%	93.3%
3	Preferred 3	Preferred 3	80.0%	95.6%	16/20	12/15	44/46	80.0%	80.0%
4	Preferred 1	Preferred 1	97.1%	100.0%	20/20	14/15	46/46	100.0%	93.3%
5	Preferred 1	Preferred 2	97.1%	97.8%	20/20	14/15	45/46	100.0%	93.3%
6	Preferred 1	Preferred 3	88.6%	95.6%	19/20	12/15	44/46	95.0%	80.0%
7	Preferred 2	Preferred 1	97.1%	97.8%	19/20	15/15	45/46	95.0%	100.0%
8	Preferred 2	Preferred 2	94.3%	97.8%	19/20	14/15	45/46	95.0%	93.3%
9	Preferred 2	Preferred 3	85.7%	95.6%	19/20	11/15	44/46	95.0%	73.3%

[00212] It can be seen that: (1) Combination Nos. 4, 5 and 7 are superior to other six combinations. Considering all tested primers and probes are preferred ones and are superior to other primers and probes (see Example 3), the three particular combinations are superior to any known combinations of BMP3 and NDRG4 primers/probes. (2) The sensitivity and specificity of the kit for CRC and AA detection comprising BMP3, NDRG4, KRAS genes and fecal hemoglobin detection is significantly better than those of BMP3 or NDRG4 single gene methylation detection. (3) The sensitivity and specificity of the kit for CRC detection in the Asian population (e.g., the Chinese population) is obviously superior to existing similar products, such as Cologuard®. (4) The sensitivity and specificity of the kit for AA detection in the Asian population (e.g., the Chinese population) is significantly better than the existing similar products, such as Cologuard®.

[00213] The disclosures, including the claims, figures and/or drawings, of each and every patent, patent application, and publication cited herein are hereby incorporated herein by reference in their entireties.

[00214] Unless defined otherwise, all technical and scientific terms herein have the same meaning as commonly understood by one of ordinary skill in the art to which this invention belongs. Although any methods and materials, similar or equivalent to those described herein, can be used in the practice or testing of the present invention, the preferred methods and materials are described herein. All publications, patents, and patent publications cited are incorporated by reference herein in their entirety for all purposes.

[00215] The publications discussed herein are provided solely for their disclosure prior to the filing date of the present application. Nothing herein is to be construed as an admission that the present invention is not entitled to antedate such publication by virtue of prior invention.

[00216] While the invention has been described in connection with specific embodiments thereof, it will be understood that it is capable of further modifications and this application is intended to cover any variations, uses, or adaptations of the invention following, in general, the principles of the invention and including such departures from the present disclosure as come within known or customary practice within the art to which the invention pertains and as may be applied to the essential features hereinbefore set forth and as follows in the scope of the appended claims.

CLAIMS

1. A kit for detecting the presence or the absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof, comprising:

- a) a first pair of primers and a first probe for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in a biological sample obtained from the patient, wherein the first pair of primers and first probe, each of which comprises a contiguous sequence of at least 16 nucleotides that is identical to, complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 1,
- b) a second pair of primers and a second probe for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in a biological sample obtained from the patient, wherein the second pair of primers and second probe, each of which comprises a contiguous sequence of at least 16 nucleotides that is identical to, complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 2,

2. The kit of claim 1,

wherein the first pair of primers and the first probe are selected from the group consisting of:

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5;
- ii) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11; and
- iii) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17;

and,

wherein the second first pair of primers and the second probe are selected from the group consisting of:

- iv) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8;
- v) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14; and
- vi) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20;

3. The kit of claim 1, wherein the kit comprises:

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and
- ii) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

4. The kit of claim 1, wherein the kit comprises:

- i) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and
- ii) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

5. The kit of claim 1, wherein the kit comprises:

- i) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and
- ii) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

6. The kit of any one of claims 1 to 5, wherein both the first probe and the second probe comprise a fluorescent donor and an acceptor fluorophore.

7. The kit of claim 6, wherein first probe and the second probe are TAQMAN® probes.

8. The kit of any one of claims 1 to 7, wherein the kit further comprises:

- (1) means for detecting the presence or absence of at least one mutation in the KRAS gene in the patient; and
- (2) means for detecting the presence or absence of hemoglobin in a biological sample obtained from the patient.

9. The kit of claim 8, wherein the means for detecting the presence or absence of at least one mutation in the KRAS gene in the patient comprises at least one pair of primers capable of amplifying the Exon 12 and/or Exon 13 region of the KRAS gene in a polymerase chain reaction (PCR).
10. The kit of claim 8, wherein the means for detecting the presence or absence of hemoglobin in the biological sample comprises an anti-hemoglobin antibody.
11. The kit of claim 9, wherein the primers are capable of amplifying a KRAS gene region comprising at least one KRAS mutation selected from the group consisting G12D, G12V, G12C, G13D, G12A, G12R, G12S, and G13C.
12. The kit of claim 10, wherein the antibody is a colloidal gold-conjugated antibody.
13. The kit of any one of claims 1 to 12, wherein the kit further comprises means for amplifying an a reference gene for quantification.
14. The kit of any one of claims 1 to 12, wherein the kit further comprises instructions for use and/or interpretation of a test result obtained by using the kit.
15. The kit of claim 10, wherein the kit further comprises means to detect a complex formed by the antibody and the hemoglobin in the biological sample.
16. The kit of any one of claims 1 to 15, wherein the biological sample obtained from the patient is a fecal sample.
17. The kit of any one of claims 1 to 16, wherein the kit further comprises a bisulfite reagent, and a container suitable for mixing the bisulfite reagent and the biological sample of the patient, or polynucleotides obtained from the biological sample.
18. The kit of any one of claims 1 to 17, wherein the kit further comprises a methylation sensitive restriction enzyme reagent.
19. The kit of any one of claims 1 to 18, wherein the kit further comprises

- (1) a positive standard and a negative standard for detecting BMP3 methylation in the biological sample, and
- (2) a positive standard and a negative standard for detecting NDRG4 methylation in the biological sample.

20. The kit of claim 19, wherein

the positive standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTCGGGTGTTTTTAAAAATAAAGCGAGGAGGGAAGGTATAGATAGATTTTG
AAAATATTCGGGTATATACGTCGCGATTTATAGTTTTTTTTTAGCGTTGGAGTGGAGACG
GCGTTTCGTAGCGTTTTGCGCGGGTGAGGTTTCGCGTAGTTGTTGGGGAAGAGTTTATTTGT
TAGGTTGCGTTGGGTAGCGTAGTAAGTGGGGTTGGTCGTTATTTTCGTTGTATTTCGGTCGC
GTTTCGGGTTTTTCGTGCGTTTTTCGTTTTAG (SEQ ID NO: 67);

the negative standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTTGGGTGTTTTTAAAAATAAAGTGAGGAGGGAAGGTATAGATAGATTTTG
AAAATATTTGGGTATATATGTTGTGATTTATAGTTTTTTTTTAGTGTTGGAGTGGAGATGG
TGTTTGTAGTGTTTTGTGTGGGTGAGGTTTGTGTAGTTGTTGGGGAAGAGTTTATTTGTTA
GGTTGTGTTGGGTAGTGTAGTAAGTGGGGTTGGTTGTTATTTTGTGTTATTGGTTGTGT
TTTGGGTTTTGTGTGTTTTTGTTTTAG (SEQ ID NO: 68);

the positive standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTCGGCGGGGGCGCGGATCGATCGGGGTGTTTTTTAGGTTTCGCGTCGCGGTTT
TCGTTTCGTTTTTTCGTTTCGTTTATCGGGTATTTTAGTCGCGTAGAAGGCGGAAGTTACGCG
CGAGGGATCGCGGTTTCGTTTCGGGATTAGTTTTAGGTTTCGGTATCGTTTCGCGGGTCGAGC
GTTTATATTCGTTAAATTTACGCGGGTACGTTTTTCGCGGCGTATCGTTTTTAGTT (SEQ ID
NO.: 69); and

the negative standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTTGGTGGGGGTGTGGATTGATTGGGGTGTTTTTTAGGTTTTGTGTTGTGGTTT
TTGTTTGTTTTTTTGTTTGTATTGGGTATTTTAGTTGTGTAGAAGGTGGAAGTTATGTGT
GAGGGATTGTGGTTTGTGTTGGGATTAGTTTTAGGTTTGGTATTGTTTTGTGGGTGAGTGT
TTATATTTGTTAAATTTATGTGGGTATGTTTTTGTGGTGTATTGTTTTTAGTT (SEQ ID NO.:
70).

21. The kit of claim 13, the means for amplifying an internal control gene comprises comprising primers for amplifying a positive control gene and/or a negative control gene.

22. A method for detecting the presence or absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof, comprising:

- a) obtaining genomic DNA from a biological sample of the patient;
- b) treating the genomic DNA of a), or a fragment thereof, with one or more reagents to convert cytosine bases that are unmethylated thereof to uracil or another base that is detectably dissimilar to cytosine in terms of hybridization properties;
- c) contacting the treated genomic DNA, or the treated fragment thereof, with a first pair of primers for detecting the presence or absence of methylation sites of a gene encoding bone morphogenetic protein 3 (BMP3) in the patient, and a second pair of primers for detecting the presence or absence of methylation sites of a gene encoding NDRG family member 4 protein (NDRG4) in the patient, wherein the first pair of primers comprise a contiguous sequence of at least 9 nucleotides that is identical to, complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 1, and wherein the second pair of primers comprise a contiguous sequence of at least 9 nucleotides that is complementary to, or hybridizes under stringent hybridization conditions to SEQ ID NO.: 2, wherein the treated genomic DNA or the fragment thereof is either amplified to produce at least one amplicate by the first pair of primers or the second pair of primers, or is not amplified; and
- d) determining the presence or absence of CRC or AA in the patient, based on a presence or absence of said amplicate, the methylation state or level of at least one CpG dinucleotide of the BMP3 gene and the NDRG4 gene in the patient.

23. The method of claim 22, wherein the first pair of primers and the first probe are selected from the group consisting of:

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5;
 - ii) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11; and
 - iii) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17;
- and,

wherein the second first pair of primers and the second probe are selected from the group consisting of:

- iv) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8;

v) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14; and

vi) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20;

24. The method of claim 22, wherein the method comprises using

i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

25. The method of claim 22, wherein the method comprises using

i) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

26. The method of claim 22, wherein the method comprises using

i) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the biological sample obtained from the patient, and

ii) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the biological sample obtained from the patient.

27. The method of any one of claims 22 to 26, wherein both the first probe and the second probe comprise a fluorescent donor and an acceptor fluorophore.

28. The method of any one of claims 22 to 26, wherein first probe and the second probe are TAQMAN® probes.

29. The method of any one of claims 22 to 28, wherein the method further comprises a step of detecting the presence or absence of at least one mutation in the KRAS gene in a biological sample obtained from the patient, and a step of detecting the presence or absence of hemoglobin in a biological sample obtained from the patient.

30. The method of any one of claim 29, wherein the step of detecting the presence or absence of at least one mutation in the KRAS gene in the patient comprises using at least one pair of primers capable of amplifying the Exon 12 and/or Exon 13 region of the KRAS gene in a polymerase chain reaction (PCR).

31. The method of any one of claim 29, wherein the step of detecting the presence or absence of hemoglobin in the biological sample comprises using an anti-hemoglobin antibody.

32. The method of claim 30, wherein the primers are capable of amplifying a KRAS gene region comprising at least one KRAS mutation selected from the group consisting G12D, G12V, G12C, G13D, G12A, G12R, G12S, and G13C.

33. The method of claim 31, wherein the antibody is a colloidal gold-conjugated antibody.

34. The method of any one of claims 22 to 33, wherein the amplification of BMP3 gene is performed in a quantitative PCR (qPCR), and the method further comprises amplifying a first reference gene to determining the Ct value of the BMP3 amplification as $\Delta Ct1$.

35. The method of any one of claims 22 to 33, wherein the amplification of NDRG4 gene is performed in a quantitative PCR (qPCR), and the method further comprises amplifying a second reference gene to determining the Ct value of the NDRG4 amplification as $\Delta Ct2$.

36. The method of any one of claims 30 to 33, wherein the amplification of mutant KRAS gene is performed in a quantitative PCR (qPCR), and the method further comprises amplifying a third reference gene to determining the Ct value of the mutant KRAS amplification as $\Delta Ct3$.

37. The method of claim 34 or 35, wherein the first and the second reference genes are the same.

38. The method of claim 37, wherein the same reference gene is a B2M gene.

39. The method of claim 36, wherein the mutant KRAS gene comprises a mutation selected from the group consisting of G12D, G13D, G12V, G12C, G12S, G12A, and G13R.

40. The method of claim 39, wherein the mutant KRAS gene is amplified by one or more pairs of primers selected from the group consisting of:

(1) a forward primer G12D-F comprising SEQ ID NO.: 35, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(2) a forward primer G13D-F comprising SEQ ID NO.: 36, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(3) a forward primer G12V-F comprising SEQ ID NO.: 37, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(4) a forward primer G12C-F comprising SEQ ID NO.: 38, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(5) a forward primer G12S-F comprising SEQ ID NO.: 39, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(6) a forward primer G12A-F comprising SEQ ID NO.: 40, and a reverse primer Kras-R comprising SEQ ID NO.: 42; and

(7) a forward primer G12R-F comprising SEQ ID NO.: 41, and a reverse primer Kras-R comprising SEQ ID NO.: 42,

and wherein the KRAS probe for the qPCR comprises SEQ ID NO.: 46.

41. The method of claim 36, wherein the third reference gene is an ACTB gene.

42. The method of claim 41, wherein qPCR primers for amplifying ACTB gene comprise SEQ ID NOs.: 43 and 44, and the probe comprise SEQ ID NO.: 46.

43. The method of claims 22 to 42, wherein the method comprises using

(1) a positive standard and a negative standard for detecting BMP3 methylation in the sample, and

(2) a positive standard and a negative standard for detecting NDRG4 methylation in the sample.

44. The method of claim 43, wherein

the positive standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTCGGGTGTTTTTAAAAATAAAGCGAGGAGGGAAGGTATAGATAGATTTTG
 AAAATATTCGGGTATATACGTCGCGATTTATAGTTTTTTTTTAGCGTTGGAGTGGAGACG
 GCGTTCGTAGCGTTTTGCGCGGGTGAGGTTTCGCGTAGTTGTTGGGGAAGAGTTTATTTGT
 TAGGTTGCGTTGGGTAGCGTAGTAAGTGGGGTTGGTCGTTATTTTCGTTGTATTCGGTCGC
 GTTTCGGGTTTCGTGCGTTTTTCGTTTTAG (SEQ ID NO: 67);

the negative standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTTGGGTGTTTTTAAAAATAAAGTGAGGAGGGAAGGTATAGATAGATTTTG
 AAAATATTTGGGTATATATGTTGTGATTATAGTTTTTTTTTAGTGTTGGAGTGGAGATGG
 TGTTTGTAGTGTTTTGTGTGGGTGAGGTTTGTGTAGTTGTTGGGGAAGAGTTTATTTGTTA
 GGTTGTGTTGGGTAGTGTAAGTGGGGTTGGTTGTTATTTGTTGTATTTGGTTGTGT
 TTTGGGTTTTGTGTGTTTTTGTTTTAG (SEQ ID NO: 68);

the positive standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTCGGCGGGGGCGCGGATCGATCGGGGTGTTTTTTAGGTTTCGCGTCGCGGTTT
 TCGTTCGTTTTTTTCGTTTCGTTTATCGGGTATTTTAGTCGCGTAGAAGGCGGAAGTTACGCG
 CGAGGGATCGCGGTTTCGTTTCGGGATTAGTTTTAGGTTTCGGTATCGTTTCGCGGGTCGAGC
 GTTTATATTCGTTAAATTTACGCGGGTACGTTTTTCGCGGCGTATCGTTTTTAGTT (SEQ ID
 NO.: 69); and

the negative standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTTGGTGGGGGTGTGGATTGATTGGGGTGTTTTTTAGGTTTGTGTTGTGGTTT
 TTGTTTGTTTTTTTGTTTGTATTGGGTATTTTAGTTGTGTAGAAGGTGGAAGTTATGTGT
 GAGGGATTGTGGTTTGTGTTGGGATTAGTTTTAGGTTTGGTATTGTTTTGTGGGTGAGTGT
 TTATATTTGTTAAATTTATGTGGGTATGTTTTGTGGTGTATTGTTTTTAGTT (SEQ ID NO.:
 70).

45. The method of any one of claims 22 to 44, wherein the method comprises amplifying a quality control standard.

46. A method for detecting the presence or absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof, comprising using a kit of any one of claims 1 to 21.

47. A method for detecting the presence or absence of colorectal cancer (CRC) or advanced adenoma (AA) in a patient in need thereof, comprising:

a) obtaining an untreated genomic DNA from a fecal sample of the patient;

- b) treating the genomic DNA of a), or a fragment thereof, with one or more reagents to convert cytosine bases that are unmethylated thereof to uracil or another base that is detectably dissimilar to cytosine in terms of hybridization properties;
- c) performing a quantitative PCR (qPCR) using the treated genomic DNA of b) as a template, and determining the Ct value of BMP3 gene in the patient as $\Delta Ct1$;
- d) performing a qPCR using the treated genomic DNA of b) as a template, and determining the Ct value of NDRG4 gene in the patient as $\Delta Ct2$;
- e) performing a qPCR using the untreated genomic DNA as a template, and determining the Ct value of a mutant KRAS gene in the patient as $\Delta Ct3$;
- f) performing a fecal immunochemical test of hemoglobin protein in the fecal sample and determining a score as FIT;
- g) determining the value of K, wherein $K = a * \Delta Ct1 + b * \Delta Ct2 + c * \Delta Ct3 + d * FIT + X$, wherein a, b, c, d, X are clinical constants; and
- h) determining the value of a comprehensive index P, wherein $P = e^K / (1 + e^K)$, wherein e is the natural constant, wherein when P is equal or more than a predetermined threshold, the patient is determined to have CRC and/or AA, and when P is less than the threshold, the patient is determined to be health.

48. The method of claim 47, wherein the qPCR for amplifying BMP3 gene comprises a first pair of primers and a first probe, wherein the first pair of primers and the first probe are selected from the group consisting of:

- i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5;
- ii) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11; and
- iii) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17;

and,

wherein the qPCR for amplifying NDRG4 gene comprises a second pair of primers and a second probe, wherein the second pair of primers and the second probe are selected from the group consisting of:

- iv) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8;

v) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14; and

vi) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20;

49. The method of claim 47, wherein the method comprises using

i) a forward primer comprising SEQ ID NO.: 3, a reverse primer comprising SEQ ID NO.: 4, and a probe comprising SEQ ID NO.: 5, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the sample, and

ii) a forward primer comprising SEQ ID NO.: 6, a reverse primer comprising SEQ ID NO.: 7, and a probe comprising SEQ ID NO.: 8, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the sample.

50. The method of claim 47, wherein the method comprises using

i) a forward primer comprising SEQ ID NO.: 9, a reverse primer comprising SEQ ID NO.: 10, and a probe comprising SEQ ID NO.: 11, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the sample, and

ii) a forward primer comprising SEQ ID NO.: 12, a reverse primer comprising SEQ ID NO.: 13, and a probe comprising SEQ ID NO.: 14, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the sample.

51. The method of claim 47, wherein the method comprises using

i) a forward primer comprising SEQ ID NO.: 15, a reverse primer comprising SEQ ID NO.: 16, and a probe comprising SEQ ID NO.: 17, for detecting the methylation state or level of at least one CpG dinucleotide of the BMP3 gene in the sample, and

ii) a forward primer comprising SEQ ID NO.: 18, a reverse primer comprising SEQ ID NO.: 19, and a probe comprising SEQ ID NO.: 20, for detecting the methylation state or level of at least one CpG dinucleotide of the NDRG4 gene in the sample.

52. The method of any one of claims 47 to 51, wherein both the first probe and the second probe comprise a fluorescent donor and an acceptor fluorophore.

53. The method of any one of claims 47 to 52, wherein first probe and the second probe are TAQMAN® probes.

54. The method of any one of claims 47 to 53, wherein the mutant KRAS gene comprising at least one KRAS mutation selected from the group consisting G12D, G12V, G12C, G13D, G12A, G12R, G12S, and G13C.

55. The method of any one of claims 47 to 54, wherein the fecal immunochemical test comprises a colloidal gold-conjugated antibody.

56. the method of any one of claims 47 to 55, wherein step c) and step d) comprises using B2M gene as a reference gene.

57. The method of claim 54, wherein the mutant KRAS gene is amplified by one or more pairs of primers selected from the group consisting of:

(1) a forward primer G12D-F comprising SEQ ID NO.: 35, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(2) a forward primer G13D-F comprising SEQ ID NO.: 36, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(3) a forward primer G12V-F comprising SEQ ID NO.: 37, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(4) a forward primer G12C-F comprising SEQ ID NO.: 38, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(5) a forward primer G12S-F comprising SEQ ID NO.: 39, and a reverse primer Kras-R comprising SEQ ID NO.: 42;

(6) a forward primer G12A-F comprising SEQ ID NO.: 40, and a reverse primer Kras-R comprising SEQ ID NO.: 42; and

(7) a forward primer G12R-F comprising SEQ ID NO.: 41, and a reverse primer Kras-R comprising SEQ ID NO.: 42,

and wherein the KRAS probe for the qPCR comprises SEQ ID NO.: 46.

58. The method of claim 57, wherein ACTB gene is used as a reference gene in the qPCR for amplifying the mutant KRAS gene.

59. The method of claim 58, wherein the qPCR primers for amplifying ACTB gene comprise SEQ ID NOs.: 43 and 44, and the qPCR probe for ACTB gene comprises SEQ ID NO.: 46.

60. The method of claims 47 to 59, wherein the method comprises using

- (1) a positive standard and a negative standard for detecting BMP3 methylation in the sample, and
- (2) a positive standard and a negative standard for detecting NDRG4 methylation in the sample.

61. The method of claim 60, wherein

the positive standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTCGGGTGTTTTTAAAAATAAAGCGAGGAGGGAAGGTATAGATAGATTTTG
AAAATATTCGGGTATATACGTCGCGATTTATAGTTTTTTTTTAGCGTTGGAGTGGAGACG
GCGTTCGTAGCGTTTTGCGCGGGTGAGGTTTCGCGTAGTTGTTGGGGAAGAGTTTATTTGT
TAGGTTGCGTTGGGTTAGCGTAGTAAGTGGGGTTGGTCGTTATTTTCGTTGTATTCGGTCGC
GTTTCGGGTTTCGTGCGTTTTTCGTTTTAG (SEQ ID NO: 67);

the negative standard for detecting BMP3 methylation comprises a polynucleotide sequence of

GTTAGTTTGGTTGGGTGTTTTTAAAAATAAAGTGAGGAGGGAAGGTATAGATAGATTTTG
AAAATATTTGGGTATATATGTTGTGATTTATAGTTTTTTTTTAGTGTTGGAGTGGAGATGG
TGTTTGTAGTGTTTTGTGTGGGTGAGGTTTGTGTAGTTGTTGGGGAAGAGTTTATTTGT
GGTTGTGTTGGGTTAGTGTAGTAAGTGGGGTTGGTTGTTATTTTGTGTTATTTGGTTGTGT
TTTGGGTTTTGTGTGTTTTTGTTTTAG (SEQ ID NO: 68);

the positive standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTCGGCGGGGGCGCGGATCGATCGGGGTGTTTTTTAGGTTTCGCGTCGCGGTTT
TCGTTTCGTTTTTTCGTTTCGTTTATCGGGTATTTTAGTCGCGTAGAAGGCGGAAGTTACGCG
CGAGGGATCGCGGTTTCGTTTCGGGATTAGTTTTAGGTTTCGGTATCGTTTCGCGGGTTCGAGC
GTTTATATTCGTTAAATTTACGCGGGTACGTTTTTCGCGGCGTATCGTTTTTAGTT (SEQ ID
NO.: 69); and

the negative standard for detecting NDRG4 methylation comprises a polynucleotide sequence of

TGAGAAGTTGGTGGGGGTGTGGATTGATTGGGGTGTTTTTTAGGTTTTGTGTTGTGGTTT
TTGTTTGTTTTTTTGTGTTTATTGGGTATTTTAGTTGTGTAGAAGGTGGAAGTTATGTGT
GAGGGATTGTGGTTTGTGTTGGGATTAGTTTTAGGTTTGGTATTGTTTTGTGGGTTGAGTGT
TTATATTTGTTAAATTTATGTGGGTATGTTTTTGTGGTGTATTGTTTTTAGTT (SEQ ID NO.:
70).

62. The method of any one of claims 47 to 61, wherein the method comprises amplifying a quality control standard in the step c) and the step d).

63. A method for diagnosing and treating a colorectal cancer (CRC) and/or advanced adenoma (AA) in a patient in need thereof, comprising determining the presence or absence of CRC and/or AA in the patient by using a kit of any one of claims 1 to 21, and treating the patient depends on the presence or absence of CRC and/or AA in the patient.

64. A method for diagnosing and treating a colorectal cancer (CRC) and/or advanced adenoma (AA) in a patient in need thereof, comprising determining the presence or absence of CRC and/or AA in the patient by using a method of any one of claims 22 to 62, and treating the patient depends on the presence or absence of CRC and/or AA in the patient.