

COMBINATION RESPIRATORY THERAPY DEVICE, SYSTEM AND METHOD

BACKGROUND

[0001] Patients with neuromuscular weakness as a result of strokes, spinal cord injuries, head trauma, or diseases such as muscular dystrophy and amyotrophic lateral sclerosis (ALS, or Lou Gehrig's disease) have an increased risk of morbidity and mortality due to a weak cough and shallow breathing (hypoventilation). The number of chronic illnesses that cause weak cough and impaired lung ventilation is large and expanding.

[0002] When a patient has an ineffective cough, chest secretions are retained in the respiratory system, causing pneumonia, lung collapse, or, where the mucus fills the windpipe, fatal respiratory arrest. Additionally, shallow breathing causes low oxygen levels and high carbon dioxide levels in the patient's bloodstream, resulting in a medically fragile state of chronic respiratory failure in which even a common cold can result in severe respiratory illness. For these reasons, pulmonary complications are viewed as a major cause of morbidity and death in patients that have neuromuscular weakness.

[0003] As a patient's condition worsens, it becomes more likely that the patient will need both cough assistance and assisted ventilation. Respiratory therapies for addressing a weak cough often involve devices that provide assisted coughing via mechanical insufflation/exsufflation, while shallow breathing is commonly addressed by a separate, mechanical ventilation device.

SUMMARY

[0004] In accordance with one embodiment, there is provided a combination respiratory therapy device that includes a blower for providing negatively pressurized air to a mouthpiece coupled to a patient's airway, and an air pulse generator configured for delivering air pulses to at least one of a garment worn by the patient or a nasal interface worn by the patient. The device further includes a network interface in communication with an associated computer network, and a controller including a processor in communication with memory storing instructions which are executable by the processor to execute a combined respiratory therapy prescription, the combined respiratory therapy prescription defining a plurality of different therapy sessions to be performed by the

combination respiratory therapy device over a period of time, each of the plurality of different therapy sessions comprising a mucus extraction therapy. The combination therapy device further includes a display in communication with the controller and configured to display a graphical user interface associated with at least one operation of the combination respiratory therapy device.

[0005] According to another embodiment, there is provided a system for problem-first device control of at least one combination respiratory therapy device. The system includes at least one combination respiratory therapy device and at least one clinician device. The at least one combination respiratory therapy device includes a network interface in communication with an associated computer network, and a controller including a processor in communication with memory storing instructions which are executable by the processor to execute a combined respiratory therapy prescription, the combined respiratory therapy prescription defining a plurality of different therapy sessions to be performed by the combination respiratory therapy device over a period of time. The at least one combination respiratory therapy device further includes a display in communication with the controller and configured to display a graphical user interface associated with at least one operation of the combination respiratory therapy device. The at least one clinician computing device is in communication with the at least one combination respiratory therapy device via the associated computer network and is configured to control the at least one operation of the at least one combined respiratory therapy device.

[0006] In another embodiment, there is provided a method for remotely controlling at least one combination respiratory therapy device by a clinician device. The method includes receiving, at the clinician device in data communication with the at least one combination respiratory therapy device, patient data representative of at least one physiological parameter associated with a patient, and generating on an associated display of the clinician device, a graphical representation of the received patient data. The method further includes receiving, via the associated display, selection data corresponding to a selected therapy adjustment, and communicating the selected therapy adjustment to the at least one combination therapy device. The clinician device includes

a processor in communication with memory storing instructions which are executed by the processor, causing the processor to perform the method.

[0007] In yet another embodiment, there is provided a clinician device for remotely controlling at least one combination respiratory therapy device by a clinician device. The clinician device includes a processor in communication with memory, a network interface in communication with the processor and configured to communicate with at least one combination respiratory therapy device via an associated computer network, and a display in communication with the processor and configured to display a graphical user interface associated with at least one operation of the at least one combination respiratory therapy device. The memory stores instructions which are executed by the processor, causing the processor to receive, via the associated network, patient data representative of at least one physiological parameter associated with a patient of the at least one combination respiratory therapy device, and generate, on the display, a graphical representation of the received patient data. The memory further stores instructions to receive, via the display, selection data corresponding to a selected therapy adjustment, and communicate the selected therapy adjustment to the at least one combination respiratory therapy device via the associated computer network.

BRIEF DESCRIPTION OF THE DRAWINGS

[0008] This disclosure is illustrated by way of example and not by way of limitation in the accompanying figures. The figures may, alone or in combination, illustrate one or more embodiments of the disclosure. For simplicity and clarity of illustration, elements illustrated in the figures are not necessarily drawn to scale. For example, the dimensions of some elements may be exaggerated relative to other elements for clarity. Further, where considered appropriate, reference labels may be repeated among the figures to indicate corresponding or analogous elements.

[0009] FIG. 1 is a simplified module diagram of at least one embodiment of a combination respiratory therapy management system;

[0010] FIG. 2 is a simplified flow diagram of at least one embodiment of a method for configuring a combination respiratory therapy device;

[0011] FIG. 3 is a simplified schematic diagram for at least one embodiment of a combination respiratory therapy management system;

[0012] FIG. 4 is a simplified elevational view of at least one embodiment of a user interface for a prescription creator module of at least one embodiment of the system of FIG. 1;

[0013] FIG. 5 is a simplified elevational view of at least one embodiment of another user interface for a patient interface module of at least one embodiment of the system of FIG. 1;

[0014] FIG. 6 is a simplified flow diagram of at least one embodiment of a method for configuring respiratory therapy using the system of FIG. 1;

[0015] FIG. 7 is a simplified flow diagram of at least one embodiment of another method for configuring respiratory therapy using the system of FIG. 1;

[0016] FIG. 8 is a simplified flow diagram of at least one embodiment of another method for configuring respiratory therapy using the system of FIG. 1;

[0017] FIG. 9 is a simplified flow diagram of at least one embodiment of another method for configuring respiratory therapy using the system of FIG. 1;

[0018] FIG. 10 is a simplified flow diagram of at least one embodiment of a method for performing respiratory therapy using the system of FIGS. 1; and

[0019] FIG. 11 is a simplified block diagram of an exemplary computing environment in connection with which the system of FIG. 1 may be implemented.

[0020] FIG. 12 is a graphical user interface representation of a main screen for operating the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0021] FIG. 13 is a graphical user interface representation of an emergency screen of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0022] FIG. 14 is a graphical user interface representation of an intensify therapy screen of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0023] FIG. 15 is a graphical user interface representation of a customize your therapy screen of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0024] FIG. 16A is a graphical user interface representation of a cough settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0025] FIG. 16B is a graphical user interface representation of a cough inspiratory pressure settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0026] FIG. 16C is a graphical user interface representation of a cough expiratory (suction) settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0027] FIG. 16D is a graphical user interface representation of a cough cycle duration adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0028] FIG. 16E is a graphical user interface representation of a cough sensitivity settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0029] FIG. 16F is a graphical user interface representation of a ventilation settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0030] FIG. 16G is a graphical user interface representation of a ventilation inspiratory pressure settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0031] FIG. 16H is a graphical user interface representation of a ventilation sensitivity settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0032] FIG. 16I is a graphical user interface representation of a duration of ventilation after coughs settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0033] FIG. 16J is a graphical user interface representation of an initial oscillation (mucus mobilization) settings adjustment screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0034] FIG. 16K is a graphical user interface representation of an oscillation (mucus mobilization) adjustment screen following FIG. 16J for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0035] FIG. 16L is a graphical user interface representation of setup completion screen for customizing a therapy of the combination respiratory therapy device in accordance with one embodiment of the subject application.

[0036] FIG. 17 is a flowchart illustrating an example methodology for remotely accessing and controlling a combination respiratory therapy device in accordance with one embodiment of the subject application.

DETAILED DESCRIPTION

[0037] While the concepts of the present disclosure are susceptible to various modifications and alternative forms, specific embodiments thereof are shown by way of example in the drawings and are described in detail below. It should be understood, however, that there is no intent to limit the concepts of the present disclosure to the particular forms disclosed, but on the contrary, the intention is to cover all modifications, equivalents, and alternatives consistent with the present disclosure and the appended claims.

[0038] Referring to FIG. 1, a system **100** for managing combined respiratory therapy provided to a person using a combination respiratory therapy device **110** includes a number of different pieces of computerized functionality, which for ease of discussion are represented herein as modules. Illustratively, these modules include a configurable user interface module **114** and a combined respiratory therapy device control module **120**. In various embodiments of the system **100**, the configurable user interface module **114** may include one or more other modules **112**, **122**, **130**, **132**, **134**, **140**, **142**, **144**. Each of the modules **114**, **120**, **112**, **122**, **130**, **132**, **134**, **140**, **142**, **144** may be implemented as

computer software, firmware, hardware, or a combination thereof, according to the requirements of a particular design or implementation of the system **100**. Further, any or all of the modules **114**, **120**, **112**, **122**, **130**, **132**, **134**, **140**, **142**, **144** may be implemented as part of the combination respiratory therapy device **110** in some embodiments (e.g., as a "standalone" system), while in other embodiments, some of the modules **114**, **120**, **112**, **122**, **130**, **132**, **134**, **140**, **142**, **144** may be implemented on one or more other computing devices (e.g., mobile devices or more traditional desktop- or laptop-style computers). As described in more detail below, the configurable user interface modules **112**, **122**, **130**, **132**, **134**, **140**, **142**, **144** can be selectively activated at the device **110** or at one or more other computing devices to facilitate interaction with and control of the combination device **110** by an authorized person or by different types of authorized persons. Such authorized persons may include clinicians (e.g., physicians), respiratory therapists, nurses, family members and other caregivers, as well as the patient receiving the combined respiratory therapy. For ease of discussion, any of such persons may be referred to herein as a "user" or "users" of the system **100**.

[0039] The prescription creator module **112** enables a user (e.g., a clinician) to create a combined respiratory therapy prescription in a user-friendly way using, for example, graphical, audio, and/or video features. As such, the prescription creator module **112** alleviates the need for users to spend time learning how to use non-intuitive buttons, dials, or the like on a piece of therapy equipment that often has limited real estate for more sophisticated user interface features. A combined respiratory therapy prescription may include one or multiple integrated combinations of different respiratory therapies that are scheduled to occur at various times of the day. For example, in some embodiments, a combined respiratory therapy prescription integrates repeated cycles of cough assistance (e.g., mucus mobilization and/or mucus extraction therapy) with lung ventilation therapy and/or lung volume recruitment therapy. The combined respiratory therapy prescription can be translated into machine-readable instructions by the prescription translator module **122**. Such instructions can be read and executed by hardware components (e.g., a microprocessor) of the combination device **110** to control the provision of combined respiratory therapy to the patient. Once created, the combined respiratory therapy prescription and/or the machine-readable version thereof may be stored in computer

memory, for example, in a computerized data structure such as a combined respiratory therapy prescription database **118**.

[0040] The stored combined respiratory therapy prescription or portions thereof can be accessed and viewed by the person who created the prescription or by other users (e.g., respiratory therapists, other clinicians, caregivers, or the patient) using the data sharing module **136**. As such, the data sharing module **136** presents the combined respiratory therapy prescription in a format that is easy to understand and customizable by or for the particular user. In some embodiments, the data sharing module **136** may allow certain categories of authorized persons to view certain data or make certain types of changes to the combined respiratory therapy prescription. For example, a clinician may be permitted to view and change any aspect of the combined respiratory therapy prescription, while a therapist may be permitted to view the prescription but only change certain parts of the prescription (e.g., cough inspiratory or cough expiratory (suction) pressure within a limited defined range of pressures, or the therapy start time) but not others (e.g., the therapy duration). Similarly, limitations may be placed on a patient or family member's ability to view and change the combined respiratory therapy prescription.

[0041] In some embodiments, the user interface module **114** includes a problem-first device control module **140**, which allows a user to make adjustments to the patient's therapy regimen while the therapy is in progress, whether in response to a change in the patient's condition observed by the user, or in response to some other triggering condition. For example, in some embodiments, the prescription translator module **122** may receive data from the combination device **110** while a therapy is in progress, such as the current status of the therapy or information about the patient's condition. Such data obtained from the device **110** can be displayed to the user in a human-understandable form using the data sharing module **136**. Changes to the patient's therapy may be input by the user in response to the data obtained from the combination device **110**, or in response to changes in the patient's condition observed by the caregiver, for example, using the problem-first device control module **140**.

[0042] Further, in some embodiments, the user interface module **114** includes an audio interface module **142**, and may also include a persona configuration module **144**.

The modules **142**, **144** are designed to further enhance the human-device interaction. For example, the audio interface **142** may play pre-recorded human voice messages explaining how to use the device **110** or explaining a therapy that is about to begin. The persona configuration module **144** may allow the user to ascribe a "personality" to the device **110**. For example, the persona configuration module **144** may allow the user to select a graphical or animated character and/or a particular tone of voice or accent to be used by the device **110** when communicating with the user. Thus, using the system **100**, a combination respiratory therapy device **110** can be used to provide multiple integrated and coordinated respiratory therapies to a respiratory patient over a period of time, in a manner that is both intuitive and non-threatening to the user or the patient. Moreover, the various treatments specified by the respiratory therapy prescription can be customized and adjusted, even down to the breath by breath level, according to the patient's needs as they might change or progress. Such customizations and adjustments can be implemented responsively by the combination respiratory device **110**.

[0043] The illustrative combination respiratory therapy device **110** is one of a family of combination respiratory therapy devices **110**, each of which can be used to provide multiple integrated respiratory therapies. The combination device **110** includes a combined respiratory therapy control module **120**, a positive-pressure air flow patient interface **124**, a negative-pressure air flow patient interface **126**, and, in some embodiments, an air pulse patient interface **128**. The positive-pressure air flow patient interface **124** is designed to supply positive (e.g., inspiratory) pressure to a patient with whom the interface **124** is engaged. The negative-pressure air flow patient interface **126** is designed to supply negative (e.g., expiratory) pressure to a patient with whom the interface **126** is engaged. The air pulse patient interface **128** is designed to supply air pulses to a patient's airway, lungs or chest area to provide, for example, Continuous High Frequency Oscillation (CHFO), Continuous Positive Expiratory Pressure (CPEP), and/or High Frequency Chest Wall Oscillation (HFCWO).

[0044] The illustrative control module **120** interfaces with the prescription creator module **112** to obtain the user-defined and/or user-modified combined respiratory therapy prescription. A version of the prescription translator module **122** may be provided at the user interface level (e.g., as part of the configurable user interface module **114**) and/or at

the device level (e.g., as part of the control module **120**) in various embodiments of the system **100**. Whether at the user interface level or at the device level, the prescription translator module **122** converts the prescription into a machine-executable form that can be used to control the application of air flow to the patient via the interfaces **124**, **126**, **128** as needed. For example, where a prescription specifies a treatment session to include a cough cycle followed by lung ventilation, the prescription translator module **122** generates the device settings needed for the combination device **110** to perform the therapy session. Those device settings may include, for instance, specific air pressure levels, an indication of whether the air pressure is to be positive (e.g., inflation) or negative (e.g., suction), a duration of time to provide the airflow, a number of times to repeat the application of air pressure, etc. As an example, the device settings for a combined respiratory treatment sequence may include "cough: interface **124** on at +25 cm water, interface **126** on at -30 cm water; vent: interface **124** on at +15 cm water, interface **124** on at +4 cm water, time=2 minutes." The translator module **122** may further translate these device settings to specific "valve open" and "valve close" control signals that can be received and acted on directly by specific electromechanical components of the device **110**.

[0045] Using the combination respiratory therapy device **110**, the lung volume recruitment and lung ventilation therapies can be integrated with assisted cough cycles such that assisted ventilation or lung volume recruitment can be automatically coordinated (e.g., alternated) with assisted coughing on a breath to breath basis if needed. In this way, the mix of cough assistance and lung ventilation or lung volume recruitment therapy can be customized to each patient's needs. Further, the combination device **110** eliminates the need to apply two separate devices sequentially. Sequential therapy is difficult for sick or weak patients to tolerate and may cause the patient to become clinically unstable. Additionally, some embodiments of the combination device **110** are portable, such that they can be mounted or stored on a wheelchair, thereby increasing the patient's quality of life. Still further, in some embodiments, the inspiratory and expiratory air circuits of the device **110** are physically separated so that the positive pressure circuit remains clean and unobstructed. Further details of the illustrative combination device **110** are described below in connection with FIG. 3.

[0046] In more detail, the illustrative prescription creator module **112** utilizes a framework of standardized terminology to describe the combined respiratory therapy prescription. As such, the framework provides a vehicle by which combined respiratory therapy prescriptions can be easily created, understood, and shared by the various healthcare practitioners that may be involved in the patient's care. In that regard, the illustrative combined respiratory therapy prescription creator module **112** includes a layering module **130**, a sequencing module **132**, and a patterning module **134**. The layering module **130** allows the caregiver to select an appropriate combination respiratory therapy device **110** simply by specifying (at a graphical user interface of the configurable interface module **114**, for example) the different types or "layers" of therapy that the patient needs.

[0047] The layering module **130** automatically maps the therapy layers selected by the caregiver to one or more combination respiratory therapy devices **110** that are capable of providing those therapies. In the illustrated embodiments, the therapy layers include mucus mobilization, mucus extraction, lung volume recruitment, and lung ventilation therapy layers. Generally speaking, mucus mobilization refers to respiratory therapy that is intended to loosen chest secretions (e.g., mucus) so that the chest secretions may be extracted from the lungs by a normal or assisted cough. Mucus mobilization therapy often involves the mechanical application of air pulses, vibrations, or oscillations to the patient's airway, lungs, chest and/or back by a device such as THE VEST, or the METANEB device, both of which are available from the Hill-Rom Company, Inc.

[0048] Mucus extraction refers to therapy that mechanically assists the patient's natural ability to cough, or which mechanically removes secretions from the lungs for the patient, if the patient is unable to cough on his or her own. To perform mucus extraction therapy, the control module **120** may configure the combination device **110** to begin the therapy by providing an extra-large expanded breath (e.g., "deep lung insufflation"). This may involve the device **110** delivering an inspiratory pressure in the range of about 30%-50% above the "chronic" inspiratory pressure that the device **110** would normally use for ventilation or lung volume recruitment therapy.

[0049] The deep breaths provided by the combination device **110** during mucus extraction therapy are intended to maximize lung recoil (faster flow during exhalation) and

may help expand any collapsed portions of the lung. But because they are big breaths, a patient can normally tolerate mucus extraction therapy only in limited "doses." During mucus extraction therapy, the control module **120** synchronizes the assisted inspiratory breath provided by the device **110** to the patient's inspiratory effort, so that the patient does not exhale during the inspiratory phase. The "deep insufflation" breath is followed by "active" exhalation or "exsufflation," during which suction (negative) pressure removes secretions from the airway. The control module **120** likewise synchronizes the suction (negative) pressure to the patient's expiratory phase of breathing. Accordingly, some forms of mucus extraction therapy that can be provided by the device **110** may be referred to as "mechanical insufflation/exsufflation." For mucus extraction, the inspiratory pressure may be set in the range of about +25 cm water and the expiratory (suction) pressure in the range of about -30 cm water, for example. Thus, a mucus extraction portion of the combined respiratory therapy prescription may be written as "+25/-30."

[0050] Lung ventilation refers to therapy that is intended to mechanically assist the patient with his or her normal breathing pattern, or to mechanically breathe for the patient, if the patient is unable to breathe on his or her own. Thus, lung ventilation therapy is generally applied in a continuous manner for a period of time (rather than in "cycles" like mucus extraction therapy). In the illustrated embodiments, in which they are integrated within a single combination respiratory therapy device **110**, both mucus extraction and lung ventilation therapy include software control algorithms that synchronize the machine-generated breathing or mucus extraction with the patient's natural breathing pattern. Lung ventilation may be achieved by using the positive-pressure air flow patient interface **124** to deliver a higher positive pressure during the patient's inhalation and a lower positive pressure during exhalation. The difference between the higher, inspiratory positive pressure level and the lower, expiratory positive pressure level (the "pressure span") is what ventilates the lungs. For lung ventilation, the inspiratory pressure may be set at +15 cm water and the expiratory pressure at +4 cm water. Thus, a lung ventilation portion of the combined respiratory therapy prescription may be written as "+15/+4" by which the lungs will be ventilated using a pressure span of 11 cm water ($15 - 4 = +11$).

[0051] To better achieve effective lung ventilation, the two positive pressure levels can be synchronized to the patient's breathing pattern. To do this, the combination respiratory

therapy device **110** may detect a small, patient-generated negative "sniff pressure" or negative (inspiratory) flow to sense that the patient is starting inhalation ("wants a breath") and then the device may support that breath with the commanded inspiratory pressure. When the device **110** detects that the patient's inspiratory flow/pressure/effort tapers off, the device **110** can conclude that the patient is ready for exhalation and switch to the lower (expiratory) positive pressure setting (which may be known as the PEEP or positive end expiratory pressure).

[0052] In some embodiments, the control module **120** includes software that synchronizes the mechanically-assisted breathing provided by the device **110** to the patient's spontaneous breathing. For example, in some embodiments, the control module **120** may keep track of the patient's respiratory pattern over time, and use a prior respiratory pattern to predict a future respiratory pattern, in terms of rate of breathing, duration of inhalation, etc. In some embodiments, the flow waveform or "(chest) rise time" generated by the device **110** (how fast the device **110** achieves the set inspiratory pressure) during lung ventilation therapy may be adjustable according to the patient's needs or preferences (a slower rise time is gentler, but too slow may leave the patient breathless). The combination respiratory therapy device **110** may use the same or similar techniques as described above to synchronize mucus extraction (assisted coughing) to the patient's breathing pattern, in order to maximize the efficacy and comfort of mucus extraction therapy or for other reasons.

[0053] Additionally, when the combination device **110** is used for lung ventilation, the control module **120** can instruct the device **110** to provide a "backup" rate of breathing if the patient stops breathing on their own or is sedated, in which case the device **110** will deliver automatic breathing on a timer. Further, when the device **110** is used for lung ventilation, the control module **120** can initiate an alarm to alert caregivers if the interface **124** becomes disconnected from the patient (e.g., tubing is pulled away, nasal mask falls off, etc.). As such, the combination device **110** can selectively provide a number of different features depending on the type of therapy for which it is being used, where some features (such as backup breathing and alarms) may be applicable to some therapies but not others. The combination respiratory therapy device **110** may use the same or similar techniques as described above to automatically provide mucus extraction (assisted

coughing) to a patient who is asleep, unconscious or sedated; for example, by providing mucus extraction automatically, on a timer, while a patient is asleep. Likewise, the control module **120** can initiate an alarm to alert caregivers if the negative airflow interface **126** becomes disconnected from the patient (e.g., falls out of the patient's mouth).

[0054] Lung volume recruitment refers to therapy that mechanically inflates the lungs episodically rather than continuously. When the combination device **110** applies positive pressure air flow for lung volume recruitment, the positive pressure inflates the lungs, and thus prevents them from collapsing. This may be done, for example, after each cough cycle of a daytime therapy session, to help the patient recover his or her breath before the next cough. Lung volume recruitment therapy may also be applied after an assisted cough therapy session has concluded, to "solidify" the lung volume improvements made during the assisted cough session and to help the patient recover from the assisted coughing. Thus, lung volume recruitment therapy may be used when a patient is awake and cooperative.

[0055] As compared to lung ventilation therapy, lung volume recruitment is generally less sophisticated. When performing lung volume recruitment therapy, the combination device **110** may synchronize the inspiratory pressure to the patient's breathing pattern, but there may not be a need for a positive expiratory pressure (PEEP). Since the patient is typically awake and assisting with the process during lung volume recruitment therapy, more advanced software algorithms to track the patient's breathing pattern or to provide a backup rate or alarms are generally not needed, as they would be for lung ventilation. For lung volume recruitment therapy, the inspiratory pressure may be set in the range of about +15 cm water and the expiratory pressure may be set in the range of about 0 cm water (e.g., exhaling to no added pressure). Thus, a lung volume recruitment portion of the combined respiratory therapy prescription may be written as "+15/0."

[0056] The sequencing module **132** allows the caregiver to define, customize, and modify the details of each respiratory therapy treatment session according to the patient's needs, where a "treatment session" generally refers to an instance or occurrence of a coordinated combination of respiratory therapies. A treatment session may include a number of sequentially-executed therapies, or treatment sequences that involve the repetition of one or more types of therapy (such as a treatment sequence made up of

cough assistance followed by lung ventilation therapy). For example, the sequencing module **132** may be used to define or specify the positive and negative pressure settings that comprise each cough cycle; to define or specify the number of cough cycles that will be applied sequentially at the start of each treatment sequence; to define or specify the duration of lung volume recruitment or lung ventilation therapy to be applied substantially immediately after the specified number of cough cycles (e.g., without interruption) to complete each treatment sequence; and/or to define or specify the number of treatment sequences that, applied sequentially, may comprise a treatment session.

[0057] The patterning module **134** allows the caregiver to define an entire pattern of respiratory care to be applied to the patient over a defined period of time (e.g., a 24-hour period), in real clock time. For example, a therapy pattern may include a treatment session comprised of a specific number of treatment sequences, as described above, which starts and ends in the morning; a similar treatment session that starts and ends in the evening, but is immediately preceded by a time period of mucus mobilization therapy; and a lung ventilation therapy session that starts at night, when the patient falls asleep, and ends in the morning, when the patient awakens. As such, the patterning module **134** allows the caregiver to specify that the combination device **110** begins or performs certain therapies at certain times of the day. In other words, the patterning module **134** can be used to associate specific start and stop times with the various coordinated therapies. The patterning module **134** can be used to illustrate a patient's specific layers of therapy and their inter-relationships in clock time. As such, the patterning module **134** can provide a useful tool for creating, modifying and sharing the combination respiratory therapy prescription among interested parties. Thus, in some embodiments, the patterning module **134** presents the pattern of respiratory therapy visually, e.g., as a graphical timeline or other visual representation of the patient's 24-hour respiratory care plan, at a user interface of the system **100**. An illustrative example of one such visual representation is shown in FIG. 4, described below.

[0058] The configurable user interface module **114** includes a software-based user interface to the various modules **112**, **122**, **136**, **140**, **142**, **144**. For example, in embodiments that include the prescription creator module **112**, the user interface module **114** includes a software-based user interface that allows a physician or other qualified

health professional to create the combined respiratory therapy prescription and store the prescription (e.g., in the database **118**) at a computing device of the system **100**. In other embodiments, the user interface may alternatively or in addition, provide the user with access to the data sharing module **136**, the problem-first device control module **140**, and/or other modules of the system **100**.

[0059] As noted above, the data sharing module **136** provides a communication interface by which a combined respiratory therapy prescription, or portions thereof, can be presented at a computing device of the system **100** used to create the combined respiratory therapy prescription to the combination device **110**, or at another computing device, such as a mobile device used by a respiratory therapist responsible for the patient's care, or a computing device located near the patient or used by the patient (such as a personal computer or mobile computing device located in the patient's hospital room or in the patient's home). In some embodiments, the prescription is shared between or among devices using an electronic communication interface accessed by the data sharing module **136**, which includes, for example, one or more input/output modules for data communication via a standard wired or wireless network interface (e.g., WIFI, cellular, Ethernet, etc.), one or more hard-wired communication ports (e.g., a Universal Standard Bus port or other port by which a flash drive or cable may be connected), or a combination thereof.

[0060] Whereas the prescription creation module **112** allows the combined respiratory therapy prescriptions to be created using a "user-friendly" interface and then transferred to the combination device **110** for execution, the illustrative data sharing module **136** enables specialized respiratory clinicians who are not at the patient's bedside to remotely change or update the patient's respiratory prescription, in order to respond in a timely manner to clinical changes in the patient's condition or for other reasons.

[0061] One illustrative example of a scenario in which the data sharing module **136** may be used is as follows. Using a computer on which the configurable user interface module **114** is configured to include the prescription creator module **112**, e.g., the prescription creator module **112** is either installed or accessible via a network (e.g., the "cloud"), such as a hospital computer used to manage electronic medical records, a physician creates a combined respiratory therapy prescription. The layering module **130**

automatically selects a combination device **110** to be used to implement the prescription. The patterning module **134** displays the chronological pattern of prescribed therapy in clock time. The sequencing module **132** generates the specific details of the treatment sequences that form the treatment sessions, including all of the requisite device settings for the selected combination device **110**. As a result, the physician can, via the configurable user interface module **114**, examine and adjust or modify the patient's combined respiratory prescription quickly and accurately, using, e.g., the hospital computer. Another version of the configurable user interface module **114**, including the data sharing module **136**, may be installed on another computer used by the physician, to allow the physician to view the combined respiratory prescriptions that he or she has created from a remote location, even from home.

[0062] Yet another version of the configurable user interface module **114**, including the data sharing module **136**, may be installed on a computer used by a respiratory therapist. With the data sharing module **136**, the respiratory therapist reviews the combined respiratory therapy prescription (e.g., at a hospital computer) previously created by the physician. The computer may be a handheld device or the therapist may transfer the prescription to another computing device that is a handheld or "mobile" device (such as a smart phone, tablet computer, or personal digital assistant) on which another instance of the configurable interface module **114** or a simplified version thereof including the data sharing module **136** is installed.

[0063] Using the data sharing module **136**, which may be installed directly on the handheld computing device or accessed via a network (e.g., the "cloud"), the therapist is able to plan his or her shift schedule by viewing the 24-hour prescription timelines (e.g., treatment patterns), and/or other details of the combined respiratory therapy prescriptions, for all of the patients under his or her care. Using the handheld device, the therapist can review the treatment plans at any time, whether the therapist is at the patient's bedside or at a distant location. With the data sharing module **136**, the therapist's handheld device can be programmed to issue reminders that can help keep the therapist on schedule by generating an audio and/or visual alert at or in advance of the start time for the patient's next scheduled therapy. Also using the data sharing module **136**, the therapist's handheld device may link to an electronic communication network (e.g., a

hospital paging system, computer network, or telecommunications network), to alert the therapist when new combined respiratory therapy prescriptions are created or existing prescriptions are modified by other caregivers, for example. Further, using the data sharing module **136**, new or updated combined respiratory therapy prescriptions can be transferred electronically (e.g., via a direct, hard-wired connection or over a wired or wireless network) to the therapist's handheld device (wherever it may be located) in real time or downloaded to the handheld device from, e.g., a hospital computer. When the therapist signs out from a shift, he or she can transfer the combined respiratory therapy prescriptions for the patients in his or her care to a handheld device used by the next therapist that is coming on duty, by linking the handheld devices using a wired or wireless (e.g., WIFI or Near Field Communication (NFC)) data communication connection.

[0064] When a patient is ready to change venues, e.g., to move to a different hospital or nursing home, or to return to the patient's own home, the data sharing module **136** can be used to send the patient's combined respiratory therapy prescription to the patient's new venue via wired or wireless data communication as described above. As the combined respiratory therapy prescription specifies the combination device **110** to be used to perform the combined respiratory therapy prescription, as well as the particular schedule and combination of respiratory therapies to be performed for the patient (including the device settings), the transfer of the prescription to a computing device located at the patient's new venue should enable the patient's respiratory care to be continued at the new location relatively seamlessly.

[0065] The illustrative problem-first device control module **140** interfaces with the prescription translation module **122** to directly implement changes to a patient's combined respiratory therapy prescription "on the fly," e.g., in real-time during the patient's therapy, according to the patient's preferences or as the patient's health condition changes. The problem-first device control module **140** includes computer logic and data (e.g., look-up tables or the like) that map various device settings of the combination device **110** to different clinical conditions. For example, the problem-first module **140** may derive triggering conditions and desired therapeutic changes from evidence-based guidelines or from the patient's own treatment history data (which may indicate therapies that have been successful or unsuccessful for the patient in the past). As such, the problem-first

module **140** allows caregivers and others to modify the patient's combined respiratory therapy prescription simply by indicating the clinical change to the problem-first module **140**. For example, a therapist may notice that his or her patient is currently unable to cough up secretions without assistance. In this scenario, the caregiver can input "unable to cough up secretions" to the problem-first module **140**, using a graphical user interface provided by the configurable user interface module **114**, for example. The prescription translation module **122** translates the clinical change into the appropriate device setting changes for the combination device **110**. In embodiments where the problem-first module **140** is not integrated with the combination device **110**, the data sharing module **136** transfers the device setting changes directly to the combination device **110**, which implements the device setting changes. In other words, the problem-first device control module **140** enables direct, automatic prescription revision, without requiring the user to view or compose an entire prescription.

[0066] The illustrative problem-first device control module **140** can also allow the user to respond to data transmitted by the combination device **110**, via the data sharing module **136**. Based on changes in the data received from the combination device **110**, which the caregiver may view using the data sharing module **136**, for example, the caregiver may determine that adjustments to the patient's combined respiratory therapy prescription are needed, and implement those adjustments using the problem-first device control module **140** as described above.

[0067] One illustrative example of a scenario in which the problem-first device control module **140** may be used is as follows. Suppose a home care company receives a new patient who had been discharged from the hospital. Using the data sharing module **136**, the patient's combined respiratory therapy prescription has been electronically sent to the home care company's computer. As a result, the company knows which combination device **110** it needs to provide for the patient and also knows the patient's specific respiratory therapy treatment plan.

[0068] Using the data sharing module **136** as implemented on his or her handheld device, a home care therapist employed by the home care company downloads the patient's combined respiratory therapy prescription to his or her handheld device and brings it to the patient's home. The therapist may then use the data sharing module **136**

to transfer the patient's combined respiratory therapy prescription directly to the control module **120** of the combination device **110** (e.g., by connecting his or her handheld device to the combination device **110**). If a prescription translation module **122** is installed on the therapist's handheld device, the prescription may be translated to machine-readable instructions at the hand-held device and then implemented directly by the combination device **110**. Alternatively, the prescription may be translated by the prescription translation module **122** of the control module **120**.

[0069] During or after a therapy session, the combination device **110** can send data about the therapy session or data relating to the patient's condition or preferences to the clinician's handheld device (e.g., in the form of a notification message). Based on this notification and, perhaps, a telephone follow-up with the patient, the clinician can use the problem-first module **140** to change the patient's combined respiratory prescription and send the new prescription to the device **110** using the data sharing module **136**. Using the data sharing module **136**, the new or changed prescription can be made available at a display of the device **110** or at other electronic devices, for viewing by the home care company, other caregivers, the patient, and/or family members of the patient, for example.

[0070] The audio interface module **142** includes a software-based user interface to the prescription creator module **112**, which allows the patient and/or family members or others associated with the patient to view the patient's combined respiratory therapy prescription and configure reminders, alerts, and other messages relating to the patient's therapy prescription. In some embodiments, the audio interface module **142** provides a software-driven, human-voice natural language interface to the combination device **110**. The audio interface **142** maps pre-recorded human voice messages (or computer-synthesized spoken natural language messages) to various aspects of the patient's combined respiratory therapy prescription, as may be desired by the patient or configured by a caregiver or family member. For example, the recorded messages may provide instructions on how to use the combination device **110** or adjust its settings. These instructional-type messages may be timed to be played prior to the start of a therapy session or upon the user's request. For example, the patient may input a coded question

such as "how do I turn this device on?" and in response, the audio interface **142** may provide the requested instructions.

[0071] In some embodiments, the audio interface **142** may be programmed to play recorded messages of an inspirational or reassuring nature at appropriate times, prior to, during, or after a therapy session. In some cases, the content and timing of these types of messages is based on the patient's preferences, focus groups, and/or research relating to the psychology of people who have chronic illnesses. For example, patients with chronic illnesses may have psychological profiles that could impair their compliance with the use of respiratory therapy devices and other medical devices, such as chronic anxiety or social or behavioral disorders. Additionally, pediatric patients can be especially fearful of mechanical devices. Thus, the use of respiratory care devices and other medical devices can be viewed as a burden rather than a benefit by many patients, resulting in poor compliance and limited device efficacy. The audio interface **142** is therefore designed to implement the concept of anthropomorphism (the attribution of human characteristics to non-living things) to improve patient compliance with his or her combined respiratory therapy prescription by making the device **110** more appealing and enjoyable to use.

[0072] The audio interface **142** also allows the patient or other user to choose (e.g., from a list of choices presented on a touchscreen display of the configurable user interface module **114**) the timing of the desired messages. For example, the patient may specify that a message is to be played in the morning, on the patient's awakening, before the start of a treatment session, during a treatment session, after a treatment session, and/or at bedtime. In other words, a particular recorded message may be linked to one or more portions of the patient's combined respiratory therapy prescription (e.g., treatment sessions and/or treatment sequences) over the course of a treatment pattern (e.g., a 24-hour timeline).

[0073] In the illustrative embodiments, the audio interface module **142** interfaces with a persona configuration module **144**, which allows the patient to select an identity or personality to be ascribed to the combination device **110**. An example of a user interface that may be implemented by the configurable user interface module **114** in connection with the persona configuration module **144** is shown in FIG. 5. In FIG. 5, an illustrative

display screen **500** includes a number of selectable choices **510**, each of which embodies a different persona or "mascot" that may be ascribed to the combination device **110**. The selection of a choice **510** automatically configures the character of the voice (e.g., pace, tone) and the content of the pre-recorded audio messages to correspond to the selected persona/mascot. In some embodiments, a graphical or animated depiction of the selected character may be displayed to the user, as well.

[0074] The illustrative display screen **500** is a touch-sensitive screen such that selection of a choice **510** can be accomplished by the patient or another user simply touching the desired choice on the screen with a hand, finger, stylus, or the like. Once a choice **510** is selected, the persona configuration module **144** configures the content, timing, voice, and intonation of the pre-recorded audio messages to correspond to the selected persona. To do so, the persona configuration module **144** may select and download recorded messages from a pre-recorded message database that are tagged or otherwise associated with the selected persona. For example, the pre-recorded message database may associate the characteristics of "reliable," "tenacious," and "determined" with the "Winston the Bulldog" persona, and so on. In this way, the combination device **110** can be customized to project a set of human characteristics that appeal to the patient, so that the patient may perceive the device **110** as an ally and companion rather than as an inanimate object that is a threat or a burden, and thereby facilitate interactions between the patient and the combination device **110**.

[0075] One illustrative example of how the audio interface module **142** and the persona configuration module **144** can be used to customize the message output provided by the combination device **110** is as follows. Suppose a patient would like her combination device **110** to exhibit the human qualities of tenacity and reliability. The patient can use the touchscreen display to input these desired characteristics to the audio interface module **142** (by, for example, selecting them from a drop-down list). The persona configuration module **144** maps the patient's choices to one or more pre-defined personas, which it displays in the list of selectable choices **510** on the display screen **500**. From the list of choices **510**, the patient selects choice D, "power girl." As a result, the audio interface **142** plays a reassuring message each night, just before the combination device **110** begins a therapy session to assist the patient's breathing during sleep, such

as, "I will help you breathe comfortably all night. If you need a cough, just let me know by pressing your thumb switch. I won't let you down. See you in the morning!"

[0076] In some embodiments, aspects of the audio interface **142** are adapted for use by clinicians, therapists, or other users, alternatively or in addition to its use in connection with the patient. For example, the audio interface **142** may be configured to output instructions for preparing a combined respiratory therapy prescription or instructions for using the combination device **110**, in a spoken natural language form, to a caregiver. As another example, the audio interface **142** may be configured as an interface to the problem-first module **140** of the device control module **140**, described above. Additional graphical user interface screens for display on the touchscreen display **500** are discussed with respect to FIGS. 12-16L, below.

[0077] Referring now to FIG. 2, an illustrative method **200** executable as computerized programs, routines, logic and/or instructions by the combined respiratory therapy prescription creator module **112** and/or one or more of the other modules of the system **100** to create a combined respiratory therapy prescription for a patient is shown. The method **200** may be viewed as one example of how the combined respiratory prescription creator module **112** could work. In that case, the primary user would likely be a physician using the method **200** to create a prescription for one of his or her patients. At block **210**, the method **200** receives information about a patient's respiratory condition. Such information may include a symptom or a recent change in the patient's clinical condition, which may be input by a clinician, such as a physician, or by a caregiver, or even by the patient or a family member, for example (using, e.g., the problem-first device control module **140**). At block **212**, the method **200** determines which of the therapeutic layers (e.g., mucus mobilization **214**, mucus extraction **216**, lung volume recruitment **218**, lung ventilation **220**) are associated with the information about the patient's condition received at block **210**. For example, if the input at block **210** indicates that the patient is having trouble generating a productive cough on his or her own, the method **200** maps that information to one or more of the therapy layers **214**, **216**, **218**, **220** that are directed to providing cough assistance. In this case, those therapy layers include the mucus extraction and either lung ventilation or lung volume recruitment therapy layers. Alternatively, the method **200** may receive information about the patient's respiratory

condition, as mentioned above, from the clinician creating the combined respiratory therapy prescription, and the clinician may select and determine the therapeutic layers by himself or herself (e.g., manually), using, for example, the combined respiratory therapy prescription creator module **112**.

[0078] The determination of whether to select lung ventilation **220** or lung volume recruitment **218** can be based on additional inputs received at one or more of the modules **114**, **120**, or stored information about the patient's current health condition or clinical history. For example, at block **212**, the method **200** may access electronic medical records associated with the patient and thereby determine that the patient has responded well in the past to lung ventilation therapy provided after mucus extraction therapy. As another example, the method **200** may access date and time information that is automatically kept by the system **100**, determine therefrom that the patient is likely to be awake and able to participate at the time of the therapy session, and select the lung volume recruitment therapy layer **218**, rather than the lung ventilation therapy layer **220**, as a result.

[0079] Once the therapy layers associated with the patient's condition have been determined, at block **222**, the method **200** automatically selects an appropriate combination device **110** from the family of combination respiratory therapy devices, based on the therapy layers determined at block **212**. Illustratively, the family of combination devices includes four combination devices **224**, **226**, **228**, **230**. The mapping of therapy layers **214**, **216**, **218**, **220** that can be provided by each of the devices **224**, **226**, **228**, **230** is shown illustratively by lines connecting the various therapies with the corresponding combination devices, in block **222**. For example, the combination device **224** can be used to provide both mucus extraction and lung ventilation therapies. The combination device **226** can provide both mucus extraction and lung volume recruitment. The combination device **228** can provide all three of mucus mobilization, mucus extraction, and lung ventilation therapies. The combination device **230** can provide all three of mucus mobilization, mucus extraction, and lung volume recruitment therapies. The method **200** determines which of the combination devices **224**, **226**, **228**, **230** to select based on each device's capabilities in relation to the therapy layers **214**, **216**, **218**, **220** determined at block **212** to be needed by the patient. For example, if the patient needs mucus

mobilization, the method **200** may select either device **228** or device **230**, but not device **224** or device **226**. If the patient needs cough assistance therapy but can otherwise breathe on his or her own, the method **200** may select device **226** or device **230**, but not device **224** or device **228**.

[0080] Once a combination device has been selected, at block **232**, the method **200** obtains the information it needs from the user to prepare a combined respiratory therapy prescription for the patient using the combination device **110** selected at block **222**. To do this, the method **200** interfaces with the user to define one or more treatment sessions **234**, define one or more treatment sequences **236** for each treatment session, and define the treatment patterning **238** over a period of time during which the patient is to receive respiratory care. In the illustrated examples, a treatment session is made up of a number of treatment sequences that are applied consecutively, e.g., repeated successively a defined number of times, where each treatment sequence includes one or more assisted cough cycles followed substantially immediately by lung ventilation or lung volume recruitment therapy. Thus, the duration of a treatment session can depend on the number of treatment sequences to be provided during the treatment session. So, the process of defining the treatment session or sessions **234** involves the method **200** interfacing with the user to specify, over the respiratory care period or "pattern," the number of treatment sessions to be performed, the start time for each treatment session, and the number of treatment sequences to be performed in each treatment session.

[0081] Next, the method **200** interfaces with the user to define the details of each of the treatment sequences **236** to be performed in each of the treatment sessions defined at block **234**. To do this, the method **200** interfaces with the user to specify the number of assisted cough cycles in each treatment sequence, the inspiratory and expiratory (suction) pressures for each cough cycle (e.g., +25 cm water inspiratory pressure, -30 cm water expiratory (suction) pressure), the amount (duration of time) of assisted ventilation to follow the cough cycles (e.g., 2 minutes), and the positive inspiratory pressure and positive expiratory pressure levels for the assisted ventilation therapy (e.g., +15 cm water inspiratory pressure, +4 cm water expiratory pressure).

[0082] At block **238**, the method **200** may interface with the user to define additional respiratory therapies that may be applied to the patient during the respiratory care period

or pattern (e.g., a 24-hour period). For example, the caregiver may wish to schedule one or more mucus mobilization therapies to occur prior to a treatment session, or add an additional daytime or night-time lung ventilation or lung volume recruitment therapy. Accordingly, at block **238**, the method **200** interfaces with the user to specify the start times, stop times, and device settings for each of the additional desired therapies. Additionally, at block **232**, the method **200** may interface with the user to receive other details relating to the combined respiratory therapy prescription. For example, the user may wish to specify that certain portions of the patient's prescription can be modified by the patient or a family member, while other portions can only be modified by the user or an authorized physician, or that some portions can be modified by the patient or family member with the user or physician's authorization. Once the combined respiratory therapy prescription is complete (to the satisfaction of the user), the method **200** electronically communicates the combined respiratory therapy prescription to the combination device **110** for execution by the device **110**. As noted above, this can be done using a wired or wireless data communication method.

[0083] One illustrative example of a user interface **400** that may be provided in connection with the process of creating a combined respiratory therapy prescription using the system **100** is shown in FIG. 4. The user interface **400** displays a timeline **410**, a legend **412** describing the key abbreviations used in the timeline **410**, various details of the respiratory care pattern in their sequential order of occurrence along the timeline **410** (described below), and a simulation feature **440** which allows the user to see an animated simulation **442** (e.g., an animated graphic or video clip) of respiratory therapy as it would be applied to the patient's lungs. To view a simulation, the user may select or highlight one of the therapies or treatment sessions displayed on the timeline **410** and then select the view simulation button. The system **100** then locates and accesses a stored simulation that corresponds to the selected therapy or treatment session (where, for example, such simulations may be indexed or tagged according to their associated therapies or treatment sessions and stored in a database).

[0084] The illustrative timeline **410** includes two treatment sessions, TS-1 and TS-2, as well as a mucus mobilization therapy session M-2 and a night-time nasal lung ventilation therapy session NIV. Each of these sessions has an associated start time. For

example, the treatment session TS-1 has a start time of 8:00 a.m., the mucus mobilization therapy session has a start time of 7:50 p.m., the treatment session TS-2 has a start time of 8:00 p.m., and the night-time nasal lung ventilation therapy session has a start time of 10:00 pm. The start time and/or end time and duration of each of these sessions can be varied using selectable markers **418, 420, 430, 432, 434, 436, 438**. For example, the caregiver may select or "click" on a marker and drag or slide it horizontally to the right or left to change the patient's therapy pattern or schedule. Moving a start-time marker (e.g., markers **418, 438, 430, 434**) to the left causes the therapy to have an earlier start time, while moving the start-time marker to the right will cause the therapy to start later in the day. Moving an end-time marker (e.g., markers **420, 430, 432, 436**) to the left will decrease the duration of the therapy, while moving the end-time marker to the right will increase the therapy duration. In the illustrated example, the marker **430** is both a start-time marker (for the treatment session TS-2) and an end-time marker (for the mucus mobilization therapy MM). This indicates to the system **100** that the treatment session TS-2 is to begin substantially immediately upon the completion of the mucus mobilization session MM. In this way, dependencies between the various therapies can be created so that their performance can be coordinated automatically.

[0085] Each of the therapy sessions (TS-1, MM, TS-2, NIV) can have an expand/contract button (e.g., **422, 426, 428**) associated with it. The expand/contract buttons **422, 426, 428** can be selected to show or hide further details about the therapy session, such as the number of treatment sequences, the device settings, etc.). In the illustrated example, the expand/contract button **422** has been selected to show further details of the treatment session TS-1. Those details are displayed in a window **414**. The window **414** shows that the treatment session TS-1 is made up of 5 treatment sequences. Each of the treatment sequences (SEQ 1, SEQ 2, SEQ 3, SEQ 4, SEQ 5) has its own start and end markers **444, 448, 450, 452, 454, 456**, which the user can slide back and forth horizontally to adjust the duration of the treatment sequence (e.g., to adjust the number of cough cycles in the treatment sequence or the duration of the lung ventilation). The markers **448, 450, 452, 454** act as both start and end markers, thereby making the start time for the beginning of the next treatment sequence (e.g., SEQ 2) dependent on the completion of the previous treatment sequence (e.g., SEQ 1) rather than on a specific

clock time. Each of the treatment sequences also has an expand/contract button **424**, **458**, **460**, **462**, **464** associated with it. As such, the caregiver can view and/or modify the details of a particular treatment sequence by selecting the corresponding expand/contract button **424**, **458**, **460**, **462**, **464**. In the illustrated example, the button **424** has been selected to show the details for the treatment sequence SEQ 1. As shown in the window **416**, these details include four assisted cough cycles with +25 inspiratory pressure/-30 expiratory (suction) pressure followed by two minutes of lung ventilation at +15 inspiratory pressure/+4 expiratory pressure. The details shown in the windows **414**, **416** can be hidden by selecting the corresponding expand/contract button (e.g., **422**, **424**) again. Similarly, the simulation **442** can be hidden by selecting the button **440** a second time. In some embodiments, the window **416** is interactive (e.g., it contains one or more text boxes) so that the details shown therein can be edited directly by the caregiver.

[0086] Referring now to FIG. 3, an illustrative control unit **300** for the combination device **110** is shown in greater detail. The control unit **300** is embodied as a housing (e.g., plastic or metal), which contains or supports, as the case may be, the electronic and mechanical components shown inside the dashed lines in FIG. 3. In some embodiments, the housing is sized and designed so that the control unit **300** is relatively lightweight and portable. For example, some embodiments of the control unit **300** are configured so that they may be mounted to a patient support apparatus, such as a wheelchair, stretcher, lift, hospital bed, or other patient transport device.

[0087] The housing has defined therein a number of ports to which a number of patient interfaces **322**, **324**, **326** can connect to provide various forms of respiratory therapy to the patient. The positive airflow patient interface **322** is one exemplary embodiment of the positive-pressure airflow patient interface **124** shown in FIG. 1. The patient interface **322** is embodied as a nasally mounted device that contains a pair of air delivery conduits, each of which is configured to engage one of the patient's nostrils. As such, the patient interface **322** is configured to supply positive-pressure airflow to the patient via the patient's nose.

[0088] The patient interface **324** is one exemplary embodiment of the negative airflow patient interface **126** of FIG. 1. The patient interface **324** is embodied as mask that is designed to engage with the patient's mouth area to supply airflow through the patient's

mouth. In the embodiment of FIG. 3, the patient interface **322** is configured to supply only positive airflow and the patient interface **324** is configured to supply only negative-pressure airflow to the patient. In other words, the illustrative patient interface **324** is only used for the negative-pressure portions of assisted cough cycles, and is not used to provide lung ventilation or lung volume recruitment therapy. Similarly, the patient interface **322** is only used for lung ventilation, lung volume recruitment, and the inspiratory pressure portion of assisted cough cycles, and is not used during the negative-pressure portions of the assisted cough cycles. This configuration of the interfaces **322**, **324** keeps the positive and negative air flow circuits separate, and thus free of contamination. In other embodiments, however, the interfaces **322**, **324** may be combined or integrated as a single patient interface; for example, as a single patient interface having separate positive and negative airflow circuits.

[0089] The air pulse patient interface **326** is one exemplary embodiment of the patient interface **128** shown in FIG. 1. In some embodiments, the interface **326** is embodied as a wearable element to which tubing can be connected to supply air pulses to the patient's chest region when worn. One example of such a device is THE VEST, available from the Hill-Rom Company, Inc. Alternatively or in addition, some embodiments may provide certain forms of air pulse therapy using the negative airflow patient interface **324**. For example, a METANEB device or other type of continuous high frequency oscillation (CHFO) device **364** may be connected to the negative airflow patient interface **324**. Other devices that provide various forms of air pulse therapy may also be used in a similar fashion.

[0090] In the embodiment of FIG. 3, all of the computer programs and other components that provide the functionality of the system **100** reside in the control unit **300**. That is, all of the various features of the system **100** provided by the various modules described above can be accessed and used directly at the control unit **300**. Accordingly, the illustrative control unit **300** includes therein a controller **310**, which may be embodied as one or more microprocessors, microcontrollers, digital signal processors, or the like. The controller **310** communicates electronically with many other elements of the control unit **300** via a data communication link or bus **316** (e.g., a Controller Area Network bus or the like). In some embodiments, the control module **120**, the configurable user interface

module **114**, the prescription database **118**, and/or any of their respective submodules, described above, are embodied as software that is stored in e.g., disk storage, and then loaded into memory **312** (e.g., random-access memory (RAM)) at runtime as needed. In some embodiments, portions of the data **118** and/or the modules **114**, **120** may be embodied as firmware residing in non-volatile memory. Further, in some embodiments, the memory **312** may be integrated with the controller **310**. Accordingly, the simplified illustration of FIG. 3 showing the configurable user interface module **114**, the prescription database **118**, and the control module **120** embodied in the memory **312**, which is accessible to the controller **310**, is intended to cover all of the various possible embodiments of the database **118** and modules **114**, **120**, whether implemented as software, firmware, hardware or a combination thereof.

[0091] The control unit **300** includes a control panel **318**, which, as indicated by the schematic of FIG. 3, may have its own power supply. The illustrative control panel **318** includes a display screen **320**, which may be embodied as a touchscreen display supported by the housing of the control unit **300**. During operation, information about the combination device **110**, such as the current device settings, as well as therapy being performed, may be displayed on the display screen **320**. Features provided by the modules **114**, **120**, and/or data accessed from the database **118**, may be made available through the control panel **318** and/or the display screen **320**, or may be provided at other computing devices as described above. In other words, any of the features of the system **100** described above may be accessible to users via the control panel **318** and/or the display screen **320**, or through other computing devices as described herein, in various embodiments of the system **100**.

[0092] The control unit **300** includes an audio circuit illustratively made up of an audio interface **372**, an audio driver **374**, an amplifier **376** and an audio controller **378**. The audio circuit is configured to allow the system **100** to process audio inputs and output audio through speakers as auditory sound, in order to implement the features of the audio interface **142** described above. While the embodiment of FIG. 3 shows the audio circuit as being part of the control unit **300**, it should be understood that portions of the audio interface **142** may be implemented at a computing device (such as a user's local computing device) using similar components. As such, the patient or another user may

interact with the combination device **110** either via the control panel **318** or via another computing device, in various embodiments of the system **100**.

[0093] The control unit **300** includes a data management module **382**, a network connector **384**, and a power management module **386**. The data management module **382** manages the communication of data (e.g., portions of the patient's combined respiratory prescription, data generated by the device **110** during operation, etc.) from the device **110** to other devices and vice versa, using the network connector **384**. The power management module **382** interfaces with a power supply (e.g., a battery or a wall socket) to supply electrical power to the various components of the control unit **300**. The network connector **384** may include a wireless network interface, Ethernet adapter, and/or other components as may be needed or desired to enable the control unit **300** to electronically communicate with other devices through either a wired or wireless network connection.

[0094] A finger switch **380** is also provided at the control unit **300** and is in electronic communication with the controller **310**. Portions of the finger switch (e.g., a lever, dial, button or toggle) are mounted to the housing of the control unit **300** to be easily accessible to the patient. The controller **310** is configured to turn the execution of the patient's combined respiratory therapy prescription on or off in response to signals received from the finger switch **380**. For instance, in some embodiments, the controller **310** is responsive to the finger switch **380** to activate or deactivate an assisted cough therapy. That is, if the patient feels congested and needs to cough, the patient may activate the finger switch to initiate an assisted cough therapy. Similarly, if a therapy is in progress and the patient becomes uncomfortable, the patient may press the finger switch to discontinue or temporarily stop the therapy.

[0095] The remaining components of the control unit **300** shown in FIG. 3 include mechanical and electromechanical components to effectuate the various aspects of the patient's combined respiratory therapy prescription through the patient interface(s) **322**, **324**, **326**. In operation, the controller **310** executes the combined respiratory therapy prescription by sending control signals to the various components at the appropriate times, via the bus **316** and a number of servo control modules **350**, **354**, **356**, **358**, **368**. The servo control modules **350**, **358** operate control circuits to control motors **386**, **360**, respectively, which operate manifolds **332**, **330**, respectively, to control the flow of air

generated by an air supply **328** (e.g., a blower) to the patient interfaces **322**, **324**, **326**. The servo control module **354** operates a control circuit to control the generation of air pulses by the air pulse generator **352** based on airflow received from the air supply **328** through the manifolds **330**, **332** and the valve **348**. The servo control module **356** controls the operation of the air supply **328** based on parameters supplied by the controller **310** (e.g., on/off, positive/negative airflow, air pressure) in accordance with the combined respiratory therapy prescription. The servo control module **368** operates a control circuit to control the operation of a air supply **362** (e.g., a compressor), which, in some embodiments, may provide airway clearance therapy such as intrapulmonary percussive ventilator (IPV) through a continuous high-frequency oscillation (CHFO) device **364**, valve **366**, and moisture generator **345** (e.g., a nebulizer) to the nasal patient interface **322**.

[0096] The patient interfaces **322**, **324** are connected to the air supply **328** via the manifolds **332**, **330**, respectively. The air circuit **340** for the positive-pressure airflow patient interface **322** also includes an airflow sensor **334**, a filter **336**, and a moisture generator **338**, to ensure that air supplied to the patient via the nose is clean, at the correct pressure vis a vis the combined respiratory therapy prescription, and somewhat moist so as to avoid overdrying the patient's nasal passages. Similarly, the air circuit for the negative-pressure airflow patient interface **324** includes an airflow sensor **342** and a filter **344**. The airflow sensors **334**, **342** and a pressure sensor **346** sense airflow and air pressure, respectively, in their corresponding circuits and provide airflow and air pressure data to a safety monitoring module **370**. The safety monitoring module **370** monitors the air circuits for the occurrence of any malfunctions and to ensure that the respiratory therapy is being provided in accordance with the patient's combined respiratory therapy prescription. In some embodiments, the sensors **334**, **342**, **346** are used to synchronize the operation of the device **110** (e.g., the timing of the application of positive or negative pressure) with the patient's normal breathing pattern as described above. For example, the sensor **346** may detect initiation of a breath by the patient based on the change in air pressure in the air circuit **340**, and initiate the inspiratory phase of a cough cycle, lung volume recruitment therapy, or lung ventilation therapy in response. As noted above, in the illustrated embodiment, the air circuits that supply air to the positive-pressure airflow patient interface **322** and the negative-pressure airflow patient interface **324**, including

the tubing connecting the interfaces **322**, **324** to their respective air manifolds **332**, **330**, are separated from one another.

[0097] FIGS. 6-9 illustrate device control algorithms that can be implemented by the system **100** to change various aspects of the patient's combined respiratory therapy prescription and its execution by the combination device **110**, in real time (e.g., while the patient is receiving respiratory care). Using the problem-first module **140**, for example, the user (e.g., a clinician, caregiver, patient, or family member, as the case may be) can input information about the patient's current respiratory condition, such as "patient can't cough up secretions." The device control module **140** automatically implements therapy adjustments based on the input according to defined algorithms, and then queries the user to determine whether the adjustments were effective. If the user answers that the adjustment didn't help the patient's condition, the system **100** will proceed to the next step in the algorithm as described below with reference to FIGS. 6-9. If the user responds that the adjustment was effective, the system **100** will continue providing therapy according to the current settings, without making any additional changes.

[0098] Referring now to FIG. 6, an illustrative method **600** executable as computerized programs, routines, logic and/or instructions by the device control module **140** and/or one or more of the other modules of the system **100** to, in real time, adjust a patterning-spacing aspect of the patient's combined respiratory therapy prescription, either automatically or in response to user input, is shown. In this example, the system **100** detects (either automatically based on sensor data or through analysis of user input) that the "patient can't cough up secretions." The method **600** starts operating at the patterning (e.g., 24-hour timeline) level at block **610**, and, at block **612**, determines (e.g., automatically or by issuing a query to the user) the spacing of the patient's various therapies over the course of the therapy timeline. For instance, at block **612**, the method **600** determines how much time currently elapses between the patient's cough assistance treatment sessions while the patient is awake (e.g., are there more than four hours between each session?). If the answer is yes, the method **600** updates the patient's prescription to change the spacing of the cough assistance treatment sessions during the patient's waking hours to occur every four hours (e.g., to increase the frequency of the treatment sessions), at block **614**. If the answer is no that means that the patient is already

receiving cough assistance therapy at least every four hours while awake. At block **616**, the method **600** determines whether the spacing between treatment sessions is more than two hours but less than or equal to four hours, while the patient is awake. In other words, is the interval between treatment sessions more than two hours but not more than four hours? If the answer is yes, the method **600** updates the patient's combined respiratory therapy prescription so that the cough assistance treatment sessions occur every two hours, at block **618** (e.g., to increase the frequency of the treatment sessions). If the answer is no, then the method **600** maintains the current frequency of cough assistance treatment sessions (deduced as being less than or equal to two hours), at block **620**.

[0099] Referring now to FIG. 7, an illustrative method **700** executable as computerized programs, routines, logic and/or instructions by the device control module **140** and/or one or more of the other modules of the system **100** to, in real time, adjust a therapy layering aspect of the patient's combined respiratory therapy prescription, is shown. Whereas the method **600** can be used to adjust the frequency of respiratory treatment sessions over the course of a period of time, the method **700** is directed to determining whether the particular combination device **110** being used by the patient should be changed. At block **710**, the layering changes algorithm starts, either automatically or in response to input from a user indicating, for example, a clinical change in the patient's health condition.

[00100] At block **712**, the method **700** determines (again, either automatically or based on user input) which features of the combination device **110** the patient is presently using, or which combination device **110** from the family of combination devices described above is currently in use. For example, the method **700** may determine whether the patient is already receiving mucus mobilization therapy in addition to mucus extraction therapy and either lung volume recruitment therapy (via the combination device **224**, for example) or lung ventilation (via the combination device **226**, for example). If the answer is yes, then the method **700** continues providing the current therapy without any changes, at block **714**. If the answer is no, then at block **716** the method **700** automatically adds mucus mobilization therapy to the patient's combined respiratory therapy prescription, or instructs the user to do so. This may be accomplished by, for example, activating the air

pulse patient interface **128** of the patient's existing combination device **110** or switching the patient to a different combination device (e.g. device **228** or device **230**).

[00101] Further at block **716**, the method **700** updates the patient's combined respiratory therapy prescription to add a ten-minute mucus mobilization therapy session before each scheduled cough assistance treatment session. At block **718**, the method **700** determines if the patient is already receiving mucus mobilization therapy using a combination device **110** that is also providing lung volume recruitment therapy to the patient (e.g., device **230**). If the patient is already receiving mucus mobilization therapy with lung volume recruitment therapy, then the method **700** displays a message suggesting that the user switch the patient to a device **110** that can provide lung ventilation therapy (in place of the lung volume recruitment therapy) as well as the mucus mobilization therapy (e.g., device **228**), at block **720**. If the answer is no (meaning that the patient is already receiving both mucus mobilization and lung ventilation therapy), then at block **722** the method **700** displays a message suggesting that the patient continue using the same device **110** without any changes.

[0100] Referring now to FIG. 8, an illustrative method **800** executable as computerized programs, routines, logic and/or instructions by the device control module **140** and/or one or more of the other modules of the system **100** to, in real time, adjust a patterning-duration aspect of the patient's combined respiratory therapy prescription, is shown. Whereas the method **600** may be used to adjust the spacing or time interval between treatment sessions over the course of a respiratory care period (e.g., timeline), the method **800** can be used to, automatically or in response to user input, adjust the length or duration of individual treatment sessions in the patient's prescription. At block **810**, the patterning duration change algorithm begins, in response to the system **100** determining, for example, a clinical change in the patient's condition. At block **812**, the method **800** determines whether the treatment sessions currently defined in the patient's combined respiratory therapy prescription are shorter in duration than fifteen minutes. If so, the method **700** adjusts the patient's combined respiratory therapy prescription to increase the duration of the treatment sessions to fifteen minutes, at block **814**. If not, the method **700** does not make any changes to the existing specifications for the duration of the patient's treatment sessions, leaving them at the current duration of fifteen minutes or

longer. At block **818**, the method **800** determines whether the patient is receiving either lung ventilation or lung volume recruitment therapy after each treatment session for at least fifteen minutes. If so, the method **800** maintains the current settings for lung ventilation or lung volume recruitment, as the case may be. If not, the method **700** adjusts the patient's combined respiratory therapy prescription to increase the duration of the patient's lung ventilation or lung volume recruitment therapy to fifteen minutes.

[0101] Referring now to FIG. 9, an illustrative method **900** executable as computerized programs, routines, logic and/or instructions by the device control module **140** and/or one or more of the other modules of the system **100** to, in real time, adjust a sequencing aspect of the patient's combined respiratory therapy prescription, is shown. That is, the method **900** is directed to adjusting specific details of cough assistance treatment sequences in response to, for example, clinical changes in the patient's condition, which the system **100** may detect automatically or receive via user input. The method **900** begins at block **910**, in response to determining that the patient is having difficulty clearing chest secretions on his or her own. At block **912**, the method **900** determines whether the treatment sequences in the patient's current combined respiratory therapy prescription already include four assisted cough cycles followed by at least two minutes of assisted ventilation. If yes, then at block **914** the method **900** continues the therapy according to the existing prescription without making any changes. If no, then at block **916** the method **900** adjusts the patient's existing respiratory prescription to include four cough cycles followed by two minutes of assisted ventilation.

[0102] At block **918**, the method **900** determines whether the inspiratory pressure used in the cough cycle is less than 30 centimeters (cm) water (e.g., delivered via the nasal patient interface **322**). If no (meaning that the inspiratory pressure is already at least 30 cm water), then the method **900** continues the therapy according to the existing prescription without making any changes, at block **920**. If yes, the method **900** increases the inspiratory pressure by 1 cm water for each of the next three treatment sessions (e.g., to a maximum of 33 cm water), at block **922**. At block **924**, the method **900** checks to see if the expiratory (suction) pressure (in this case, to the mouth) during the assisted cough cycles is less than -40 cm water (suction), i.e. if the negative pressure is less forceful than 40 cm of water. If not (meaning the expiratory (suction) pressure is already at least -40

cm water), the method **900** continues the therapy according to the existing prescription without making any changes, at block **926**. If yes, the method **900** updates the patient's prescription to increase the expiratory (suction) pressure by two cm water for each of the next three treatment sessions (e.g., to a maximum of -40 cm water). For example, if the negative pressure is currently set at -30 cm water), the suction pressure will increase by 2 cm of water for each of the next 3 treatment sessions, i.e., to a maximum of -36 cm of water. At block **930**, the method **900** determines whether the patient's treatment sessions are currently defined to include at least five treatment sequences. If yes, the method **900** continues the therapy according to the existing prescription without making any changes. If no, the method **900** updates the patient's prescription to increase the number of consecutive treatment sequences in each treatment session to five.

[0103] Referring now to FIG. 10, an illustrative method **1000** executable as computerized programs, routines, logic and/or instructions by the control module **120** and/or one or more of the other modules of the system **100** to, in real time, operate the combination device **110** to provide the appropriate respiratory therapy to the patient at the appropriate times, is shown. At block **1010**, the method **1000** monitors the clock time to determine whether to start a portion of the patient's therapy, according to the patient's combined respiratory therapy prescription. For example, if the patient's prescription indicates that a treatment session is to begin at 8:00 a.m., the method **1000** compares the current clock time to the start time, 8:00 a.m., and when the comparison is successful, proceeds to block **1012**. If no therapy is scheduled to begin at the current clock time, according to the patient's prescription, then the method **1000** simply continues to monitor the clock time at block **1010**.

[0104] At block **1012**, the method **1000** determines the type of therapy that it needs to initiate using the combination device **110**. If the therapy type is mucus extraction, then at blocks **1014** and **1018**, the method **1000** initiates a treatment sequence, which illustratively includes a number of consecutively executed cough cycles. If the therapy type is something other than mucus extraction, then the method **1000** initiates the provision of therapy at block **1016** by configuring the settings of the combination device **110** for the therapy in accordance with the patient's prescription (e.g., therapy pressures, duration, etc.), and performs the therapy at block **1036**, for the prescribed period of time.

While the therapy is in progress, at block **1036**, the system **100** may receive inputs from the user (e.g., a clinician, a caregiver, the patient, or a family member), at block **1038** and adjust the device settings according to the input, at block **1040**. For example, the patient may wish to reduce the inspiratory or expiratory pressure and signal the device **110** to do so using the finger switch **380** described above. The method **1000** monitors the time elapsed during the performance of the therapy at block **1036**, and at block **1042** determines whether it is time for the therapy to conclude (based on the therapy duration specified in the patient's respiratory prescription). If the requisite amount of time has elapsed, the method proceeds to block **1052** and ends the therapy session. If not, the method returns to block **1036** and continues the current therapy.

[0105] Returning to block **1014**, the difference between the treatment sequences and other types of therapy is that the duration of a treatment sequence is based at least in part on the number of repetitions of the prescribed cough cycle, rather than on clock time. Thus, for treatment sequences, the method **1000** keeps track of the number of treatment sequences that have already been performed in the current treatment session. So, at block **1018**, the treatment sequence counter is initially set to zero. Once a treatment sequence is initiated at block **1018**, the method proceeds to configure the combination device **110** for the mucus extraction therapy and for the number of cough cycles specified in the patient's prescription, at block **1020**. At block **1022**, the method **1000** begins performing the mucus extraction therapy (e.g., by providing the number of cough cycles specified in the prescription). As the treatment sequence specifies that the cough cycles are followed substantially immediately by a short period of either lung volume recruitment or lung ventilation therapy, the method **1000** configures the combination device **110** to provide the lung volume recruitment therapy or lung ventilation therapy upon the completion of the cough cycle, at block **1024**, and performs the lung volume recruitment or lung ventilation therapy, at block **1026**. To do so, the method **1000** changes the pressure settings from those used for mucus extraction therapy to those that are appropriate for lung volume recruitment or lung ventilation therapy.

[0106] Upon completion of the lung volume recruitment therapy or lung ventilation therapy (e.g., upon expiration of the time for providing that therapy, or therapy duration), the method **1000** marks the end of a completed cough assistance treatment sequence,

at block **1028**, and increments the number of treatment sequences, at block **1030**. At block **1032**, the method **1000** compares the current number of treatment sequences (e.g., the treatment sequence counter value) to the total number of treatment sequences to be performed during the therapy session, as specified in the patient's combined respiratory therapy prescription. If the value of the treatment sequence counter equals the total number of treatment sequences to be performed during the therapy session, then the therapy session had been completed and the method **1000** proceeds to block **1034** where it resets the treatment sequence counter to zero, and then ends the therapy session, at block **1052**. If the value of the treatment sequence counter is less than the total number to be performed during the therapy session, then the method returns to block **1022** to perform another treatment sequence. As with the other forms of therapy, the treatment sequences may be interrupted and modified by user input in real time. This is illustrated by the loops **1044**, **1046** and **1048**, **1050**, each of which operates in a similar manner to the loop **1038**, **1040** described above. As such, the description will not be repeated here. At block **1052**, the method **1000** returns back to the beginning, block **1010**, to continue monitoring the clock time for the start of the next therapy that is to occur according to the patient's combined respiratory therapy prescription.

[0107] In the foregoing description of the methods **600**, **700**, **800**, **900**, and **1000**, references may be made to the method or the system **100** "determining," "checking," "asking the user," etc. It should be understood that whenever a method or another aspect of the system **100** is executing computer logic, the required inputs may be received from a user, calculated automatically, or accessed from a storage location in computer memory. For example, if an illustrative method described herein indicates that the method asks the user for input, it should be understood that other embodiments may not require such user input, and instead may obtain the needed information from, for example, calculations or by accessing a stored database or lookup table. Likewise, illustrative methods described herein as "determining" certain things may do so by obtaining user input, accessing stored information, or performing calculations, as needed. Further, in the illustrative methods **600**, **700**, **800**, **900**, and **1000**, and in other examples described herein, specific values are mentioned (e.g., air pressures, time durations, etc.). It should

be understood that such values are provided for illustration purposes only, and that this disclosure is not limited thereby.

[0108] Referring now to FIG. 11, an exemplary computing environment **1100** in which the system **100** may be implemented is shown. Whereas FIG. 3 illustrates an embodiment in which all of the features of the system **100** may be accessible directly at the combination device **110**, FIG. 11 illustrates an embodiment in which some of the features of the system **100** may be provided on other devices. Even so, while the computing environment **1100** is shown as involving multiple components and devices, it should be understood that in some embodiments, the computing environment **1100** may constitute a single computing device (e.g., a hospital computer, or a mobile computing device) in combination with the device **110** and/or other devices. In other words, as used herein, the terms "system" and "environment" may refer to a single computing device or a combination of computing devices and/or other components.

[0109] The illustrative computing environment **1100** includes a physician computing device **1110**, a therapist computing device **1130**, a patient computing device **1150**, and one or more other computing devices **1170**, which are in electronic communication with each other, with other computing devices or systems **1170**, and with the combination respiratory therapy device **110**, via one or more electronic communication networks and/or telecommunications networks **1180**. Each of the devices **1110**, **1130**, **1150** is configured to use a variation of the system **100** that is appropriate for the type of user. For example, in some embodiments, various permissions and access controls may be selected for each type of user when the system **100** is initially set up or as new users are added.

[0110] Illustratively, the prescription creator module **112** resides on the physician computing device **1110**, and portions **118A**, **118B**, **118C** of the combined respiratory therapy prescription database **118** are stored on each of the computing devices **1110**, **1130**, **1150**, respectively. The different portions **118A**, **118B**, **118C** of the database **118** may each include subsets of the database **118**. For example, the portions **118A** and **118B** may include the prescriptions for only those patients under the care of the particular physician or therapist using the devices **1110**, **1130**, and the portion **118C** may contain only the prescription for the particular patient using the device **1150**. Similarly, portions

136A, 136B, 136C of the data sharing module **136**, portions **140A, 140B, 140C** of the device control module **140**, and portions **122A, 122B, and 122C** of the prescription translator module **122** may be configured specifically for the user of the corresponding computing device **1110, 1130, 1150**. For example, the data sharing portion **136A** and the device control portion **140A** may include an extended set of features and capabilities, while the data sharing portions **136B, 136C** and the device control portions **140B, 140C** may include more limited functionality, based on the intended users of the respective computing devices **1110, 1130, 1150**. The prescription translation portions **122A, 122B, 122C** may each have the same or similar functionality, or, in some embodiments, the portion **122A** may have greater prescription translation capabilities than the portions **122B** or **122C**, for example. As shown, the audio interface **142** and the persona configuration module **144** reside on the patient computing device **1150**. However, as discussed above, portions or variations of these modules **142, 144** may be adapted for use by other users, such as physicians or therapists, and those portions or alternative versions may reside on one or both of the physician computing device **1110** and the therapist computing device **1130**, respectively.

[0111] In some embodiments, the computerized modules of the system **100** are embodied as a downloadable software application or "app," which can be obtained from a centralized storage location on a network (such as a private hospital or home care company "app store" or "app market"). In these embodiments, there may be a single app that is downloadable by all types of users, which is then configured for the particular user once installed on the user's local computing device. Alternatively, an app store may provide different downloadable apps for different user types, so that the user may select and download the app that contains the functionality needed by that user. For example, one app may contain the prescription creator module **112** while another app may contain the audio interface and persona configuration modules **142, 144** but not the prescription creator module **112**. Of course, permissions and access controls for downloading the apps may be set by an authorized person such as a hospital system administrator.

[0112] Each of the illustrative computing devices **1110, 1130, 1150** includes at least one processor **1112, 1132, 1152** (e.g. a microprocessor, microcontroller, digital signal processor, etc.), memory **1114, 1134, 1154**, and an input/output (I/O) subsystem **1116**,

1136, 1156. The computing devices **1110, 1130, 1150** may be embodied as any type of computing device such as server, an enterprise computer system, a network of computers, a combination of computers and other electronic devices, a personal electronic device such as a mobile, portable, or handheld computing device, smart phone, personal digital assistant, laptop computer, tablet computer, or desktop computer.

[0113] Although not specifically shown, it should be understood that the I/O subsystems **1116, 1136, 1156** typically include, among other things, an I/O controller, a memory controller, and one or more I/O ports. The processors **1112, 1132, 1152** and the I/O subsystems **1116, 1136, 1156** are communicatively coupled to the memory **1114, 1134, 1154**. The memory **1114, 1134, 1154** may be embodied as any type of suitable computer memory device (e.g., volatile memory such as various forms of random access memory). In the illustrative environment **1100**, the I/O subsystems **1116, 1136, 1156** are communicatively coupled to a number of hardware components including various input devices **1118, 1140, 1158** (e.g., a touchscreen, microphone, physical keyboard or keypad, button, or hard panel control), at least one data storage **1126, 1146, 1166**, various output devices **1120, 1140, 1160** (e.g., an LED, display screen, speaker), one or more other peripheral devices **1122, 1142, 1162** (e.g., sound, graphics or media adaptors), and one or more network interfaces **1124, 1144, 1164**.

[0114] The data storage **1126, 1146, 1166** may include one or more hard drives or other suitable data storage devices (e.g., flash memory, memory cards, memory sticks, and/or others). In some embodiments, portions of the prescription database **118A, 118B, 118C** reside at least temporarily in the data storage **1126, 1146, 1166**. Portions of the prescription database **118A, 118B, 118C** may be copied to the memory **1114, 1134, 1154** during operation, for faster processing or other reasons. Further, in some embodiments, portions of any of the software modules of the system **100** may be stored in the data storage **1126, 1146, 1166** and loaded to the memory at runtime.

[0115] The network interfaces **1124, 1144, 1164** may communicatively couple the computing devices **1110, 1130, 1150** to one or more networks **1180**. Such other networks may include a local area network, wide area network, enterprise cloud, and/or the Internet, for example. Accordingly, the network interfaces **1124, 1144, 1164** may include a wired or wireless Ethernet, mobile/cell network, WI-FI, BLUETOOTH, VPN (Virtual

Private Network), or NFC (Near Field Communication) device or adapter as may be needed, pursuant to the specifications and/or design of the particular network **1180**. Thus, the skilled artisan will appreciate that the network interfaces **1124**, **1144**, and **1164** enable bi-directional communication amongst the respective devices **1110**, **1130**, and **1150**.

[0116] Each of the other computing devices/systems **1170** may be embodied as any suitable type of computing device such as, for example, a server, an enterprise computer system, a network of computers, a combination of computers and other electronic devices, a mobile device, any of the aforementioned types of electronic devices, or other electronic devices. For example, in some embodiments, the other computing devices **1170** may include other computers or computer systems of a hospital or other healthcare facility, which run enterprise-type software applications such as electronic medical records (EMR) systems **1172**, admission, discharge, and transfer (ADT) systems **1174**, and healthcare communication systems (e.g., nurse call systems) **1176**. Thus, in some embodiments, the system **100** may communicate with one or more of the systems **1172**, **1174**, **1176**. For example, if a patient undergoing combined respiratory therapy using the device **110** has a clinical change in his or her health condition that requires medical attention, the system **100** may communicate an alert to the responsible nurse or therapist via the healthcare communication system **1176**. As another example, the system **100** may obtain data about the patient's medical history or history of previous treatments from the electronic medical records system **1172** and use that information to configure or adjust the patient's current therapy prescription. Additionally, the system **100** may interface with a healthcare facility's admission, transfer and discharge system **1174** to, for example, automatically send the patient's combined respiratory prescription to a remote computing device upon the patient's discharge from the facility.

[0117] The computing environment **1100** may include other components, sub-components, and devices not illustrated in FIG. 11 for clarity of the description. In general, the components of the computing environment **1100** are communicatively coupled as shown in FIG. 11 by signal paths, which may be embodied as any type of wired or wireless signal paths capable of facilitating communication between the respective devices and components.

[0118] Turning now to FIG. 12, there is shown a main screen **1200** depicting a graphical user interface for interaction with the combined respiratory therapy devices **224**, **226**, **228**, and **230** depicted in FIG. 2. It will be appreciated that while the implementation shown in FIG. 2 depicts four combined respiratory therapy devices **224-230**, any number of such devices may be utilized in accordance with the systems and methods described herein, and the use of four devices **224-230** is intended solely as one example embodiment. As shown in FIG. 12, the main screen **1200** is suitably displayed via the touchscreen **500** described above. The main screen **1200** includes four button icons, preferably color-coded, allowing for interaction with the devices **224-230** by the patient, therapist, clinician, etc., for the control of the devices **224-230**, as explained below.

[0119] Further, the skilled artisan will appreciate that the main screen **1200** may be adapted in accordance with the underlying device **224**, **226**, **228**, or **230** to which the display **500** is associated. Thus, for example, device 1 **224** (mucus extraction and lung ventilation) and device 3 **228** (mucus extraction, lung ventilation and mucus mobilization (oscillations)) may illustrate all four icons, whereas device 2 **226** (mucus extraction and lung volume recruitment) and device 4 **230** (mucus extraction, lung volume recruitment, and mucus mobilization) may only illustrate three of the icons, as discussed in greater detail below.

[0120] FIG. 12 illustrates an “Emergency” button icon **1202**, an “Intensify Therapy” button icon **1206**, a “Customize Your Therapy” button icon **1208**, and a “Cough” button icon **1204**. According to one embodiment, the emergency button icon **1202** is illustrated with a “problem-first” phrase, e.g., “I Can’t Breathe!” to provide the user, such as a patient, an indication as to what problem is related to the underlying button icon **1202-1208**. In varying embodiments, the icons **1202-1208** are color-coded, thereby providing a further distinguishing characteristic should the patient be confused, in distress, etc. According to the illustration of FIG. 12, the emergency icon **1202** is depicted in red on the touchscreen display **500** associated with the combination respiratory therapy devices 1 and 3, respectively noted as **224** and **228** in FIG. 2.

[0121] Upon depressing the button (i.e. touching the icon **1202**) corresponding to an emergency, a corresponding emergency algorithm, stored in memory **312** and executed by the controller **310** (shown in FIG. 3), is activated in accordance with one embodiment

of the subject application. The skilled artisan will appreciate that such activation of the emergency button icon **1202** enables the device **224** and/or **228** to respond immediately to acute respiratory distress and shortness of breath via the enactment of continuous lung ventilation on the patient's usual settings. FIG. 13 provides an illustrative graphical user interface screen **1300** displayed on the touchscreen **500** upon selection of the emergency button icon **1202**.

[0122] As depicted in FIG. 13, the screen **1300** provides an indication as to the mode selected, in this instance, the emergency button icon **1202** is presented, flashing or having another visual indicia associated therewith, indicating the emergency algorithm has been activated. The cough button icon **1204** remains displayed on the screen **1300**, and a new "Escape" button icon **1302** is presented on the screen **1300**, including illustration of a problem-first explanation, e.g., "return to usual settings." When the patient depresses, i.e. selects, the cough button icon **1204**, the device **224** or **228** provides three cough cycles, on demand.

[0123] According to one implementation contemplated herein, when the emergency button icon **1202** has been selected and emergency operations begun, the device **224** and/or **228** may switch automatically to battery power if power is lost, e.g. power failure, unplugging of the device, or the like, such as when an ambulance arrives and the patient must be transported. It will be appreciated by those skilled in the art that the "usual settings" refers to the usual treatment plan, i.e. a prescription for combination respiratory therapy as discussed and described in U.S. Patent No. 9,795,752, by layering (the particular device **224-230** chosen for its capabilities, such as mucus extraction and lung ventilation (device **1 224**)), patterning (the 24-hour timeline indicating the timing of all treatments, such as overnight lung ventilation and cough treatment every four hours); and sequencing (the specific number of cough cycles and settings that comprise a treatment session. For example, that each cough cycle is applied at pressures of +30/-30 cm water pressure, that after every 5 cough cycles the patient gets one minute of assisted ventilation, and that each treatment session consists of five sets of cough cycles followed by ventilation).

[0124] Returning to the main screen **1200** displayed to the patient/caregiver, selection of the "Intensify Therapy" button icon **1206** prompts the touchscreen **500** to display the

interface screen **1400** of FIG. 14. As shown in FIG. 14, the intensify therapy and customize your therapy button icons **1206-1208** are removed and two new options are presented, an “Intensify Level 1 12hr Treatment” button icon **1402** and an “Intensify Level 2 25hr Treatment” button icon **1404**. The intensify treatment screen **1400** further includes the cough icon button **1204**, the escape button icon **1302** and the emergency button icon **1202**, each accessible to the patient/caregiver from the intensify treatment screen **1400**.

[0125] Upon selection of the Intensify Level 1 12hr Treatment button icon **1402** or the Intensify Level 2 25hr Treatment button icon **1404**, the selected icon begins flashing or provides some other suitable visual or auditory indicia relating to its selection. It will be appreciated that the patient/caregiver may present issues relating to the inability, for example, to cough up phlegm during the usual treatment session. The patient/caregiver may then select the intensify therapy button icon **1206** to intensify therapy automatically when patients have retained respiratory secretions that they cannot cough up due to weakness, e.g., with colds. By intensifying therapy immediately, at home, using automated algorithms, patients can avoid hospitalizations due to lung collapse, pneumonia and respiratory failure. Further, the skilled artisan will appreciate that the buttons **1202-1208** and corresponding algorithms associated therewith may also be used by clinicians/caregivers in home care, hospital, and rehabilitation settings as they simplify management and allow immediate responsiveness to changes in the patient’s clinical status.

[0126] The intensify level button icons **1402-1404** allow for the increase in the frequency of cough treatment sessions, and for devices that provide lung ventilation, also increases the duration of lung ventilation in each 24 hour time period. In some embodiments, the device **224**, **226**, **228**, or **230** may continue therapy at the algorithm’s final level of intensity (e.g., treatment sessions every four hours around the clock) until the escape button icon **1302** is activated, which returns treatments to the patient’s usual settings/prescription(s). Selection of the escape button icon **1302** on the screen **500** causes the device **224-230** to cease continuous ventilation and return to the usual treatment plan.

[0127] With respect to Device 2 **226** (mucus extraction and lung volume recruitment), intensity level 1 provides treatment sessions every one hour for four hours, then every

two hours for eight hours, then every four hours around the clock. During this operation, all sets of cough cycles are followed by lung volume recruitment on the patient's usual pressure setting. For example, one minute of lung volume recruitment after every set of five cough cycles. Intensity level 2 for Device 2 **226** provides treatment sessions every twenty minutes for one hour, then every one hour for eight hours, then every two hours for sixteen hours, then every four hours around the clock. During this operation, all sets of cough cycles are followed by lung volume recruitment on the patient's usual pressure setting. For example, one minute of lung volume recruitment after every set of five cough cycles. These treatments continue until such time as the escape button icon **1302** is selected for a return to the patient's usual treatment plan.

[0128] For Device 4 **230** (mucus extraction, lung volume recruitment, and mucus mobilization), intensity level 1 provides treatment sessions every one hour for four hours, then every two hours for eight hours, then every four hours around the clock. During this operation, each session is preceded by oscillatory mucus mobilization and followed by lung volume recruitment on the patient's usual pressure setting. Intensity level 2 for Device 4 **230** provides treatment sessions every twenty minutes for one hour, then every one hour for eight hours, then every two hours for sixteen hours, then every four hours around the clock. As with intensity level 1, each session is preceded by oscillatory mucus mobilization and followed by lung volume recruitment. These treatments continue until such time as the escape button icon **1302** is selected for a return to the patient's usual treatment plan.

[0129] Regarding Device 1 **224** (mucus extraction, lung ventilation), intensity level 1 provides treatment sessions every one hour for four hours, with continuously assisted ventilation after the treatment sessions. Thereafter, every two hours for eight hours, with one hour of assisted ventilation after each treatment session. Next, treatment is provided every four hours around the clock with assisted ventilation for two hours after each treatment session. According to one embodiment, an exception may be provided regarding assisted ventilation throughout the patient's sleep time, if indicated in the patient's usual treatment prescription. Intensity level 2 for Device 1 **224** provides treatment sessions every twenty minutes for one hour, then every one hour for eight hours, then every two hours for sixteen hours, then every four hours around the clock.

Each session is followed by continuously assisted ventilation. These treatments continue until such time as the escape button icon **1302** is selected for a return to the patient's usual treatment plan.

[0130] With respect to Device 3 **228** (mucus extraction, lung ventilation, and mucus mobilization (oscillations)), intensity level 1 provides treatment sessions every one hour for four hours, with continuously assisted ventilation after the treatment sessions. Thereafter, every two hours for eight hours, with one hour of assisted ventilation after each treatment session. Next, treatment is provided every four hours around the clock with assisted ventilation for two hours after each treatment session. Each session is preceded by oscillatory mucus mobilization. According to one embodiment, an exception may be provided regarding assisted ventilation throughout the patient's sleep time, if indicated in the patient's usual treatment prescription. Intensity level 2 for Device 3 **228** provides treatment sessions every twenty minutes for one hour, then every one hour for eight hours, then every two hours for sixteen hours, then every four hours around the clock. Each session is preceded by oscillatory mucus mobilization and followed by continuously assisted ventilation. These treatments continue until such time as the escape button icon **1302** is selected for a return to the patient's usual treatment plan.

[0131] The customize your therapy button icon **1208** may be selected by the patient/caregiver to customize the therapy applied to the patient. It will be appreciated that the customize your therapy button icon **1208** is visible on the touchscreen **500** of each Device 1-4 (**224-230**) in accordance with one embodiment of the subject application. Upon selection of the customization icon **1208**, a series of sequential subsequent displays are presented to the patient/caregiver to enable the user to fine tune the device **224-230** settings based on the algorithms presented herein, as well as based upon each patient's preferences, optimizing comfort and effectiveness.

[0132] Turning now to FIG. 15, there is shown a graphical user interface display screen **1500** associated with the customization initiated by selection of the customization button icon **1208** by the patient/caregiver. FIGS. 16A-16L provide additional illustrations of the various selection components displayed in the selection box **1502** for customizing settings of the Device 1-4 (**224-230**) in accordance with one embodiment of the subject application. It will be appreciated by those skilled in the art that the order in which FIGS.

16A-16L are shown are intended to provide one possible implementation of the interaction between the Device **224-230** and the user (patient/caregiver) during customization, and other sequences of interactions are contemplated herein. Accordingly, the sequential series of graphical user interfaces depicted in FIGS. 16A-16L used for customization of the combined respiratory therapy device **224**, **226**, **228**, or **230** may occur in the sequence illustrated, or may be modified in accordance with the capabilities of the device **224-230** or existing therap(ies) of the patient.

[0133] Beginning at FIG. 16A, the user is presented with a customization option to adjust cough settings **1504** on the touchscreen display **500** of any of Devices 1-4 (**224-230**). However, the skilled artisan will appreciate that for Devices 2 and 4 (**226** and **230**), the emergency button icon **1202** is not displayed. The user is presented with an option to adjust the cough settings by selecting the yes button icon **1506** or declining the option and skipping to ventilation settings (discussed in FIG. 16F). Upon selection of the yes icon **1506**, the display **500** is advanced to FIG. 16B, whereupon the user is presented with three choices for cough inspiratory pressure: low **1512**, normal **1514**, and high **1516**.

[0134] After selection of the inspiratory pressure, the display **500** advances to FIG. 16C, whereupon the user is presented with three options to adjust the cough expiratory (suction) pressure **1518**: low **1520**, normal **1522**, and high **1524**. As used herein, inspiratory and expiratory (suction) pressure levels of the cough cycle: controls the level of inspiratory pressure at the start of the cough cycle and the level of expiratory (suction) pressure that completes the cough cycle and extracts mucus. In an exemplary embodiment, the pressures (low, normal, high) for inspiratory and expiratory (suction) pressures may incorporate the following: “Low” is 2 cm water below the set pressure level for inspiratory pressure and is 2 cm water less negative than the set pressure level for expiratory (suction) pressure, “Normal” for both the inspiratory and expiratory (suction) pressure is the set pressure level, and “High” is 2 cm above the set pressure level for inspiratory pressure and 2 cm water more negative than the set pressure level for expiratory suction pressure.

[0135] After selection of the expiratory (suction) pressure **1518**, the customization screen **1500** advances on the touchscreen display **500** to adjustment of the cough cycle duration **1526** in FIG. 16D. In FIG. 16D, the user is presented with three options for the

duration of the cough cycle: short **1528**, normal **1530**, and long **1532**. In accordance with one embodiment, the duration of the cough cycle **1526** controls the length of the inspiratory pressure, expiratory suction pressure, and the pause. In one embodiment, the following ranges may be used for the presented options **1528-1532**: Short = 1 second duration for inspiratory pressure, 1 second for expiratory suction pressure, then a 1 second pause; Normal = 2 seconds duration for inspiratory pressure, 2 seconds for expiratory suction pressure, then a 2 second pause; and Long = 3 seconds duration for inspiratory pressure, 3 seconds for expiratory suction pressure, then a 3 second pause.

[0136] Thereafter, the customization screen **1500** advances to enable the user to the adjust cough sensitivity **1534** in FIG. 16E. As used herein, the cough sensitivity controls the inspiratory pressure level that triggers the inspiratory phase of the cough cycle. The user may select easy **1536**, normal **1538**, or hard **1540** options corresponding to the sensitivity by which the device **224-230** assists coughing. In one exemplary embodiment, easy **1536** triggers with minimal inspiratory pressure (e.g., 1 lpm flow or 1 cm water pressure), normal **1538** triggers with average pressure (e.g. 2 lpm flow or 1.5 cm water pressure), and hard **1540** triggers with above average pressure (e.g., 3 lpm flow or 2 cm water pressure).

[0137] After adjusting the cough sensitivity **1534**, or directly progressing from the skip option **1508** presented in FIG. 16A, the customization screen **1500** progresses to FIG. 16F, whereupon the user is able to adjust ventilation settings **1542**. It will be appreciated that the option for ventilation settings **1542** applies to Devices 1 and 3 (**224** and **228**) and Devices 2 and 4 (**226** and **230**) skip the steps (screens) relating thereto shown in FIGS. 16F-16I. As shown in FIG. 16F, the user is presented with the option to adjust ventilation settings **1542** for Devices 1 and 3 (**224** and **228**) via selection of the yes button icon **1544** or the no, skip to oscillations button icon **1546** (for Devices 3 and 4 **228-230**). After selection of the yes icon **1544**, the customization display screen **1500** presents the ability to adjust the ventilation inspiratory pressure **1548** as shown in FIG. 16G. The user is presented, in FIG. 16G, with three options: Low **1550** = 2 cm water below the set inspiratory pressure; Normal **1552** = the set pressure level; and High **1554** = 2 cm water above the set pressure level. The customization display screen **1500** thereafter advances

to present the user with the ability adjust the ventilation sensitivity **1556**, as shown in FIG. 16H.

[0138] FIG. 16H illustrates three options for the user (i.e. patient/caregiver) to customize the ventilation sensitivity of the device **224** or **228**, easy **1558**, normal **1560**, and hard **1562**. As contemplated herein, easy **1558** triggers with minimal inspiratory pressure or flow (e.g., 1 lpm flow or 1 cm water pressure), normal **1560** triggers with average pressure or flow (e.g. 2 lpm flow or 1.5 cm water pressure), and hard **1562** triggers with above average pressure or flow (e.g., 3 lpm flow or 2 cm water pressure). Once the sensitivity **1556** has been adjusted by the user, the customization display screen **1500** advances to present the user with the ability to adjust the duration of ventilation after coughs **1564**.

[0139] FIG. 16I illustrates the adjust duration of ventilation after coughs **1564** customization option of the device **224** or **228**. Duration of lung ventilation after the cough cycle **1564** controls the length of time the patient “recovers” from assisted coughing with lung ventilation. For example, after each group of 5 cough cycles within a Treatment Session, the device provides lung ventilation for the following period of time: short **1566** = 30 seconds; normal **1568** = 1 minute; and long **1570** = 3 minutes. After setting the ventilation duration after coughs for device 3 **228**, or after bypassing FIGS. 16H-16I for device 4 (**230**) as discussed above, the customization screen **1500** progresses as shown in FIG. 16J to adjust oscillation settings **1572**, where the user is presented options to proceed to adjustment via the yes button icon **1574** to FIG. 16K or the no button icon **1576**, whereupon the touchscreen display **500** reverts to the main screen **1200** shown in FIG. 12.

[0140] When the user selects the yes button icon **1574**, the customization screen **1500** progresses as shown in FIG. 16K to the adjust oscillations **1578** setting. Via FIG. 16K, the user is able to control the rate and depth of the oscillations provided by Devices 3 and 4 (**228** and **230**), allowing for slow/hard (e.g. 4 Hz) **1580**, normal (e.g. 8 Hz) **1582**, and fast/shallow (e.g. 12 Hz) **1584** oscillations. In accordance with one embodiment, the slow/hard **1580** option provides slower, deeper oscillations, the normal **1582** option provides average frequency and depth of oscillations, and the fast/shallow **1584** option provides faster, shallower oscillations.

[0141] Once the user has selected the desired value for oscillation adjustment, the customization screen **1500** progresses to FIG. 16L, which depicts the setup complete **1586** indication, presenting two options to the user. The user may select the start button icon **1588** which starts the new settings and escapes/returns to the main screen **1200**, or the return button icon **1590** which returns to FIG. 16A to make further adjustments. In varying embodiments, additional customization may be provided, e.g., expanding the scales, where using the touch screen **500**, advanced users can expand the scales to program intermediate settings for each parameter. For example, for the cough cycle duration: Short ||| Normal ||| Long, in which the ||| symbol represents multiple intermediate levels between the named settings.

[0142] Referring now to FIG. 17, in conjunction with the features and structure disclosed and described above with respect to FIGS. 11-16L, the skilled artisan will appreciate the capabilities provided herein for a physician (via physician computing device **1110**) or a therapist (via therapist computing device **1130**) to remotely control the combination respiratory therapy device **224-230**. The capabilities and functionalities with respect to the components set forth in FIG. 11 above are capable of being expanded upon via interaction of the graphical user interfaces depicted in FIGS. 12-16L such that user (patient) selection of an icon (as described below) can initiate the algorithms (e.g. the actions associated with the icon selection as described above with respect to FIGS. 12-16L) to permit remote control of a combination respiratory therapy device 1-4 (**224-230**) by the physician computing device **1110**, the therapist computing device **1130**, or the like.

[0143] FIG. 17 provides an illustrative methodology for the bi-directional interactions between a clinician device (e.g. devices **1110**, **1130** of FIG. 11) and any one or more of the combination respiratory therapy devices **224-230** in accordance with one embodiment of the subject application. The exemplary method begins at **1702**, whereupon the clinician device **1110**, **1130** establishes a bi-directional communication link with one or more combination respiratory therapy devices **224-230**. It will be appreciated that the bi-directional communications may be established via the network interfaces **1124**, **1144**, **1164** as described above with respect to FIG. 11 between the various devices **1110**, **1130**, **224-130**. For example purposes herein, reference will be made to the clinician device as the physician device **1110**, i.e. a tablet, smartphone, laptop, etc., thus the skilled artisan

will appreciate that the term “clinician device” is representative of either the physician device **1110** and/or the therapist device **1130**.

[0144] After establishing the bi-directional communication between the clinician device **1110**, **1130** and the one or more combination respiratory therapy devices **224-230**, operations proceed to **1704**. At **1704**, patient data is received by the clinician device **1110**, **1130** from the one or more combination respiratory therapy device **224-230** via the computer network **1180**. According to one example embodiment, the device **224-230** transmits data (i.e. physiological patient data) to the managing clinician (e.g. the clinician device **1110**, **1130**). Examples of data that the clinician device **1110**, **1130** may receive from the device **224-230** include, without limitation, vital signs such as heart rate, respiratory rate, body temperature, oxyhemoglobin saturation (which may be captured by ancillary devices, e.g. a nasal thermistor for respiratory rate, a pulse oximeter for oxyhemoglobin saturation), physiological data, captured by the combination device **224-230** itself, such as (for devices capable of providing assisted ventilation) the patient's minute ventilation (tidal volume x respiratory rate), tidal volume (volume of gas in each breath), percent of assisted breaths that are patient-triggered, leak flows, inspiratory and expiratory pressures that are actually delivered (compared with set pressures), and the patient's compliance with the prescribed therapies (as discussed above), such as hours per day of device **224-230** use.

[0145] At **1706**, the clinician device **1110**, **1130** generates a graphical display of the received patient data. It will further be appreciated that any or all of the above referenced vital signs or physiological data (i.e. parameters) may be presented on the clinician device **1110**, **1130** in a graphical representation or montage on the associated output device **1120**, **1140**. According to one embodiment, the values of the associated parameters that are normal may be displayed in green, abnormal in red, and the like. In such an embodiment, the ranges of normal values relevant to the particular parameter may be displayed with the patient's measured/actual value shown within or outside the expected range. A particular parameter may be selected at **1708** via clinician interaction with the clinician device **1110**, **1130** for additional information or access to the therapy device **224-230**. The clinician device **1110**, **1130** then sends a query to a particular therapy device **224-230** in accordance with the selected parameter at **1710**.

[0146] The data associated with the selected parameter is then received by the clinician device **1110**, **1130** at **1712** from the particular therapy device **224-230**, whereupon the graphical display on the clinician device **1110**, **1130** is updated to reflect the selected parameter and the received data at **1714**. A determination is then made at **1716** whether the clinician requires a voice/video connection with the patient associated with the therapy device **224-230**. Upon a positive determination, a real-time, bi-directional audio and/or video connection is established between the therapy device **224-230** and the clinician device **1110**, **1130** at **1718**. After the connection is established at **1718**, or upon a determination that no patient connection is required at **1716**, operations progress to **1720**, whereupon a determination is then made at whether a therapy adjustment has been selected. In accordance with one embodiment, the graphical user interface displayed on the clinician device **1110**, **1130** illustrates a button icon similar to that which is described above with respect to FIGS. 16A-16L. That is, the clinician may be presented with one of the button icons described above on the clinician device **1110**, **1130**, remotely from the patient and associated therapy device **224-230**, and may select and activate the various algorithms associated therewith, i.e. provide three coughs, emergency, intensify therapy, etc. In the event that no therapy adjustment is selected, e.g., the clinician merely desires to speak with the patient, direct repositioning of a sensor or mask, etc., operations proceed from **1720** to **1732**, as discussed in greater detail below.

[0147] After a determination is made that a therapy adjustment is desired, operations proceed to **1722**, whereupon the graphical user interface displayed on the clinician device **1110**, **1130** generates a display of the patient parameters and expected ranges associated with the selected therapy being adjusted. Operations then proceed to **1724**, whereupon a therapy adjustment associated with the selected parameter is received at the clinician device **1110**, **1130**. The therapy adjustment is then communicated from the clinician device **1110**, **1130** to the particular therapy device **224-230** at **1726**. As indicated above, the therapy adjustment may occur in the manner described above with respect to the direct interaction at the graphical user interface **1200** of the therapy device **224-230** in FIGS. 12-16L.

[0148] Operations thereafter proceed to **1728**, wherein the clinician device **1110**, **1130** receives update data from the therapy device **224-230** responsive to the administration

of the adjusted therapy. In accordance with one exemplary embodiment, after the therapy device **224-230** provides the adjusted therapy to the patient, vital signs, physiological parameters, etc. (as referenced above at **1706** and **1722**) is provided on the clinician device **1110**, **1130** to enable the clinician to ascertain whether the adjusted therapy is effective/successful. A determination is then made at **1730** whether any additional adjustment to the therapy is needed, e.g., additional coughs to be provided, increase in inspiratory/expiratory pressure, etc. Upon a positive determination at **1730**, operations return to **1722**, whereupon a display of patient parameters and expected ranges associated with the selected therapy adjustment is generated on the clinician device **1110**, **1130**. Operations thereafter proceed through **1724-1730** as discussed above. Upon a determination at **1730** that no additional adjustments to the selected therapy are needed, operations progress to **1732**.

[0149] At **1732**, whereupon a determination is made whether another parameter is selected at the clinician device **1110**, **1130**. When another parameter is to be viewed and/or adjusted by the remote clinician device **1110**, **1130**, operations return to **1708** and operations of the methodology **1700** continue as described above for the next selected parameter. In the event that no additional therapy adjustment is required, operations with respect to FIG. 17 terminate thereafter.

[0150] The preceding methodology **1700** may best be understood in conjunction with one example interaction between the clinician device **1110**, **1130** and the therapy device **224-230**. In this example, the display of the clinician device **1110**, **1130** includes a button icon associated with Tidal Volume adjustment (for devices **224-230** that are capable of providing lung ventilation). Thereafter, the display of the clinician device **1110**, **1130** generates a graphical interface showing the expected ranges of the physiological parameters and the value of the patient's Tidal Volume parameter. For example, if Tidal Volume is expected to be 4 liters +/- 0.5 liter, the graphical display might register a value of 4.25 liters in green, indicating Normal, and a value of 3 liters in red, indicating Abnormal. As indicated above at **1722**, the clinician via the clinician device **1110**, **1130** may contact the patient (via device **224-230**) directly via electronic communications (e.g. audio and/or video) through the respective devices **1110**, **1130**, **224-230** or via telephone. When the Tidal Volume button icon is selected at the clinician device **1110**, **1130** (e.g. steps **1716-**

1726 of FIG. 17), the clinician device **1110**, **1130** modifies the displayed graphical user interface to present new button icons, e.g., a “Tidal Volume Too Low” button icon and a “Tidal Volume Too High” button icon.

[0151] As with the interactions described above with respect to FIGS. 12-16L, control of the combined respiratory therapy device **224-230** is effectuated by the clinician device **1110**, **1130** via the selection of a button icon which triggers the corresponding control algorithm directing operations of the corresponding therapy device **224-230**. Thus, if the clinician touches “Tidal Volume Too Low” button icon displayed on the clinician device **1110**, **1130**, it activates an algorithm (stored in memory **1114**, **1126** and executed by the processor **1112**) that increases the inspiratory ventilation pressure by 2 cm water and transmits the new setting to the device controller **120** of the therapy device **224-230**. Similarly, if the clinician touches the “Tidal Volume Too High” button icon displayed on the clinician device **1110**, **1130**, an algorithm (similarly stored in memory **1114** or **1126** and executed by the processor **1112**) is activated that decreases the inspiratory ventilation pressure by 2 cm water and transmits the new setting to the device controller **120** of the corresponding therapy device **224-230**.

[0152] The clinician may then observe the change in Tidal Volume that occurs with the new settings and determine whether the Tidal Volume now falls within the acceptable range, i.e. in the green zone. When the Tidal Volume has been restored to the green zone, the clinician may select an Escape button icon displayed on the clinician device **1110**, **1130** (as described above with respect to the therapy device **224-230** in FIGS. 12-16L), and the new settings will continue at the corresponding therapy device **224-230**. When the Tidal Volume is still abnormal, the clinician can re-touch the appropriate Button—i.e., “D” if the tidal volume is still too low, or “E” if the tidal volume is still too high. According to one particular embodiment of the methodology **1700** of FIG. 17, as an added safety feature, after 4 adjustments, the clinician can no longer use the Buttons to change the inspiratory ventilation pressure settings.

[0153] It will be appreciated that not all of the parameters transmitted from the therapy device **224-230** to the clinician device **1110**, **1130** are appropriate for button icons associated with device control algorithms. For example, if the therapy device **224-230** is transmitting a body temperature that is too high, the patient has a fever, and the clinician

may choose to contact the patient or caregiver via the audio/video connection using the clinician device **1110**, **1130** or alternatively via telephone. If the therapy device **224-230** is transmitting that there is a high Leak Flow, the patient's nasal mask may not fit properly, and the clinician may choose to troubleshoot this over the telephone or electronically with the patient and caregiver utilizing the clinician device **1110**, **1130**.

[0154] Likewise, at times the clinician may choose to change the combination respiratory device prescription manually such as by entering the Customize Therapy Button (as described in greater detail above) to adjust individual device settings semi-quantitatively with the relevant algorithm, or quantitatively, using the scales to program specific settings for each parameter (as discussed above). In such circumstances, the bi-directional connection between the clinician device **1110**, **1130** and the therapy device **224-230** enables remotely effectuating these changes, thereby allowing the patient to remain in isolation and limit the clinician's exposure to potential pathogens.

[0155] In addition to the foregoing, varying embodiments contemplated herein provide for the incorporation of the graphical representation of the Vital Signs and physiological parameters into the patient's electronic medical record. As will be appreciated by those skilled in the art, a clinician who is reviewing a patient's electronic chart or writing a medical note can pull up the real-time data and see if any adjustments or communications are needed. In a further embodiment, the methodology **1700** described herein enables the use of the clinician device **1110**, **1130** memory (or the device **224-230** memory) to store various parameters recorded over specific periods of time. The skilled artisan will appreciate that the ability to see the trend in the various parameters over specified periods of time; for example, 6 hours, 24 hours, 3 days and 1 week, provides the clinician with valuable insight into the patient's response to treatment, prognosis, and effect of changes.

[0156] It will be appreciated that such an implementation, for example from a central location, allows specialized management of a plurality of combination respiratory therapy devices at distant sites, or alternatively, simply from outside a hospitalized patient's room, so the clinician does not have to enter the room and be exposed to the patient's exhaled pathogens.

[0157] In such implementations, the physician computing device **1110** or therapist computing device **1130** may be implemented as a tablet, smartphone app, personal

computer, or other suitable electronic computing device located outside of the patient's room and in data communication with the corresponding combination respiratory therapy device which is attending to the patient. The skilled artisan will appreciate that the display generated on such a computing device may correspond to the display illustrated in FIGS. 12-17, or alternatively may comprise a specialized view that adds additional configuration options solely for use by the clinician. For example, a remotely observing clinician, e.g. via video or audio monitoring of a patient, may determine that the patient is under extreme duress. Rather than entering the patient's room, the clinician can remotely select the "Emergency" button icon **1202**, directing the combination respiratory therapy device 1 or 3 (**224** or **228**) to perform the emergency algorithms described above.

[0158] In accordance with another example implementation, the remotely observing clinician may observe that the patient associated with a combination respiratory therapy device is experiencing some discomfort that can be alleviated via a particular operation of the combination respiratory therapy device. The clinician, via the physician computing device **1110** or the therapist computing device **1130**, may remotely direct the combination respiratory therapy device 1, 2, 3, or 4 (**224**, **226**, **228**, or **230**) to increase the amount of therapy provided to the patient, i.e. via remote selection of the "Intensify Therapy" button icon **1206**, triggering the algorithms relating to increased therapy described above with respect to FIGS. 12-17.

[0159] In accordance with the systems and methods described above, the skilled artisan will appreciate that 1) the patient or caregiver can adjust settings using the Algorithms via a touchscreen on the device itself or perhaps with a tablet or other remote control. This offers patient/caregiver autonomy, therapy customization, and rapid response to clinical changes. 2) a clinician can remotely manage one patient's device via the algorithms and activate the "professional algorithms" that depend on 2-way communication that connects the clinician device **1110**, **1130** to the device controller **120** of the corresponding device **224-230**. As discussed in greater detail above, a home care therapist or a physician in his office can adjust settings remotely in response to clinical changes via the clinician device **1110**, **1130** either 1) reported by the patient/caregiver or 2) triggered by data received from the device as it communicates with the managing clinical center. Accordingly, as will be appreciated by those skilled in the art, the

preceding embodiments enable a clinician to manage the device **224-230** remotely without entering a patient's clinical environment, avoiding exposure to aerosolized pathogens like Covid-19. In some circumstances, the clinician device **1110**, **1130** may be a tablet or smartphone in electronic communication with the device **224-230**, just outside the patient's room enabling close, but still remote management of the device **224-230**. Still further, the skilled artisan will appreciate that a group of clinicians may manage many patients remotely, allowing a limited number of experts to manage many patients remotely, preventing a clinician from perhaps inadvertently exposing one patient to a pathogen acquired from another patient, reducing the amount of personal protective equipment used per patient, and the like. It will be understood that this is especially relevant in the case of a pandemic, such as Covid-19, in which care is impaired by the shortage of skilled clinicians with the knowledge to manage devices that deliver assisted ventilation and mucus clearance.

[0160] In the foregoing description, numerous specific details, examples, and scenarios are set forth in order to provide a more thorough understanding of the present disclosure. It will be appreciated, however, that embodiments of the disclosure may be practiced without such specific details. Further, such examples and scenarios are provided for illustration, and are not intended to limit the disclosure in any way. Those of ordinary skill in the art, with the included descriptions, should be able to implement appropriate functionality without undue experimentation.

[0161] References in the specification to "an embodiment," etc., indicate that the embodiment described may include a particular feature, structure, or characteristic, but every embodiment may not necessarily include the particular feature, structure, or characteristic. Such phrases are not necessarily referring to the same embodiment. Further, when a particular feature, structure, or characteristic is described in connection with an embodiment, it is believed to be within the knowledge of one skilled in the art to effect such feature, structure, or characteristic in connection with other embodiments whether or not explicitly indicated.

[0162] Embodiments in accordance with the disclosure may be implemented in hardware, firmware, software, or any combination thereof. Embodiments may also be implemented as instructions stored using one or more machine-readable media, which

may be read and executed by one or more processors. A machine-readable medium may include any mechanism for storing or transmitting information in a form readable by a machine (e.g., a computing device or a "virtual machine" running on one or more computing devices). For example, a machine-readable medium may include any suitable form of volatile or non-volatile memory.

[0163] In the drawings, specific arrangements or orderings of schematic elements may be shown for ease of description. However, the specific ordering or arrangement of such elements is not meant to imply that a particular order or sequence of processing, or separation of processes, is required in all embodiments.

[0164] In general, schematic elements used to represent instruction blocks or modules may be implemented using any suitable form of machine-readable instruction, and each such instruction may be implemented using any suitable programming language, library, application programming interface (API), and/or other software development tools or frameworks. Similarly, schematic elements used to represent data or information may be implemented using any suitable electronic arrangement or data structure. Further, some connections, relationships or associations between elements may be simplified or not shown in the drawings so as not to obscure the disclosure.

[0165] This disclosure is to be considered as exemplary and not restrictive in character, and all changes and modifications that come within the spirit of the disclosure are desired to be protected. For example, while aspects of the present disclosure may be described in connection with particular types and features of respiratory therapy devices, it should be understood that the various aspects are applicable to other types and features of such devices.

CLAIMS:

1. A combination respiratory therapy device comprising:
 - a blower for providing negatively pressurized air to a mouthpiece coupled to a patient's airway;
 - an air pulse generator configured for delivering air pulses to at least one of a garment worn by the patient or a nasal interface worn by the patient;
 - a network interface in communication with an associated computer network;
 - a controller including a processor in communication with memory storing instructions which are executable by the processor to execute a combined respiratory therapy prescription, the combined respiratory therapy prescription defining a plurality of different therapy sessions to be performed by the combination respiratory therapy device over a period of time, each of the plurality of different therapy sessions comprising a mucus extraction therapy; and
 - a display in communication with the controller and configured to display a graphical user interface associated with at least one operation of the combination respiratory therapy device.
2. The combination respiratory therapy device of claim 1, wherein the graphical user interface further comprises a cough icon associated with a cough on demand operation of the combination respiratory therapy device and, wherein responsive to a selection of the cough icon, the controller is configured to operate the air pulse generator to provide inspiratory airflow to the patient and the blower to provide expiratory suction pressure to the patient.
3. The combination respiratory therapy device of claim 1, wherein the graphical user interface, further comprises an emergency icon associated with an emergency operation of the combination respiratory device, and wherein responsive to a selection of the emergency icon, the controller is configured to operate the blower and the air pulse generator to provide continuous lung ventilation to the patient.
4. The combination respiratory therapy device of claim 1, wherein the graphical user interface further comprises an intensify therapy icon associated with an intensify therapy

operation of the combination respiratory therapy device, and wherein responsive to a selection of the intensify therapy icon, the controller is configured to operate the combination respiratory therapy device at a preselected level.

5. The combination respiratory therapy device of claim 4, wherein the graphical user interface responsive to the selection of the intensify therapy icon displays a first level intensify icon and a second level intensify icon, wherein each of the first and second level intensify icons correspond to respective first and second therapy intensification operations of the combination respiratory therapy device.

6. The combination respiratory therapy device of claim 1, wherein the graphical user interface further comprises a customization icon associated with a customization of a therapy operation of the combination respiratory therapy device, and wherein responsive to a selection of the customization icon, the controller is configured to sequentially generate a series of screens associated with operations of the combination respiratory therapy device.

7. The combination respiratory therapy device of claim 6, wherein the graphical user interface, responsive to the selection of the customization icon, displays a cough inspiratory pressure adjustment screen including icons for low pressure, normal pressure, and high pressure.

8. The combination respiratory therapy device of claim 7, wherein the graphical user interface, responsive to the selection of at least one of the low, normal or high icons of the cough inspiratory pressure adjustment screen, displays a cough expiratory (suction) pressure adjustment screen including icons for low pressure, normal pressure, and high pressure.

9. The combination respiratory therapy device of claim 8, wherein the graphical user interface, responsive to the selection of at least one of the low, normal or high icons of the cough expiratory (suction) pressure adjustment screen, displays a cough cycle duration adjustment screen including icons for short duration, normal duration, and long duration.

10. The combination respiratory therapy device of claim 9, wherein the graphical user interface, responsive to the selection of at least one of the short, normal or long icons of the cough cycle duration adjustment screen, displays a cough sensitivity adjustment screen including icons for easy, normal, and hard.

11. The combination respiratory therapy device of claim 10, wherein the graphical user interface, responsive to the selection of at least one of the easy, normal or hard icons of the cough sensitivity adjustment screen, displays a ventilation inspiratory pressure adjustment screen including icons for low pressure, normal pressure, and high pressure.

12. The combination respiratory therapy device of claim 11, wherein the graphical user interface, responsive to the selection of at least one of the low, normal or high pressure icons of the ventilation inspiratory adjustment screen, displays a ventilation sensitivity adjustment screen including icons for easy, normal and hard.

13. The combination respiratory therapy device of claim 8, wherein the graphical user interface, responsive to the selection of at least one of the easy, normal or hard icons of the ventilation sensitivity adjustment screen, displays a ventilation duration after coughs adjustment screen including icons for short duration, normal duration, and long duration.

14. The combination respiratory therapy device of claim 13, wherein the graphical user interface, responsive to the selection of at least one of the short, normal or long icons of the ventilation duration adjustment screen, displays an oscillations adjustment screen including icons for slow/hard, normal, and fast/shallow oscillations.

15. The combination respiratory therapy device of claim 14, wherein the controller of the combination respiratory therapy stores each selection received via the graphical user interface in associated memory.

16. The combination respiratory therapy device of claim 1, wherein the graphical user interface includes a plurality of icons associated with a corresponding plurality of operations of the combination respiratory therapy device, and wherein the controller is

configured to receive a selection of at least one of the plurality of icons from at least one of a physician computing device or a therapist computing device via the network interface.

17. The combination respiratory therapy device of claim 16, wherein the at least one of the physician computing device or therapist computing device is remotely located from the combination respiratory therapy device.

18. The combination respiratory therapy device of claim 1, wherein the blower and the air pulse generator are resident on physically separate circuits.

19. The combination respiratory therapy device of claim 18, wherein each of the plurality of different therapy sessions comprises a lung ventilation therapy or a lung volume recruitment therapy substantially immediately following the mucus extraction therapy.

20. A system for problem-first device control of at least one combination respiratory therapy device, comprising:

at least one combination respiratory therapy device comprising:

a network interface in communication with an associated computer network,

a controller including a processor in communication with memory storing instructions which are executable by the processor to execute a combined respiratory therapy prescription, the combined respiratory therapy prescription defining a plurality of different therapy sessions to be performed by the combination respiratory therapy device over a period of time, and

a display in communication with the controller and configured to display a graphical user interface associated with at least one operation of the combination respiratory therapy device; and

at least one clinician computing device in communication with the at least one combination respiratory therapy device via the associated computer network and configured to control the at least one operation of the at least one combined respiratory therapy device.

21. The system of claim 20, wherein the graphical user interface further comprises a plurality of icons corresponding to a plurality of operations of the at least one combined respiratory therapy device, and wherein selection of one of the plurality of icons controls the at least one combined respiratory therapy device to perform the corresponding one of the plurality of operations.

22. The system of claim 21, wherein the selection of one of the plurality of icons is received via the display.

23. The system of claim 21, wherein the at least one clinician computing device further comprises a graphical user interface corresponding to the graphical user interface of the at least one combined respiratory therapy device.

24. The system of claim 23, wherein the selection of one of the plurality of icons is received via the graphical user interface of the at least one clinician computing device by the at least one combined respiratory therapy device over the associated computer network.

25. A method for remotely controlling at least one combination respiratory therapy device by a clinician device, the method comprising:

- receiving, at a clinician device in data communication with the at least one combination respiratory therapy device, patient data representative of at least one physiological parameter associated with a patient;

- generating on an associated display of the clinician device, a graphical representation of the received patient data;

- receiving, via the associated display, selection data corresponding to a selected therapy adjustment; and

- communicating the selected therapy adjustment to the at least one combination therapy device,

- wherein the clinician device includes a processor in communication with memory

storing instructions which are executed by the processor, causing the processor to perform the method.

26. The method of claim 25, further comprising establishing at least one of a video or an audio communication between the clinician device and the at least one combination respiratory therapy device.

27. The method of claim 25, further comprising generating, on the associated display of the clinician device, a graphical user interface associated with at least one operation of the at least one combination respiratory therapy device.

28. The method of claim 25, wherein the graphical user interface further comprises a plurality of icons corresponding to a plurality of operations of the at least one combined respiratory therapy device, and wherein selection of one of the plurality of icons controls the at least one combined respiratory therapy device to perform the corresponding one of the plurality of operations.

29. A clinician device for remotely controlling at least one combination respiratory therapy device by a clinician device, the clinician device comprising:

- a processor in communication with memory;

- a network interface in communication with the processor and configured to communicate with at least one combination respiratory therapy device via an associated computer network; and

- a display in communication with the processor and configured to display a graphical user interface associated with at least one operation of the at least one combination respiratory therapy device,

wherein the memory stores instructions which are executed by the processor, causing the processor to:

- receive, via the associated network, patient data representative of at least one physiological parameter associated with a patient of the at least one combination respiratory therapy device;

generate, on the display, a graphical representation of the received patient data;
receive, via the display, selection data corresponding to a selected therapy adjustment; and

communicate the selected therapy adjustment to the at least one combination respiratory therapy device via the associated computer network.

30. The clinician device of claim 29, wherein the memory further stores instructions for establishing at least one of a video or an audio communication between the clinician device and the at least one combination respiratory therapy device.

31. The clinician device of claim 29, wherein the graphical user interface further comprises a plurality of icons corresponding to a plurality of operations of the at least one combined respiratory therapy device, and wherein selection of one of the plurality of icons controls the at least one combined respiratory therapy device to perform the corresponding one of the plurality of operations.