

JPOXMLDOC01-seq1
SEQUENCE LISTING

<110> Sumitomo Dainippon Pharma Co., Ltd.
HEALIOS K.K.

<120> Method of producing retinal pigment epithelial cells

<130> 092513

<150> JP 2015-176896
<151> 2015-09-08

<160> 8

<170> PatentIn version 3.5

<210> 1
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> Forward primer for RPE65

<400> 1
tccccaatac aactgccact 20

<210> 2
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> Reverse primer for RPE65

<400> 2
ccttggcatt cagaatcagg 20

<210> 3
<211> 19
<212> DNA
<213> Artificial Sequence

<220>
<223> Forward primer for BEST1

<400> 3
tagaaccatc agcgccgtc 19

<210> 4
 <211> 21
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> Reverse primer for BEST1

 <400> 4
 tgagtgtagt gtgtatgttg g 21

 <210> 5
 <211> 20
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> Forward primer for CRALBP

 <400> 5
 gagggtgcaa gagaaggaca 20

 <210> 6
 <211> 20
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> Reverse primer for CRALBP

 <400> 6
 tgcagaagcc attgatttga 20

 <210> 7
 <211> 20
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> Forward primer for GAPDH

 <400> 7
 accacagtcc atgcatcac 20

 <210> 8
 <211> 20
 <212> DNA
 <213> Artificial Sequence

<220>

<223> Reverse primer for GAPDH

<400> 8

tccaccaccc tgttgctgta

20