

eolf-seql (1)
SEQUENCE LISTING

<110> ORION DIAGNOSTICA OY
<120> DETECTION OF NUCLEIC ACIDS
<130> N400592W0
<140> tbc
<141> tbc
<150> EP13194118.9
<151> 2013-11-22
<160> 48
<170> PatentIn version 3.5
<210> 1
<211> 16
<212> DNA
<213> Artificial Sequence
<220>
<223> DNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (16)..(16)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 1
agccgatgac taatgc

<210> 2
<211> 21
<212> DNA
<213> Artificial Sequence
<220>
<223> DNA PROBE

16

eolf-seql (1)

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 2
acgaaagccg atgactaatg c

21

<210> 3
<211> 16
<212> RNA
<213> Artificial Sequence

<220>
<223> RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (16)..(16)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 3
agccgaugac uaaugc

16

<210> 4
<211> 21
<212> RNA
<213> Artificial Sequence

<220>
<223> RNA PROBE

eolf-seql (1)

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 5' TETRAMETHYLRHODAMINE TAG

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 4
 acgaaagccg augacuaaug c

21

<210> 5
 <211> 16
 <212> RNA
 <213> Artificial Sequence

<220>
 <223> RNA PROBE

<220>
 <221> misc_feature
 <222> (1)..(16)
 <223> 2'-O-METHYL RNA

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 5' TETRAMETHYLRHODAMINE TAG

<220>
 <221> misc_feature
 <222> (16)..(16)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 5
 agccgaugac uaaugc

16

<210> 6
 <211> 21
 <212> RNA
 <213> Artificial Sequence

<220>

eolf-seql (1)

<223> RNA PROBE

<220>

<221> misc_feature

<222> (1)..(21)

<223> 2'-O-METHYL RNA

<220>

<221> misc_feature

<222> (1)..(1)

<223> 5' TETRAMETHYLRHODAMINE TAG

<220>

<221> misc_feature

<222> (21)..(21)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 6

acgaaagccg augacuaaug c

21

<210> 7

<211> 16

<212> RNA

<213> Artificial Sequence

<220>

<223> RNA PROBE

<220>

<221> misc_feature

<222> (1)..(16)

<223> 2'-FLUORO RNA

<220>

<221> misc_feature

<222> (1)..(1)

<223> 5' TETRAMETHYLRHODAMINE TAG

<220>

<221> misc_feature

<222> (16)..(16)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 7

agccgaugac uaaugc

16

eolf-seql (1)

<210> 8
 <211> 21
 <212> RNA
 <213> Artificial Sequence

<220>
 <223> RNA PROBE SEQUENCE

<220>
 <221> misc_feature
 <222> (1)..(21)
 <223> 2'-FLUORO RNA

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 5' TETRAMETHYLRHODAMINE TAG

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 8
 acgaaagccg augacuaaug c

21

<210> 9
 <211> 12
 <212> RNA
 <213> Artificial Sequence

<220>
 <223> LNA PROBE

<220>
 <221> misc_feature
 <222> (1)..(12)
 <223> LOCKED NUCLEIC ACID

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 5' TETRAMETHYLRHODAMINE TAG

eolf-seql (1)

<220>
 <221> misc_feature
 <222> (12)..(12)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 9
 aguugaugug ua 12

<210> 10
 <211> 14
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> PNA PROBE

<220>
 <221> misc_feature
 <222> (1)..(12)
 <223> PEPTIDE NUCLEIC ACID

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> TETRAMETHYLRHODAMINE TAG

<220>
 <221> misc_feature
 <222> (12)..(12)
 <223> 3' T CONJUGATED TO LYSINE AND BLACK HOLE QUENCHER 2 TAG

<400> 10
 agttgatgtg tact 14

<210> 11
 <211> 21
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> DNA PRIMER SEQUENCE

<400> 11
 aacaagaagg cgtactcgac c 21

eolf-seql (1)

<210> 12
 <211> 20
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> DNA PRIMER SEQUENCE

<400> 12
 agttgatgtg tactgagatc 20

<210> 13
 <211> 24
 <212> RNA
 <213> Artificial Sequence

<220>
 <223> RNA PRIMER

<220>
 <221> misc_feature
 <222> (1)..(24)
 <223> 2'-FLUORO RNA

<400> 13
 aguugaugug uacugacuga gauc 24

<210> 14
 <211> 24
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> MIXED DNA/RNA PRIMER

<220>
 <221> misc_feature
 <222> (1)..(20)
 <223> 2'-FLUORO RNA

<220>
 <221> misc_feature

eolf-seql (1)

<222> (21)..(24)

<223> DNA

<400> 14

aguugaugug uacugacuga gatc

24

<210> 15

<211> 61

<212> DNA

<213> Artificial Sequence

<220>

<223> STRAND INVASION OLIGONUCLEOTIDE

<220>

<221> misc_feature

<222> (48)..(60)

<223> 2'-O-METHYL RNA

<220>

<221> misc_feature

<222> (61)..(61)

<223> INVERTED dTTP

<400> 15

ttgtccatag actgctcgac ctgatacacg ttatcggtcca tacggatucg ggaucaucaua

60

t

61

<210> 16

<211> 16

<212> DNA

<213> Artificial Sequence

<220>

<223> DNA PROBE SEQUENCE

<220>

<221> misc_feature

<222> (1)..(1)

<223> 6-CARBOXYL-X-RHODAMINE

<220>

<221> misc_feature

eolf-seql (1)

<222> (16)..(16)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 16

agttgatgtg tactca

16

<210> 17

<211> 16

<212> RNA

<213> Artificial Sequence

<220>

<223> RNA PROBE SEQUENCE

<220>

<221> misc_feature

<222> (1)..(1)

<223> 6-CARBOXYL-X-RHODAMINE

<220>

<221> misc_feature

<222> (1)..(16)

<223> 2'-FLUORO RNA

<220>

<221> misc_feature

<222> (16)..(16)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 17

aguugaugug uacuga

16

<210> 18

<211> 12

<212> RNA

<213> Artificial Sequence

<220>

<223> LNA PROBE SEQUENCE

<220>

<221> misc_feature

<222> (1)..(12)

<223> LOCKED NUCLEIC ACID

eolf-seql (1)

<220>
<221> misc_feature
<222> (1)..(1)
<223> 6-CARBOXYL-X-RHODAMINE

<220>
<221> misc_feature
<222> (12)..(12)
<223> DABCYL

<400> 18
aguugaugug ua

12

<210> 19
<211> 16
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE/PRIMER

<220>
<221> misc_feature
<222> (1)..(10)
<223> 2'-FLUORO RNA

<220>
<221> misc_feature
<222> (1)..(1)
<223> IOWA BLACK FQ

<220>
<221> misc_feature
<222> (11)..(11)
<223> FLUORESCCEIN LINKED TO THE 5' POSITION OF THYMIDINE

<220>
<221> misc_feature
<222> (12)..(12)
<223> 2'-FLUORO RNA

<220>
<221> misc_feature
<222> (13)..(16)
<223> DNA

eolf-seql (1)

<400> 19
aguugaugug tactga 16

<210> 20
<211> 72
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA TARGET NUCLEIC ACID SEQUENCE

<400> 20
aacaagaagg cgtactcgac ctgatacacg ttatcgtcca tacggattcg ggatctcagt 60
acacatcaac tg 72

<210> 21
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA PRIMER SEQUENCE

<400> 21
aaccaaagtg gagtgttaca a 21

<210> 22
<211> 16
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA PRIMER SEQUENCE

<400> 22
tggggcaaaa tattta 16

<210> 23
<211> 57
<212> DNA
<213> Artificial Sequence

eolf-seql (1)

<220>

<223> STRAND INVASION OLIGONUCLEOTIDE

<220>

<221> misc_feature

<222> (44)..(56)

<223> 2'-O-METHYL RNA

<220>

<221> misc_feature

<222> (57)..(57)

<223> INVERTED dTTP

<400> 23

tcctcctgta cctcgttaca aacaggtgta tttagtagacag aagauggauu uaaaauat

57

<210> 24

<211> 16

<212> RNA

<213> Artificial Sequence

<220>

<223> RNA PROBE SEQUENCE

<220>

<221> misc_feature

<222> (1)..(1)

<223> 6-CARBOXYL-X-RHODAMINE

<220>

<221> misc_feature

<222> (1)..(16)

<223> 2'-FLUORO RNA

<220>

<221> misc_feature

<222> (16)..(16)

<223> IOWA BLACK FQ

<400> 24

uggggcaaaa uauuuua

16

<210> 25

<211> 15

eolf-seql (1)

<212> RNA
<213> Artificial Sequence

<220>
<223> RNA PROBE SEQUENCE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' CYANINE 5 TAG

<220>
<221> misc_feature
<222> (1)..(15)
<223> 2'-FLUORO RNA

<220>
<221> misc_feature
<222> (15)..(15)
<223> 3' IOWA BLACK RQ TAG

<400> 25
uggggcaaaa uauuu 15

<210> 26
<211> 66
<212> DNA
<213> Artificial Sequence

<220>
<223> TARGET SEQUENCE

<400> 26
aaccaaagtg gagtggtaca aacaggtgta tttagtagag aagatggatt taaatatttt 60

gccccca 66

<210> 27
<211> 18
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA PRIMER SEQUENCE

eolf-seql (1)

<400> 27
ctggcgatat tgggtgttt 18

<210> 28
<211> 16
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA PRIMER SEQUENCE

<400> 28
accgcaggaa acgttg 16

<210> 29
<211> 63
<212> DNA
<213> Artificial Sequence

<220>
<223> STRAND INVASION OLIGONUCLEOTIDE

<220>
<221> misc_feature
<222> (47)..(62)
<223> 2'-O-METHYL RNA

<220>
<221> misc_feature
<222> (63)..(63)
<223> INVERTED dTTP

<400> 29
tcctcctctt cctttgttta tggggtcgtt ctacattgac agaatccuca guuuuucaac 60

gat 63

<210> 30
<211> 13
<212> RNA
<213> Artificial Sequence

<220>
<223> RNA PROBE SEQUENCE

eolf-seql (1)

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 6-CARBOXYL-X-RHODAMINE TAG

<220>
 <221> misc_feature
 <222> (1)..(13)
 <223> 2'-FLUORO RNA

<220>
 <221> misc_feature
 <222> (13)..(13)
 <223> IOWA BLACK FQ TAG

<400> 30
 accgcaggaa acg 13

<210> 31
 <211> 86
 <212> DNA
 <213> Salmonella typhimurium

<400> 31
 tgaatatcgt actggcgata ttggtgttta tggggtcggt ctacattgac agaatcctca 60
 gtttttcaac gtttcctgcg gtactg 86

<210> 32
 <211> 16
 <212> RNA
 <213> Artificial Sequence

<220>
 <223> RNA PROBE SEQUENCE

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 6-CARBOXYL-X-RHODAMINE

<220>
 <221> misc_feature

eolf-seql (1)

<222> (1)..(16)
 <223> 2'-O-METHYL RNA

<220>
 <221> misc_feature
 <222> (16)..(16)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 32
 aguugaugug uacuca 16

<210> 33
 <211> 21
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> DNA PRIMER SEQUENCE

<400> 33
 gtgtacagag catttaagat t 21

<210> 34
 <211> 21
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> TARGET DNA SEQUENCE

<400> 34
 gcattagtca tcggctttcg t 21

<210> 35
 <211> 21
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> DNA/RNA PROBE

<220>
 <221> misc_feature
 <222> (1)..(1)

eolf-seql (1)

<223> 5' TETRAMETHYLRHODAMINE TAG

<220>

<221> misc_feature

<222> (1)..(1)

<223> RNA

<220>

<221> misc_feature

<222> (7)..(7)

<223> RNA

<220>

<221> misc_feature

<222> (11)..(11)

<223> RNA

<220>

<221> misc_feature

<222> (15)..(15)

<223> RNA

<220>

<221> misc_feature

<222> (19)..(19)

<223> RNA

<220>

<221> misc_feature

<222> (21)..(21)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 35

acgaaagccg atgactaaug c

21

<210> 36

<211> 22

<212> DNA

<213> Artificial Sequence

<220>

<223> DNA/RNA PROBE

<220>

<221> misc_feature

<222> (1)..(1)

eolf-seql (1)

<223> 5' TETRAMETHYLRHODAMINE TAG

<220>

<221> misc_feature

<222> (1)..(1)

<223> RNA

<220>

<221> misc_feature

<222> (3)..(3)

<223> RNA

<220>

<221> misc_feature

<222> (5)..(5)

<223> RNA

<220>

<221> misc_feature

<222> (7)..(7)

<223> RNA

<220>

<221> misc_feature

<222> (9)..(9)

<223> RNA

<220>

<221> misc_feature

<222> (11)..(11)

<223> RNA

<220>

<221> misc_feature

<222> (13)..(13)

<223> RNA

<220>

<221> misc_feature

<222> (15)..(15)

<223> RNA

<220>

<221> misc_feature

<222> (17)..(17)

<223> RNA

<220>

eolf-seql (1)

<221> misc_feature
<222> (19)..(19)
<223> RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 36
tacgaaagcc gatgactaau gc

22

<210> 37
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (1)..(3)
<223> RNA

<220>
<221> misc_feature
<222> (5)..(7)
<223> RNA

<220>
<221> misc_feature
<222> (9)..(11)
<223> RNA

<220>

eolf-seql (1)

<221> misc_feature
<222> (13)..(15)
<223> RNA

<220>
<221> misc_feature
<222> (17)..(19)
<223> RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 37
acgaaagccg atgactaaug c

21

<210> 38
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (7)..(7)
<223> 2' FLUORO RNA

<220>

eolf-seql (1)

<221> misc_feature
<222> (11)..(11)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (15)..(15)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (19)..(19)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 38
acgaaagccg atgactaaug c

21

<210> 39
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (1)..(1)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (3)..(3)
<223> 2' FLUORO RNA

<220>

eolf-seql (1)

<221> misc_feature
<222> (5)..(5)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (7)..(7)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (9)..(9)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (11)..(11)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (13)..(13)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (15)..(15)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (17)..(17)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (19)..(19)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

eolf-seql (1)

<400> 39
acgaaagccg atgactaaug c

21

<210> 40
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (1)..(3)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (5)..(7)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (9)..(11)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (13)..(15)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (17)..(19)
<223> 2' FLUORO RNA

<220>
<221> misc_feature
<222> (21)..(21)
<223> 2' FLUORO RNA

eolf-seql (1)

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 40
acgaaagccg atgactaaug c

21

<210> 41
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' TETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (1)..(1)
<223> 2' O METHYL RNA

<220>
<221> misc_feature
<222> (7)..(7)
<223> 2' O METHYL RNA

<220>
<221> misc_feature
<222> (11)..(11)
<223> 2' O METHYL RNA

<220>
<221> misc_feature
<222> (15)..(15)
<223> 2' O METHYL RNA

<220>
<221> misc_feature
<222> (19)..(19)
<223> 2' O METHYL RNA

eolf-seql (1)

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 41
 acgaaagccg atgactaaug c

21

<210> 42
 <211> 21
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> RNA PROBE

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 5' TETRAMETHYLRHODAMINE TAG

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (3)..(3)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (5)..(5)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (7)..(7)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (9)..(9)
 <223> 2' O METHYL RNA

eolf-seql (1)

<220>
 <221> misc_feature
 <222> (11)..(11)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (13)..(13)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (15)..(15)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (17)..(17)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (19)..(19)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 42
 acgaaagccg atgactaaug c

21

<210> 43
 <211> 21
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> DNA/RNA PROBE

eolf-seql (1)

<220>
 <221> misc_feature
 <222> (1)..(1)
 <223> 5' TETRAMETHYLRHODAMINE TAG

<220>
 <221> misc_feature
 <222> (1)..(3)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (5)..(7)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (9)..(11)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (13)..(15)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (17)..(19)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 2' O METHYL RNA

<220>
 <221> misc_feature
 <222> (21)..(21)
 <223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 43
 acgaaagccg atgactaaug c

21

<210> 44
 <211> 20
 <212> DNA

<213> Artificial Sequence

<220>

<223> TARGET DNA SEQUENCE

<400> 44

gatctcagta cacatcaact

20

<210> 45

<211> 21

<212> DNA

<213> Artificial Sequence

<220>

<223> DNA/RNA PROBE

<220>

<221> misc_feature

<222> (1)..(1)

<223> 5' CARBOYTETRAMETHYLRHODAMINE TAG

<220>

<221> misc_feature

<222> (7)..(7)

<223> RNA BASE

<220>

<221> misc_feature

<222> (21)..(21)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 45

acgaaagccg atgactaatg c

21

<210> 46

<211> 21

<212> DNA

<213> Artificial Sequence

<220>

<223> DNA/RNA PROBE

<220>

<221> misc_feature

eolf-seql (1)

<222> (1)..(1)
<223> 56-TAMN TAG

<220>
<221> misc_feature
<222> (7)..(7)
<223> RNA BASE

<220>
<221> misc_feature
<222> (19)..(19)
<223> RNA BASE

<220>
<221> misc_feature
<222> (21)..(21)
<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 46
acgaaagccg atgactaaug c

21

<210> 47
<211> 21
<212> DNA
<213> Artificial Sequence

<220>
<223> DNA/RNA PROBE

<220>
<221> misc_feature
<222> (1)..(1)
<223> 5' CARBOXYTETRAMETHYLRHODAMINE TAG

<220>
<221> misc_feature
<222> (7)..(7)
<223> RNA BASE

<220>
<221> misc_feature
<222> (15)..(15)
<223> RNA BASE

<220>
<221> misc_feature

eolf-seql (1)

<222> (19)..(19)

<223> RNA BASE

<220>

<221> misc_feature

<222> (21)..(21)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 47

acgaaagccg atgactaaug c

21

<210> 48

<211> 21

<212> DNA

<213> Artificial Sequence

<220>

<223> DNA/RNA PROBE

<220>

<221> misc_feature

<222> (1)..(1)

<223> 5' CARBOXYTETRAMETHYLRHODAMINE TAG

<220>

<221> misc_feature

<222> (1)..(1)

<223> RNA BASE

<220>

<221> misc_feature

<222> (7)..(7)

<223> RNA BASE

<220>

<221> misc_feature

<222> (15)..(15)

<223> RNA BASE

<220>

<221> misc_feature

<222> (19)..(19)

<223> RNA BASE

<220>

<221> misc_feature

eolf-seql (1)

<222> (21)..(21)

<223> 3' BLACK HOLE QUENCHER 2 TAG

<400> 48

acgaaagccg atgactaaug c

21