



(51) International Patent Classification:

A61K 39/23 (2006.01) C12N 7/02 (2006.01)
B01D 15/38 (2006.01)

(21) International Application Number:

PCT/US2018/067627

(22) International Filing Date:

27 December 2018 (27.12.2018)

(25) Filing Language:

English

(26) Publication Language:

English

(30) Priority Data:

62/611,709 29 December 2017 (29.12.2017) US

(71) Applicants: **BAXALTA INCORPORATED** [US/US];
1200 Lakeside Drive, Bannockburn, IL 60015 (US). **BAX-
ALTA GMBH** [CH/CH]; Zaehlerweg 4, 6300 Zug (CH).

(72) Inventors: **FIEDLER, Christian**; Raaberbahngasse 18/36,
1100 Vienna (AT). **HASSLACHER, Meinhard**; Vor-
gartenstrasse 221/1/7, 1020 Vienna (AT). **KOEHN, Jad-
ranka**; Industriestrasse 8, 88471 Laupheim (DE).

(74) Agent: **RUBE, Daniel, A.**; Troutman Sanders LLP, 875
Third Avenue, New York, NY 10022 (US).

(81) Designated States (unless otherwise indicated, for every
kind of national protection available): AE, AG, AL, AM,
AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ,
CA, CH, CL, CN, CO, CR, CU, CZ, DE, DJ, DK, DM, DO,
DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN,
HR, HU, ID, IL, IN, IR, IS, JO, JP, KE, KG, KH, KN, KP,
KR, KW, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME,
MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ,
OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA,
SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN,
TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every
kind of regional protection available): ARIPO (BW, GH,
GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ,
UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ,
TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK,
EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV,
MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM,
TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW,
KM, ML, MR, NE, SN, TD, TG).

(54) Title: ADENO-ASSOCIATED VIRUS PURIFICATION METHODS




FIG. 1

(57) Abstract: Provided herein are methods of producing an adeno-associated virus (AAV) product and methods of purifying adeno-associated virus. AAV is loaded onto an affinity resin, wash steps are undertaken, and AAV is eluted from the affinity resin. Various buffers are disclosed for use in the wash steps and elution.

[Continued on next page]

Published:

— *with international search report (Art. 21(3))*