

TRANSGENIC NON-HUMAN ANIMAL AND USES THEREOF

5 FIELD OF THE INVENTION

The present invention relates generally to transgene constructs, transgenic non-human animals comprising transgene constructs, methods of making and methods of using the transgenic non-human animals comprising transgene constructs. An
10 embodiment of the invention relates to methods of assaying for GPCR ligands non-invasively in whole animals, tissue slices, or in native cells using a transgenic model containing a bioluminescent transgene reporter system that is responsive to pathway modulation following ligand binding to GPCR receptors.

15 BACKGROUND OF THE INVENTION

In drug development, attrition rates are high with only one in five compounds making it through development to Food and Drug Administration approval (FDA) (DiMasi, JA, et al, J Health Econ 22,151-185, 2003). Moreover, despite dramatically increased
20 investment, the rate of introduction of novel drugs has remained relatively constant over the past 30 years, with only two to three advances in new drug classes per year eventually making it to market (Lindsay MA, Nature Rev Drug Discovery, 2, 831-838, 2003).

25 Molecular and functional imaging applied to the initial stages of drug development can provide evidence of biological activity and confirm the putative drug having an effect on its intended target. Thus, there is considerable expectation that investment in molecular imaging technology will enhance drug development (Rudin M, Progress in Drug Res vol 62). The advantage of molecular imaging techniques over more
30 conventional readouts is that they can be performed in the intact organism with sufficient spatial and temporal resolution for studying biological processes in vivo. Furthermore, it allows a repetitive, non-invasive, uniform and relatively automated study of the same biological model at different time points, thus increasing the

statistical power of longitudinal studies plus reducing the number of animals required and thereby reducing cost of drug development.

Molecular Imaging

5

Molecular imaging refers to the convergence of approaches from various disciplines (cell and molecular biology, chemistry, medicine, pharmacology, physics, bioinformatics and engineering) to exploit and integrate imaging techniques in the evaluation of specific molecular processes at the cellular and sub-cellular levels in living organism. (Massoud T.F., Genes Dev. 17:545-580, 2003)

10

The advent of genetic engineering has brought about major changes to applied science, including for example the drug discovery pipeline. In the same way, the development and exploitation of animal imaging procedures is providing new means for pre-clinical studies (Maggie A. and Ciana P., Nat.Rev.Drug Discov. 4, 249-255, 2005). Animal models traditionally have been cumbersome because of the difficulty in quantifying physiological events in real-time. Over the years new imaging methods have been developed to overcome this difficulty, such as magnetic resonance imaging (MRI) and positron emission tomography (PET). More recently bioluminescence imaging based on in vivo expression of luciferase, the light-emitting enzyme of the firefly, has been used for non-invasive detection.

15

20

Molecular Imaging: Bioluminescence

25

In vivo bioluminescent imaging (BLI) is a sensitive tool that is based on detection of light emission from cells or tissues. The utility of reporter gene technology makes it possible to analyze specific cellular and biological processes in a living animal through in vivo imaging methods. Bioluminescence, the enzymatic generation of visible light by a living organism, is a naturally occurring phenomenon in many non-mammalian species (Contag, C.H., et al, Mol. Microbiol. 18:593-603, 1995).

30

Luciferases are enzymes that catalyze the oxidation of a substrate to release photons of light (Greer L.F., III, Luminescence 17:43-74, 2002). Bioluminescence from the North American firefly is the most widely studied. The firefly luciferase gene (luc)

expression produces the enzyme luciferase which converts the substrate D-luciferin to non-reactive oxyluciferin, resulting in green light emission at 562 nm. Because mammalian tissues do not naturally emit bioluminescence, in vivo BLI has considerable appeal because images can be generated with very little background signal.

BLI requires genetic engineering of cells or tissues with an expression cassette consisting of the bioluminescent reporter gene under the control of a selected gene promoter constitutively driving the light reporter (Figure 3). In order to induce light production, the substrates such as luciferin are administered by intracerebroventricular (icv), intravenous (iv), intraperitoneal (ip) or subcutaneous (sq) injection.

The light emitted by luciferase is able to penetrate tissue depths of several millimeters to centimeters; however photon intensity decreases 10 fold for each centimeter of tissue depth (Contag, C.H., et al, Mol. Microbiol. 18:593-603, 1995). Sensitive light-detecting instruments must be used to detect bioluminescence in vivo. The detectors measure the number of photons emitted per unit area. Low levels of light at wavelengths between 400 and 1000 nm can be detected with charge coupled device cameras that convert the light photons that strike silicon wafers into electrons (Spibey CP et al electrophoresis 22:829-836, 2001). The software is able to convert electron signals into a two-dimensional image. The software is also able to quantify the intensity of the emitted light (number of emitted photons striking the detectors) and convert these numerical values into a pseudocolor graphic or grayscale image (Figure 2A and 2B). The actual data is measured in photons, but the pseudocolor graphic enables rapid visual interpretation. Quantitative measurements within a region of interest may be necessary for more subtle differences. The use of cooled charge coupled device (CCD) cameras reduces the thermal noise and a light-tight box allows luciferase-produced light to be optimally visualized and quantified (Contag C.H. and Bachmann, M.H., Annu. Rev. Biomed. Eng. 4:235-260, 2002).

It is useful to have the luciferase image superimposed on another type of image such as an autograph or radiograph for anatomical location of the emission signal (Figure 2B). The software superimposes images for visualization and interpretation.

5 By combining animal engineering with molecular imaging techniques, it has become possible to conduct dynamic studies on specific molecular processes in living animals. This approach could potentially impact on pre-clinical protocols thus widely changing all aspects of medicine (Maggie A. Trends Pharmacolo. Sci 25, 337-342, 2004)

10 **G-protein coupled receptors (GPCRs) – GPCRs as drug targets**

GPCRs constitute a large super family of cell surface receptors that are classified into more than 100 subfamilies on the basis of their shared topological structure; GPCRs
15 are also referred to as seven transmembrane (7TM) receptors. GPCRs are the most frequently addressed drug targets in the pharmaceutical industry. Approximately 30% of all marketed prescription drugs target GPCRs, which makes this protein family pharmaceutically the most successful target class (Jacoby, E; Chem. Med. Chem., 1: 761-782, 2006).

20 The interaction between GPCRs and their extracellular ligands has proven to be an attractive point of interference for therapeutic agents. For this reason, the pharmaceutical industry has developed biochemical drug discovery assays to investigate these ligand-GPCR interactions. Interaction of an activated GPCR with a
25 heterotrimeric G-protein catalyzes the exchange of guanosine diphosphate (GDP) by guanosine triphosphate (GTP) enabling the interaction with several downstream effectors (Cabrera-Vera T.M., Endocr. Rev. 24:765-781, 2003). Signaling downstream is dependent on the G-alpha isoform that is preferred by the GPCR of interest. Proteins of the G-alpha_{q/11} family stimulate phospholipase C (PLC), while
30 representatives of the G-alpha_{i/o} and G-alpha_s families mostly modulate adenylate cyclase (AC) activity. If the GPCR of interest signals via PLC, then the most broadly applied reporter based technique to measure GPCR activation is a calcium (Ca⁺²) release assay, either measured in a fluorescent format using Ca⁺²-sensitive

fluorophores (Sullivan E, *Methods Mol. Biol.* 114:125-133, 1999) or in a luminescent format using aequorin and a chemiluminescent substrate (Dupriez V.J., *Receptors Channels* 8: 319-330, 2002). If the GPCR of interest signals via AC, then cytosolic cyclic adenosine monophosphate (cAMP) content may be determined using various
5 detection technologies (Gabriel D. *Assay Drug Dev. Technol.* 1:291-303, 2003)

GPCR reporter based assays have been extensively used in current drug discovery programs. Typically, GPCR reporters have been introduced into cell based systems to support in vitro high-throughput screening (HTS) of large pharmaceutical libraries
10 to identify ligands or compounds that activate or modulate the specific GPCR.

Secondary and follow-up cell based assays confirm and refine any "hits" identified in HTS against a specific GPCR; but again, these assays rely on recombinant DNA methods to introduce a cloned GPCR into a transformed cell type. While transformed cell types have excellent proliferative capacity to support large screening programs,
15 they often display aberrant genetic and functional characteristics and consequently significant attrition of putative "hits" from HTS is encountered using this paradigm.

For several years, bioluminescence-based reporter gene assays have been employed to measure functional activity of GPCRs (Hill, S.J. *Curr. Opin. Pharmacol.* 1: 526-532, 2001). This assay format is very sensitive owing to the low signal
20 background of the bioluminescent readout and the signal amplifications steps between GPCR activation and the cumulative reporter gene expression.

A cAMP response element (CRE) in the promoter of the reporter gene enables the
25 specific monitoring of G protein dependent signaling. When a ligand binds to the GPCR it causes a conformational change in the GPCR which allows it to activate an associated G-protein. The enzyme adenylylase is a cellular protein that can be regulated by G-proteins. Adenylylase activity is either activated or inhibited when it binds to a subunit of the activated G protein. Signal transduction depends on
30 the type of G protein. Adenylylase acts to either increase or decrease cAMP production in the cell. The cAMP produced is a second messenger in cellular metabolism and is an allosteric activator to protein kinase A (PKA). When there is no cAMP, the PKA complex is inactive. When cAMP binds to the regulatory subunits of

PKA, their conformation is altered, causing the dissociation of the regulatory subunits, which activates protein kinase A and allows further biological effects. PKA then phosphorylates and activates the transcription factor CREB. CREB binds to certain DNA sequences called cAMP response elements (CRE) and thereby increases or decreases transcription, and thus the expression, of certain genes, such as the luciferase reporter gene.

The CreLuc transgene is designed to assay activation of all three major GPCRs either directly through the cAMP intracellular signaling pathway or indirectly through signaling via PLC. Because any one cell type contains many different types of GPCRs on their cell surface, (thus any cell would have GPCRs signaling via G- $\alpha_{q/11}$, G- $\alpha_{i/o}$ and G- α_s occurring simultaneously within a cell) conventional wisdom would suggest that it would be improbable that a transgene such as CreLuc would be specific enough to discriminate any one specific GPCR ligand. However, we demonstrate here the CreLuc transgene is able to discriminate GPCR ligands. We predict that the bioluminescent signal for the luciferase reporter in cells, tissues slices, or the whole animal will be increased with forskolin and be modulated by ligands for Gs, Gq or Gi receptors. Table 1 shows the anticipated effect that GPCR activation/inhibition will have on the CreLuc reporter system upon binding to a GPCR ligand. Further, we show data that our novel CreLuc reporter system can discriminate different classes of GPCR ligands and that such a reporter system is applicable for identifying novel GPCR ligands when used in cells, tissue slices and the whole animal.

Table 1. Predicted change in bioluminescent signal from the CreLuc reporter upon ligand binding to specific GPCRs

Receptor Type	Agonist	Antagonist; Inverse Agonist
Gs; Gq	Increase	Decrease
Gi	Decrease	Increase

A GPCR Bioimaging Reporter Transgenic Model

Significant attrition of potential drug candidates in the current drug discovery paradigm is encountered in the phase transition from cell-based reporter assays to in vivo models. Numerous in vivo models are available that recapitulate either all or part of a particular human disease. Demonstrating lead compound activity in these models is a significant milestone for progression of new chemical GPCR drugs. Animal disease models typically require a large number of animals and time to allow for the development of their phenotype and an accurate assay of the candidate compound's impact on altering the disease outcome. Following in vitro testing, the next level of testing a drug candidate in a complex system is using in vivo testing or in vivo models of disease states which are mechanistic based. Failures to alter the induced disease outcomes are poorly understood but yet result in the large attrition rates of candidate compounds in the drug development pipeline.

A transgenic model containing a GPCR ligand binding and activation reporter assay would be a significant improvement in the current drug discovery paradigm for GPCRs. For instance, an embodiment of this invention describes a transgene containing the cAMP reporter assay based on a luciferase reporter (CreLuc) that is combined with molecular imaging in whole animals, tissues, or cells which would significantly accelerate GPCR ligand drug discovery (Bhaumik, S. and Gambhir, S.S., Proc. Natl. Acad. Sci. USA, 99:377-382 2002; Hasan M.T., et al., Genesis 29:116-122, 2001). As described herein, embodiments of the transgenic non-human animal of the instant invention offer the following non-limiting advantages:

1. Tissues or cell based assays have the same reporter system as in the transgenic in vivo model assay thus reducing the number of unknowns in complex intact biological systems.
2. Non-invasive imaging allows quantitative analysis of ligand or compound activity in a time-course assay in the same animal.
3. Non-invasive imaging reduces the number of animals per study and leads to greater statistical power by each animal being its own control wherein the control would be the animal assayed at time zero.

4. The transgenic animal would be a source of cells and tissues to support parallel assays done in vitro or ex vivo.
5. The assay of the transgenic animal would support the assay of ligand activity in native cell types which leads to a more realistic profile of ligand: receptor interaction.
6. The transgenic animal allows for simultaneous assessment of pharmacodynamics and pharmacokinetics of GPCR ligands.
7. The transgenic animal allows for simultaneously identification of tissue and cell-type specificity at either the organ or whole animal level.
8. The transgenic animal allows for cross breeding with other genetically altered models to reveal novel signaling pathways and their response to specific ligands.

Many transgenic animals engineered with different reporters are being employed in the study of molecular processes such as drug metabolism (Zhang W., et al. Drug Metab. Dispos. 31:1054-1064, 2003), genotoxicity (Gossen J.A., et al., Proc. Natl. Acad. Sci. USA 86:7971-7975, 1989) and the effects of toxic compounds (Sacco M.G. et al., Nat. Biotechnol. 15:1392-1397, 1997). To achieve their design goals, a GPCR reporter animal suitable for molecular imaging studies has to incorporate several elements arranged to allow both high levels of reporter expression to support a large window of bioluminescent detection as well as expression in every cell type to support broad acute in vivo assays on biodistribution of the ligand or compound under study.

The complexity and diversity of the mechanisms involved in gene expression will never allow researchers to construct genes capable in all cases of being expressed in transgenic animals in a fully predictable manner (Pinkert, C. A. (ed.) 1994. Transgenic animal technology: A laboratory handbook. Academic Press, Inc., San Diego, Calif.; Monastersky G. M. and Robl, J. M. (ed.) (1995) Strategies in Transgenic Animal Science. ASM Press. Washington D.C). Only through extensive trial and error can unique combinations of transgene structures be arrived at to deliver model design goals as required for bioimaging of GPCR reporters.

Utility of a Transgenic GPCR Reporter over Recombinant Cell Assays

As screening technology advances to the point of understanding the behaviors of individual GPCRs, it is clear that rather than being on-off switches, these receptors are acting more as microprocessors of information. This has introduced the phenomenon of functional selectivity, whereby certain ligands initiate only portions of the signaling mechanism mediated by a given receptor, which has opened new horizons for drug discovery. The need to discover new GPCR ligand relationships and quantify the effect of the drug on these complex systems to guide medicinal chemistry puts significantly higher demands on any pharmacological reporter assay. This concept drives the return to whole-system assays from the reductionist recombinant cell based screening systems. Profiling a ligand's activity with a specific GPCR or set of GPCRs in a native cellular environment is expected to improve the success rate of identifying new drugs against a key class of pharmaceutically important receptors (Kenakin TP, Nat. Rev. Drug Discov. 8,617-625, 2009). An animal model containing a bioluminescent GPCR reporter transgene is a highly desirable molecular imaging strategy to define GPCR ligand activity in an intact biologically complex system with the goal of improving drug discovery to fight human diseases.

Because activation of CRE/CREB is involved many varied biological processes, there has been considerable interest in studying the activation of CRE by using a CRE/CREB reporter expression system. Cyclic adenosine monophosphate (cAMP) is a second messenger in intracellular signal transduction following receptor activation and subsequent activation of protein kinase, thereby being involved in the regulation of many biological processes. CREB (cAMP responsive element binding protein), phosphorylated by kinase activated by cAMP, binds to the cAMP responsive element (CRE) in the promoter region of many genes and activates transcription (Shaywitz and Greenberg, Annul. Rev. Biochem., 68:821-861, 1999). Transgenic mice carrying six tandem CREs with a minimal herpes simplex virus (HSV) promoter driving beta-galactosidase expression were used to study CRE-mediated gene expression in brain slices in response to chronic antidepressant treatment (Thome J., et al., J. Neurosci. 20:4030-4036, 2000). Similarly, transgenic mice carrying four copies of rat

somatostatin gene promoter CRE fused to a thymidine kinase promoter and the luciferase gene have been used to study CRE activation in histological brain slices or homogenates (Boer et al, PloS One, May 9; 2(5):e431, 2007). However, studies to date have been hampered by the need to screen large numbers of transgenic lines to find a suitable animal model. Further, after the appropriate line has been identified, relatively low reporter expression levels require the transgenic animal be euthanized in order to measure the reporter gene, thus requiring large number of animals be used to in a single experimental paradigm.

An embodiment of the invention is the development of a transgene comprising insulator elements, response elements, promoter elements, reporter genes, and functional elements. The transgene can be quickly introduced into non-human animals because of its high rate of integration and high level of reporter gene expression, thus transgenic animals can be easily developed as models to study regulatory element activation *in vivo* (i.e., in the living animal), *in situ* (e.g., brain slices, intact whole organ) or *in vitro* (e.g., primary cells cultured from the transgenic animal, tissue homogenates).

An embodiment of the invention is a transgene comprising a CRE Luc reporter system used in transgenic non-human animals as models to quantify GPCR ligand activities through the regulation of intracellular cAMP levels *in vivo*. As a non-limiting example, we have demonstrated changes in the luciferase reporter via bioluminescence in isolated primary cells and in whole animals using general cAMP regulators. In another embodiment, activation of the reporter has been assayed and confirmed in tissue extracts using luciferase assays *ex vivo*. The response of the CRE Luc transgene has been documented in multiple mouse lines and exhibits either single or multiple tissue activation profiles. Furthermore, as non-limiting examples, we demonstrate that specific GPCR ligands activated the CRE Luc transgene in whole animals, tissue slices, and primary cells.

BRIEF SUMMARY OF THE INVENTION

In general, the invention provides transgene constructs, transgenic non-human animals comprising transgene constructs, methods of making and methods of using the transgenic non-human animals comprising transgene constructs. An embodiment of the invention provides a transgene construct comprising the CRE Luc reporter system. An embodiment of the invention is the introduction of the transgene construct comprising the CRE Luc reporter system into a non-human animal.

Since cAMP modulation is a key activation pathway for GPCRs, the invention serves as a platform for quantifying in whole animals, tissue slices, or cells the activation of a GPCR by a ligand or compound through the activation of a reporter gene wherein the reporter gene provides for a measurable bioluminescent signal, for example, metabolism of luciferin by luciferase. This invention supplies tools to improve the transition of new drug discovery entities such as ligands or compounds from cell based assays to whole animals. An embodiment of the invention uses the same reporter system in native cells which will reduce the attrition rate for new GPCR ligands while simultaneously supplying bioavailability data.

An embodiment of the invention is a transgenic non-human animal having a genome comprising a transgene comprising a first insulator element, a response element, a promoter, a bioluminescent reporter, a functional element and a second insular element.

An embodiment of the invention is a transgenic non-human animal wherein the first insulator element is selected from the group consisting of matrix attachment regions (MAR), DNase I-hypersensitive site (HS4) and inverted terminal repeats (ITR). A further embodiment of the invention is a transgenic non-human animal wherein the second insulator element is selected from the group consisting of matrix attachment element (MAR), HS4 and ITR. A further embodiment of the invention is a transgenic non-human animal wherein the first insulator element is the same as the second insulator element.

An embodiment of the invention encompasses a transgenic non-human animal wherein the response element is selected from the group consisting of cAMP response element (CRE), activator protein 1 (ASP1), glucocorticoid response element (GRE), heat shock response element (HSE), serum response element (SRE), thyroid response element (TRE) and estrogen response element (ERE). A further
5 embodiment of the invention is a transgenic non-human animal wherein the response element is repeated in tandem two to twenty-four times. A further embodiment of the invention is a transgenic non-human animal wherein the response element is repeated in tandem six times. A further embodiment of the invention is a transgenic
10 non-human animal wherein the response element is CRE, further wherein the CRE response element may be a single element or repeated two to twenty-four times.

An embodiment of the invention is a transgenic non-human animal wherein the promoter is herpes simplex virus thymidine kinase minimal (HSV TK min).

15 An embodiment of the invention is a transgenic non-human animal wherein the bioluminescent reporter is selected from the group consisting of luciferase, chloramphenicol acetyltransferase (CAT), beta-galactosidase, secreted alkaline phosphatase (SEAP), human growth hormone (HGH) and green fluorescent protein
20 (GFP).

An embodiment of the invention is a transgenic non-human wherein the functional element is human growth hormone (hGH) gene.

25 An embodiment of the invention is a transgenic non-human wherein the transgene comprises SEQ ID NO: 18.

An embodiment of the invention is a transgenic non-human wherein the transgene comprises SEQ ID NO: 19.

30 An embodiment of the invention is a cell isolated from the transgenic non-human animal or a tissue slice isolated from the transgenic non-human animal of claim 1.

An embodiment of the invention is a method of identifying a G protein-coupled receptor (GPCR) ligand, the method comprising (a) measuring an amount of bioluminescence in the transgenic non-human animal disclosed herein; (b) administering a test agent to the transgenic non-human animal; (c) measuring an amount of bioluminescence of the transgenic non-human animal at one or more time points following administration of the test agent; and (d) comparing the amount of bioluminescence measured in (a) to the amount of bioluminescence measured in (c) wherein a difference in the amount of bioluminescence in (a) compared to (c) identifies the test agent as a GPCR ligand.

An embodiment of the invention is a method of identifying a G protein-coupled receptor (GPCR) ligand, the method comprising (a) preparing a tissue slice from the transgenic non-human animal disclosed herein; (b) measuring an amount of bioluminescence in the tissue slice; (c) administering a test agent to the tissue slice; (d) measuring an amount of bioluminescence of the tissue slice at one or more time points following administration of the test agent; and (e) comparing the amount of bioluminescence measured in (b) to the amount of bioluminescence measured in (d) wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the test agent as a GPCR ligand.

An embodiment of the invention is a method of identifying a G protein-coupled receptor (GPCR) ligand, the method comprising (a) preparing a cell isolated from the transgenic non-human animal disclosed herein; (b) measuring an amount of bioluminescence in the cell; (c) administering a test agent to the cell; (d) measuring an amount of bioluminescence in the cell at one or more time points following administration of the test agent; and (e) comparing the amount of bioluminescence measured in (b) to the amount of bioluminescence measured in (d) wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the test agent as a GPCR ligand.

An embodiment of the invention is a method of monitoring GPCR function in a non-human animal, the method comprising (a) transgenically modifying a non-human animal to express a transgene comprising a first insulator element, a response

element, a promoter, a bioluminescent reporter, a functional element and a second insular element; (b) monitoring bioluminescence from the non-human animal; and (c) correlating said bioluminescence to GPCR function.

- 5 An embodiment of the invention is a method of monitoring GPCR function in a non-human animal, the method comprising (a) transgenically modifying a non-human animal to express a transgene comprising a first insulator element, a response element, a promoter, a bioluminescent reporter, a functional element and a second insular element; (b) monitoring luciferase from the non-human animal; and (c)
10 correlating said bioluminescence to GPCR function.

- An embodiment of the invention is a method of monitoring GPCR function in a non-human animal, the method comprising (a) transgenically modifying a non-human animal to express a transgene comprising a first insulator element, a response
15 element, a promoter, a bioluminescent reporter, a functional element and a second insular element; (b) manipulating the non-human animal to mimic an aspect of a disease state; (c) monitoring bioluminescence from the non-human animal; and (d) correlating said bioluminescence to GPCR function.

- 20 An embodiment of the invention is a method of monitoring GPCR function in a non-human animal, the method comprising (a) transgenically modifying a non-human animal to express a transgene comprising a first insulator element, a response element, a promoter, a bioluminescent reporter, a functional element and a second insular element; (b) manipulating the non-human animal to mimic an aspect of a
25 disease state; (c) monitoring luciferase from the non-human animal; and (d) correlating said bioluminescence to GPCR function.

- An embodiment of the invention is a method of making a non-human transgenic animal for use in monitoring GPCR function, the method comprising (a) transgenically
30 modifying a non-human animal to express a transgene comprising a first insulator element, a response element, a promoter, a bioluminescent reporter, a functional element and a second insular element; (b) measuring an amount of bioluminescence in the transgenic non-human animal of (a); (c) administering a GPCR ligand to the

transgenic non-human animal; (d) measuring an amount of bioluminescence of the transgenic non-human animal at one or more time points following administration of the GPCR ligand; and (e) comparing the amount of bioluminescence measured in (b) to the amount of bioluminescence measured in (d) wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the non-human transgenic animal for use in monitoring GPCR function.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) preparing a cell isolated from the transgenic non-human animal disclosed herein; (b) measuring an amount of bioluminescence in the cell; (c) administering a test agent to the cell; (d) measuring an amount of bioluminescence in the cell at one or more time points following administration of the test agent; and (e) comparing the amount of bioluminescence measured in (b) to the amount of bioluminescence measured in (d) wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the test agent as a GPCR ligand.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing cells isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a control to one or more receptacles; (c) administering a test agent to one or more receptacles; and (d) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle(s) comprising control compared to the amount of luciferase in the receptacle(s) comprising test agent indicates the compound as modulating a GPCR.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing cells isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a general cAMP modulator to one or more receptacles; (c) administering a test agent to one or more receptacles; and (d) measuring an amount of luciferase in the receptacles, wherein a difference in the

amount of luciferase measured in the receptacle(s) comprising only the general cAMP modulator is compared to the amount of luciferase in the receptacle(s) comprising the general cAMP modulator and the test agent indicates the compound as modulating a GPCR.

5

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing tissue slices isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a control to one or more
10 receptacles; (c) administering a test agent to one or more receptacles; and (d) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle(s) comprising control compared to the amount of luciferase in the receptacle(s) comprising test agent indicates the compound as modulating a GPCR.

15

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing tissue slices isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a general cAMP modulator to
20 one or more receptacles; (c) administering a test agent to one or more receptacles; and (d) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle(s) comprising only the general cAMP modulator is compared to the amount of luciferase in the receptacle(s) comprising the general cAMP modulator and the test agent indicates the compound
25 as modulating a GPCR.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing cells isolated from the transgenic non-human animal disclosed herein into
30 one or more receptacles; (b) administering a cell stimulator to one or more receptacles; (c) administering a test agent to one or more receptacles; and (e) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle(s) comprising cell stimulator

compared to the amount of luciferase in the receptacle(s) comprising test agent and cell stimulator indicates the compound as modulating a GPCR.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing cells isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a cell stimulator to one or more receptacles; (c) administering a general cAMP modulator to one or more receptacles; (d) administering a test agent to one or more receptacles; and (e) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle comprising the cell stimulator and general cAMP modulator is compared to the amount of luciferase in the receptacle(s) comprising the cell stimulator and general cAMP modulator and the test agent indicates the compound as modulating a GPCR.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing tissue slices isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a cell stimulator to one or more receptacles; (c) administering a test agent to one or more receptacles; and (d) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle(s) comprising cell stimulator compared to the amount of luciferase in the receptacle(s) comprising test agent and cell stimulator indicates the compound as modulating a GPCR.

An embodiment of the invention is a method of identifying a compound that modulates a G protein-coupled receptor (GPCR), the method comprising (a) providing tissue slices isolated from the transgenic non-human animal disclosed herein into one or more receptacles; (b) administering a cell stimulator to one or more receptacles; (c) administering a general cAMP modulator to one or more receptacles; (d) administering a test agent to one or more receptacles; and (e) measuring an amount of luciferase in the receptacles, wherein a difference in the amount of luciferase measured in the receptacle comprising the cell stimulator and general

cAMP modulator is compared to the amount of luciferase in the receptacle(s) comprising the cell stimulator and general cAMP modulator and the test agent indicates the compound as modulating a GPCR.

5

BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 shows the CreLuc bioimaging mouse model. This figure illustrates the intracellular activation of the CRE Luc reporter transgene by all three types of GPCRs (Gi, Gs, and Gq) either directly through the cAMP pathway or indirectly through the PLC pathway (panel A). The change in the bioluminescence of the luciferase reporter in response to forskolin induction is illustrated for the three types of GPCRs (panel B). Forskolin will increase Gs and Gq signaling and thus CreLuc bioluminescence will increase, while Gi induction will decrease the signal from the reporter. G α s activates the cAMP dependent pathway by direct stimulation of AC, G α s inhibits the production of cAMP and G α q stimulates PLC resulting in the generation of the two second messengers IP3 and DAG. Abbreviations: α , α -subunit of the G protein; β , β -subunit of the G protein; γ , γ -subunit of the G-protein; AC, adenylate cyclase; PLC, phospholipase C; PKA, protein kinase A; PKC, protein kinase C; DAG, diacylglycerol; IP3, inositol triphosphate; Ca⁺², calcium; CaMK, calcium/calmodulin protein kinase; cAMP, cyclic adenosine monophosphate CRE, cAMP response element; CREB, cAMP responsive element binding protein,

Figure 2A shows real-time in vivo bioluminescence imaging and describes the benefits of using the system. The IVIS100 (Xenogen) bioimaging instrument with a computer analyzer system allows real-time in vivo imaging utilizing the light emitted by a bioluminescent reporter gene, for instance luciferase, expressed in vivo. The software supports quantification of the signal non-invasively and longitudinally. Real-time in vivo imaging has many advantages over traditional in vivo compound testing. Traditional animal studies require individual mice at multiple treatment points while studies utilizing bioimaging models allow the same animals to be sampled at multiple time points and reused for multiple treatments. As shown in this figure, a time course of 0 hours, 2 hours, 4 hours and 8 hours would require 24 animals (n=6 per time

point) using current methodology whereas only 6 animals would be required using bioimaging technology. This results in several benefits which include: higher throughput because fewer test animals are required allowing more compounds to be tested for efficacy; greater data content and quality since temporal and spatial data can be collected from the same animal; and decreases in statistical error which improves the quality of decisions made about individual compounds.

Figure 2B shows a typical visual image of compound induction in a CreLuc transgenic mouse. Administration of isoproterenol (right panel) increases spinal cord expression of the CRE Luc reporter compared to basal levels (left panel). The bioluminescent detection is represented visually on a white-light image of the animal as a pseudocolor representation in grayscale.

Figure 3 shows a schematic representation of a transgene structure comprising multiple DNA elements to enhance expression. The schematic transgene structure comprises the following elements: an insulator element shown in this figure as matrix attachment regions (MAR) to generate position independent expression; a response element represented by CRE-cAMP repeated six times (6X CRE); a promoter element shown as a herpes simplex virus thymidine kinase minimal promoter (HSV TK min); a reporter element which is represented by a luciferase gene optimized for mammalian expression (LUC2); and a functional element depicted by human growth hormone gene with poly A tail (hGH poly A) to enhance transgene expression.

Figure 4: *In vitro* validation of CreLuc transgene vectors. Hybrid or synthetic (synth) vectors at 30, 100 or 199ng DNA were transfected into CHO cells (along with a renilla luciferase positive control vector) with Lipofectamine 2000 (Invitrogen, Carlsbad, CA, cat #11668-019). Two days later, the cells were stimulated with 3 μ M forskolin (Sigma, St Louis, MO, cat #F6886) for 4 hours and then luciferase activity was measured with Dual Glo Luciferase Assay System (Promega, Madison, WI, cat# E2920). Results show a dose dependant increase in the luciferase signal with both vectors. Higher levels of induction were achieved with the hybrid vector.

Figure 5A shows the effects of PDE inhibitors on cAMP levels in normal mice. (reproduced from Cheng JB, JPET, 280, 621-626). Balb/c mice were dosed with either vehicle or drugs as listed on the x-axis, the blood was harvested after 20 minutes and assayed for cAMP by cAMP radioimmunoassay. Both CP-80,633 and rolipram significantly increase plasma cAMP levels at 10mg/kg.

Figure 5B shows *in vivo* stimulation of cAMP in plasma. FVB/Tac females were dosed with either vehicle (1% DMSO) or drugs *i.p.*, then 30 minutes later, blood samples were collected and assayed for cAMP by ELISA (Assay Designs, Ann Arbor, MI, cat#900-163). The drugs used are as follows: 5mg/kg forskolin (F) (Sigma F6886), 5mg/kg water soluble forskolin (H20F) (Calbiochem 344273), 10mg/kg rolipram (R) (Sigma R6520) or combinations of either forskolin plus rolipram (F/R). Statistically significant increases, 14 fold, were observed with treatment of water soluble forskolin in combination with rolipram as determined by t-test. Forskolin increases cAMP levels by activating adenylylate cyclase while inhibitors of PDE4, such as rolipram, raise plasma cAMP by preventing hydrolysis of cAMP. A combination of rolipram and water soluble forskolin increases cAMP levels *in vivo* 14 fold. This combination was used to provide a large window of induction for founder screening by bioimaging, a representative study is shown in Figure 6.

Figure 6 shows the results of the initial founder induction and line selection for the CreLuc reporter mouse model. Multiple transgenic lines were screened for luciferase induction with forskolin and rolipram *in vivo* and then tissues were isolated and assayed for luciferase enzymes. Transgenic mice were bioimaged pre-dosing (basal expression levels), then the same mouse was dosed *i.p.* with 10mg/kg rolipram and 5mg/kg water soluble forskolin and bioimaged 4 hours post-dosing (induced expression). A (subline #90): forskolin/rolipram administration increased basal expression of CreLuc reporter transgene in the lungs and other tissues; B (subline #219): induction of basal expression is mainly in the gut; C(subline #44): undetectable basal expression and reporter is induced in brain plus other tissues; D (subline #28): undetectable basal expression that is increased in thymus and liver; E (subline #187): undetectable basal expression that is induced in the brain and spinal cord. As expected for a randomly integrated transgene, there was variation between lines in

basal expression, tissue distribution, and response to induction. Twenty lines were identified to have greater than 5X induction in one or more tissues. Variation in tissue profile demonstrates that single tissue (i.e. lung, liver, brain) allows imaging devoid of background tissue response while multiple tissues allows a broad compound response profile to be generated.

Figure 7 shows a general schematic of CreLuc screening assays in vivo or ex vivo.

Figure 8A shows the effects of isoproterenol (ISO) and AMN082 (AMN) on luminescence in whole brain slices from CreLuc mice (line 187).

Figure 8B shows the effects of forskolin on luminescence in whole brain slices from CreLuc mice (line 44) over time. Time is represented on the X-axis in minutes. Forskolin at 50uM or vehicle (DMSO) was added at time = 2880 marked by an arrow in the bottom panel.

Figure 9 shows a schematic representing the isolation and compound treatment of primary neuronal cells from the CreLuc mice.

Figure 10 shows Gs modulation via β -adrenergic receptor (AD β R) activation and D1 dopamine receptor (DRD1) activation. Neurons were isolated from the cortices of line 187 E18 embryos. On day three in culture, test compounds were added Forskolin 5 μ M (F), rolipram at 10 μ M (R), forskoline and rolipram in combination (F/R) isoproterenol at 10 μ M, isoproterenol and rolipram in combination (I/R); SKF82958 at 10 μ M, and SKF82958 and rolipram in combination (S/R). Data is shown as counts per second (cps)

Figure 11 shows the effects of prokineticin 2 (PROK2) peptide on luciferase expression in primary cortical neurons. Primary cortical neurons were harvested from line 187 (inducible luciferase in brain and spinal cord) on E18. The assay was run on day 3 in culture for 4 hours or 8 hours. The PROK2 peptide is added as an aqueous solution at 1 nM and 100nM. Data is shown as counts per second (cps)

Figure 12 shows the effects of prokineticin 2 (PROK2) peptide on luciferase expression in primary cortical neurons from different CreLuc lines. Primary cortical neurons were harvested from four different CreLuc lines at E18. The assays were run in triplicate at day three in culture with either 1nM or 100nM PROK2 peptide at two timepoints, 4 hours and 24 hours. BrightGlo was used for the assay and read on a TopCount. Data is shown as counts per second (cps)

Figure 13A shows the effects of the mGluR7 agonist, AMN082 on luciferase expression in primary cortical neurons. Cortical neurons were harvested from E18 embryos (line 187). The assay was run at day 3 in culture. Forskolin was used at 10µM. The agonist, AMN082 was used in combination with forskolin at 1nM, 10nM, 100nM and 1µM. The assay was read on a TopCount with Bright Glo (Promega) at 4 hours, and 8 hours. Data is shown as counts per second (cps)

Figure 13B shows the results of screening unknown compounds for the ability to modulate Gi activity in primary cortical neurons. Cortical neurons were harvested from E18 embryos (line 187). The assay was run at day 3 in culture. Forskolin was used at 10µM. AMN082 or unknown compounds A, B or C was tested in combination with forskolin at different concentration and EC50 values were calculated. The assay was read on a TopCount with Bright Glo (Promega) at 4 hours. Data is shown as counts per second (cps)

Figure 14 shows the effects of the mGluR7 agonist, AMN082 on luciferase expression in primary cortical neurons. Primary cortical neurons were isolated from E18 embryos from line 187. The assay, in triplicates, was run at day 7 in culture for 6 hours. A concentration curve for AMN082 was run in combination with 50µM forskolin and 10µM rolipram. The assay was read on a TopCount luminometer with Bright Glo substrate (Promega). Data is shown as counts per second (cps)

Figure 15A shows Gi modulation of luciferase expression in primary cortical neurons from different CreLuc lines by the CB1 agonist, CP 55,940. Primary cortical neurons were harvested from four different CreLuc lines at E18. The assays were run on day three in culture. The CB1 agonist was used at 10µM, forskolin at 5µM and rolipram at

10 μ M. Two timepoints were run, four hours and twenty-four hours. Bright Glo luciferase assay substrate was then added, and the assay read on a Topcount luminometer. Data shown is the average of the triplicates. Data is shown as counts per second (cps).

5

Figure 15B shows Gi modulation of luciferase in primary cortical neurons from CreLuc mice by the CB1 agonist, CP 55,940. Cortical neurons were isolated from E18 embryos (line 187). The assay was run on day 3 in culture. Forskolin (F) and rolipram (R) were used at 10 μ M. The agonist was added at concentrations of 10 μ M, 1
10 μ M and 100nM. The assay was read on a TopCount with BrightGlo (Promega) at 8 hours. Data is shown as counts per second (cps).

Figure 16 shows induced luciferase expression in CreLuc striatal neurons by forskolin and rolipram, and Gs agonists DRD1 and AD β R. Striatum neurons were
15 isolated from E14 embryos (line 187). The assays were run at day 4 in culture. Forskolin (F) was used at 5 μ M, rolipram (R) at 10 μ M. The Gs agonists isoproterenol (iso), dopamine (dopa) and SKF82958 (chloro) were used at 10 μ M, 3 μ M and 1 μ M. The assay was read at 5 hours with a TopCount luminometer and Bright Glo luciferase reagent (Promega). Data is shown as counts per second (cps).

20

Figure 17 shows the effects of general cAMP inducers such as forskolin (F) and rolipram (R), as well as Gs agonists on luciferase expression in whole splenocyte preps isolated from CreLuc mice. Line 64 splenocytes were stimulated for 24 hours with anti-CD3 antibody (CD), the other half were untreated (unstim). At 24 hours,
25 compounds were added to the plates for an additional 4 hours. The co-treatment of forskolin and rolipram (F/R) was 5uM forskolin and 10um rolipram. The Gs agonists used are: EX00000173A (173A) as an EP2 agonist, BW245C as a DP1 agonist and isoproterenol as an AD β R agonist. All Gs agonists were used at 10uM. The assay is run as triplicates. After 4 hours, 100ul of BrightGlo was added and the assay was
30 read on a TopCount luminometer. Data is shown as counts per second (cps).

Figure 18 shows the effects of general cAMP activation by rolipram and forskolin in T cells isolated from five different sublines of CreLuc mice. The cells were stimulated

with anti CD3 antibodies (1ug/ml). After 18 hours, 10µM rolipram and 5µM forskolin were added to the plates for an additional 4 hours. BrightGlo was added and the assay was read on the TopCount. Data is shown as luminescence (counts per second) in the top panel, and as fold increase over media only controls in the bottom panel.

Figure 19 shows the effects of Gs agonists on luciferase levels in anti CD3 stimulated CD4+ T cells isolated from CreLuc mice (line 64). The cells, 1.5×10^5 per well, were plated on 96 well white opaque plates and then stimulated with anti CD3 antibodies (1ug/ml). After 24 hours, compounds were added for an additional 4 hours. Gs agonists BW245C, EX00000173A (1734A) and isoproterenol (iso) were all used at 10µM. Forskolin (F) was added at 5uM and rolipram (R) at 10uM. BrightGlo was added and the assay was read on the TopCount. Data is shown as counts per second (cps).

Figure 20 shows the effects of general cAMP activation by rolipram and forskolin in B cells isolated from two different sublines of CreLuc mice. Cells were plated 2.0×10^5 per well on 96 well white opaque plates and stimulated with 10ng/ml lipopolysaccharide (LPS). After 18 hours, 10µM rolipram and 5µM forskolin were added to the plates for an additional 4 hours. BrightGlo was added and the assay was read on the TopCount. Data is shown as luminescence (counts per second), and as fold increase over media only controls.

Figure 21 shows the effects of Gs agonists on luciferase levels in LPS stimulated B220+ B cells isolated from CreLuc mice. The cells were plated, 2.0×10^5 per well on 96 well white opaque plates and then stimulated with 10ng/ml lipopolysaccharide (LPS). After 24 hours, compounds were added for an additional 4 hours. Gs agonists BW245C, EX00000173A (1734A) and isoproterenol (iso) were all used at 10µM. Forskolin (F) was added at 5uM and rolipram (R) at 10uM. BrightGlo (Promega, Madison, WI, cat#E2610) was added and the assay was read on the TopCount. Data is shown as counts per second (cps).

Figure 22 shows the Induced luciferase expression in isolated microglia (line 64) by the general cAMP activators, forskolin (F) and rolipram (R), and an agonist for the DP receptor, BW245C. Primary microglia were isolated from the cortices from P2 mice and plated in 96 well format on Poly-D-Lysine-coated plates. Cells were either left untreated or stimulated for 2 hours with 100 ng/ml LPS. Compounds were then added for an additional 4 hours before the Bright Glo assay was run. The compounds used were 5 μ M forskolin, 10 μ M rolipram or the combination of the two, or the Gs agonist for the DP1 receptor, BW245C at 10 μ M. Data is shown as counts per second (cps).

Figure 23 shows the effects of intrathecally injected forskolin (F) and rolipram (R) on the induction of luciferase expression in the brain and spinal cord of CreLuc mice (line 187). N=3-4 mice per group, 3 month old males. Group A: DMSO control, Group B: 1ug forskolin/ 10ug rolipram, Group C: 10ug forskolin/ 10ug rolipram, Group D: 40ug forskolin/ 10ug rolipram. The animals were dosed via intrathecal injection, lumbar region, and volume of 5 μ l per mouse. They were imaged at 4 hours post dosing. The data for both spinal cord and brain is shown as the average peak radiance, photons per second per cm².

Figure 24 shows the effects of the EP2 agonist, EX00000173A on luciferase expression in the brain and spinal cord of CreLuc mice. Mice (line 187) were injected i.p. with either vehicle (5% DMSO, 0.05% tween 80, PBS) or 10 mg/kg EX00000173A. Animals were bioimaged at 4 hours post doing. Data shown as photons per second per cm².

Figure 25 shows the effects of the EP2 agonist, EX00000173A on luciferase expression in CreLuc mice. Mice were dosed with either vehicle control or varying doses of the EP2 agonist EX00000173A. Mice were dosed by intrathecal injection (5 μ l per mouse) and were bioimaged 4 hours later on the IVIS bioimager. Data is shown as the mean of the five mice, average peak radiance, photons per second per cm².

Figure 26 shows induction of luciferase in different tissues by the adrenoceptor beta3 (Adrb3) agonist, CL316,243 (1 mg/kg, ip) in CRE-Luc mice. The luciferase assay was performed in tissue homogenates.

Figure 27A shows induction of luciferase by the Adrb3 agonist CL316,243 (1 mg/kg, ip) in lines 11 (n=2) and 115 (n=3) of CRE-Luc mice. BLIs were taken before and 4-5 hours after treatment. Luciferase activities in tissue homogenates is shown below the pictures.

Figure 27B shows induction of luciferase by the Adrb3 agonist CL316,243 (1 mg/kg, ip) in lines 31 (n=2) and 175 (n=3) of CRE-Luc mice. BLIs were taken before and 4-5 hours after treatment. Luciferase activities in tissue homogenates is shown below the pictures.

Figure 28 shows induction of luciferase reporter by the glucagon-like peptide 1 receptor (GLP-1R) agonist, AVE0010, in three independent lines of CRE-Luc mice. Baseline images were acquired on day 1. On day 2, mice were treated with AVE0010 (0.1 mg/kg, sc) and imaged after 4 hours. Fold induction over baseline at indicated at the bottom.

Figure 29 shows induction of luciferase reporter by the glucagon-like peptide 1 receptor (GLP-1R) agonist, AVE0010, in three independent lines of CRE-Luc mice. Mice were treated with AVE0010 (0.1 mg/kg, sc) for 4 hours. Luciferase activities in 8 different tissues were measured.

Figure 30 shows the effects of the beta-cell toxin streptozotocin (STZ) on the induction of CRE-Luc by AVE0010. Male CRE-luc mice (line 11) were imaged before ("uninduced"; top panel) and after AVE0010 was given at 0.1 mg/kg, sc ("induction by AVE0010", middle panel). All mice were responsive to AVE0010 (middle panel). Then, the animals were treated with vehicle (control) or STZ (200 mpk, ip). Four days later, they were imaged again after AVE0010 treatment (bottom panel).

Figure 31 shows that the induction of CRE-Luc by AVE0010 is likely beta-cell-specific. Animals were treated as described in Figure 30. Blood glucose levels were measured by tail vein nicking on unfasted mice. Glucose levels were read on a Bayer glucometer. Glucose levels are shown as mg glucose/ml. Fold induction is the

luciferase bioimaging levels of AVE10 dosing versus the baseline signals. Blood glucose levels (BG) were increased by STZ (upper left panel). Non-fasting BG levels were reduced by AVE0010 (0.1 mg/kg, sc). BLI data shown in Figure 30 were quantified.

5

Figure 32 shows CreLuc bone marrow engraftments into NOD scid gamma (NSG) mice. Bone marrow cells were harvested from lines 44 heterozygotes and line 64 homozygotes. The cells were then engrafted via tail vein injections of cells into irradiated NSG mice at 1 million or 5 million cells per mouse. For line 44: mouse 1 and 2 received 5 million cells, mouse 3 and 4 received 1 million cells; for line 64: mouse 1 received 5 million cells, mouse 2, 3, and 4 received 1 million cells. (4 NSG mice per CreLuc line). The animals were bioimaged at 4 weeks (data not shown) and then again at 8 weeks (data shown). Prior to imaging, the line 64 mice were induced for 5 hours with 5mg/kg forskolin and 10mg/kg rolipram.

15

Figure 33 shows the effects of forskolin, rolipram and isoproterenol on luciferase expression in mouse embryonic fibroblasts. Mouse embryonic fibroblasts were cultured from E12 embryos from six independent CreLuc lines and plated at 20,000 cells per well. Compounds tested include 10 μ M forskolin (F), 5 μ M rolipram (R) and 10 μ M isoproterenol (iso). Data is shown as counts per second (cps).

20

Figure 34 show the effects of zymosan treatment on luciferase levels in CreLuc mice (line 187). Animals in the treated group were injected s.c in both rear paws with zymosan (zymo) to induce a pain response. The animals were then bioimaged daily for 4 days (denoted as d1, d2, d3 and d4).

25

Figure 35 shows the effects of forskolin and rolipram and isoproterenol on luciferase levels in cardiomyocytes. Cardiomyocytes were isolated from P3 pups (line 229). The cells were cultured in a 96 well plate. Compounds tested include 10 μ M forskolin (F), 5 μ M rolipram (R) and 10 μ M isoproterenol (iso). Data is shown as counts per second (cps).

30

DETAILED DESCRIPTION OF THE INVENTION

Unless otherwise defined, all technical and scientific terms used herein have the same meaning as commonly understood by one of ordinary skill in the art to which
5 this invention belongs.

Each publication, patent application, patent, and other reference cited herein is incorporated by reference in its entirety to the extent that it is not inconsistent with this present disclosure.

10 It is noted here that, as used in this specification and the appended claims, the singular forms "a," "an," and "the" include plural reference unless the context clearly dictates otherwise.

15 Furthermore, in accordance with the present invention there may be employed conventional molecular biology, microbiology, and recombinant DNA techniques within the skill of the art. Such techniques are explained fully in the literature. See, e.g., Sambrook, Fritsch & Maniatis, *Molecular Cloning: A Laboratory Manual*, Second Edition (1989) Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York
20 (herein "Sambrook et al., 1989"); *DNA Cloning: A Practical Approach*, Volumes I and II (D.N. Glover ed. 1985); *Oligonucleotide Synthesis* (M.J. Gait ed. 1984); *Nucleic Acid Hybridization* [B.D. Hames & S.J. Higgins eds. (1985)]; *Transcription And Translation* [B.D. Hames & S.J. Higgins, eds. (1984)]; *Animal Cell Culture* [R.I. Freshney, ed. (1986)]; *Immobilized Cells And Enzymes* [IRL Press, (1986)]; B.
25 Perbal, *A Practical Guide To Molecular Cloning* (1984); F.M. Ausubel et al. (eds.), *Current Protocols in Molecular Biology*, John Wiley & Sons, Inc. (1994).

A "test agent" is interpreted broadly to include any material such as a compound or chemical compounds, e.g., organic chemical entities, inorganic chemical entities,
30 biologic compounds or biological materials, e.g., antibodies and antigen recognizing fragments and constructs, nucleic acids, e.g., RNAi, etc. A test agent encompasses a single agent or multiple agents applied together.

As used herein, a “transgenic animal” is a non-human animal, a non-limiting example being a mammal, in that one or more of the cells of the animal includes a genetic modification as defined herein. Further non-limiting examples includes rodents such as a rat or mouse. Other examples of transgenic animals include non-human
5 primates, sheep, dogs, cows, goats, chickens, amphibians etc. The choice of transgenic animal is only limited by the ability of light generated from the reporter to cross tissues and reach the surface where detection can occur.

As used herein, a “genetic modification” is one or more alterations in the non-human
10 animal’s gene sequences. A non-limiting example is insertion of a transgene into the genome of the transgenic animal.

As used herein, the term “transgene” refers to exogenous DNA containing a promoter, reporter gene, poly adenylation signal and other elements to enhance
15 expression (insulators, introns). This exogenous DNA integrates into the genome of a 1-cell embryo from which a transgenic animal develops and the transgene remains in the genome of the mature animal. The integrated transgene DNA can occur at single or multiple places in the genome of the egg or mouse and also single to multiple (several hundred) tandem copies of the transgene can integrate at each
20 genomic location.

The term “general cAMP modulator” refers to chemical compounds, e.g., organic chemical entities, inorganic chemical entities, biologic compounds or biological materials, e.g., antibodies and antigen recognizing fragments and constructs, nucleic
25 acids, e.g., RNAi, etc capable of increasing or maintaining cAMP levels. Non-limiting examples include forskolin and rolipram. A general cAMP modulator encompasses a single cAMP modulator or multiple cAMP modulators applied together.

The term “cell stimulator” refers to chemical compounds, e.g., organic chemical
30 entities, inorganic chemical entities, biologic compounds or biological materials, e.g., antibodies and antigen recognizing fragments and constructs, nucleic acids, e.g., RNAi, etc capable of activating the cell or causing the cell to be in a more activated state. Non-limiting examples include lipopolysaccharide and anti CD3.

An embodiment of the invention uses a control. A control is a term of art well understood by skilled artisans. An appropriate control may be dependent on the assay parameters utilized or the experimental question under investigation. Typically, a control is a vehicle control in which the control is the same buffer or solvent used to dissolve test agent or compounds. A non-limiting example is if phosphate-buffered saline is used to dissolve compound then the vehicle control would be phosphate buffered saline. Similarly, if DMSO is used to dissolve test agents, then the control is DMSO. Often, more than one control must be used per experiment or assay because more than one diluent is used for the compounds tested.

As used herein, "luciferase" refers not only to luciferase enzyme activity but also to actual amounts of luciferase protein.

In accordance with the present invention there may be employed conventional techniques known to those skilled in the art to generate transgenic non-human animals. For instance, Pinkert, C. A. (ed.) 1994. Transgenic animal technology: A laboratory handbook. Academic Press, Inc., San Diego, Calif.; Monastersky G. M. and Robl, J. M. (ed.) (1995) Strategies in transgenic animal science. ASM Press. Washington D.C.. and Nagy A, Gertsenstein, M, Vintersten, K, Behringer R 2003. Manipulating the Mouse Embryo; A laboratory Manual 3rd edition. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York.

Transgene Elements

An embodiment of the invention relates to a transgene. The transgene may comprise insulator elements, response elements, promoter elements, reporter elements, and functional elements.

Response elements are palindromic DNA sequences that respond to cellular signals such as hormones, enzymes, or other key signaling proteins within a cell. Non-limiting examples of response elements include CRE (cAMP response element), estrogen response elements and others listed in Table 2. Response elements may

be incorporated into the transgene as a single DNA sequence or in tandem repeats. For instance, CRE response elements repeated four times or six times have been used in transgene construction and validated in vitro (Deutsch P.J., et al., J. Biol. Chem., 263;18466-18472, 1988; Oetjen E JBC 269;27036-27044, 1994). Cre response elements have also been compared in vivo and the increase in multimers has correlated with an increase transcriptional response to cAMP pathway activators (Montoliu, L. et al., Proc. Natl. Acad. Sci. USA 92;4244, 1995; Boer et al, PloS One, May 9; 2(5):e431, 2007). An embodiment of the invention may utilize any known response element, either as a single sequence or in multiple tandem repeats, for instance, tandem six repeat of CRE (6X CRE).

Table 2: Cis acting response elements

Cis acting enhancer element	Transcription Factor(s)	Signal transduction pathway
Activator protein 1 (ASP1)	c-jun/c-fos	JNK
cAMP response element (CRE)	ATF2/CREB	JNK/p38, PKA
Estrogen response element (ERE)	Estrogen receptor	Estrogen receptor
Glucocorticoid response element (GRE)	GR	Glucocorticoid/ HSP90
Heat shock response element (HSE)	HSF	Heat shock response
Serum response element (SRE)	Lek-1/SRF	MAPK/JNK
Thyroid response element (TRE)	Thyroid receptor	Thyroid hormone receptor

DNA promoter elements are regions of DNA that facilitate the transcription of a particular gene. Promoters are typically located near the genes they regulate, on the same strand and upstream (towards the 5' region of the sense strand). Promoters

contain specific DNA sequences and response elements which provide a binding site for RNA polymerase and for proteins called transcription factors that recruit RNA polymerase. DNA promoters are highly variable in their size and internal substructures that contribute to the regulation of a particular gene's expression in time and space. In a non-limiting example of promoter elements, the herpes simplex virus thymidine kinase minimal promoter (HSV TK min) is designed to allow expression of the reporter gene in every cell type. Only its core expression elements are retained to impart ubiquitous expression either in vitro or in vivo (Park, J., et al., DNA Cell Bio., 12:1147-1149, 1994). An embodiment of the invention may utilize any known promoter element. As a non-limiting example, any known promoter element may be used such that the promoter element when combined with the CRE cis-activating response element allows gene expression of the reporter to be responsive to cAMP pathway modulation.

Since the response elements such as CRE and the promoter element HSV TK min are small in size, a transgene promoter containing these elements regulating the transcription of a reporter element will be very sensitive to position effects and lead to poor expression responses to ligands especially at low ligand concentrations in vivo. Thus, an embodiment of the invention is to include to additional elements to be added to the transgene to achieve high levels of expression and wide distribution of a functioning transgene throughout all the cellular compartments. These elements include insulator elements and functional elements (Sun F.L and Elgin S.C, Cell 99:459-462, 1999).

An embodiment of the invention utilizes functional elements or functional enhancer elements within the transgene. A nonlimiting example of a functional element is the human growth hormone (hGH) gene. The hGH sequence used in the CreLuc transgene contains several design elements that contribute and interact with the insulator element to achieve the bioimaging model design goals. The hGH sequence contains all the hGH genomic structure and thus supplies several important elements, but is not transcribed or translated into a protein. The critical influence of the hGH sequence to improve the production of a functioning transgene has been demonstrated since 1990 for several transgenic models Erickson LA , Nature 346:

74-76, 1990. While there is not a comparative analysis of its importance, the hGH structure does contain several important and critical DNA elements:

- a. Intron splicing: Initially transcribed mRNA contains both intronic and exon sequences which then are exported from the nucleus and further processed to remove intronic sequence resulting in a mature mRNA containing only exon sequences. This transport and trimming process connect the maturing mRNA strand to additional translational machinery for the final production of a protein. Including introns in transgene cDNA structure has been shown to improve the expression level of the transgene cDNA (Palmiter, R.D., et al., Proc. Natl. Acad. Sci. USA 88:478-482, 1988)
- b. Genomic structure with intact 3'UTR: In an embodiment of the invention, the hGH sequence in the transgene contains an intact 3'UTR which imparts a high degree of mRNA stability and thus higher levels of transgene expression
- c. Genomic structure with poly A (PA+) structures The hGH sequence contains its native PA+ structure imbedded in the nature 3' UTR. Typically PA+ signals are from viral sequences (SV40, RSV etc.) and are minimal structures added to the end of unrelated 3'UTRs. In an embodiment of the invention, the entire 3'UTR structure with the nature PA+ signal is preserved in the even wider genomic context of the full hGH gene.

Insulator elements are sequences of DNA that generate position independent expression (Giraldo et al, Transgenic Research 12: 751–755, 2003). Insulator sequences were described in the 1980s for the globulin locus (Sun F.L. and Elgin, S.C., Cell 99:459-462, 1999) and were reported to increase the chance of obtaining correct and responsive transgenic expression in selected tissues to support the model design goals for bioimaging. (Pinkert, C. A. (ed.) 1994. Transgenic animal technology: A laboratory handbook. Academic Press, Inc., San Diego, Calif.; Monastersky G. M. and Robl, J. M. (ed.) (1995) Strategies in transgenic animal science. ASM Press. Washington D.C) Insulators are DNA elements that create open chromatin domains permissive to gene expression and constitute a barrier against the influence of distal silencer/enhancer sequences and against acetylation and methylation events. They should significantly increase the number of independent

transgenic founder lines that have the reporter gene expressed at detectable levels for bioimaging. Insulator elements have been shown to increase the number of luciferase expressing clones in a transient transfection assay from 40 to 70% thus enhancing the inducibility of luciferase expression in an ERE-luci transgene.

(Ottobrini L., Mol Cell Endo 246, 69-75). However a full review of the application of insulators to transgene reporter expression leads to the conclusion that in practical terms, it remains difficult to utilize insulators. Their mechanism of action is only partly known and their effect is not fully predictable. Non-limiting examples of insulator elements are included in Table 3.

Table 3: Examples of known insulator elements to improve transgene DNA expression

Insulator	Gene or origin	References
DNase I-hypersensitive site (HS4)	Human beta-globin	Chung et al 1993
Matrix attachment regions (MAR)	Chicken lysozyme	Stief et al 1989
Inverted terminal repeats (ITR)	Adeno-associated virus	Fu et al, 1998

In an embodiment of the Cre-Luc transgene, the inclusion of insulator elements significantly increased the frequency of generating lines with a functional reporter as detected by bioimaging. For an analysis of the contribution of insulator elements to luciferase expression in our CreLuc transgenic mouse lines, see section VI of "Examples" below.

An embodiment of the invention relates to a transgene comprising a reporter element or gene. A reporter gene includes any gene that expresses a detectable gene product, which may be an RNA or a protein. Many reporter genes are known in the art, including, but not limited to beta-galactosidase and alkaline phosphatase. In another embodiment, the transgene comprises a bioluminescent reporter gene. Many bioluminescent reporter genes are known in the art, including, but not limited to

luciferase. There are many sources of luciferase, nonlimiting examples include firefly luciferase and bacterial luciferase. An embodiment of the invention may utilize any known bioluminescent reporter, for instance, luciferase. Other non-limiting examples of reporter elements are shown in Table 5.

5

TABLE 4. Genetic reporter systems

Reporter Gene	In vitro assay	In vivo assay
Chloramphenicol acetyltransferase (CAT)	Chromatography, differential extraction, fluorescence or immunoassay	RNA or enzyme assay of tissue extracts
Luciferase - firefly	Bioluminescence with luminometer or scintillation counter	Assay in live cells with luciferin esters as substrate
Luciferase – renilla	Bioluminescence with luminometer or scintillation counter	Assay in live cells with luciferin esters as substrate
Beta- galactosidase	Colorimetric, fluorescence, chemiluminescence	Histochemical staining with X-gal substrate: bioluminescence assay in live cells with fluorescein di-beta –D- galactopyransodie (FDG)
Secreted alkaline phosphatase (SEAP)	Colorimetric bioluminescence, or chemiluminescence	RNA or enzyme assay from tissue extracts
Human growth hormone (hGH)	Radioimmunoassay	RNA or protein assay from tissue extracts
Green fluorescent protein (GFP)	Fluorescence (UV light box or fluorescence imaging device)	Fluorescence (UV light box, fluorescence microscopy or FACS

Other embodiments of the invention can incorporate modified versions of the luciferase enzyme, luciferase enzyme from different species or any other protein that can produce light able to cross animal tissues or any enzyme that can emit light able to cross animal tissues when provided with a suitable substrate. The reporter protein of the present invention is only limited by the fact that signal attenuation depends on the wavelength of the light being emitted and the tissue properties surrounding the emitting cells. Generally, blue-green light (400 590 nm) is strongly attenuated while red to near-infrared light (590 800 nm) suffers much less attenuation. Most types of luciferase have peak emission at blue to yellow-green wavelengths, the emission spectrum is broad such that there is significant emission at red wavelengths (>600 nm) that penetrate quite deeply into tissue. For small rodents such as mice, this allows detection of signals throughout the entire animal.

The limits of light detection in vivo depend on the type of bioluminescent reporter, the surrounding physiology of the animal and on the source depth. Typically, bioluminescent cells in animals can be observed from 1 3 cm deep with sensitive CCD cameras, depending on the number and location of the cells. Scattering of photons as they propagate through tissue limits the spatial resolution of images detected on the animal surface. In general, spot size or resolution on the surface is approximately equal to the depth of the source below the surface. Using physics based diffusion models, improvements in spatial resolution approaching the millimeter level can be achieved. Using cooled scientific grade CCD arrays, the limit in signal detection is determined by the read noise associated with reading CCD pixels after an image is taken, which is on the order of a few photons per pixel (Honigman et al., Mol. Ther. 4:239-249, 2001). There may be additional background light coming from the animal due to phosphorescence of the fur, skin, or perhaps contaminants on the animal. Typically, background light is at a low level and only has a deleterious effect on images of deep low-level bioluminescent sources. However, background light can be eliminated by using use of an appropriate optical filter.

An embodiment of the invention utilizes CCD cameras such as the IVIS (Xenogen Corporation, 860 Atlantic Avenue, Alameda, Calif. 94501, USA). The IVIS.RTM. Imaging System includes a sensitive CCD camera, a dark imaging chamber to

minimize incident light, and specialized software to quantify and analyze the results. IVIS is a registered trademark of Xenogen Corporation. However, any such bioluminescence imaging system can be applied to the instant invention.

5 Real-time in vivo imaging allows the quantification of the bioluminescent reporter gene non-invasively, i.e., the animal does not need to be euthanized, and longitudinally, i.e., the measurements can continuous or repeated over a prolonged time course. Real-time in vivo imaging requires fewer test animals (e.g., because the same animal can be used over a specified time period) and less time (e.g., because
10 fewer animals need to be handled) than conventional protocols allowing more test compounds to be tested for efficacy. real-time in vivo imaging provides for a higher data content and higher data quality for many reasons. For instance, temporal and spatial data can be collected from the same animal and data can be collected without need for time-consuming histological assessment. Higher data quality decreases
15 statistical error and improves the quality of test compound assessment and decision making.

An embodiment of the invention is use in high throughput screening (HTS) methods. HTS is the automated, simultaneous testing of thousands of distinct chemical
20 compounds in assays designed to model biological mechanisms or aspects of disease pathologies. More than one compound, e.g., a plurality of compounds, can be tested simultaneously, e.g., in one batch. In one embodiment, the term HTS screening method refers to assays which test the ability of one compound or a plurality of compounds to influence the readout of choice.

25 Liquid handling systems, analytical equipment such as fluorescence readers or scintillation counters and robotics for cell culture and sample manipulation are well known in the art. Mechanical systems such as robotic arms or "cherry-picking" devices are available to the skilled artisan. Commercial plate readers are available to
30 analyze conventional 96- well or 384-well plates. Single sample, multiple sample or plate sample readers are available that analyze predetermined wells and generate raw data reports. The raw data can be transformed and presented in a variety of ways.

An embodiment of the invention comprises an array of receptacles that can receive cells, tissue slices and other materials such as culture media. An array of receptacles can be any number of receptacles from at least one or more than one receptacle
5 suitable for holding cells or tissue slices within the scope of the invention. Examples include but are not limited to flasks, culture dishes, tubes such as 1.5 ml tubes, 12 well plates, 96 well plates, 384 well plates and miniaturized microtiter plates with perhaps 4000 receptacles (U.S. Patent Application 20050255580). The array of receptacles may be amendable to the addition of a protective covering thus
10 preventing against entry of contaminants or evaporation of contents.

A further characteristic of the receptacles is that the receptacle may allow for analysis, non-limiting examples include, spectrophotometric analysis, scintillation counting and fluorescence measurements. However, this is not a limitation to
15 receptacles that can be used within the scope of the invention given that samples can be transferred to a suitable container amendable for further analysis. A non limiting example is to modify the method such that the method further comprises providing a second array of receptacles wherein the step of lysing the cells further comprises separating supernatant from cell debris and the next step further comprises adding a
20 detectable compound capable of intercalating into DNA fragments to at least one receptacle of said second array of receptacles containing a sample of said separated supernatant.

Additional aspects and details of the invention will be apparent from the following
25 examples, which are intended to be illustrative and are not meant to limit the scope of the invention in any way.

All experimental work involving animals was performed in accordance with federal guidelines and protocols were prior reviewed and approved by the sanofi-aventis site
30 Institutional Animal Care and Use Committee (IACUC).

EXAMPLES

I. Vector Backbone for the MultiSite Gateway Pro Plus Cloning System (Fig. 3)

5 A destination vector, pDest2XMARS was designed and cloned specifically for use with the MultiSite Gateway Pro Plus Cloning System (Invitrogen Carlsbad, CA, cat #12537-100). All elements inserted into the vector were either PCR (polymerase chain reaction) cloned or synthetically generated. For all PCR cloning steps, the fragments were amplified from the specified vector through the use of PCR SuperMlx
10 Hi Fidelity (Invitrogen, Carlsbad, CA, cat # 10790-020) and fragment specific primers. The PCR products were then subcloned into pCR2.1 vectors by TOPO cloning. TOPO Cloning is a molecular biology technique in which DNA fragments amplified by either Taq or Pfu polymerases are cloned into specific vectors without the requirement for DNA ligases.

15 First, insulator elements were cloned by polymerases chain reaction (PCR) from the vector pCpG-LacZ (InvivoGen San Diego, CA, cat #pcpg-lacz), sequence verified and then subcloned into restriction sites of pNEB193 (New England Biolabs Ipswich, MA, cat#N3051S). The insulator elements were used to reduce variability in
20 expression of the transgene due to integration site dependant position effects.

The PCR primers used to clone the fragments are the following:

a. Human IFN- β MAR primer sequences (primers have EcoRV sites for subcloning)

25 SEQ ID NO: 1, IFN forward primer:

5'-GGGGGATATCAGTCAATATGTTACCCCA-3'

SEQ ID NO: 2, IFN reverse primer:

30 5'-GGGGGATATCCTACTGTTTTAATTAAGC-3'

b. Human β -globin MAR primer sequencesSEQ ID NO: 3, β -globin forward primer:

5'-AAGGATCCTTAATTAATAATTATCTCTAAGGC-3'

5

SEQ ID NO: 4, β -globin reverse primer:

5'-GGATCCCTGCAGGAATTCCTTTTAAT-3'

10 The β -globin PCR fragment was topo cloned using pCR2.1 (Invitrogen, Carlsbad, CA, cat#K2030-01). After sequence confirmation, the fragment was cut out with BamHI and then subcloned into the BamHI site of pNEB193 (β gloMAR-pNEB193). The IFN- β PCR fragment was also topo cloned with pCR2.1, sequenced, and then cut out with EcoRV and subcloned into the HincII site of β gloMAR-pNEB193 (2XMARS-
15 pNEB193).

For the final Gateway destination vector, a linker containing XbaI ends and an internal EcoRV site was subcloned into XbaI site of 2XMARS-pNEB193. The blunt ended Gateway conversion cassette, RfA was then inserted into the EcoRV site of 2XMARS-
20 pNEB193 (2XMARSpDest). The final vector, 2XMARSpDest was the final destination vector used to generate the transgenes.

SEQ ID NO: 5

25 TCGCGCGTTTCGGTGATGACGGTGAAAACCTCTGACACATGCAGCTCCCGGAG
ACGGTCACAGCTTGTCTGTAAGCGGATGCCGGGAGCAGACAAGCCCGTCAGGG
CGCGTCAGCGGGTGTGTTGGCGGGTGTGCGGGGCTGGCTTAAGTATGCGGCATCAG
AGCAGATTGTACTGAGAGTGCACCATATGCGGTGTGAAATACCGCACAGATGCG
TAAGGAGAAAATACCGCATCAGGCGCCATTCGCCATTCAGGCTGCGCAACTGTT
30 GGGAAGGGCGATCGGTGCGGGCCTCTTCGCTATTACGCCAGCTGGCGAAAGG
GGGATGTGCTGCAAGGCGATTAAGTTGGGTAACGCCAGGGTTTTCCAGTCAC
GACGTTGTAAAACGACGGCCAGTGAATTCGAGCTCGGTACCCGGGGGCGCGCC
GGGATCCTTAATTAATAATTATCTCTAAGGCATGTGAACTGGCTGTCTTGGTTTTTC

ATCTGTACTTCATCTGCTACCTCTGTGACCTGAAACATATTTATAATTCCATTAAG
CTGTGCATATGATAGATTTATCATATGTATTTTCCTTAAAGGATTTTTGTAAGAAC
TAATTGAATTGATACCTGTAAAGTCTTTATCACACTACCCAATAAATAATAAATCT
CTTTGTTTCAGCTCTCTGTTTCTATAAATATGTACCAGTTTTATTGTTTTTAGTGGA
5 GTGATTTTTATTCTCTTTCTATATATATACACACACATGTGTGCATTCATAAATATAT
ACAATTTTTATGAATAAAAAATTATTAGCAATCAATATTGAAAACCACTGATTTTTG
TTTATGTGAGCAAACAGCAGATTAAAAGGAATTCCTGCAGGATCCTTAATTAAGT
TCTAGATCACAAGTTTGTACAAAAAGCTGAACGAGAAACGTAAAATGATATAAA
TATCAATATATTAATTAGATTTTGCATAAAAAACAGACTACATAATACTGTAAAAC
10 ACAACATATCCAGTCACTATGGCGGCCGCATTAGGCACCCCAGGCTTTACACTT
TATGCTTCCGGCTCGTATAATGTGTGGATTTTGAGTTAGGATCCGTCGAGATTTT
CAGGAGCTAAGGAAGCTAAAATGGAGAAAAAATCACTGGATATACCACCGTTG
ATATATCCCAATGGCATCGTAAAGAACATTTTGAGGCATTTTCAGTCAGTTGCTCA
ATGTACCTATAACCAGACCGTTCAGCTGGATATTACGGCCTTTTTAAAGACCGTA
15 AAGAAAAATAAGCACAAAGTTTTATCCGGCCTTTATTCACATTCTTGCCCGCCTGA
TGAATGCTCATCCGGAATTCCGTATGGCAATGAAAGACGGTGAGCTGGTGATAT
GGGATAGTGTTACCCCTTGTTACACCGTTTTCCATGAGCAAACCTGAAACGTTTTTC
ATCGCTCTGGAGTGAATACCACGACGATTTCCGGCAGTTTCTACACATATATTG
CAAGATGTGGCGTGTTACGGTGAAAACCTGGCCTATTTCCCTAAAGGGTTTATT
20 GAGAATATGTTTTTCGTCTCAGCCAATCCCTGGGTGAGTTTCACCAGTTTTGATT
TAAACGTGGCCAATATGGACAACCTTCTTCGCCCCCGTTTTACCATGGGCAAATA
TTATACGCAAGGCGACAAGGTGCTGATGCCGCTGGCGATTCAGGTTTCATCATGC
CGTTTGTGATGGCTTCCATGTCGGCAGAATGCTTAATGAATTACAACAGTACTGC
GATGAGTGGCAGGGCGGGGCGTAAACGCGTGGATCCGGCTTACTAAAAGCCAG
25 ATAACAGTATGCGTATTTGCGCGCTGATTTTTGCGGTATAAGAATATATACTGAT
ATGTATACCCGAAGTATGTCAAAAAGAGGTATGCTATGAAGCAGCGTATTACAGT
GACAGTTGACAGCGACAGCTATCAGTTGCTCAAGGCATATATGATGTCAATATCT
CCGGTCTGGTAAGCACAAACCATGCAGAATGAAGCCCGTCGTCTGCGTGCCGAA
CGCTGGAAAGCGGAAAATCAGGAAGGGATGGCTGAGGTCGCCCCGGTTTATTGA
30 AATGAACGGCTCTTTTGTGACGAGAACAGGGGCTGGTGAAATGCAGTTTAAGG
TTTACACCTATAAAAGAGAGAGCCGTTATCGTCTGTTTGTGGATGTACAGAGTGA
TATTATTGACACGCCCGGGCGACGGATGGTGATCCCCCTGGCCAGTGCACGTC
TGCTGTCAGATAAAGTCTCCCGTGAACCTTACCCGGTGGTGATATCGGGGATG

AAAGCTGGCGCATGATGACCACCGATATGGCCAGTGTGCCGGTCTCCGTTATC
GGGGAAGAAGTGGCTGATCTCAGCCACCGCGAAAATGACATCAAAAACGCCATT
AACCTGATGTTCTGGGGAATATAAATGTCAGGCTCCCTTATACACAGCCAGTCTG
CAGGTCGACCATAGTGACTGGATATGTTGTGTTTTACAGTATTATGTAGTCTGTT
5 TTTTATGCAAAATCTAATTTAATATATTGATATTTATATCATTTTTACGTTTCTCGTTC
AGCTTTCTTGTACAAAGTGGTGTCTAGACTAGAGTCATCAGTCAATATGTTTAC
CCCAAAAAGCTGTTTGTAACTTGTCAACCTCATTCTAAAATGTATATAGAAGCC
CAAAAGACAATAACAAAAATATTCTTGTAGAACAAAATGGGAAAGAATGTTCCAC
TAAATATCAAGATTTAGAGCAAAGCATGAGATGTGTGGGGATAGACAGTGAGGC
10 TGATAAAATAGAGTAGAGCTCAGAAACAGACCCATTGATATATGTAAGTGACCTA
TGAAAAAATATGGCATTTTACAATGGGAAAATGATGATCTTTTTCTTTTTTAGAA
AAACAGGGAAATATATTTATATGTAAAAAATAAAAGGGAACCCATATGTCATACCA
TACACACAAAAAAATTCCAGTGAATTATAAGTCTAAATGGAGAAGGCAAAACCTTT
AAATCTTTTAGAAAATAATATAGAAGCATGCCATCAAGACTTCAGTGTAGAGAAA
15 AATTTCTTATGACTCAAAGTCCTAACCACAAAGAAAAGATTGTTAATTAGATTGCA
TGAATATTAAGACTTATTTTTAAAATTA AAAAACCATTAAGAAAAGTCAGGCCATA
GAATGACAGAAAATATTTGCAACACCCCAGTAAAGAGAATTGTAATATGCAGATT
ATAAAAAGAAGTCTTACAAATCAGTAAAAAATAAAACTAGACAAAAATTTGAACAG
ATGAAAGAGAAACTCTAAATAATCATTACACATGAGAAACTCAATCTCAGAAATCA
20 GAGAACTATCATTGCATATACACTAAATTAGAGAAATATTA AAAAGGCTAAGTAACA
TCTGTGGCTTAATTA AACACAGTAGGATGACTGTTTAAACCTGCAGGCATGCAAGC
TTGGCGTAATCATGGTCATAGCTGTTTCCTGTGTGAAATTGTTATCCGCTCACAA
TTCCACACAACATACGAGCCGGAAGCATAAAGTGTAAGCCTGGGGTGCCCTAAT
GAGTGAGCTAACTCACATTAATTGCGTTGCGCTCACTGCCCCGCTTTCCAGTCGG
25 GAAACCTGTCGTGCCAGCTGCATTAATGAATCGGCCAACGCGCGGGGAGAGGC
GGTTTGCGTATTGGGCGCTCTTCCGCTTCCTCGCTCACTGACTCGCTGCGCTCG
GTCGTTCCGGCTGCGGCGAGCGGTATCAGCTCACTCAAAGGCGGTAATACGGTT
ATCCACAGAATCAGGGGATAACGCAGGAAAGAACATGTGAGCAAAAGGCCAGC
AAAAGGCCAGGAACCGTAAAAAGGCCGCGTTGCTGGCGTTTTTCCATAGGCTCC
30 GCCCCCTGACGAGCATCACAAAAATCGACGCTCAAGTCAGAGGTGGCGAAAC
CCGACAGGACTATAAAGATACCAGGCGTTTCCCCCTGGAAGCTCCCTCGTGCG
CTCTCCTGTTCCGACCCTGCCGCTTACCGGATACCTGTCCGCCTTTCTCCCTTC
GGGAAGCGTGGCGCTTTCTCATAGCTCACGCTGTAGGTATCTCAGTTCGGTGTA

GGTCGTTGCTCCAAGCTGGGCTGTGTGCACGAACCCCCCGTTCAGCCCGACC
GCTGCGCCTTATCCGGTAACTATCGTCTTGAGTCCAACCCGGTAAGACACGACT
TATCGCCACTGGCAGCAGCCACTGGTAACAGGATTAGCAGAGCGAGGTATGTA
GGCGGTGCTACAGAGTTCTTGAAGTGGTGGCCTAACTACGGCTACACTAGAAGA
5 ACAGTATTTGGTATCTGCGCTCTGCTGAAGCCAGTTACCTTCGGAAAAAGAGTT
GGTAGCTCTTGATCCGGCAAACAAACCACCGCTGGTAGCGGTGGTTTTTTTGT
TGCAAGCAGCAGATTACGCGCAGAAAAAAGGATCTCAAGAAGATCCTTTGATC
TTTTCTACGGGGTCTGACGCTCAGTGGAACGAAAACCTCACGTTAAGGGATTTTG
GTCATGAGATTATCAAAAAGGATCTTCACCTAGATCCTTTTAAATTAATAATGAAG
10 TTTTAAATCAATCTAAAGTATATATGAGTAACTTGGTCTGACAGTTACCAATGCT
TAATCAGTGAGGCACCTATCTCAGCGATCTGTCTATTTTCGTTTCATCCATAGTTGC
CTGACTCCCCGTCGTGTAGATAACTACGATACGGGAGGGCTTACCATCTGGCCC
CAGTGCTGCAATGATACCGCGAGACCCACGCTCACCGGCTCCAGATTTATCAGC
AATAAACCCAGCCAGCCGGAAGGGCCGAGCGCAGAAGTGGTCCTGCAACTTTAT
15 CCGCCTCCATCCAGTCTATTAATTGTTGCCGGGAAGCTAGAGTAAGTAGTTTCGC
CAGTTAATAGTTTGCGCAACGTTGTTGCCATTGCTACAGGCATCGTGGTGTCAC
GCTCGTCGTTTGGTATGGCTTCATTACGCTCCGGTTCCCAACGATCAAGGCGAG
TTACATGATCCCCCATGTTGTGCAAAAAAGCGGTTAGCTCCTTCGGTCCTCCGAT
CGTTGTCAGAAGTAAGTTGGCCGCGAGTGTTATCACTCATGGTTATGGCAGCACT
20 GCATAATTCTCTTACTGTCATGCCATCCGTAAGATGCTTTTCTGTGACTGGTGAG
TACTCAACCAAGTCATTCTGAGAATAGTGTATGCGGCGACCGAGTTGCTCTTGC
CCGGCGTCAATACGGGATAATACCGCGCCACATAGCAGAACTTTAAAAGTGCTC
ATCATTGGAAAACGTTCTTCGGGGCGAAAACCTCTCAAGGATCTTACCGCTGTTG
AGATCCAGTTCGATGTAACCCACTCGTGACCCAACTGATCTTCAGCATCTTTTA
25 CTTTCACCAGCGTTTCTGGGTGAGCAAAAACAGGAAGGCAAAATGCCGCAAAAA
AGGGAATAAGGGCGACACGGAAATGTTGAATACTCATACTCTTCCTTTTTCAATA
TTATTGAAGCATTTATCAGGGTTATTGTCTCATGAGCGGATACATATTTGAATGTA
TTAGAAAAATAAACAAATAGGGGTTCGCGGCACATTTCCCCGAAAAGTGCCAC
CTGACGTCTAAGAAACCATTATTATCATGACATTAACCTATAAAAATAGGCGTATC
30 ACGAGGCCCTTTCGTC

II. Transgene components

A. CREs:

- 5 Two different 6XCRE elements were used in the generation of the transgenes. The first is a synthetic CRE element (noted as “synthetic”) which was synthetically generated at DNA2.0, Menlo Park, CA. The synthetic sequence has attL1 and attR5 sites for 4-fragment Gateway cloning.
- 10 Multiple attempts were made to PCR clone out the 6X CRE from different vectors; however this was not possible due to a hairpin structure in the middle of the CREs. This fragment was then cloned synthetically. The second 6X CRE element used is a hybrid version generated by PCR cloning the CRE elements from a Clontech vector (Mountainview, CA, cat#PT3336-5). Clontech claims that there are 2x CREs in the
- 15 vector but sequence analysis indicated that there were actually three present. Two different PCR reactions were used to clone the fragment. Att sites for Gateway cloning as well as an EcoR1 site were introduced in the primers to enable the two fragments to be “glued” together and then recombined into the Invitrogen Gateway pDONR P1-P5r vector (Carlsbad, CA, cat #12537-100).

20

1. Synthetic CRE element

SEQ ID NO: 6:

25

5'AAATAATGATTTTATTTTGACTGATAGTGACCTGTTCGTTGCAACAAATT
GATGAGCAATGCTTTTTTATAATGCCAACTTTGTACAAAAAAGCAGGCTTA
CTGTGACAATTGCGTCATACTGTGACGTCTTTCAGACACCCCATGACG
TCAATGGGATTGACGTCAATGGGGTGTCTGAAAGACGTCACAGTATGACC
CGGGCTCGAGCCTCCTTGGCTGACGTCAGAGAGAGAGGGCCGGCCCCCTT
ACGTCAGAGGCGAGAATTGACAACCTTTGTATACAAAAGTTGAACGAGAA
30 ACGTAAAATGATATAAATATCAATATATTAAATTAGATTTTGCATAAAAAAC
AGACTACATAATACTGTAAAACACAACATATCCAGTCACTATG3'

2. Hybrid CRE: Standard cloning, PCR amplified the CRE element from the Clontech pCreLuc vector which contains 3XCREs, PCR primers have an EcoR1 restriction site at one end (CRE3X-B, CRE3X-C), and an att site, either attB1 or attB5r for Gateway cloning at the other end (CRE3X-A, CRE3X-D).
5 The two fragments were then combined to create one fragment with 6XCREs. The fragment gets recombined into the Gateway pDONR P1-P5r vector.

SEQ ID NO: 7: CRE forward primer A

5'-GGGGACAAGTTTGTACAAAAAAGCAGGCTTAGCACCAGACAGTGA-3'

SEQ ID NO: 8: CRE reverse primer B

5'-GGGAATTCGTTCTCCCATGACGTCA-3'

SEQ ID NO: 9: CRE forward primer C

5'-GGGAATTCGCACCAGACAGTGACGTC-3'

SEQ ID NO: 10: CRE reverse primer D

5'-GGGGACAACTTTTGTATACAAAGTTGTGTTCTCCCATGACGTCA-3'

SEQ ID NO: 11: Hybrid CRE sequence:

GCTTAGCACCAGACAGTGACGTCAGCTGCCAGATCCCATGGCCGTCATA
CTGTGACGTCTTTCAGACACCCCATTGACGTCAATGGGAGAACGAATTCG
CACCAGACAGTGACGTCAGCTGCCAGATCCCATGGCCGTCATACTGTGA
CGTCTTTCAGACACCCCATTGACGTCAATGGGAGAACCA

B. Human Growth Hormone Poly A Tail:

Human growth hormone poly A tail was PCR cloned from the vector pOGH (Nichols Institute Diagnostics (San Juan Capistrano, CA, Cat # 40-2205). Primer sequences have an attB3 site on the forward primer and an attB2 site on the reverse primer for recombination into the Gateway pDONR P3-P2 vector.

SEQ ID NO: 12: hGH forward primer:

5'-GGGGACAACCTTTGTATAATAAAGTTGGATCCCAAGGCCCAACTCC-3'

5 SEQ ID NO: 13: hGH reverse primer

5'-GGGGACCACTTTGTACAAGAAAGCTGGGTACAACAGGCATCTACT-3'

C. HSV TK Minimal Promoter:

10

HSV TK Minimal Promoter was PCR cloned from an in house vector called "661 CreLuc". Primer sequences have an attB5 site on the forward primer and an attB4 site on the reverse primer for recombination into the Gateway pDONR P5-P4 vector.

15 SEQ ID NO: 14: TK forward primer

5'-GGGGACAACCTTTGTATACAAAAGTTGTGGAACACGCAGATGCAGT-3'

SEQ ID NO: 15: TK reverse primer

5'-GGGGACAACCTTTGTATAGAAAAGTTGGGTGGATCTGCGGCACGCT-3'

20

D. Luciferase cDNA:

25 Luciferase cDNA was PCR cloned from the vector pGL4.10 (Promega Madison, WI Cat #E6651). Luc primer sequences have an attB4r site on the forward primer and an attB3r site on the reverse primer for recombination into the Gateway pDONR P4r-P3r vector.

SEQ ID NO: 16: luci forward primer

30

5'-GGGGACAACCTTTTCTATACAAAGTTGATGGAAGATGCCAAAAACA3'

SEQ ID NO: 17: luci reverse primer

5'-GGGGACAACCTTTATTATACAAAGTTGTTTACACGGCGATCTTGCC3'

III. Final Vector Construction

The four components of the transgene, CRE (either synthetic or hybrid version), HSVTK min, luciferase cDNA, and the hGH poly A tail in the pDONR vectors were recombined with the pDest2XMARs destination vector according to the standard Invitrogen protocol. The two transgenes were then sequenced, and transfected into CHOK1 to test for function. Both transgenes are functional in vitro as shown in Figure 5. Figure 5 shows the in vitro analysis of CreLuc vectors by forskolin induction.

IV. Transgenic Mouse Generation

A. Transgene preparation:

The quality of the transgene constructs were verified and confirmed by running aliquots on an agarose gel. No trace of degradation was observed. Finally, restriction analysis using diagnostic XhoI and PvuII and PstI restriction enzymes yielded the expected restriction profiles. The transgene plasmids were then digested by Acc65I and PmeI and the fragments containing the transgenes, were isolated from the vector backbones by running the digests on an agarose gel. The transgene fragments were then cut out of the gel, purified with a Qiaquick gel extraction kit (Qiagen, Valencia, CA cat#28706) and diluted before injection into fertilized oocytes. The purity and concentration of the isolated transgene were verified by agarose gel electrophoresis

SEQ ID NO: 18: CreLuc synthetic transgene (Acc65I/PmeI digest, 5550bp)

```
5'GTACCCGGGGGCGCGCCGGGATCCTTAATTAAAATTATCTCTAAGGCATGTGA
ACTGGCTGTCTTGGTTTTTCATCTGTACTTCATCTGCTACCTCTGTGACCTGAAAC
ATATTTATAATTCCATTAAGCTGTGCATATGATAGATTTATCATATGTATTTTCCTT
AAAGGATTTTTGTAAGAACTAATTGAATTGATACCTGTAAAGTCTTTATCACACTA
CCCAATAAATAATAAATCTCTTTGTTTCAGCTCTCTGTTTCTATAAATATGTACCAG
TTTTATTGTTTTTAGTGGTAGTGATTTTATTCTCTTTCTATATATACACACACAT
GTGTGCATTCATAAATATATACAATTTTTATGAATAAAAAATTATTAGCAATCAATA
TTGAAAACCACTGATTTTTGTTTATGTGAGCAAACAGCAGATTAAAAGGAATTCCT
GCAGGATCCTTAATTAAAGTTCTAGATCCAAGTTTGTACAAAAAAGCAGGCTTACT
GTCGACAATTGCGTCATACTGTGACGTCTTTCAGACACCCCATTTGACGTCAATG
GGATTGACGTCAATGGGGTGTCTGAAAGACGTACAGTATGACCCGGGGCTCGA
GCCTCCTTGGCTGACGTCAGAGAGAGAGAGGCCGGCCCTTACGTCAGAGGCGAG
```

AATTCGACAACTTTGTATACAAAAGTTGTGGAACACGCAGATGCAGTCGGGGCG
GCGCGGTCCCAGGTCCACTTCGCATATTAAGGTGACGCGTGTGGCCTCGAACA
CCGAGCGACCCTGCAGCGACCCGCTTAACAGCGTCAACAGCGTGCCGCAGATC
CACCCAACTTTTCTATACAAAAGTTGCTATGGAAGATGCCAAAAACATTAAGAAGG
5 GCCCAGCGCCATTCTACCCACTCGAAGACGGGACCGCCGGCGAGCAGCTGCA
CAAAGCCATGAAGCGCTACGCCCTGGTGCCCGGCACCATCGCCTTTACCGACG
CACATATCGAGGTGGACATTACCTACGCCGAGTACTTCGAGATGAGCGTTCCGGC
TGGCAGAAGCTATGAAGCGCTATGGGCTGAATACAAACCATCGGATCGTGGTGT
GCAGCGAGAATAGCTTGCACTTCTTCATGCCCGTGTTGGGTGCCCTGTTTCATCG
10 GTGTGGCTGTGGCCCCAGCTAACGACATCTACAACGAGCGCGAGCTGCTGAAC
AGCATGGGCATCAGCCAGCCCACCGTCGTATTTCGTGAGCAAGAAAGGGCTGCA
AAAGATCCTCAACGTGCAAAAGAAGCTACCGATCATACAAAAGATCATCATCATG
GATAGCAAGACCGACTACCGAGGGCTTCCAAAGCATGTACACCTTCGTGACTTCC
CATTTGCCACCCGGCTTCAACGAGTACGACTTCGTGCCCGAGAGCTTCGACCG
15 GGACAAAACCATCGCCCTGATCATGAACAGTAGTGGCAGTACCGGATTGCCCAA
GGGCGTAGCCCTACCGCACCGCACCGCTTGTGTCCGATTTCAGTCATGCCCGCG
ACCCCATCTTCGGCAACCAGATCATCCCCGACACCGCTATCCTCAGCGTGGTGC
CATTTACCCACGGCTTCGGCATGTTACCCACGCTGGGCTACTTGATCTGCGGCT
TTCGGGTCTGTGCTCATGTACCGCTTCGAGGAGGAGCTATTCTTGCGCAGCTTGC
20 AAGACTATAAGATTCAATCTGCCCTGCTGGTGCCCACTATTTAGCTTCTTCGC
TAAGAGCACTCTCATCGACAAGTACGACCTAAGCAACTTGCACGAGATCGCCAG
CGGCGGGGGCGCCGCTCAGCAAGGAGGTAGGTGAGGCCGTGGCCAAACGCTTC
CACCTACCAGGCATCCGCCAGGGCTACGGCCTGACAGAAACAACCAGCGCCAT
TCTGATCACCCCCGAAGGGGACGACAAGCCTGGCGCAGTAGGCAAGGTGGTGC
25 CTTTCTTCGAGGCTAAGGTGGTGGACTTGGACACCGGTAAGACACTGGGTGTG
AACCAGCGCGGCGAGCTGTGCGTCCGTGGCCCCATGATCATGAGCGGCTACGT
TAACAACCCCGAGGCTACAAACGCTCTCATCGACAAGGACGGCTGGCTGCACA
GCGGCGACATCGCCTACTGGGACGAGGACGAGCACTTCTTCATCGTGGACCGG
CTGAAGAGCCTGATCAAATACAAGGGCTACCAGGTAGCCCCAGCCGAACCTGGA
30 GAGCATCCTGCTGCAACACCCCAACATCTTCGACGCCGGGGTTCGCCGGCCTGC
CCGACGACGATGCCGGCGAGCTGCCCGCCGCGAGTCGTCTGTGCTGGAACACGG
TAAAACCATGACCGAGAAGGAGATCGTGGACTATGTGGCCAGCCAGGTTACAAC
CGCCAAGAAGCTGCGCGGTGGTGTGTTGTGTTCTGTGGACGAGGTGCCTAAAGGAC
TGACCGGCAAGTTGGACGCCCCGAAGATCCGCGAGATTCTCATTAAGGCCAAG
35 AAGGGCGGCAAGATCGCCGTGTAAACAACTTTGTATAATAAAGTTGCTGATCCC
AAGGCCCAACTCCCCGAACCACTCAGGGTCCTGTGGACAGCTCACCTAGCTGC
AATGGCTACAGGTAAGCGCCCCATAAATCCCTTTGGGCACAATGTGTCCTGAGG
GGAGAGGCAGCGACCTGTAGATGGGACGGGGGCACTAACCCTCAGGTTTGGG
GCTTCTGAATGTGAGTATCGCCATGTAAGCCCAGTATTTGGCCAATCTCAGAAA
40 GCTCCTGGTCCCTGGAGGGATGGAGAGAGAAAAACAACAGCTCCTGGAGCAG
GGAGAGTGCTGGCCTCTTGCTCTCCGGCTCCCTCTGTTGCCCTCTGGTTTCTCC
CCAGGCTCCCGGACGTCCCTGCTCCTGGCTTTTGGCCTGCTCTGCCTGCCCTG
GCTTCAAGAGGGCAGTGCTTCCCAACCATTCCTTATCCAGGCTTTTTGACAA
CGCTATGCTCCGCGCCCATCGTCTGCACCAGCTGGCCTTTGACACCTACCAGG
45 AGTTTGTAAAGCTCTTGGGGAATGGGTGCGCATCAGGGGTGGCAGGAAGGGGTG
ACTTTCCCCCGCTGGGAAATAAGAGGAGGAGACTAAGGAGCTCAGGGTTTTTCC
CGAAGCGAAAATGCAGGCAGATGAGCACACGCTGAGTGAGGTTCCAGAAAAG
TAACAATGGGAGCTGGTCTCCAGCGTAGACCTTGGTGGGCGGTCTTCTCCTAG
GAAGAAGCCTATATCCCAAAGGAACAGAAGTATTCATTCTGCAGAACCCCCAG
50 ACCTCCCTCTGTTTCTCAGAGTCTATTCCGACACCCTCCAACAGGGAGGAAACA

CAACAGAAATCCGTGAGTGGATGCCTTCTCCCCAGGCGGGGATGGGGGAGACC
TGTAGTCAGAGCCCCCGGGCAGCACAGCCAATGCCCGTCCTTCCCCTGCAGAA
CCTAGAGCTGCTCCGCATCTCCCTGCTGCTCATCCAGTCGTGGCTGGAGCCCG
TGCAGTTCCTCAGGAGTGTCTTCGCCAACAGCCTGGTGTACGGCGCCTCTGACA
5 GCAACGTCTATGACCTCCTAAAGGACCTAGAGGAAGGCATCCAAACGCTGATGG
GGGTGAGGGTGGCGCCAGGGGTCCCCAATCCTGGAGCCCCACTGACTTTGAGA
GCTGTGTTAGAGAAACACTGCTGCCCTCTTTTAGCAGTCAGGCCCTGACCCAA
GAGAACTCACCTTATTCTTCATTTCCCCTCGTGAATCCTCCAGGCCTTTCTCTAC
ACCCTGAAGGGGAGGGAGGAAAATGAATGAATGAGAAAGGGAGGGGAACAGTAC
10 CCAAGCGCTTGGCCTCTCCTTCTCTTCTTCACTTTGCAGAGGCTGGAAGATGG
CAGCCCCCGGACTGGGCAGATCTTCAAGCAGACCTACAGCAAGTTCGACACAA
ACTCACACAACGATGACGCACTACTCAAGAACTACGGGCTGCTCTACTGCTTCA
GGAAGGACATGGACAAGGTGAGACATTCTGCGCATCGTGCAGTGCCGCTCT
GTGGAGGGCAGCTGTGGCTTCTAGCTGCCCGGGTGGCATCCCTGTGACCCCTC
15 CCCAGTGCCTCTCCTGGCCCTGGAAGTTGCCACTCCAGTGCCCAACCAGCCTTG
TCCTAATAAAATTAAGTTGCATCATTTTGTCTGACTAGGCGTCCTTCTATAATATT
ATGGGGTGGAGGGGGGTGGTATGGAGCAAGGGGCAAGTTGGGAAGACAACCT
GTAGGGCCTGCGGGGTCTATTGGGAACCAAGCTGGAGTGCAGTGGCACAATCT
TGGCTCACTGCAATCTCCGCCTCCTGGGTTCAAGCGATTCTCCTGCCTCAGCCT
20 CCCGAGTTGTTGGGATTCCAGGCATGCATGACCAGGCTCAGCTAATTTTTGTTTT
TTTGGTAGAGACGGGGTTTCACCATATTGGCCAGGCTGGTCTCCAACCTCCTAAT
CTCAGGTGATCTACCCACCTTGGCCTCCCAAATTGCTGGGATTACAGGCGTGAA
CCACTGCTCCCTTCCCTGTCCTTCTGATTTTAAATAACTATAACCAGCAGGAGGA
CGTCCAGACACAGCATAGGCTACCTGGCCATGCCCAACCGGTGGGACATTTGA
25 GTTGTGTTGCTTGGCACTGTCCTCTCATGCGTTGGGTCCACTCAGTAGATGCCTG
TTGTACCCAGCTTTCTTGTACAAAGTGGGATCTAGACTAGAGTCATCAGTCAATA
TGTTACCCCCAAAAAAGCTGTTTGTTAACTTGTCAACCTCATTCTAAAATGTATAT
AGAAGCCCCAAAAGACAATAACAAAAATATTCTTGTAGAACAAAATGGGAAGAAT
GTTCCACTAAATATCAAGATTTAGAGCAAAGCATGAGATGTGTGGGGATAGACA
30 GTGAGGCTGATAAAATAGAGTAGAGCTCAGAAACAGACCCATTGATATATGTAA
GTGACCTATGAAAAAATATGGCATTTTACAATGGGAAAATGATGATCTTTTTCTT
TTTTAGAAAAACAGGGAAATATATTTATATGTAAAAAATAAAAGGGAACCCATATG
TCATACCATACACACAAAAAAATTCCAGTGAATTATAAGTCTAAATGGAGAAGGC
AAAACTTTAAATCTTTTAGAAAAATAATATAGAAGCATGCCATCAAGACTTCAGTGT
35 AGAGAAAAATTTCTTATGACTCAAAGTCCTAACCACAAAGAAAAGATTGTTAATTA
GATTGCATGAATATTAAGACTTATTTTTTAAAATTAATAAACCATTAAGAAAAGTCA
GGCCATAGAATGACAGAAAATATTTGCAACACCCCAGTAAAGAGAATTGTAATAT
GCAGATTATAAAAAAGAAGTCTTACAAATCAGTAAAAAATAAACTAGACAAAAATT
TGAACAGATGAAAGAGAAACTCTAAATAATCATTACACATGAGAACTCAATCTC
40 AGAAATCAGAGAACTATCATTGCATATACACTAAATTAGAGAAATATTAAGGCT
AAGTAACATCTGTGGCTTAATTAACAGTAGGATGACTGTTT 3'

SEQ ID NO: 19: CreLuc hybrid transgene sequence (Acc65I/PmeI digest, 5562bp)

45 5'GTACCCGGGGGCGCGCCGGGATCCTTAATTAAATTATCTCTAAGGCATGTGA
ACTGGCTGTCTTGGTTTTTCATCTGTACTTCATCTGCTACCTCTGTGACCTGAAAC
ATATTTATAATTCCATTAAGCTGTGCATATGATAGATTTATCATATGTATTTTCCTT
AAAGGATTTTGTGAAGAACTAATTGAATTGATACCTGTAAAGTCTTTATCACACTA

CCCAATAAATAATAAATCTCTTTGTTTCAGCTCTCTGTTTCTATAAATATGTACCAG
TTTTATTGTTTTTAGTGGTAGTGATTTTATTCTCTTTCTATATATATACACACACAT
GTGTGCATTTCATAAATATATACAATTTTTATGAATAAAAAATTATTAGCAATCAATA
TTGAAAACCACTGATTTTTGTTTATGTGAGCAAACAGCAGATTAAAAGGAATTCCT
5 GCAGGATCCTTAATTAAGTTCTAGATCCAAGTTTGTACAAAAAAGCAGGCTTAGC
ACCAGACAGTGACGTCAGCTGCCAGATCCCATGGCCGTCATACTGTGACGTCTT
TCAGACACCCCATTGACGTCAATGGGAGAACGAATTCGCACCAGACAGTGACGT
CAGCTGCCAGATCCCATGGCCGTCATACTGTGACGTCTTTCAGACACCCCATTG
ACGTCAATGGGAGAACACAACCTTTGTATACAAAAGTTGTGGAACACGCAGATGC
10 AGTCGGGGCGGCGCGGTCCCAGGTCCACTTCGCATATTAAGGTGACGCGTGTG
GCCTCGAACACCGAGCGACCCTGCAGCGACCCGCTTAACAGCGTCAACAGCGT
GCCGCAGATCCACCCAACTTTTCTATACAAAGTTGCTATGGAAGATGCCAAAAAC
ATTAAGAAGGGGCCAGCGCCATTCTACCCACTCGAAGACGGGACCGCCGGCGA
GCAGCTGCACAAAGCCATGAAGCGCTACGCCCTGGTGCCCGGCACCATCGCCT
15 TTACCGACGCACATATCGAGGTGGACATTACCTACGCCGAGTACTTCGAGATGA
GCGTTCGGCTGGCAGAAGCTATGAAGCGCTATGGGCTGAATACAAACCATCGG
ATCGTGGTGTGCAGCGAGAATAGCTTGCAGTTCTTCATGCCCGTGTTGGGTGCC
CTGTTTCATCGGTGTGGCTGTGGCCCCAGCTAACGACATCTACAACGAGCGCGA
GCTGCTGAACAGCATGGGCATCAGCCAGCCACCGTCGTATTCGTGAGCAAGA
20 AAGGGCTGCAAAAGATCCTCAACGTGCAAAAGAAGCTACCGATCATACAAAAGA
TCATCATCATGGATAGCAAGACCGACTACCAGGGCTTCCAAAGCATGTACACCT
TCGTGACTTCCCATTTGCCACCCGGCTTCAACGAGTACGACTTCGTGCCCGAGA
GCTTCGACCGGGACAAAACCATCGCCCTGATCATGAACAGTAGTGGCAGTACC
GGATTGCCCAAGGGCGTAGCCCTACCGCACCGCACCGCTTGTGTCCGATTGAG
25 TCATGCCCGCGACCCCATCTTCGGCAACCAGATCATCCCCGACACCGCTATCCT
CAGCGTGGTGCCATTTACACACGGCTTCGGCATGTTACACACGCTGGGCTACTT
GATCTGCGGCTTTCGGGTTCGTGCTCATGTACCGCTTCGAGGAGGAGCTATTCTT
GCGCAGCTTGCAAGACTATAAGATTCAATCTGCCCTGCTGGTGCCACACTATT
TAGCTTCTTCGCTAAGAGCACTCTCATCGACAAGTACGACCTAAGCAACTTGAC
30 GAGATCGCCAGCGGCGGGGCGCCGCTCAGCAAGGAGGTAGGTGAGGCCGTGG
CCAAACGCTTCCACCTACCAGGCATCCGCCAGGGCTACGGCCTGACAGAAACA
ACCAGCGCCATTCTGATCACCCCCGAAGGGGACGACAAGCCTGGCGCAGTAGG
CAAGGTGGTGCCCTTCTTCGAGGCTAAGGTGGTGGACTTGGACACCGGTAAGA
CACTGGGTGTGAACCAGCGCGGCGAGCTGTGCGTCCGTGGCCCCATGATCATG
35 AGCGGCTACGTTAACAACCCCGAGGCTACAAACGCTCTCATCGACAAGGACGG
CTGGCTGCACAGCGGCGACATCGCCTACTGGGACGAGGACGAGCACTTCTTCA
TCGTGGACCGGCTGAAGAGCCTGATCAAATACAAGGGCTACCAGGTAGCCCCA
GCCGAACCTGGAGAGCATCCTGCTGCAACACCCCAACATCTTCGACGCCGGGGT
CGCCGGCCTGCCCGACGACGATGCCGGCGAGCTGCCCGCCGAGTCGTGCTG
40 CTGGAACACGGTAAAACCATGACCGAGAAGGAGATCGTGGACTATGTGGCCAG
CCAGGTTACAACCGCCAAGAAGCTGCGCGGTGGTGTGTTGTGTTTCGTGGACGAGG
TGCCTAAAGGACTGACCGGCAAGTTGGACGCCCCGCAAGATCCGCGAGATTCTC
ATTAAGGCCAAGAAGGGCGGCAAGATCGCCGTGTAAACAACCTTTGTATAATAAA
GTTGCTGATCCCAAGGCCCAACTCCCCGAACCACTCAGGGTCCTGTGGACAGC
45 TCACCTAGCTGCAATGGCTACAGGTAAGCGCCCCTAAAATCCCTTTGGGCACAA
TGTGTCCTGAGGGGAGAGGCAGCGACCTGTAGATGGGACGGGGGCACTAACC
CTCAGGTTTGGGGCTTCTGAATGTGAGTATCGCCATGTAAGCCCAGTATTTGGC
CAATCTCAGAAAGCTCCTGGTCCCTGGAGGGATGGAGAGAGAAAAACAAACAG
CTCCTGGAGCAGGGAGAGTGCTGGCCTCTTGCTCTCCGGCTCCCTCTGTTGCC
50 CTCTGGTTTCTCCCCAGGCTCCCGGACGTCCCTGCTCCTGGCTTTTGGCCTGCT

CTGCCTGCCCTGGCTTCAAGAGGGCAGTGCCTTCCCAACCATTCCCTTATCCAG
GCTTTTTTGACAACGCTATGCTCCGCGCCCATCGTCTGCACCAGCTGGCCTTTGA
CACCTACCAGGAGTTTGTAAAGCTCTTGGGGAATGGGTGCGCATCAGGGGTGGC
AGGAAGGGGTGACTTTCCCCCGCTGGGAAATAAGAGGAGGAGACTAAGGAGCT
5 CAGGGTTTTTCCCGAAGCGAAAATGCAGGCAGATGAGCACACGCTGAGTGAGG
TTCCCAGAAAAGTAACAATGGGAGCTGGTCTCCAGCGTAGACCTTGGTGGGCG
GTCCTTCTCCTAGGAAGAAGCCTATATCCCAAAGGAACAGAAGTATTCATTCTG
CAGAACCCCCAGACCTCCCTCTGTTTCTCAGAGTCTATTCCGACACCCTCCAAC
AGGGAGGAAACACAACAGAAATCCGTGAGTGGATGCCTTCTCCCCAGGCGGGG
10 ATGGGGGAGACCTGTAGTCAGAGCCCCCGGGCAGCACAGCCAATGCCCGTCCT
TCCCCTGCAGAACCTAGAGCTGCTCCGCATCTCCCTGCTGCTCATCCAGTCGTG
GCTGGAGCCCGTGACAGTTCCTCAGGAGTGTCTTCGCCAACAGCCTGGTGTACG
GCGCCTCTGACAGCAACGTCTATGACCTCCTAAAGGACCTAGAGGAAGGCATCC
AAACGCTGATGGGGGTGAGGGTGGCGCCAGGGGTCCCCAATCCTGGAGCCCC
15 ACTGACTTTGAGAGCTGTGTTAGAGAAACACTGCTGCCCTCTTTTTAGCAGTCAG
GCCCTGACCCAAGAGAACTCACCTTATTCTTCATTTCCCCTCGTGAATCCTCCAG
GCCTTTCTCTACACCCTGAAGGGGAGGGAGGAAAATGAATGAATGAGAAAGGG
AGGGAACAGTACCCAAGCGCTTGGCCTCTCCTTCTCTTCCTTCACTTTGCAGAG
GCTGGAAGATGGCAGCCCCCGGACTGGGCAGATCTTCAAGCAGACCTACAGCA
20 AGTTCGACACAACTCACACAACGATGACGCACTACTCAAGAACTACGGGCTGC
TCTACTGCTTCAGGAAGGACATGGACAAGGTGAGACATTCTGCGCATCGTGC
AGTGCCGCTCTGTGGAGGGCAGCTGTGGCTTCTAGCTGCCCGGGTGGCATCCC
TGTGACCCCTCCCCAGTGCCTCTCCTGGCCCTGGAAGTTGCCACTCCAGTGCC
CACCAGCCTTGTCTAATAAAATTAAGTTGCATCATTTTGTCTGACTAGGCGTCC
25 TTCTATAATATTATGGGGTGGAGGGGGGTGGTATGGAGCAAGGGGCAAGTTGG
GAAGACAACCTGTAGGGCCTGCGGGGTCTATTGGGAACCAAGCTGGAGTGCAG
TGGCACAATCTTGGCTCACTGCAATCTCCGCCTCCTGGGTTCAGCGATTCTCC
TGCCTCAGCCTCCCGAGTTGTTGGGATTCCAGGCATGCATGACCAGGCTCAGCT
AATTTTTGTTTTTTGGTAGAGACGGGGTTTCACCATATTGGCCAGGCTGGTCTC
30 CAACTCCTAATCTCAGGTGATCTACCCACCTTGGCCTCCCAAATTGCTGGGATTA
CAGGCGTGAACCACTGCTCCCTTCCCTGTCCTTCTGATTTTAAAATAACTATACC
AGCAGGAGGACGTCCAGACACAGCATAGGCTACCTGGCCATGCCCAACCGGTG
GGACATTTGAGTTGTTTGCTTGGCACTGTCCTCTCATGCGTTGGGTCCACTCAG
TAGATGCCTGTTGTACCCAGCTTTCTTGTACAAAGTGGGATCTAGACTAGAGTCA
35 TCAGTCAATATGTTACCCCCAAAAAGCTGTTTGTTAACTTGTCAACCTCATTCTA
AAATGTATATAGAAGCCCCAAAAGACAATAACAAAAATATTCTTGTAGAACAAAATG
GGAAAGAATGTTCCACTAAATATCAAGATTTAGAGCAAAGCATGAGATGTGTGG
GGATAGACAGTGAGGCTGATAAAATAGAGTAGAGCTCAGAAACAGACCCATTGA
TATATGTAAGTGACCTATGAAAAAATATGGCATTTTACAATGGGAAAATGATGAT
40 CTTTTCTTTTTTAGAAAAACAGGGAAATATATTTATATGTAAAAAATAAAAGGGA
ACCCATATGTCATACCATAACACAAAAAAATTCCAGTGAATTATAAGTCTAAATG
GAGAAGGCAAACTTTAAATCTTTTAGAAAATAATATAGAAGCATGCCATCAAGA
CTTCAGTGTAGAGAAAAATTTCTTATGACTCAAAGTCCTAACCACAAAGAAAAGA
TTGTTAATTAGATTGCATGAATATTAAGACTTATTTTTAAAATTA AAAAACCATTAA
45 GAAAAGTCAGGCCATAGAATGACAGAAAATATTTGCAACACCCCAAGTAAAGAGA
ATTGTAATATGCAGATTATAAAAAGAAGTCTTACAAATCAGTAAAAAATAAACTA
GACAAAAATTTGAACAGATGAAAGAGAACTCTAAATAATCATTACACATGAGAA
ACTCAATCTCAGAAATCAGAGAACTATCATTGCATATACACTAAATTAGAGAAATA
TTAAAAGGCTAAGTAACATCTGTGGCTTAATTAAAACAGTAGGATGACTGTTT

B. Founder generation:

The transgenes were separately microinjected into the pronuclei of FVB/Taconic one-celled embryos. General strategies for generating transgenic (Tg) animals has been well described and are known to the skilled artisan.

Using a PCR genotyping strategy, transgenic positive founder mice were identified shortly after birth through tail biopsies and a PCR for genotype. The selected primer pair allows the amplification of a short DNA sequence within the transgene sequence, yielding a specific 405-bp PCR product. PCR reaction mix is 1X Qiagen PCR buffer with MgCl₂ (Qiagen, Valencia, CA cat #201205), 0.2mM dNTPs, 0.4uM primers, 0.5 units of Taq Polymerase. Cycling conditions are as follows: 1 cycle at 94°C; 5min, 35 cycles at 94°C; 45sec/ 57°C; 45sec/72°C; 1 min , 1 cycle at 72°C; 5min

The genotyping primers are as follows:

SEQ ID NO: 20: Luc2-forward primer:

5'-GAAGATGCCAAAAACATTAAGAAG-3'

SEQ ID NO: 21: Luc2-reverse primer:

5'-GATCTTTTGCAGCCCTTTCT -3'

For the synthetic CreLuc transgene, 39 transgene positive founders were produced (out of 181 born), and for the hybrid CreLuc transgene, 73 transgene positives (out of 244 born). The transgenic mouse lines are abbreviated as CreLuc.

V. Founder expression analysis:

Figure 5A is a comparison of the effects of different phosphodiesterase (PDE) inhibitors on plasma cAMP levels in normal mice (reproduced from Cheng et al, JPET, 280(2):621-626, 1997). The mice were dosed p.o. with compound or vehicle(VEH) in groups of four animals. Twenty minutes later, blood was collected

and processed for cAMP measurements. Based on this published experiment, rolipram was chosen as a general inducer along with forskolin for the CreLuc founder screen for its ability and its availability to increase cAMP levels. The additive effect of the two compounds when dosed together, provide a larger window of induction when
5 examining luciferase levels in the animals.

Figure 5B is a cAMP assay on plasma levels in wild-type mice following various treatment with forskolin (F), rolipram (R) or a combination of forskolin and rolipram (F/R). Two month old FVB/Tac females were dosed i.p. with one of the following:
10 Group A: vehicle control, 1% dimethyl sulfoxide (DMSO) in phosphate buffered saline (PBS), Group B: 5mg/kg forskolin (Sigma, St Louis MO, cat#F6886), and Group C: 10mg/kg rolipram (Sigma St Louis, MO, cat# R6520), Group D: 5mg/kg forskolin plus 10mg/kg rolipram, Group E: 5mg/kg water soluble forskolin (Calbiochem Gibbstown, NJ cat#344273), and Group F: 5mg/kg water soluble forskolin plus 10mg/kg rolipram.
15 Group G received no treatment. Compared to vehicle alone, rolipram and/or forskolin were statistically significance in their elevation of circulating cAMP levels. Forskolin alone did not significantly increase circulating cAMP over the vehicle control. The combination of water soluble forskolin (H2OF) and rolipram together had the greatest elevation in cAMP plasma levels and therefore indicated this treatment would be
20 optimal for screening CreLuc founders for expression of the transgene in vivo.

All founders (founders listed in Tables 5-8 and representative founders shown in Figure 6) were tested for transgene expression by dosing the mice i.p. with a combination of 10mg/kg rolipram (Sigma St Louis, MO, cat# R6520) and 5mg/kg
25 water soluble forskolin (Calbiochem Gibbstown, NJ cat#344273). Four hours later, the mice were bioimaged using the Xenogen IVIS Lumina system. The mice were anesthetized with isoflurane, and 250mg/kg luciferin was injected s.c. Eight minutes after the luciferin injection, mice were imaged. Of the 39 transgene positive synthetic CreLuc transgene, 34 or 87% were expression positive. For the hybrid transgene, 55
30 out of 71 or 78% of the mice tested were expression positive. Expression was considered to be positive when luciferase levels of any tissue are above the background levels for the system.

Various forskolin rolipram response profiles assayed by bioimaging of different independent lines of CreLuc mice is illustrated in Figure 6. In the uninduced state we observed either detectable basal bioimaging signals as illustrated in line 90 and 219 or no basal activity as illustrated in lines 44, 28, and 187. After forskolin rolipram induction the bioimaging pattern of CreLuc reporter induction occurred in single or multiple tissues and was unique in each independent CreLuc line. Many of the CreLuc lines had one or more tissues expressing detectable luciferase bioimages in therapeutically relevant areas.

Following a bioimaging screen of CreLuc expression in founders with both basal and forskolin rolipram induction, CreLuc lines that met a filter window of 5x bioimaging induction in one or more tissues (2x in the brain) were analyzed in more detail by assaying the induced levels of luciferase enzyme in tissue extracts.

For tissue expression analysis, animals were either left untreated for a baseline signal, or were induced for 4 hours. For induction, animals were dosed i.p. with a combination of 5mg/kg Forskolin (Calbiochem, Gibbstown, NJ cat#344273) and 10mg/kg Rolipram (Sigma St Louis, MO, cat# R6520) with 1%DMSO (Sigma, St Louis, MO, cat#D2650) in Dulbecco's PBS (Invitrogen, Carlsbad, CA, cat#14040). Four hours later, the animals were sacrificed by CO₂ and the various organs were removed and frozen on dry ice. For the luciferase assay, the Luciferase Assay System (Promega, Madison, WI, cat#E1500) was used. The tissues were homogenized in 1ml of Cell Culture Lysis Buffer (Promega, Madison, WI, cat#E1531), incubated on ice for 5 minutes, and then spun in a centrifuge for 5 minutes. 20µl the supernatant was used in the assay as per manufacturers instructions. Protein concentrations of the samples were measured with the DC protein assay kit (BioRad, Hercules, CA, cat#500-0111). Data is shown in Tables 5 – 8 and is presented as fold induction (induced RLU/ug protein divided by basal RLU/ug protein). RLU induced is the RLU/ug protein of the induced tissue sample.

5 **TABLE 5: Luciferase Tissue Expression in Brain, Spleen and Kidney**

Line	Tissue					
	Brain		Spleen		Kidney	
	Fold Induction	RLU Induced	Fold Induction	RLU Induced	Fold Induction	RLU Induced
561-24	1.6	1.55	66.6	71.3	2	60.9
561-44	1.8	1.2	7.3	98.7	0.2	5.0
561-90	4.5	1.0	20.6	12.1	21.5	132.5
562-11	3.9	0.6	11.1	4.7	17.4	1.4
562-16	0.9	2.7	18.2	5.4	147.5	82.0
562-31	1.0	48.9	3.0	3.1	3.0	119.3
562-33	1.0	3.3	1.0	2.1	80.6	53.8
562-64	1.0	3.7	39.5	565.2	45.9	47.7
562-69	2.5	1.7	52.2	44.8	30.6	36.1
562-175	2.1	81.3	0.5	0.4	3.6	12.9
562-180	1.2	0.2	38.9	6.9	131.9	23.8
562-184	1.0	2.5	11.0	26.0	6.0	71.2
562-187	3.1	3.0	9.6	6.0	17.7	4.2
562-219	2.0	4.8	2.0	3.8	2.0	103.0
562-229	1.4	3.1	76.1	111.9	302.6	154.3

5

Table 6: Luciferase Tissue Expression in Liver, Thymus and Pancreas

Line	Tissue					
	Liver		Thymus		Pancreas	
	Fold Induction	RLU Induced	Fold Induction	RLU Induced	Fold Induction	RLU Induced
561-24	29.4	25	10.6	93.4	152	108
561-44	16.2	3.5	2.3	3.2	4.5	2.0
561-90	14.5	19.2	1.3	11.8	117.8	78.6
562-11	1.5	0.5	1.5	3.6	28.7	26.4
562-16	89.9	12.9	2.7	9.7	753.7	442.6
562-31	768.0	132.6	2.0	11.4	204.0	199.6
562-33	53.6	40.3	3.0	4.9	0.6	0.9
562-64	62.8	104.3	4.5	505.5	321.7	900.7
562-69	115.7	108.6	0.6	18.5	508.7	197.3
562-175	17.9	2.5	3.2	3.5	0.4	0.5
562-180	259.7	21.3	7.3	4.9	889.0	130.5
562-184	9.0	7.7	1.0	16.0	118.0	76.7
562-187	18.7	105.3	4.9	23.6	496.8	155.6
562-219	1.0	2.2	1.0	3.1	15.0	20.3
562-229	337.7	206.0	33.5	157.7	258.0	503.1

Table 7: Luciferase Tissue Expression in Heart, Lungs and Spinal Cord

Line	Tissue					
	Heart		Lungs		Spinal Cord	
	Fold Induction	RLU Induced	Fold Induction	RLU Induced	Fold Induction	RLU Induced
561-24	2.1	2.4	6.2	10.9	1.1	14.5
561-44	2.4	0.7	0.6	1.9	1.2	8.7
561-90	3.0	1.5	6.8	7.7		
562-11	11.1	2.9	2.5	7.4		
562-16	1.1	1.9	3.2	7.3		
562-31	4.0	2.2	26.0	17.4		
562-33	2.1	1.9	8.1	10.0		
562-64	3.9	8.4	3.2	36.7		
562-69	0.6	0.7	0.5	1.1		
562-175	1.4	0.6	12.5	228.0		
562-180	1.5	0.6	1.8	0.9		
562-184	1.0	0.6	1.0	2.1		
562-187	1.3	1.4	0.4	5.0	0.9	2.9
562-219	2.0	3.6	27.0	318.4		
562-229	5.2	7.4	17.3	40.8	0.3	7.6

Table 8: Luciferase Tissue Expression in Skin and Intestine

Line	Tissue			
	Skin		Intestine	
	Fold Induction	RLU Induced	Fold Induction	RLU Induced
561-24	11.1	59	2.6	16
561-44				
561-90				
562-11	4.5	77.2		
562-16	7	114	0.2	7.7
562-31	1	3.2	34	80.5
562-33	27	142	35	31.2
562-64	1.1	69	10.6	49.9
562-69	1.3	37.7	0.5	1.1
562-175	0.3	15.7	0	1.4
562-180	1	14.4		
562-184	1	6.6	1	
562-187				
562-219	4	51	6	4.2
562-229	12.1	144	1.5	43.3

5

This detailed level of CreLuc expression in each line allowed the selection of optimal lines with the highest expression levels in target tissues and the selection of lines with unique single or multiple patterns of expression. The enzyme assay also insured that CreLuc expression was clearly associated with a particular tissue, In general, as designed in the transgene we observed luciferase expression in nearly all tissues across the majority of the lines. The top 3- 5 lines from key pharmaceutically relevant target tissues (brain, pancreas, lung, spleen) were selected for further reference compound response profiling using GPCR ligands.

10

VI. Analysis of the Effects of Insulator Elements on Expression Levels

In Table 9, the influence of adding insulator elements is compared across 11 different transgenes and including 225 independent transgenic lines. The production of functional transgenic lines ranges typically from 10-30% (historical data taken from our laboratory). As can be seen in Table 9, 10 of 11 transgenes and 190 lines obtained expression levels at the high end of this range (mean number of functional lines is 37 or 32.7%) due to the inclusion of the hGH minigene or functional enhancer. Significantly, the inclusion of insulator elements more than doubled this expression frequency and for 112 lines 78% express a functional reporter. The sample sizes for both comparisons and the diversity of the transgenes are large enough to overcome normal biological variation due to random integration of transgenes in the genome.

Generating transgenic bioimaging reporter lines is a low frequency, high attrition rate process. A large number of independent founder lines (50-100) are needed to obtain lines with the desired detection level for the reporter and to have the reporter expressed in a wide selection of target tissues appropriate for drug candidate screening. Typically 5% or less of lines are sufficient for drug screening applications and therefore the transgene containing the reporter has to include many elements to maximize the expression in each transgenic line generated. The insulator elements coupled to other expression enhancers (hGH minigene) has been demonstrated to yield a majority of transgenic lines that express the transgene and are key to finding transgenic lines with an expression pattern that is applicable to drug discovery (especially for various GPCR signaling pathways to detect a reduction in signal upon ligand to receptor binding).

Table 9. Comparison of the frequency of generating functional transgenic lines generated with either the hGH minigene and insulator elements

Transgene	Promoter	Functional Enhancer	Founder Lines			Insulator
			Total #	# Positive	% Positive	
CD11ArtT A x tetEGFP **	CD11A / tet p	hGH minigene	10	5	50%	none
CD11ArtT A X tetIL13 **	CD11A / tetp	hGH minigene	9	3	33%	none
CC10rtTA x tetEGFP **	CC10 / tetp	hGH minigene	11	7	63%	none
hCXCR5-CD11a	CD11a	hGH minigene	10	5	50%	none
hCXCR5-CD11b	CD11b	hGH minigene	16	5	50%	none
MBP luci	MBP	hGH minigene	35	6	17%	none
GP1 CRE	GP1	hGH minigene	7	2	29%	none
GP1 GFP	GP1	hGH minigene	7	4	57%	none
PF4 GFP	PF4	hGH minigene	8	0	0	none
Totals			113	37	Ave. = 32.7%	
CreLuc Synthetic	HSV TK + 6X CRE	hGH minigene	39	33	85%	yes *
CreLuc hybrid	HSV TK + 6X CRE	hGH minigene	73	55	75%	yes *
Totals			112	88	Ave. = 78%	

* beta globin / IFN, ** binary model

Legend: CD11 ArtTA, CD11a promoter expressing the reverse activator for tetracycline transactivator; TetEGFP is the tetracycline promoter expressing the EGFP reporter; TetIL13 is the tetracycline promoter expressing IL13; CC10rtTA; the CC10 lung specific promoter expressing the reverse activator for the tetracycline transactivator; hCXCR5- CD11a or b is the human CXCR5 promoter expressing the CD11a or CD11b coding cDNA; MBP luci is the myelin basic protein (MBP) promoter expressing luciferase (luci); GP1 and PF4- GFP; GP1 and PF4 are platelet specific promoters expressing the GFP reporter.

Deposit of Material:

Preparation of Mouse Embryonic Fibroblasts: Day 12 embryos were removed from pregnant females. The embryos were dissected out of the uterus in a 10cm dish containing phosphate-buffered saline (PBS) (Invitrogen, Carlsbad, CA, cat#14040). They were then transferred to a new dish of PBS. The embryos were teased apart with forceps and the heart and liver were removed. The tissue pieces were then incubated in 2-3 mls of cold 0.25% trypsin (Invitrogen, Carlsbad, CA, cat#25200) overnight at 4°C. The next day, most of the trypsin was aspirated off and the tissue pellet was incubated for 30 minutes at 37°C. 2-3 ml of media, DMEM (Invitrogen, Carlsbad, CA, cat#11965), 5% heat inactivated FBS (Invitrogen, Carlsbad, CA, cat#16140), 1% pen strep (Invitrogen, Carlsbad, CA, cat# 15140) was added and pipetted up and down to break up the remaining clumps of tissue. The cell suspension was then plated onto one T75 tissue culture flasks. Cells were cultured for 3 passages and then frozen down in 1 ml vials, 1-2 million cells per vial, in freezing media, 15% heat inactivated FBS (Invitrogen, Carlsbad, CA, cat#16140), 5% DMSO (Sigma, St Louis, MO, cat#D2650), DMEM (Invitrogen, Carlsbad, CA, cat#11965).

The following materials have been deposited with the American Type Culture Collection, 10801 University Blvd., Manassas, Va. 20110-2209, USA (ATCC):

Table 10. Deposit of Material

Line	ATCC Deposit #	Identification Reference by Depositer
562-11	PTA-10536	Mouse embryonic fibroblast:FVB/Creluc11
562-187	PTA-10537	Mouse embryonic fibroblast:FVB/Creluc187
562-229	PTA-10538	Mouse embryonic fibroblast:FVB/Creluc229

The deposit was made under the provisions of the Budapest Treaty on the International Recognition of the Deposit of Microorganisms for the Purpose of Patent Procedure and the Regulations thereunder (Budapest Treaty). This assures

maintenance of a viable culture of the deposit for 30 years from the date of deposit. The deposit will be made available by ATCC under the terms of the Budapest Treaty, and subject to an agreement between sanofi-aventis US Inc. and ATCC, which assures permanent and unrestricted availability of the progeny of the culture of the deposit to the public upon issuance of the pertinent U.S. patent or upon laying open to the public of any U.S. or foreign patent application, whichever comes first, and assures availability of the progeny to one determined by the U.S. Commissioner of Patents and Trademarks to be entitled thereto.

The assignee of the present application has agreed that if a culture of the materials on deposit should die or be lost or destroyed when cultivated under suitable conditions, the materials will be promptly replaced on notification with another of the same. Availability of the deposited material is not to be construed as a license to practice the invention in contravention of the rights granted under the authority of any government in accordance with its patent laws.

VII. In Vivo and Ex Vivo Screening Assays

Figure 7 shows embodiments of the invention by representing the flow of CreLuc screening assays either in vivo or ex vivo. For in vivo assays, the animals are dosed with luciferin only and then bioimaged on an IVIS bioimager to get baseline signal. The animals are then dosed with compound, and then bioimaged again at a set timepoint to examine the effects of the compound on the luciferase signal. Data can be analyzed as the fold induction of the induced signal versus the baseline signal. For the ex vivo assays, mice are dosed with either compound or vehicle, and then at a set timepoint, tissues are harvested and frozen. The tissues can then be homogenized and a luciferase assay such as Promega Luciferase assay system (Promega cat #E1500) run on the tissue extracts and the assay read on a luminometer such as the Topcount. Data can be analyzed as the relative light units per μg protein, comparing vehicle versus compound groups.

Additional aspects and details of the invention will be apparent from the following examples, which are intended to be illustrative and are not meant to limit the scope of the invention in any way.

5

A. Whole Brain Slices Experiments

Whole brain slices were used to show the effects of compounds on luminescence in CreLuc mice (Figure 8A, 8B). In general, brains from CreLuc mice can be harvested and then either cut into slices on a microtome or various subregions can also be dissected out and then sliced on a microtome. The slices can then be incubated with different compounds such as forskolin and rolipram, or GPCR specific compounds. Luciferase levels can then be measured with an IVIS bioimager or a photomultiplier with a carousel to take quantitative timecourse measurements.

15

Brain slices from the hippocampal region of CreLuc mice (line # 187) were incubated with 1uM isoproterenol (Sigma, cat#5627) for 2 days. Isoproterenol (ISO) is a Gs agonist that activates the cAMP pathway (see Figure 1 for pathway schematic) and induces luciferase expression levels in CreLuc mice. (Figure 8A, bottom panel). Thus, isoproterenol was used to increase the luciferase signal window which allows for screening GPCR specific compounds such as Gi agonists. AMN082 (AMN, Tocris Bioscience, Cat. No. 2385) is an mGluR7 agonist and mGluR7 is a receptor that is coupled to Gi. When brain slices were treated with 1uM ISO in combination with 1uM AMN082, the luciferase signal was diminished in comparison to ISO treatment alone.

25

Compound activity can also be quantitated over time in a viable brain slice and luciferase expression can be visualized by a bioimager. In Figure 8B, brain slices from CreLuc mice (line 44) were responsive to forskolin. Forskolin is a general activator of cAMP. Treatment of the brain slice with 50uM forskolin (arrow marks the time of forskolin was added to the slice) caused an induction in time dependent manner compared to vehicle.

30

B. In Vitro Cell Culture Experiments

1. Primary Cortical Neurons

5 A general schematic of the workflow for the isolation and compound treatment of primary neuronal cells from the CreLuc model is depicted in Figure 9. Cortices or other brain subregions are removed from embryos, usually embryonic (E) day 17 or 18, but also as early as E14. The individual neurons are then isolated and plated onto 96 well assay plates. Assays are run either in triplicate or quadruplicate. The
10 assays are run from day 3 in culture up until day 7 or day 10 in culture depending on the expression levels of the target GPCR. Compounds are added to the plates then at a set timepoint, Bright Glo (Promega, Cat # E2610) is added to the plates and the assay is then run on a TopCount or other luminometer.

15 a) Gs Modulation via β Adrenergic Receptors ($\text{AD}\beta\text{R}$) and D1 Dopamine Receptors (DRD1)

Primary cell cultures from the CreLuc mice may be used to screen compounds that
20 are capable of modulating Gs coupled receptors. Neurons were isolated from the cortices of line 187 E18 embryos. Line 187 has been shown to have inducible levels of luciferase in both the brain and spinal cord by whole tissue luciferase assays (see text below; Figure 24). The assay was run in triplicates on day three in culture. Forskolin was used at 5 μM , rolipram at 10 μM and the Gs agonists, $\text{AD}\beta\text{R}$ -
25 isoproterenol; and Dopamine- SKF82958 (MFG, Cat # C130) at 10 μM . At 4 hours, Bright Glo was added and the assays run on a Topcount. Forskolin is a general activator of cAMP. Rolipram does not activate cAMP, instead, rolipram blocks the breakdown of cAMP, thereby stabilizing the cAMP levels. Forskoline and rolipram act synergistically to increase luciferase expression. An induction of over 100 fold was
30 observed with a combination of forskolin and rolipram. Increased luci levels were seen with $\text{AD}\beta\text{R}$ agonist isoproterenol (11 fold) vs. DMSO control. While the D1 DR agonist, SKF82958, did not significantly increase luciferase levels, a 2-fold induction was observed when SKF82958 was administered in combination with rolipram when

compared to rolipram alone. Thus, line 187 is a useful model for screening compounds to determine if the compound is able to modulate Gs coupled receptors.

5 **b) Gq Modulation via Prokineticin 2 Receptor (PROK2R)**

Primary cell cultures from the CreLuc mice may be used to screen compounds that are capable of modulating Gs coupled receptors. The PROK2R is Gq coupled, therefore, prokineticin 2 (PROK2) which is an agonist for the PROK2R was used to
10 see whether CreLuc mice would be responsive to Gq modulation. Gq bypasses cAMP but utilizes the PLC pathway to activate CREB (see pathway schematic of Figure 1). Primary cortical neurons were harvested from line 187 (which was shown to have inducible luciferase in brain and spinal cord- see text below and Figure 24) on E18. The assay was run on day 3 in culture for 4 hours, or 8 hours. The PROK2
15 peptide (Peprotech, Cat#100-46) is added as an aqueous solution at 1 nM and 100nM. Bright Glo was used for the luciferase signal, and the assay was run on a TopCount luminometer. Highly significant increases in luciferase levels vs. the media only control are observed at both timepoints, for both concentrations of PROK2. Thus, the CreLuc mice can be used to screen for compounds that modulate Gq
20 coupled receptors.

The effects of PROK2 peptide on luciferase expression in primary cortical neurons from different CreLuc lines was examined. Primary cortical neurons were harvested from four different CreLuc lines at E18. Lines were selected based on earlier
25 experiments which determined inducible luciferase levels in whole brain extracts (data not shown). The assays were run in triplicate at day three in culture with either 1nM or 100nM PROK2 peptide at two timepoints, 4 hours and 24 hours. BrightGlo was used for the assay and read on a TopCount. Statistically significant increases in luciferase levels were observed at both timepoints and at both concentrations for
30 lines 219 and 175. Even though all lines responded to the PROK2 ligand to some degree, differences in luciferase expression levels may be due to differences in transgene integration.

c) Gi Modulation via mGluR7 receptors

Primary cell cultures from the CreLuc mice may be used to screen compounds that are capable of modulating Gi coupled receptors. Figure 13A shows the effects of the mGluR7 agonist, AMN082 on luciferase expression in primary cortical neurons.

Cortical neurons were harvested from E18 embryos (line 187). The assay was run at day 3 in culture. Forskolin was used at 10 μ M to give a larger signal window to be used to detect Gi activity. The mGluR7 agonist, AMN082 was used in combination with forskolin at 1nM, 10nM, 100nM and 1 μ M. The assay was read on a TopCount with Bright Glo (Promega) at 4 hours, and 8 hours. Significant decreases in luciferase levels, versus forskolin only, are observed at 100nM and 1 μ M AMN082 for both timepoints. The CreLuc mice can be used to screen for compounds that modulate Gq coupled receptors. For example, "unknown" compounds designated "A", "B" and "C" were screened in primary cortical neuronal cultures using AMN082 as an internal control to determine if any compounds exhibited Gi modulating ability (Figure 13B). Compound A would be considered an active compound with a calculated EC₅₀ = 1.529e-06 compared to the known mGluR7 agonist AMN082 which had an EC₅₀ = 1.735e-07 while compound B would be considered a non-specific compound and compound C is inactive.

Screening compounds for the ability to modulate Gi coupled receptors using primary cell cultures from the CreLuc mice may be done in the presence of a cAMP stimulator such as forskolin as demonstrated above. However, to increase the signal window even wider, rolipram (blocks the breakdown of cAMP) was used in combination with forskolin. Primary cortical neurons were isolated from E18 embryos from line 187 which was previously determined to have inducible brain and spinal cord expression. The assay, in triplicates, was run at day 7 in culture for 6 hours. A dose curve, 10 μ M down to 1nM, for AMN082 was run in combination with 50 μ M forskolin and 10 μ M rolipram. The assay was read on a TopCount luminometer with Bright Glo substrate (Promega). Statistically significant decreases, $p < 0.0005$, in luciferase levels (59-61% vs. forskolin and rolipram controls) are seen at 1 μ M, 100nM and 1nM.

Gi modulation of luciferase expression in primary cortical neurons from different CreLuc lines by the CB1 agonist, CP 55,940 was performed in the presence of forskolin and rolipram (Figure 15). Primary cortical neurons were harvested from lines 69, 187, 175 and 219 at E18. All lines used were previously determined to have inducible luciferase expression levels in whole brain tissue extracts (data not shown). The assays were run on day three in culture. The CB1 agonist was used at 10 μ M, forskolin at 5 μ M and rolipram at 10 μ M. Two timepoints were run, four hours and twenty-four hours. Bright Glo luciferase assay substrate was then added, and the assay read on a Topcount luminometer. The assay was run in triplicate. Data shown is the average of the triplicates. Data is shown as counts per second (cps). Decreased signals in luciferase are observed in all four lines at both timepoints with the addition of the CB1 agonist. Slight differences in response levels may be due to transgene integration, but all lines were responsive to Gi modulation and would be amendable to screening compounds for Gi activity.

Gi modulation of luciferase in primary cortical neurons from CreLuc mice by the CB1 agonist, CP 55,940 was in a dose-dependent manner (Figure 15B). Cortical neurons were isolated from E18 embryos. The assay was run on day 3 in culture. Forskolin and rolipram were used at 10 μ M. The agonist was added at concentrations of 10 μ M, 1 μ M and 100nM. The assay was read on a TopCount with BrightGlo (Promega) at 8 hours. Significant decreases in luciferase levels, (agonist plus forskolin and rolipram versus forskolin and rolipram alone) are observed at all 3 concentrations.

2. Primary Striatal Neurons

Any cell culture derived from the CreLuc mice may be used to screen compounds for the ability to modulate GPCRs. Luciferase expression was induced in CreLuc striatal neurons by forskolin and rolipram, and Gs agonists for DRD1 and AD β R (Figure 16). Striatum neurons were isolated from E14 embryos (line 187). The assays were run at day 4 in culture. Forskolin was used at 5 μ M, rolipram at 10 μ M. The Gs agonists (isoproterenol is an agonist at AD β R, dopamine and SKF82958 are agonists at D1DR) were used at 10 μ M, 3 μ M and 1 μ M. The assay was read at 5 hours with a

TopCount luminometer and Bright Glo luciferase reagent (Promega). A highly significant increase of 28 fold is observed in cells treated with the forskolin and rolipram combination. Significant increases are also observed with 10 μ M dopamine, 2.7 fold, and all 3 concentrations of isoproterenol. Thus, the transgene is functional in striatal as well cortical neurons (see Figure 10).

3. Whole Splenocytes

Whole splenocyte preps isolated from CreLuc mice (line 64) were used to show the effects of general cAMP inducers such as forskolin and rolipram, as well as Gs agonists on luciferase expression (Figure 16). Spleens were harvested from the animals in 1X HBSS (Invitrogen, Carlsbad, CA, cat#14025). The cells were then isolated by mechanically disrupting the spleen capsule using the end of a 5ml syringe in 5mls of D-PBS (Invitrogen, Carlsbad, CA, cat#14190). The cell suspension was then passed through a 70 μ m strainer into a 50ml conical tube. The cells were then spun down at 800rpm, then resuspended and incubated for 6 minutes at room temperature in 5ml of 1X Pharm Lyse solution (BD BioSciences, cat#555899). Cells were then washed by adding 28 mls of media consisting of RPMI 1640 (Invitrogen, Carlsbad, CA, cat#11875), 10% FCS (Invitrogen, Carlsbad, CA, cat#16000), 1% pen/strep (Invitrogen, Carlsbad, CA, cat#15070) and 0.1% β -mercaptoethanol (Invitrogen, Carlsbad, CA, cat#21985). Cells were diluted to 2x10⁶/ml, and 100 μ l of the cells were plated per well in a 96 well white opaque plate. Half of the cells were stimulated for 24 hours with anti-CD3 antibody (BD Pharmingen, cat#553058), the other half were untreated. At 24 hours, compounds were added to the plates for an additional 4 hours. The general inducers include rolipram (Sigma R6520), forskolin (Sigma F6886). The Gs agonists used are: EX00000173A (173A; in house synthesized) is an agonist at the Gs coupled prostaglandin E2 (EP2) receptor; BW245C (Sigma B9305) is an agonist at the Gs coupled prostaglandin D2 receptor (DP1) and isoproterenol (Sigma I5627) is an agonist at the Gs coupled β -adrenergic receptor (AD β R). The assay was run in triplicate. After 4 hours, 100 μ l of BrightGlo (Promega, Madison, WI, cat#E2610) and the assay was read on a TopCount luminometer. An increase of 14 fold versus DMSO is observed in CD3 stimulated

cells in the presence of rolipram and forskolin. Statistically significant increases (by t-test) are also observed with all three of the Gs agonists versus the DMSO only control. This experiment shows that the transgene is functional in whole splenocyte populations. The preparation used in the experiment was a whole splenocyte preparation which is a mixed population of cells. Experiments described below look at luciferase expression levels in subpopulations such as T cells and B cells.

4. Isolated T Cells

The effects of general cAMP activation by rolipram and forskolin in T cells isolated from five different sublines of CreLuc mice was studied (Figure 18). Whole splenocyte populations were prepared by mechanical disruption as described above. CD4⁺ cells were then isolated using MACS Magnetic Separation with Positive Selection Column (Miltenyi Biotec cat#130-049-201). The cells, 1.5×10^5 per well, were then plated on 96 well white opaque plates and then stimulated with anti CD3 antibodies (BD Pharmingen cat#553058). After 18 hours, 10 μ M rolipram (Sigma R6520) and 5 μ M forskolin (Sigma F6886) were added to the plates for an additional 4 hours. BrightGlo (Promega, Madison, WI, cat#E2610) was added and the assay was read on the TopCount. Data is shown as luminescence (counts per second), and as fold increase over media only controls. Increases in expression of luciferase were observed in all lines tested with line 64 giving the highest levels of induction demonstrating that the cAMP pathway activated by general modulators in the CreLuc mice.

The effects of different Gs agonists on luciferase levels in anti CD3 stimulated CD4⁺ T cells isolated from CreLuc mice (line 64) was studied (Figure 19). Whole splenocyte populations were prepared by mechanical disruption. CD4⁺ cells were then isolated by MACS Magnetic Separation with Positive Selection (Miltenyi Biotec cat#130-049-201). The cells, 1.5×10^5 per well, were then plated on 96 well white opaque plates and then stimulated with anti CD3 antibodies (BD Pharmingen cat#553058). After 24 hours, compounds were added for an additional 4 hours. Gs agonists for DP (BW245C), EP2 (EX00000173A) and AD β R (isoproterenol) were all used at 10 μ M. Forskolin 5 μ M, and rolipram 10 μ M. BrightGlo (Promega, Madison, WI,

cat#E2610) was added and the assay was read on the TopCount. A 25-fold increase is observed in cells treated with forskolin and rolipram. Highly significant increases are observed with all three Gs agonists, 3-fold for BW245C, 2-fold for 173A, and 5 fold for isoproterenol. The transgene was responsive to specific Gs agonists as well as general modulators of the cAMP pathway.

5. Isolated B Cells

The effects of general cAMP activation by rolipram and forskolin in B cells isolated from two different sublines of CreLuc mice was examined (Figure 20). Whole splenocyte populations were prepared by mechanical disruption as described above. B220+ cells were then isolated using MACS Magnetic Separation with Positive Selection Column (Miltenyi Biotec, cat# 130-049-501). The cells were then plated, 2.0×10^5 per well on 96 well white opaque plates and then stimulated with 10ng/ml lipopolysaccharide (LPS) (Sigma L-2630). After 18 hours, 10 μ M rolipram (Sigma R6520) and 5 μ M forskolin (Sigma F6886) were added to the plates for an additional 4 hours. BrightGlo (Promega, Madison, WI, cat#E2610) was added and the assay was read on the TopCount. Data is shown as luminescence (counts per second), and as fold increase over media only controls. Increases in luciferase expression are observed in line 64 but not line 229.

6. Microglia

Induction of luciferase expression in microglia isolated from CreLuc mice by the general cAMP activators, forskolin and rolipram, and an agonist for the DP receptor, BW245C was examined (Figure 22). Primary microglia were isolated from the cortices from P2 mice (line 64) in media consisting of DMEM (GIBCO Cat# 11995), 10% FBS (GIBCO Cat# 16140) and 1% Penicillin-Streptomycin 100X (GIBCO Cat# 15140). The cells were plated in 96 well format on Poly-D-Lysine-coated plates. Cells were either left untreated or stimulated for 2 hours with 100 ng/ml LPS. Compounds were then added for an additional 4 hours before the Bright Glo assay was run. The compounds used were 5 μ M forskolin, 10 μ M rolipram or the combination of the two, or the Gs agonist for the DP1 receptor, BW245C at 10 μ M. In unstimulated conditions,

microglia are unresponsive to general cAMP modulators and the specific Gs agonist for the DP receptor. However, when microglia were stimulated with LPS, the cells became responsive to forskolin and BW245C.

7. Mouse Embryonic Fibroblasts

The effects of forskolin, rolipram and isoproterenol on luciferase expression in mouse embryonic fibroblasts was investigated (Figure 33). Mouse embryonic fibroblasts were cultured from E12 embryos from six independent CreLuc lines. Cells were plated at 20,000 cells per well. Compounds tested include 10 μ M forskolin, 5 μ M rolipram and 10 μ M isoproterenol (AD β R agonist). Significant increases were observed in all lines in response to the combination of forskolin and rolipram. Significant increases were also observed in three of the lines in response to isoproterenol.

7. Cardiomyocytes

The effects of forskolin and rolipram and isoproterenol on luciferase levels in cardiomyocytes was studied (Figure 35). Cardiomyocytes were isolated from P3 pups from line 229. The cells were cultured in a 96 well plate. Compounds tested include 10 μ M forskolin, 5 μ M rolipram and 10 μ M isoproterenol (AD β R agonist). Significant increase in luciferase levels were observed with the combination of rolipram and forskolin (25 fold) and a significant increase was also observed with the agonist isoproterenol (2 fold).

B. In Vivo and Ex Vivo Experiments

1. Effects of General cAMP Modulators

The effects of intrathecally injected forskolin and rolipram on the induction of luciferase expression was studied in the brain and spinal cord of line 187 CreLuc mice (Figure 23). The line of mice chosen was previously determined to have

expression of the transgene in both brain and spinal cord (data not shown). The line of mice selected for this assay (line187) has inducible levels of luciferase in both the brain and spinal cord. N=3-4 mice per group, four treatment groups, 3 month old males. Group A: DMSO control, Group B: 1µg forskolin/ 10µg rolipram, Group C: 10µg forskolin/ 10µg rolipram, Group D: 40µg forskolin/ 10µg rolipram. The animals were dosed via intrathecal injection, lumbar region, and volume of 5µl per mouse. They were imaged at 4 hours post dosing. The data for both spinal cord and brain is shown as the average peak radiance, photons per second per cm². Statistically significant increase is observed in the spinal cord, a significant increase in luciferase signal in brain is only seen in the highest concentration of forskolin. When general modulators or cAMP are injected into the spinal cord, there is a localized increased expression of the transgene in the spinal cord with a lower response observed in the brain.

2. Effects of Gs Agonists

The effects of the EP2 agonist, EX00000173A on luciferase expression in the brain and spinal cord of CreLuc mice was studied (Figure 24). Male mice, 5 months old, from line 187, n=5, were injected i.p. with either vehicle (5% DMSO, 0.05% tween 80, PBS) or 10 mg/kg EX00000173A (in house synthesis). Animals were bioimaged at 4 hours post doing. Data shown is the average of the 5 mice for both brain and spinal cord, as photons per second per cm². Statistically significant increases are seen with agonist treatment in both the brain and spinal cord. Thus, a Gs agonist administered i.p. activated the transgene.

Further, the effects of the EP2 agonist, EX00000173A on luciferase expression in CreLuc mice was dose-responsive (Figure 25). Line 187 was selected because this line has high inducible expression in both brain and spinal cord. The assay consists of five groups of six week old mice, with an n=5. Mice were dosed with either vehicle control (D-PBS; Dulbeccos phosphate buffered saline, Invitrogen, cat#14040) or varying doses of the EP2 agonist EX00000173A (in house synthesis) 1.5nmol, 5nmol, 15nmol and 50 nmol. Mice were dosed by intrathecal injection (5µl per mouse) and were bioimaged 4 hours later on the IVIS bioimager. Data for both the brain and

spinal cord is shown as the mean of the five mice, average peak radiance, photons per second per cm². Statistically significant (t-test) increases with agonist treatment are observed at three concentrations in the spinal cord, but not in the brain. An increase (not significant) in the brain is only observed at the highest concentration as expected with an intrathecal injection.

Figure 26 shows induction of luciferase in different tissues by the adrenoceptor beta3 (Adrb3) agonist, CL316,243 (1 mg/kg, ip) in CRE-Luc mice. The luciferase assay was performed in tissue homogenates. Among 12 independent transgenic lines screened, 7 lines showed greater than 10 fold of induction in the adipose tissue and lung. 4 lines showed visible induction via bioimaging (see Figure 27A and 27B) .

3. Specific Effects of the Gs Agonist AVE0010

Transgene activation co-localizes with tissue specific receptor activation as exemplified by the induction of luciferase reporter by the glucagon-like peptide 1 receptor (GLP-1R) agonist, AVE0010, which was studied in three independent lines of CRE-Luc mice (Figure 28). The GLP-1R is Gs coupled. Baseline images were acquired on day 1. On day 2, mice were treated with AVE0010 (0.1 mg/kg, sc) and imaged after 4 hours. Folds of induction over baseline were indicated at the bottom. As expected since the GLP-1R receptor is mainly expressed in pancreatic tissue, there was a strong induction of transgene activation observed in the pancreas. Luciferase activities in 8 different tissues were measured for lines 11, 16 and 90 (Figure 29) following treatment with AVE0010 (0.1 mg/kg, s.c.) for 4 hours. The luciferase activity in tissue homogenates confirmed the pancreas-specific induction of luciferase. The activity of AVE0010 is limited to the pancreas even though GLP-1R is found in different tissue types.

The induction of Cre-Luc by AVE0010 is likely beta-cell mediated. The effects of the beta-cell toxin streptozotocin (STZ) on the induction of CRE-Luc by AVE0010 was investigated (Figure 30). Male CRE-luc mice (line 11) were imaged before and after AVE0010 (0.1 mg/kg, sc) treatment. The data indicated that all mice were responsive to AVE0010 (see middle panel of Figure 30). Then, the animals were treated with

vehicle or STZ (200 mpk, ip). Four days later, they were imaged again after AVE0010 treatment. Compared to vehicle group, the STZ group had decreased luciferase induction.

5 The induction of CRE-Luc by AVE0010 is likely beta-cell-specific (Figure 31). Animals were treated as described in Figure 30. Blood glucose levels were measured by tail vein nicking on unfasted mice. Glucose levels were read on a Bayer glucometer. Glucose levels are shown as mg glucose/ml. Fold induction is the luciferase bioimaging levels of AVE10 dosing versus the baseline signals. Blood glucose levels
10 (BG) were increased by STZ (upper left panel). Non-fasting BG levels were reduced by AVE0010 (0.1 mg/kg, sc). BLI data shown in Figure 30 were quantified.

4. Bone Marrow Engraftments

15 CreLuc bone marrow engraftments were performed using NOD scid gamma (NSG) mice (Figure 32). NSD mice are immunocompromised mice lacking mature T and B cells, functional natural killer cells and are deficient in cytokine signaling allowing for engraftment of hematopoietic cells. Bone marrow cells were harvested from lines 44
20 heterozygotes (having high basal luciferase levels; data not shown) and line 64 homozygotes (having inducible luciferase levels). The cells were then engrafted via tail vein injections of cells into irradiated NSG mice at 1 million or 5 million. For line 44, mouse 1 and 2 received 5 million cells while mouse 3 and 4 received 1 million cells. For line 64, mouse 1 received 5 million cells, mouse 2, 3, and 4 received 1
25 million cells per mouse. (4 NSG mice per CreLuc line). The animals were bioimaged at 4 weeks (data not shown) and then again at 8 weeks. Prior to imaging, the line 64 mice were induced for 5 hours with 5mg/kg forskolin and 10mg/kg rolipram. Bioimaging pictures are shown for the 8 week timepoint. The luciferase levels in the NSG mice engrafted with line 44 bone marrow cells mimics the image seen with the
30 CreLuc line 44 (data not shown), with expression observed in the joints, spinal cord, head, and breastbone. Inducible luciferase expression is observed in the spleens of the NSG mice engrafted with line 64 bone marrow cells.

4. Animal Models

The CreLuc mice can be used to study animal models of disease or aspects of disease states. Furthermore, the CreLuc mice can be used to screen for compounds that are able to modulate the disease or aspects of the disease that was induced in the CreLuc mice. For example, the effects of zymosan treatment on luciferase levels in CreLuc line 187 was studied (Figure 34). CreLuc mice (line 187) in the treated group were injected s.c in both rear paws with zymosan to induce a pain response. The animals were then bioimaged daily for 4 days. Statistically significant Increases in luciferase expression are observed in the paws of the animals, in response to zymosan at all time points. Zymosan is a yeast cell wall component that strongly activates an inflammatory response. Thus, the CreLuc mice are a tool in which the inflammatory response can be monitored over time in the same animal. Also, the zymosan treated animals can be used a screening tool to assess the abilities of test compounds to reverse or exacerbate the inflammatory induced by zymosan.

CLAIMS:

1. A transgenic non-human animal having a genome comprising a transgene
5 comprising a first insulator element, a response element, a promoter, a
bioluminescent reporter, a functional element and a second insular element.
2. The transgenic non-human animal of claim 1, wherein the first insulator element is
10 selected from the group consisting of matrix attachment regions (MAR), DNase
I-hypersensitive site (HS4) and inverted terminal repeats (ITR).
3. The transgenic non-human animal of claim 1, wherein the second insulator
15 element is selected from the group consisting of matrix attachment element
(MAR), HS4 and ITR.
4. The transgenic non-human animal of claim 1, wherein the first insulator element is
the same as the second insulator element.
5. The transgenic non-human animal of claim 1, wherein the response element is
20 selected from the group consisting of cAMP response element (CRE), activator
protein 1 (ASP1), glucocorticoid response element (GRE), heat shock response
element (HSE), serum response element (SRE), thyroid response element
(TRE) and estrogen response element (ERE).
- 25 6. The transgenic non-human animal of claim 5, wherein the response element is
repeated in tandem two to twenty-four times.
7. The transgenic non-human animal of claim 6, wherein the response element is
30 repeated in tandem six times.
8. The transgenic non-human animal of claim 7, wherein the response element is
CRE.

9. The transgenic non-human animal of claim 1, wherein the promoter is herpes simplex virus thymidine kinase minimal (HSV TK min).

10. The transgenic non-human animal of claim 1, wherein the bioluminescent reporter
5 is selected from the group consisting of luciferase, chloramphenicol acetyltransferase (CAT), beta-galactosidase, secreted alkaline phosphatase (SEAP), human growth hormone (HGH) and green fluorescent protein (GFP).

11. The transgenic non-human animal of claim 1, wherein the functional element is
10 human growth hormone (hGH) gene.

12. The transgenic non-human animal of claim 1, wherein the transgene comprises
SEQ ID NO: 18.

13. The transgenic non-human animal of claim 11, wherein the transgene comprises
15 SEQ ID NO: 19.

14. A cell isolated from the transgenic non-human animal of claim 1.

15. A tissue slice isolated from the transgenic non-human animal of claim 1.
20

16. A method of identifying a G protein-coupled receptor (GPCR) ligand, the method comprising:

(a) measuring an amount of bioluminescence in the transgenic non-human
25 animal of claim 1;

(b) administering a test agent to the transgenic non-human animal;

(c) measuring an amount of bioluminescence of the transgenic non-human
animal at one or more time points following administration of the test
agent; and

(d) comparing the amount of bioluminescence measured in (a) to the amount
30 of bioluminescence measured in (c)

wherein a difference in the amount of bioluminescence in (a) compared to (c)
identifies the test agent as a GPCR ligand.

17. A method of identifying a G protein-coupled receptor (GPCR) ligand, the method comprising:

- (a) preparing a tissue slice from the transgenic non-human animal of claim 1;
 - 5 (b) measuring an amount of bioluminescence in the tissue slice;
 - (c) administering a test agent to the tissue slice;
 - (d) measuring an amount of bioluminescence of the tissue slice at one or more time points following administration of the test agent; and
 - (e) comparing the amount of bioluminescence measured in (b) to the amount
 - 10 of bioluminescence measured in (d)
- wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the test agent as a GPCR ligand.

18. A method of identifying a G protein-coupled receptor (GPCR) ligand, the method comprising:

- (a) preparing a cell isolated from the transgenic non-human animal of claim 1;
 - (b) measuring an amount of bioluminescence in the cell;
 - (c) administering a test agent to the cell;
 - (d) measuring an amount of bioluminescence in the cell at one or more time
 - 20 points following administration of the test agent; and
 - (e) comparing the amount of bioluminescence measured in (b) to the amount of bioluminescence measured in (d)
- wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the test agent as a GPCR ligand.

19. A method of monitoring GPCR function in a non-human animal, the method comprising:

- (a) transgenically modifying a non-human animal to express a transgene comprising a first insulator element, a response element, a promoter, a
- 30 bioluminescent reporter, a functional element and a second insular element;
- (b) monitoring bioluminescence from the non-human animal; and
- (c) correlating said bioluminescence to GPCR function.

20. A method of making a non-human transgenic animal for use in monitoring GPCR function, the method comprising:

5 (a) transgenically modifying a non-human animal to express a transgene comprising a first insulator element, a response element, a promoter, a bioluminescent reporter, a functional element and a second insular element;

(b) measuring an amount of bioluminescence in the transgenic non-human animal of (a);

10 (c) administering a GPCR ligand to the transgenic non-human animal;

(d) measuring an amount of bioluminescence of the transgenic non-human animal at one or more time points following administration of the GPCR ligand; and

15 (e) comparing the amount of bioluminescence measured in (b) to the amount of bioluminescence measured in (d)

wherein a difference in the amount of bioluminescence in (b) compared to (d) identifies the non-human transgenic animal for use in monitoring GPCR function.