

**PCT**WORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau

INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁶: A61L 2/18, G02C 13/00	A1	(11) International Publication Number: WO 97/26923 (43) International Publication Date: 31 July 1997 (31.07.97)
(21) International Application Number: PCT/US96/18837 (22) International Filing Date: 6 December 1996 (06.12.96) (30) Priority Data: 60/010,402 22 January 1996 (22.01.96) US (71) Applicant: BAUSCH & LOMB INCORPORATED [US/US]; One Bausch & Lomb Place, Rochester, NY 14604-2701 (US). (72) Inventors: ROGALSKYJ, Jill, S.; 96 Main Street, Livonia, NY 14487 (US). HEILER, David, J.; 173 Wadsworth Avenue, Avon, NY 14414 (US). DOBIE, Alyce, K.; 6932 Tuckahoe Road, Williamson, NY 14589 (US). (74) Agents: BLACK, Edward, W. et al.; Bausch & Lomb Incorporated, One Bausch & Lomb Place, Rochester, NY 14604-2701 (US).		(81) Designated States: AL, AM, AT, AU, AZ, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, IS, JP, KE, KG, KP, KR, KZ, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, TJ, TM, TR, TT, UA, UG, UZ, VN, ARIPO patent (KE, LS, MW, SD, SZ, UG), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG). Published <i>With international search report.</i> <i>Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i>
(54) Title: DUAL NEUTRALIZATION SYSTEM FOR IODINE TREATMENT OF CONTACT LENSES (57) Abstract <p>A method and system for disinfecting contact lenses with iodine using a neutralizing solution having a first, rapid neutralizing component which completes its reaction with available iodine in less than five seconds and a second, slower neutralizing component which completes its reaction with available iodine in about 3 to 30 minutes.</p>		