

A method for detection of human papillomavirus (HPV) type

FIELD OF THE INVENTION

The present invention relates to detection, typing and quantitation of human papillomaviruses (HPV) in clinical samples. More specifically, it relates to
5 oligonucleotide sequences which can be used in amplification reactions and used for the simultaneously detection, typing and quantification of a plurality of HPV genomes and a diagnostic kit thereof.

BACKGROUND OF THE INVENTION

10 Human papillomavirus (HPV) infection was reported as a cause of cervical cancer in the 1980's. The relationship between cancer malignancy and HPV genotype was considered to be of interest. More than 120 HPV types have since then been identified. Thirteen of the HPV types are recognised as being oncogenic in mucosal epithelia. These types include HPV 16, 18, 31, 33, 35, 39, 45, 51, 52, 56,
15 58, 59 and 66.

Human papillomavirus infections are confined to the epithelium, which limits their contact with the immune apparatus and they are poorly immunogenic. Thus serology is not a useful test for HPV. Similarly, culture is not a useful diagnostic technique. Although it is possible to reproduce the virus life cycle in vitro, this is a
20 highly sophisticated tissue culture technique requiring in vitro differentiating epithelial cell cultures. Furthermore, high grade cervical neoplasias and cancers do not produce virus. Therefore, HPV testing relies on the detection of viral nucleic acids.

HPV testing has been shown to have a higher sensitivity than the PAP test
25 (Papanicolaou test) for high-grade cervical dysplasia, although its specificity is lower due to the high prevalence of benign, transient HPV infections.

Many countries have adopted HPV testing as a secondary screening (triage) test for low-grade and equivocal dysplasia and its use as in primary screening has also been suggested.

The dominant HPV test is at present the Digene HCII test, whose utility is supported by a large body of scientific literature. In this test HPV DNA in the sample is hybridised to complementary RNA and the resulting DNA-RNA hybrids are captured by antibodies on a solid surface where they subsequently bind
5 antibody-enzyme conjugates that are detected by enzyme-catalysed chemiluminescence. The LCII test does not distinguish between the different oncogenic HPV types.

It is also possible to detect HPV in fixed cells using in-situ hybridisation with probes labelled either with fluorescent moieties (fluorescence in-situ hybridisation,
10 FISH) or with haptens that can bind to an enzyme conjugate which generates an insoluble coloured precipitate in situ. In situ hybridisation can provide valuable information on the location and physical state of HPV, but type specific detection of multiples of HPV types in the same sample is impracticable.

Many screening authorities express a preference for genotyping tests in secondary
15 screening since this allows serial transient infections which are benign, to be distinguished from type-specific persistent infections which are not benign. In addition, some HPV types, particularly HPV16, are markedly more oncovirulent.

Finally, the use of a genotyping test provides useful epidemiological information which will be of particular value in monitoring the population effects of vaccination
20 against HPV16 and HPV18.

The available HPV genotyping tests are mostly based on PCR followed by reverse hybridisation in some format. The PCR step targets the L1 or the E1 genes. PCR amplicons are labelled under amplification and hybridised to immobilised probes on a microarray, a membrane strip or in the wells of a microtiter plate. The bound
25 labelled amplicon is then detected, either directly by fluorimetry if the amplicon is labelled with a fluorescent moiety or using enzyme conjugates that bind to the label and which are detected using colorimetric or luminometric techniques. Such tests involve handling of PCR amplicons in the open laboratory which necessitates stringent hygienic precautions. This, combined with the high price of such tests
30 places the available tests beyond the reach of most mass-screening programs.

A genotyping test that detects mRNA from the viral oncogenes E6 and E7 using SDA (sequence dependent amplification) is also available. However, this method detects only five of the thirteen oncogenic types and only when there is active oncogene expression.

5 WO 2007115582 A1 describes sets of PCR primers and probes for detection of multiple HPV types by real-time PCR. These sets involves a complex combination of primers. In the method according to present invention, only one forward primer and four reverse primers are required. This simplifies quality control and manufacture and reduces cost. Further, the genes that are targeted in WO
10 2007115582 A1 are subject to intratypic sequence variation. This variation is not taken into account. Thus false negative results may occur in clinical samples that contain HPV sequence variants other than the prototype sequences. According to the present invention known sequence variation within the target region of the L1 gene has been taken into account prior to primer and probe design and all variant
15 positions are either avoided or accounted for by modification of primer sequences.

Further, in WO 2007115582 A1 it is stated that sensitivities down to 100 copies for all targeted HPV types are achieved. However, the experiments described do not include the use of carrier DNA. This gives an unrealistically optimistic estimate of
20 sensitivity as the human DNA present in clinical samples suppresses PCR sensitivity. The sensitivity of the method according to the present invention is measured in the presence of 500ng of human DNA, which correspond to $7,5 \times 10^4$ human cells and thus accurately mimics the conditions in a cell-rich clinical sample. Thus the present invention provides enhanced sensitivity under realistic
25 conditions.

Further, WO 2007115582 A1 describes the use of a single tube 13-plex multiplex PCR test. However, according to the state of the art there does not seem to be an instrument capable of performing such a test. The method according to the
30 present invention may be operated on existing obtainable instruments in the format proposed or a variant or an analog thereof.

In WO 2009011472 A1 the target gene which has been used for the PCR tests is E1. However, this gene is a hotspot for deletions and insertions during cancer development. The probability of a false negative result increases thus with the severity of disease. The method according to the present invention targets the L1
5 gene, which is much less prone to deletion.

Further, in WO 2009011472 A1 the methods for detection and typing of HPV involve handling of PCR products in the open laboratory. This involves the risk of carry-over contamination. The method according to the present invention is
10 performed in a closed-system which eliminates this risk, thus ensuring more reliable results and eliminating the extra labour required for decontamination of equipment and facilities after open-laboratory post-PCR work.

Based on the above, there is thus a need in the field for new methods, primers, probes and kits for detection and typing of HPV which are highly sensitive,
15 reproducible and less expensive than the existing tests and which thus are appropriate for mass-screening programs.

SUMMARY OF THE INVENTION

The present invention describes the development of a quantitative realtime PCR
20 method using consensus primers and multiplexed type-specific TaqMan probes that allows sensitive and specific detection and typing of the thirteen oncogenic HPV types HPV16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 and 66. Further, the method according to the present invention detects HPV types HPV6,11, 26, 40, 42, 43, 44, 53, 54, 61, 68, 70, 72, 73, 82.

25 More specifically, the present invention relates to a method for detecting and/or typing and/or quantitation of the presence of a human papillomavirus (HPV) type in a biological sample, the test comprising the following steps:

- (i) amplification of a nucleic acid fragment of HPV in the presence of a DNA polymerase, forward PCR primer with an oligonucleotide
30 sequence of SEQ ID NO:1 or an analog thereof, reverse PCR primer with an oligonucleotide sequence selected from the group consisting

- of: SEQ ID NO:2, SEQ ID NO:3, SEQ ID NO:4, SEQ ID NO:5 or
 SEQ ID NO:6 or an analog thereof, a probe with a oligonucleotide
 sequence selected from the group consisting of: SEQ ID NO:7, SEQ
 ID NO:8, SEQ ID NO:9, SEQ ID NO:10, SEQ ID NO:11, SEQ ID
 NO:12, SEQ ID NO:13, SEQ ID NO:14, SEQ ID NO:15, SEQ ID
 NO:16, SEQ ID NO:17, SEQ ID NO:18, SEQ ID NO:19, SEQ ID
 NO:20, SEQ ID NO:21, SEQ ID NO:22, SEQ ID NO:23, SEQ ID
 NO:25, SEQ ID NO:26, SEQ ID NO:27, SEQ ID NO:30, SEQ ID
 NO:31, SEQ ID NO:32, SEQ ID NO:34, SEQ ID NO:35, SEQ ID
 NO:37, SEQ ID NO:38, SEQ ID NO:39, SEQ ID NO:41, SEQ ID
 NO:44, SEQ ID NO:45, SEQ ID NO:46, SEQ ID NO:47, SEQ ID
 NO:48, SEQ ID NO:49, SEQ ID NO:50, SEQ ID NO:51, SEQ ID
 NO:52, SEQ ID NO:55, SEQ ID NO:56, SEQ ID NO:57, SEQ ID
 NO:58, SEQ ID NO:59, SEQ ID NO:61, SEQ ID NO:62, SEQ ID
 NO:63, SEQ ID NO:64, SEQ ID NO:65, SEQ ID NO:66, SEQ ID
 NO:68, SEQ ID NO:70, SEQ ID NO:71, SEQ ID NO:72, SEQ ID
 NO:73, SEQ ID NO:74, SEQ ID NO:75, SEQ ID NO:77, SEQ ID
 NO:78, SEQ ID NO:79 and SEQ ID NO:80 or an analog thereof,
 labelled with a fluorescent dye and a quencher molecule, and an
 internal control;
- (ii) detection of a change of fluorescence and determining HPV type.

Further, the method according to the present invention refers to detecting and/or
 typing and/or quantitation simultaneously the presence of HPV 18, 52, 59, 39, 51,
 56, 66, 16, 45, 58, 31, 33 and 35 in a biological sample, the method comprising
 the following four parallel reactions:

Reaction 1- a) amplification of the nucleic acid fragments of HPV 18, 52
 and 59 in the presence of

- (i) a nucleic acid polymerase
- (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID
 NO: 1 or an analog thereof,

- (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 2 or an analog thereof,
- (iv) a probe with an oligonucleotide sequence of SEQ ID NO: 7,8 and 9 or an analog thereof, labelled with a fluorophore and a quencher molecule and

b) detection of a change of fluorescence the wavelength thereof determining the HPV type 18, 52 and 59 ;

Reaction 2-a) amplification of the nucleic acid fragments of HPV 39,51,56 and 66 in the presence of

- (i) a nucleic acid polymerase
- (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
- (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 3 or an analog thereof,
- (iv) a probe with an oligonucleotide sequence of SEQ ID NO: 10, 11, 12 and 13 or an analog thereof, labelled with a fluorophore and a quencher molecule and

b) detection of a change of fluorescence the wavelength thereof determining the HPV type 39, 51, 56 and 66;

Reaction 3- a) amplification of the nucleic acid fragments of HPV 16,45 and 58 in the presence of

- (i) a nucleic acid polymerase
- (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
- (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 4 or an analog thereof,

(iv) a probe with an oligonucleotide sequence of SEQ ID NO: 14, 15 and 16 or an analog thereof, labelled with a fluorophore and a quencher molecule and

5 b) detection of a change of fluorescence the wavelength thereof determining the HPV type 16, 45 and 58;

Reaction 4- a) amplification of the nucleic acid fragments of HPV 31, 33,35 and internal control in the presence of

- 10 (i) a nucleic acid polymerase
- (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
- (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 5 or an analog thereof,
- 15 (iv) a probe with an oligonucleotide sequence of SEQ ID NO: 17, 18, 19 and 20 or an analog thereof, labelled with a fluorophore and a quencher molecule and
- (v) internal control and

20 b) detection of a change of fluorescence the wavelength thereof determining the HPV type 31, 33, 35 and internal control.

The present invention also refers to a primer for HPV amplification selected from the group consisting of: SEQ ID NO: 1, SEQ ID NO:2, SEQ ID NO: 3, SEQ ID NO:4, SEQ ID NO: 5 and SEQ ID NO:6.

25

Further, the invention refers to a probe for HPV detection and/or typing and/or quantitation selected from the group consisting of: SEQ ID NO:7, SEQ ID NO:8, SEQ ID NO:9, SEQ ID NO:10, SEQ ID NO:11, SEQ ID NO:12, SEQ ID NO:13, SEQ ID NO:14, SEQ ID NO:15, SEQ ID NO:16, SEQ ID NO:17, SEQ ID NO:18, SEQ ID NO:19, SEQ ID NO:20, SEQ ID NO:21, SEQ ID NO:22, SEQ ID NO:23,

30

SEQ ID NO:25, SEQ ID NO:26, SEQ ID NO:27, SEQ ID NO:30, SEQ ID NO:31,
SEQ ID NO:32, SEQ ID NO:34, SEQ ID NO:35, SEQ ID NO:37, SEQ ID NO:38,
SEQ ID NO:39, SEQ ID NO:41, SEQ ID NO:44, SEQ ID NO:45, SEQ ID NO:46,
SEQ ID NO:47, SEQ ID NO:48, SEQ ID NO:49, SEQ ID NO:50, SEQ ID NO:51,
5 SEQ ID NO:52, SEQ ID NO:55, SEQ ID NO:56, SEQ ID NO:57, SEQ ID NO:58,
SEQ ID NO:59, SEQ ID NO:61, SEQ ID NO:62, SEQ ID NO:63, SEQ ID NO:64,
SEQ ID NO:65, SEQ ID NO:66, SEQ ID NO:68, SEQ ID NO:70, SEQ ID NO:71,
SEQ ID NO:72, SEQ ID NO:73, SEQ ID NO:74, SEQ ID NO:75, SEQ ID NO:77,
SEQ ID NO:78, SEQ ID NO:79 and SEQ ID NO:80 or an analog thereof.

10

The present invention also refers to a diagnostic kit for detection and/or typing
and/or quantitation of HPV in a biological sample, comprising forward and reverse
primers, probes for HPV11, 16, 18, 26, 31, 33, 35, 39, 40, 42, 43, 44, 45, 51, 52,
53, 54, 56, 59, 59, 61, 66, 68, 70, 72, 73 and 82, internal positive control and DNA
15 polymerase.

BRIEF DESCRIPTION OF THE FIGURES AND THE TABLES

Figures 1-13: shows analysis of samples and amplification curves for consensus
multiplex realtime PCR type for the specific detection of HPV types 16, 18, 31, 33, 35,
20 39, 45, 51, 52, 56, 58, 59 and 66 wherein:

Figure 1: shows amplification of HPV18 in multiplex group 1.

Figure 2 : shows amplification of HPV52 in multiplex group 1.

Figure 3: shows amplification of HPV59 in multiplex group 1.

Figure 4: shows amplification of HPV39 in multiplex group 2.

25 Figure 5: shows amplification of HPV51 in multiplex group 2.

Figure 6: shows amplification of HPV56 in multiplex group 2.

Figure 7: shows amplification of HPV66 in multiplex group 2.

Figure 8: shows amplification of HPV16 in multiplex group 3.

Figure 9: shows amplification of HPV45 in multiplex group 3.

Figure 10: shows amplification of HPV58 in multiplex group 3.

Figure 11: shows amplification of HPV31 in multiplex group 4.

Figure 12: shows amplification of HPV33 in multiplex group 4

5 Figure 13: shows amplification of HPV35 in multiplex group 4.

Figure 14: shows amplification curves for the internal control.

Figures 15-20: show amplification of dilution series of HPV16, HPV18, HPV31, HPV33, HPV35 and HPV56 with primers PTf and PTr2.

Table 1: Refers to the list of designed primers.

10 Table 2: Refers to the list of designed probes.

Table 3: Refers to the composition of the Multiplex PCR of Multiplex Group 1, Multiplex Group 2, Multiplex Group 3 and Multiplex Group 4.

Table 4: Refers to evaluation of consensus realtime method on clinical samples.

Table 5: Refers to HPV Group 1, mastermix formulation.

15 Table 6: Refers to HPV Group 2, mastermix formulation.

Table 7: Refers to HPV Group 3, mastermix formulation.

Table 8: Refers to HPV Group 4, mastermix formulation

Table 9: Refers to the analytical sensitivity for the target HPV types.

20 DETAILED DESCRIPTION OF THE INVENTION

The quantitative realtime PCR method according to the present invention detects all the recognised oncogenic mucosal HPV types down to levels of 400 GU per sample or less. The method is intrinsically semi-quantitative since results are read off as C_T values which are proportional to the initial concentration of target DNA. It

is also possible to run the reactions in a fully-quantitative mode by including dilution series of the target type(s) in order to construct standard curves.

The quantitative nature of the method has several advantages. In epidemiological studies it would be possible to establish a constant cutoff value for all types, thus
5 eliminating the biased prevalence estimates that result from type-specific differences in analytical sensitivity. It has been suggested that only high-level HPV infections of greater 50 000 GU per sample are clinically significant, although this remains a matter of debate. In the consensus multiplex method it is possible to
10 adjust the cutoff to any quantitative value required. Lastly, use of a quantitative method may be of use in following the progress of infections and the effects of treatment.

The analytical sensitivity of the method is 4-400 GU for all types in the presence of 500ng of human DNA , which corresponds to 0.75×10^5 cells. Thus the presence of a single virus copy in 0.5% of the cells would give a positive result. At lower
15 human DNA concentrations (50ng), corresponding to 0.75×10^4 cells, the sensitivity improves by a factor of 10. This means that false positive results due to cell-poor samples are unlikely to be a problem. However, it also implies that operating the test in a fully quantitative mode will require standardisation of the cellular DNA concentration.

20 Compared with the hybrid capture (HCII) test, HPV was detected in more clinical samples, as is expected considering that the detection sensitivity of the consensus multiplex method according to the present invention was at least an order of magnitude higher.

25 The method according to the present invention thus covers all confirmed oncogenic HPV types in one test. The method has been specifically designed to give a four-tube format. This gives a method that optimally uses the 8x12 (96-well) and 16x24 (384 well) formats that are standard in most realtime PCR instruments. However, it is obvious that the method according to the present invention may be
30 used in a different format without departing from the scope of the invention.

Further, in the method according to present invention, only one forward primer and four reverse primers are required. This simplifies quality control and manufacture and reduces cost.

- 5 According to the present invention known sequence variation within the target region of the L1 gene has been taken into account prior to primer and probe design and all variant positions are either avoided or accounted for by modification of primer sequences.
- 10 The sensitivity of the method according to the present invention is measured in the presence of 500ng of human DNA, which correspond to $7,5 \times 10^4$ human cells and thus accurately mimics the conditions in a cell-rich clinical sample. This provides for an enhanced sensitivity for all targeted HPV types.
- 15 The method according to the present invention targets the L1 gene, which is much less prone to deletion.

The primers according to the present invention could with further modifications be used for the detection of related HPV types such as inter alia HPV11, 26, 40, 42, 43, 44, 53, 54, 61, 68, 70, 72, 73, 82.

- 20 In the real-time PCR instrument used in the method according to the present invention there are two vacant colour channels which allows the possibility of adding further oncogenic HPV types to the test.

- The multiplex grouping described in the present application is one of a number of possible groupings. For example the primers used in group 3 are known to be able
25 to amplify HPV types 18, 31, 33, 39, 51, 58, 59 and 66 and the internal control HPV6, and thus any of these types might therefore be assigned to this group.

- In the present invention a real-time PCR instrument is used having four colour channels. However with real-time PCR instruments having five, six, seven or more colour channels alternative variants of the test would be possible. The primers and
30 probes could then be combined in three or two parallel reactions for detection of the same HPV types.

In the present test format 20 samples can be analysed in one working day, making the test most suitable for low- to moderate-throughput laboratories. Larger batches of samples could be analysed using a 384-well format realtime PCR instrument, although the capacity of available automatic DNA extraction instruments is at present a more important limit on throughput than the PCR analysis time.

Realtime PCR is a closed-system test. It is thus not necessary to handle amplified samples in the open laboratory. Problems of PCR carry-over contamination are thus greatly reduced, and this in turn reduces the resources that must otherwise be devoted to the stringent contamination control associated with open-laboratory post-PCR analysis.

The method according to the present invention is inter alia of value in screening programs where HPV genotyping is preferred.

Preferred embodiments of the invention are described. However it is to be understood that the invention is not limited to those precise embodiments, and that various changes and modifications may be effected therein by one skilled in the art without departing from the scope or spirit of the invention as defined in the appended claims. The following examples illustrate, but do not limit the invention.

EXAMPLES

Example 1- Plasmids

Plasmids containing the target region of the L1 genes of HPV6, 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 and 66 were transformed into *Escherichia coli* DH5-alpha and isolated by caesium chloride density gradient centrifugation (Fritsch et al) or Qiagen HiSpeed Midi Kit (Qiagen GmbH, Hilden, Germany). DNA concentrations were determined by ultraviolet spectrophotometry.

Example 2- Clinical samples

Clinical samples were cervical cytobrush samples sent to our laboratory for routine HPV analysis for secondary screening of low grade dysplasia (LSIL) or equivocal cytology (ASCUS or inadequate sample). Samples were collected and transported
5 in CYTYC thin prep transport (medium CYTYC, Crawley, UK) according to the manufacturers recommendations. 10 ml of sample was centrifuged at 3000 rpm for 10 min and the pellet was resuspended in 100µl phosphate buffered saline before extraction of DNA using the MagNAPure automatic DNA extraction instrument with the MagNAPure LC DNA extraction kit (Roche Diagnostics, Penzberg, Germany).
10 DNA was eluted in 100µl. Clinical specimens that had previously been tested with the HCII test and typed using a consensus PCR and reverse line blot were used to evaluate the test.

Example 3- Hybrid Capture Test

15 The Digene Hybrid Capture Test (Digene, Gaithersburg, MD) was performed using 25 ml of sample according to the manufacturers' instructions.

Example 4- Primer Design

In order to identify suitable targets for consensus PCR primers, the sequences of
20 the L1 genes of HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 and 66, plus HPV6 were aligned using CLUSTALW. In order to generate input sequences for the CLUSTAL alignments the prototype sequence of each HPV type was used in a BLAST search in order to identify intratypic variant sequences. The sequence was then edited at variant position by the substitution of the appropriate IUPAC
25 ambiguity code. The regions corresponding to nt 6348 – 6374 and 6579-6557 on the HPV16 genome were identified as highly conserved and were chosen for the placing of the forward and reverse primers respectively. This led to the construction of Primers PTf and PTr. As this primer pair did not amplify all target types satisfactorily, a phenomenon which could be attributed to low stability of
30 duplexes with PTr, less degenerate version of PTr were sought.

The sequences were grouped into four homology groups on the basis of a distance tree of the concatenated forward and reverse primer target sequences constructed in CLUSTALW (<http://www.ebi.ac.uk/Tools/clustalw2/index.html>). The groups were HPV18, HPV52 and HPV59 (Group 1), HPV39, HPV51, HPV56 and HPV66 (Group 2), HPV16, HPV45 and HPV58 (Group 3) and HPV31, HPV33 and HPV35 plus the internal control HPV6 (Group 4). For each group a new CLUSTAL alignment was constructed and group-specific reverse primers were designed in order to reduce primer degeneracy and improve sensitivity. Variable positions were accommodated by using inosine as a wildcard base, thymine as a universal pyrimidine, guanine as a universal purine (exploiting the stability of the noncanonical G:T basepair) or incorporating mixed nucleotides. This led to the design of primers PTrGr1, PTrGr2 and PTrGr4. The original PTr sequence was retained for group 3.

An alternative target region corresponding to positions 6644 – 6620 on HPV16 was also investigated for placing of the reverse primer. The primer PTr2, was designed to target this region. PTr2 is suitable for the amplification of HPV types 16, 18, 31, 33, 35, 45, 52, 56, 58, 58 and 59.

Tabell 1. List of primers

Primer	Sequence	SEQ ID NO
PTf	5-TG(CT)AAA TAT CC(AT) GAT TAT (AT)TI IAA ATG-3	1
PTrGr1b	5-TGT AGC CAG TAT GG(CT) TTA TTG AA-3	2
PTrGr2	5-TG(GT) AGC CAA TAA GGT TTA TTA AA-3	3
PTr957	5- TGI AIC CAA TAI GG(CT) TTA TTG AA -3	4

PTrGr4	5- TGT AIC CAA TAT GGt TTA TTA AA-3	5
PTr2	5- C(GT)G T(AG)G TAT CIC CAI CAA GTA ACA AA-3	6

Example 5- Probe Design

TaqMan MGB probes were designed using Primer Express software (Applied Biosystems, Foster City, CA) using the target amplicon as input sequence and specifying an amplicon length of ≥ 200 bp. The software was set to find primers with a T_m of 60°C and probes with a T_m 10°C higher (70°C). Where necessary, in order to compensate for the software's inability to locate probes without simultaneously locating primers, the sequence was bracketed by dummy primer sequences. TaqMan MGB probes were supplied by Applied Biosystems and were labelled with FAM, VIC or NED which are the recommended labels for the three standard colour channels of the AB Prism 7000 series of realtime PCR instruments. MGB probes incorporate a minor groove binder, which raises probe T_m and allows shorter probes to be used, and a dark (non-fluorescent quencher). A fourth channel, usually used for a passive reference, is available on ABI 7000 series instruments. This channel ideally uses the dye ROX. ROX-labelled MGB probes are not obtainable. In order to achieve the required T_m , two approaches were employed. In the first approach probes were first designed as TaqMan MGB probes, then redesigned as LNA probes on the basis of the suggestions of the suppliers (Sigma Genosys,) or using the website <http://lna-tm.com/>. LNA probes incorporate a modified deoxyribose in the sugar-phosphate backbone that stabilises the helical structure and increases T_m . Alternatively Primer Express was used to design conventional (non-MGB) TaqMan probes using the program module for this purpose and probes <22 bp were further evaluated. All probes were designed with a T_m of $70 \pm 2^\circ\text{C}$. ROX-labelled probes were fluorophore-labelled at the 5' end and quenched with the dark quencher BHQ2 at the 3' end.

Candidate probes were compared with CLUSTALW multiple sequence alignments of the amplicon region of the thirteen target HPV types and HPV6 with variant positions marked. Probe sequences targetting intra-type invariant but inter-type specific sequences were selected for further testing.

5

Example 6- PCR Program

All PCR reactions were run on the ABI Prism 7000 or 7300 genetic analyzer (Applied Biosystems, Foster City, CA) with the ROX passive reference function disabled in order to allow the ROX colour channel to be used for detection. For all
10 experiments the PCR program was: 95°C, 12 min followed by 45 cycles of 95°C, 15s; 45°C, 30s; 60°C, 60s. For experiments using SYBR green, dissociation analysis was appended to the end of the program. The reaction volume was 50µl and 5µl of sample DNA was used in all reactions. PCR mastermixes were used according to the manufacturers instructions.

15

Example 7- Testing of Primers

Primers were used at a concentration of 600nM. Analytical sensitivity was assessed for all target HPV types by amplifying 10-fold dilution series of HPV plasmid DNA in 100 ng/µl human DNA (Sigma Aldrich, St Louis, MI). Amplified
20 DNA was detected using SYBR green with the Finnzymes DyNAmo SYBR green qPCR kit F400 (Finnzymes, Espoo, Finland) using the PCR program described combined with dissociation analysis. Analytical sensitivity was defined as the lowest target concentration that gave an unequivocal type-specific dissociation peak. The criterion for primer acceptability at this stage of the analysis was an
25 analytical sensitivity of 4000 GU or better for all target HPV types.

Example 8- Probe assessment and titration

Probes were initially used at a concentration of 200 nM. Assessment of probes was done by TaqMan PCR assay using Finnzymes DyNAmo probe PCR kit F450

(Finnzymes, Espoo, Finland) or, in later experiments, PerfeCta qPCR FastMix, UNG Quanta Mastermix (Quanta Biosciences, Gaithersburg, MD) using a 10-fold serial dilution of the target HPV type in 100 ng/μl human DNA and the cognate primers. A positive amplification signal from 400 GU of target DNA was the
5 minimum criterion for approval of the probe for further testing.

When used at high concentrations in multiplex reactions, probes tend to interfere with one another. Therefore probe titration was performed in order to determine the minimum probe concentration that could be used without loss of sensitivity. Probe concentrations of 200, 150, 100 and 50 nM were tested against 4000, 400,
10 40 and 4 GU of the target type. The lowest concentration of probe that could be used without loss of analytical sensitivity was used in multiplex reactions.

Example 9- Testing of Probes in Multiplex

For multiplex testing, the probes were combined at the chosen concentration
15 together with the cognate primers and tested against dilution series of the three or four target HPV types in PCR reactions as described. Criteria for acceptance of a multiplex combination were analytical sensitivity 400 GU or better and no cross-reaction with the other types in the group. Where unacceptable loss of analytical sensitivity or cross-reaction were detected in a multiplex group, reactions were
20 repeated with fewer probes in all combinations in order to identify the interfering probe.

Successful groups of probes were tested against $4 \cdot 10^4$ GU of each of the fourteen HPV types in the test in order to detect cross-reaction. Where cross-reaction was observed, the cross-reacting probe was tested against the cross-reacting HPV
25 type(s) in a simplex reaction in order to determine whether the cross-reaction was due to a direct interaction between the probe and the cross-reacting HPV type, or due to secondary effects caused by multiplexing.

Where probes were rejected due to interference or cross-reaction, new probes were designed and tested as described above.

30 The probes which were assessed are shown in table 2.

Table 2. List of probes

Grp	HPV	Sequence	SEQ ID NO	Label	Type	Target
1	PTp18(7)	ATA CAC ACA GCT GCC AGG T	21	FAM	TaqMan MGB	HPV18
1	PTp18(8)	CTTCACCTGGCAGCTG	22	FAM	TaqMan MGB	HPV18
1	PTp18(130)	CACTGTGCCTCAATC	23	FAM	TaqMan MGB	HPV18
1	PTp18(131)	ACTGTGCCTCAATCC	7	FAM	TaqMan MGB	HPV18
1	PTp18(177)	ACCTGGCAGCTGTGT	25	FAM	TaqMan MGB	HPV18
1	PTp 52	TCT AAC TCT GGC AAT A	26	ROX/BHQ2	TaqMan MGB	HPV52
1	PTp52(4)	tc+Taa+Ctc+Tggc+Aa+Tac+Tgc	8	ROX/BHQ2	LNA	HPV52
1	PTp 59(1)	TAC TGA CAT ACG TGC C	9	NED	TaqMan MGB	HPV59
1	PTp59(77q)	AGGCGAGAACAGGTTTTTGCCAGACA	29	NED	TaqMan MGB	HPV59
1	PTp59(132)	ATCAACTTCCTGAATCACT	30	NED	TaqMan MGB	HPV59
1	PTp59(159)	AAGGTACTGACATACGTGCC	31	NED	TaqMan MGB	HPV59
2	PTp 39(1)	CAC AGA TAT ACG TGC AAA	32	ROX/BHQ2	TaqMan MGB	HPV39
2	PTp39(4)	ca+Ca+Gatat+Acg+Tg+Caa+Acc	10	ROX/BHQ2	LNA	HPV39
2	PTp39(5)	caca+Gatat+Acg+Tg+Caa+Acc	34	ROX/BHQ2	LNA	HPV39
2	PTp 51	AGT GGT AAT GGC CGT G	35	FAM	TaqMan MGB	HPV51
2	PTp51(156)	AGTGGTAATGGCCGTGAC	11	FAM	TaqMan MGB	HPV51

2	PTp51(158)	TGGTAATGGCCGTGACC	37	FAM	TaqMan MGB	HPV51
2	PTp51(168)	CGTGACCCTATAGAAAG	38	FAM	TaqMan MGB	HPV51
2	PTp 56	TAG CAA TGG TAG AGA AC	39	VIC	TaqMan MGB	HPV56
2	PTp56(2)	CAATACCTGCAGAGTTATA	12	VIC	TaqMan MGB	HPV56
2	PTp56(3)	ATATTAAAGGGTAGCAATGGT	41	VIC	TaqMan MGB	HPV56
2	PTp 66(1)	TTG GAA GGG TGG CAA T	13	NED	TaqMan MGB	HPV66
3	PTp 16(6)	ATT TGC AGT AGA CCC AGA G	14	FAM	TaqMan MGB	HPV16
3	PTp16(94)	CTGTTGGTGAAAATGTACC	44	FAM	TaqMan MGB	HPV16
3	PTp16(101)	TGAAAATGTACCAGACGATT	45	FAM	TaqMan MGB	HPV16
3	PTp16(155s)	AGCCAGTTCAAATTAT	46	FAM	TaqMan MGB	HPV16
3/4	PTp 33	TCT ATT CAA AGC AGT GC	47	VIC	TaqMan MGB	HPV33
3/4	PTp33(2)	ctac+Tgc+Ctc+Ta+Tt+Ca+Aag+Cag	48	ROX/BHQ2	LNA	HPV33
3/4	PTp33(8)	TGCCTCTATTCAAAGCAGTGC	49	NED, VIC	TaqMan MGB	HPV33
3/4	PTp33(10)	GGAAC TACTGCCTCTATTCAAAGCAGT GCTGC	50	VIC	TaqMan	HPV33
3 /4	PTp33(141)	AGAGGCTGTTCCCGATGACCTG	51	ROX	TaqMan	HPV33
3/4	PTp33(149)	TTCCCGATGACCTGTACATTAAAGGT	52	ROX	TaqMan	HPV33
3/4	PTp33(9)	tg+Cctctatt+Caaag+Cagtgc	19	ROX/BHQ2	LNA	HPV33

3	PTp 58	TCC GGT AAT ACT GCA G	16	NED	TaqMan MGB	HPV58
3	PTp58(36)	CCAGTGAACCTTATGGG	55	NED	TaqMan MGB	HPV58
3	PTp58(118 m)	AAAACTTGGCGAGGCTG	56	NED	TaqMan MGB	HPV58
3	PTp58(156)	AAGGGTCCGTAATAC	57	NED	TaqMan MGB	HPV58
3	PTp58(170)	ACTGCAGTTATCCAAAGTA	58	NED	TaqMan MGB	HPV58
3/1	PTp 45(3)	ATA CAC ACA ACT GCC A	59	VIC	TaqMan MGB	HPV45
3/1	PTp45(6)	GTGAAACCCCTGGCAGT	15	VIC	TaqMan MGB	HPV45
3/1	PTP45(4)	aa+Tac+Aca+Caa+Ctg+Ccagg	61	ROX/BHQ2	LNA	HPV45
3/1	PTp45(5)	ct+Agc+Gct+Aat+Atg+Cgtgaa	62	ROX/BHQ2	LNA	HPV45
3/1	PTp45(172 Q)	TGCGTGAAACCCCTGGCAGTTG	63	VIC	TaqMan	HPV45
3/1	PTp45(127)	ACACAGTACCTACGGACCTA	64	VIC	TaqMan MGB	HPV45
3/1	PTp45(164)	CGCTAATATGCGTGAAAC	65	VIC	TaqMan MGB	HPV45
4	PTp 35	TAC CAC TGG CAC ATT G	66	FAM	TaqMan MGB	HPV35
4	PTp35(4)	GTACCACTGGCACATTG	20	FAM	TaqMan MGB	HPV35
4	pTP35(5)	TACCACTGGCACATTG	68	FAM	TaqMan MGB	HPV35
4	PTp 31	TCC GGT TCA ACA GCT AC	18	VIC	TaqMan MGB	HPV31
4	PTp31 (115g)	TCAGGCGCGTTGGTGAATCGGT	70	VIC	TaqMan	HPV31

4	PTp31(168)	TTCAACAGCTACTTTAGC	71	VIC	TaqMan MGB	HPV31
4	PTp31(162)	CTCCGGTTCAACAGCTA	72	VIC	TaqMan MGB	HPV31
4	PTp6(2)	TGCCTGATACTCTTATAATTAA	73	NED	TaqMan MGB	HPV6
4	PTp6(197)	ACCCCGAGCGGCT	74	NED	TaqMan MGB	HPV6
4	PTp6(190)	TGTTAACACCCCGAGCGG	75	NED	TaqMan MGB	HPV6
4	PTp6(154)	TAGTGGAATCGCACGTCT	17	NED	TaqMan MGB	HPV6
4	PTp6(101)	AACAGGGCTGGCGAGGTGGG	77	ROX/BHQ2	TaqMan	HPV6
4	PTp6(98)	TTTAACAGGGCTGGCGAGGTGGG	78	ROX/BHQ2	TaqMan	HPV6
4	PTp6(122)	ACAGGGCTGGCGAGGTGGG	79	ROX/BHQ2	TaqMan	HPV6
4	PTp6(116)	TTTTTAACAGGGCTGGCGAGGT	80	ROX/BHQ2	TaqMan	HPV6

The multiplex group to which the types are currently assigned is given in the first column. Locked nucleic acid (LNA™) sequences are written with lower case letters for normal residues and uppercase letters prefixed with '+' for LNA residues.

- 5 Probes in use in the multiplex mixes shown below are in bold text. Grey text indicates that the probe has been replaced. The fluorophors used (FAM, VIC, NED and ROX) are those used in the present method. Alternative fluorophors and combinations of fluorophors may be used as known by the skilled man.
- 10 After 1 – 3 rounds of probe/primer redesign, the multiplex formulation in Table 3 was arrived at.

Table 3. Composition of the Multiplex PCR mixes**Multiplex Group 1. HPV 18,52,59.**

Primer	Sequence	SEQ ID NO	Concentration
PTf	5-TG(CT)AAA TAT CC(AT) GAT TAT (AT)TI IAA ATG-3	1	600nM
PTGr1b	5-TGT AGC CAG TAT GG(CT) TTA TTG AA-3	2	600nM

Probe	Sequence	SEQ ID NO	Type	Concentration
PTp18(131)	5- FAM ACT GTG CCT CAA TCC-3	7	TaqMan MGB	150nM
PTp52(4)	5- ROX TC+TAA+CTC+TGGC+AA+ TAC+TGC-3	8	LNA	50nM
PTp59(1)	5- NED TAC TGA CAT ACG TGC C-3	9	TaqMan MGB	100nM

5 **Multiplex Group 2. HPV 39,51,56,66**

Primer	Sequence	SEQ ID NO	Concentration
PTf	5-TG(CT)AAA TAT CC(AT) GAT TAT (AT)TI IAA ATG-3	1	600nM
PTGr2	5-TG(GT) AGC CAA TAA GGT TTA TTA AA-3	3	600nM

Probe	Sequence	SEQ ID NO	Type	Concentration
PTp 39(4)	5-ROX CA+CA+GATAT+ACG+TG+CAA+ACC-3	10	LNA	50nM
PTp 51(156)	5- FAM AGT GGT AAT GGC CGT GAC-3	11	TaqMan MGB	100nM
PTp 56(2)	5- VIC CAA TAC CTG CAG AGT TAT A-3	12	TaqMan MGB	50nM
PTp 66 (1)	5- NED TTG GAA GGG TGG CAA T-3	13	TaqMan MGB	100nM

Multiplex Group 3. HPV 16,45,58

Primer	Sequence	SEQ ID NO	Concentration
PTf	5-TG(CT)AAA TAT CC(AT) GAT TAT (AT)TI IAA ATG-3	1	600Nm
PTr957	5- TGI AIC CAA TAI GG(CT) TTA TTG AA -3	4	600Nm

Probe	Sekvens	SEQ ID NO	Type	Concentration
PTp16(6)	5- FAM ATT TGC AGT AGA CCC AGA G-3	14	TaqMan MGB	150Nm
PTp45(6)	5- VIC GTG AAA CCC CTG GCA GT-3	15	TaqMan MGB	150Nm
PTp58	5- NED TCC GGT AAT ACT GCA G-3	16	TaqMan MGB	50Nm

Multiplex Group 4. HPV 6,31,33,35

Primer	Sequence	SEQ ID NO	Concentration
Oli1+	5- TG(CT)AAA TAT CC(AT) _ATT AT (AT)TI IAA ATG -3	1	600nM
PTrGr4	5- TGT AIC CAA TAT GGt TTA TTA AA-3	5	600nM

Probe	Sekvens	SEQ ID NO	Kons.i Quantamix	C.off
PTp 6(154)	5-TAG TGG AAA TCG CAC GTC T- 3 NED	17	50nM	6-7
PTp 31	5-TCC GGT TCA ACA GCT AC-3 VIC	18	50nM	31-9
PTp 33(9)	5-ROX tg(+C)ctctatt(+C)aaag(+C)agtcg-BHQ2	19	150nM	33-7
PTp 35(4)	5- GTA CCA CTG GCA CAT TG-3 FAM	20	100nM	35-7

Example 10- Internal Control

- 5 A plasmid containing the genome of the non-oncogenic HPV type HPV6 was used as an internal control. Design and testing of the HPV6 internal control probe was as for the other probes. In order to determine a suitable concentration of internal control, HPV6 was added to group 4 PCR reactions at 10x serially diluted concentrations from $4 \cdot 10^3$ GU to 40 GU per reaction. The multiplex mix containing
- 10 internal control was used in amplification test using HPV types 31, 33 and 35 10x dilution series in order to determine which internal control concentrations, if any, reduced analytical sensitivity for the target types.

Example 11- Evaluation against clinical samples

- 15 Once the final formulation of the test was arrived at based on the results of the above tests, the test was evaluated on 37 clinical specimens and the results were

compared with the previous results of the hybrid capture HCII test and previous typing results using a reverse line-blot test.

Table 4 compares the results of consensus multiplex PCR for 37 clinical samples, compared with the hybrid capture (HCII) test and the results of previous typing using reversed line blot. No false negative results relative to the previous tests were encountered. One sample was false negative relative to the HCII test and in one sample the consensus multiplex test failed to detect a type found by the previous typing test.

The results thus shows that in 5 samples the consensus multiplex test according to the present invention detected additional types not detected in the previous analyses.

Table 4. Evaluation of Consensus Realtime Test on Clinical Samples

Sample	Result HCII	Previous Typing Result	Result Consensus Multiplex
1	+	HPV16	HPV16
2	-	-	-
3	-	-	-
4	+	HPV35	HPV16, 35, 52
5	-	-	-
6	+	HPV31	HPV31, 52
7	-	-	-
8	+	HPV51	HPV51, 56
9	+	HPV51	HPV51
10	-	-	-
11	-	-	-
12	-	-	HPV51, 66
13	-	-	-
14	-	-	-
15	-	-	-

16	-	-	-
17	-	-	-
18	+	HPV39	HPV39
19	-	-	-
20	-	HPV51	HPV51
21	-	-	HPV16, 45
22	-	-	-
23	-	-	-
24	-	-	-
25	-	-	-
26	+	-	-
27	-	-	-
28	-	-	-
29	-	-	-
30	±	HPV45	HPV45
31	-	-	-
32	+	HPV16	HPV16
33	-	-	-
34	-	HPV31	HPV31
35	+	HPV51, 52, 59	HPV51, 52
36	-	-	-
37	-	-	-

Table 4 thus compares the results of consensus multiplex PCR for 37 clinical samples, compared with the HCII test and the results of previous typing using reversed line blot. There was complete agreement of results for 28 samples, of which five were positive and 23 were negative. For two samples HCII was negative while reverse line blot and consensus multiplex PCR gave identical results. In 5 samples the consensus multiplex test detected extra types not detected by reverse line blot. These types were HPV66 (one sample), which is not included in the reverse line blot assay, HPV52 (two samples) and HPV56 (one sample) for which the reverse line blot assay is known to have reduced sensitivity

and HPV35, HPV16 and HPV45 (one case each). In one sample the consensus multiplex test failed to detect a type (HPV59) detected by reverse line blot. One sample was positive only by the HCII test.

5 Example 12- Consensus realtime PCR setup

The following procedure was used to analyse 21 clinical samples. The samples were DNA from cervical brush samples extracted using the MagNAPure automatic DNA extractor.

- 10 1. PCR mixes for multiplex groups 1,2,3 and 4 were prepared in separate tubes according to the tables below.

Table 5

HPV Group 1, mastermix formulation: 50µl reaction.

Primer concentration 600nM.

Probe concentration 150 nM for HPV18, 50nM for HPV52 and 100nM for HPV59.

Component	Volume Per sample	Volume 26 samples and controls
Rx-mix QuantaFast Cycling qPCR	25µl	650µl
PTf 100 pmol/µl	0.3 µl	7.8µl
PtrGr1 100 pmol/µl	0.3 µl	7.8µl
Probes 18(131), 52(4), 59(1) 5 pmol/µl	1.5+0.5+1 .0µl	39+13+ 26µl

Table 6**HPV Group 2, mastermix formulation: 50µl reaction.**

Primer concentration 600nM.

Probe concentration 50nM for HPV39, 100nM for HPV51, 50nM for HPV56 and

5 100nM for HPV66

Component	Volume Per sample	Volume 26 samples and controls
Rx-mix QuantaFast Cycling qPCR	25µl	650µl
PTf 100 pmol/µl	0.3 µl	7.8µl
PTrGr2 100pmol/µl	0.3 µl	7.8µl
Probes 39(4), 51(156), 56(2)66(1) 5 pmol/µl	0.5+1.0+0 .5+ 1.0µl	13+26+ 13+26µl
H₂O	16.4µl	426µl

Table 7**HPV Group 3, mastermix formulation: 50µl reaction.**

Primer concentration 600nM.

Probe concentration 150 nM for HPV16, 150nM for HPV45 and 50nM for HPV58

5

Component	Volume Per sample	Volume 26 samples and controls
Rx-mix QuantaFast Cycling qPCR	25µl	650µl
PTf 100 pmol/µl	0.3 µl	7.8µl
PTr 100 pmol/µl	0.3 µl	7.8µl
Probes 16(6), 45(6), 58 5 pmol/µl	1.5+1.5+0 .5µl	39+39+ 13µl
H₂O	15.9µl	413µl

Table 8**HPV Group 4, mastermix formulation: 50µl reaction.**

Primer concentration 600nM.

Probe concentration 50 nM for HPV6, 50nM for HPV31, 150nM for HPV33 and
100nM for HPV35

Component	Volume Per sample	Volume 26 samples and controls
Rx-mix QuantaFast Cycling qPCR	25µl	650µl
PTf 100 pmol/µl	0.3 µl	7.8µl
PtrGr 4 100 pmol/µl	0.3 µl	7.8µl
Probes 6(154), 31, 33(9), 35(4) 5 pmol/µl	0.5+0.5+1 .5+ 1.0µl	13+13+ 39+26µl
IK	0.5µl	13µl
H₂O	15.9µl	413µl

2. PCR program: 95°C, 12 min followed by 45 cycles of 95°C, 15s;
45°C, 30s; 60°C, 60s. AB7300 realtime PCR machine. The following 96
well plate setup is used:

Rows 1, 2, 3: detectors for Group 1,

Rows 4, 5, 6: detectors for Group 2,

Rows 7, 8, 9 detectors for Group 3 and

Row 10, 11, 12 detectors for Group 4.

Sample number for each of the 21 samples is also entered.

Each sample is run in all groups. In addition, group specific control mixes
5 containing all the HPV types in the group and negative controls are
included.

Thus the following controls are included:

Positive controls

13-type control, containing 4000 GU of HPV16, 18, 31,33, 35, 39, 45, 51,
10 52, 56, 58, 59 and 66 per 5µl aliquot.

Group positive controls, containing 4000 GU of HPV18, 52 and 59 (group
1), HPV39, 51, 59 and 66 (group 2), HPV16, 45 and 58 (group 3) and HPV
31, 33 and 35 (group 4) were used with their respective groups.

For each multiplex group a 13-type multiplex control, a group control and a
15 negative control (no DNA added) were run.

3. PCR mixes as described in (1) are aliquoted , 45µl to a well,
following the plate setup.

5 µl sample DNA, positive control or negativ control are added according to
20 the plate setup.

The PCR plate is sealed and placed in the PCR machine. The PCR
program takes approximately 2.5 hours.

Logistic considerations

The time required for completion of the analysis of 20 samples is 8 hours – 4.5
25 hours for sample processing, 15 minutes for PCR setup, 2.5 hours PCR run time
and 1 hour for analysis. Of this, 2.5 hours is hands-on time.

Analysis of samples and amplification curves for consensus multiplex realtime PCR type for the specific detection of HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 and 66 are shown in Figures 1-13. The vertical axis is log (fluorescence), background subtracted, for the appropriate colour channel. The horizontal axis is PCR cycle. The horizontal line marked with an arrowhead on the right hand vertical axis is the threshold line which is placed manually after correction for background fluorescence. Amplification curves crossing the threshold line are considered positive, and the fractional PCR cycle where the amplification curve crosses the threshold line is the threshold cycle (C_T) given in the accompanying report. Figure 4 shows amplification of the internal control, HPV6.

Results for Multiplex Group 1. HPV18, 52 and 59.

Figure 1 shows amplification of HPV18 in multiplex group 1. The two positive curves are positive controls. There are no HPV18-positive samples.

Figure 2 shows amplification of HPV52 in multiplex group 1. 4 positive curves are visible. The orange and grey curves to the right are positive controls. The two blue curves are positive samples (samples four and 6) containing HPV52.

Figure 3 shows amplification of HPV59 in multiplex group 1. The two positive curves are positive controls. There are no HPV59-positive samples.

Results for Multiplex Group 2. HPV39, 51, 56 and 66.

Figure 4 shows amplification of HPV39 in multiplex group 2. The green curve is a positive control containing the target types of group 2. In a second positive control containing all fourteen target types HPV39 failed to amplify. The black curve is a strongly positive sample (sample 18) containing HPV39.

Figure 5 shows amplification of HPV51 in multiplex group 2. Six positive curves are visible. The green and brown curves crossing the threshold line at cycle 29 are positive controls. The four remaining curves are HPV51-positive samples, two of which (samples 8 and 9 are very strongly positive). The black and blue curves for

these two samples are almost coincident. Sample 20 (green curve) crosses the threshold line close to the positive controls, indicating that it contains similar amounts of HPV51 (4000 GU). Sample 12 (red curve) is weakly positive.

Figure 6 shows amplification of HPV56 in multiplex group 2. 3 positive curves are visible. Green and brown curves are positive controls. The blue curve is a positive sample (sample 8). The flattened curve is due to the competitive effect of the very high concentrations of HPV51 in this sample.

Figure 7 shows amplification of HPV66 in multiplex group 2. The green and brown positive curves at the extreme right are positive controls. The red curve is a positive sample.

Results for Multiplex Group 3. HPV16, 45 and 58.

Figure 8 shows amplification of HPV16 in multiplex group 3. 4 positive curves are visible. The brown and red curves crossing the threshold line at cycle 34 are positive controls. The brown curve is noticeably flatter than the red curve, probably because the positive control represented by the brown curve contains all fourteen target types, producing a greater potential for competition. The black and green curves are HPV16-positive samples (samples 1 and 21) one of which is very weakly positive.

Figure 9 shows amplification of HPV45 in multiplex group 3. 3 positive curves are visible. The brown and red curves are positive controls. The green curve is an HPV45 positive sample (sample 21).

Figure 10 shows amplification of HPV58 in multiplex group 3. The two positive curves are positive controls. There are no HPV58-positive samples.

Results for Multiplex Group 4. HPV31, 33, 35 and Internal Control (HPV6)

Figure 11 shows amplification of HPV31 in multiplex group 4. 3 positive curves are visible. The two similar green curves are positive controls. The blue curve is a positive sample (sample 6) containing large amounts of HPV31.

Figure 12 shows amplification of HPV33 in multiplex group 4. The two similar green curves are positive controls which are the only positive signals. There are no HPV33-positive samples.

Figure 13 shows amplification of HPV35 in multiplex group 4. 3 positive curves are visible. The two similar green curves are positive controls. The maroon curve is a positive sample (sample 4) containing HPV35.

Figure 14 shows internal control amplification curves. The samples all contain the same amount of HPV6 and the curves are clustered, indicating absence of inhibition. The curves to the right of the main cluster of curves are indicative of partial inhibition or competition. Five such curves are seen. In one curve (sample 4) the amplification of the internal control is nearly completely inhibited. This is caused by competition with a strong HPV35 amplification signal in the same tube. A second slightly inhibitory sample, sample 9, contains HPV51 at unusually high concentrations. Although HPV51 is not included in multiplex group 4, the similarity of the primer sequences is probably sufficient to generate a competitive amplification in such a hot sample. The remaining three samples (3, 8 and 17) are negative and are probably partially inhibitory.

The internal control, HPV6 could thus be added to the group 4 multiplex reaction at levels of 400 - 4000 GU per reaction without compromising the analytical sensitivity for HPV31, HPV33 or HPV35. 4000 GU per reaction was chosen as standard concentration for the internal control.

Table 9. Analytical sensitivity

HPV type	Analytical sensitivity (GU)	Analytical sensitivity (GU)
	(100 ng/μl human DNA)	(10 ng/μl human DNA)
HPV16	40	4
HPV18	40	4
HPV31	4	4
HPV33	400	40
HPV35	400	40
HPV39	40	4
HPV45	400	4

HPV51	4	4
HPV52	400	40
HPV56	40	4
HPV58	400	40
HPV59	40	4
HPV66	400	40

This table refers to the analytical sensitivity for multiplex PCR reactions. Reactions were performed using the multiplex formulations described in tables 5-8. Target DNA was a 10 fold serial dilution of HPV plasmid in 100 ng/μl or 10ng/μl human DNA, corresponding to 500 or 50ng respectively of human DNA in the sample. The analytical sensitivity is the lowest number of copies (GU) that gives an amplification curve that can be clearly distinguished from background fluorescence.

Example 13- Use of primer PTr2

Amplification was performed with dilution series of HPV16, HPV18, HPV31, HPV33, HPV35 and HPV56 with primers PTf and PTr2.

The amplification was thus performed in a 100μl reaction using 500 ng of each primer, 1,5u Amplitaq Gold (Applied Biosystems), 1.5 mM MgCl₂, and 200 μM dNTPs in Amplitaq Gold buffer. The PCR program was 93°C, 30s; 40°C, 30s; 72°C, 90s; 40 cycles. 10μl of the amplified material was applied to a 1.5% agarose gel and subjected to electrophoresis at 130 V for 90 minutes, stained with SYBR green according to the manufacturers instructions and photographed under 300 nm UV illumination. The results are shown in Figures 15-22.

PTr2 is thus suitable for the amplification of HPV types.

Example 14- Diagnostic kit

The diagnostic kit according to the present invention may have the following construction, but it is obvious that various changes and modifications may be

effected by one skilled in the art without departing from the scope or spirit of the invention.

The kit may contain four rows of six tubes (Tubes 1 – 6) containing all the reagents necessary for detection of HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59
5 and 66.

Tube 1 contains 500 µl of a solution containing forward primer PTf and reverse primer PTrGr1, plus TaqMan MGB probes for HPV18 (FAM label), HPV59 (NED label) and an LNA probe for HPV52 (ROX label).

10 **Tube 2** contains 500 µl of a solution containing forward primer PTf and reverse primer PTrGr2, plus TaqMan MGB probes for HPV51 (FAM label), HPV56 (VIC label) HPV66 (NED label) and an LNA probe for HPV39 (ROX label).

Tube 3 contains 500 µl of a solution containing forward primer PTf and reverse primer PTr, plus TaqMan MGB probes for HPV16 (FAM label), HPV45 (VIC label) and HPV58 (NED label)

15 **Tube 4** contains 500 µl of a solution containing forward primer PTf and reverse primer PTrGr4, plus TaqMan MGB probes for HPV35 (FAM label), HPV31 VIC label, and an LNA probe for HPV52 (ROX label). Also included is a TaqMan probe for the HPV6 internal control and HPV6 internal control target at a concentration of 800 copies per µl.

20 **Tube 5** contains 20µl of positive control containing 4000 copies of HPV types 16, 18, 31, 33, 35, 35, 39, 45, 51, 52, 56, 58, 59 and 66.

Tube 6 contains 2500 µl of 2 x PCR mastermix containing thermostable DNA polymerase and deoxynucleotide triphosphates and magnesium chloride in a proprietary buffer solution.

25 According to the following procedure 23 samples may be analysed simultaneously. Three of these samples may be negative controls. Further, DNA is extracted from the samples using a suitable method.

1. Each of tubes 1 to 6 are allowed to thaw at room temperature.

2. Group 1 mastermix of example 12 is prepared by adding 600 µl of PCR mastermix (tube 6) to tube 1 and mixed by inversion.
3. Group 2 mastermix of example 12 is prepared by adding 600 µl of PCR mastermix (tube 6) to tube 2 and mixed by inversion.
- 5 4. Group 3 mastermix of example 12 is prepared by adding 600 µl of PCR mastermix (tube 6) to tube 3 and mixed by inversion .
5. Group 4 mastermix of example 12 is prepared by adding 600 µl of PCR mastermix (tube 6) to tube 4 and mixed by inversion.
- 10 6. 45µl of group 1 mastermix is transferred to each well of rows A and E of a 96-well 500 µl optical PCR plate.
7. 45µl of group 2 mastermix is transferred to each well of rows B and F.
8. 45µl of group 3 mastermix is transferred to each well of rows C and G.
9. 45µl of group 4 mastermix to each well of rows D and H.
10. Add 5 µl of sample DNA 1 to wells A1, B1, C1 and D1.
- 15 11. Add 5 µl of sample DNA 2 to wells A2, B2, C2 and D2.
12. Continue adding samples in the same pattern to wells A3 – D3 through A12 – D12.
13. Continue adding samples in the same pattern to wells D1 – H1 through D11 – H11.
- 20 14. Add 5 µl of positive control to each of wells D12 – H12.
15. Seal the plate, place in the realtime PCR machine and start the PCR program.
16. HPV types 18, 52 and 59 will be detected in rows A and E.
17. HPV types 39, 51, 56 and 66 will be detected in rows B and F.
- 25 18. HPV types 16, 45 and 58 will be detected in rows C and G.
19. HPV types 31, 33, 35 and the internal control will be detected in rows D and H.

The PCR program should be:

Polymerase activation and sample denaturation: 10 mins at 95°C.

30 Followed by 45 cycles of:

Denaturation 15 seconds at 95°C

Annealing 30 seconds at 45°C

Extension 60 seconds at 60°C.

The realtime PCR machine should be programmed to collect fluorescence data in four colour channels optimized for detection of the fluorescent dyes FAM, VIC,
5 NED and ROX.

The above Examples are only meant as illustrating the invention. It must be understood that a person skilled in the art can modify the method, kit, probes and primers herein described without deviating from the concept and scope of this
10 invention set forth in the claims.

15

20

25

Patent claims

5

1. A method for detecting and/or typing and/or quantitation of the presence of a human papillomavirus (HPV) type in a biological sample, the test comprising the following steps:

- (i) amplification of a nucleic acid fragment of HPV in the presence of a DNA polymerase, forward PCR primer with an oligonucleotide sequence of SEQ ID NO:1 or an analog thereof, reverse PCR primer with an oligonucleotide sequence selected from the group consisting of: SEQ ID NO:2, SEQ ID NO:3, SEQ ID NO:4, SEQ ID NO:5 or SEQ ID NO:6 or an analog thereof, a probe with an oligonucleotide sequence selected from the group consisting of SEQ ID NO:7, SEQ ID NO:8, SEQ ID NO:9, SEQ ID NO:10, SEQ ID NO:11, SEQ ID NO:12, SEQ ID NO:13, SEQ ID NO:14, SEQ ID NO:15, SEQ ID NO:16, SEQ ID NO:17, SEQ ID NO:18, SEQ ID NO:19, SEQ ID NO:20, SEQ ID NO:21, SEQ ID NO:22, SEQ ID NO:23, SEQ ID NO:25, SEQ ID NO:26, SEQ ID NO:27, SEQ ID NO:30, SEQ ID NO:31, SEQ ID NO:32, SEQ ID NO:34, SEQ ID NO:35, SEQ ID NO:37, SEQ ID NO:38, SEQ ID NO:39, SEQ ID NO:41, SEQ ID NO:44, SEQ ID NO:45, SEQ ID NO:46, SEQ ID NO:47, SEQ ID NO:48, SEQ ID NO:49, SEQ ID NO:50, SEQ ID NO:51, SEQ ID NO:52, SEQ ID NO:55, SEQ ID NO:56, SEQ ID NO:57, SEQ ID NO:58, SEQ ID NO:59, SEQ ID NO:61, SEQ ID NO:62, SEQ ID NO:63, SEQ ID NO:64, SEQ ID NO:65, SEQ ID NO:66, SEQ ID NO:68, SEQ ID NO:70, SEQ ID NO:71, SEQ ID NO:72, SEQ ID NO:73, SEQ ID NO:74, SEQ ID NO:75, SEQ ID NO:77, SEQ ID NO:78, SEQ ID NO:79 and SEQ ID NO:80 or an analog thereof, labelled with a fluorescent dye and a quencher molecule, and an internal control;
- (ii) detection of a change of fluorescence and determining HPV type.

30

2. A method according to claim 1, wherein the oligonucleotide specifically hybridizes to a HPV amplicon wherein the HPV types are selected from HPV6,11,
5 16, 18, 26, 31, 33, 35, 39, 40, 42, 43, 44, 45, 51, 52, 53, 54, 56, 59, 59, 61, 66, 68, 70, 72, 73 and 82.
3. A method according to claims 1-2, wherein said HPV types are selected from HPV16,18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 and 66.
10
4. A method according to claims 1-3, wherein said probes are selected from the group consisting of: SEQ ID NO:7, SEQ ID NO:8, SEQ ID NO:9, SEQ ID NO:10, SEQ ID NO:11, SEQ ID NO:12, SEQ ID NO:13, SEQ ID NO:14, SEQ ID NO:15, SEQ ID NO:16, SEQ ID NO:17, SEQ ID NO:18, SEQ ID NO:19 and SEQ
15 ID NO:20 or an analog thereof.
5. A method according to claims 1-4, wherein two or more parallel reactions are run according to claim 1 for simultaneously detection of said HPV types.
- 20 6. A method according to claims 1-5, for detecting and/or typing and/or quantitation simultaneously the presence of HPV 18, 52, 59, 39, 51, 56, 66, 16, 45, 58, 31, 33 and 35 in a biological sample, the method comprising the following four parallel reactions:
- 25 Reaction 1- a) amplification of the nucleic acid fragments of HPV 18, 52 and 59 in the presence of
- (i) a nucleic acid polymerase
 - (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
 - 30 (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 2 or an analog thereof,

(iv) a probe with an oligonucleotide sequence of SEQ ID NO: 7,8 and 9 or an analog thereof, labelled with a fluorophore and a quencher molecule and

5 b) detection of a change of fluorescence the wavelength thereof determining the HPV type 18, 52 and 59 ;

Reaction 2- a) amplification of the nucleic acid fragments of HPV 39,51,56 and 66 in the presence of

- 10 (i) a nucleic acid polymerase
- (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
- (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 3 or an analog thereof,
- 15 (iv) a probe with an oligonucleotide sequence of SEQ ID NO: 10, 11, 12 and 13 or an analog thereof, labelled with a fluorophore and a quencher molecule and

 b) detection of a change of fluorescence the wavelength thereof
20 determining the HPV type 39, 51, 56 and 66;

Reaction 3- a) amplification of the nucleic acid fragments of HPV 16,45 and 58 in the presence of

- 25 (i) a nucleic acid polymerase
- (ii) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
- (iii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 4 or an analog thereof,
- 30 (iv) a probe with an oligonucleotide sequence of SEQ ID NO: 14, 15 and 16 or an analog thereof, labelled with a fluorophore and a quencher molecule and

b) detection of a change of fluorescence the wavelength thereof determining the HPV type 16, 45 and 58;

Reaction 4- a) amplification of the nucleic acid fragments of HPV
31, 33, 35 and internal control in the presence of

- (i) a nucleic acid polymerase
- (i) a forward PCR primer with an oligonucleotide sequence of SEQ ID NO: 1 or an analog thereof,
- (ii) a reverse PCR primer with an oligonucleotide sequence of SEQ ID NO: 5 or an analog thereof,
- (iii) a probe with an oligonucleotide sequence of SEQ ID NO: 17, 18, 19 and 20 or an analog thereof, labelled with a fluorophore and a quencher molecule and
- (iv) internal control and

b) detection of a change of fluorescence the wavelength thereof determining the HPV type 31, 33, 35 and internal control.

7. A method according to claim 6, wherein said reactions are performed in parallel and in separate vials.

8. A method according to claims 1-7, wherein the four fluorophores are FAM, VIC, NED and ROX.

9. A method according to claims 1-7, wherein the internal control is HPV6.

10. A primer for HPV amplification selected from the group consisting of: SEQ ID NO: 1, SEQ ID NO:2, SEQ ID NO: 3, SEQ ID NO:4, SEQ ID NO: 5 and SEQ ID NO:6.

11. A primer according to claim 10, wherein said primer is especially adapted to real-time multiplex HPV amplification.

12. A primer according to claims 1-6, 10 and 11, for use in the detection and/or typing and/or quantitation of HPV types present in a biological sample.

5 13. A probe for HPV detection and/or typing and/or quantitation selected from the group consisting of: SEQ ID NO:7, SEQ ID NO:8, SEQ ID NO:9, SEQ ID NO:10, SEQ ID NO:11, SEQ ID NO:12, SEQ ID NO:13, SEQ ID NO:14, SEQ ID NO:15, SEQ ID NO:16, SEQ ID NO:17, SEQ ID NO:18, SEQ ID NO:19, SEQ ID NO:20, SEQ ID NO:21, SEQ ID NO:22, SEQ ID NO:23, SEQ ID NO:25, SEQ ID
10 NO:26, SEQ ID NO:27, SEQ ID NO:30, SEQ ID NO:31, SEQ ID NO:32, SEQ ID NO:34, SEQ ID NO:35, SEQ ID NO:37, SEQ ID NO:38, SEQ ID NO:39, SEQ ID NO:41, SEQ ID NO:44, SEQ ID NO:45, SEQ ID NO:46, SEQ ID NO:47, SEQ ID NO:48, SEQ ID NO:49, SEQ ID NO:50, SEQ ID NO:51, SEQ ID NO:52, SEQ ID NO:55, SEQ ID NO:56, SEQ ID NO:57, SEQ ID NO:58, SEQ ID NO:59, SEQ ID
15 NO:61, SEQ ID NO:62, SEQ ID NO:63, SEQ ID NO:64, SEQ ID NO:65, SEQ ID NO:66, SEQ ID NO:68, SEQ ID NO:70, SEQ ID NO:71, SEQ ID NO:72, SEQ ID NO:73, SEQ ID NO:74, SEQ ID NO:75, SEQ ID NO:77, SEQ ID NO:78, SEQ ID NO:79 and SEQ ID NO:80 or an analog thereof.

20 14. A probe for HPV detection and/or typing and/or quantitation according to claim 13 selected from the group consisting of: SEQ ID NO:7, SEQ ID NO:8, SEQ ID NO:9, SEQ ID NO:10, SEQ ID NO:11, SEQ ID NO:12, SEQ ID NO:13, SEQ ID NO:14, SEQ ID NO:15, SEQ ID NO:16, SEQ ID NO:17, SEQ ID NO:18, SEQ ID NO:19 and SEQ ID NO:20.

25 15. A probe according to claims 1-6 and 13, for use in the detection and/or typing and/or quantitation of HPV present in a biological sample.

30 16. The use of the primers according to claims 1 and 10 for detecting and/or typing and/or quantitation of the presence of HPV.

17. The use of the probes according to claims 1 and 13 for detecting and/or typing and/or quantitation of the presence of HPV.

18. A diagnostic kit for detection and/or typing and/or quantitation of HPV in a biological sample, comprising forward and reverse primers, probes for HPV6,11, 16, 18, 26, 31, 33, 35, 39, 40, 42, 43, 44, 45, 51, 52, 53, 54, 56, 59, 59, 61, 66, 68,
5 70, 72, 73 and 82, internal positive control and DNA polymerase.
19. A diagnostic kit according to claim 18, wherein said primers are according to claim10 and said probes are according to claim 13.
- 10 20. A diagnostic kit according to claims 18 and 19, wherein said probes are for HPV16,18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 and 66.