

Optical Channel: FAM (Flu A)





**B** Optical Channel: HEX (SwH1)

FIGURE 1

C Optical Channel: TexRD (Flu B)**D Optical Channel: CY5 (H3N2)****FIGURE 1 (CONTINUED)**




FIGURE 2

SEQUENCE LISTING

<110> Agency for Science, Technology and Research

<120> INFLUENZA DETECTION METHOD AND KIT THEREFOR

<130> FP5107

<150> 200904777-0

<151> 2009-07-13

<160> 19

<170> PatentIn version 3.4

<210> 1

<211> 22

<212> DNA

<213> Artificial

<220>

<223> Forward Primer of Flu A (490:A145F)

<400> 1

ttctaaccga ggtcgaaacg ta

22

<210> 2

<211> 22

<212> DNA

<213> Artificial

<220>

<223> Reverse Primer of Flu A (488:A145R)

<400> 2

attggctcttg tcttttagcca yt

22

<210> 3

<211> 20

<212> DNA

<213> Artificial

<220>

<223> Probe of Flu A (500:A145U)

<400> 3

cccctcaaag ccgagatcgc

20

<210> 4

<211> 21

<212> DNA

<213> Artificial

<220>

<223> Forward Primer of Set A of Swine-Origin H1N1 (470:H1s138F)

<400> 4

ggctgctttg aattttacca c

21

<210> 5

<211> 20

<212> DNA

<213> Artificial

<220>

<223> Reverse Primer of Set A of Swine-Origin H1N1(471:H1s138R)
 <400> 5
 tgttgattcc agctttaccc 20

 <210> 6
 <211> 23
 <212> DNA
 <213> Artificial

 <220>
 <223> Probe of Set A of Swine-Origin H1N1 (475:138P)
 <400> 6
 gcgataacac gtgcatggaa agt 23

 <210> 7
 <211> 18
 <212> DNA
 <213> Artificial

 <220>
 <223> Forward Primer of Set B of Swine-Origin H1N1 (476:H1s142F)
 <400> 7
 aactttggac taccacga 18

 <210> 8
 <211> 19
 <212> DNA
 <213> Artificial

 <220>
 <223> Reverse Primer of Set B of Swine-Origin H1N1 (477:H1s142R)
 <400> 8
 ttgacactt tccatgcac 19

 <210> 9
 <211> 23
 <212> DNA
 <213> Artificial

 <220>
 <223> Probe of Set B of Swine-Origin H1N1 (478:142U)
 <400> 9
 acggctgctt tgaattttac cac 23

 <210> 10
 <211> 21
 <212> DNA
 <213> Artificial

 <220>
 <223> Forward Primer of Seasonal H1N1 (468:H1h218F)
 <400> 10
 cctgttatgc acatttacag c 21

 <210> 11
 <211> 19
 <212> DNA

<213> Artificial
 <220>
 <223> Reverse Primer of Seasonal H1N1 (469:H1h218R)
 <400> 11
 ctaagatcca cccggcaac 19
 <210> 12
 <211> 21
 <212> DNA
 <213> Artificial
 <220>
 <223> Probe of Seasonal H1N1 (72:H1h218P)
 <400> 12
 ctaacaactc gaccgacact g 21
 <210> 13
 <211> 20
 <212> DNA
 <213> Artificial
 <220>
 <223> Forward Primer of Seasonal H3N2 (482:H3,116F)
 <400> 13
 aatgacaaca gcacggcaac 20
 <210> 14
 <211> 20
 <212> DNA
 <213> Artificial
 <220>
 <223> Reverse Primer of Seasonal H3N2 (483:H3,116R)
 <400> 14
 tgaaccagct cagtagcatt 20
 <210> 15
 <211> 23
 <212> DNA
 <213> Artificial
 <220>
 <223> Probe of Seasonal H3N2 (484:H3,116U)
 <400> 15
 ccatgcagta ccaaacggaa cga 23
 <210> 16
 <211> 20
 <212> DNA
 <213> Artificial
 <220>
 <223> Forward Primer of Flu B (494:B188F)
 <400> 16
 ttggagaca caattgccta 20

<210> 17
<211> 18
<212> DNA
<213> Artificial

<220>
<223> Reverse Primer 1 of Flu B (495:B188R1)

<400> 17
aagcakatag aggcacca

18

<210> 18
<211> 18
<212> DNA
<213> Artificial

<220>
<223> Reverse Primer 2 of Flu B (496:B188R2)

<400> 18
aagcagatag aggcgcca

18

<210> 19
<211> 21
<212> DNA
<213> Artificial

<220>
<223> Probe of Flu B (497:B188L)

<400> 19
tcaaattctt tcccaccraa c

21