

Sequence Listing

<110> Chang Gung Memorial Hospital, Linkou

<120> METHODS FOR ASSESSING THE RISK OF DEVELOPING SEVERE CUTANEOUS ADVERSE DRUG REACTIONS INDUCED BY DISEASE-MODIFYING ANTIRHEUMATIC DRUGS, DETECTION KIT THEREOF AND USES THEREOF

<130> CGMH-LK-19002-PCT

<160> 16

<170> PatentIn version 3.5

<210> 1

<211> 13

<212> DNA

<213> Artificial Sequence

<220>

<223> HLA-B*1502 forward primer

<400> 1

atggcgcccc ggg

13

<210> 2

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> HLA-B*1502 reverse primer

<400> 2

tagtagccgc gcaggttcc

19

<210> 3

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> HLA-B*1502 probe 1

<400> 3

aacacacaga tctacaagg

19

<210> 4

<211> 19

<212> DNA

<213> Artificial Sequence

<220>
 <223> HLA-B*1502 probe 2

 <400> 4
 aacacacaga tctccaaga 19

<210> 5
 <211> 17
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> HLA-B*3802 forward primer

 <400> 5
 gccgcgagtc cgagaga 17

<210> 6
 <211> 20
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> HLA-B*3802 reverse primer

 <400> 6
 gtgcgcaggt tctctcggta 20

<210> 7
 <211> 15
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> HLA-B*3802 probe 1

 <400> 7
 ccggagtatt gggac 15

<210> 8
 <211> 15
 <212> DNA
 <213> Artificial Sequence

<220>
 <223> HLA-B*3802 probe 2

 <400> 8
 ccggaatatt gggac 15

<210> 9
 <211> 17
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> HLA-B*1301 forward primer

 <400> 9
 agccccgctt catcacc 17

<210> 10
 <211> 20
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> HLA-B*1301 reverse primer

 <400> 10
 tccttgccgt cgtaggctaa 20

<210> 11
 <211> 17
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> HLA-B*1301 probe 1

 <400> 11
 cacatcatcc agaggat 17

<210> 12
 <211> 16
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> HLA-B*1301 probe 2

 <400> 12
 acacttggca gacgat 16

<210> 13
 <211> 19
 <212> DNA
 <213> Artificial Sequence

 <220>
 <223> HLA-B*3901 forward primer

<400> 13
gcgagtccga gagaggagc 19

<210> 14
<211> 19
<212> DNA
<213> Artificial Sequence

<220>
<223> HLA-B*3901 reverse primer

<400> 14
tagtagccgc gcaggttcc 19

<210> 15
<211> 17
<212> DNA
<213> Artificial Sequence

<220>
<223> HLA-B*3901 probe 1

<400> 15
tccaattcac agactga 17

<210> 16
<211> 15
<212> DNA
<213> Artificial Sequence

<220>
<223> HLA-B*3901 probe 2

<400> 16
caacacacag actga 15