

-1-

DECONVOLUTION AND DETECTION OF RARE DNA IN PLASMA
CROSS REFERENCE TO RELATED APPLICATION

[0001] This application claims the priority benefit of U.S. Provisional Patent Application No. 62/438,401, filed December 21, 2016, the entire contents of which is
5 incorporated herein by reference.

STATEMENT REGARDING FEDERALLY SPONSORED R&D

[0002] This invention was made with government support under Grant Number R01GM097253 awarded by the National Institutes of Health. The government has certain rights in the invention.

10 BACKGROUND OF THE INVENTION

Field of the Invention

[0003] Some embodiments described herein relate to compositions and methods for detecting a target nucleic acid in a sample. For example, some embodiments may be used for non-invasive detection of tumors or organ damage. Other embodiments may be used, for
15 example, for detection of fetal aneuploidy at a very early stage of pregnancy.

Description of Related Art

[0004] CpG methylation in mammalian genomes is a relatively stable epigenetic modification, which can be transmitted across cell division (Wigler et al. (1981)) through DNMT1, and dynamically established, or removed by DNMT3 A/B and TET proteins. Due to
20 the processivity of some of these enzymes, physically adjacent CpG sites on the same DNA molecules can share similar methylation status, although discordant CpG methylation has also been observed, especially in cancer cells. The theoretical framework of linkage disequilibrium (Slatkin (2008)), which was developed to model the coordinated segregation of adjacent genetic variants on human chromosomes among human populations, can be applied to the
25 analysis of CpG co-methylation in cell populations. A number of studies related to the concepts of methylation haplotypes, epi-alleles, or epi-haplotypes have been reported, albeit at small numbers of genomic regions or limited numbers of cell/tissue types. Recent data production efforts, especially by large consortia such as the NIH RoadMap Epigenomics project (Bernstein et al. (2010)) and the EU Blueprint Epigenome project (Jones et al. (2005)) have produced a

-2-

large number of whole-genome, base-resolution bisulfite sequencing data sets for many tissue and cell types. These public data sets, in combination with additional WGBS data generated, allowed full-genome characterization of local coupled CpG methylation across the largest set of human tissue types available to date to be performed, and annotate these blocks of co-methylated CpGs as a distinct set of genomic features.

[0005] DNA methylation is cell-type specific, and the pattern can be harnessed for deconvoluting the relative cell composition of heterogeneous samples, such as different white blood cells in whole blood (Houseman et al. (2016)), fetal components in maternal cell-free DNA (Sun et al. (2015)), or circulating tumor DNA in plasma (Sun et al. (2015)). Most of these recent efforts rely on the methylation level of individual CpG sites, and are fundamentally limited by the technical noise and sensitivity in measuring single CpG methylation. Very recently, Lehmann-Werman et al. demonstrated a superior sensitivity with multi-CpG haplotypes in detecting tissue-specific signatures in circulating DNA (Lehmann-Werman et al. (2016)). The markers in that study were discovered from Infinium 450k methylation array data, which represent only a very limited fraction of the human genome.

SUMMARY OF THE INVENTION

[0006] Some embodiments described herein provide accurate tissue-of-origin mapping based on comparing the patterns and abundance of methylation haplotypes against a reference set of human reference tissues and provide accurate quantitative estimation of the cancer DNA fraction.

[0007] Some embodiments are described in the following numbered paragraphs:

1. A method for detecting the presence of one or more nucleic acids indicative of a health condition, tissue of origin, germ layer of origin or organ of origin in a mixture of nucleic acids comprising:
 - performing methylation analysis on a sample comprising a plurality of nucleic acids; and
 - determining whether said sample includes a plurality of methylation haplotype blocks indicative of the presence one or more nucleic acids indicative of a health condition, tissue of origin, germ layer of origin, organ of origin or any combination thereof, wherein said methylation haplotype blocks comprise a plurality of methylation sites for which the methylation status is coordinated.

-3-

2. The method of Paragraph 1, wherein said methylation analysis is performed on cell-free DNA.

3. The method of any one of Paragraphs 1 and 2, wherein said methylation analysis is performed on cell-free DNA in a blood sample.

5 4. The method of any one of Paragraphs 1-3 wherein said plurality of methylation haplotype blocks comprises at least 2, at least 3, at least 4, at least 5, at least 10, at least 20, at least 40, at least 50, at least 100, at least 200, at least 300, at least 400, at least 500 or more than 500 methylation haplotype blocks.

10 5. The method of any one of Paragraphs 1-4, wherein said health condition is a tumor.

6. The method of Paragraph 5, further comprising determining whether said sample includes a plurality of methylation haplotype blocks indicative of the presence of one or more nucleic acids indicative of a normal tissue or normal organ corresponding to the tissue or organ of origin of said tumor.

15 7. The method of any one of Paragraphs 1-6, wherein said health condition is fetal aneuploidy.

8. The method of any one of Paragraphs 1-7, wherein said sample is a blood sample.

20 9. The method of any one of Paragraphs 1-8, further comprising quantitating the level of said one or more nucleic acids indicative of a health condition, tissue of origin, organ of origin or any combination thereof in said sample.

25 10. The method of any one of Paragraphs 1-9, wherein said methylation analysis is performed using a technique selected from the group consisting of bisulfite methylation analysis, reduced representation bisulfite sequencing, WGBS, BSPP, micro-droplet PCR, selector probe based methods, and MeDiP.

30 11. The method of any one of Paragraphs 1-10, further comprising determining a methylated haplotype load or unmethylated haplotype load for each methylation haplotype block, wherein said methylation haplotype load comprises the normalized fraction of methylated haplotypes at different lengths and the unmethylated haplotype load comprises the normalized fraction of unmethylated haplotypes at different lengths.

-4-

12. The method of any one of Paragraphs 1-11, wherein said methylation haplotype blocks have an average size of 95bp.

13. The method of any one of Paragraphs 1-12, wherein said methylation haplotype blocks have a minimum of 3 CpGs per block.

5 14. The method of any one of Paragraphs 1-13 further comprising quantifying the level of said plurality of methylation haplotype blocks indicative of the presence of one or more nucleic acids indicative of a health condition, tissue of origin, germ layer of origin, organ of origin or any combination thereof in said sample.

10 15. The method of Paragraph 6 further comprising quantifying the level of said plurality of methylation haplotype blocks indicative of the presence of a tumor in said sample and quantifying the level of said plurality of methylation haplotype blocks indicative of the presence of one or more nucleic acids indicative of a normal tissue or normal organ corresponding to the tissue or organ of origin of said tumor in said sample.

15 16. A method of identifying methylation haplotype blocks comprising:
determining methylation haplotypes in a plurality of nucleic acid segments;
combining the methylation haplotypes and calculating methylation linkage disequilibrium on the combined methylation haplotypes; and
partitioning each segment into a plurality of methylation haplotype blocks,
wherein said methylation haplotype blocks comprise a plurality of methylation sites for
20 which the methylation status is coordinated.

17. The method of Paragraph 16, wherein said methylation haplotypes are determined for a portion of a genome.

18. The method of any one of Paragraphs 16 and 17, wherein said methylation haplotypes are determined across a whole genome.

25 19. The method of any one of Paragraphs 16-18, wherein said methylation haplotype blocks are defined as the genomic region in which the r^2 value of two adjacent CpG sites is no less than 0.5.

20 20. The method of any one of Paragraphs 16-19, wherein said methylation haplotype blocks are identified in nucleic acids from a tumor tissue.

30 21. The method of any one of Paragraphs 16-20, wherein said methylation haplotype blocks are identified in nucleic acids from a known type of tissue.

-5-

22. The method of any one of Paragraphs 16-21, wherein said methylation haplotype blocks are identified in nucleic acids from a fetus.

23. The method of any one of Paragraphs 16-22, wherein said methylation haplotype blocks are identified in nucleic acids from an embryonic stem cell.

5 24. The method of any one of Paragraphs 16-23, wherein said methylation haplotype blocks are identified in nucleic acids from a known germ layer.

25. The method of any one of Paragraphs 16-24, wherein said methylation haplotype blocks have an average size of 95bp.

10 26. The method of any one of Paragraphs 16-25, wherein said methylation haplotype blocks have a minimum of 3 CpGs per block.

27. The method of any one of Paragraphs 16-26, wherein said methylation haplotype blocks are identified in nucleic acid regions from a Whole Genome Bisulfite Sequencing analysis.

15 28. The method of any one of Paragraphs 16-27, wherein said methylation haplotype blocks are identified in data sets from methylation analysis of ENCODE cell lines or tissue samples.

29. The method of any one of Paragraphs 16-28, wherein said methylation haplotype blocks are identified in data sets from methylation analysis of Infinium HumanMethylation450 BeadChip (HM450K).

20 30. The method of any one of Paragraphs 16-29, further comprising calculating the pairwise correlation coefficient of adjacent CpG methylation levels across different sample sets for block partitioning.

31. The method of any one of Paragraphs 16-30, further comprising determining methylation haplotype load for each methylation haplotype block, wherein
25 said methylation haplotype load comprises the normalized fraction of methylated haplotypes at different lengths.

32. The method of any one of Paragraphs 16-30, further comprising determining unmethylated haplotype load for each methylation haplotype block, wherein said unmethylated haplotype load comprises the normalized fraction of
30 unmethylated haplotypes at different lengths.

BRIEF DESCRIPTION OF THE DRAWINGS

[0008] Figures 1(a)-1(e). Identification and characterization of human methylation haplotype blocks (MHBs). Fig. 1(a) Schematic overview of data collection, generation, and analysis. Fig. 1(b) An example of MHB at the promoter of the gene APC. Fig. 1(c) Smooth scatterplots of methylation linkage disequilibrium decay of adjacent CpG sites over larger
5 distances. 500,000 adjacent CpG loci in MHB regions were randomly sampled and their corresponding r^2 values within different sets of cell types were plotted to demonstrate the differential decay characteristics. The dotted lines at high linkage ($r^2 = 0.9$) is where stem and progenitor cells (10 WGBS samples), normal adult tissue cells (49 WGBS samples), and primary tumor (6 WGBS samples) cells were found to have 94.8%, 91.2% and 87.8% of CpGs
10 respectively. Fig. 1(d) Co-localization analysis of MHBs with known genomic features. Genome distribution (left) and CpG-island relationships (right). Fig. 1(e). Enrichment of MHBs in known genomic features. Bootstrap random sampling regions with same size for 10,000 times to estimate empirical statistical significance and enrichment factor (fold-change).

[0009] Figure 2. Comparison of methylation haplotype load with four metrics used in
15 the literature. Five patterns of methylation haplotype combinations are used to illustrate the difference between methylation frequency, methylation entropy, epi-polymorphism and methylation haplotype load. Methylation haplotype load can discriminate all the five patterns while other metrics cannot.

[0010] Figures 3(a)-3(b). Heatmaps of tissue specific MHBs. Fig. 3(a) Heatmap of
20 MHL values for tissue specific MHBs selected by MHL. Fig. 3(b) Heatmap of uMHL values for tissue specific MHBs selected by uMHL.

[0011] Figure 4. Comparison of signal to noise ratio between average methylation frequency (AMF) and methylation haplotype load (MHL) metrics at tissue specific differentially methylated loci (Lokk et al. (2014)). For most loci, the MHL metric had a much
25 higher signal to noise ratio than the AMF metric even though these loci were selected using AMF for another data set.

[0012] Figure 5. Heatmap of the MHL in plasma samples in test sets and primary cancer tissues for colon cancer markers (left) and lung cancer markers (right). Cancer markers were identified by comparison of the respective cancer tissues against background (normal
30 plasma set aside for features selection).

[0013] Figure 6. Boxplots of average MHL in plasma samples in test sets and cancer tissues for colon cancer markers (left) and lung cancer markers (right). The cancer markers

-7-

were significantly (Two Sample one-sided T-Test) different between the test sets of normal plasma (NP) and cancer plasma (CCP = colon cancer plasma, LCP = lung cancer plasma). The marker regions were highly methylated in tumor tissues (CCT = colon cancer tissue, LCT = lung cancer tissue).

5 **[0014]** Figure 7. Simulated standard curves for estimation of tumor fraction. The mean for each cancer fraction from 20 simulations plotted with standard deviations as error bars. A fitted linear model on the standard curve for colon cancer had an adjusted r^2 of 0.9621 and the fitted linear model for lung cancer had an adjusted r^2 of 0.9573.

10 **[0015]** Figure 8. The tumor fraction between test sets of cancer and normal plasma samples. Cancer plasma (n = 30 colon cancer, n= 29 lung cancer) had significantly elevated tumor fractions (Two Sample t-test with unequal variance) when compared with normal plasma (n=23) using either colon cancer markers (left) or lung cancer markers (right).

15 **[0016]** Figures 9(a)-9(c). Results of MARS- based binary classification for test plasma data sets. The set of tissue-specific markers derived from an independent set of normal tissue data was able to segregate each class of plasma samples from the other classes in a binary fashion. Fig. 9(a) Normal vs other. Fig. 9(b) Colon vs other. Fig. 9(c) Lung vs other.

20 **[0017]** Figure 10. Confusion matrix of MARS-based feature selection and PLSDA prediction of tissue-of-origin in the test plasma data set. The set of lung and colon tissue specific markers identified from an independent data set was able to correctly classify 74% of the test plasma samples.

25 **[0018]** Figures 11(a)-11(c). Characteristics of MHB in human genome. Fig. 11(a) Distribution of MHB sizes. Fig. 11(b) Distribution of CpG density (CpGs/bp) in MHB regions. Fig. 11(c) Co-localization of MHB with known genomics features breaking down based on CpG density. MHBs were placed in quartiles where the CpGs/bp of MHBs within quartiles were as follows (0, 0.46), (0.046, 0.096), (0.096, 0.155), and (0.155, 0.6). The 1st quartile (MHBs with lowest CpG density) was mostly CGI shelves or shores, and was enriched for LAD, LOCK, and enhancers.

30 **[0019]** Figures 12(a)-12(b). Validation of MHB with Illumina 450k methylation array and RRBS data. Fig. 12(a) Pearson correlation coefficient (r^2) versus absolute LD r^2 . Fig. 12(b) The Pearson correlation coefficient r^2 in RRBS and HM450K were significantly higher in overlapped MHBs with WGBS compared with the MHBs without overlapping with

-8-

WGBS MHBs. IN: denotes RRBS or HM450K regions within MHB. OUT: denotes RRBS or HM450K regions beyond MHB regions.

[0020] Figure 13. The distribution of incidence of cancer-associated high-methylated haplotypes (caHMH) in colon cancer (CRC) and lung cancer (LC) plasma samples. Y-axis denotes the frequency of caHMH and x-axis denotes the incidence (number of samples) of the caHMH in cancer plasmas. A majority of caHMH are patient specific.

DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENT

[0021] Unless defined otherwise, all technical and scientific terms used herein have the same meaning as commonly understood by one of ordinary skill in the art pertinent to the methods and systems described. As used herein, the following terms and phrases have the meanings ascribed to them unless specified otherwise. All patents, applications, published applications and other publications referred to herein are incorporated by reference in their entirety.

[0022] “A,” “an,” and “the” include plural referents, unless the context clearly indicates otherwise. For example, “a nucleic acid” as used herein is understood to represent one or more nucleic acids. As such, the terms “a” (or “an”), “one or more,” and “at least one” can be used interchangeably herein.

[0023] “About” is used herein to mean approximately, roughly, around, or in the region of. When the term “about” is used in conjunction with a numerical range, it modifies that range by extending the boundaries above and below the values set forth. In general, the term “about” is used herein to modify a numerical value above and below the stated value by a deviation of $\pm 10\%$ and preferably $\pm 5\%$.

[0024] “Comprise,” “comprises,” “comprising,” “include,” “includes,” and “including” are interchangeable and not intended to be limiting. Furthermore, where the description of one or more embodiments uses the term “comprising,” those skilled in the art would understand that, in some specific instances, the embodiment or embodiments can be alternatively described using the language “consisting essentially of” and/or “consisting of.”

[0025] The following abbreviations are used throughout the application: Methylation Haplotype Block (MHB); Methylation Haplotype Load (MHL); Unmethylated Haplotype Load (uMHL); Group Specific Index (GSI); Circulating cell-free DNA (cfDNA); Reduced Representation Bisulfite Sequencing (RRBS); single-cell Reduced Representation Bisulfite

Sequencing (scRRBS); Whole Genome Bisulfite Sequencing (WGBS); The Cancer Genome Atlas project (TCGA); The Encyclopedia of DNA Elements (ENCODE); Gene Expression Omnibus (GEO); Lung Cancer (LC); Colorectal or colon cancer (CRC); cancer associated High Methylation Haplotype (caHMH); tissue-specific Methylation Haplotype Block regions
5 (tsMHB); Colorectal or colon cancer tissue (CCT); colorectal or colon cancer plasma (CCP); lung cancer tissue (LCT); lung cancer plasma (LCP); normal plasma (NP).

[0026] "Amplification" refers to any known procedure for obtaining multiple copies of a target nucleic acid or its complement, or fragments thereof. The multiple copies may be referred to as amplicons or amplification products. Amplification, in the context of fragments,
10 refers to production of an amplified nucleic acid that contains less than the complete target nucleic acid or its complement, e.g., produced by using an amplification oligonucleotide that hybridizes to, and initiates polymerization from, an internal position of the target nucleic acid. Known amplification methods include, for example, replicase-mediated amplification, polymerase chain reaction (PCR), reverse transcription polymerase chain reaction (RT-PCR),
15 ligase chain reaction (LCR), strand-displacement amplification (SDA), and transcription-mediated or transcription-associated amplification.

[0027] "Complementary" means that a contiguous nucleic acid base sequence is capable of hybridizing to another base sequence by standard base pairing (hydrogen bonding) between a series of complementary bases. Complementary sequences may be completely
20 complementary (*i.e.* no mismatches in the nucleic acid duplex) at each position in an oligomer sequence relative to its target sequence by using standard base pairing (e.g., G:C, A:T or A:U pairing) or sequences may contain one or more positions that are not complementary by base pairing (e.g., there exists at least one mismatch or unmatched base in the nucleic acid duplex), but such sequences are sufficiently complementary because the entire oligomer sequence is
25 capable of specifically hybridizing with its target sequence in appropriate hybridization conditions (*i.e.* partially complementary). Contiguous bases in an oligomer are typically at least 80%, preferably at least 90%, and more preferably completely complementary to the intended target sequence.

[0028] "Configured to" denotes an actual arrangement of a nucleic acid sequence
30 configuration of a referenced oligonucleotide. For example, a primer that is configured to generate a specified amplicon from a target nucleic acid has a nucleic acid sequence that hybridizes to the target nucleic acid or a region thereof and can be used in an amplification

reaction to generate the amplicon. Also as an example, an oligonucleotide that is configured to specifically hybridize to a target nucleic acid or a region thereof has a nucleic acid sequence that specifically hybridizes to the referenced sequence under stringent hybridization conditions.

[0029] “Configured to specifically hybridize to” means that an oligonucleotide is designed to have a nucleic acid sequence that can hybridize with a target nucleic acid or region thereof. The oligonucleotide is designed to function as a component of an assay for amplification and detection of a target nucleic acid (or a region thereof) in a sample, and therefore is designed to hybridize with a target nucleic acid (or a region thereof) in the presence of other nucleic acids that may be found in testing samples.

10 [0030] “Fragment” refers to a piece of contiguous nucleic acid that contains fewer nucleotides than the complete nucleic acid.

[0031] “Hybridization” or “annealing” refer to the base-pairing interaction of one nucleic acid with another nucleic acid (typically an antiparallel nucleic acid) that results in formation of a duplex or other higher-ordered structure (*i.e.* a hybridization complex). The primary interaction between the antiparallel nucleic acid molecules is typically base specific, e.g., A/T and G/C. It is not a requirement that two nucleic acids have 100% complementarity over their full length to achieve hybridization. Nucleic acids hybridize due to a variety of well characterized physio-chemical forces, such as hydrogen bonding, solvent exclusion, base stacking and the like. An extensive guide to the hybridization of nucleic acids is found in
15 Tijssen (1993) Laboratory Techniques in Biochemistry and Molecular Biology--Hybridization with Nucleic Acid Probes part I chapter 2, “Overview of principles of hybridization and the strategy of nucleic acid probe assays,” (Elsevier, New York), as well as in Ausubel (Ed.) Current Protocols in Molecular Biology, Volumes I, II, and III, 1997, which is incorporated by reference.

25 [0032] “Nucleic acid” or “nucleic acid molecule” refers to a multimeric compound comprising two or more covalently bonded nucleosides or nucleoside analogs having nitrogenous heterocyclic bases, or base analogs, where the nucleosides are linked together by phosphodiester bonds or other linkages to form a polynucleotide. Nucleic acids include RNA, DNA, or chimeric DNA-RNA polymers or oligonucleotides, and analogs thereof. A nucleic acid backbone can be made up of a variety of linkages, including one or more of sugar-phosphodiester linkages, peptide-nucleic acid bonds, phosphorothioate linkages, methylphosphonate linkages, or combinations thereof. Sugar moieties of the nucleic acid can
30

-11-

be ribose, deoxyribose, or similar compounds having known substitutions (e.g. 2'-methoxy substitutions and 2'-halide substitutions). Nitrogenous bases can be conventional bases (A, G, C, T, U) or analogs thereof (e.g., inosine, 5-methylisocytosine, isoguanine). A nucleic acid can comprise only conventional sugars, bases, and linkages as found in RNA and DNA, or can include conventional components and substitutions (e.g., conventional bases linked by a 2'-methoxy backbone, or a nucleic acid including a mixture of conventional bases and one or more base analogs). Nucleic acids can include "locked nucleic acids" (LNA), in which one or more nucleotide monomers have a bicyclic furanose unit locked in an RNA mimicking sugar conformation, which enhances hybridization affinity toward complementary sequences in single-stranded RNA (ssRNA), single-stranded DNA (ssDNA), or double-stranded DNA (dsDNA). Nucleic acids can include modified bases to alter the function or behavior of the nucleic acid (e.g., addition of a 3'-terminal dideoxynucleotide to block additional nucleotides from being added to the nucleic acid). Synthetic methods for making nucleic acids in vitro are well known in the art although nucleic acids can be purified from natural sources using routine techniques. Nucleic acids can be single-stranded or double-stranded.

[0033] "Primer" refers to an enzymatically extendable oligonucleotide, generally with a defined sequence that is designed to hybridize in an antiparallel manner with a complementary, primer-specific portion of a target nucleic acid. A primer can initiate the polymerization of nucleotides in a template-dependent manner to yield a nucleic acid that is complementary to the target nucleic acid when placed under suitable nucleic acid synthesis conditions (e.g. a primer annealed to a target can be extended in the presence of nucleotides and a DNA/RNA polymerase at a suitable temperature and pH). Suitable reaction conditions and reagents are known to those of ordinary skill in the art. A primer is typically single stranded for maximum efficiency in amplification, but may alternatively be double stranded. If double stranded, the primer is generally first treated to separate its strands before being used to prepare extension products. The primer generally is sufficiently long to prime the synthesis of extension products in the presence of the inducing agent (e.g. polymerase). Specific length and sequence will be dependent on the complexity of the required DNA or RNA targets, as well as on the conditions of primer use such as temperature and ionic strength. Preferably, the primer is about 5-100 nucleotides. Thus, a primer can be, e.g., 5, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95 or 100 nucleotides in length. A primer does not need to have 100%

-12-

complementarity with its template for primer elongation to occur; primers with less than 100% complementarity can be sufficient for hybridization and polymerase elongation to occur. A primer can be labeled if desired. The label used on a primer can be any suitable label, and can be detected by, for example, spectroscopic, photochemical, biochemical, immunochemical, chemical, or other detection means. A labeled primer therefore refers to an oligomer that hybridizes specifically to a target sequence in a nucleic acid, or in an amplified nucleic acid, under conditions that promote hybridization to allow selective detection of the target sequence.

[0034] "Sample preparation" refers to any steps or methods that prepare a sample for subsequent sequencing, amplification, and/or detection of target nucleic acids present in the sample. Sample preparation may include any known method of concentrating components, such as nucleic acids, from a larger sample volume. Sample preparation may include physical disruption and/or chemical lysis of cellular components to release intracellular components into a substantially aqueous or organic phase and removal of debris, such as by using filtration, centrifugation or adsorption. Sample preparation may include use of a nucleic acid oligonucleotide that selectively or non-specifically captures a target nucleic acid and separates it from other sample components.

[0035] "Sequencing" refers to any known procedure, method, or technology for determining the precise order of the nucleosides or nucleoside analogs of a target nucleic acid molecule, or its complement, or fragments thereof. Sequencing, in the context of fragments, refers to determining the precise order of nucleosides or nucleotides within a nucleic acid molecule that contains less bases than the complete target nucleic acid molecule e.g., determined by sequencing amplicons produced by using an amplification oligonucleotide that hybridizes to, and initiates polymerization from, an internal position of the target nucleic acid. Known sequencing methods include, for example, whole-genome sequencing as well as targeted sequencing wherein only subset of genes or regions of the genome are isolated and sequenced.

[0036] Some embodiments described herein allow a quantitative deconvolution of biological samples (for example, human plasma) that contain mixed DNA molecules, based on comparing the patterns and abundance of methylation haplotypes against a reference set of human tissues. Detection and quantification of low-abundant cancer DNA in plasma, and simultaneous mapping to the tissue of origin has been successfully demonstrated. With the high sensitivity and the capability of detecting rare species in heterogeneous samples. Some

embodiments described herein are, for example, suitable for non-invasive detection of tumor or other organ damage (diabetes, stroke etc.) in a subject's plasma, urine, stool, or cerebrospinal fluid, and detection of fetal aneuploidy in maternal blood at a very early stage of pregnancy.

5 **[0037]** Some embodiments described herein provide methods of creating a database of patient specific methylation haplotypes and unmethylated haplotypes. Some embodiments described herein provide methods of preparing data of comparative methylation haplotypes and unmethylated haplotypes. Some embodiments described herein include further steps for methods of treatment of diagnosed diseases or conditions associated with methylation haplotypes and unmethylated haplotypes.

10 **[0038]** Some embodiments described herein utilize the co-methylation status of multiple adjacent CpG sites in single DNA molecules rather than the methylation status of individual loci and or their local average. Some embodiments described herein suppress the stochastic noise from the low coverage methylation data and therefore have high sensitivity and specificity in both detecting circulating tumor DNA and mapping the tissue-of-origin of the tumor. Additionally, some embodiments described herein provide a set of MHL based biomarkers and a statistical approach to the detection.

15 **[0039]** For the quantitative estimation of the level of cancer DNA fractions, a strategy based on a tissue-specific MHL sampling technique which is distinctly different from existing approaches, such as quadratic programming, was developed. Due to the low level of cancer and fetal DNA fragments in the blood (plasma), some embodiments are applicable to even low coverage data based on methylation sequencing, such as RRBS and WGBS, whereas other methods would fail due to the difficulties in dealing with missing values.

25 **[0040]** For non-invasive cancer diagnosis, current mature commercial assays in the market use bronchial fluid (EPI PROLUNG BL Reflex assay), or can only be used for specific cancer type detection (EPI PROCOLON blood assay for colon cancer) without origin-tissue prediction. Some embodiments described herein provide accurate tissue-of-origin mapping based on comparing the patterns and abundance of methylation haplotypes against a reference set of human reference tissues and provide accurate quantitative estimation of the cancer DNA fraction.

30 **[0041]** Some embodiments described herein relate to a set of markers for cancer detection and tissue-of-origin mapping based on an exhaustive search across the whole genome for methylation haplotype blocks (MHB) that have a higher level of methylation haplotype load

(MHL) or higher level of unmethylated haplotype load (uMHL). Some embodiments described herein relate to a statistical approach that takes the information of these markers for the tissue-of-origin mapping and cancer detection. Compared with existing prediction models (such as tree-based: random forest), some embodiments described herein are suitable for sparse methylation sequencing data, which is typical in the clinical setting due to the low amounts of cell-free DNA that can be extracted from patients. Some embodiments described herein are also applicable to detection of fetal aneuploidy in maternal blood. Currently there is a commercial TISSUE OF ORIGIN test from Cancer Genetics Inc. which is an invasive assay that uses tissue biopsies from patients, as compared to non-invasive liquid biopsies that allow far less discomfort for the patients.

[0042] Human peripheral blood contains low levels of DNA molecules from other tissues or cell types, such as circulating cancer stem cells, cell-free DNA from apoptotic cancer cells in cancer patients, or fetal DNA in pregnant women. DNA methylation signals along with the DNA molecules are released into the blood simultaneously and these methylation signals are tissue-specific and can be applied to identify the tissue source of the DNA fragment. Methylation haplotypes provide high sensitivity detection of DNA molecules from DNA samples (WO2015/116837 and PCT/US2015/013562 incorporated herein by reference).

[0043] In some embodiments, DNA molecules are extracted from plasma samples, from for example cancer patients or healthy individuals, and the DNA methylation status of cell-free DNA molecules are assayed by bisulfite methylation sequencing (reduced representation bisulfite sequencing (RRBS), Meissner et al. (2005)). Note that alternative technologies, such as BSPP (Diep et al. (2012)), micro-droplet PCR (Komori et al. (2011)), Selector probes (Johansson et al. (2011)), or MeDiP (Papageorgiou et al. (2011)) can also potentially be used with some differences in the requirement of input materials and/or cost. Other methods for determining the methylation status of cell-free DNA may also be used.

[0044] In some embodiments, regardless of the specific sample preparation methods or sequencing platforms used, bisulfite sequencing reads (single-ends or paired-ends) are used as the input for the proposed analytical framework for detection of cancer and their tissue-of-origin. Methylation haplotypes and their abundance may be derived from the raw sequencing reads. Each haplotype represents the combination of binary methylation status (methylated or unmethylated) at multiple CpG sites of one sequencing read. A computational pipeline (implemented with Perl) was developed to deal with methylation haplotype and further derived

metrics, such as methylation haplotype (successive methylation allele combination), methylation haplotype block (linkage disequilibrium regions), highly methylated haplotype (HMH), and methylation/unmethylated haplotype load (MHL/uMHL) (weighted methylation status by haplotype length).

5 **[0045]** To identify the most informative marker for detecting cancer or tissue mapping, the term “methylation haplotype block” (MHB) was defined to describe genomic regions in which the methylation status of CpG pairs within are in linkage disequilibrium and could be taken as co-methylation regions. A set of Whole Genome Bisulfite Sequencing (WGBS) data was compiled, such as from stem cell, cancer cell, and normal adult tissues, in order to
10 exhaustively search the entire human genome. A total of 147,888 MHBs in the human genome were identified.

[0046] A metric was defined, called methylation haplotype load (MHL) for each MHB (as described in WO2015/116837), which is the normalized fraction of methylated haplotypes at different lengths to indicate the methylation status and complexity of the methylation
15 population. MHL is superior to other metrics (such as methylation level, methylation entropy (Xie, H. et al, (2011)), and epi-polymorphism (Landan et al, (2012)) for distinguishing different methylation patterns. In addition, MHL is bounded between 0 and 1, which allows for direct comparison of different regions across many data sets without normalization.

[0047] A metric was defined, called unmethylated haplotype load (uMHL) for each
20 MHB, which is the normalized fraction of unmethylated haplotypes at different lengths to indicate the methylation status and complexity for the unmethylated population.

[0048] A MHL matrix for a set of human tissues (61 published samples/data sets) was built encompassing all the 147,888 identified MHBs in the entire human genome. From this matrix, tissue-specific MHBs were identified. After including additional training data subsets
25 comprising plasma from cancer patients and healthy controls, sets of 154 MHBs for binary prediction of tissue-of-origin between colon, lung, and normal (no solid tissue) and 295 MHBs for cancer types classification (colon versus lung) non-invasively in the blood were identified (Table 5(a)-(b)). The invention demonstrated binary classification AUCs of 0.856, 0.725, and 0.751 in independent test sets from normal, colon cancer, and lung cancer plasma respectively
30 (Fig. 9(a)-9(c)). Distinguishing between colon cancer and lung cancer plasma achieved classification accuracies of 70% and 78% for test sets comprising colon cancer and lung cancer plasma respectively using cancer type classification (Fig. 10).

[0049] The disclosed method was also applicable to tracking metastasized secondary tumors. Significantly higher original tissue markers (46 and 79 for two metastasized cancer samples) in the metastasized cancers (primary colon to liver and primary breast to lung) were identified with MHL measurement and a maximum marker count approach to infer its primary cancer source.

[0050] Adjacent CpG sites in mammalian genomes tend to be co-methylated due to the processivity of enzymes responsible for adding or removing the methyl group. Yet discordant methylation patterns have also been observed, and found to be related to stochastic or uncoordinated molecular processes. The invention focused on a systematic search and investigation of regions in the human genome that exhibit highly coordinated methylation.

[0051] By examining the co-methylation patterns of multiple adjacent CpG sites, termed methylation haplotypes, in single bisulfite sequencing reads, a greedy-searching strategy was applied to define blocks of tightly coupled CpG sites, called methylation haplotype blocks (MHBs), based on 53 sets of whole genome bisulfite sequencing (WGBS) data, including 43 published sets from human adult tissues, ESC, and in vitro differentiated cell lines, as well as 10 sets from human adult tissues generated in this invention. The MHBs were further validated with 101 sets of RRBS ENCODE data, and 1,274 sets of Illumina450k methylation array data from TCGA tumor and normal samples. Globally, MHBs are enriched in, but only partially overlap with, several well-known genomic features, including CpG islands, promoters, enhancers and VMRs.

[0052] To perform quantitative analysis of the MHBs, a metric called methylation haplotype load (MHL) was computed, which covers both average methylation level and methylation complexity, and therefore is more informative than average methylation level or Shannon entropy. Using a feature selection strategy, a set of tissue-specific MHBs was identified that cluster by developmental germ-layers. Interestingly, examination of these MHBs revealed two distinct mechanisms for fate commitment during development: epigenetic silencing of pluripotent genes, such as NANOG, for mesoderm induction; and epigenetic induction (or de-suppression) of lineage-specific factors for ectoderm commitment.

[0053] Furthermore, to examine MHBs in cancers and explore its clinical utility, 162 sets of RRBS data from primary tumor tissues and matched plasma from patients with lung cancer, colon cancer, or pancreatic cancer were generated, as well as plasma controls from healthy individuals. Compared with normal tissues and stem cells, primary tumor tissues

-17-

exhibit a distinct methylation pattern within MHBs, related to locally disordered methylation recently discovered in chronic lymphocytic leukemia (CLL). Importantly, a subset of blocks (Table 4) was derived that can estimate the tumor content from circulating DNA in the plasma (Fig. 8). Finally, prediction of tissue-of-origin was performed on plasma using tissue specific methylation haplotype blocks (Table 5).

5 [0054] A manuscript entitled “Identification of methylation haplotype blocks aids in deconvolution of heterogeneous tissue samples and tumor tissue-of-origin mapping from plasma DNA” by Guo et al. (Nat Genet 2017 Apr;49(4):635-642), the disclosure of which is incorporated herein by reference in its entirety, including all figures, tables, supplementary
10 figures and supplementary tables therein or related thereto.

Examples

[0055] An exhaustive search of tissue-specific methylation haplotype blocks across the entire human genome was performed, and block-level metrics including methylation haplotype load (MHL) and unmethylated haplotype load (uMHL) were proposed, for a systematic
15 discovery of informative markers. Applying the disclosed method’s analytical framework and identified markers, accurate determination of tissue origin as well as estimation of tumor load in clinical plasma samples from patients of lung cancer (LC) and colorectal cancer (CRC) (Fig. 1(a)) was demonstrated.

[0056] Identification and characterization of methylation haplotype blocks. To
20 investigate the co-methylation status of adjacent CpG sites along single DNA molecules, the concept of genetic linkage disequilibrium was extended (Slatkin (2008); Shoemaker et al. (2010)) and the r^2 metric determined to quantify the degree of coupled CpG methylation among different DNA molecules of the same samples. CpG methylation status of multiple CpG sites in single- or paired-end Illumina sequencing reads were extracted to form methylation
25 haplotypes, and pairwise “linkage disequilibrium” of CpG methylation r^2 was calculated from the abundance of different methylation haplotypes (see Methods). The full human genome was partitioned into blocks of tightly coupled CpG methylation sites, called methylation haplotype blocks (MHBs) (Fig. 1(b)), using a r^2 cutoff of 0.5. Using slightly different cutoff values, such as 0.3 or 0.7, for partitioning of the human genome into genetic haplotype blocks, resulted in
30 only minor quantitative differences in the block size and number without affecting the global pattern (data not shown).

[0057] To characterize the global pattern and distribution of MHBs, starting with 51 sets of published Whole Genome Bisulfite Sequencing (WGBS) data from human primary tissues (Schultz et al. (2015); Heyn et al. (2012)), as well as the H1 human embryonic stem cells, *in vitro* derived progenitors (Xie, W. et al. (2013)), and human cancer cell lines (Blattler et al. (2014); Heyn et al. (2016)). An additional WGBS data set from 10 adult tissues from a single human donor was also generated and included. Across this set of 61 samples (>2000x combined genome coverage) a total of ~55 billion methylation haplotype informative reads covering 58.2% of autosomal CpGs were identified. The uncovered CpG sites were either in regions with low mappability, or CpG sparse regions where there are too few sites within Illumina read pairs for deriving informative haplotypes. 147,888 MHBs with an average size of 95bp and minimum 3 CpGs per block were identified, representing ~0.5% of the human genome that tends to be tightly co-regulated on the epigenetic status at the level of single DNA molecules (Fig. 11(a)-11(b)). The majority of CpG sites within the same MHBs are near perfectly coupled ($r^2 \sim 1.0$) regardless of the sample type. Methylation LD was found to extend further along the DNA in stem cells and progenitors, compared with normal adult tissue, both in the fraction of tightly coupled CpG pairs (94.8% versus 91.2%, $P\text{-value} < 2.6 \times 10^{-16}$), and the over-representation of partially coupled CpG pairs that are over 100 bp apart while the linkage was slightly decayed in primary cancer data sets (87.8%, mixture of CRC and LC), which was validated by another independent WGBS data set from kidney cancer (Chen et al. (2016)) (Fig. 1(c)). Gene Ontology analysis show cancer loss of linkage regions was significantly associated with number of cancer related pathways and functions. This is consistent with previous observations on a smaller BSPP data set comprising 2,020 CpG islands (Shoemaker et al. (2010)) for culture cell lines and another previous report (Shao et al. (2014)). Interestingly, in tumor samples, a reduction of perfectly coupled CpG pairs was observed, which could be related to the pattern of discordant methylation recently reported in variable methylation regions (VMR) (Landau et al. (2014); Hansen et al. (2011)).

[0058] While WGBS data allowed MHBs across the entire genome to unbiasedly be identified, the 61 sets of data did not represent the full diversity of human cell/tissue types. To validate the presence of MHBs in a wider range of human tissues and culture cells, 101 published reduced representation bisulfite sequencing (RRBS) data sets from ENCODE cell lines and tissue samples were examined, as well as 637 sets of Infinium HumanMethylation450 BeadChip (HM450K) data including 11 normal human tissues from the TCGA project. The

ENCODE RRBS data sets were generated with short (36bp) Illumina sequencing reads, greatly limiting the length of methylation haplotypes that can be identified. Similarly, Illumina methylation arrays only report average CpG methylation of all DNA molecules in a sample, preventing a methylation linkage disequilibrium analysis. Therefore, the invention calculated the pairwise correlation coefficient of adjacent CpG methylation levels across different sample sets for block partitioning. Note that the presence of such correlated methylation blocks is a necessary but not sufficient condition for MHBs (Fig. 12(a)). Nonetheless, the absence of correlated methylation blocks in these data sets would invalidate the pattern of MHBs. 23,517 and 2,212 correlated methylation blocks from ENCODE RRBS and TCGA HM450K array data, respectively, were identified, among which 8,920 and 1,258 have significant overlaps with WGBS-defined MHBs. Additionally, significantly higher correlation among the CpGs within the MHB regions compared to CpG loci outside MHBs in the HM450K and RRBS data sets was observed, further supporting the block-like organization of local CpG co-methylation across a wide variety of cells and tissues (Fig. 12(b)). Taken together, the MHBs identified represent a distinct class of genomic features where local CpG methylation is established or removed in a highly coordinated manner at the level of single DNA molecules, presumably due to the processive activities of the related enzymes coupled with the local density of CpG dinucleotides.

[0059] Co-localization of methylation haplotype blocks with known regulatory elements. The MHBs established by 61 sets of WGBS data appear to represent a distinct type of genomic feature that partially overlaps with multiple well-documented genomic elements (Fig. 1(d), Fig. 11). Among all the methylation blocks, 60,828 (41.1%) were located in intergenic regions and 87,060 (58.9%) regions in transcribed regions. These MHBs were significantly ($p\text{-value} < 10^{-6}$) enriched in enhancers (enrichment factor=7.6), super-enhancers (enrichment factor=2.3), promoter regions (enrichment factor=14.5), CpG islands (enrichment factor=70.4), and imprinted genes (enrichment factor=54.6). In addition, modest depletion in LAD was observed (Guelen et al. (2008)) and LOCK regions (Wen et al. (2009)) (46% and 37% of the expected values), and modest enrichment in TAD (Dixon et al. (2012)). Importantly, a very strong (26-fold) enrichment in variable methylation regions (VMR) was observed (Hansen et al. (2011)) (Fig. 1(e)), suggesting that increased epigenetic variability in a cell population or tissue can be coordinated locally among hundreds of thousands of genomic regions (Pujadas et al. (2012)). A subset of MHBs that do not overlap with CpG islands was

examined, and a consistent enrichment pattern (Fig. 1(e)) was observed, suggesting that local CpG density alone does not account for the enrichment.

[0060] Previous studies on mouse and human samples (Irizarry et al. (2009); Ziller et al. (2013)) demonstrated that dynamically methylated regions were associated with regulatory regions such as enhancer-like regions marked by H3K27ac and transcription factor binding sites. In human, 21.8% of autosomal CpGs were found to be differentially methylated across 30 human cell and tissue types (Hansen et al. (2011)). These CpGs were enriched at low to intermediate CpG density promoters. Using publicly available histone mapping data for human adult tissues, co-localization of methylation haplotype blocks with marks for active promoters (H3K4me3 with H3K27ac) was found, but not for active enhancers (Leung et al. (2015)) (no peak for H3K4me1). Meanwhile, enhancers were found to tend to overlap with CpG sparse MHBs, whereas the overlap with super-enhancers were independent of CpG density (Fig. 11(c)). Therefore, MHBs likely capture the local coherent epigenetic signatures that are directly or indirectly coupled with transcriptional regulation.

[0061] Block-level analysis of human normal tissues and stem cell lines with methylation haplotype load. To enable quantitative analysis of the methylation patterns within individual MHBs across many samples, a single metric to define the methylated pattern of multiple CpG sites within each block is needed. Ideally this metric is not only a function of average methylation level for all the CpG sites in the block, but can also capture the pattern of co-methylation on single DNA molecules. For this purpose, methylation haplotype load (MHL) and unmethylated haplotype load (uMHL) were defined, where the first is a weighted mean of the fraction of fully methylated haplotypes and substrings at different lengths (i.e. all possible substrings) and the latter is a weighted mean of the fraction of fully unmethylated haplotypes and substrings at different lengths. Compared with other metrics used in the literature (methylation level, methylation entropy, epi-polymorphism, and haplotypes counts), MHL is capable of distinguishing blocks that have the same average methylation but various degrees of coordinated methylation (Fig. 2). In addition, MHL and uMHL are bounded between 0 and 1, which allows for direct comparison of different regions across many data sets without normalization.

[0062] The invention addressed whether treating MHBs as individual genomic elements and performing quantitative analysis based on MHL and uMHL would provide an advantage over previous approaches using the weighted average methylation in genomic

windows. To this end, the invention sought to identify tissue specific MHBs from a collection of human solid tissues WGBS datasets based on the MHL and uMHL. A group specific index for each MHB (see Methods) was computed using either MHL or uMHL for every tissue type and then each MHB was assigned to one or more tissues based on their respective GSI. If no other tissue type have a GSI that is at least 80% of the maximum GSI, then only the tissue type with the maximum GSI is reported. From the top 500 MHBs assigned to each tissue, the top tissue specific MHL or uMHL regions were selected using average MHL or uMHL values in blood cells of less than 0.05, and average non-blood cell values greater than 0.4. Using a maximum GSI of greater than 0.6 as a threshold resulted in the identification of 1,290 and 15,377 tissue specific MHL and uMHL regions respectively. To identify the most informative MHL or uMHL marker regions for cancer detection or tissue mapping the top 10% of the identified MHL or uMHL regions determined by GSI value are reported in Tables 1(a) and 1(b). The tissue specific regions using MHL and uMHL were visualized using heatmaps (Fig. 3(a)-3(b)). To demonstrate that MHL was a better metric than AMF, a set of MHBs overlapping with published tissue specific methylated regions (Lokk et al. (2014)) were identified. Using this set of MHBs, the performance between MHL and average methylation fraction in the MHL regions (AMF) were compared. Both MHL and AMF were able to identify tissue specific MHBs, but MHL has better signal (average value from within the tissue specific group) to noise (average value from other groups) than average methylation (Fig. 4).

[0063] The human adult tissues used in this invention have various degrees of similarity amongst each other. The invention hypothesized that this is primarily defined by their developmental lineage, and that the related MHBs might reveal epigenetic insights related to germ layer speciation. All the data sets based on the three germ layers were grouped, and searched for MHBs that have differential MHL. In total 114 ectoderm-specific MHBs (99 hyper- and 15 hypo-methylated), 75 endoderm specific MHBs (58 hyper- and 17 hypo-methylated) and 31 mesoderm specific MHBs (9 hyper- and 22 hypo-methylated) were identified. The invention speculated that some of these MHBs might capture binding events of transcription factors (TF) specific to developmental germ-layers. Compared with ENCODE TFBS data (The ENCODE Project Consortium (2012)), distinctive patterns of TFs binding to layer specific MHBs were observed. For layer specific MHBs with hypo-methylation MHL, which tends to represent activation signals, 53 TF binding events in mesoderm specific MHBs, 71 in endoderm specific MHB, and 2 in ectoderm specific MHBs were identified. Gene

ontology analysis showed TFs binding to mesoderm exhibit negative regulator activity, while TFs binding to endoderm exhibited positive regulator activity. For layer specific MHBs with hyper-methylation MHL, which tend to represent repressive signals, 38 TF binding events in mesoderm specific MHBs, 102 in endoderm specific MHB, and 145 in ectoderm specific MHBs were identified. Interestingly, ectoderm and endoderm shared few bounded TFs, while mesoderm tissues share multiple groups of TFs with ectoderm and endoderm. Two endoderm specific hyper-MHL regions were identified, which are related to *ESRRA* and *NANOG*. This is consistent with a previous finding that mouse ES cells differentiated spontaneously into visceral/parietal endoderm upon NANOG knock-out (Mitsui et al. (2003)). Gene ontology analysis showed that hypo-MHL regions shared by mesoderm and endoderm might have regulatory functions in the fate commitment towards multiple tissues, whereas ectoderm specific hyper-MHL regions might induce the ectoderm development by suppressing the path towards the immune lineage. These observations are indicative of two distinctive “push” and “pull” mechanisms in the transition of cell states that have been harnessed for the induction of pluripotency by over-expressing lineage specifiers (Shu et al. (2013)).

[0064] Methylation-haplotype based analysis of circulating cell-free DNA in cancer patients and healthy donors. A unique aspect of methylation haplotype analysis is that the pattern of co-methylation, especially within MHBs, is robust for capturing low-frequency alleles among a heterogeneous population of molecules or cells, in the presence of biological noise or technical variability (i.e. incomplete bisulfite conversion or sequencing errors). To explore the clinical potential, the invention focused on the methylation haplotype analysis of cell-free DNA from healthy donors and cancer patients, of which various low fractions of DNA molecules were released from tumor cells and potentially carry epigenetic signatures different from blood. 4-122ng (average 20ng) of cell-free DNA from an average of 866μL human plasma from 75 normal individuals and 59 cancer patients were isolated, except for four with unusually high yield due to cell lysis. Due to the limited DNA availability, scRRBS was performed (Guo et al. (2013)) on 1 to 10 ng of cfDNA from 134 plasma samples and obtained an average of 13 million paired-end 150bp reads per sample. On average, 57.7% WGBS-defined MHBs were covered in the RRBS data set from the clinical samples.

[0065] The invention sought to detect the presence of tumor specific signatures in the plasma samples, using the methylation haplotypes identified in the reference tumor tissues and in normal samples used as the negative controls. For five lung cancer plasma samples and five

colorectal cancer plasma samples, matched primary tumor tissues were also obtained, and generated RRBS data (30 million reads per sample) from 100ng of tumor genomic DNA. The invention focused on MHBs with low MHL (i.e. genomic regions that have low or no methylation) in the blood, and determined whether cancer-associated highly methylated haplotypes (caHMH) can be detected. The invention provides in some embodiments that such haplotypes are present only in the tumor tissues and the matched plasma from the same patient, but not in whole blood or any other non-cancer samples. These highly confident tumor signatures in circulating DNA were considered. caHMH in all cancer patient plasma samples (Average=36, interquartile range (IQR)=17) were detected. These caHMHs were associated with 320 genes, some of which are known to be aberrantly methylated in human cancers such as *WDR37*, *VAX1*, *SMPD1* (Table 2). The 49 additional cancer plasma samples with no matched tumor samples were analyzed, using 65 normal plasma samples as background and negative control. On average 60 (IQR=31) caHMH were identified for each cancer plasma sample. Interestingly, a significant fraction (35%) of caHMH identified with matched tumor-plasma pairs were also detected in the expanded set of cancer patient plasma samples. A majority of caHMHs were found to be individual specific while few caHMHs were present in at least 53% (16/30) and 62% (18/29) cancer plasma samples for CRC and LC (Fig. 13). Improving the sampling depth, by either using more input cfDNA or reducing sample loss during analysis, will likely increase the number of caHMHs commonly observed in multiple patients.

[0066] The tumor load in cancer plasma samples was quantified, using non-negative decomposition with quadratic programming, on the RRBS data from primary cancer biopsies (LC and CRC) and the WGBS data from 10 normal tissues. (Table 3) The invention estimated that a predominant fraction, 72.0% (IQR=40%) in cancer and normal plasma samples were contributed by white blood cells, which is consistent with the levels recently reported based on shallow whole genome bisulfite sequencing (69.4%) (Sun et al. (2015)). Primary tumor and normal tissue-of-origin contributed at the similar level of 2.3% (IQR=3.7%) and 3.0% (IQR=4.4%). In contrast, the similar analysis applied to normal plasma only found residual tumor contributions (0.17%, IQR=2.9% for CRC and 1.0%, IQR=3.1% LC) to normal plasma, which were significantly lower ($P=3.4 \times 10^{-5}$ and 5.2×10^{-10} for CRC and LC, respectively) than for cancer plasma. 76.7% plasma samples from CRC patients and 89.6% from LC patients were found to have a detectable contribution from tumor tissues while only 13% and 26%

normal plasmas have certain (low) tumor contribution. Therefore, circulating cell-free DNA contains a relatively stable fraction of molecules released from various normal tissues, whereas tumor cells released DNA molecules that can be more abundant in tumor cells than in normal tissues (Tables 3(a)-3(d)). The fractions of white blood cells observed are lower than what was
5 reported previously (Sun et al. (2015)), most likely due to the inclusion of 10 normal tissue types in the deconvolution analysis.

[0067] The invention sought to use the information from normal human tissues, primary tumor biopsies, and cancer cell lines to improve the detection of cfDNA. The invention selected a subset of MHBs that show high MHL in primary cancer biopsies and low MHL in
10 normal control plasma. A subset of MHBs that have high MHL in cancer tissues and low MHLs in normal plasma was identified for each cancer type (Table 4(a)-(b)). Cancer plasma showed significantly higher MHL in these regions than independent normal plasma ($P=1.2 \times 10^{-13}$ and 2.2×10^{-16} for CRC and LC, respectively) (Fig. 5 and Fig. 6). By computationally
15 mixing the sequencing reads from cancer tissues and whole blood samples (WBC), synthetic admixtures at various levels of tumor fraction were created to calibrate the relationship between tumor load and the MHL values in these regions (Fig. 7). A fitted linear model on the standard curve for colon cancer markers had an adjusted R-squared value of 0.9621 and for lung cancer
20 markers the adjusted R-squared value was 0.9573. Note that these MHBs were selected without using any information from the cancer plasma samples, and hence they should be applicable to other cancer plasma samples. Tumor load estimation was performed on test sets of cancer and
normal plasma samples. Cancer plasma was found to have significantly higher tumor load than normal plasma (Two Sample t-test with unequal variance, $P=0.005074$ for colon cancer plasma versus normal and $P=0.0002783$ for lung cancer plasma versus normal plasma) (Fig. 8).

[0068] Recent studies (Sun et al. (2015); Lehmann-Werman et al. (2016); Snyder et al.
25 (2016)) have demonstrated that epigenetic information imbedded in cfDNA has the potential for predicting tumors tissue-of-origin. Here the invention addresses whether a MHL-based framework and a set of targets derived from whole genome data would allow for the prediction of tissue-of-origin with quantifiable sensitivity and specificity, which is crucial for future
clinical applications. Training WGBS tissue data from the following tissues: colon, lung, neural, heart, liver, lung, pancreas, and stomach were compiled and a set of 15,000 MHBs
30 selected using their Group Specific Index (GSI) was created. Note that the MHBs with a high GSI score tended to be methylated in fewer tissue types. 20 colon cancer plasma, 20 lung cancer

plasma, and 30 normal plasma samples were randomly sampled to create training data sets from the RRBS plasma data sets. The remaining (10 colon, 10 lung, and 39 normal) were held out as test data sets. An ensemble MARS (Multivariate adaptive regression splines) model was generated using the training data sets (Friedman, (1991)). The invention also employed K-means clustering of the MHBs using WGBS data to generate 50 clusters and each MARS model in the ensemble only saw a feature set consisting of 3 features sampled from a cluster (3 features x 150 clusters = 450 features for each model). To minimize overfitting of training data sets each model selected at most 2 features to use during classification. The first classifier was an ensemble model with 100 MARS models with 154 MHB features (Table 5(a)). In a manner analogous to a random forest classifier, the scores from the resulting ensemble are averaged to compute a prediction score. By utilizing these prediction scores as a binary classifier, the invention was able to obtain a colon vs non-colon AUC of 0.725, a lung vs non-lung AUC of 0.751, and a normal vs abnormal AUC of 0.856 (Fig. 9(a)-9(c)). Note that these AUC values were obtained using an initial marker set identified from an independent training set of normal WGBS tissue samples only; no cancer tissue samples were used in the initial feature selection, suggesting that plasma samples can be segregated using only tissue specific markers.

[0069] To obtain a prediction of tissue-of-origin, the invention focused on colon and lung cancer plasma. First, the top MHBs features that were assigned to either colon, lung, or both tissues ranked by GSI were selected. An ensemble of 500 MARS models identifying 295 unique tissue specific features was made (Table 5(b)). PLSDA (partial least squares discriminant analysis) implemented in the caret R package was used to perform the classification. The classification using this model was assessed using independent test data sets and the accuracy for classifying colon cancer plasma to colon tissue was found to be 70% (7/10) and the accuracy for classifying lung cancer plasma to lung tissue was 78% (7/9) (Fig. 10).

[0070] In this invention a well-established concept in population genetics, linkage disequilibrium (LD), was extended to the analysis of co-methylated CpG patterns. While the mathematical representations are identical, there are two key differences. First, traditional linkage disequilibrium (LD) was defined for human individuals in a population, whereas in this invention the analysis was performed on the diploid genome of individual cells in a heterogeneous cell population. Second, linkage disequilibrium in human populations depend on the mutation rate, frequency of meiotic recombination, effective population size and demographic history. The LD level typically decays over the range of hundreds of kilobases

-26-

to megabases. In contrast, CpG co-methylation depends on DNA methyltransferases and demethylases, which tend to have lower processivity, and, in the case of hemimethyltransferases, much lower fidelity compared with DNA polymerases (Williams et al. (2011)). Therefore, methylation LD decays over much shorter distances ranging from tens to
5 hundreds of bases, with the exception of imprinting regions. Even if longer-read sequencing methods were used no radical change of the block-like pattern presented in this work is expected, which is supported by a recent study (Saito et al. (2015)). Nonetheless, these short and punctuated blocks capture discrete entities of epigenetic regulation in individual cells widespread in the human genome. Such a phenomenon can be harnessed to improve the
10 robustness and sensitivity of DNA methylation analysis, such as the deconvolution of data from heterogeneous samples including circulating cell-free DNA.

[0071] Epigenetic abnormalities tend to be more widespread across the genome (compared with somatic mutations), and hence enabling the integration of the sparse coverage across many loci to achieve very accurate prediction by direct counting of methylated
15 haplotypes within the appropriate tissue-specific features.

Methods

[0072] Normal and cancer samples. Ten human primary tissues were purchased from BioChain. Cancer tissue and plasma samples were collected from UCSD Moores Cancer Center and normal plasma samples were obtained from UCSD Shiley Eye center under IRB
20 protocols approved by UCSD Human Research Protections Program (HRPP).

[0073] Generation of DNA libraries for sequencing. Extracted genomic DNA were prepared for bisulfite sequencing using published protocols. For whole genome bisulfite (WGBS) and reduced representation bisulfite sequencing (RRBS), the DNA fragments were adapted to barcoded methylated adaptors (Illumina). For WGBS, the adapted DNA were
25 converted using the EZ DNA Methylation Lightning kit (Zymo Research) and amplified for 10 cycles using iQ SYBR Green Supermix (BioRad). For RRBS, the adapted DNA were converted using the MethylCode™ Bisulfite Conversion kit (Thermo Fisher Scientific) and amplified using the PfuTurboC_x polymerase (Agilent) for 12-14 cycles. Libraries were pooled and size selected using 6% TBE polyacrylamide gels. Libraries were sequenced using the
30 Illumina HiSeq platform for paired-end 100-111 cycles, the Illumina MiSeq platform for paired-end 75 cycles, and the GAIIx (WGBS only) for single-end 36 cycles.

-27-

[0074] Methylation haplotype blocks (MHB). The human genome was separated into non-overlapping “sequenceable and mappable” segments using a set of generated WGBS data from 10 tissues from a 25-year adult male individual. Mapped reads from WGBS data sets were converted into methylation haplotypes in each segment. Methylation linkage disequilibrium was calculated on the combined methylation haplotypes. Each segment was partitioned into methylation haplotype blocks (MHBs). MHBs were defined as the genomic region in which the r^2 value of two adjacent CpG sites is no less than 0.5. MHB regions inferred by WGBS data sets were also validated by bulk data of methylation level. Takai and Jones's sliding-window algorithm (Takai et al. (2002)) was applied for methylation high linkage regions in the HM450K (TCGA) and the RRBS (ENCODE) data set. Finally, simulation analysis to investigate the relationship between LD and correlation of average 5mC of two CpG loci were conducted based on random sampling of different methylation haplotypes with 1000 individuals, and each individual sampling 10 methylation haplotypes.

[0075] Methylation haplotype load (MHL). A methylated haplotype load (MHL) for each candidate region was defined, which is the normalized fraction of methylated haplotypes (MH) at different lengths:

$$MHL = \frac{\sum_{i=1}^l w_i \times P(MH_i)}{\sum_{i=1}^l w_i}$$

$$w_i = i$$

[0076] Unmethylated haplotype load (uMHL). An unmethylated haplotype load (uMHL) for each candidate region was defined, which is the normalized fraction of unmethylated haplotypes (UMH) at different lengths:

$$MHL = \frac{\sum_{i=1}^l w_i \times P(UMH_i)}{\sum_{i=1}^l w_i}$$

$$w_i = i$$

[0077] Where l is the length of haplotypes, $P(MH_i)$ or $P(UMH_i)$ is the fraction of fully successive methylated or unmethylated haplotype with i loci respectively. For a haplotype of length L , all the sub-strings with length from 1 to L were considered in this calculation. w_i is the weight for i -locus haplotype. $w_i = i$ or $w_i = i^2$ was typically used to favor the contribution of longer haplotypes. In the present invention, $w_i = i$ was applied. Quantile normalization,

standardization (scale) as well as the batch effect elimination (Johnson et al. (2007)) were applied.

[0078] Developmental germ layers and tissue specific MHB regions. In order to investigate the layer and tissue specific MHB regions, group specific index (see below) was applied. An empirical threshold of 0.6 was selected to filter out layer and tissue specific MHB regions. Layer specific MHB regions were selected to show the ability to distinguish between the different development layers. Tissue specific MHB regions were further used for tissue mapping and cancer diagnosis.

$$GSI = \frac{\sum_{j=1}^n 1 - \frac{10^{MHL_j}}{10^{MHL_{max}}}}{n - 1}$$

10 n indicates the number of the groups. $MHL(j)$ denotes the average of MHL of j^{th} group. MHL_{max} denotes the average of MHL of highest methylated group.

[0079] Deconvolution analysis. The deconvolution references were constructed from normal human solid tissues, WBC, colorectal cancer tissues (CCT), and lung cancer tissues (LCT). Tissue specific MHB regions for normal human tissues were selected for brain, colon, esophagus, heart, intestine, kidney, liver, lung, and stomach using candidate features for deconvolution based on non-negative decomposition with quadratic programming (Sun et al. (2015); Houseman et al. (2012); Gong et al. (2013)). Raw MHL signals were logit transformed before deconvolution analysis on plasma samples. Samples with less than 30% whole blood content from deconvolution analysis were considered to have failed due to poor library complexity.

[0080] Tumor load estimation using tumor specific MHBs. Pruning and K-nearest neighbors (KNN) imputation was performed on the MHL matrix with only the RRBS plasma and RRBS tumor tissue samples removing samples with low coverage and imputed missing values. Thus 30 colon cancer plasma, 29 lung cancer plasma, 69 normal plasma, 4 colon cancer tissues, and 5 lung cancer tissues remained. The 69 healthy normal plasma samples were split into “training” and “test” sets; with 46 samples set aside for feature selection and training while the remaining 23 samples were used as a completely independent data set to test the quantitation. Tumor specific methylation haplotype blocks were identified by 2-tailed t-test with a False Discovery Rate (FDR) of 0.001 and a minimum difference cutoff of 0.3. Two sets of markers were identified separately for colon cancer and lung cancer. To calibrate the

-29-

- relationship between tumor load and MHL values, 20 sets of simulations in which mixed sequencing reads from cancer tissue samples and normal plasma samples at a 1:5, 1:10, 1:20, 1:100, and 0:1 ratio were generated (totaling 100 simulated data sets for each cancer tissue). The average MHL value for these regions was computed for each region and a linear regression
- 5 model was generated using these values and the known cancer tissue proportions for each cancer tissue. The model was applied to the average MHL value in these regions for test sets which included 30 colon cancer plasma, 29 lung cancer plasma, and 23 normal plasma samples. The estimated tumor loads for normal versus colon cancer and for normal versus lung cancer were compared using the colon cancer and lung cancer markers respectively.
- 10 **[0081]** Cancer plasma classification. Training WGBS data were collected to generate an MHL matrix and the top 15,000 MHBs by GSI were selected. These MHBs were clustered using K-means into 50 marker groups with the WGBS data matrix. Each cluster should have similar methylation patterns across each tissue type. In order to avoid overfitting, 20 colon cancer plasma, 20 lung cancer plasma, and 30 normal plasma samples were randomly selected
- 15 to create a training data set from the RRBS plasma samples. The remaining samples (10 colon, 10 lung, and 39 normal) were held out as a test data set. For building the classifiers, an ensemble MARS (Multivariate Adaptive Regression Splines) model (Friedman, (1991)) implemented in the Earth R package (<https://cran.r-project.org/web/packages/earth/index.html>) was trained on the training data set using features from the top 15,000 MHBs ranked by the Group Specific
- 20 Index (GSI). The ensemble MARS model was used to perform classification on the test plasma data set using the binary classifier. For cancer type classification, the features selected using MARS were used in a PLSDA (partial least squares discriminant analysis) to distinguish colon versus lung.

-30-

Data Availability

[0082] WGBS and RRBS data are available at the Gene Expression Omnibus (GEO) under accession GSE79279.

Tables:

- 5 **[0083]** Table 1. Top 10% tissue specific MHL and uMHL markers identified by GSI.
- [0084]** Table 2. Complete list of high methylated haplotype shared between matched primary tumor tissues and plasma for colon cancer (CRC) and lung cancer (LC) patients.
- [0085]** Table 3. Deconvolution of colon cancer (CRC), lung cancer (LC) and normal plasma into ten tissues using non-negative decomposition with quadratic programming.
- 10 **[0086]** Table 4. Differentially methylated MHB regions between cancer tissues and normal plasma.
- [0087]** Table 5. The sets of cancer specific and tissue specific markers derived from MARS based features selection on training data sets.

15

Table 1a:
Top 10% tissue specific uMHL markers identified by GSI

uMHL Markers			
Region	Group	GSI	refMax
chr19:5894163:5894242	vessel	8.75E-01	9.94E-01
chr16:521904:521925	neural	8.73E-01	9.60E-01
chr2:43295338:43295363	neural	8.71E-01	9.72E-01
chr2:69345847:69345875	vessel	8.71E-01	9.63E-01
chr2:8360246:8360318	vessel	8.69E-01	9.60E-01
chr2:110103840:110103879	vessel	8.69E-01	9.44E-01
chr10:73767213:73767231	neural	8.67E-01	1.00E+00
chr21:39450802:39450857	vessel	8.64E-01	9.63E-01
chr8:23201701:23201725	vessel	8.63E-01	9.29E-01
chr8:96706051:96706130	vessel	8.63E-01	9.62E-01
chr10:17281034:17281085	vessel	8.63E-01	9.76E-01
chr18:9535925:9535962	vessel	8.62E-01	9.64E-01
chr1:196373497:196373569	vessel	8.61E-01	9.74E-01
chr7:73314135:73314205	vessel	8.61E-01	9.57E-01
chr5:168192470:168192555	vessel	8.59E-01	9.41E-01
chr14:83966994:83967047	vessel	8.59E-01	9.62E-01
chr10:14012644:14012740	vessel	8.59E-01	9.93E-01
chr16:4420959:4421041	vessel	8.59E-01	9.93E-01
chr2:8360333:8360384	vessel	8.58E-01	9.41E-01
chr1:2899575:2899616	vessel	8.57E-01	9.97E-01
chr19:41932380:41932387	intestine	8.57E-01	9.37E-01
chr2:232087015:232087102	vessel	8.56E-01	9.89E-01
chr17:738931:738960	vessel	8.56E-01	9.80E-01
chr17:48243247:48243305	vessel	8.55E-01	9.83E-01
chr19:768642:768715	vessel	8.54E-01	9.58E-01
chr2:38460795:38460935	vessel	8.54E-01	9.42E-01
chr6:57123002:57123073	vessel	8.53E-01	9.82E-01
chr1:243368706:243368788	vessel	8.53E-01	9.90E-01
chr16:73086441:73086558	vessel	8.53E-01	9.53E-01
chr6:1702288:1702366	vessel	8.52E-01	9.38E-01
chr10:45676961:45677042	vessel	8.52E-01	9.71E-01
chr12:20700689:20700721	vessel	8.52E-01	9.80E-01
chr12:20254130:20254176	vessel	8.51E-01	9.35E-01
chr2:3496911:3496968	vessel	8.51E-01	9.23E-01
chr11:69235366:69235449	vessel	8.51E-01	9.82E-01
chr17:38605979:38605996	liver	8.51E-01	1.00E+00
chr8:10001048:10001097	neural	8.50E-01	8.60E-01
chr11:44994687:44994725	vessel	8.50E-01	9.53E-01
chr2:110861027:110861197	intestine	8.50E-01	8.87E-01
chr3:64702129:64702144	vessel	8.49E-01	9.43E-01
chr7:73389579:73389642	vessel	8.49E-01	9.67E-01
chr10:14013737:14013765	vessel	8.49E-01	9.82E-01
chr3:4458679:4458863	vessel	8.49E-01	9.59E-01
chr7:703821:703897	vessel	8.49E-01	9.78E-01
chr13:101302763:101302821	neural	8.49E-01	8.35E-01
chr19:16178379:16178427	vessel	8.48E-01	9.83E-01

uMHL Markers			
Region	Group	GSI	refMax
chr7:73389660:73389669	vessel	8.48E-01	9.48E-01
chr2:43492333:43492418	vessel	8.48E-01	9.47E-01
chr14:93113577:93113776	vessel	8.48E-01	9.56E-01
chr22:40845219:40845265	vessel	8.48E-01	9.59E-01
chr1:115610338:115610366	vessel	8.47E-01	9.66E-01
chr6:131312703:131312868	vessel	8.47E-01	9.70E-01
chr16:81520175:81520330	vessel	8.47E-01	9.83E-01
chr20:19474573:19474667	vessel	8.46E-01	9.81E-01
chr1:244217148:244217329	neural	8.46E-01	8.57E-01
chr16:49822649:49822660	vessel	8.46E-01	9.62E-01
chr1:202170004:202170040	vessel	8.46E-01	9.99E-01
chr3:8562700:8562922	vessel	8.46E-01	9.60E-01
chr8:1187134:1187185	neural	8.45E-01	1.00E+00
chr15:54832956:54833002	vessel	8.45E-01	9.57E-01
chr11:3168353:3168372	vessel	8.45E-01	9.83E-01
chr12:116864174:116864293	vessel	8.45E-01	9.69E-01
chr17:31128334:31128406	vessel	8.45E-01	9.83E-01
chr2:72162546:72162581	vessel	8.44E-01	9.72E-01
chr9:137553753:137553885	vessel	8.44E-01	9.59E-01
chr16:66957496:66957553	vessel	8.44E-01	9.88E-01
chr6:169568133:169568353	vessel	8.44E-01	9.36E-01
chr8:133466185:133466297	vessel	8.43E-01	9.54E-01
chr12:124774360:124774380	vessel	8.43E-01	1.00E+00
chr19:3466975:3467064	vessel	8.43E-01	9.82E-01
chr9:98829531:98829605	vessel	8.43E-01	9.58E-01
chr17:37279963:37280015	vessel	8.43E-01	9.39E-01
chr9:116247804:116247934	vessel	8.43E-01	9.72E-01
chr1:87223214:87223344	vessel	8.43E-01	9.52E-01
chr11:66138099:66138125	intestine	8.42E-01	9.22E-01
chr5:142533336:142533503	vessel	8.42E-01	9.86E-01
chr7:4065672:4065679	vessel	8.42E-01	9.66E-01
chr17:73831565:73831633	vessel	8.42E-01	9.67E-01
chr16:1373394:1373461	neural	8.42E-01	8.62E-01
chr17:40477211:40477303	vessel	8.42E-01	9.78E-01
chr6:165341974:165342035	neural	8.42E-01	8.81E-01
chr18:58648459:58648474	vessel	8.42E-01	9.33E-01
chr10:3928761:3928829	vessel	8.41E-01	9.68E-01
chr15:67457875:67458134	vessel	8.41E-01	9.44E-01
chr19:10233053:10233111	vessel	8.41E-01	9.42E-01
chr18:74171483:74171505	vessel	8.41E-01	9.92E-01
chr22:49409034:49409082	vessel	8.41E-01	9.19E-01
chr21:40047317:40047326	vessel	8.41E-01	9.33E-01
chr1:34451152:34451164	vessel	8.41E-01	9.86E-01
chr9:136357330:136357347	vessel	8.41E-01	1.00E+00
chr9:116681630:116681838	vessel	8.41E-01	9.21E-01
chr7:5011476:5011523	esophagus	8.41E-01	9.33E-01
chr16:87261081:87261131	vessel	8.39E-01	9.76E-01
chr5:172194371:172194450	vessel	8.39E-01	9.66E-01

uMHL Markers			
Region	Group	GSI	refMax
chr19:32450501:32450625	pancreas	8.39E-01	8.40E-01
chr15:79052333:79052347	vessel	8.39E-01	1.00E+00
chr12:116756805:116756874	vessel	8.39E-01	9.58E-01
chr12:2457684:2457778	vessel	8.39E-01	9.96E-01
chr13:36273480:36273646	vessel	8.39E-01	9.52E-01
chr2:145764662:145764766	vessel	8.39E-01	9.62E-01
chr4:140968580:140968766	vessel	8.38E-01	9.44E-01
chr7:158890050:158890132	vessel	8.38E-01	9.15E-01
chr3:71586325:71586633	vessel	8.38E-01	9.82E-01
chr3:125819901:125819917	vessel	8.38E-01	9.20E-01
chr18:76551153:76551172	vessel	8.38E-01	9.48E-01
chr22:29347978:29348075	vessel	8.38E-01	9.88E-01
chr2:11526450:11526505	vessel	8.38E-01	9.49E-01
chr14:91765021:91765059	vessel	8.38E-01	9.94E-01
chr2:10544929:10545012	vessel	8.38E-01	9.64E-01
chr2:217839781:217839805	vessel	8.37E-01	9.71E-01
chr12:109179149:109179194	vessel	8.37E-01	9.33E-01
chr1:226128702:226128727	vessel	8.37E-01	9.38E-01
chr7:40240562:40240640	vessel	8.37E-01	9.49E-01
chr14:75039798:75039894	vessel	8.37E-01	9.50E-01
chr20:56721581:56721650	vessel	8.37E-01	9.64E-01
chr8:6652013:6652077	vessel	8.36E-01	9.56E-01
chr16:1146322:1146337	vessel	8.36E-01	9.60E-01
chr10:15667379:15667446	vessel	8.36E-01	9.69E-01
chr15:89560186:89560238	vessel	8.36E-01	9.35E-01
chr7:4707844:4707866	vessel	8.36E-01	9.43E-01
chr5:151082200:151082218	vessel	8.36E-01	8.84E-01
chr3:193715471:193715560	vessel	8.36E-01	9.58E-01
chr9:93727370:93727392	vessel	8.36E-01	9.47E-01
chr8:97596962:97597040	vessel	8.36E-01	9.38E-01
chr21:44484510:44485019	vessel	8.36E-01	9.71E-01
chr3:8279761:8279814	vessel	8.36E-01	9.59E-01
chr16:73454353:73454372	vessel	8.36E-01	9.45E-01
chr3:14279072:14279273	vessel	8.35E-01	9.78E-01
chr6:168498836:168498871	pancreas	8.35E-01	9.74E-01
chr4:173973114:173973166	vessel	8.35E-01	9.77E-01
chr9:74431770:74431840	vessel	8.35E-01	9.82E-01
chr12:2396495:2396507	vessel	8.35E-01	9.51E-01
chr13:109807772:109807837	vessel	8.34E-01	9.61E-01
chr8:145019179:145019191	vessel	8.34E-01	9.41E-01
chr12:109182282:109182346	vessel	8.34E-01	9.87E-01
chr8:23201626:23201644	vessel	8.34E-01	8.96E-01
chr1:171326156:171326309	vessel	8.34E-01	9.06E-01
chr3:30538064:30538174	vessel	8.34E-01	9.46E-01
chr16:68857417:68857469	vessel	8.34E-01	9.67E-01
chr17:60774267:60774330	vessel	8.34E-01	9.67E-01
chr7:64020677:64020788	vessel	8.34E-01	9.43E-01
chr1:201748589:201748658	vessel	8.34E-01	9.79E-01

uMHL Markers			
Region	Group	GSI	refMax
chr19:18783073:18783151	vessel	8.34E-01	8.95E-01
chr13:112161714:112161827	vessel	8.34E-01	9.56E-01
chr11:117070485:117070493	vessel	8.34E-01	9.57E-01
chr7:128468433:128468462	vessel	8.34E-01	9.36E-01
chr1:2188830:2188947	neural	8.33E-01	8.29E-01
chr8:143695065:143695081	neural	8.33E-01	8.86E-01
chr6:71790637:71790777	vessel	8.33E-01	9.84E-01
chr7:73406729:73406811	vessel	8.33E-01	9.84E-01
chr12:1584268:1584438	vessel	8.33E-01	9.45E-01
chr2:1657462:1657556	vessel	8.33E-01	9.28E-01
chr2:121677850:121677924	vessel	8.33E-01	9.24E-01
chr13:99104177:99104223	vessel	8.33E-01	9.82E-01
chr6:23004460:23004497	vessel	8.33E-01	9.08E-01
chr7:33758956:33759144	vessel	8.33E-01	9.59E-01
chr11:65326863:65326909	vessel	8.33E-01	9.57E-01
chr3:134052418:134052557	vessel	8.33E-01	9.33E-01
chr16:1038028:1038079	vessel	8.33E-01	9.95E-01
chr7:25810506:25810570	vessel	8.32E-01	9.44E-01
chr11:19601607:19601696	vessel	8.32E-01	9.91E-01
chr3:54746711:54746836	vessel	8.32E-01	9.89E-01
chr3:126720805:126720813	vessel	8.32E-01	9.64E-01
chr21:46453929:46453954	vessel	8.32E-01	9.61E-01
chr4:6773182:6773272	vessel	8.32E-01	9.45E-01
chr5:149980906:149980972	vessel	8.32E-01	9.59E-01
chr5:77830138:77830190	vessel	8.32E-01	8.86E-01
chr10:80339365:80339515	vessel	8.32E-01	9.23E-01
chr17:12989679:12989753	vessel	8.32E-01	9.26E-01
chr10:114344929:114344989	vessel	8.32E-01	9.41E-01
chr2:36719487:36719574	vessel	8.32E-01	9.65E-01
chr14:75446146:75446254	vessel	8.32E-01	8.61E-01
chr11:133895052:133895366	vessel	8.32E-01	9.55E-01
chr11:62370671:62370715	vessel	8.31E-01	9.97E-01
chr5:52897443:52897563	vessel	8.31E-01	8.89E-01
chr11:130545528:130545623	vessel	8.31E-01	9.84E-01
chr2:217405801:217405854	vessel	8.31E-01	9.64E-01
chr6:169018100:169018147	vessel	8.31E-01	9.28E-01
chr17:75598600:75598659	vessel	8.31E-01	9.48E-01
chr10:45675475:45675491	vessel	8.31E-01	9.64E-01
chr18:53074420:53074508	neural	8.31E-01	7.93E-01
chr2:54199164:54199243	vessel	8.31E-01	9.28E-01
chr5:168451596:168451712	vessel	8.31E-01	9.28E-01
chr2:241185226:241185267	vessel	8.31E-01	9.35E-01
chr7:70120035:70120076	vessel	8.30E-01	9.39E-01
chr10:30109892:30109931	vessel	8.30E-01	9.74E-01
chr7:2951480:2951552	vessel	8.30E-01	9.50E-01
chr16:79166842:79167002	vessel	8.30E-01	9.67E-01
chr1:3464719:3464876	vessel	8.30E-01	9.07E-01
chr13:110775448:110775493	vessel	8.30E-01	9.19E-01

uMHL Markers			
Region	Group	GSI	refMax
chr20:19272195:19272440	vessel	8.30E-01	9.84E-01
chr1:9386783:9386877	vessel	8.30E-01	9.70E-01
chr6:41549109:41549179	vessel	8.29E-01	9.60E-01
chr11:63826186:63826257	neural	8.29E-01	8.91E-01
chr1:3286118:3286163	vessel	8.29E-01	9.32E-01
chr16:88279603:88279618	vessel	8.29E-01	9.79E-01
chr3:187769246:187769343	vessel	8.29E-01	9.64E-01
chr2:121298318:121298336	vessel	8.29E-01	9.61E-01
chr10:87796214:87796244	vessel	8.29E-01	9.35E-01
chr11:111784374:111784426	vessel	8.29E-01	9.77E-01
chr16:49686422:49686568	vessel	8.29E-01	9.39E-01
chr10:711379:711694	vessel	8.28E-01	9.68E-01
chr4:71941315:71941425	vessel	8.28E-01	9.57E-01
chr13:52304379:52304567	vessel	8.28E-01	9.81E-01
chr10:10336914:10336945	neural	8.28E-01	8.62E-01
chr15:99498725:99498811	vessel	8.28E-01	9.59E-01
chr15:88153699:88153769	vessel	8.28E-01	9.76E-01
chr5:149980623:149980743	vessel	8.28E-01	9.73E-01
chr15:47765811:47765900	vessel	8.28E-01	9.64E-01
chr15:36894077:36894154	vessel	8.28E-01	9.72E-01
chr3:73414820:73414881	vessel	8.28E-01	9.26E-01
chr10:44197849:44197884	pancreas	8.28E-01	8.67E-01
chr18:60387698:60387879	neural	8.28E-01	8.60E-01
chr20:17822381:17822534	vessel	8.28E-01	9.58E-01
chr13:31470063:31470114	vessel	8.28E-01	9.67E-01
chr2:66074251:66074308	vessel	8.28E-01	9.59E-01
chr9:131398646:131398694	pancreas	8.28E-01	8.48E-01
chr1:234594452:234594635	liver	8.27E-01	8.36E-01
chr11:19778332:19778414	vessel	8.27E-01	9.35E-01
chr6:155593701:155593747	neural	8.27E-01	8.79E-01
chr5:1542556:1542597	vessel	8.27E-01	9.14E-01
chr17:2362969:2363062	vessel	8.27E-01	8.98E-01
chr10:131488146:131488308	vessel	8.27E-01	9.41E-01
chr16:86461393:86461510	vessel	8.27E-01	9.20E-01
chr13:111024091:111024104	vessel	8.27E-01	9.35E-01
chr10:105879872:105879934	vessel	8.27E-01	9.57E-01
chr2:216299024:216299143	vessel	8.27E-01	9.23E-01
chr2:206620645:206620696	vessel	8.27E-01	9.67E-01
chr6:90708046:90708215	vessel	8.27E-01	9.21E-01
chr1:3460222:3460269	vessel	8.27E-01	9.80E-01
chr16:88096861:88096956	vessel	8.27E-01	9.48E-01
chr7:159004650:159004678	vessel	8.26E-01	9.08E-01
chr8:22443383:22443437	vessel	8.26E-01	9.93E-01
chr20:35169012:35169055	vessel	8.26E-01	9.93E-01
chr19:39889822:39889941	liver	8.26E-01	8.82E-01
chr10:129721730:129721745	vessel	8.26E-01	9.38E-01
chr4:120636770:120636877	intestine	8.26E-01	8.81E-01
chr5:141628521:141628586	vessel	8.26E-01	9.23E-01

uMHL Markers			
Region	Group	GSI	refMax
chr5:59792001:59792053	vessel	8.26E-01	9.75E-01
chr6:14925566:14925731	vessel	8.26E-01	9.37E-01
chr3:148729585:148729652	vessel	8.26E-01	9.84E-01
chr5:117889340:117889421	vessel	8.26E-01	9.03E-01
chr18:74918812:74918906	vessel	8.26E-01	9.13E-01
chr19:2523639:2523852	vessel	8.26E-01	9.87E-01
chr2:218692718:218692749	vessel	8.26E-01	9.67E-01
chr7:70133955:70134155	vessel	8.26E-01	9.19E-01
chr14:37504457:37504520	vessel	8.26E-01	9.28E-01
chr2:121371664:121371708	vessel	8.26E-01	9.93E-01
chr1:203491536:203491617	vessel	8.26E-01	9.47E-01
chr11:114237690:114237718	vessel	8.26E-01	9.20E-01
chr13:111000221:111000344	vessel	8.26E-01	9.10E-01
chr18:12646025:12646103	vessel	8.26E-01	9.00E-01
chr11:64631168:64631193	vessel	8.26E-01	9.61E-01
chr7:1135564:1135633	vessel	8.25E-01	9.72E-01
chr5:158945587:158945820	vessel	8.25E-01	9.58E-01
chr18:46353742:46353814	vessel	8.25E-01	9.54E-01
chr2:137330457:137330466	vessel	8.25E-01	9.43E-01
chr5:77844053:77844099	vessel	8.25E-01	9.26E-01
chr5:177694377:177694398	vessel	8.25E-01	9.18E-01
chr20:35168933:35168988	vessel	8.25E-01	9.46E-01
chr3:70048580:70048701	vessel	8.25E-01	9.41E-01
chr11:12097899:12098112	vessel	8.25E-01	9.52E-01
chr2:20055744:20055772	vessel	8.25E-01	9.46E-01
chr1:10614315:10614356	liver	8.25E-01	8.10E-01
chr10:14012441:14012468	vessel	8.25E-01	9.80E-01
chr12:111707943:111708004	vessel	8.25E-01	9.64E-01
chr1:3460137:3460182	vessel	8.25E-01	9.98E-01
chr20:17830926:17830968	vessel	8.25E-01	8.96E-01
chr9:133895273:133895401	vessel	8.25E-01	9.81E-01
chr16:68386106:68386174	intestine	8.25E-01	9.04E-01
chr7:2719029:2719102	neural	8.25E-01	9.63E-01
chr2:174060123:174060220	vessel	8.25E-01	9.47E-01
chr11:12136324:12136406	vessel	8.25E-01	9.86E-01
chr2:42226193:42226234	vessel	8.25E-01	9.50E-01
chr22:33313719:33313805	vessel	8.25E-01	9.73E-01
chr11:126033028:126033048	vessel	8.25E-01	9.75E-01
chr11:1418719:1418773	neural	8.25E-01	8.11E-01
chr22:31489110:31489197	vessel	8.24E-01	9.79E-01
chr11:45260506:45260585	vessel	8.24E-01	9.07E-01
chr2:47040392:47040486	vessel	8.24E-01	9.90E-01
chr3:186819136:186819178	vessel	8.24E-01	9.62E-01
chr3:9153801:9153977	vessel	8.24E-01	9.38E-01
chr8:73559260:73559345	vessel	8.24E-01	9.72E-01
chr16:29142419:29142508	liver	8.24E-01	9.36E-01
chr1:9410453:9410468	vessel	8.24E-01	9.66E-01
chr21:47512644:47512690	vessel	8.24E-01	9.36E-01

uMHL Markers			
Region	Group	GSI	refMax
chr9:16247410:16247417	vessel	8.24E-01	9.76E-01
chr18:46309140:46309233	vessel	8.24E-01	9.79E-01
chr3:66413019:66413233	vessel	8.23E-01	9.14E-01
chr20:10641514:10641548	vessel	8.23E-01	9.54E-01
chr10:17280749:17280812	vessel	8.23E-01	9.77E-01
chr7:5592713:5592836	intestine	8.23E-01	8.14E-01
chr14:23318835:23318857	vessel	8.23E-01	9.86E-01
chr7:42779397:42779584	vessel	8.23E-01	9.61E-01
chr1:61914656:61914841	vessel	8.23E-01	9.79E-01
chr19:13962429:13962502	vessel	8.23E-01	9.20E-01
chr4:76270190:76270273	vessel	8.23E-01	9.13E-01
chr10:81915105:81915176	vessel	8.23E-01	9.79E-01
chr11:94573359:94573509	vessel	8.23E-01	9.41E-01
chr11:44729756:44729789	vessel	8.23E-01	9.53E-01
chr7:73424958:73425002	vessel	8.23E-01	9.17E-01
chr8:141960637:141960714	vessel	8.23E-01	8.98E-01
chr2:3452386:3452438	fat	8.23E-01	8.23E-01
chr2:101402810:101403044	vessel	8.23E-01	9.79E-01
chr8:77316559:77316617	neural	8.23E-01	8.38E-01
chr2:240733175:240733324	vessel	8.23E-01	9.25E-01
chr15:26135892:26136025	vessel	8.23E-01	9.09E-01
chr2:107193220:107193348	vessel	8.23E-01	9.14E-01
chr20:39392677:39392771	vessel	8.22E-01	9.35E-01
chr1:3314113:3314144	vessel	8.22E-01	8.57E-01
chr8:141057820:141057872	vessel	8.22E-01	9.57E-01
chr10:5660858:5660911	vessel	8.22E-01	9.61E-01
chr5:8696414:8696485	vessel	8.22E-01	9.70E-01
chr13:45289725:45289758	vessel	8.22E-01	9.17E-01
chr7:37039869:37039998	vessel	8.22E-01	9.52E-01
chr12:32293438:32293717	neural	8.22E-01	9.04E-01
chr5:2515778:2515802	neural	8.22E-01	8.27E-01
chr10:10336989:10337032	neural	8.22E-01	8.49E-01
chr9:72012687:72012707	vessel	8.22E-01	9.64E-01
chr16:88976391:88976404	vessel	8.22E-01	9.77E-01
chr17:30001481:30001529	vessel	8.22E-01	9.53E-01
chr16:80752918:80753020	vessel	8.22E-01	9.47E-01
chr11:130283367:130283389	vessel	8.22E-01	9.10E-01
chr16:72928691:72928785	vessel	8.22E-01	9.75E-01
chr11:10714431:10714466	vessel	8.22E-01	9.51E-01
chr19:18783185:18783216	vessel	8.22E-01	9.40E-01
chr16:61245717:61245888	vessel	8.22E-01	9.57E-01
chr9:129517613:129517699	intestine	8.21E-01	7.94E-01
chr9:124615315:124615387	vessel	8.21E-01	8.99E-01
chr6:15455978:15456143	vessel	8.21E-01	8.72E-01
chr5:137727491:137727661	neural	8.21E-01	8.51E-01
chr3:13560010:13560122	vessel	8.21E-01	9.49E-01
chr1:16346068:16346121	vessel	8.21E-01	9.30E-01
chr10:711944:712088	vessel	8.21E-01	9.18E-01

uMHL Markers			
Region	Group	GSI	refMax
chr22:49412171:49412287	vessel	8.21E-01	9.49E-01
chr2:40355055:40355244	vessel	8.21E-01	9.16E-01
chr22:29146497:29146625	vessel	8.21E-01	9.23E-01
chr16:14418702:14418748	vessel	8.21E-01	9.57E-01
chr21:42496716:42496736	vessel	8.21E-01	9.44E-01
chr2:1820299:1820400	vessel	8.21E-01	9.69E-01
chr13:111922451:111922625	vessel	8.21E-01	9.42E-01
chr7:158207285:158207423	neural	8.21E-01	8.86E-01
chr14:75086815:75086860	vessel	8.21E-01	9.32E-01
chr6:45784561:45784755	vessel	8.21E-01	8.71E-01
chr12:98869058:98869255	vessel	8.21E-01	9.13E-01
chr17:39680190:39680212	neural	8.21E-01	9.27E-01
chr2:3473947:3474025	vessel	8.21E-01	9.20E-01
chr6:16725500:16725574	vessel	8.21E-01	9.80E-01
chr3:123419686:123419782	vessel	8.21E-01	1.00E+00
chr5:78099270:78099340	vessel	8.21E-01	9.53E-01
chr4:82057849:82058045	vessel	8.20E-01	9.35E-01
chr7:73465618:73465668	vessel	8.20E-01	9.54E-01
chr7:30818374:30818439	vessel	8.20E-01	9.14E-01
chr4:88461671:88461843	vessel	8.20E-01	9.81E-01
chr1:3459727:3459820	vessel	8.20E-01	9.70E-01
chr12:109179337:109179529	vessel	8.20E-01	9.65E-01
chr14:51237067:51237132	vessel	8.20E-01	9.13E-01
chr3:141151136:141151212	vessel	8.20E-01	9.68E-01
chr1:203541608:203541617	vessel	8.20E-01	9.70E-01
chr18:3622294:3622340	vessel	8.20E-01	9.36E-01
chr1:3459875:3459893	vessel	8.20E-01	9.34E-01
chr12:52938381:52938412	vessel	8.20E-01	9.74E-01
chr4:5736721:5736775	vessel	8.20E-01	9.30E-01
chr1:39933015:39933081	vessel	8.20E-01	9.80E-01
chr19:16178735:16178766	vessel	8.20E-01	9.03E-01
chr20:57473842:57473881	vessel	8.20E-01	9.78E-01
chr12:123738593:123738662	vessel	8.19E-01	9.66E-01
chr16:88677685:88677878	neural	8.19E-01	7.95E-01
chr1:3142374:3142411	vessel	8.19E-01	8.72E-01
chr17:57565015:57565070	vessel	8.19E-01	9.36E-01
chr10:134527781:134527839	neural	8.19E-01	8.61E-01
chr10:30083908:30083982	vessel	8.19E-01	9.79E-01
chr2:72366403:72366433	vessel	8.19E-01	9.41E-01
chr4:24732279:24732305	vessel	8.19E-01	9.72E-01
chr17:75861540:75861581	vessel	8.19E-01	8.97E-01
chr7:66935160:66935181	vessel	8.19E-01	9.33E-01
chr14:69410822:69410995	vessel	8.19E-01	9.15E-01
chr5:175119757:175119816	vessel	8.19E-01	9.77E-01
chr16:86965611:86965651	vessel	8.19E-01	9.94E-01
chr5:52658535:52658625	vessel	8.19E-01	9.62E-01
chr5:109775450:109775504	vessel	8.19E-01	9.29E-01
chr5:73234564:73234688	vessel	8.19E-01	9.39E-01

uMHL Markers			
Region	Group	GSI	refMax
chr2:161777367:161777468	vessel	8.19E-01	9.72E-01
chr13:33922144:33922332	vessel	8.19E-01	9.79E-01
chr7:120164933:120164963	vessel	8.19E-01	9.30E-01
chr13:110918379:110918481	vessel	8.19E-01	9.16E-01
chr17:2011086:2011131	vessel	8.19E-01	9.69E-01
chr9:7984790:7984898	vessel	8.19E-01	9.40E-01
chr18:20682264:20682338	intestine	8.19E-01	8.40E-01
chr7:104978120:104978212	vessel	8.19E-01	9.63E-01
chr5:166634253:166634335	vessel	8.19E-01	9.03E-01
chr18:56534875:56534925	vessel	8.19E-01	9.70E-01
chr3:187988060:187988137	vessel	8.19E-01	9.76E-01
chr10:121438397:121438569	vessel	8.19E-01	9.68E-01
chr19:39182926:39182945	vessel	8.19E-01	9.35E-01
chr13:99667976:99668059	neural	8.18E-01	8.29E-01
chr7:127584205:127584527	neural	8.18E-01	8.66E-01
chr11:35651590:35651804	vessel	8.18E-01	9.19E-01
chr11:103407665:103407717	vessel	8.18E-01	8.98E-01
chr7:135433499:135433541	vessel	8.18E-01	9.64E-01
chr1:23098725:23098744	vessel	8.18E-01	9.51E-01
chr1:3494166:3494235	vessel	8.18E-01	9.52E-01
chr2:2107292:2107346	neural	8.18E-01	7.72E-01
chr8:101427557:101427651	vessel	8.18E-01	9.11E-01
chr10:14016066:14016097	vessel	8.18E-01	9.41E-01
chr11:10647667:10647725	vessel	8.18E-01	9.56E-01
chr16:75278952:75279001	vessel	8.18E-01	9.60E-01
chr9:124089044:124089101	vessel	8.18E-01	9.70E-01
chr17:21185776:21185831	vessel	8.18E-01	9.89E-01
chr3:11178498:11179102	vessel	8.17E-01	9.40E-01
chr17:57925734:57925805	intestine	8.17E-01	8.84E-01
chr20:16106391:16106635	vessel	8.17E-01	9.07E-01
chr2:241536131:241536192	vessel	8.17E-01	9.91E-01
chr5:76714624:76714674	vessel	8.17E-01	9.08E-01
chr8:42521965:42522125	vessel	8.17E-01	9.75E-01
chr11:36003464:36003500	vessel	8.17E-01	9.59E-01
chr20:23251275:23251362	vessel	8.17E-01	9.59E-01
chr1:168578476:168578512	vessel	8.17E-01	9.42E-01
chr20:50186917:50186939	vessel	8.17E-01	8.95E-01
chr3:49556866:49557077	neural	8.17E-01	8.96E-01
chr10:131324082:131324193	vessel	8.17E-01	9.67E-01
chr1:205780103:205780116	vessel	8.17E-01	9.51E-01
chr1:42003200:42003221	vessel	8.17E-01	9.57E-01
chr11:2730324:2730441	vessel	8.17E-01	9.80E-01
chr5:156988163:156988221	vessel	8.16E-01	8.84E-01
chr19:6744481:6744634	vessel	8.16E-01	8.99E-01
chr20:56472461:56472592	liver	8.16E-01	8.74E-01
chr6:1624318:1624399	vessel	8.16E-01	9.29E-01
chr12:106517592:106517795	vessel	8.16E-01	9.85E-01
chr6:44195262:44195347	vessel	8.16E-01	9.04E-01

uMHL Markers			
Region	Group	GSI	refMax
chr19:39182585:39182856	vessel	8.16E-01	9.65E-01
chr10:126104092:126104263	vessel	8.16E-01	9.75E-01
chr2:66528170:66528286	vessel	8.16E-01	9.69E-01
chr22:19971912:19971946	vessel	8.16E-01	9.42E-01
chr21:47403346:47403426	vessel	8.16E-01	9.03E-01
chr16:87234734:87234896	vessel	8.16E-01	9.64E-01
chr2:38611005:38611050	vessel	8.16E-01	9.82E-01
chr2:12677328:12677416	vessel	8.16E-01	9.22E-01
chr9:101749596:101749621	vessel	8.16E-01	9.13E-01
chr9:92300357:92300439	vessel	8.16E-01	9.28E-01
chr14:76461726:76461806	vessel	8.16E-01	9.80E-01
chr6:159624504:159624617	vessel	8.16E-01	9.82E-01
chr4:57909186:57909239	vessel	8.16E-01	9.74E-01
chr21:40195044:40195130	vessel	8.16E-01	9.57E-01
chr1:2899714:2899770	neural,vessel	8.16E-01	9.66E-01
chr1:9384754:9384813	vessel	8.16E-01	9.67E-01
chr16:57149767:57149848	neural	8.16E-01	8.76E-01
chr2:1656844:1657056	vessel	8.16E-01	9.72E-01
chr15:68777293:68777315	vessel	8.16E-01	9.89E-01
chr6:164506404:164506437	vessel	8.16E-01	9.35E-01
chr2:218722879:218722891	vessel	8.15E-01	9.67E-01
chr9:137536001:137536060	vessel	8.15E-01	9.40E-01
chr13:106439974:106440058	vessel	8.15E-01	9.44E-01
chr3:123132072:123132158	vessel	8.15E-01	9.64E-01
chr14:105802382:105802561	neural	8.15E-01	7.82E-01
chr18:7607703:7607797	vessel	8.15E-01	8.85E-01
chr1:2868537:2868570	vessel	8.15E-01	8.53E-01
chr19:1081908:1081925	intestine	8.15E-01	9.06E-01
chr7:3322236:3322467	vessel	8.15E-01	9.44E-01
chr19:10464320:10464378	liver	8.15E-01	9.20E-01
chr12:125106220:125106281	vessel	8.15E-01	9.49E-01
chr2:242668092:242668152	neural	8.15E-01	8.47E-01
chr15:26050364:26050411	vessel	8.15E-01	9.66E-01
chr16:1033828:1033849	neural	8.15E-01	7.58E-01
chr15:92573244:92573323	vessel	8.15E-01	9.82E-01
chr16:29206978:29207023	vessel	8.15E-01	9.69E-01
chr10:50155412:50155484	vessel	8.15E-01	9.06E-01
chr4:7844482:7844728	vessel	8.15E-01	9.33E-01
chr5:177799176:177799310	vessel	8.15E-01	9.86E-01
chr7:71563965:71564049	vessel	8.15E-01	9.22E-01
chr18:77398817:77398865	vessel	8.15E-01	9.82E-01
chr19:15361966:15362035	vessel	8.15E-01	9.85E-01
chr1:156111362:156111379	vessel	8.15E-01	9.21E-01
chr4:141207852:141207918	neural	8.15E-01	8.40E-01
chr18:46316977:46317049	vessel	8.15E-01	9.60E-01
chr19:7580118:7580276	vessel	8.15E-01	9.76E-01
chr2:238322615:238322681	vessel	8.15E-01	9.55E-01
chr11:4726406:4726437	vessel	8.15E-01	9.63E-01

uMHL Markers			
Region	Group	GSI	refMax
chr16:1243486:1243578	liver	8.14E-01	8.16E-01
chr9:101749649:101749733	vessel	8.14E-01	9.54E-01
chr11:12184939:12185109	vessel	8.14E-01	9.59E-01
chr15:26047396:26047413	vessel	8.14E-01	9.12E-01
chr1:204401170:204401223	vessel	8.14E-01	9.58E-01
chr1:59986591:59986617	vessel	8.14E-01	9.06E-01
chr17:43188128:43188157	vessel	8.14E-01	9.70E-01
chr13:24658482:24658536	neural	8.14E-01	7.67E-01
chr5:168168241:168168258	vessel	8.14E-01	8.15E-01
chr1:55909399:55909463	vessel	8.14E-01	9.66E-01
chr2:20001248:20001401	vessel	8.14E-01	9.92E-01
chr22:45907965:45908083	vessel	8.14E-01	9.42E-01
chr8:21499839:21499917	vessel	8.14E-01	9.58E-01
chr1:203526848:203526898	vessel	8.14E-01	9.59E-01
chr6:148561802:148561903	vessel	8.14E-01	9.39E-01
chr2:21676719:21676761	vessel	8.14E-01	9.26E-01
chr21:30502731:30502809	vessel	8.14E-01	9.58E-01
chr15:71621595:71621649	vessel	8.14E-01	9.32E-01
chr2:25078720:25078801	vessel	8.14E-01	9.30E-01
chr6:2355526:2355693	vessel	8.13E-01	9.42E-01
chr19:1136134:1136202	intestine	8.13E-01	8.14E-01
chr5:141628696:141628776	vessel	8.13E-01	9.55E-01
chr14:89895046:89895080	vessel	8.13E-01	8.91E-01
chr19:38485266:38485425	pancreas	8.13E-01	7.91E-01
chr21:47530618:47530659	vessel	8.13E-01	9.26E-01
chr17:42875273:42875297	vessel	8.13E-01	9.35E-01
chr20:50368876:50368956	vessel	8.13E-01	9.89E-01
chr7:75595987:75596208	vessel	8.13E-01	9.61E-01
chr9:80384641:80384724	neural	8.13E-01	8.16E-01
chr10:123901799:123901843	vessel	8.13E-01	9.40E-01
chr13:98082290:98082421	neural	8.13E-01	8.78E-01
chr19:33674931:33674961	vessel	8.13E-01	9.23E-01
chr16:88195380:88195408	vessel	8.13E-01	9.09E-01
chr17:8584339:8584382	vessel	8.13E-01	9.72E-01
chr11:44729854:44729898	vessel	8.13E-01	9.74E-01
chr3:129899423:129899521	heart,vessel	8.13E-01	9.71E-01
chr17:14047988:14048092	vessel	8.13E-01	9.35E-01
chr22:32976039:32976101	liver	8.13E-01	8.36E-01
chr11:57184788:57184846	vessel	8.13E-01	9.39E-01
chr22:43079771:43079865	vessel	8.13E-01	9.98E-01
chr1:32045230:32045272	vessel	8.13E-01	9.39E-01
chr2:242775475:242775494	vessel	8.13E-01	9.83E-01
chr1:2899810:2899955	vessel	8.13E-01	9.03E-01
chr22:47331564:47331659	liver	8.13E-01	8.83E-01
chr6:37046101:37046167	vessel	8.12E-01	9.11E-01
chr4:36056662:36056826	vessel	8.12E-01	9.21E-01
chr2:1566219:1566233	vessel	8.12E-01	8.87E-01
chr1:3282806:3282902	vessel	8.12E-01	9.80E-01

uMHL Markers			
Region	Group	GSI	refMax
chr12:2749661:2749704	vessel	8.12E-01	9.48E-01
chr7:73119356:73119403	vessel	8.12E-01	8.94E-01
chr21:44690743:44690783	neural	8.12E-01	8.09E-01
chr14:104018772:104018857	vessel	8.12E-01	9.82E-01
chr2:46893583:46893684	vessel	8.12E-01	9.49E-01
chr14:95875180:95875232	vessel	8.12E-01	9.74E-01
chr15:27573948:27574019	neural	8.12E-01	8.87E-01
chr20:46756899:46756983	vessel	8.12E-01	9.89E-01
chr2:46531085:46531128	vessel	8.12E-01	9.51E-01
chr1:217203173:217203214	vessel	8.12E-01	9.66E-01
chr18:71442341:71442413	vessel	8.12E-01	8.94E-01
chr7:151133987:151134097	vessel	8.12E-01	9.16E-01
chr2:144282564:144282630	vessel	8.12E-01	9.76E-01
chr3:188624320:188624396	vessel	8.12E-01	9.38E-01
chr16:4092830:4092845	vessel	8.12E-01	9.07E-01
chr9:136437866:136437930	vessel	8.12E-01	9.56E-01
chr1:61657038:61657079	vessel	8.12E-01	9.58E-01
chr8:142420181:142420195	vessel	8.12E-01	9.40E-01
chr8:69783211:69783254	vessel	8.12E-01	9.73E-01
chr17:37833823:37833848	vessel	8.12E-01	9.77E-01
chr11:12255071:12255245	vessel	8.12E-01	9.29E-01
chr15:100708207:100708246	vessel	8.12E-01	9.37E-01
chr16:51868437:51868468	vessel	8.12E-01	9.06E-01
chr17:46343933:46343974	intestine	8.12E-01	8.75E-01
chr10:30110183:30110189	vessel	8.12E-01	8.90E-01
chr8:70846160:70846211	vessel	8.12E-01	9.12E-01
chr7:132037792:132037854	vessel	8.11E-01	9.53E-01
chr16:49686581:49686695	vessel	8.11E-01	9.46E-01
chr17:30821748:30821864	neural	8.11E-01	7.47E-01
chr20:24631061:24631099	vessel	8.11E-01	9.29E-01
chr18:74172210:74172240	vessel	8.11E-01	9.23E-01
chr6:41419353:41419617	vessel	8.11E-01	9.58E-01
chr17:78002601:78002654	vessel	8.11E-01	8.84E-01
chr2:99427338:99427414	vessel	8.11E-01	9.86E-01
chr1:156271781:156271977	liver	8.11E-01	8.87E-01
chr15:56208850:56209086	vessel	8.11E-01	9.61E-01
chr19:16178499:16178673	vessel	8.11E-01	8.57E-01
chr9:126081436:126081458	vessel	8.11E-01	9.68E-01
chr5:156988956:156989239	vessel	8.11E-01	8.74E-01
chr8:23201483:23201545	vessel	8.11E-01	9.40E-01
chr17:12554165:12554199	vessel	8.11E-01	9.45E-01
chr7:2500197:2500212	neural	8.11E-01	7.78E-01
chr11:114167495:114167693	vessel	8.11E-01	9.80E-01
chr12:123579604:123579635	vessel	8.11E-01	9.18E-01
chr5:36569776:36569815	vessel	8.11E-01	9.62E-01
chr9:101721017:101721046	vessel	8.11E-01	9.44E-01
chr11:7507122:7507225	vessel	8.11E-01	9.34E-01
chr16:66737688:66737762	vessel	8.11E-01	9.54E-01

uMHL Markers			
Region	Group	GSI	refMax
chr19:46286747:46286876	vessel	8.11E-01	9.81E-01
chr2:107381544:107381755	vessel	8.11E-01	8.53E-01
chr17:74736513:74736621	vessel	8.10E-01	9.65E-01
chr13:39205072:39205118	vessel	8.10E-01	9.40E-01
chr9:92277644:92277734	neural	8.10E-01	7.89E-01
chr7:465066:465102	vessel	8.10E-01	9.47E-01
chr4:41157591:41157674	vessel	8.10E-01	9.20E-01
chr10:131694575:131694682	neural	8.10E-01	8.35E-01
chr2:109788545:109788604	vessel	8.10E-01	9.53E-01
chr9:285924:285974	vessel	8.10E-01	9.13E-01
chr1:9431472:9431636	vessel	8.10E-01	9.65E-01
chr18:77246357:77246388	vessel	8.10E-01	9.48E-01
chr16:66652537:66652773	vessel	8.10E-01	9.71E-01
chr19:3466761:3466827	vessel	8.10E-01	9.39E-01
chr2:1565960:1566095	vessel	8.10E-01	8.86E-01
chr12:123963634:123963661	vessel	8.10E-01	9.71E-01
chr13:21572143:21572171	vessel	8.10E-01	9.65E-01
chr14:52373819:52373891	vessel	8.10E-01	9.29E-01
chr1:9436026:9436361	vessel	8.10E-01	9.81E-01
chr15:100851337:100851552	vessel	8.10E-01	9.56E-01
chr10:13801406:13801461	vessel	8.10E-01	8.38E-01
chr12:47491451:47491535	vessel	8.10E-01	9.70E-01
chr12:57197704:57197748	vessel	8.10E-01	9.71E-01
chr17:77460860:77460869	neural	8.10E-01	8.69E-01
chr7:143169791:143169807	intestine	8.10E-01	7.65E-01
chr15:92462388:92462506	vessel	8.10E-01	9.12E-01
chr10:11733400:11733428	vessel	8.10E-01	9.38E-01
chr6:1867570:1867643	pancreas	8.09E-01	8.40E-01
chr13:113698994:113699050	neural	8.09E-01	7.75E-01
chr12:77624317:77624349	vessel	8.09E-01	9.63E-01
chr2:12461342:12461469	vessel	8.09E-01	9.53E-01
chr11:12352080:12352113	neural	8.09E-01	8.04E-01
chr18:72063318:72063459	vessel	8.09E-01	9.14E-01
chr21:44162140:44162220	vessel	8.09E-01	9.31E-01
chr13:99023307:99023384	vessel	8.09E-01	9.91E-01
chr1:179322741:179322768	vessel	8.09E-01	9.68E-01
chr9:126927295:126927337	vessel	8.09E-01	9.43E-01
chr6:85592864:85592982	vessel	8.09E-01	9.46E-01
chr8:106155638:106155841	neural	8.09E-01	8.65E-01
chr13:23951893:23951910	vessel	8.09E-01	9.70E-01
chr22:27501587:27501648	vessel	8.09E-01	8.66E-01
chr1:7728694:7728716	neural	8.09E-01	7.96E-01
chr9:137554083:137554352	vessel	8.09E-01	9.53E-01
chr18:53074062:53074313	neural	8.09E-01	8.70E-01
chr16:72892761:72892801	vessel	8.09E-01	9.49E-01
chr10:44339266:44339451	vessel	8.09E-01	9.17E-01
chr10:126724009:126724099	vessel	8.09E-01	9.27E-01
chr9:73204687:73204768	vessel	8.09E-01	9.63E-01

uMHL Markers			
Region	Group	GSI	refMax
chr22:31440272:31440362	vessel	8.09E-01	9.83E-01
chr20:44457872:44457929	vessel	8.09E-01	8.28E-01
chr16:79751767:79751800	vessel	8.09E-01	8.38E-01
chr1:22658835:22658882	vessel	8.09E-01	8.65E-01
chr2:8389473:8389584	vessel	8.09E-01	9.80E-01
chr7:98719280:98719357	vessel	8.09E-01	9.45E-01
chr1:203526708:203526778	vessel	8.09E-01	9.86E-01
chr5:53742361:53742402	vessel	8.09E-01	9.19E-01
chr2:1820467:1820494	vessel	8.08E-01	8.77E-01
chr15:69841259:69841298	vessel	8.08E-01	9.39E-01
chr20:61013114:61013193	vessel	8.08E-01	9.26E-01
chr15:76470078:76470143	vessel	8.08E-01	9.68E-01
chr7:70480615:70480687	vessel	8.08E-01	8.88E-01
chr9:116362067:116362205	vessel	8.08E-01	8.40E-01
chr3:48732676:48732782	neural	8.08E-01	8.43E-01
chr10:103406210:103406384	vessel	8.08E-01	9.45E-01
chr7:30885459:30885500	vessel	8.08E-01	9.82E-01
chr15:92509445:92509609	vessel	8.08E-01	9.34E-01
chr15:84361392:84361425	vessel	8.08E-01	9.40E-01
chr2:238322190:238322513	vessel	8.08E-01	9.58E-01
chr3:196910630:196910640	neural	8.08E-01	8.13E-01
chr20:50187210:50187265	vessel	8.08E-01	9.67E-01
chr2:48710712:48710871	vessel	8.08E-01	9.51E-01
chr18:34366845:34366901	vessel	8.08E-01	9.40E-01
chr21:43147530:43147588	vessel	8.08E-01	9.70E-01
chr11:36037604:36037621	vessel	8.08E-01	9.24E-01
chr6:160257580:160257599	vessel	8.08E-01	9.58E-01
chr15:66914008:66914032	vessel	8.08E-01	9.67E-01
chr10:134881126:134881168	vessel	8.08E-01	9.23E-01
chr7:7032037:7032089	intestine	8.08E-01	9.17E-01
chr10:13881755:13881845	vessel	8.08E-01	9.46E-01
chr6:11214781:11215022	vessel	8.08E-01	9.81E-01
chr14:76446554:76446750	vessel	8.08E-01	9.56E-01
chr12:124004371:124004416	vessel	8.08E-01	9.48E-01
chr1:112400108:112400129	neural	8.08E-01	8.33E-01
chr15:74496556:74496593	vessel	8.07E-01	9.46E-01
chr10:126749370:126749418	vessel	8.07E-01	9.55E-01
chr13:49528296:49528490	vessel	8.07E-01	9.66E-01
chr7:135440322:135440342	vessel	8.07E-01	8.04E-01
chr10:16869686:16869746	vessel	8.07E-01	9.52E-01
chr21:32532775:32532824	neural	8.07E-01	8.82E-01
chr9:97682657:97682800	vessel	8.07E-01	9.53E-01
chr20:49982779:49982992	vessel	8.07E-01	9.52E-01
chr16:19872075:19872241	vessel	8.07E-01	9.51E-01
chr9:92300116:92300178	vessel	8.07E-01	9.54E-01
chr6:37045638:37045671	vessel	8.07E-01	9.71E-01
chr1:225838270:225838425	vessel	8.07E-01	9.04E-01
chr2:5701195:5701296	vessel	8.07E-01	9.68E-01

uMHL Markers			
Region	Group	GSI	refMax
chr1:19740620:19740711	vessel	8.07E-01	9.07E-01
chr21:44161774:44161828	vessel	8.07E-01	9.53E-01
chr6:143434808:143435073	liver	8.07E-01	8.80E-01
chr10:121040981:121041031	vessel	8.07E-01	9.32E-01
chr6:154795868:154796097	vessel	8.07E-01	8.55E-01
chr2:233105800:233105832	kidney	8.07E-01	7.86E-01
chr1:2006023:2006038	neural	8.07E-01	7.73E-01
chr1:241912764:241912999	vessel	8.07E-01	8.86E-01
chr1:80251396:80251438	vessel	8.07E-01	9.04E-01
chr8:89270337:89270376	vessel	8.07E-01	9.38E-01
chr8:22458405:22458469	vessel	8.07E-01	8.77E-01
chr5:107754120:107754132	vessel	8.07E-01	9.57E-01
chr3:194493573:194493711	vessel	8.07E-01	9.55E-01
chr5:158224572:158224690	vessel	8.07E-01	9.82E-01
chr2:20262560:20262730	vessel	8.07E-01	8.69E-01
chr19:3460058:3460132	vessel	8.07E-01	9.41E-01
chr9:92503363:92503383	vessel	8.07E-01	9.51E-01
chr18:46375226:46375356	vessel	8.07E-01	9.01E-01
chr20:17850931:17851022	vessel	8.06E-01	9.04E-01
chr19:2186543:2186661	vessel	8.06E-01	9.61E-01
chr9:97682842:97683006	vessel	8.06E-01	9.71E-01
chr6:158460945:158461070	vessel	8.06E-01	9.43E-01
chr9:2816879:2817023	vessel	8.06E-01	9.76E-01
chr1:183187558:183187637	vessel	8.06E-01	9.47E-01
chr12:121344919:121344964	vessel	8.06E-01	9.67E-01
chr20:35169120:35169198	vessel	8.06E-01	9.00E-01
chr15:86185647:86185814	vessel	8.06E-01	9.46E-01
chr16:84617788:84617843	vessel	8.06E-01	9.38E-01
chr17:43926509:43926525	vessel	8.06E-01	9.65E-01
chr8:30420219:30420335	vessel	8.06E-01	9.57E-01
chr6:16541787:16541933	vessel	8.06E-01	9.64E-01
chr17:78999585:78999726	neural	8.06E-01	8.25E-01
chr10:16759844:16759898	vessel	8.06E-01	9.28E-01
chr14:95875024:95875083	vessel	8.06E-01	9.93E-01
chr21:34192080:34192208	vessel	8.06E-01	9.26E-01
chr21:46454122:46454165	vessel	8.06E-01	9.39E-01
chr18:13436078:13436235	vessel	8.06E-01	9.80E-01
chr10:130094981:130095001	vessel	8.06E-01	9.48E-01
chr6:169178590:169178687	vessel	8.06E-01	9.89E-01
chr15:58815128:58815195	vessel	8.06E-01	9.20E-01
chr2:67487963:67487997	vessel	8.06E-01	8.92E-01
chr18:72913659:72913731	vessel	8.06E-01	9.16E-01
chr22:23624534:23624645	vessel	8.06E-01	9.44E-01
chr12:2494075:2494199	vessel	8.06E-01	9.18E-01
chr17:55663160:55663251	vessel	8.06E-01	8.71E-01
chr2:65543343:65543411	vessel	8.06E-01	9.81E-01
chr20:56305784:56305798	vessel	8.06E-01	9.53E-01
chr18:43097228:43097251	vessel	8.06E-01	9.86E-01

uMHL Markers			
Region	Group	GSI	refMax
chr9:98829385:98829459	vessel	8.06E-01	9.00E-01
chr9:98812048:98812090	vessel	8.06E-01	9.27E-01
chr12:93521229:93521358	vessel	8.06E-01	9.68E-01
chr10:3146921:3146961	vessel	8.06E-01	9.36E-01
chr19:28609319:28609354	vessel	8.06E-01	8.70E-01
chr13:113698901:113698952	neural	8.06E-01	7.96E-01
chr2:216298546:216298643	vessel	8.05E-01	9.78E-01
chr1:2991642:2991733	vessel	8.05E-01	9.61E-01
chr17:76858155:76858303	vessel	8.05E-01	9.78E-01
chr17:19390397:19390458	vessel	8.05E-01	9.75E-01
chr4:1205974:1206025	vessel	8.05E-01	9.28E-01
chr11:11606063:11606090	vessel	8.05E-01	9.75E-01
chr7:22757575:22757624	vessel	8.05E-01	9.70E-01
chr11:19752027:19752292	vessel	8.05E-01	9.37E-01
chr9:97587850:97588035	vessel	8.05E-01	9.29E-01
chr2:375931:375958	vessel	8.05E-01	9.16E-01
chr17:71773282:71773397	vessel	8.05E-01	9.61E-01
chr12:128807935:128807959	vessel	8.05E-01	9.65E-01
chr13:23732081:23732140	vessel	8.05E-01	9.48E-01
chr3:59003065:59003100	liver	8.05E-01	8.20E-01
chr1:157984204:157984269	vessel	8.05E-01	9.40E-01
chr7:55112569:55112635	vessel	8.05E-01	9.85E-01
chr20:17830714:17830792	vessel	8.05E-01	9.31E-01
chr4:95331866:95332202	vessel	8.05E-01	9.51E-01
chr3:122692837:122692868	vessel	8.05E-01	9.03E-01
chr16:2177303:2177329	neural,vessel	8.05E-01	1.00E+00
chr15:85826138:85826201	vessel	8.05E-01	8.76E-01
chr14:76461910:76462138	vessel	8.05E-01	9.54E-01
chr5:155755991:155756019	vessel	8.05E-01	8.44E-01
chr4:57922136:57922184	vessel	8.05E-01	9.57E-01
chr9:15401524:15401532	vessel	8.05E-01	9.68E-01
chr19:3603394:3603412	vessel	8.05E-01	9.17E-01
chr19:47316636:47316653	vessel	8.05E-01	8.82E-01
chr11:76360139:76360318	vessel	8.05E-01	8.86E-01
chr10:729913:729968	vessel	8.05E-01	9.23E-01
chr11:114003718:114003831	vessel	8.05E-01	9.42E-01
chr17:75775848:75775881	vessel	8.05E-01	9.53E-01
chr12:116737894:116737988	vessel	8.05E-01	9.32E-01
chr21:44694448:44694461	liver	8.04E-01	8.95E-01
chr5:167372652:167372700	vessel	8.04E-01	9.18E-01
chr10:61408493:61408669	vessel	8.04E-01	9.36E-01
chr5:156989353:156989478	vessel	8.04E-01	9.34E-01
chr8:97419822:97419867	vessel	8.04E-01	9.02E-01
chr8:22516858:22516919	vessel	8.04E-01	9.54E-01
chr6:38237680:38237861	vessel	8.04E-01	9.53E-01
chr1:37773917:37773943	neural	8.04E-01	7.70E-01
chr7:95818174:95818283	vessel	8.04E-01	9.66E-01
chr7:73428567:73428618	vessel	8.04E-01	9.58E-01

-47-

uMHL Markers			
Region	Group	GSI	refMax
chr13:114632959:114633003	vessel	8.04E-01	9.20E-01
chr1:59912018:59912049	vessel	8.04E-01	9.35E-01
chr2:109941286:109941359	vessel	8.04E-01	8.91E-01
chr1:168054694:168054754	vessel	8.04E-01	9.53E-01
chr3:69061666:69061765	vessel	8.04E-01	9.66E-01
chr6:2382649:2382767	vessel	8.04E-01	9.59E-01
chr12:63124047:63124109	vessel	8.04E-01	9.21E-01
chr14:73025645:73025673	vessel	8.04E-01	9.50E-01
chr15:25813912:25814131	vessel	8.04E-01	9.46E-01
chr21:40044046:40044198	vessel	8.04E-01	9.52E-01
chr19:38526875:38526924	vessel	8.04E-01	9.34E-01
chr2:67488351:67488376	vessel	8.04E-01	9.53E-01
chr17:117230:117455	vessel	8.04E-01	9.76E-01
chr4:3703023:3703043	vessel	8.04E-01	9.47E-01
chr18:73882379:73882586	intestine	8.04E-01	7.71E-01
chr16:77465342:77465471	vessel	8.04E-01	8.69E-01
chr14:91750605:91750669	vessel	8.04E-01	9.33E-01
chr11:44028454:44028472	vessel	8.04E-01	9.13E-01
chr2:9377096:9377111	vessel	8.04E-01	9.67E-01
chr20:62166877:62166882	liver	8.03E-01	9.93E-01
chr15:81384367:81384382	kidney	8.03E-01	7.88E-01
chr12:106246647:106246707	vessel	8.03E-01	9.89E-01
chr2:79807509:79807533	vessel	8.03E-01	8.49E-01
chr7:39793986:39794185	vessel	8.03E-01	9.14E-01
chr17:12553963:12554077	vessel	8.03E-01	9.40E-01
chr3:194961215:194961310	vessel	8.03E-01	9.57E-01
chr5:176884349:176884370	vessel	8.03E-01	9.77E-01
chr4:24732341:24732420	vessel	8.03E-01	9.63E-01
chr22:20186181:20186409	intestine	8.03E-01	7.68E-01
chr3:123397265:123397302	vessel	8.03E-01	9.79E-01
chr1:212734867:212734967	vessel	8.03E-01	9.75E-01
chr9:89882438:89882551	vessel	8.03E-01	9.42E-01
chr6:150848318:150848392	vessel	8.03E-01	9.59E-01
chr15:61073748:61073811	vessel	8.03E-01	8.65E-01
chr3:126720886:126721050	vessel	8.03E-01	9.77E-01
chr21:25397216:25397285	vessel	8.03E-01	9.59E-01
chr2:8143877:8143939	vessel	8.03E-01	9.97E-01
chr16:82336047:82336082	vessel	8.03E-01	9.42E-01
chr9:738321:738429	vessel	8.03E-01	9.02E-01
chr5:77843468:77843772	vessel	8.03E-01	9.25E-01
chr5:65165845:65165992	intestine	8.03E-01	8.30E-01
chr16:1037729:1037750	vessel	8.03E-01	8.66E-01
chr9:36566815:36566887	vessel	8.03E-01	9.76E-01
chr1:2900163:2900177	vessel	8.03E-01	9.39E-01
chr2:218805147:218805216	vessel	8.03E-01	9.53E-01
chr16:88976144:88976169	vessel	8.03E-01	9.71E-01
chr1:32220068:32220206	neural	8.03E-01	8.26E-01
chr20:49983123:49983144	vessel	8.03E-01	9.11E-01

uMHL Markers			
Region	Group	GSI	refMax
chr7:30947529:30947552	vessel	8.03E-01	9.48E-01
chr2:97572877:97572933	vessel	8.03E-01	9.59E-01
chr5:111890374:111890446	vessel	8.03E-01	9.57E-01
chr10:50346525:50346634	vessel	8.02E-01	8.45E-01
chr10:100151545:100151593	vessel	8.02E-01	9.76E-01
chr20:48651408:48651539	vessel	8.02E-01	9.53E-01
chr5:138440901:138441072	vessel	8.02E-01	9.50E-01
chr13:40090279:40090501	vessel	8.02E-01	9.20E-01
chr8:22516977:22517039	vessel	8.02E-01	9.26E-01
chr12:125332183:125332213	neural	8.02E-01	8.43E-01
chr8:30658350:30658451	vessel	8.02E-01	8.37E-01
chr10:10312643:10312797	vessel	8.02E-01	8.86E-01
chr22:44484948:44485088	vessel	8.02E-01	9.42E-01
chr2:138721303:138721318	liver	8.02E-01	9.18E-01
chr16:16161126:16161243	vessel	8.02E-01	9.40E-01
chr2:97019386:97019473	vessel	8.02E-01	9.83E-01
chr20:24568297:24568342	neural	8.02E-01	7.78E-01
chr19:15360079:15360255	vessel	8.02E-01	9.06E-01
chr22:44756580:44756669	vessel	8.02E-01	9.47E-01
chr22:29726826:29726886	kidney	8.02E-01	9.21E-01
chr18:59685565:59685607	vessel	8.02E-01	9.01E-01
chr5:168477679:168477730	vessel	8.02E-01	9.07E-01
chr20:50187002:50187064	vessel	8.02E-01	9.17E-01
chr18:46468631:46468643	vessel	8.02E-01	9.55E-01
chr4:159011359:159011402	vessel	8.02E-01	9.82E-01
chr1:3446617:3446738	vessel	8.01E-01	9.87E-01
chr5:170910865:170910881	vessel	8.01E-01	9.44E-01
chr5:101493402:101493492	vessel	8.01E-01	8.57E-01
chr19:18782908:18783006	vessel	8.01E-01	8.00E-01
chr4:169765695:169765832	vessel	8.01E-01	9.79E-01
chr11:12182856:12182943	vessel	8.01E-01	9.38E-01
chr3:98988954:98989145	vessel	8.01E-01	9.72E-01
chr2:28958706:28958762	vessel	8.01E-01	9.29E-01
chr17:77069737:77069835	vessel	8.01E-01	9.47E-01
chr15:101912443:101912517	vessel	8.01E-01	9.09E-01
chr12:124959289:124959504	vessel	8.01E-01	9.17E-01
chr11:122205905:122205985	vessel	8.01E-01	8.54E-01
chr21:39903689:39903803	vessel	8.01E-01	9.42E-01
chr11:75269142:75269197	vessel	8.01E-01	9.60E-01
chr20:10616067:10616113	vessel	8.01E-01	9.57E-01
chr9:73539583:73539651	vessel	8.01E-01	9.47E-01
chr6:12349194:12349298	vessel	8.01E-01	9.30E-01
chr7:6291930:6292069	vessel	8.01E-01	9.02E-01
chr2:1820169:1820243	vessel	8.01E-01	8.40E-01
chr8:25936821:25936871	vessel	8.01E-01	9.56E-01
chr11:10715379:10715450	vessel	8.01E-01	9.66E-01
chr20:35169310:35169362	vessel	8.01E-01	9.68E-01
chr9:127047244:127047277	vessel	8.01E-01	9.61E-01

uMHL Markers			
Region	Group	GSI	refMax
chr12:124925719:124925860	vessel	8.01E-01	9.17E-01
chr6:157039978:157040028	vessel	8.01E-01	8.63E-01
chr12:54071769:54071918	neural	8.01E-01	8.47E-01
chr16:4001865:4002028	vessel	8.00E-01	9.70E-01
chr18:13542218:13542473	liver	8.00E-01	8.54E-01
chr3:71213389:71213444	vessel	8.00E-01	8.33E-01
chr14:91622695:91622754	liver	8.00E-01	8.73E-01
chr5:78115075:78115276	vessel	8.00E-01	9.36E-01
chr22:30737410:30737507	vessel	8.00E-01	9.69E-01
chr2:109788226:109788430	vessel	8.00E-01	9.26E-01
chr9:112841248:112841282	vessel	8.00E-01	9.45E-01
chr5:95194167:95194199	vessel	8.00E-01	9.88E-01
chr11:113513456:113513577	vessel	8.00E-01	8.11E-01
chr7:41956427:41956582	vessel	8.00E-01	9.28E-01
chr10:31108592:31108620	vessel	8.00E-01	9.61E-01
chr17:1038460:1038480	vessel	8.00E-01	9.74E-01
chr5:158230112:158230269	vessel	8.00E-01	9.67E-01
chr11:75444671:75444686	vessel	8.00E-01	8.50E-01
chr12:116746846:116746866	vessel	8.00E-01	9.52E-01
chr9:97535350:97535437	vessel	8.00E-01	9.56E-01
chr18:10376914:10377078	vessel	8.00E-01	9.25E-01
chr11:1552850:1552939	vessel	8.00E-01	9.73E-01
chr18:76551122:76551148	vessel	8.00E-01	8.19E-01
chr9:133734023:133734044	vessel	8.00E-01	9.43E-01
chr5:136955094:136955232	vessel	8.00E-01	9.90E-01
chr19:42528857:42528937	vessel	8.00E-01	9.74E-01
chr17:77993088:77993180	liver	8.00E-01	8.60E-01
chr9:139466646:139466672	vessel	8.00E-01	9.51E-01
chr4:129440877:129440895	vessel	8.00E-01	9.07E-01
chr2:135118406:135118696	vessel	8.00E-01	9.57E-01
chr8:1895423:1895593	vessel	8.00E-01	8.96E-01
chr16:89900228:89900472	intestine	8.00E-01	8.66E-01
chr18:53533397:53533431	vessel	8.00E-01	8.91E-01
chr8:1734515:1734583	vessel	8.00E-01	9.27E-01
chr5:146783414:146783535	vessel	8.00E-01	9.21E-01
chr10:129721163:129721395	vessel	8.00E-01	9.26E-01
chr16:88976512:88976545	vessel	8.00E-01	8.82E-01
chr17:12534611:12534640	vessel	8.00E-01	8.99E-01
chr8:124552202:124552287	vessel	8.00E-01	9.68E-01
chr5:4983978:4984033	vessel	8.00E-01	8.65E-01
chr4:2589345:2589376	vessel	8.00E-01	9.62E-01
chr14:23318984:23319004	vessel	8.00E-01	9.39E-01
chr3:125978692:125978717	vessel	8.00E-01	9.38E-01
chr7:137668500:137668714	vessel	7.99E-01	9.45E-01
chr22:29707734:29707797	neural	7.99E-01	8.65E-01
chr6:168811258:168811288	vessel	7.99E-01	8.85E-01
chr11:48086058:48086140	vessel	7.99E-01	9.48E-01
chr1:32581854:32581901	vessel	7.99E-01	9.52E-01

uMHL Markers			
Region	Group	GSI	refMax
chr10:63817039:63817096	vessel	7.99E-01	9.83E-01
chr10:104679620:104679676	vessel	7.99E-01	9.47E-01
chr18:72913820:72914028	vessel	7.99E-01	8.76E-01
chr18:8381961:8382108	vessel	7.99E-01	9.32E-01
chr2:43691730:43692049	vessel	7.99E-01	9.77E-01
chr16:66260739:66260808	vessel	7.99E-01	8.07E-01
chr1:3314034:3314111	vessel	7.99E-01	8.05E-01
chr15:47732211:47732323	vessel	7.99E-01	9.53E-01
chr16:88228206:88228258	intestine	7.99E-01	8.48E-01
chr7:75595743:75595885	vessel	7.99E-01	9.21E-01
chr17:77302405:77302421	vessel	7.99E-01	9.62E-01
chr7:133043576:133043707	vessel	7.99E-01	9.54E-01
chr22:31640170:31640177	vessel	7.99E-01	9.54E-01
chr1:25916675:25916756	vessel	7.99E-01	9.21E-01
chr1:184120477:184120508	vessel	7.99E-01	9.43E-01
chr5:131713182:131713241	vessel	7.99E-01	9.30E-01
chr14:75019600:75019659	vessel	7.99E-01	9.45E-01
chr16:16088347:16088528	vessel	7.99E-01	9.66E-01
chr7:1135687:1135764	vessel	7.99E-01	9.33E-01
chr2:98775072:98775082	vessel	7.99E-01	9.37E-01
chr10:5605504:5605530	vessel	7.99E-01	9.39E-01
chr14:92391559:92391578	vessel	7.99E-01	9.16E-01
chr8:20834957:20835064	vessel	7.99E-01	9.36E-01
chr7:73414105:73414180	vessel	7.99E-01	1.00E+00
chr3:114225764:114225848	vessel	7.99E-01	9.48E-01
chr19:47242589:47242709	vessel	7.99E-01	9.22E-01
chr3:41042033:41042081	vessel	7.99E-01	9.30E-01
chr19:57429326:57429374	vessel	7.98E-01	9.44E-01
chr7:41735421:41735459	vessel	7.98E-01	9.70E-01
chr11:10886639:10886702	vessel	7.98E-01	9.58E-01
chr10:664325:664473	vessel	7.98E-01	9.17E-01
chr5:38469431:38469488	vessel	7.98E-01	9.51E-01
chr16:15237718:15237812	vessel	7.98E-01	9.86E-01
chr19:4571622:4571663	vessel	7.98E-01	9.20E-01
chr7:139458285:139458300	liver	7.98E-01	8.65E-01
chr9:116967656:116967726	vessel	7.98E-01	9.49E-01
chr16:19656482:19656576	vessel	7.98E-01	9.59E-01
chr16:12865231:12865387	vessel	7.98E-01	9.77E-01
chr6:3766566:3766739	vessel	7.98E-01	9.28E-01
chr15:60705703:60705745	vessel	7.98E-01	9.76E-01
chr1:85070440:85070481	vessel	7.98E-01	8.60E-01
chr13:99622769:99622853	heart,vessel	7.98E-01	9.55E-01
chr2:197962437:197962594	vessel	7.98E-01	9.66E-01
chr8:70180168:70180286	vessel	7.98E-01	8.72E-01
chr7:141768483:141768503	vessel	7.98E-01	9.25E-01
chr10:102642706:102642787	vessel	7.98E-01	9.89E-01
chr11:19837610:19837706	vessel	7.98E-01	8.96E-01
chr16:2176909:2176936	vessel	7.98E-01	9.19E-01

uMHL Markers			
Region	Group	GSI	refMax
chr15:99414209:99414271	vessel	7.98E-01	8.64E-01
chr1:3459917:3459985	vessel	7.98E-01	9.31E-01
chr17:60830341:60830394	vessel	7.98E-01	8.74E-01
chr13:53608851:53608955	vessel	7.98E-01	9.19E-01
chr7:143504126:143504224	vessel	7.97E-01	9.07E-01
chr14:93530866:93530975	vessel	7.97E-01	9.50E-01
chr2:9375605:9375853	vessel	7.97E-01	8.80E-01
chr17:2870919:2870940	vessel	7.97E-01	9.13E-01
chr5:158018853:158018897	vessel	7.97E-01	9.54E-01
chr13:112271255:112271317	vessel	7.97E-01	9.54E-01
chr11:3183785:3183815	neural,vessel	7.97E-01	9.50E-01
chr12:131430753:131430967	vessel	7.97E-01	9.64E-01
chr14:80052309:80052328	vessel	7.97E-01	9.26E-01
chr16:79943780:79943836	vessel	7.97E-01	9.17E-01
chr17:55498598:55498641	stomach	7.97E-01	8.66E-01
chr4:39530807:39530904	liver	7.97E-01	8.09E-01
chr19:35629398:35629459	vessel	7.97E-01	9.85E-01
chr10:30743560:30743638	vessel	7.97E-01	9.73E-01
chr2:1656489:1656558	vessel	7.97E-01	9.44E-01
chr10:65831088:65831186	vessel	7.97E-01	9.41E-01
chr21:35718410:35718483	vessel	7.97E-01	9.56E-01
chr1:3474045:3474120	vessel	7.97E-01	9.15E-01
chr18:77246010:77246208	vessel	7.97E-01	9.92E-01
chr2:234390268:234390333	vessel	7.97E-01	8.79E-01
chr7:73309532:73309575	vessel	7.97E-01	7.92E-01
chr1:183013562:183013626	vessel	7.97E-01	9.03E-01
chr16:70725685:70726085	vessel	7.97E-01	9.58E-01
chr1:154401741:154401879	neural	7.97E-01	7.74E-01
chr4:36166725:36166900	neural	7.97E-01	8.20E-01
chr9:96366610:96366664	vessel	7.97E-01	7.96E-01
chr9:95371219:95371367	vessel	7.97E-01	9.70E-01
chr10:13999112:13999167	vessel	7.97E-01	9.68E-01
chr9:137535735:137535830	vessel	7.97E-01	9.05E-01
chr7:94025970:94026143	vessel	7.96E-01	9.60E-01
chr20:46475609:46475683	vessel	7.96E-01	9.58E-01
chr12:94580545:94580701	neural	7.96E-01	8.90E-01
chr1:160877880:160877935	vessel	7.96E-01	9.71E-01
chr10:104380407:104380589	vessel	7.96E-01	9.58E-01
chr22:28191297:28191348	vessel	7.96E-01	8.51E-01
chr4:99981639:99981676	vessel	7.96E-01	9.38E-01
chr12:58287209:58287324	vessel	7.96E-01	9.50E-01
chr6:168975596:168975609	vessel	7.96E-01	9.74E-01
chr14:104178824:104178840	vessel	7.96E-01	9.66E-01
chr9:87051327:87051398	vessel	7.96E-01	8.95E-01
chr11:4212243:4212284	vessel	7.96E-01	9.59E-01
chr5:125720663:125720715	vessel	7.96E-01	8.48E-01
chr21:33570835:33570867	vessel	7.96E-01	8.36E-01
chr7:1578135:1578254	vessel	7.96E-01	9.53E-01

-52-

uMHL Markers			
Region	Group	GSI	refMax
chr16:15239066:15239131	vessel	7.96E-01	9.12E-01
chr4:119711755:119711864	vessel	7.96E-01	9.04E-01
chr21:40459773:40459854	vessel	7.96E-01	9.31E-01
chr11:61655771:61655827	vessel	7.96E-01	9.60E-01
chr1:2286708:2286716	vessel	7.96E-01	8.17E-01
chr6:3749643:3749718	vessel	7.96E-01	8.99E-01
chr5:71805380:71805534	vessel	7.96E-01	9.34E-01
chr8:88927532:88927575	vessel	7.96E-01	9.48E-01
chr1:24437962:24437971	vessel	7.96E-01	9.54E-01
chr5:158232314:158232393	vessel	7.96E-01	9.37E-01
chr12:52938468:52938505	vessel	7.96E-01	9.35E-01
chr3:101230854:101230960	intestine	7.96E-01	7.53E-01
chr7:156160942:156161041	liver	7.96E-01	8.01E-01
chr1:32054262:32054281	vessel	7.96E-01	9.52E-01
chr16:3843446:3843552	vessel	7.96E-01	9.69E-01
chr15:41316602:41316676	vessel	7.96E-01	9.72E-01
chr16:82668011:82668094	vessel	7.96E-01	9.75E-01
chr9:82195114:82195183	pancreas	7.96E-01	7.64E-01
chr10:31108828:31108849	vessel	7.96E-01	9.33E-01
chr10:121006638:121006679	vessel	7.96E-01	8.31E-01
chr2:240421723:240421775	vessel	7.96E-01	9.41E-01
chr10:3616894:3616970	heart	7.96E-01	8.22E-01
chr10:17275675:17275857	vessel	7.96E-01	9.73E-01
chr3:105448599:105448786	vessel	7.96E-01	9.54E-01
chr2:12597916:12598117	vessel	7.96E-01	9.03E-01
chr8:49059786:49059859	vessel	7.95E-01	9.46E-01
chr11:19952720:19952784	vessel	7.95E-01	8.63E-01
chr2:129936106:129936152	vessel	7.95E-01	9.60E-01
chr20:10620456:10620662	neural,vessel	7.95E-01	9.38E-01
chr5:176515821:176515829	kidney	7.95E-01	9.03E-01
chr11:93869233:93869339	vessel	7.95E-01	9.31E-01
chr7:30951682:30951741	vessel	7.95E-01	9.44E-01
chr7:70096103:70096122	vessel	7.95E-01	9.06E-01
chr9:92716323:92716482	vessel	7.95E-01	9.39E-01
chr7:111599315:111599369	vessel	7.95E-01	8.85E-01
chr2:27136976:27137029	vessel	7.95E-01	9.26E-01
chr7:608326:608527	vessel	7.95E-01	9.51E-01
chr9:97684443:97684562	vessel	7.95E-01	9.63E-01
chr5:66944023:66944179	vessel	7.95E-01	9.38E-01
chr2:190701403:190701607	vessel	7.95E-01	8.86E-01
chr22:40673658:40673728	vessel	7.95E-01	9.55E-01
chr14:92366261:92366344	vessel	7.95E-01	9.22E-01
chr7:75640653:75640744	vessel	7.95E-01	9.11E-01
chr16:73059129:73059169	vessel	7.95E-01	9.09E-01
chr9:129155111:129155328	vessel	7.95E-01	9.61E-01
chr2:183292642:183292828	vessel	7.95E-01	9.59E-01
chr5:172123630:172123730	vessel	7.95E-01	9.93E-01
chr3:47613082:47613088	neural	7.95E-01	8.23E-01

uMHL Markers			
Region	Group	GSI	refMax
chr10:105408583:105408600	vessel	7.95E-01	9.13E-01
chr5:2387936:2388031	vessel	7.95E-01	8.59E-01
chr5:2729492:2729508	vessel	7.95E-01	9.36E-01
chr20:41766979:41766991	vessel	7.95E-01	9.07E-01
chr3:184409648:184409750	vessel	7.95E-01	9.11E-01
chr10:88493041:88493164	vessel	7.95E-01	9.25E-01
chr6:132270440:132270682	vessel	7.95E-01	9.80E-01
chr22:44756989:44757220	vessel	7.94E-01	9.13E-01
chr14:68831611:68831749	vessel	7.94E-01	8.91E-01
chr14:92367927:92368071	vessel	7.94E-01	9.11E-01
chr12:2411416:2411598	vessel	7.94E-01	9.24E-01
chr22:18049263:18049365	liver	7.94E-01	8.26E-01
chr10:14065428:14065473	vessel	7.94E-01	8.83E-01
chr7:111977285:111977345	vessel	7.94E-01	9.65E-01
chr1:27445356:27445391	vessel	7.94E-01	9.69E-01
chr1:7074131:7074170	vessel	7.94E-01	9.38E-01
chr17:15313380:15313443	vessel	7.94E-01	9.62E-01
chr12:102820870:102820961	vessel	7.94E-01	9.04E-01
chr17:71936239:71936248	vessel	7.94E-01	8.45E-01
chr13:22665786:22665867	vessel	7.94E-01	9.02E-01
chr16:85187520:85187614	vessel	7.94E-01	9.75E-01
chr8:6442958:6443017	vessel	7.94E-01	9.11E-01
chr14:91107941:91108028	vessel	7.94E-01	9.90E-01
chr8:25806426:25806443	vessel	7.94E-01	9.64E-01
chr16:85620451:85620532	vessel	7.94E-01	8.04E-01
chr21:43951304:43951331	neural	7.94E-01	7.97E-01
chr12:94581206:94581313	neural	7.94E-01	8.38E-01
chr10:134527606:134527777	neural	7.94E-01	7.96E-01
chr22:26143790:26143817	vessel	7.94E-01	8.69E-01
chr7:4296052:4296381	vessel	7.94E-01	8.87E-01
chr3:134034372:134034636	vessel	7.94E-01	9.45E-01
chr10:28111467:28111604	vessel	7.94E-01	8.67E-01
chr17:79391325:79391366	vessel	7.94E-01	9.35E-01
chr14:96520629:96520665	vessel	7.94E-01	9.12E-01
chr2:18699763:18699825	vessel	7.94E-01	9.71E-01
chr17:925744:925782	vessel	7.94E-01	8.75E-01
chr21:17442467:17442664	neural	7.94E-01	8.12E-01
chr14:90042589:90042653	vessel	7.94E-01	9.12E-01
chr6:158509133:158509238	vessel	7.94E-01	9.58E-01
chr8:141629654:141629773	kidney	7.94E-01	8.44E-01
chr16:85196321:85196364	vessel	7.94E-01	9.71E-01
chr1:203525849:203525916	vessel	7.94E-01	9.61E-01
chr16:521619:521899	neural	7.94E-01	8.59E-01
chr17:38612725:38612772	vessel	7.94E-01	9.17E-01
chr15:81605401:81605471	vessel	7.94E-01	9.10E-01
chr9:135994555:135994742	neural	7.93E-01	9.64E-01
chr17:79543666:79544106	vessel	7.93E-01	9.50E-01
chr17:79012377:79012405	vessel	7.93E-01	9.70E-01

uMHL Markers			
Region	Group	GSI	refMax
chr21:45968533:45968547	vessel	7.93E-01	9.44E-01
chr21:42619547:42619640	vessel	7.93E-01	9.69E-01
chr20:19830589:19830718	vessel	7.93E-01	9.76E-01
chr10:14712906:14712935	neural	7.93E-01	8.81E-01
chr2:145465753:145465804	vessel	7.93E-01	8.25E-01
chr5:59504641:59504726	vessel	7.93E-01	9.03E-01
chr6:4613067:4613243	vessel	7.93E-01	9.23E-01
chr18:11269160:11269181	vessel	7.93E-01	9.54E-01
chr9:84113892:84113906	vessel	7.93E-01	9.48E-01
chr16:85418413:85418592	intestine	7.93E-01	8.38E-01
chr20:17817173:17817268	vessel	7.93E-01	9.43E-01
chr18:77398686:77398709	vessel	7.93E-01	1.00E+00
chr1:2005083:2005710	neural	7.93E-01	9.11E-01
chr2:44502707:44502779	intestine	7.93E-01	7.91E-01
chr9:132743907:132744043	liver	7.93E-01	8.91E-01
chr8:30343080:30343196	vessel	7.93E-01	9.10E-01
chr14:56056039:56056048	vessel	7.93E-01	9.31E-01
chr1:6445478:6445510	vessel	7.93E-01	9.52E-01
chr1:10168307:10168493	vessel	7.93E-01	9.08E-01
chr3:11645230:11645368	vessel	7.93E-01	9.65E-01
chr4:13415653:13415732	vessel	7.93E-01	9.51E-01
chr7:33732757:33732843	vessel	7.93E-01	9.04E-01
chr3:9289465:9289502	neural	7.93E-01	7.70E-01
chr14:76445973:76446021	vessel	7.93E-01	9.92E-01
chr3:72640590:72640767	vessel	7.93E-01	9.38E-01
chr7:67319367:67319425	neural	7.93E-01	7.89E-01
chr2:159950921:159951154	vessel	7.92E-01	9.46E-01
chr2:10545557:10545645	vessel	7.92E-01	9.72E-01
chr13:110874293:110874360	vessel	7.92E-01	9.51E-01
chr10:74677129:74677157	vessel	7.92E-01	9.23E-01
chr12:124948015:124948049	neural	7.92E-01	7.85E-01
chr17:72445156:72445231	vessel	7.92E-01	9.75E-01
chr15:26076118:26076279	vessel	7.92E-01	9.71E-01
chr7:100761027:100761057	vessel	7.92E-01	9.50E-01
chr10:5587760:5587894	vessel	7.92E-01	9.11E-01
chr2:217544422:217544603	vessel	7.92E-01	9.73E-01
chr2:1743862:1743888	vessel	7.92E-01	8.09E-01
chr19:2513246:2513267	vessel	7.92E-01	9.54E-01
chr9:110501886:110502051	vessel	7.92E-01	9.73E-01
chr12:66025456:66025636	vessel	7.92E-01	9.75E-01
chr19:45323868:45323940	vessel	7.92E-01	9.27E-01
chr5:141225844:141225865	neural	7.92E-01	9.00E-01
chr17:79026865:79026951	vessel	7.92E-01	9.85E-01
chr1:16625787:16625849	vessel	7.92E-01	9.80E-01
chr2:238321665:238321816	vessel	7.92E-01	9.20E-01
chr5:148345506:148345687	vessel	7.92E-01	9.73E-01
chr1:172113662:172113921	vessel	7.92E-01	9.71E-01
chr4:7369242:7369265	vessel	7.92E-01	9.10E-01

uMHL Markers			
Region	Group	GSI	refMax
chr1:11979065:11979214	vessel	7.92E-01	9.70E-01
chr16:585886:585937	liver	7.92E-01	8.80E-01
chr4:124707773:124707972	vessel	7.92E-01	9.42E-01
chr8:109355301:109355329	vessel	7.92E-01	8.33E-01
chr12:107799287:107799331	vessel	7.92E-01	8.97E-01
chr2:218722958:218723051	vessel	7.92E-01	8.87E-01
chr9:132476798:132476839	vessel	7.92E-01	9.13E-01
chr19:13108902:13109094	vessel	7.92E-01	9.59E-01
chr6:155398687:155398778	vessel	7.92E-01	9.07E-01
chr10:35256271:35256338	vessel	7.92E-01	9.64E-01
chr4:15429515:15429573	vessel	7.92E-01	9.29E-01
chr8:67405522:67405594	neural,vessel	7.92E-01	8.89E-01
chr16:1533917:1533944	neural	7.92E-01	7.09E-01
chr10:63797212:63797266	vessel	7.92E-01	8.82E-01
chr12:125242958:125242996	vessel	7.92E-01	9.37E-01
chr21:33835784:33836018	vessel	7.92E-01	7.98E-01
chr8:134194805:134194894	vessel	7.92E-01	9.23E-01
chr1:53992643:53992661	vessel	7.92E-01	9.97E-01
chr2:72058564:72058590	vessel	7.92E-01	9.25E-01
chr17:47768524:47768612	vessel	7.91E-01	9.59E-01
chr7:70868724:70868738	vessel	7.91E-01	9.27E-01
chr7:1135787:1135880	vessel	7.91E-01	9.19E-01
chr11:69981466:69981483	vessel	7.91E-01	9.04E-01
chr10:44406319:44406371	vessel	7.91E-01	9.72E-01
chr17:62321657:62321769	vessel	7.91E-01	8.43E-01
chr9:16276328:16276511	vessel	7.91E-01	9.60E-01
chr22:26150622:26150652	vessel	7.91E-01	9.03E-01
chr13:111173782:111173824	vessel	7.91E-01	9.64E-01
chr9:132465990:132466044	vessel	7.91E-01	8.89E-01
chr9:137601605:137601641	vessel	7.91E-01	9.03E-01
chr12:2378226:2378306	vessel	7.91E-01	9.42E-01
chr6:3868913:3868926	vessel	7.91E-01	8.56E-01
chr13:110956979:110957046	vessel	7.91E-01	9.19E-01
chr3:10619418:10619432	vessel	7.91E-01	9.31E-01
chr16:73454450:73454527	vessel	7.91E-01	7.99E-01
chr15:63305915:63305934	vessel	7.91E-01	9.61E-01
chr7:4734422:4734577	vessel	7.91E-01	9.52E-01
chr10:121006525:121006629	vessel	7.91E-01	9.38E-01
chr8:144951293:144951312	pancreas	7.91E-01	9.25E-01
chr2:223917511:223917627	vessel	7.91E-01	9.00E-01
chr16:72892847:72892919	vessel	7.91E-01	8.71E-01
chr15:88525288:88525358	vessel	7.91E-01	8.77E-01
chr18:76680559:76680740	pancreas	7.91E-01	8.21E-01
chr16:79785861:79785893	vessel	7.91E-01	9.69E-01
chr13:113257964:113257997	vessel	7.91E-01	9.67E-01
chr15:62851448:62851528	vessel	7.90E-01	9.67E-01
chr17:47651712:47651734	vessel	7.90E-01	9.79E-01
chr15:94039831:94039861	vessel	7.90E-01	8.86E-01

uMHL Markers			
Region	Group	GSI	refMax
chr12:108788025:108788069	vessel	7.90E-01	9.26E-01
chr11:63823800:63823819	neural	7.90E-01	8.00E-01
chr5:106553748:106553861	vessel	7.90E-01	9.36E-01
chr12:43135893:43136097	vessel	7.90E-01	8.94E-01
chr1:20081194:20081286	vessel	7.90E-01	9.44E-01
chr3:73560709:73560905	vessel	7.90E-01	9.25E-01
chr3:123836445:123836611	vessel	7.90E-01	9.04E-01
chr1:7603622:7603720	vessel	7.90E-01	8.99E-01
chr6:158464875:158465035	vessel	7.90E-01	9.46E-01
chr6:131958144:131958342	vessel	7.90E-01	8.98E-01
chr15:67458370:67458404	vessel	7.90E-01	9.35E-01
chr8:22458256:22458378	vessel	7.90E-01	9.49E-01
chr7:94025305:94025363	vessel	7.90E-01	9.64E-01
chr6:16431273:16431391	vessel	7.90E-01	9.09E-01
chr20:56523021:56523183	vessel	7.90E-01	9.04E-01
chr5:50219227:50219301	vessel	7.90E-01	7.47E-01
chr21:47403142:47403214	vessel	7.90E-01	8.68E-01
chr9:136801147:136801304	intestine	7.90E-01	8.17E-01
chr18:9535708:9535902	vessel	7.90E-01	9.25E-01
chr18:9809125:9809169	vessel	7.90E-01	9.09E-01
chr4:101089727:101089795	vessel	7.90E-01	9.35E-01
chr20:60833682:60833702	vessel	7.90E-01	9.20E-01
chr2:66528452:66528529	vessel	7.90E-01	9.80E-01
chr20:22707146:22707191	vessel	7.90E-01	8.44E-01
chr10:44776635:44776761	vessel	7.90E-01	8.66E-01
chr3:52869083:52869162	vessel	7.90E-01	9.16E-01
chr9:129155547:129155601	vessel	7.90E-01	8.97E-01
chr1:11039869:11039880	neural	7.90E-01	8.00E-01
chr16:4461722:4462177	neural	7.90E-01	8.30E-01
chr5:173252378:173252404	vessel	7.90E-01	8.73E-01
chr9:124042382:124042419	vessel	7.90E-01	9.78E-01
chr18:56719116:56719243	vessel	7.90E-01	9.24E-01
chr11:11425708:11425747	vessel	7.89E-01	8.33E-01
chr12:52056721:52056773	liver	7.89E-01	8.09E-01
chr1:39847630:39847734	vessel	7.89E-01	9.26E-01
chr12:114327643:114328085	vessel	7.89E-01	9.52E-01
chr16:27256643:27256871	vessel	7.89E-01	9.53E-01
chr13:111830279:111830297	vessel	7.89E-01	9.46E-01
chr19:49804030:49804037	vessel	7.89E-01	8.62E-01
chr11:133371858:133371874	vessel	7.89E-01	9.54E-01
chr10:105378386:105378436	vessel	7.89E-01	9.49E-01
chr2:21500056:21500255	vessel	7.89E-01	8.37E-01
chr14:69419282:69419313	vessel	7.89E-01	9.48E-01
chr10:80913511:80913602	vessel	7.89E-01	8.73E-01
chr14:69203310:69203388	vessel	7.89E-01	8.57E-01
chr16:29706207:29706243	vessel	7.89E-01	9.41E-01
chr8:49493773:49493912	vessel	7.89E-01	8.98E-01
chr10:13869492:13869597	vessel	7.89E-01	9.58E-01

-57-

uMHL Markers			
Region	Group	GSI	refMax
chr17:75724157:75724205	vessel	7.89E-01	8.95E-01
chr7:105134734:105134810	vessel	7.89E-01	9.09E-01
chr10:78157310:78157339	vessel	7.89E-01	9.38E-01
chr1:213155773:213155827	vessel	7.89E-01	9.55E-01
chr17:12453277:12453304	vessel	7.89E-01	9.07E-01
chr11:130358485:130358572	vessel	7.89E-01	9.83E-01
chr19:3536938:3536990	kidney	7.89E-01	7.85E-01
chr5:172030681:172030738	vessel	7.89E-01	9.62E-01
chr1:24429671:24429689	vessel	7.89E-01	9.74E-01
chr3:128098054:128098098	vessel	7.89E-01	8.66E-01
chr15:67142607:67142698	vessel	7.89E-01	9.30E-01
chr17:64536087:64536138	neural	7.89E-01	8.10E-01
chr4:77598869:77598976	vessel	7.89E-01	8.99E-01
chr10:10965036:10965181	pancreas	7.89E-01	7.53E-01
chr2:223917632:223917948	vessel	7.89E-01	8.31E-01
chr1:61742965:61742988	vessel	7.89E-01	8.91E-01
chr14:21292114:21292310	neural	7.89E-01	7.29E-01
chr2:46564649:46564698	vessel	7.89E-01	8.81E-01
chr9:109561663:109561825	vessel	7.89E-01	9.15E-01
chr10:112395034:112395067	vessel	7.89E-01	7.88E-01
chr7:73406549:73406623	vessel	7.89E-01	9.58E-01
chr19:13176607:13176640	vessel	7.89E-01	9.29E-01
chr4:166300973:166301086	vessel	7.89E-01	9.67E-01
chr16:1559867:1559928	vessel	7.89E-01	9.81E-01
chr1:85069466:85069782	vessel	7.89E-01	9.30E-01
chr2:8360458:8360518	vessel	7.89E-01	7.94E-01
chr12:116008258:116008282	neural	7.89E-01	8.08E-01
chr9:132456706:132456716	vessel	7.89E-01	9.51E-01
chr10:10686890:10686959	vessel	7.89E-01	8.74E-01
chr11:46740703:46740795	liver	7.88E-01	8.46E-01
chr22:37419994:37420244	liver	7.88E-01	7.79E-01
chr21:47512230:47512291	vessel	7.88E-01	8.70E-01
chr15:36169148:36169197	vessel	7.88E-01	9.59E-01
chr14:65147082:65147109	vessel	7.88E-01	9.19E-01
chr16:49867825:49867955	vessel	7.88E-01	9.73E-01
chr9:116464320:116464372	vessel	7.88E-01	9.07E-01
chr13:110953976:110954153	vessel	7.88E-01	8.15E-01
chr15:68236450:68236481	vessel	7.88E-01	9.46E-01
chr3:122803555:122803707	liver	7.88E-01	9.15E-01
chr11:130735452:130735480	vessel	7.88E-01	9.24E-01
chr10:73608417:73608452	neural	7.88E-01	7.89E-01
chr13:111013814:111013913	vessel	7.88E-01	9.63E-01
chr12:92290560:92290631	vessel	7.88E-01	9.07E-01
chr16:14503474:14503520	vessel	7.88E-01	9.32E-01
chr16:49615634:49615896	vessel	7.88E-01	9.41E-01
chr8:130499865:130499936	vessel	7.88E-01	9.17E-01
chr10:50216911:50217019	vessel	7.88E-01	9.07E-01
chr12:20503508:20503619	vessel	7.88E-01	9.37E-01

-58-

uMHL Markers			
Region	Group	GSI	refMax
chr4:7287509:7287534	vessel	7.88E-01	9.35E-01
chr13:80511183:80511242	vessel	7.88E-01	8.71E-01
chr11:12166090:12166124	vessel	7.88E-01	8.65E-01
chr22:38864829:38864862	vessel	7.88E-01	8.37E-01
chr2:39945836:39946013	vessel	7.88E-01	9.44E-01
chr3:71289392:71289422	vessel	7.88E-01	9.31E-01
chr2:1689778:1689882	vessel	7.88E-01	8.04E-01
chr15:67471125:67471172	vessel	7.88E-01	9.53E-01
chr13:114076193:114076253	vessel	7.88E-01	9.24E-01
chr3:4468948:4469000	vessel	7.88E-01	9.07E-01
chr13:114064797:114064981	vessel	7.88E-01	8.95E-01
chr16:75520052:75520187	intestine	7.88E-01	8.45E-01
chr2:86452555:86452614	vessel	7.88E-01	8.71E-01
chr16:65105398:65105960	vessel	7.88E-01	9.54E-01
chr2:235286886:235286969	vessel	7.88E-01	9.20E-01
chr10:3182566:3182624	vessel	7.88E-01	9.04E-01
chr4:129445552:129445744	vessel	7.87E-01	9.32E-01
chr16:49750543:49750635	vessel	7.87E-01	9.44E-01
chr1:150164178:150164217	vessel	7.87E-01	8.23E-01
chr13:114797349:114797404	vessel	7.87E-01	9.06E-01
chr13:47824530:47824646	vessel	7.87E-01	9.27E-01
chr11:115882971:115883000	vessel	7.87E-01	9.58E-01
chr18:77398729:77398773	vessel	7.87E-01	9.37E-01
chr1:3282991:3283074	vessel	7.87E-01	9.49E-01
chr17:15416725:15416832	vessel	7.87E-01	9.02E-01
chr11:12098290:12098375	vessel	7.87E-01	9.37E-01
chr9:13444156:13444186	vessel	7.87E-01	8.91E-01
chr7:135430783:135430838	vessel	7.87E-01	9.31E-01
chr3:25927064:25927249	vessel	7.87E-01	8.72E-01
chr5:65124029:65124117	vessel	7.87E-01	9.78E-01
chr11:12455168:12455287	vessel	7.87E-01	9.57E-01
chr2:33295587:33295754	vessel	7.87E-01	8.15E-01
chr2:240386966:240387067	vessel	7.87E-01	9.92E-01
chr12:2411116:24111340	vessel	7.87E-01	9.61E-01
chr22:46302183:46302264	vessel	7.87E-01	9.38E-01
chr5:148851569:148851626	intestine	7.87E-01	8.54E-01
chr6:170491419:170491447	kidney	7.87E-01	8.52E-01
chr2:116891317:116891376	vessel	7.87E-01	8.23E-01
chr3:20108380:20108465	vessel	7.87E-01	9.06E-01
chr9:132156886:132156911	neural	7.87E-01	7.58E-01
chr11:76290294:76290405	vessel	7.87E-01	9.16E-01
chr22:45945164:45945215	vessel	7.87E-01	9.08E-01
chr2:42279054:42279101	vessel	7.87E-01	9.56E-01
chr3:58151255:58151409	vessel	7.87E-01	9.97E-01
chr8:30073333:30073369	vessel	7.87E-01	9.47E-01
chr9:15250121:15250232	vessel	7.87E-01	9.80E-01
chr17:17780294:17780370	vessel	7.87E-01	9.60E-01
chr13:39969110:39969133	vessel	7.87E-01	9.51E-01

uMHL Markers			
Region	Group	GSI	refMax
chr22:31502874:31502905	vessel	7.87E-01	9.74E-01
chr19:31125579:31125600	pancreas	7.87E-01	8.04E-01
chr14:84557314:84557472	vessel	7.87E-01	8.52E-01
chr21:43514198:43514493	vessel	7.87E-01	9.64E-01
chr18:22307046:22307231	vessel	7.87E-01	9.57E-01
chr8:21916771:21916803	vessel	7.87E-01	9.14E-01
chr8:22498729:22498809	vessel	7.86E-01	9.68E-01
chr20:30318764:30318781	vessel	7.86E-01	8.10E-01
chr6:136257638:136257678	vessel	7.86E-01	8.92E-01
chr8:25250313:25250341	neural,vessel	7.86E-01	9.44E-01
chr1:41272962:41273105	vessel	7.86E-01	9.17E-01
chr9:97771885:97771906	vessel	7.86E-01	9.85E-01
chr13:113695160:113695394	vessel	7.86E-01	8.81E-01
chr2:20380163:20380252	liver	7.86E-01	8.58E-01
chr16:19863013:19863115	neural	7.86E-01	9.00E-01
chr8:125677517:125677565	vessel	7.86E-01	9.55E-01
chr7:127670761:127670841	neural	7.86E-01	7.69E-01
chr6:168195164:168195233	pancreas	7.86E-01	7.64E-01
chr15:85493580:85493683	vessel	7.86E-01	9.53E-01
chr22:43031651:43031853	vessel	7.86E-01	9.73E-01
chr13:113496626:113496879	vessel	7.86E-01	8.70E-01
chr16:14489990:14490064	neural	7.86E-01	7.96E-01
chr6:57042660:57042734	vessel	7.86E-01	8.96E-01
chr3:46941526:46941607	vessel	7.86E-01	9.27E-01
chr20:35987305:35987352	vessel	7.86E-01	8.54E-01
chr16:14401759:14401838	vessel	7.86E-01	9.52E-01
chr1:14469003:14469067	vessel	7.86E-01	9.33E-01
chr20:56534156:56534196	vessel	7.86E-01	9.66E-01
chr9:132236150:132236184	vessel	7.86E-01	9.52E-01
chr1:3050597:3050633	vessel	7.86E-01	9.60E-01
chr7:37524088:37524305	vessel	7.86E-01	8.47E-01
chr13:111024166:111024198	vessel	7.86E-01	9.20E-01
chr12:6202364:6202394	vessel	7.86E-01	8.65E-01
chr17:2116086:2116211	vessel	7.86E-01	9.50E-01
chr7:4647117:4647169	neural	7.86E-01	8.68E-01
chr5:55821161:55821181	vessel	7.86E-01	9.22E-01
chr7:2666687:2666748	vessel	7.86E-01	8.39E-01
chr22:45914756:45914798	vessel	7.86E-01	8.64E-01
chr16:20624541:20624773	vessel	7.86E-01	8.82E-01
chr2:237551381:237551403	vessel	7.86E-01	9.84E-01
chr3:31363948:31364089	vessel	7.86E-01	9.52E-01
chr2:71927938:71927991	vessel	7.86E-01	9.43E-01
chr16:80972085:80972168	vessel	7.85E-01	9.86E-01
chr19:49673692:49673775	intestine,stomach	7.85E-01	8.60E-01
chr2:236843445:236843509	neural	7.85E-01	7.69E-01
chr1:15751369:15751446	intestine	7.85E-01	8.46E-01
chr5:73045778:73045803	kidney	7.85E-01	9.14E-01
chr6:128738839:128739002	vessel	7.85E-01	9.17E-01

-60-

uMHL Markers			
Region	Group	GSI	refMax
chr4:85725549:85725639	vessel	7.85E-01	9.38E-01
chr3:176622338:176622418	vessel	7.85E-01	8.82E-01
chr5:172926847:172926905	vessel	7.85E-01	9.46E-01
chr16:85240886:85240936	vessel	7.85E-01	8.75E-01
chr19:15324808:15324820	vessel	7.85E-01	9.51E-01
chr10:16636258:16636309	vessel	7.85E-01	9.18E-01
chr17:79394818:79394831	vessel	7.85E-01	9.02E-01
chr19:4373886:4374073	vessel	7.85E-01	9.41E-01
chr22:27641211:27641248	vessel	7.85E-01	8.96E-01
chr11:1418459:1418562	neural	7.85E-01	7.35E-01
chr7:151511780:151511833	neural	7.85E-01	7.76E-01
chr15:68677027:68677109	vessel	7.85E-01	9.37E-01
chr6:157470211:157470320	vessel	7.85E-01	9.24E-01
chr19:18567964:18568032	vessel	7.85E-01	9.50E-01
chr4:183152531:183152684	vessel	7.85E-01	9.72E-01
chr2:218800243:218800316	vessel	7.85E-01	9.83E-01
chr2:239367426:239367504	vessel	7.85E-01	9.60E-01
chr13:33163751:33163849	vessel	7.85E-01	9.22E-01
chr6:5068481:5068499	liver	7.85E-01	8.40E-01
chr1:42091868:42091905	vessel	7.85E-01	8.60E-01
chr7:134374786:134374800	vessel	7.85E-01	9.22E-01
chr16:15916856:15916909	vessel	7.85E-01	9.74E-01
chr8:41061579:41061651	vessel	7.85E-01	8.97E-01
chr10:30400645:30400733	vessel	7.85E-01	9.62E-01
chr8:30421070:30421152	vessel	7.85E-01	9.36E-01
chr8:139613769:139613920	neural	7.85E-01	7.22E-01
chr10:114760378:114760459	vessel	7.85E-01	9.89E-01
chr5:167361642:167361742	vessel	7.85E-01	9.28E-01
chr2:220634112:220634129	vessel	7.85E-01	9.17E-01
chr1:3157480:3157571	vessel	7.85E-01	9.26E-01
chr3:177545779:177546015	liver	7.85E-01	8.04E-01
chr8:54866663:54866734	vessel	7.85E-01	9.50E-01
chr2:65804282:65804447	vessel	7.84E-01	9.02E-01
chr20:35943041:35943225	vessel	7.84E-01	9.24E-01
chr19:1254191:1254317	neural	7.84E-01	8.08E-01
chr8:8699113:8699285	vessel	7.84E-01	9.32E-01
chr16:51610883:51611016	vessel	7.84E-01	8.73E-01
chr10:131744697:131744726	vessel	7.84E-01	8.56E-01
chr13:110147363:110147475	vessel	7.84E-01	9.55E-01
chr5:134583058:134583072	vessel	7.84E-01	8.62E-01
chr1:27891699:27891725	vessel	7.84E-01	9.55E-01
chr2:47234673:47234721	vessel	7.84E-01	9.27E-01
chr6:74965618:74965651	esophagus	7.84E-01	8.32E-01
chr4:7911612:7911672	vessel	7.84E-01	9.21E-01
chr19:5059366:5059423	neural	7.84E-01	9.24E-01
chr13:24192625:24192669	vessel	7.84E-01	9.16E-01
chr11:27536140:27536237	vessel	7.84E-01	9.69E-01
chr5:116066604:116066649	vessel	7.84E-01	9.30E-01

uMHL Markers			
Region	Group	GSI	refMax
chr11:7508112:7508197	vessel	7.84E-01	8.82E-01
chr6:166879764:166879786	vessel	7.84E-01	9.59E-01
chr6:168068881:168068949	liver	7.84E-01	8.89E-01
chr11:113424086:113424149	vessel	7.84E-01	9.27E-01
chr5:151059157:151059183	vessel	7.84E-01	9.50E-01
chr22:43821384:43821428	vessel	7.84E-01	9.18E-01
chr3:53807276:53807287	vessel	7.84E-01	9.49E-01
chr17:62167638:62167728	vessel	7.84E-01	9.69E-01
chr8:123860971:123861059	vessel	7.84E-01	9.26E-01
chr11:68894912:68894969	vessel	7.84E-01	9.47E-01
chr16:16098416:16098605	vessel	7.84E-01	8.77E-01
chr5:55193804:55193902	vessel	7.84E-01	8.07E-01
chr10:84118577:84118605	heart	7.84E-01	7.93E-01
chr1:226852160:226852233	vessel	7.84E-01	9.89E-01
chr19:36180754:36180817	vessel	7.84E-01	9.35E-01
chr20:31144652:31144702	pancreas	7.84E-01	7.25E-01
chr16:79296157:79296186	vessel	7.84E-01	8.19E-01
chr22:19863472:19863495	intestine,vessel	7.83E-01	9.22E-01
chr12:105021929:105021960	vessel	7.83E-01	9.50E-01
chr15:56209215:56209303	vessel	7.83E-01	8.85E-01
chr15:29261980:29262074	neural	7.83E-01	7.43E-01
chr7:132087715:132087737	vessel	7.83E-01	9.36E-01
chr2:151333241:151333340	vessel	7.83E-01	9.54E-01
chr19:6485754:6485800	intestine	7.83E-01	8.25E-01
chr13:31272143:31272303	vessel	7.83E-01	8.99E-01
chr4:166301164:166301272	vessel	7.83E-01	9.20E-01
chr2:121603185:121603207	vessel	7.83E-01	9.02E-01
chr6:45530885:45530911	vessel	7.83E-01	9.34E-01
chr14:94447303:94447348	vessel	7.83E-01	9.37E-01
chr1:242948252:242948360	vessel	7.83E-01	8.99E-01
chr19:49481383:49481409	vessel	7.83E-01	9.30E-01
chr10:52905224:52905373	vessel	7.83E-01	9.02E-01
chr8:25867146:25867194	vessel	7.83E-01	9.14E-01
chr15:101912644:101912745	vessel	7.83E-01	9.07E-01
chr4:40847393:40847513	vessel	7.83E-01	9.08E-01
chr17:70442942:70443022	vessel	7.83E-01	9.65E-01
chr16:86060972:86061154	intestine	7.83E-01	8.22E-01
chr17:17448903:17448931	vessel	7.83E-01	9.41E-01
chr5:134577993:134578042	vessel	7.83E-01	9.07E-01
chr10:126644018:126644097	vessel	7.83E-01	9.57E-01
chr9:132248466:132248663	vessel	7.83E-01	9.33E-01
chr3:149310847:149311097	vessel	7.82E-01	8.96E-01
chr21:46549535:46549556	vessel	7.82E-01	9.41E-01
chr8:97399458:97399485	vessel	7.82E-01	9.50E-01
chr10:134976229:134976248	neural	7.82E-01	8.67E-01
chr17:40672856:40672864	vessel	7.82E-01	9.66E-01
chr16:29706041:29706152	vessel	7.82E-01	9.84E-01
chr12:93708655:93708836	intestine	7.82E-01	7.67E-01

uMHL Markers			
Region	Group	GSI	refMax
chr17:4231044:4231126	vessel	7.82E-01	8.75E-01
chr10:104883494:104883599	vessel	7.82E-01	9.68E-01
chr3:134115902:134115909	vessel	7.82E-01	9.55E-01
chr10:81182412:81182519	vessel	7.82E-01	8.66E-01
chr12:43129874:43129987	vessel	7.82E-01	9.02E-01
chr16:50619336:50619567	vessel	7.82E-01	9.37E-01
chr11:133917437:133917528	vessel	7.82E-01	9.22E-01
chr12:23568298:23568404	vessel	7.82E-01	9.75E-01
chr9:13485157:13485208	vessel	7.82E-01	8.87E-01
chr1:170677495:170677579	vessel	7.82E-01	8.77E-01
chr17:79027073:79027118	neural	7.82E-01	8.77E-01
chr6:4351916:4351953	vessel	7.82E-01	9.58E-01
chr8:13371792:13371850	vessel	7.82E-01	9.81E-01
chr10:3366176:3366407	vessel	7.82E-01	9.36E-01
chr7:137669463:137669576	vessel	7.82E-01	9.20E-01
chr9:124535135:124535408	vessel	7.82E-01	9.16E-01
chr6:122273295:122273448	vessel	7.82E-01	8.98E-01
chr10:30932923:30932972	vessel	7.82E-01	9.29E-01
chr8:102457944:102457961	vessel	7.82E-01	8.08E-01
chr6:89871857:89871883	vessel	7.82E-01	9.05E-01
chr11:40301875:40301907	neural	7.82E-01	7.81E-01
chr2:100241347:100241372	vessel	7.82E-01	9.52E-01
chr6:43464262:43464281	vessel	7.82E-01	9.57E-01
chr11:2871840:2871864	neural	7.82E-01	8.29E-01
chr13:53617306:53617318	vessel	7.82E-01	9.45E-01
chr12:2482078:2482201	vessel	7.82E-01	7.57E-01
chr1:172114039:172114125	vessel	7.82E-01	9.42E-01
chr2:216299434:216299546	vessel	7.82E-01	9.67E-01
chr13:27254470:27254486	vessel	7.82E-01	7.89E-01
chr2:38387148:38387190	vessel	7.82E-01	8.78E-01
chr9:137426232:137426270	neural	7.82E-01	8.22E-01
chr1:19778768:19779073	vessel	7.82E-01	9.21E-01
chr2:20834577:20834625	vessel	7.82E-01	9.27E-01
chr19:12732552:12732570	neural	7.82E-01	9.70E-01
chr5:134605567:134605740	vessel	7.82E-01	8.57E-01
chr11:1929407:1929484	vessel	7.82E-01	9.66E-01
chr7:140048950:140049086	neural	7.82E-01	8.47E-01
chr16:73086674:73086828	vessel	7.82E-01	8.89E-01
chr1:62353776:62353847	vessel	7.82E-01	9.58E-01
chr3:16707220:16707260	intestine	7.82E-01	7.96E-01
chr15:99050082:99050108	vessel	7.82E-01	9.28E-01
chr17:60706523:60706581	vessel	7.81E-01	8.84E-01
chr15:67142417:67142575	vessel	7.81E-01	9.28E-01
chr20:56576504:56576598	vessel	7.81E-01	9.24E-01
chr19:5036492:5036562	neural	7.81E-01	8.58E-01
chr20:25030492:25030618	vessel	7.81E-01	9.49E-01
chr6:1766006:1766029	vessel	7.81E-01	9.67E-01
chr18:72839832:72839844	intestine	7.81E-01	7.43E-01

-63-

uMHL Markers			
Region	Group	GSI	refMax
chr15:99995054:99995113	vessel	7.81E-01	9.19E-01
chr12:130612543:130612596	vessel	7.81E-01	9.07E-01
chr6:37648587:37648654	vessel	7.81E-01	9.27E-01
chr12:124984979:124985033	vessel	7.81E-01	8.60E-01

Table 1b:
Top 10% tissue specific MHL markers identified by GSI

MHL Markers			
Region	Group	GSI	refMax
chr5:122422637:122422689	vessel	8.44E-01	9.80E-01
chr16:88293071:88293119	vessel	8.35E-01	9.71E-01
chr14:91790551:91790559	intestine	8.29E-01	9.76E-01
chr7:560607:560650	vessel	8.29E-01	9.12E-01
chr18:34823918:34823977	intestine	8.28E-01	8.70E-01
chr5:122422972:122423004	vessel	8.27E-01	9.45E-01
chr18:52613464:52613527	neural	8.25E-01	7.82E-01
chr4:174440618:174440652	vessel	8.20E-01	8.88E-01
chr6:6000351:6000385	vessel	8.20E-01	9.56E-01
chr6:85476611:85476636	vessel	8.20E-01	9.53E-01
chr2:66810606:66810640	vessel	8.20E-01	9.31E-01
chr16:73098764:73098784	pancreas	8.18E-01	9.15E-01
chr1:3310806:3310814	pancreas	8.17E-01	8.36E-01
chr8:99961711:99961827	vessel	8.16E-01	9.06E-01
chr7:154720589:154720672	neural	8.15E-01	8.17E-01
chr5:122423184:122423204	vessel	8.10E-01	8.88E-01
chr8:97165649:97165702	vessel	8.03E-01	9.14E-01
chr2:45231186:45231227	vessel	8.03E-01	9.24E-01
chr17:59532637:59532659	vessel	8.02E-01	9.46E-01
chr5:122422007:122422064	vessel	8.00E-01	9.31E-01
chr11:2114366:2114498	kidney	7.97E-01	8.05E-01
chr6:6000008:6000046	vessel	7.96E-01	8.05E-01
chr5:72731939:72731975	vessel	7.96E-01	8.41E-01
chr8:71287244:71287415	neural	7.95E-01	9.26E-01
chr6:39966947:39967010	neural	7.93E-01	8.73E-01
chr15:96909551:96909595	vessel	7.93E-01	9.30E-01
chr7:560864:560875	vessel	7.92E-01	8.38E-01
chr15:53087377:53087420	pancreas	7.91E-01	7.83E-01
chr7:560685:560695	vessel	7.91E-01	8.47E-01
chr6:159655065:159655082	fat	7.90E-01	8.71E-01
chr5:122422161:122422184	vessel	7.90E-01	8.74E-01
chr1:27848198:27848263	liver	7.88E-01	8.51E-01
chr7:561160:561184	pancreas,vessel	7.88E-01	1.00E+00
chr8:99959774:99959900	vessel	7.87E-01	7.89E-01
chr16:73098734:73098753	pancreas	7.87E-01	8.02E-01
chr3:73620624:73620814	neural	7.86E-01	9.13E-01
chr4:24801971:24802055	vessel	7.85E-01	8.70E-01
chr19:13209981:13210007	vessel	7.84E-01	7.81E-01

MHL Markers			
Region	Group	GSI	refMax
chr2:177004040:177004111	vessel	7.84E-01	7.65E-01
chr17:55520634:55520788	liver	7.83E-01	8.43E-01
chr7:560782:560796	vessel	7.82E-01	7.80E-01
chr5:158532062:158532093	vessel	7.80E-01	8.87E-01
chr7:5336543:5336571	liver	7.80E-01	8.23E-01
chr2:45240094:45240129	vessel	7.79E-01	9.15E-01
chr3:23653540:23653753	liver	7.79E-01	7.28E-01
chr2:66810478:66810502	vessel	7.77E-01	8.71E-01
chr11:64509762:64509801	vessel	7.77E-01	8.83E-01
chr5:72676020:72676058	vessel	7.76E-01	9.25E-01
chr7:25892505:25892545	vessel	7.76E-01	8.16E-01
chr15:96897012:96897055	vessel	7.75E-01	7.68E-01
chr10:126407964:126408032	vessel	7.75E-01	9.42E-01
chr1:47909185:47909226	vessel	7.74E-01	9.58E-01
chr7:157479397:157479515	neural	7.74E-01	8.12E-01
chr5:72597625:72597716	vessel	7.74E-01	7.40E-01
chr9:97431689:97431863	neural	7.73E-01	7.47E-01
chr8:99961971:99962204	vessel	7.72E-01	7.37E-01
chr5:81652981:81653357	neural	7.72E-01	7.79E-01
chr12:54088920:54088960	vessel	7.72E-01	8.92E-01
chr8:99951082:99951152	vessel	7.71E-01	8.81E-01
chr4:15412389:15412548	neural	7.70E-01	8.16E-01
chr15:41217790:41217801	vessel	7.70E-01	8.59E-01
chr17:48206061:48206096	vessel	7.70E-01	9.27E-01
chr2:45231539:45231586	vessel	7.69E-01	7.61E-01
chr9:98785107:98785147	vessel	7.69E-01	8.36E-01
chr4:1194418:1194505	liver	7.69E-01	8.79E-01
chr15:74426581:74426618	fat	7.68E-01	7.56E-01
chr17:17628486:17628493	intestine	7.68E-01	8.62E-01
chr12:30976279:30976313	neural	7.68E-01	8.54E-01
chr16:86968517:86968564	vessel	7.67E-01	9.11E-01
chr16:86959129:86959213	neural	7.66E-01	8.41E-01
chr11:110581348:110581376	vessel	7.66E-01	7.83E-01
chr17:17628638:17628722	liver	7.64E-01	8.20E-01
chr16:51669186:51669276	neural	7.64E-01	8.14E-01
chr10:119972986:119973016	neural	7.63E-01	7.86E-01
chr13:28000468:28000558	vessel	7.63E-01	7.92E-01
chr8:99954669:99954686	vessel	7.63E-01	7.38E-01
chr18:60173469:60173666	neural	7.62E-01	7.69E-01
chr10:105452540:105452576	liver	7.62E-01	7.91E-01
chr16:89005381:89005398	vessel	7.60E-01	9.62E-01
chr6:168533689:168533704	vessel	7.60E-01	7.51E-01
chr8:80740825:80740880	vessel	7.59E-01	8.20E-01
chr19:2240088:2240161	liver	7.59E-01	7.52E-01
chr12:115251376:115251453	kidney	7.58E-01	8.58E-01
chr18:12287438:12287452	vessel	7.57E-01	7.15E-01
chr20:39320779:39320848	vessel	7.57E-01	8.55E-01
chr12:54807242:54807336	neural	7.57E-01	7.35E-01

MHL Markers			
Region	Group	GSI	refMax
chr2:177037468:177037632	vessel	7.57E-01	8.72E-01
chr10:5489561:5489690	neural	7.57E-01	8.28E-01
chr12:54400441:54400545	vessel	7.56E-01	8.44E-01
chr21:36901623:36901692	neural	7.56E-01	7.36E-01
chr2:241395207:241395347	vessel	7.55E-01	8.71E-01
chr1:145440313:145440506	vessel	7.55E-01	7.93E-01
chr14:105126561:105126572	pancreas	7.55E-01	8.19E-01
chr8:102506473:102506597	vessel	7.55E-01	8.05E-01
chr13:95354223:95354262	neural	7.54E-01	7.52E-01
chr22:43659543:43659642	neural	7.53E-01	7.43E-01
chr2:10231429:10231487	pancreas	7.53E-01	7.77E-01
chr5:72676905:72676924	vessel	7.53E-01	8.17E-01
chr19:19571777:19571806	liver	7.53E-01	7.92E-01
chr17:32705890:32705981	neural	7.52E-01	7.92E-01
chr15:53098366:53098407	pancreas	7.52E-01	7.64E-01
chr10:123355655:123355837	vessel	7.52E-01	8.06E-01
chr20:3053093:3053103	liver	7.51E-01	8.58E-01
chr1:170635953:170635964	vessel	7.51E-01	8.93E-01
chr5:20041646:20041789	neural	7.51E-01	8.13E-01
chr20:39319423:39319447	esophagus	7.51E-01	7.25E-01
chr7:27194521:27194570	vessel	7.51E-01	7.12E-01
chr10:119300620:119300704	vessel	7.51E-01	8.79E-01
chr8:99962451:99962688	vessel	7.51E-01	7.67E-01
chr15:53087876:53087895	pancreas	7.50E-01	7.66E-01
chr9:139740765:139740775	pancreas	7.50E-01	8.86E-01
chr14:105944278:105944310	vessel	7.50E-01	9.04E-01
chr9:102587733:102587805	vessel	7.49E-01	8.72E-01
chr16:81030772:81030844	neural	7.49E-01	8.30E-01
chr2:89166335:89166578	neural	7.49E-01	8.89E-01
chr8:99951364:99951421	vessel	7.48E-01	7.27E-01
chr8:76316319:76316353	vessel	7.48E-01	8.59E-01
chr12:54408685:54408713	kidney	7.47E-01	7.53E-01
chr9:136567964:136567975	liver	7.47E-01	7.47E-01
chr5:72526768:72526786	vessel	7.46E-01	7.66E-01
chr2:66809772:66809852	vessel	7.46E-01	8.14E-01
chr2:96814448:96814515	liver	7.46E-01	8.86E-01
chr5:38368165:38368317	neural	7.46E-01	7.62E-01
chr1:110610678:110610715	vessel	7.45E-01	7.56E-01
chr12:54440712:54440753	vessel	7.45E-01	8.38E-01
chr5:122421804:122421820	vessel	7.44E-01	9.30E-01
chr7:19149999:19150182	kidney,vessel	7.44E-01	8.83E-01
chr19:13124365:13124391	pancreas	7.44E-01	7.18E-01
chr12:118314033:118314075	pancreas	7.43E-01	6.85E-01

Table 2:
Complete list of high methylated haplotype shared between matched primary tumor tissues
and plasma for colon cancer (CRC) and lung cancer (LC) patients.

Chr	Start	End	Coordinate	Symbol	Annotation
1	3527738	3527759	chr1:3527738-3527759	MEGF6	Promoter
1	13839806	13839815	chr1:13839806-13839815	n/a	n/a
1	14219586	14219639	chr1:14219586-14219639	n/a	n/a
1	17019874	17019889	chr1:17019874-17019889	ESPNP	UTR3
1	21836158	21836216	chr1:21836158-21836216	ALPL	UTR5
1	22889799	22889812	chr1:22889799-22889812	EPHA8	Promoter
1	23279723	23279744	chr1:23279723-23279744	n/a	n/a
1	23280192	23280211	chr1:23280192-23280211	n/a	n/a
1	27676453	27676616	chr1:27676453-27676616	SYTL1	Intron
1	29563577	29563710	chr1:29563577-29563710	PTPRU	Promoter
1	39981227	39981247	chr1:39981227-39981247	BMP8A	Intron
1	44883591	44883606	chr1:44883591-44883606	RNF220	Intron
1	46859961	46859974	chr1:46859961-46859974	FAAH	UTR5
1	48190891	48190924	chr1:48190891-48190924	n/a	n/a
1	50884393	50884411	chr1:50884393-50884411	DMRTA2	Exon
1	50884419	50884430	chr1:50884419-50884430	DMRTA2	Exon
1	50884472	50884632	chr1:50884472-50884632	DMRTA2	Exon
1	59280289	59280357	chr1:59280289-59280357	n/a	n/a
1	59280369	59280455	chr1:59280369-59280455	n/a	n/a
1	61517642	61517933	chr1:61517642-61517933	n/a	n/a
1	74663749	74663776	chr1:74663749-74663776	LRRIQ3	UTR5
1	92946665	92946767	chr1:92946665-92946767	GFI1	Intron
1	108508530	108508549	chr1:108508530-108508549	VAV3	Promoter
1	119543104	119543127	chr1:119543104-119543127	TBX15	Enhancer
1	150293718	150293852	chr1:150293718-150293852	PRPF3	Promoter
1	156828865	156828914	chr1:156828865-156828914	NTRK1	Promoter
1	158151057	158151116	chr1:158151057-158151116	CD1D	Intron
1	171810397	171810513	chr1:171810397-171810513	DNM3	Promoter
1	203096889	203096934	chr1:203096889-203096934	ADORA1	UTR5
1	203598610	203598622	chr1:203598610-203598622	ATP2B4	UTR5
1	214153443	214153464	chr1:214153443-214153464	PROX1	Enhancer
1	215256127	215256195	chr1:215256127-215256195	KCNK2	Promoter
1	234040873	234041006	chr1:234040873-234041006	SLC35F3	Promoter
1	234350457	234350469	chr1:234350457-234350469	SLC35F3	Intron
1	236558377	236558653	chr1:236558377-236558653	EDARADD	Promoter
1	242687600	242687610	chr1:242687600-242687610	PLD5	UTR5
1	244894219	244894231	chr1:244894219-244894231	n/a	n/a
1	246952304	246952348	chr1:246952304-246952348	LOC149134	Promoter
2	3751318	3751336	chr2:3751318-3751336	ALLC	Downstream
2	8833584	8833597	chr2:8833584-8833597	n/a	n/a
2	11809996	11810041	chr2:11809996-11810041	NTSR2	Promoter
2	26521758	26521880	chr2:26521758-26521880	n/a	n/a
2	31360797	31360816	chr2:31360797-31360816	GALNT14	Promoter
2	39102723	39102771	chr2:39102723-39102771	DHX57	UTR5
2	47596444	47596455	chr2:47596444-47596455	EPCAM	UTR5
2	47596483	47596505	chr2:47596483-47596505	EPCAM	UTR5

Chr	Start	End	Coordinate	Symbol	Annotation
2	70056907	70056960	chr2:70056907-70056960	GMCL1	UTR5
2	70994743	70994753	chr2:70994743-70994753	ADD2	UTR5
2	73429837	73429858	chr2:73429837-73429858	NOTO	Promoter
2	88470034	88470056	chr2:88470034-88470056	THNSL2	UTR5
2	105459218	105459234	chr2:105459218-105459234	POU3F3	Enhancer
2	110372264	110372284	chr2:110372264-110372284	ANKRD57	Promoter
2	113956590	113956654	chr2:113956590-113956654	LOC440839	UTR3
2	119604037	119604049	chr2:119604037-119604049	EN1	Exon
2	119916314	119916319	chr2:119916314-119916319	C1QL2	UTR5
2	127783309	127783371	chr2:127783309-127783371	n/a	n/a
2	175594922	175594966	chr2:175594922-175594966	n/a	n/a
2	176945353	176945373	chr2:176945353-176945373	EVX2	Exon
2	176971839	176971857	chr2:176971839-176971857	HOXD11	Promoter
2	176972805	176972814	chr2:176972805-176972814	HOXD11	Promoter
2	207307674	207307712	chr2:207307674-207307712	ADAM23	Promoter
2	220313255	220313271	chr2:220313255-220313271	SPEG	Exon
2	233351443	233351467	chr2:233351443-233351467	ECEL1	UTR5
2	239755874	239755895	chr2:239755874-239755895	TWIST2	Promoter
3	9988661	9989203	chr3:9988661-9989203	PRRT3	Exon
3	10749398	10749432	chr3:10749398-10749432	n/a	n/a
3	10749434	10749467	chr3:10749434-10749467	n/a	n/a
3	12046446	12046504	chr3:12046446-12046504	SYN2	Promoter
3	13324006	13324119	chr3:13324006-13324119	n/a	n/a
3	16554839	16555196	chr3:16554839-16555196	RFTN1	UTR5
3	42304944	42304981	chr3:42304944-42304981	CCK	Exon
3	54156903	54156978	chr3:54156903-54156978	CACNA2D3	Promoter
3	54156985	54157005	chr3:54156985-54157005	CACNA2D3	Promoter
3	69591310	69591445	chr3:69591310-69591445	n/a	n/a
3	119528983	119529218	chr3:119528983-119529218	NR1I2	Intron
3	122641209	122641233	chr3:122641209-122641233	SEMA5B	Exon
3	130646225	130646274	chr3:130646225-130646274	ATP2C1	Intron
3	151178854	151178938	chr3:151178854-151178938	IGSF10	Promoter
3	187676563	187676643	chr3:187676563-187676643	n/a	n/a
4	124515	124758	chr4:124515-124758	ZNF718	Intron
4	467665	467684	chr4:467665-467684	ZNF721	UTR5
4	658004	658032	chr4:658004-658032	PDE6B	Intron
4	3873182	3873271	chr4:3873182-3873271	n/a	n/a
4	5713317	5713393	chr4:5713317-5713393	EVC2	Promoter
4	11370394	11370433	chr4:11370394-11370433	MIR572	Promoter
4	11370452	11370520	chr4:11370452-11370520	MIR572	UTR5
4	30719538	30719763	chr4:30719538-30719763	PCDH7	Promoter
4	42153515	42153591	chr4:42153515-42153591	BEND4	Intron
4	44449839	44449861	chr4:44449839-44449861	KCTD8	Promoter
4	55098179	55098209	chr4:55098179-55098209	PDGFRA	UTR5
4	55991626	55991684	chr4:55991626-55991684	KDR	UTR5
4	62067517	62067536	chr4:62067517-62067536	n/a	n/a
4	77819160	77819271	chr4:77819160-77819271	ANKRD56	Promoter
4	103940573	103940885	chr4:103940573-103940885	NHEDC1	UTR5
4	103940890	103941101	chr4:103940890-103941101	NHEDC1	Promoter

Chr	Start	End	Coordinate	Symbol	Annotation
4	126236257	126236931	chr4:126236257-126236931	FAT4	Promoter
4	144833142	144833212	chr4:144833142-144833212	GYPE	Enhancer
4	151000141	151000246	chr4:151000141-151000246	DCLK2	UTR5
5	191803	191818	chr5:191803-191818	LRRC14B	Promoter
5	5139924	5139935	chr5:5139924-5139935	ADAMTS16	Promoter
5	8457467	8457735	chr5:8457467-8457735	n/a	n/a
5	43603785	43604084	chr5:43603785-43604084	NNT	UTR5
5	43604116	43604165	chr5:43604116-43604165	NNT	UTR5
5	55117727	55117748	chr5:55117727-55117748	n/a	n/a
5	76476761	76476791	chr5:76476761-76476791	n/a	n/a
5	76507004	76507082	chr5:76507004-76507082	PDE8B	Promoter
5	77147764	77147911	chr5:77147764-77147911	n/a	n/a
5	113391874	113391904	chr5:113391874-113391904	n/a	n/a
5	115298740	115298778	chr5:115298740-115298778	LVRN	Promoter
5	150536722	150536746	chr5:150536722-150536746	ANXA6	UTR5
5	155108288	155108355	chr5:155108288-155108355	n/a	n/a
5	174151522	174151577	chr5:174151522-174151577	MSX2	UTR5
5	179780701	179780801	chr5:179780701-179780801	GFPT2	Promoter
5	180486476	180486537	chr5:180486476-180486537	BTNL9	Exon
6	1390421	1390427	chr6:1390421-1390427	FOXF2	Promoter
6	1555484	1555571	chr6:1555484-1555571	n/a	n/a
6	3849234	3849541	chr6:3849234-3849541	FAM50B	Promoter
6	11242082	11242134	chr6:11242082-11242134	NEDD9	Intron
6	18122712	18122718	chr6:18122712-18122718	NHLRC1	Promoter
6	29760303	29760314	chr6:29760303-29760314	HCG4	UTR3
6	35888708	35888855	chr6:35888708-35888855	SRPK1	UTR5
6	42072328	42072372	chr6:42072328-42072372	C6orf132	Exon
6	137818835	137818915	chr6:137818835-137818915	OLIG3	Promoter
6	152129664	152129700	chr6:152129664-152129700	ESR1	Promoter
6	166582820	166582835	chr6:166582820-166582835	T	Promoter
6	170581003	170581095	chr6:170581003-170581095	LOC154449	Enhancer
7	27138133	27138172	chr7:27138133-27138172	HOXA1	Promoter
7	27146154	27146541	chr7:27146154-27146541	HOXA3	UTR3
7	27147986	27148068	chr7:27147986-27148068	HOXA3	Exon
7	27162225	27162404	chr7:27162225-27162404	HOXA3	UTR5
7	27182613	27183574	chr7:27182613-27183574	HOXA5	UTR5
7	27196517	27196529	chr7:27196517-27196529	HOXA7	Promoter
7	28997143	28997166	chr7:28997143-28997166	TRIL	Promoter
7	32110174	32110180	chr7:32110174-32110180	PDE1C	UTR5
7	37392844	37393124	chr7:37392844-37393124	ELMO1	UTR5
7	43798072	43798080	chr7:43798072-43798080	BLVRA	Promoter
7	44349388	44349523	chr7:44349388-44349523	CAMK2B	Intron
7	45002209	45002518	chr7:45002209-45002518	MYO1G	UTR3
7	45002526	45002634	chr7:45002526-45002634	MYO1G	Intron
7	45614570	45614580	chr7:45614570-45614580	ADCY1	Promoter
7	45615005	45615102	chr7:45615005-45615102	ADCY1	Promoter
7	50850278	50850659	chr7:50850278-50850659	GRB10	UTR5
7	56355680	56355715	chr7:56355680-56355715	n/a	n/a
7	64022988	64023250	chr7:64022988-64023250	ZNF680	Promoter

Chr	Start	End	Coordinate	Symbol	Annotation
7	65971098	65971186	chr7:65971098-65971186	n/a	n/a
7	82792010	82792128	chr7:82792010-82792128	PCLO	UTR5
7	100203335	100203366	chr7:100203335-100203366	PCOLCE	Exon
7	103629981	103630092	chr7:103629981-103630092	RELN	Promoter
7	113727611	113727622	chr7:113727611-113727622	n/a	n/a
7	127672078	127672112	chr7:127672078-127672112	LRR4	Promoter
7	127743728	127743766	chr7:127743728-127743766	n/a	n/a
7	130130739	130131267	chr7:130130739-130131267	MESTIT1	UTR5
7	130131358	130131518	chr7:130131358-130131518	MEST	UTR5
7	149746006	149746019	chr7:149746006-149746019	n/a	n/a
7	150812726	150812750	chr7:150812726-150812750	AGAP3	Intron
7	151216757	151216773	chr7:151216757-151216773	RHEB	UTR5
7	153583591	153583622	chr7:153583591-153583622	DPP6	Promoter
7	155247552	155247562	chr7:155247552-155247562	EN2	Promoter
7	155595896	155595952	chr7:155595896-155595952	SHH	Exon
7	156400470	156400500	chr7:156400470-156400500	n/a	n/a
7	157486226	157486275	chr7:157486226-157486275	PTPRN2	Intron
7	158938119	158938146	chr7:158938119-158938146	VIPR2	Promoter
8	2585693	2585757	chr8:2585693-2585757	n/a	n/a
8	4851492	4851508	chr8:4851492-4851508	CSMD1	Promoter
8	9009136	9009388	chr8:9009136-9009388	PPP1R3B	Promoter
8	23260741	23260788	chr8:23260741-23260788	LOXL2	UTR5
8	31497559	31497576	chr8:31497559-31497576	NRG1	Promoter
8	37552122	37552160	chr8:37552122-37552160	ZNF703	Promoter
8	37699481	37699558	chr8:37699481-37699558	GPR124	Exon
8	38034584	38034641	chr8:38034584-38034641	LSM1	Promoter
8	41166680	41166708	chr8:41166680-41166708	SFRP1	UTR5
8	49426959	49427414	chr8:49426959-49427414	n/a	n/a
8	54163561	54163585	chr8:54163561-54163585	OPRK1	Promoter
8	54163604	54163694	chr8:54163604-54163694	OPRK1	UTR5
8	55380019	55380033	chr8:55380019-55380033	n/a	n/a
8	58055200	58055257	chr8:58055200-58055257	n/a	n/a
8	65493709	65493763	chr8:65493709-65493763	BHLHE22	Promoter
8	67874080	67874104	chr8:67874080-67874104	n/a	n/a
8	72756057	72756082	chr8:72756057-72756082	MSC	Promoter
8	98290011	98290080	chr8:98290011-98290080	TSPYL5	UTR5
8	103750881	103750903	chr8:103750881-103750903	n/a	n/a
8	127568853	127569069	chr8:127568853-127569069	FAM84B	Exon
8	141108442	141109280	chr8:141108442-141109280	TRAPPC9	Intron
8	144511400	144511448	chr8:144511400-144511448	MAFA	Downstream
8	145104394	145104454	chr8:145104394-145104454	OPLAH	Downstream
9	113865	113881	chr9:113865-113881	n/a	n/a
9	113884	113897	chr9:113884-113897	n/a	n/a
9	19788900	19788911	chr9:19788900-19788911	SLC24A2	Promoter
9	25677605	25677627	chr9:25677605-25677627	TUSC1	UTR3
9	35689643	35689690	chr9:35689643-35689690	TPM2	Promoter
9	38424066	38424081	chr9:38424066-38424081	IGFBPL1	Promoter
9	89560709	89560739	chr9:89560709-89560739	GAS1	Exon
9	95572080	95572086	chr9:95572080-95572086	ANKRD19	UTR5

Chr	Start	End	Coordinate	Symbol	Annotation
9	99983989	99984041	chr9:99983989-99984041	KIAA1529	Enhancer
9	101471709	101471724	chr9:101471709-101471724	GABBR2	Promoter
9	101706293	101706314	chr9:101706293-101706314	COL15A1	UTR5
9	120507462	120507562	chr9:120507462-120507562	n/a	n/a
9	123656794	123657026	chr9:123656794-123657026	PHF19	Enhancer
9	123657048	123657162	chr9:123657048-123657162	PHF19	Enhancer
9	124888893	124889126	chr9:124888893-124889126	n/a	n/a
9	132382398	132382811	chr9:132382398-132382811	C9orf50	Promoter
9	133536491	133536515	chr9:133536491-133536515	PRDM12	Promoter
9	133536616	133536699	chr9:133536616-133536699	PRDM12	Promoter
9	135462555	135462589	chr9:135462555-135462589	BARHL1	Intron
9	137979579	137979590	chr9:137979579-137979590	OLFM1	Intron
9	139964715	139964731	chr9:139964715-139964731	C9orf140	Promoter
10	7708553	7708596	chr10:7708553-7708596	ITIH5	Promoter
10	20104705	20104724	chr10:20104705-20104724	PLXDC2	Promoter
10	25465355	25465408	chr10:25465355-25465408	LOC100128811	Promoter
10	29698362	29698585	chr10:29698362-29698585	LOC387647	UTR5
10	71626579	71626666	chr10:71626579-71626666	COL13A1	Intron
10	77158757	77158887	chr10:77158757-77158887	C10orf41	Promoter
10	80898884	80899123	chr10:80898884-80899123	ZMIZ1	UTR5
10	81163306	81163339	chr10:81163306-81163339	ZCCHC24	Intron
10	83634362	83634433	chr10:83634362-83634433	NRG3	Promoter
10	88123205	88123264	chr10:88123205-88123264	GRID1	Intron
10	101282028	101282143	chr10:101282028-101282143	NKX2-3	Enhancer
10	102495446	102495452	chr10:102495446-102495452	PAX2	Enhancer
10	102507681	102507717	chr10:102507681-102507717	PAX2	Intron
10	105344583	105344617	chr10:105344583-105344617	NEURL	Exon
10	105452852	105452884	chr10:105452852-105452884	SH3PXD2A	Intron
10	105453074	105453169	chr10:105453074-105453169	SH3PXD2A	Intron
10	111216768	111216809	chr10:111216768-111216809	n/a	n/a
10	119301950	119302046	chr10:119301950-119302046	EMX2	UTR5
10	119311968	119311994	chr10:119311968-119311994	EMX2OS	Enhancer
10	119313192	119313239	chr10:119313192-119313239	EMX2OS	Enhancer
10	125732491	125732516	chr10:125732491-125732516	n/a	n/a
10	131767467	131767523	chr10:131767467-131767523	EBF3	Enhancer
10	133999328	133999363	chr10:133999328-133999363	DPYSL4	Promoter
10	134222564	134222659	chr10:134222564-134222659	PWWP2B	UTR3
10	135090324	135090391	chr10:135090324-135090391	ADAM8	UTR5
11	397076	397141	chr11:397076-397141	PKP3	Exon
11	518994	519003	chr11:518994-519003	LRRCL2	Enhancer
11	726322	726388	chr11:726322-726388	EPS8L2	Exon
11	1874410	1874461	chr11:1874410-1874461	LSP1	Promoter
11	2021030	2021337	chr11:2021030-2021337	MIR675	Promoter
11	17297912	17298333	chr11:17297912-17298333	NUCB2	UTR5
11	17740909	17740930	chr11:17740909-17740930	MYOD1	Promoter
11	24518517	24518550	chr11:24518517-24518550	LUZP2	Promoter
11	35641254	35641291	chr11:35641254-35641291	FJX1	Exon
11	47236054	47236189	chr11:47236054-47236189	DDB2	Promoter
11	47611788	47611855	chr11:47611788-47611855	C1QTNF4	Exon

Chr	Start	End	Coordinate	Symbol	Annotation
11	60692224	60692379	chr11:60692224-60692379	TMEM132A	Promoter
11	61880088	61880140	chr11:61880088-61880140	INCENP	Enhancer
11	69589140	69589199	chr11:69589140-69589199	FGF4	Promoter
11	69924935	69924948	chr11:69924935-69924948	ANO1	Promoter
11	71951194	71951198	chr11:71951194-71951198	PHOX2A	Exon
11	72295582	72295589	chr11:72295582-72295589	PDE2A	UTR3
11	72295726	72295758	chr11:72295726-72295758	PDE2A	UTR3
11	72533078	72533338	chr11:72533078-72533338	ATG16L2	Intron
11	79148648	79148661	chr11:79148648-79148661	ODZ4	UTR5
11	82443937	82443948	chr11:82443937-82443948	FAM181B	Promoter
11	100998291	100998355	chr11:100998291-100998355	PGR	Exon
11	109963265	109963348	chr11:109963265-109963348	ZC3H12C	Promoter
11	109964113	109964167	chr11:109964113-109964167	ZC3H12C	UTR5
11	134146600	134146621	chr11:134146600-134146621	GLB1L3	UTR5
12	2163267	2163279	chr12:2163267-2163279	CACNA1C	Promoter
12	2800445	2800521	chr12:2800445-2800521	CACNA1C	UTR3
12	3309861	3309888	chr12:3309861-3309888	TSPAN9	UTR5
12	7781181	7781237	chr12:7781181-7781237	n/a	n/a
12	9217328	9217429	chr12:9217328-9217429	LOC144571	Promoter
12	29936643	29936653	chr12:29936643-29936653	TMT1	UTR5
12	54764364	54764584	chr12:54764364-54764584	ZNF385A	Intron
12	57869148	57869420	chr12:57869148-57869420	ARHGAP9	Intron
12	58025887	58025901	chr12:58025887-58025901	B4GALNT1	Exon
12	58119853	58120184	chr12:58119853-58120184	LOC100130776	UTR5
12	58131738	58132045	chr12:58131738-58132045	AGAP2	UTR5
12	122016340	122016373	chr12:122016340-122016373	KDM2B	Intron
12	132195645	132195994	chr12:132195645-132195994	SFRS8	UTR5
13	19918950	19918983	chr13:19918950-19918983	LOC100101938	UTR3
13	20139192	20139253	chr13:20139192-20139253	n/a	n/a
13	20692669	20692685	chr13:20692669-20692685	n/a	n/a
13	25115907	25115943	chr13:25115907-25115943	n/a	n/a
13	26625786	26625915	chr13:26625786-26625915	SHISA2	Promoter
13	48893192	48893246	chr13:48893192-48893246	RB1	Intron
13	52703312	52703361	chr13:52703312-52703361	NEK5	Promoter
13	100608204	100608226	chr13:100608204-100608226	n/a	n/a
13	109147798	109147937	chr13:109147798-109147937	n/a	n/a
13	109148352	109148477	chr13:109148352-109148477	n/a	n/a
13	110959180	110959184	chr13:110959180-110959184	COL4A2	Promoter
13	112723104	112723110	chr13:112723104-112723110	SOX1	Exon
13	113764991	113765210	chr13:113764991-113765210	F7	Intron
13	114462331	114462426	chr13:114462331-114462426	FAM70B	Promoter
14	28733691	28733732	chr14:28733691-28733732	n/a	n/a
14	37051685	37051713	chr14:37051685-37051713	NKX2-8	UTR5
14	38091516	38091570	chr14:38091516-38091570	n/a	n/a
14	42077408	42077482	chr14:42077408-42077482	LRFN5	UTR5
14	48143559	48143579	chr14:48143559-48143579	MDGA2	Promoter
14	77737169	77737210	chr14:77737169-77737210	NGB	Promoter
14	91720063	91720097	chr14:91720063-91720097	GPR68	UTR5
14	97499682	97499715	chr14:97499682-97499715	n/a	n/a

Chr	Start	End	Coordinate	Symbol	Annotation
14	101925421	101925446	chr14:101925421-101925446	n/a	n/a
14	105640938	105641002	chr14:105640938-105641002	NUDT14	Intron
14	105767212	105767276	chr14:105767212-105767276	BRF1	UTR5
15	23892769	23892883	chr15:23892769-23892883	MAGEL2	Promoter
15	26107989	26108171	chr15:26107989-26108171	ATP10A	Promoter
15	26108184	26108248	chr15:26108184-26108248	ATP10A	UTR5
15	27213029	27213059	chr15:27213029-27213059	GABRG3	Promoter
15	29034062	29034076	chr15:29034062-29034076	n/a	n/a
15	29034154	29034183	chr15:29034154-29034183	n/a	n/a
15	33010242	33010280	chr15:33010242-33010280	GREM1	UTR5
15	35047291	35047296	chr15:35047291-35047296	GJD2	Promoter
15	40583268	40583755	chr15:40583268-40583755	PLCB2	Intron
15	47477259	47477276	chr15:47477259-47477276	n/a	n/a
15	56035589	56035677	chr15:56035589-56035677	PRTG	Promoter
15	65360274	65360278	chr15:65360274-65360278	RASL12	UTR5
15	70354770	70354890	chr15:70354770-70354890	TLE3	Intron
15	74044709	74044798	chr15:74044709-74044798	C15orf59	Promoter
15	79576092	79576118	chr15:79576092-79576118	ANKRD34C	UTR5
15	79724688	79724701	chr15:79724688-79724701	KIAA1024	Promoter
15	89148392	89148460	chr15:89148392-89148460	MIR7-2	Enhancer
15	97491220	97491238	chr15:97491220-97491238	n/a	n/a
15	100386502	100386576	chr15:100386502-100386576	n/a	n/a
16	128297	128580	chr16:128297-128580	MPG	UTR5
16	610097	610110	chr16:610097-610110	C16orf11	Promoter
16	1203760	1203786	chr16:1203760-1203786	CACNA1H	Promoter
16	2009454	2009590	chr16:2009454-2009590	NDUFB10	UTR5
16	2317602	2317630	chr16:2317602-2317630	RNPS1	UTR5
16	28074226	28074311	chr16:28074226-28074311	GSG1L	Promoter
16	30572753	30572772	chr16:30572753-30572772	ZNF764	Promoter
16	47177577	47177628	chr16:47177577-47177628	NETO2	UTR5
16	51185385	51185428	chr16:51185385-51185428	SALL1	Promoter
16	67313433	67313443	chr16:67313433-67313443	PLEKHG4	UTR5
16	67687484	67687553	chr16:67687484-67687553	RLTPR	Intron
16	68269381	68269396	chr16:68269381-68269396	ESRP2	Promoter
16	68679063	68679166	chr16:68679063-68679166	CDH3	UTR5
16	69419782	69420122	chr16:69419782-69420122	TERF2	UTR5
16	75284386	75284475	chr16:75284386-75284475	BCAR1	Intron
16	79633394	79633613	chr16:79633394-79633613	MAF	Exon
16	86541944	86542182	chr16:86541944-86542182	FOXF1	Promoter
16	88454809	88454843	chr16:88454809-88454843	n/a	n/a
16	88600767	88600893	chr16:88600767-88600893	ZFPM1	Exon
16	90113989	90114077	chr16:90113989-90114077	LOC100130015	UTR5
17	3438996	3439015	chr17:3438996-3439015	TRPV3	Intron
17	14205002	14205182	chr17:14205002-14205182	MGC12916	Promoter
17	17399374	17399399	chr17:17399374-17399399	RASD1	Promoter
17	18061388	18061524	chr17:18061388-18061524	MYO15A	Intron
17	37761831	37761946	chr17:37761831-37761946	NEUROD2	Exon
17	43047732	43047918	chr17:43047732-43047918	C1QL1	Promoter
17	43507012	43507098	chr17:43507012-43507098	SH3D20	Exon

Chr	Start	End	Coordinate	Symbol	Annotation
17	46621993	46622021	chr17:46621993-46622021	HOXB2	Promoter
17	46641662	46641747	chr17:46641662-46641747	HOXB3	UTR5
17	46641966	46642036	chr17:46641966-46642036	HOXB3	UTR5
17	47574666	47574903	chr17:47574666-47574903	NGFR	Intron
17	56833042	56833161	chr17:56833042-56833161	PPM1E	Promoter
17	62774654	62774696	chr17:62774654-62774696	LOC146880	UTR3
17	74072736	74072749	chr17:74072736-74072749	ZACN	Promoter
17	74381072	74381119	chr17:74381072-74381119	SPHK1	UTR5
17	79455513	79455575	chr17:79455513-79455575	n/a	n/a
17	81057628	81057660	chr17:81057628-81057660	n/a	n/a
18	4455202	4455210	chr18:4455202-4455210	n/a	n/a
18	5891068	5891178	chr18:5891068-5891178	TMEM200C	Promoter
18	8706308	8706439	chr18:8706308-8706439	KIAA0802	Enhancer
18	10726392	10726415	chr18:10726392-10726415	FAM38B	Intron
18	74332372	74332448	chr18:74332372-74332448	n/a	n/a
18	76740034	76740059	chr18:76740034-76740059	SALL3	Promoter
18	77376899	77377025	chr18:77376899-77377025	n/a	n/a
19	1071232	1071353	chr19:1071232-1071353	HMHA1	Intron
19	1450118	1450129	chr19:1450118-1450129	APC2	Promoter
19	2576229	2576291	chr19:2576229-2576291	GNG7	UTR5
19	3404842	3405130	chr19:3404842-3405130	NFIC	Intron
19	3933421	3933495	chr19:3933421-3933495	ITGB1BP3	UTR5
19	6274084	6274159	chr19:6274084-6274159	MLLT1	Intron
19	6744852	6744978	chr19:6744852-6744978	TRIP10	Exon
19	10444873	10444951	chr19:10444873-10444951	RAVER1	Promoter
19	10445195	10445607	chr19:10445195-10445607	RAVER1	Promoter
19	10531599	10531608	chr19:10531599-10531608	PDE4A	Promoter
19	12306248	12306298	chr19:12306248-12306298	n/a	n/a
19	14584456	14584478	chr19:14584456-14584478	PTGER1	Exon
19	16022797	16022848	chr19:16022797-16022848	CYP4F2	Enhancer
19	16181371	16181533	chr19:16181371-16181533	TPM4	Intron
19	16437563	16437597	chr19:16437563-16437597	KLF2	Intron
19	17392927	17393042	chr19:17392927-17393042	ANKLE1	Promoter
19	18303568	18304395	chr19:18303568-18304395	MPV17L2	UTR5
19	23653826	23653836	chr19:23653826-23653836	n/a	n/a
19	30017014	30017126	chr19:30017014-30017126	VSTM2B	Promoter
19	34113352	34113366	chr19:34113352-34113366	CHST8	UTR5
19	38042272	38042323	chr19:38042272-38042323	ZNF540	UTR5
19	39261605	39261612	chr19:39261605-39261612	LGALS7	UTR3
19	39798183	39798269	chr19:39798183-39798269	LRFN1	UTR3
19	40732582	40732618	chr19:40732582-40732618	CNTD2	UTR5
19	46974787	46974808	chr19:46974787-46974808	PNMAL1	UTR5
19	48946520	48946656	chr19:48946520-48946656	GRWD1	Promoter
19	48983752	48983868	chr19:48983752-48983868	CYTH2	UTR3
19	52104699	52104735	chr19:52104699-52104735	n/a	n/a
19	54040763	54041012	chr19:54040763-54041012	ZNF331	UTR5
19	54412970	54412994	chr19:54412970-54412994	CACNG7	Promoter
19	54483370	54483483	chr19:54483370-54483483	MIR935	Promoter
19	57351416	57351756	chr19:57351416-57351756	ZIM2	UTR5

Chr	Start	End	Coordinate	Symbol	Annotation
19	58858757	58858820	chr19:58858757-58858820	A1BG	Exon
19	58867719	58867836	chr19:58867719-58867836	A1BG	Promoter
20	1784445	1784482	chr20:1784445-1784482	n/a	n/a
20	2674280	2674284	chr20:2674280-2674284	EBF4	Promoter
20	5296982	5297005	chr20:5296982-5297005	PROKR2	Promoter
20	5484944	5485022	chr20:5484944-5485022	LOC149837	UTR3
20	5485089	5485267	chr20:5485089-5485267	LOC149837	UTR5
20	10647776	10647921	chr20:10647776-10647921	JAG1	Intron
20	13201045	13201055	chr20:13201045-13201055	ISM1	Promoter
20	22549193	22549242	chr20:22549193-22549242	C20orf56	UTR3
20	22562939	22562952	chr20:22562939-22562952	C20orf56	Promoter
20	44639238	44639267	chr20:44639238-44639267	MMP9	Exon
20	55841104	55841148	chr20:55841104-55841148	BMP7	Promoter
20	55841151	55841254	chr20:55841151-55841254	BMP7	UTR5
20	57089850	57089875	chr20:57089850-57089875	APCDD1L	UTR5
20	57464110	57464262	chr20:57464110-57464262	GNAS	UTR5
20	59827735	59827770	chr20:59827735-59827770	CDH4	Promoter
20	59827789	59827831	chr20:59827789-59827831	CDH4	Promoter
20	60877547	60877638	chr20:60877547-60877638	ADRM1	Promoter
21	28216583	28216627	chr21:28216583-28216627	ADAMTS1	Exon
21	28217665	28217711	chr21:28217665-28217711	ADAMTS1	UTR5
21	28338515	28338535	chr21:28338515-28338535	ADAMTS5	Promoter
21	32930369	32930388	chr21:32930369-32930388	TIAM1	UTR5
21	32931085	32931279	chr21:32931085-32931279	TIAM1	UTR5
21	32931284	32931321	chr21:32931284-32931321	TIAM1	UTR5
21	38630603	38630727	chr21:38630603-38630727	DSCR3	Intron
21	40033433	40033449	chr21:40033433-40033449	ERG	UTR5
22	17083312	17083542	chr22:17083312-17083542	psiTPTE22	UTR5
22	19137099	19137268	chr22:19137099-19137268	GSC2	Promoter
22	19702409	19702465	chr22:19702409-19702465	SEPT5	Promoter
22	21738255	21738732	chr22:21738255-21738732	RIMBP3B	Promoter
22	38349631	38350197	chr22:38349631-38350197	POLR2F	UTR5
22	39572333	39572455	chr22:39572333-39572455	n/a	n/a
22	46423321	46423397	chr22:46423321-46423397	n/a	n/a
22	46476250	46476464	chr22:46476250-46476464	LOC400931	Enhancer
22	46484636	46484829	chr22:46484636-46484829	LOC400931	UTR5
22	50706449	50706483	chr22:50706449-50706483	MAPK11	Intron
22	51158653	51158710	chr22:51158653-51158710	SHANK3	Exon
22	51158990	51159182	chr22:51158990-51159182	SHANK3	Exon

-75-

Table 3a:
Deconvolution of plasma samples by 10 normal tissues, LCT, and CCT

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
CCP	0.027	0.015	0.019	0.036	0.030	0.035	0.031	0.030	0.145	0.046	0.044	0.543
NP	0.015	0.002	0.002	0.001	0.037	0.010	0.013	0.010	0.056	0.044	0.003	0.808
LCP	0.045	0.013	0.057	0.046	0.047	0.041	0.048	0.035	0.095	0.044	0.042	0.488

Average values from only samples with WB > 0.3

Table 3b:

Deconvolution for colon cancer plasma

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.CRC.P.001	0.117	0.027	0.062	0.053	0.049	0.062	0.168	0.062	0.168	0.050	0.075	0.108
UCSD.CRC.P.002	0.020	0.021	0.043	0.128	0.192	0.097	0.121	0.026	0.077	0.087	0.072	0.117
UCSD.CRC.P.003	0.103	0.010	0.058	0.237	0.087	0.035	0.103	0.016	0.103	0.064	0.139	0.045
UCSD.CRC.P.004	0.046	0.080	0.145	0.082	0.119	0.076	0.113	0.000	0.077	0.080	0.121	0.061
UCSD.CRC.P.005	0.085	0.020	0.061	0.143	0.112	0.111	0.130	0.000	0.042	0.086	0.113	0.096
UCSD.CRC.P.006	0.074	0.014	0.038	0.047	0.017	0.055	0.046	0.035	0.142	0.242	0.061	0.228
UCSD.CRC.P.007	0.008	0.033	0.033	0.091	0.275	0.042	0.017	0.064	0.242	0.054	0.064	0.078
UCSD.CRC.P.008	0.045	0.000	0.055	0.095	0.257	0.098	0.040	0.031	0.212	0.065	0.064	0.039
UCSD.CRC.P.009	0.037	0.211	0.007	0.000	0.000	0.031	0.061	0.245	0.134	0.015	0.203	0.057
UCSD.CRC.P.010	0.025	0.000	0.063	0.032	0.029	0.019	0.015	0.000	0.028	0.001	0.038	0.750
UCSD.CRC.P.011	0.187	0.105	0.058	0.061	0.037	0.063	0.190	0.023	0.072	0.072	0.085	0.048
UCSD.CRC.P.012	0.023	0.029	0.032	0.066	0.029	0.083	0.026	0.045	0.077	0.005	0.059	0.525
UCSD.CRC.P.013	0.054	0.109	0.047	0.073	0.168	0.099	0.060	0.042	0.062	0.089	0.133	0.064
UCSD.CRC.P.014	0.032	0.000	0.035	0.102	0.095	0.181	0.069	0.038	0.237	0.045	0.089	0.077
UCSD.CRC.P.015	0.032	0.015	0.096	0.019	0.122	0.062	0.036	0.054	0.418	0.050	0.057	0.039
UCSD.CRC.P.016	0.012	0.085	0.017	0.017	0.024	0.054	0.072	0.200	0.109	0.171	0.114	0.124
UCSD.CRC.P.017	0.103	0.032	0.014	0.046	0.043	0.105	0.068	0.077	0.068	0.093	0.155	0.197
UCSD.CRC.P.018	0.080	0.105	0.052	0.178	0.050	0.066	0.083	0.100	0.071	0.076	0.082	0.058
UCSD.CRC.P.019	0.100	0.421	0.004	0.039	0.030	0.017	0.126	0.012	0.069	0.042	0.092	0.050
UCSD.CRC.P.020	0.023	0.072	0.003	0.072	0.033	0.097	0.033	0.155	0.083	0.007	0.116	0.305
UCSD.CRC.P.021	0.133	0.108	0.055	0.142	0.017	0.045	0.045	0.019	0.137	0.053	0.119	0.127
UCSD.CRC.P.022	0.115	0.000	0.022	0.033	0.017	0.031	0.011	0.024	0.162	0.086	0.041	0.457
UCSD.CRC.P.023	0.000	0.000	0.013	0.036	0.024	0.031	0.018	0.000	0.013	0.000	0.016	0.849
UCSD.CRC.P.024	0.055	0.050	0.010	0.081	0.091	0.116	0.088	0.018	0.168	0.040	0.209	0.074

39595779.1

-76-

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.CRC.P.025	0.006	0.002	0.000	0.028	0.000	0.037	0.016	0.035	0.307	0.248	0.005	0.318
UCSD.CRC.P.026	0.000	0.000	0.006	0.000	0.008	0.000	0.009	0.000	0.026	0.000	0.024	0.926
UCSD.CRC.P.027	0.087	0.064	0.037	0.000	0.043	0.038	0.069	0.427	0.080	0.046	0.083	0.026
UCSD.CRC.P.028	0.044	0.000	0.022	0.013	0.045	0.010	0.089	0.011	0.324	0.012	0.041	0.390
UCSD.CRC.P.029	0.058	0.013	0.043	0.045	0.052	0.072	0.069	0.076	0.272	0.022	0.117	0.160
UCSD.CRC.P.030	0.006	0.028	0.006	0.041	0.085	0.009	0.063	0.000	0.288	0.051	0.059	0.363

Table 3c:
Deconvolution for lung cancer plasma

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.LC.001	0.118	0.178	0.037	0.018	0.046	0.152	0.107	0.099	0.055	0.056	0.065	0.070
UCSD.LC.002	0.234	0.084	0.049	0.037	0.087	0.053	0.049	0.077	0.151	0.081	0.068	0.029
UCSD.LC.003	0.076	0.060	0.051	0.056	0.232	0.066	0.067	0.037	0.223	0.028	0.066	0.038
UCSD.LC.004	0.023	0.054	0.030	0.077	0.079	0.028	0.131	0.000	0.212	0.102	0.074	0.190
UCSD.LC.005	0.070	0.037	0.066	0.139	0.057	0.103	0.157	0.004	0.135	0.076	0.102	0.054
UCSD.LC.006	0.115	0.017	0.073	0.149	0.066	0.064	0.065	0.022	0.112	0.104	0.088	0.126
UCSD.LC.007	0.091	0.032	0.053	0.055	0.098	0.060	0.063	0.027	0.207	0.060	0.212	0.041
UCSD.LC.008	0.052	0.022	0.059	0.117	0.098	0.083	0.139	0.015	0.080	0.078	0.089	0.167
UCSD.LC.009	0.075	0.013	0.065	0.100	0.144	0.153	0.114	0.042	0.066	0.045	0.126	0.057
UCSD.LC.010	0.051	0.058	0.220	0.066	0.098	0.016	0.062	0.089	0.063	0.045	0.091	0.141
UCSD.LC.011	0.066	0.032	0.049	0.060	0.066	0.090	0.053	0.114	0.115	0.100	0.071	0.185
UCSD.LC.012	0.045	0.037	0.028	0.010	0.022	0.022	0.087	0.020	0.180	0.027	0.064	0.458
UCSD.LC.013	0.156	0.021	0.089	0.000	0.000	0.021	0.000	0.128	0.057	0.008	0.009	0.511
UCSD.LC.014	0.040	0.235	0.037	0.129	0.037	0.055	0.195	0.000	0.118	0.045	0.079	0.029
UCSD.LC.015	0.135	0.057	0.085	0.148	0.042	0.114	0.074	0.037	0.074	0.062	0.104	0.066
UCSD.LC.016	0.239	0.055	0.039	0.109	0.070	0.055	0.125	0.049	0.085	0.094	0.024	0.055
UCSD.LC.017	0.075	0.017	0.069	0.076	0.098	0.084	0.106	0.011	0.203	0.120	0.099	0.042
UCSD.LC.018	0.051	0.028	0.145	0.187	0.088	0.091	0.149	0.013	0.039	0.053	0.130	0.026
UCSD.LC.019	0.052	0.028	0.089	0.147	0.020	0.090	0.036	0.035	0.168	0.058	0.086	0.192
UCSD.LC.020	0.010	0.010	0.041	0.052	0.138	0.005	0.061	0.015	0.036	0.058	0.016	0.558
UCSD.LC.021	0.058	0.004	0.043	0.160	0.067	0.102	0.062	0.046	0.223	0.064	0.126	0.047
UCSD.LC.022	0.064	0.024	0.041	0.039	0.096	0.042	0.066	0.000	0.086	0.049	0.072	0.420

-77-

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.LC.023	0.162	0.045	0.074	0.097	0.103	0.089	0.108	0.017	0.087	0.052	0.100	0.065
UCSD.LC.025	0.031	0.000	0.039	0.049	0.066	0.085	0.078	0.071	0.208	0.077	0.075	0.221
UCSD.LC.026	0.014	0.000	0.113	0.035	0.012	0.022	0.057	0.060	0.089	0.133	0.035	0.431
UCSD.LC.027	0.168	0.000	0.036	0.107	0.104	0.074	0.044	0.023	0.120	0.175	0.119	0.031
UCSD.LC.028	0.015	0.000	0.041	0.154	0.041	0.043	0.065	0.000	0.016	0.032	0.061	0.532
UCSD.LC.029	0.011	0.000	0.045	0.030	0.020	0.128	0.002	0.019	0.203	0.000	0.036	0.506
UCSD.LC.030	0.067	0.034	0.122	0.077	0.050	0.069	0.100	0.023	0.241	0.057	0.111	0.048

Table 3d:

Deconvolution for normal plasma

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.NP.001	0.000	0.000	0.277	0.000	0.000	0.000	0.000	0.386	0.000	0.000	0.051	0.287
UCSD.NP.002	0.000	0.000	0.191	0.213	0.105	0.036	0.000	0.000	0.355	0.000	0.055	0.046
UCSD.NP.003	0.111	0.000	0.065	0.048	0.075	0.070	0.089	0.000	0.354	0.000	0.106	0.082
UCSD.NP.004	0.000	0.000	0.047	0.107	0.071	0.059	0.067	0.000	0.491	0.000	0.095	0.062
UCSD.NP.005	0.000	0.000	0.170	0.000	0.110	0.096	0.000	0.349	0.000	0.086	0.000	0.188
UCSD.NP.006	0.000	0.000	0.000	0.000	0.023	0.000	0.000	0.000	0.000	0.005	0.000	0.973
UCSD.NP.007	0.000	0.000	0.044	0.000	0.257	0.000	0.000	0.000	0.018	0.377	0.118	0.186
UCSD.NP.008	0.000	0.000	0.016	0.000	0.230	0.000	0.000	0.000	0.359	0.322	0.073	0.000
UCSD.NP.009	0.000	0.000	0.000	0.000	0.003	0.138	0.650	0.179	0.000	0.000	0.031	0.000
UCSD.NP.010	0.000	0.011	0.000	0.000	0.000	0.122	0.000	0.039	0.000	0.000	0.000	0.828
UCSD.NP.011	0.000	0.000	0.000	0.000	0.000	0.130	0.670	0.168	0.000	0.000	0.032	0.000
UCSD.NP.012	0.000	0.000	0.081	0.027	0.301	0.000	0.000	0.000	0.061	0.389	0.141	0.000
UCSD.NP.013	0.000	0.052	0.000	0.021	0.899	0.010	0.000	0.000	0.000	0.000	0.017	0.000
UCSD.NP.014	0.000	0.000	0.000	0.000	0.021	0.048	0.751	0.124	0.004	0.000	0.052	0.000
UCSD.NP.015	0.000	0.074	0.000	0.000	0.012	0.366	0.000	0.215	0.048	0.000	0.103	0.181
UCSD.NP.016	0.000	0.000	0.000	0.000	0.022	0.049	0.759	0.112	0.006	0.000	0.053	0.000
UCSD.NP.017	0.000	0.000	0.030	0.076	0.029	0.008	0.003	0.000	0.248	0.000	0.596	0.010
UCSD.NP.018	0.000	0.000	0.028	0.037	0.240	0.000	0.000	0.000	0.012	0.573	0.109	0.000
UCSD.NP.019	0.000	0.000	0.017	0.042	0.218	0.000	0.002	0.000	0.000	0.634	0.087	0.000
UCSD.NP.020	0.000	0.000	0.033	0.056	0.081	0.036	0.036	0.000	0.036	0.000	0.692	0.030
UCSD.NP.021	0.000	0.000	0.080	0.025	0.298	0.000	0.000	0.000	0.056	0.401	0.139	0.000
UCSD.NP.022	0.000	0.069	0.000	0.000	0.029	0.386	0.000	0.236	0.067	0.000	0.114	0.099

-78-

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.NP.023	0.000	0.083	0.000	0.000	0.532	0.223	0.000	0.062	0.000	0.000	0.000	0.100
UCSD.NP.024	0.000	0.000	0.021	0.082	0.093	0.054	0.035	0.000	0.462	0.000	0.075	0.177
UCSD.NP.025	0.000	0.000	0.000	0.003	0.001	0.000	0.000	0.000	0.000	0.000	0.002	0.994
UCSD.NP.026	0.000	0.000	0.000	0.000	0.000	0.000	0.380	0.120	0.000	0.000	0.000	0.501
UCSD.NP.027	0.000	0.000	0.000	0.038	0.133	0.366	0.000	0.172	0.074	0.030	0.101	0.087
UCSD.NP.028	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.451	0.000	0.000	0.549
UCSD.NP.029	0.000	0.000	0.000	0.000	0.013	0.000	0.000	0.000	0.041	0.012	0.000	0.934
UCSD.NP.030	0.000	0.000	0.000	0.076	0.035	0.034	0.021	0.018	0.096	0.000	0.700	0.020
UCSD.NP.031	0.000	0.000	0.078	0.019	0.296	0.000	0.000	0.000	0.057	0.403	0.146	0.000
UCSD.NP.032	0.000	0.000	0.000	0.154	0.012	0.011	0.000	0.000	0.157	0.000	0.655	0.011
UCSD.NP.033	0.295	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.027	0.000	0.679
UCSD.NP.034	0.000	0.008	0.000	0.000	0.007	0.000	0.000	0.000	0.063	0.000	0.000	0.922
UCSD.NP.035	0.000	0.000	0.000	0.000	0.000	0.026	0.000	0.000	0.000	0.007	0.000	0.967
UCSD.NP.036	0.000	0.000	0.038	0.000	0.000	0.000	0.000	0.056	0.000	0.000	0.000	0.905
UCSD.NP.037	0.000	0.000	0.000	0.000	0.000	0.000	0.006	0.000	0.000	0.000	0.011	0.983
UCSD.NP.038	0.000	0.000	0.000	0.000	0.200	0.000	0.000	0.000	0.223	0.309	0.041	0.228
UCSD.NP.039	0.000	0.000	0.000	0.007	0.002	0.000	0.001	0.000	0.000	0.000	0.006	0.984
UCSD.NP.040	0.000	0.000	0.028	0.092	0.227	0.158	0.061	0.000	0.077	0.000	0.301	0.056
UCSD.NP.041	0.000	0.043	0.238	0.000	0.000	0.000	0.000	0.429	0.249	0.000	0.000	0.041
UCSD.NP.042	0.000	0.000	0.149	0.033	0.281	0.000	0.000	0.000	0.052	0.362	0.122	0.000
UCSD.NP.043	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.354	0.000	0.000	0.646
UCSD.NP.044	0.000	0.000	0.000	0.000	0.082	0.000	0.000	0.000	0.000	0.135	0.000	0.783
UCSD.NP.045	0.000	0.000	0.000	0.000	0.254	0.000	0.000	0.000	0.000	0.000	0.000	0.746
UCSD.NP.046	0.000	0.000	0.069	0.006	0.288	0.000	0.000	0.000	0.047	0.399	0.131	0.059
UCSD.NP.047	0.000	0.000	0.000	0.000	0.109	0.000	0.000	0.000	0.000	0.180	0.000	0.711
UCSD.NP.048	0.000	0.000	0.000	0.000	0.000	0.000	0.042	0.000	0.000	0.000	0.000	0.958
UCSD.NP.049	0.000	0.000	0.000	0.000	0.020	0.000	0.000	0.000	0.000	0.060	0.000	0.920
UCSD.NP.050	0.000	0.000	0.000	0.000	0.099	0.346	0.000	0.232	0.072	0.000	0.081	0.170
UCSD.NP.051	0.000	0.000	0.000	0.000	0.009	0.000	0.000	0.000	0.030	0.030	0.000	0.930
UCSD.NP.052	0.000	0.000	0.091	0.087	0.239	0.005	0.146	0.000	0.010	0.220	0.201	0.000
UCSD.NP.053	0.000	0.025	0.000	0.000	0.000	0.000	0.000	0.000	0.654	0.000	0.009	0.312
UCSD.NP.054	0.049	0.017	0.000	0.000	0.000	0.193	0.000	0.146	0.000	0.068	0.000	0.528
UCSD.NP.055	0.187	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.813

-79-

	Brain	CCT	Colon	Esophagus	Heart	Intestine	Kidney	LCT	Liver	Lung	Stomach	WB
UCSD.NP.056	0.000	0.000	0.000	0.007	0.002	0.000	0.002	0.000	0.000	0.000	0.005	0.984
UCSD.NP.057	0.842	0.000	0.000	0.000	0.000	0.000	0.061	0.000	0.000	0.087	0.010	0.000
UCSD.NP.058	0.000	0.000	0.000	0.000	0.093	0.000	0.000	0.000	0.000	0.178	0.000	0.729
UCSD.NP.059	0.000	0.000	0.008	0.000	0.194	0.000	0.000	0.000	0.000	0.299	0.046	0.454
UCSD.NP.060	0.034	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.966
UCSD.NP.061	0.024	0.000	0.494	0.000	0.108	0.048	0.015	0.000	0.106	0.000	0.172	0.032
UCSD.NP.062	0.000	0.000	0.000	0.000	0.142	0.000	0.000	0.000	0.000	0.221	0.014	0.622
UCSD.NP.063	0.000	0.000	0.000	0.000	0.000	0.026	0.005	0.009	0.000	0.000	0.000	0.959
UCSD.NP.064	0.000	0.000	0.370	0.000	0.044	0.000	0.000	0.528	0.007	0.000	0.021	0.030
UCSD.NP.065	0.000	0.000	0.000	0.000	0.182	0.000	0.000	0.000	0.000	0.000	0.000	0.818
UCSD.NP.066	0.000	0.112	0.000	0.000	0.000	0.031	0.000	0.001	0.795	0.000	0.041	0.021
UCSD.NP.067	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.026	0.974
UCSD.NP.068	0.000	0.000	0.000	0.000	0.000	0.000	0.054	0.000	0.103	0.000	0.000	0.843
UCSD.NP.069	0.000	0.000	0.000	0.000	0.066	0.000	0.000	0.000	0.031	0.166	0.000	0.737
UCSD.NP.070	0.000	0.000	0.010	0.000	0.000	0.000	0.000	0.000	0.021	0.000	0.000	0.969
UCSD.NP.071	0.000	0.000	0.082	0.027	0.300	0.000	0.000	0.000	0.064	0.387	0.140	0.000
UCSD.NP.072	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.311	0.000	0.000	0.689
UCSD.NP.073	0.000	0.000	0.000	0.000	0.121	0.000	0.000	0.000	0.000	0.200	0.003	0.675
UCSD.NP.074	0.000	0.000	0.000	0.012	0.018	0.000	0.000	0.000	0.000	0.025	0.000	0.945
UCSD.NP.075	0.000	0.000	0.000	0.000	0.021	0.000	0.000	0.000	0.000	0.031	0.000	0.948

Table 4a:
Differentially methylated MHB regions between colon cancer tissues (CCT) and normal plasma.

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr1	110693161	110693212	2.07E-06	4.11E-01	9.35E-04	9.98E-01
chr3	187676563	187676644	3.69E-16	3.66E-01	6.60E-04	9.98E-01
chr2	225266234	225266263	1.96E-04	5.49E-01	7.05E-02	8.86E-01
chr7	391332	391394	4.29E-07	5.12E-01	2.80E-02	9.48E-01
chr1	78957190	78957314	5.73E-05	3.36E-01	2.73E-03	9.92E-01
chr9	86152595	86152671	7.83E-11	6.14E-01	2.95E-04	1.00E+00
chr10	73496013	73496122	2.04E-16	3.55E-01	2.30E-02	9.39E-01
chr7	104581511	104581826	2.49E-07	5.81E-01	6.48E-02	9.00E-01
chr8	17509194	17509223	9.03E-12	5.46E-01	3.56E-02	9.39E-01
chr20	23029535	23029655	1.68E-19	4.42E-01	3.67E-03	9.92E-01
chr16	29164276	29164303	1.91E-04	6.68E-01	1.71E-01	7.96E-01
chr7	64030120	64030307	1.93E-08	3.13E-01	1.07E-03	9.97E-01
chr21	47971700	47971722	7.52E-05	6.37E-01	9.40E-02	8.71E-01
chr1	19251999	19252025	8.87E-06	7.39E-01	2.08E-01	7.80E-01
chr19	18113497	18113578	6.77E-06	5.71E-01	4.65E-02	9.25E-01
chr4	2814121	2814173	1.66E-09	3.71E-01	7.82E-03	9.79E-01
chr5	169931353	169931389	3.79E-04	3.20E-01	2.99E-03	9.91E-01
chr8	97157926	97157941	2.07E-10	3.50E-01	2.73E-04	9.99E-01
chr2	232186706	232186761	7.37E-05	7.07E-01	1.17E-01	8.58E-01
chr17	40822613	40822653	7.50E-06	5.25E-01	1.05E-01	8.34E-01
chr1	154375706	154375773	5.38E-05	4.17E-01	1.67E-02	9.61E-01
chr2	219738569	219738577	1.03E-04	3.86E-01	6.54E-02	8.55E-01
chr11	19735331	19735420	7.10E-04	3.33E-01	3.74E-03	9.89E-01
chr2	109834010	109834148	1.14E-04	4.17E-01	4.07E-02	9.11E-01
chr14	94254076	94254092	8.05E-17	5.27E-01	2.20E-02	9.60E-01
chr17	74141537	74141566	1.36E-17	4.39E-01	1.43E-02	9.68E-01
chr21	44819099	44819123	4.56E-09	4.21E-01	4.98E-02	8.94E-01
chr16	30621003	30621081	8.06E-04	3.33E-01	6.27E-05	1.00E+00
chr12	108169217	108169314	1.06E-06	3.43E-01	9.73E-03	9.72E-01
chr3	183903549	183903707	4.15E-04	5.73E-01	4.67E-02	9.25E-01
chr1	145039361	145039649	1.81E-04	3.39E-01	1.70E-03	9.95E-01

-81-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr13	43148866	43148914	5.73E-04	3.33E-01	3.62E-04	9.99E-01
chr7	49813424	49813437	4.97E-10	5.38E-01	5.44E-02	9.08E-01
chr8	38757858	38757899	5.02E-17	3.53E-01	1.47E-02	9.60E-01
chr8	142276522	142276533	8.84E-14	4.36E-01	1.94E-02	9.57E-01
chr10	114136302	114136365	1.74E-06	4.50E-01	8.45E-03	9.82E-01
chr4	122301928	122301963	8.37E-08	5.00E-01	3.37E-03	9.93E-01
chr1	12123521	12123556	9.44E-04	3.38E-01	1.85E-04	9.99E-01
chr19	30016411	30016441	7.98E-08	6.67E-01	2.05E-03	9.97E-01
chr4	17767383	17767525	1.96E-04	7.11E-01	4.31E-02	9.43E-01
chr14	74707404	74707425	2.24E-05	4.26E-01	6.81E-04	9.98E-01
chr1	225655011	225655039	1.76E-08	5.56E-01	2.22E-02	9.62E-01
chr4	111553898	111553914	3.01E-07	6.67E-01	1.81E-02	9.74E-01
chr19	3179842	3179904	4.49E-04	3.85E-01	4.09E-02	9.04E-01
chr10	123923869	123923887	1.59E-12	4.61E-01	2.60E-02	9.47E-01
chr21	46334450	46334460	4.07E-06	5.85E-01	1.06E-01	8.46E-01
chr2	70057016	70057043	3.09E-04	4.66E-01	6.67E-02	8.75E-01
chr17	40822041	40822061	1.97E-05	6.32E-01	1.43E-01	8.15E-01
chr2	25425395	25425417	5.46E-04	6.06E-01	1.72E-01	7.79E-01
chr16	68028052	68028131	2.17E-10	4.37E-01	3.45E-02	9.27E-01
chr2	177016618	177016669	1.63E-04	3.67E-01	1.44E-02	9.62E-01
chr16	88769941	88770131	3.03E-10	4.42E-01	6.08E-03	9.86E-01
chr11	123301152	123301245	1.81E-05	3.83E-01	6.85E-05	1.00E+00
chr3	192126690	192126712	9.40E-04	3.99E-01	4.97E-02	8.89E-01
chr9	133536748	133536768	5.96E-07	3.92E-01	1.82E-03	9.95E-01
chr1	9764850	9764871	2.36E-05	5.61E-01	5.82E-02	9.06E-01
chr3	48520626	48520665	6.96E-06	3.83E-01	2.73E-02	9.33E-01
chr17	13504776	13504797	1.79E-06	4.26E-01	9.32E-03	9.79E-01
chr7	32110082	32110100	2.27E-06	4.38E-01	1.28E-02	9.72E-01
chr13	109793395	109793413	8.19E-04	4.28E-01	7.93E-02	8.44E-01
chr17	36734910	36734967	2.53E-08	3.70E-01	4.52E-03	9.88E-01
chr20	37356313	37356330	9.86E-07	4.05E-01	0.00E+00	1.00E+00
chr17	59482675	59482690	1.55E-04	3.60E-01	2.63E-02	9.32E-01
chr19	40871656	40871678	2.77E-05	4.00E-01	5.21E-02	8.85E-01

39595779.1

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr10	38383097	38383309	5.65E-04	3.33E-01	0.00E+00	1.00E+00
chr6	31527706	31527920	8.63E-06	5.79E-01	5.17E-02	9.18E-01
chr7	4248707	4248866	4.90E-07	4.23E-01	4.92E-02	8.96E-01
chr17	72462820	72462912	4.98E-07	3.71E-01	1.52E-02	9.61E-01
chr5	2740503	2740527	2.56E-09	5.68E-01	2.80E-03	9.95E-01
chr8	54163561	54163586	2.32E-07	3.25E-01	4.56E-03	9.86E-01
chr1	203643083	203643220	1.35E-05	4.44E-01	1.84E-02	9.60E-01
chr15	59397205	59397420	5.93E-06	5.03E-01	8.37E-05	1.00E+00
chr11	69920045	69920065	4.11E-09	5.94E-01	4.04E-02	9.36E-01
chr19	55766296	55766320	2.26E-04	3.54E-01	3.13E-02	9.19E-01
chr8	55379351	55379448	4.82E-05	3.05E-01	3.35E-03	9.89E-01
chr4	1195837	1195893	3.91E-06	4.05E-01	5.45E-03	9.87E-01
chr8	103136001	103136031	6.36E-04	6.67E-01	7.57E-02	8.98E-01
chr19	23253315	23253335	2.38E-04	3.51E-01	3.86E-02	9.01E-01
chr11	64107477	64107619	7.34E-07	6.59E-01	5.80E-02	9.19E-01
chr3	27769863	27770001	2.45E-04	3.44E-01	5.70E-03	9.84E-01
chr11	107328048	107328380	3.75E-08	6.67E-01	3.48E-03	9.95E-01
chr7	158938119	158938147	9.53E-07	5.18E-01	8.18E-03	9.84E-01
chr8	97158006	97158018	2.13E-13	3.56E-01	9.96E-04	9.97E-01
chr12	115103838	115103858	9.69E-04	4.00E-01	1.08E-02	9.74E-01
chr1	226924888	226924924	8.06E-11	3.06E-01	2.56E-03	9.92E-01
chr4	142053481	142053783	7.48E-05	3.66E-01	3.79E-04	9.99E-01
chr7	150329304	150329321	1.42E-05	3.40E-01	9.41E-03	9.73E-01
chr19	30018006	30018051	9.22E-06	5.29E-01	1.82E-02	9.67E-01
chr1	181287633	181287642	4.86E-07	4.57E-01	3.60E-02	9.27E-01
chr14	85998084	85998139	9.36E-08	3.71E-01	1.51E-02	9.61E-01
chr14	58332805	58332827	1.66E-08	3.46E-01	1.56E-03	9.95E-01
chr20	37303000	37303007	1.89E-10	3.74E-01	1.79E-03	9.95E-01
chr10	130832367	130832381	2.29E-04	3.86E-01	6.27E-02	8.60E-01
chr7	155556566	155556600	4.72E-07	4.60E-01	5.61E-02	8.91E-01
chr8	142413012	142413121	2.98E-05	3.14E-01	1.42E-02	9.57E-01
chr2	37880456	37880503	1.37E-04	6.04E-01	8.98E-02	8.71E-01
chr16	68390618	68390663	1.88E-04	3.82E-01	5.70E-02	8.70E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr4	111554064	111554077	9.68E-05	6.94E-01	3.40E-02	9.53E-01
chr17	75525139	75525180	3.85E-07	4.23E-01	6.74E-04	9.98E-01
chr2	236574688	236574774	1.46E-06	5.08E-01	4.47E-03	9.91E-01
chr11	1779926	1779973	4.88E-06	4.04E-01	4.30E-02	9.04E-01
chr20	34189501	34189792	1.30E-06	4.64E-01	1.35E-04	1.00E+00
chr2	220313356	220313384	8.74E-10	3.57E-01	1.08E-02	9.71E-01
chr3	196367690	196367895	1.55E-23	4.37E-01	5.65E-03	9.87E-01
chr5	171605385	171605414	5.13E-04	4.87E-01	3.74E-02	9.29E-01
chr9	37903889	37903895	2.18E-07	3.37E-01	3.02E-02	9.18E-01
chr19	58446423	58446430	3.81E-08	3.27E-01	9.55E-03	9.72E-01
chr8	56852201	56852233	5.19E-04	4.61E-01	5.51E-02	8.93E-01
chr8	98290011	98290081	4.79E-09	3.28E-01	4.05E-04	9.99E-01
chr7	28449954	28449970	5.94E-04	3.33E-01	1.22E-03	9.96E-01
chr16	51184771	51184886	4.57E-05	4.04E-01	1.87E-02	9.56E-01
chr11	117186061	117186115	1.02E-06	5.21E-01	2.79E-03	9.95E-01
chr2	133403686	133403797	5.38E-05	4.17E-01	6.77E-02	8.60E-01
chr13	28498448	28498613	1.28E-08	3.31E-01	3.90E-03	9.88E-01
chr19	9517624	9517738	1.44E-05	3.18E-01	7.66E-03	9.76E-01
chr13	58205958	58206058	1.30E-08	6.25E-01	2.51E-03	9.96E-01
chr8	97157855	97157878	5.77E-04	3.25E-01	5.11E-03	9.85E-01
chr12	124950719	124950776	8.41E-05	5.91E-01	6.21E-02	9.05E-01
chr12	5997166	5997213	4.42E-05	4.69E-01	4.88E-02	9.06E-01
chr20	59827586	59827624	2.66E-08	5.32E-01	1.44E-02	9.74E-01
chr11	46367940	46367974	6.13E-04	4.06E-01	4.36E-02	9.03E-01
chr9	95821885	95821901	6.73E-14	4.20E-01	1.20E-02	9.72E-01
chr5	150326192	150326312	8.47E-05	6.67E-01	5.16E-02	9.28E-01
chr12	68043432	68043502	7.15E-04	3.33E-01	1.53E-03	9.95E-01
chr9	134465995	134466015	2.07E-05	6.18E-01	5.18E-02	9.23E-01
chr5	132158827	132158843	1.33E-04	3.64E-01	6.29E-02	8.53E-01
chr5	3606576	3606611	2.32E-06	5.51E-01	3.77E-02	9.36E-01
chr20	56247297	56247315	4.58E-06	4.67E-01	4.68E-02	9.09E-01
chr19	13125032	13125104	2.45E-08	5.10E-01	9.97E-03	9.81E-01
chr9	132891110	132891134	2.12E-04	4.00E-01	2.67E-02	9.37E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr13	95364358	95364385	4.34E-09	6.65E-01	8.16E-03	9.88E-01
chr12	43946100	43946110	3.65E-07	3.08E-01	7.33E-03	9.77E-01
chr22	2355969	23559146	4.14E-16	7.27E-01	3.67E-02	9.52E-01
chr7	156797101	156797116	3.91E-04	3.41E-01	1.50E-02	9.58E-01
chr12	124941370	124941476	1.29E-07	5.15E-01	3.26E-02	9.41E-01
chr1	62106961	62107137	2.69E-04	1.00E+00	1.56E-01	8.65E-01
chr7	70596957	70596971	3.20E-08	3.51E-01	3.16E-03	9.91E-01
chr9	21031614	21031726	3.97E-04	3.40E-01	4.44E-04	9.99E-01
chr6	32813530	32813577	2.53E-08	6.16E-01	4.71E-02	9.29E-01
chr8	140714721	140714739	1.88E-04	3.75E-01	1.97E-02	9.50E-01
chr16	68318872	68318933	7.09E-23	4.58E-01	1.58E-02	9.67E-01
chr12	12163189	12163248	1.45E-06	4.16E-01	1.72E-02	9.60E-01
chr20	31215990	31216058	3.23E-06	4.93E-01	9.67E-02	8.36E-01
chr17	78238406	78238540	4.68E-08	4.89E-01	3.60E-02	9.31E-01
chr8	22547919	22547980	8.90E-04	4.65E-01	1.03E-01	8.19E-01
chr9	2159888	2159931	3.16E-06	3.25E-01	1.23E-02	9.64E-01
chr4	85419917	85419951	5.74E-04	3.94E-01	5.63E-02	8.75E-01
chr12	108169078	108169144	6.83E-07	3.38E-01	7.19E-04	9.98E-01
chr7	45068249	45068308	3.69E-16	4.85E-01	2.03E-02	9.60E-01
chr2	177016575	177016595	1.35E-07	3.21E-01	1.31E-02	9.61E-01
chr3	9595511	9595532	1.10E-06	3.73E-01	2.02E-03	9.95E-01
chr6	90894747	90894799	3.15E-05	4.60E-01	6.88E-02	8.70E-01
chr19	41641865	41641879	5.23E-05	5.00E-01	2.98E-02	9.44E-01
chr12	66276079	66276191	2.18E-09	4.36E-01	7.66E-03	9.83E-01
chr20	61638096	61638108	3.92E-06	7.08E-01	7.34E-02	9.06E-01
chr7	155250669	155250691	1.25E-05	3.78E-01	8.29E-03	9.79E-01
chr2	39187231	39187239	3.25E-08	3.73E-01	3.18E-02	9.21E-01
chr5	116283745	116283908	4.84E-06	4.18E-01	2.02E-02	9.54E-01
chr7	152622786	152622794	1.88E-08	5.37E-01	2.33E-03	9.96E-01
chr5	170735618	170735626	1.34E-10	3.75E-01	2.41E-03	9.94E-01
chr20	41817998	41818052	9.73E-09	3.40E-01	1.10E-02	9.69E-01
chr10	31074290	31074303	1.39E-04	3.70E-01	6.01E-03	9.84E-01
chr13	112726560	112726579	2.92E-07	4.21E-01	1.11E-02	9.74E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr21	45577522	45577578	9.07E-07	4.20E-01	3.57E-02	9.22E-01
chr7	102066465	102066505	4.06E-18	4.54E-01	8.68E-03	9.81E-01
chr17	27912311	27912415	6.05E-07	4.28E-01	5.00E-02	8.95E-01
chr17	80358785	80358913	9.03E-04	5.27E-01	9.93E-02	8.42E-01
chr1	169637173	169637186	2.32E-04	8.33E-01	1.18E-01	8.76E-01
chr2	139537698	139537710	7.21E-06	4.70E-01	1.74E-03	9.96E-01
chr6	10391058	10391093	8.03E-07	5.17E-01	3.30E-02	9.40E-01
chr17	7686714	7686754	6.50E-04	4.76E-01	1.20E-01	7.99E-01
chr12	7068301	7068310	3.54E-07	6.40E-01	6.66E-02	9.06E-01
chr14	70014828	70014873	3.96E-05	3.74E-01	9.82E-04	9.97E-01
chr10	118031044	118031062	1.85E-05	3.05E-01	2.50E-03	9.92E-01
chr16	49518429	49518453	1.01E-04	4.44E-01	4.09E-02	9.16E-01
chr7	71801464	71801484	1.03E-07	3.16E-01	1.18E-02	9.64E-01
chr9	117693070	117693110	1.28E-08	6.67E-01	1.75E-02	9.74E-01
chr5	173315775	173315832	9.76E-04	3.33E-01	3.86E-04	9.99E-01
chr14	103568958	103569000	9.32E-09	4.90E-01	1.09E-02	9.78E-01
chr7	100946173	100946190	4.32E-06	4.36E-01	6.45E-02	8.71E-01
chr19	8656804	8656818	6.75E-05	3.51E-01	3.08E-02	9.19E-01
chr9	139741469	139741493	2.81E-04	5.28E-01	1.73E-01	7.53E-01
chr12	51717855	51717865	8.67E-11	3.52E-01	7.87E-03	9.78E-01
chr5	142785053	142785165	2.93E-04	3.89E-01	8.72E-03	9.78E-01
chr8	41559482	41559550	5.14E-56	9.93E-01	1.34E-03	9.99E-01
chr6	20320223	20320368	3.75E-04	8.33E-01	8.84E-02	9.04E-01
chr6	27114547	27114655	9.26E-05	3.85E-01	7.12E-02	8.44E-01
chr19	17439683	17439702	3.00E-07	3.15E-01	5.19E-03	9.84E-01
chr5	176858541	176858704	1.44E-05	4.04E-01	5.84E-02	8.74E-01
chr9	92099055	92099114	5.71E-14	5.04E-01	8.57E-03	9.83E-01
chr7	24324675	24324683	1.26E-04	3.64E-01	2.16E-02	9.44E-01
chr1	221050180	221050273	1.23E-11	3.91E-01	1.73E-02	9.58E-01
chr5	2748507	2748516	2.48E-13	5.26E-01	2.01E-02	9.63E-01
chr6	143858815	143858914	2.44E-04	3.36E-01	2.30E-02	9.36E-01
chr22	47070535	47070641	1.03E-16	5.93E-01	1.10E-02	9.82E-01
chr19	3672066	3672080	6.98E-09	3.34E-01	2.99E-02	9.18E-01

-86-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr15	90357499	90357521	6.18E-04	3.87E-01	2.77E-02	9.33E-01
chr9	127258187	127258274	1.07E-16	6.99E-01	2.86E-02	9.61E-01
chr7	2271925	2272011	9.71E-04	3.33E-01	2.25E-04	9.99E-01
chr7	71801154	71801190	5.87E-06	3.50E-01	9.54E-03	9.73E-01
chr7	156797306	156797346	1.32E-04	3.89E-01	1.19E-02	9.70E-01
chr15	70550975	70551061	1.30E-08	4.91E-01	3.19E-02	9.39E-01
chr19	49973605	49973632	1.23E-05	5.00E-01	1.69E-02	9.67E-01
chr10	26505199	26505213	2.88E-08	4.72E-01	1.34E-03	9.97E-01
chr2	129037557	129037569	2.16E-05	6.21E-01	1.42E-01	8.14E-01
chr2	46579092	46579122	7.41E-12	3.22E-01	6.32E-03	9.81E-01
chr7	1959874	1959887	2.04E-16	5.51E-01	1.27E-02	9.77E-01
chr13	24839234	24839455	2.35E-08	5.83E-01	4.57E-02	9.27E-01
chr12	22486913	22486995	8.05E-04	3.33E-01	0.00E+00	1.00E+00
chr1	112058209	112058223	1.03E-05	3.84E-01	5.20E-02	8.81E-01
chr9	134148683	134148714	2.91E-09	3.94E-01	2.23E-02	9.47E-01
chr8	132054036	132054055	2.21E-07	3.37E-01	1.44E-02	9.59E-01
chr10	133110551	133110580	3.64E-05	5.00E-01	3.23E-03	9.94E-01
chr1	207842657	207842691	1.96E-06	3.65E-01	3.29E-02	9.17E-01
chr3	195313667	195313820	3.08E-07	6.86E-01	8.41E-02	8.91E-01
chr6	43612899	43612908	5.34E-07	3.10E-01	3.73E-03	9.88E-01
chr2	219827596	219827660	6.20E-04	3.47E-01	4.28E-03	9.88E-01
chr16	85296061	85296086	7.72E-10	4.69E-01	2.30E-02	9.53E-01
chr3	65342371	65342384	6.53E-04	6.51E-01	1.73E-01	7.90E-01
chr16	50308403	50308417	2.47E-11	6.25E-01	3.88E-02	9.42E-01
chr2	182322264	182322274	3.20E-08	3.13E-01	1.08E-02	9.67E-01
chr1	145562949	145562967	5.95E-07	4.17E-01	6.64E-04	9.98E-01
chr4	6697966	6698085	3.52E-04	5.07E-01	1.53E-01	7.68E-01
chr8	128929859	128929964	3.28E-09	5.76E-01	4.10E-03	9.93E-01
chr19	16197083	16197146	6.00E-13	5.41E-01	2.90E-02	9.49E-01
chr1	27160026	27160126	2.07E-08	4.50E-01	4.34E-02	9.12E-01
chr6	127440662	127440735	4.13E-04	3.66E-01	3.15E-02	9.21E-01
chr7	81417609	81417823	9.42E-07	9.21E-01	9.84E-02	9.03E-01
chr19	55485208	55485344	4.56E-05	5.15E-01	7.33E-02	8.75E-01

39595779.1

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr5	140871224	140871238	5.56E-06	6.46E-01	8.89E-02	8.79E-01
chr6	20024224	20024319	1.83E-17	4.53E-01	1.49E-02	9.68E-01
chr11	8290108	8290126	1.97E-07	3.31E-01	1.86E-03	9.94E-01
chr7	56355663	56355674	8.20E-05	3.59E-01	5.34E-02	8.71E-01
chr11	111411737	111411788	8.26E-08	3.18E-01	6.09E-04	9.98E-01
chr11	119580472	119580588	1.30E-04	3.43E-01	3.90E-02	8.98E-01
chr16	31488816	31488849	5.30E-04	4.28E-01	5.63E-02	8.84E-01
chr12	131464602	131464618	2.56E-10	7.08E-01	2.99E-02	9.59E-01
chr10	133528489	133528587	4.09E-05	6.23E-01	8.32E-02	8.82E-01
chr19	29284698	29284703	9.80E-06	4.54E-01	2.27E-02	9.52E-01
chr19	13215304	13215327	7.35E-10	3.38E-01	3.02E-02	9.18E-01
chr4	2814877	2814968	4.30E-06	7.29E-01	6.81E-02	9.15E-01
chr10	7708862	7708891	9.37E-06	3.66E-01	2.67E-03	9.93E-01
chr8	132052589	132052738	4.46E-09	6.09E-01	8.27E-03	9.87E-01
chr8	142339677	142339715	4.70E-06	3.91E-01	2.72E-02	9.35E-01
chr9	87905314	87905325	1.01E-06	3.99E-01	4.55E-02	8.98E-01
chr13	109793272	109793287	2.49E-09	3.29E-01	1.23E-02	9.63E-01
chr1	1476095	1476154	2.66E-06	4.24E-01	2.04E-03	9.95E-01
chr1	170630065	170630080	5.31E-09	5.10E-01	3.25E-02	9.40E-01
chr10	43578285	43578402	5.83E-08	3.56E-01	2.92E-02	9.24E-01
chr21	38597974	38598066	7.35E-05	4.06E-01	4.48E-02	9.01E-01
chr2	105480382	105480390	4.87E-06	4.69E-01	3.29E-02	9.34E-01
chr13	112726834	112726848	2.66E-06	4.95E-01	9.77E-03	9.81E-01
chr17	8090935	8091039	6.78E-04	3.18E-01	1.61E-04	9.99E-01
chr13	43149103	43149114	9.10E-06	4.17E-01	2.58E-03	9.94E-01
chr16	29160257	29160326	2.13E-13	6.15E-01	2.26E-02	9.65E-01
chr4	174450772	174450785	1.01E-05	3.16E-01	5.93E-03	9.82E-01
chr11	70672834	70672890	1.16E-05	5.29E-01	7.20E-02	8.80E-01
chr12	39539256	39539268	2.15E-07	4.42E-01	2.31E-03	9.95E-01
chr14	48143617	48143641	7.10E-10	7.50E-01	2.38E-03	9.97E-01
chr8	898406	898416	1.12E-05	5.80E-01	5.73E-02	9.10E-01
chr20	22557745	22557769	8.90E-05	3.33E-01	6.63E-03	9.81E-01
chr19	3180080	3180104	2.78E-07	4.22E-01	1.17E-02	9.73E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	14739326	14739369	6.29E-05	8.66E-01	2.02E-01	8.11E-01
chr18	20888061	20888186	7.97E-04	3.33E-01	8.28E-03	9.76E-01
chr6	43613038	43613046	6.73E-06	5.18E-01	6.17E-03	9.88E-01
chr19	2085137	2085160	6.66E-08	3.72E-01	2.62E-02	9.34E-01
chr3	67706421	67706705	1.65E-12	5.16E-01	1.30E-02	9.75E-01
chr19	46997441	46997466	1.91E-07	6.67E-01	1.86E-02	9.73E-01
chr21	38069675	38069705	4.64E-05	3.52E-01	2.42E-03	9.93E-01
chr7	24324036	24324045	1.69E-08	3.17E-01	1.30E-02	9.60E-01
chr17	70112343	70112463	3.75E-13	4.69E-01	8.02E-03	9.83E-01
chr2	119602826	119602837	2.45E-06	3.89E-01	1.97E-02	9.52E-01
chr14	91580279	91580507	6.16E-08	5.00E-01	5.73E-03	9.89E-01
chr20	52825234	52825333	6.17E-05	3.61E-01	2.59E-03	9.93E-01
chr7	71800926	71800934	2.65E-07	3.50E-01	3.51E-03	9.90E-01
chr14	70756090	70756120	1.04E-10	3.79E-01	2.14E-02	9.47E-01
chr19	681393	681413	3.32E-06	3.32E-01	2.75E-02	9.23E-01
chr2	242140582	242140589	2.73E-07	6.87E-01	1.12E-01	8.59E-01
chr8	134510944	134511107	6.83E-08	4.05E-01	9.83E-03	9.76E-01
chr9	136710340	136710362	1.05E-14	4.34E-01	1.48E-02	9.67E-01
chr15	31658049	31658168	1.22E-06	5.00E-01	1.27E-03	9.97E-01
chr20	62199474	62199549	7.05E-06	4.58E-01	4.94E-02	9.03E-01
chr12	4647687	4648274	2.79E-04	3.33E-01	6.67E-04	9.98E-01
chr13	102568657	102568663	3.55E-08	3.20E-01	6.51E-03	9.80E-01
chr9	71789124	71789541	1.00E-09	3.43E-01	6.53E-04	9.98E-01
chr12	45444497	45444506	1.32E-07	4.20E-01	2.63E-03	9.94E-01
chr7	73443734	73443767	7.98E-06	5.35E-01	1.86E-02	9.66E-01
chr19	1068467	1068484	1.53E-12	4.70E-01	2.04E-02	9.58E-01
chr14	103740485	103740496	2.57E-09	3.25E-01	8.69E-03	9.74E-01
chr11	91958384	91958465	3.83E-06	3.18E-01	4.60E-03	9.86E-01
chr12	121530481	121530644	3.67E-08	6.58E-01	9.17E-02	8.78E-01
chr10	13701335	13701356	1.20E-10	4.23E-01	4.02E-02	9.13E-01
chr1	33219502	33219520	5.33E-07	4.01E-01	4.21E-02	9.05E-01
chr18	13137599	13137680	6.56E-04	3.33E-01	9.90E-03	9.71E-01
chr13	58204132	58204154	1.24E-06	5.45E-01	1.82E-02	9.68E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr8	66864145	66864243	2.33E-16	6.82E-01	1.40E-02	9.80E-01
chr7	45067930	45067978	1.82E-06	5.63E-01	2.75E-02	9.53E-01
chr16	216025	216056	2.33E-07	5.27E-01	2.38E-04	1.00E+00
chr6	166582771	166582783	2.09E-05	3.52E-01	1.44E-02	9.61E-01
chr19	54464269	54464480	3.86E-05	3.89E-01	9.87E-03	9.75E-01
chr20	10199369	10199501	4.84E-04	3.33E-01	6.32E-04	9.98E-01
chr10	96990630	96990718	7.76E-07	4.81E-01	4.94E-03	9.90E-01
chr6	151004402	151004610	6.41E-18	6.22E-01	3.29E-02	9.50E-01
chr8	57025563	57025929	1.79E-08	4.37E-01	7.51E-04	9.98E-01
chr11	65408444	65408551	3.34E-09	4.28E-01	2.10E-02	9.53E-01
chr11	130271850	130271875	3.99E-04	5.33E-01	8.38E-02	8.64E-01
chr7	644366	644393	8.83E-04	3.67E-01	5.83E-02	8.63E-01
chr4	7651988	7652061	2.29E-06	6.15E-01	3.25E-02	9.50E-01
chr1	47900111	47900152	2.16E-06	5.11E-01	5.26E-02	9.07E-01
chr13	99064107	99064184	7.02E-06	4.36E-01	1.79E-02	9.61E-01
chr2	119916553	119916571	1.67E-10	4.55E-01	3.37E-03	9.93E-01
chr16	54967319	54967361	3.49E-12	5.52E-01	1.20E-02	9.79E-01
chr20	44037518	44037624	5.37E-11	4.80E-01	5.82E-02	8.92E-01
chr4	8263785	8263806	8.13E-04	4.40E-01	1.04E-01	8.09E-01
chr16	2317671	2317694	2.75E-04	4.83E-01	7.04E-02	8.73E-01
chr10	22542557	22542582	2.43E-06	6.13E-01	4.71E-02	9.29E-01
chr12	123518638	123518773	3.48E-06	3.24E-01	2.29E-02	9.34E-01
chr6	20024500	20024569	1.60E-04	5.75E-01	8.49E-02	8.71E-01
chr13	25744912	25744924	1.93E-08	3.81E-01	7.69E-03	9.80E-01
chr2	96987068	96987098	5.17E-06	4.44E-01	9.26E-04	9.98E-01
chr2	144694722	144694809	3.36E-07	4.51E-01	1.80E-04	1.00E+00
chr19	22605150	22605168	6.83E-04	3.48E-01	2.12E-02	9.43E-01
chr8	24813563	24813599	7.79E-04	4.00E-01	6.58E-03	9.84E-01
chr13	30423680	30423763	4.02E-04	3.33E-01	5.14E-04	9.98E-01
chr11	62212143	62212195	1.18E-08	3.59E-01	7.04E-03	9.81E-01
chr13	21050140	21050156	1.56E-06	3.90E-01	1.94E-02	9.53E-01
chr15	84115968	84116104	2.79E-05	3.08E-01	1.96E-04	9.99E-01
chr14	69260675	69260712	1.03E-04	4.00E-01	1.26E-02	9.69E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	27655609	27656143	6.86E-05	4.99E-01	3.74E-02	9.30E-01
chr1	228195435	228195456	6.08E-05	7.33E-01	7.09E-02	9.12E-01
chr5	5140453	5140465	4.25E-04	3.74E-01	2.57E-02	9.36E-01
chr13	109148187	109148195	5.08E-11	3.76E-01	1.20E-03	9.97E-01
chr9	132586568	132586652	2.77E-04	3.33E-01	3.97E-04	9.99E-01
chr10	126300538	126300694	7.36E-07	4.40E-01	4.38E-02	9.09E-01
chr6	37534044	37534257	2.23E-05	3.80E-01	5.34E-03	9.86E-01
chr17	33700564	33700621	1.68E-08	6.67E-01	1.15E-02	9.83E-01
chr11	111411876	111412047	1.58E-08	5.98E-01	2.19E-04	1.00E+00
chr1	26616292	26616311	2.23E-08	7.14E-01	3.35E-02	9.55E-01
chr11	66326731	66326743	4.66E-04	3.87E-01	3.13E-02	9.25E-01
chr5	122422290	122422309	1.28E-09	4.98E-01	3.23E-02	9.39E-01
chr5	44388581	44388799	2.39E-06	5.04E-01	2.62E-03	9.95E-01
chr22	20078423	20078511	8.12E-05	5.99E-01	1.44E-01	8.06E-01
chr2	220313255	220313272	1.28E-07	3.85E-01	3.21E-02	9.23E-01
chr11	67205095	67205113	2.14E-04	5.55E-01	8.17E-02	8.72E-01
chr16	29675880	29675899	1.18E-07	5.85E-01	4.11E-02	9.34E-01
chr4	5053572	5053596	7.98E-04	3.13E-01	6.27E-03	9.80E-01
chr9	14346053	14346085	1.57E-11	4.32E-01	2.77E-02	9.40E-01
chr6	22147356	22147383	7.87E-04	5.06E-01	1.21E-01	8.06E-01
chr19	11708610	11708736	9.24E-04	3.34E-01	2.52E-04	9.99E-01
chr14	75981561	75981659	2.83E-07	3.18E-01	9.15E-03	9.72E-01
chr16	66621433	66621441	1.69E-08	5.77E-01	5.34E-02	9.15E-01
chr7	153749755	153749802	1.13E-07	3.77E-01	8.69E-03	9.77E-01
chr8	55371721	55371728	8.96E-08	3.39E-01	1.02E-02	9.71E-01
chr14	105528116	105528216	2.15E-05	4.04E-01	3.21E-02	9.26E-01
chr22	19711226	19711257	6.48E-16	3.61E-01	4.03E-03	9.89E-01
chr17	643532	643598	6.39E-07	9.52E-01	7.22E-02	9.30E-01
chr6	1625110	1625384	2.40E-11	5.84E-01	1.12E-02	9.81E-01
chr5	178422156	178422206	4.68E-04	3.69E-01	4.56E-02	8.90E-01
chr22	30662059	30662077	1.31E-07	4.30E-01	3.86E-02	9.18E-01
chr20	46044275	46044343	1.68E-06	5.60E-01	4.31E-02	9.28E-01
chr16	88767062	88767096	1.38E-12	6.54E-01	9.07E-04	9.99E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr14	60952393	60952405	1.04E-07	3.85E-01	7.13E-03	9.82E-01
chr2	8715562	8715637	4.34E-09	7.06E-01	4.17E-02	9.44E-01
chr18	2847848	2847923	4.36E-04	3.33E-01	2.20E-03	9.93E-01
chr12	45444152	45444203	2.50E-05	4.17E-01	1.83E-02	9.58E-01
chr8	145900685	145900715	3.40E-06	5.08E-01	6.76E-02	8.83E-01
chr7	107643271	107643389	7.67E-04	3.13E-01	4.56E-04	9.99E-01
chr6	137244466	137244484	1.17E-08	3.48E-01	1.89E-03	9.95E-01
chr4	37585982	37585993	2.98E-07	3.86E-01	6.66E-02	8.53E-01
chr12	2030316	2030331	5.79E-15	4.57E-01	1.73E-02	9.63E-01
chr10	134827218	134827246	6.49E-06	3.69E-01	1.80E-02	9.53E-01
chr6	26722714	26722740	3.16E-08	6.67E-01	4.21E-04	9.99E-01
chr7	644133	644190	4.24E-07	4.04E-01	3.96E-02	9.11E-01
chr5	49736943	49736950	3.36E-05	3.75E-01	2.09E-04	9.99E-01
chr4	3387340	3387383	3.25E-08	3.19E-01	1.68E-02	9.50E-01
chr5	158532301	158532330	7.76E-04	6.67E-01	6.03E-02	9.17E-01
chr1	119543104	119543128	2.91E-06	3.13E-01	2.34E-03	9.93E-01
chr22	42095558	42095578	6.66E-04	4.25E-01	4.76E-02	8.99E-01
chr22	37881408	37881495	3.45E-09	3.90E-01	1.83E-02	9.55E-01
chr19	908886	908901	7.92E-04	3.57E-01	5.16E-02	8.74E-01
chr1	224804868	224804896	3.31E-06	3.45E-01	8.41E-03	9.76E-01
chr21	27945586	27945599	1.13E-13	3.14E-01	8.73E-03	9.73E-01
chr6	150285766	150285775	3.77E-05	3.40E-01	5.25E-03	9.85E-01
chr11	47416527	47416633	3.84E-08	5.75E-01	4.61E-02	9.26E-01
chr16	55365257	55365267	3.14E-05	3.27E-01	1.45E-02	9.58E-01
chr17	79420888	79420920	2.11E-09	4.24E-01	2.06E-02	9.54E-01
chr13	37248265	37248280	3.46E-08	4.58E-01	1.50E-02	9.68E-01
chr19	1467151	1467183	1.87E-06	3.27E-01	3.97E-03	9.88E-01
chr14	106328280	106328426	9.09E-04	3.33E-01	2.49E-03	9.93E-01
chr20	55965241	55965253	5.53E-10	3.29E-01	1.04E-02	9.69E-01
chr3	16883944	16884162	4.70E-08	4.83E-01	3.09E-02	9.40E-01
chr21	45577764	45577826	1.09E-06	6.30E-01	3.40E-02	9.49E-01
chr6	170581003	170581096	1.39E-12	3.34E-01	1.64E-02	9.53E-01
chr1	1564676	1564690	2.23E-04	6.21E-01	1.11E-01	8.49E-01

-92-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr4	1397570	1397650	4.86E-06	3.16E-01	1.04E-02	9.68E-01
chr19	3179929	3179985	9.47E-15	5.51E-01	2.15E-02	9.62E-01
chr19	51951770	51951911	2.51E-06	4.51E-01	2.43E-02	9.49E-01
chr8	687582	687636	1.73E-06	4.98E-01	1.22E-03	9.98E-01
chr1	6086599	6086652	6.50E-04	3.94E-01	4.47E-02	8.98E-01
chr10	106401933	106401941	2.13E-05	4.01E-01	3.34E-02	9.23E-01
chr3	239683	239720	2.25E-05	5.71E-01	3.28E-02	9.46E-01
chr17	73839230	73839273	6.56E-06	9.52E-01	1.25E-01	8.84E-01
chr1	221050478	221050491	8.96E-04	3.91E-01	6.61E-02	8.55E-01
chr11	1332237	1332286	5.34E-04	5.06E-01	6.97E-02	8.79E-01
chr2	73928331	73928354	1.10E-06	5.11E-01	5.85E-02	8.97E-01
chr1	3077720	3077741	6.84E-05	4.70E-01	4.09E-02	9.20E-01
chr2	47206751	47206855	3.65E-06	6.23E-01	9.86E-02	8.63E-01
chr9	29214375	29214383	3.37E-05	3.61E-01	4.00E-03	9.89E-01
chr14	85997836	85997946	5.16E-07	3.48E-01	2.85E-03	9.92E-01
chr22	43627412	43627566	3.71E-05	3.52E-01	3.99E-03	9.89E-01
chr9	71788636	71788870	4.33E-08	3.07E-01	9.42E-04	9.97E-01
chr20	21492917	21492931	4.43E-05	3.88E-01	2.56E-02	9.38E-01
chr5	45696095	45696108	5.24E-10	4.60E-01	8.05E-03	9.83E-01
chr3	123242028	123242130	7.00E-05	5.00E-01	2.80E-02	9.47E-01
chr21	44819374	44819432	2.06E-06	6.48E-01	4.71E-02	9.32E-01
chr2	220313472	220313492	2.04E-16	4.88E-01	2.99E-03	9.94E-01
chr13	23733902	23733940	4.44E-05	3.91E-01	2.28E-03	9.94E-01
chr8	35093160	35093180	2.20E-08	3.10E-01	8.25E-03	9.74E-01
chr7	32338303	32338342	4.92E-07	4.92E-01	1.92E-03	9.96E-01
chr14	50829747	50830001	3.10E-04	4.57E-01	2.06E-02	9.57E-01
chr2	5813752	5813893	3.13E-05	3.24E-01	1.74E-02	9.49E-01
chr17	73839499	73839523	1.41E-09	4.75E-01	8.65E-03	9.82E-01
chr16	12997270	12997416	2.81E-04	3.47E-01	4.97E-03	9.86E-01
chr7	2774625	2774646	1.17E-05	5.44E-01	9.90E-02	8.46E-01
chr20	62200085	62200109	1.92E-05	5.81E-01	1.07E-01	8.44E-01
chr2	8723890	8723900	8.38E-04	5.80E-01	8.54E-02	8.72E-01
chr19	38886138	38886154	2.50E-05	4.26E-01	5.02E-02	8.95E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	1620159	1620204	1.76E-09	3.11E-01	3.46E-03	9.89E-01
chr2	99439256	99439395	6.54E-12	3.56E-01	7.40E-04	9.98E-01
chr6	163670668	163670704	1.36E-07	4.46E-01	3.55E-02	9.26E-01
chr6	20024352	20024422	4.34E-08	3.57E-01	2.60E-02	9.32E-01
chr3	38081466	38081537	7.06E-05	6.67E-01	3.53E-02	9.50E-01
chr10	7452356	7452391	3.84E-07	6.35E-01	2.72E-03	9.96E-01
chr9	117266875	117266928	5.04E-07	3.74E-01	3.04E-02	9.25E-01
chr12	7060734	7060909	2.53E-04	3.29E-01	1.40E-02	9.59E-01
chr19	10870440	10870451	2.31E-07	3.18E-01	1.55E-02	9.54E-01
chr20	3653597	3653644	8.22E-05	3.77E-01	5.74E-04	9.98E-01
chr3	47051337	47051366	3.72E-06	4.49E-01	1.82E-02	9.61E-01
chr16	56313147	56313258	7.19E-08	5.55E-01	8.26E-02	8.70E-01
chr8	65282679	65282693	4.30E-12	3.39E-01	6.09E-03	9.82E-01
chr7	3134380	3134520	7.69E-05	6.69E-01	1.25E-01	8.43E-01
chr12	6664477	6664942	1.86E-11	3.38E-01	1.42E-03	9.96E-01
chr8	54792994	54793073	7.78E-04	3.22E-01	8.14E-03	9.75E-01
chr1	77333197	77333207	4.46E-06	3.25E-01	7.63E-03	9.77E-01
chr20	49158144	49158301	5.89E-04	6.49E-01	1.04E-01	8.61E-01
chr15	79382228	79382307	2.92E-04	3.20E-01	8.81E-03	9.73E-01
chr5	92907787	92907881	1.16E-05	4.53E-01	1.78E-02	9.62E-01
chr1	146369387	146369494	5.36E-05	6.67E-01	8.52E-02	8.87E-01
chr11	123355219	123355301	1.34E-05	3.97E-01	6.22E-02	8.65E-01
chr1	229568006	229568012	1.28E-08	3.49E-01	8.86E-03	9.75E-01
chr19	48918187	48918247	1.48E-06	4.28E-01	1.01E-02	9.77E-01
chr19	42703701	42703778	8.19E-05	4.86E-01	4.96E-02	9.07E-01
chr14	75800718	75800727	6.11E-04	4.17E-01	2.90E-02	9.35E-01
chr4	76860974	76861103	2.44E-04	3.33E-01	1.53E-03	9.95E-01
chr13	110357693	110357853	1.33E-04	3.68E-01	1.92E-02	9.50E-01
chr20	61885738	61885787	1.29E-04	3.33E-01	3.10E-02	9.15E-01
chr20	62199644	62199673	3.69E-16	3.68E-01	1.60E-02	9.58E-01
chr9	132650733	132650800	1.81E-12	3.44E-01	2.01E-02	9.45E-01
chr2	45232449	45232493	2.46E-05	4.83E-01	2.27E-02	9.55E-01
chr4	134072180	134072187	4.01E-04	3.99E-01	3.52E-02	9.19E-01

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr6	20320098	20320141	2.23E-05	5.83E-01	9.09E-02	8.65E-01
chr13	109793301	109793322	3.09E-04	5.05E-01	5.32E-02	9.05E-01
chr8	121021099	121021235	2.99E-07	5.14E-01	2.89E-02	9.47E-01
chr19	3688173	3688181	1.21E-05	5.17E-01	7.75E-02	8.70E-01
chr17	46749238	46749330	3.19E-07	4.89E-01	1.35E-02	9.73E-01
chr10	101089443	101089462	6.37E-08	3.92E-01	1.62E-03	9.96E-01
chr7	4848865	4848900	2.46E-08	3.62E-01	2.40E-02	9.38E-01
chr17	4389757	4389831	2.21E-08	6.67E-01	8.98E-03	9.87E-01
chr10	6094775	6094827	3.21E-08	3.57E-01	4.55E-02	8.87E-01
chr22	50319330	50319373	7.12E-10	6.37E-01	3.88E-02	9.43E-01
chr8	8654900	8654936	2.75E-07	5.07E-01	4.60E-02	9.17E-01
chr11	70864168	70864186	8.09E-04	4.67E-01	3.92E-02	9.23E-01
chr3	194118678	194118738	2.87E-09	4.69E-01	2.16E-03	9.95E-01
chr4	104641038	104641088	1.12E-06	4.28E-01	1.91E-03	9.96E-01
chr22	40796311	40796475	1.68E-18	6.09E-01	1.07E-02	9.83E-01
chr1	221049959	221050007	1.23E-06	5.30E-01	4.38E-02	9.24E-01
chr8	67873654	67873665	5.20E-04	3.82E-01	5.98E-02	8.65E-01
chr11	63637374	63637410	1.66E-09	3.58E-01	1.54E-02	9.59E-01
chr9	36458494	36458734	7.27E-12	4.67E-01	2.69E-02	9.45E-01
chr8	142276277	142276297	3.08E-11	5.17E-01	4.18E-02	9.25E-01
chr7	1552786	1552837	1.45E-04	3.84E-01	4.44E-02	8.96E-01
chr3	133661417	133661556	2.91E-13	3.68E-01	2.48E-02	9.37E-01
chr2	67625058	67625155	9.95E-04	3.33E-01	7.73E-03	9.77E-01
chr15	64244503	64244684	3.63E-09	5.45E-01	7.20E-02	8.83E-01
chr1	57888963	57889003	8.31E-05	5.00E-01	2.32E-02	9.56E-01
chr6	132722899	132722939	3.55E-06	5.14E-01	5.03E-03	9.90E-01
chr8	98290100	98290333	5.19E-07	4.92E-01	4.52E-03	9.91E-01
chr2	109855158	109855264	1.22E-06	4.06E-01	2.56E-02	9.41E-01
chr11	65408370	65408410	4.85E-05	4.04E-01	5.45E-02	8.81E-01
chr19	4916881	4917023	8.52E-08	3.55E-01	2.71E-02	9.29E-01
chr2	239335984	239336031	3.79E-05	3.29E-01	1.84E-03	9.94E-01
chr17	46619163	46619183	1.00E-06	4.34E-01	1.36E-02	9.70E-01
chr17	79134306	79134315	3.45E-04	3.33E-01	6.19E-03	9.82E-01

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr8	41166324	41166334	2.39E-06	3.51E-01	4.22E-02	8.93E-01
chr8	494278	494286	1.09E-07	3.19E-01	8.84E-03	9.73E-01
chr13	25745208	25745232	3.28E-07	3.74E-01	1.56E-02	9.60E-01
chr8	687326	687341	3.48E-07	3.36E-01	3.08E-03	9.91E-01
chr14	73181172	73181203	1.06E-09	5.24E-01	3.07E-02	9.45E-01
chr6	20024006	20024141	2.60E-08	5.71E-01	5.03E-02	9.19E-01
chr7	138780361	138780396	3.52E-06	3.76E-01	7.03E-02	8.42E-01
chr7	155246087	155246111	5.14E-04	3.33E-01	8.53E-03	9.75E-01
chr1	119530576	119530606	2.47E-08	4.16E-01	9.57E-03	9.77E-01
chr16	85936479	85936599	1.14E-04	4.24E-01	4.78E-02	8.99E-01
chr7	45039481	45039607	8.08E-04	3.33E-01	1.27E-04	1.00E+00
chr22	44391923	44391940	1.76E-07	5.42E-01	7.15E-02	8.84E-01
chr5	135266182	135266207	7.63E-09	3.18E-01	3.92E-03	9.88E-01
chr6	42335118	42335231	3.72E-08	4.31E-01	3.88E-02	9.17E-01
chr12	29936617	29936632	3.31E-04	3.53E-01	5.09E-02	8.74E-01
chr7	70597943	70597953	1.18E-08	3.42E-01	2.35E-02	9.36E-01
chr6	126112458	126112707	6.93E-04	3.33E-01	4.53E-03	9.87E-01
chr8	25203033	25203231	2.57E-06	9.33E-01	6.59E-02	9.34E-01
chr2	131797613	131797636	2.63E-11	3.58E-01	2.87E-02	9.26E-01
chr13	112723104	112723111	1.22E-06	3.35E-01	1.37E-02	9.61E-01
chr1	229567768	229567822	1.50E-05	5.89E-01	4.39E-02	9.31E-01
chr10	11212536	11212603	8.90E-08	4.49E-01	2.76E-02	9.42E-01
chr15	74658177	74658190	4.08E-07	4.26E-01	1.06E-03	9.98E-01
chr2	237416280	237416327	1.85E-04	3.33E-01	3.27E-02	9.11E-01
chr12	115109842	115109855	7.17E-14	4.31E-01	1.72E-02	9.62E-01
chr22	39148239	39148256	7.94E-04	3.80E-01	3.62E-02	9.13E-01
chr19	50058236	50058275	1.06E-05	3.46E-01	4.48E-02	8.85E-01
chr6	37551119	37551243	1.16E-06	8.06E-01	8.25E-02	9.07E-01
chr7	36679543	36679660	2.97E-04	4.79E-01	9.14E-02	8.40E-01
chr1	1118430	1118437	3.47E-13	6.30E-01	3.08E-02	9.53E-01
chr16	54967023	54967077	4.46E-09	6.76E-01	1.08E-02	9.84E-01
chr17	16570355	16570361	6.01E-04	4.97E-01	1.51E-01	7.67E-01
chr6	166976776	166976796	7.60E-05	4.65E-01	8.63E-02	8.44E-01

-96-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr8	57026038	57026118	7.09E-23	8.86E-01	3.66E-03	9.96E-01
chr19	41119578	41119638	2.01E-04	3.57E-01	1.11E-03	9.97E-01
chr8	127568853	127569070	9.32E-04	4.16E-01	2.62E-02	9.41E-01
chr2	66672481	66672520	8.77E-12	3.23E-01	1.74E-02	9.49E-01
chr9	134128101	134128135	2.83E-18	4.42E-01	5.26E-03	9.88E-01
chr6	133562773	133562790	5.91E-09	6.67E-01	2.38E-03	9.96E-01
chr12	41086129	41086227	1.86E-09	6.63E-01	2.39E-03	9.96E-01
chr3	129694365	129694378	1.80E-07	3.98E-01	1.99E-03	9.95E-01
chr22	46367952	46367968	2.53E-08	3.22E-01	4.24E-03	9.87E-01
chr6	17962707	17962816	3.67E-05	3.45E-01	4.31E-02	8.89E-01
chr16	54970335	54970379	1.79E-09	3.10E-01	4.58E-03	9.85E-01
chr19	37742927	37742978	5.47E-07	4.50E-01	2.20E-02	9.53E-01
chr5	45696035	45696048	1.19E-07	3.42E-01	1.72E-02	9.52E-01
chr2	239140304	239140340	1.20E-05	3.13E-01	2.90E-03	9.91E-01
chr2	220349250	220349265	7.57E-08	6.67E-01	3.87E-04	9.99E-01
chr7	102715190	102715298	1.57E-04	3.67E-01	1.44E-02	9.62E-01
chr15	79382801	79382824	6.68E-04	3.47E-01	2.70E-03	9.92E-01
chr6	3752576	3752604	7.33E-04	3.33E-01	2.13E-03	9.94E-01
chr14	67955224	67955250	2.41E-05	4.57E-01	6.28E-02	8.79E-01
chr10	11207159	11207259	2.57E-08	4.54E-01	5.02E-04	9.99E-01
chr19	3179353	3179432	3.37E-09	5.34E-01	2.26E-02	9.59E-01
chr11	67171082	67171115	6.93E-05	6.65E-01	1.49E-01	8.17E-01
chr9	95728320	95728335	8.40E-11	3.15E-01	6.35E-03	9.80E-01
chr21	44819309	44819337	2.01E-10	4.87E-01	5.37E-02	9.01E-01
chr19	47152996	47153006	1.79E-08	4.92E-01	3.17E-02	9.39E-01
chr1	149223680	149223845	3.59E-04	3.33E-01	6.77E-03	9.80E-01
chr11	63974228	63974540	7.17E-13	5.04E-01	9.13E-03	9.82E-01
chr20	3653378	3653398	1.02E-12	4.55E-01	9.84E-03	9.79E-01
chr7	2758845	2758885	4.20E-04	6.25E-01	8.22E-02	8.84E-01
chr2	198063780	198063957	1.05E-05	4.24E-01	4.26E-02	9.09E-01
chr19	19571691	19571731	1.02E-10	3.44E-01	2.02E-02	9.45E-01
chr6	166970525	166970559	6.04E-04	3.33E-01	8.25E-03	9.76E-01
chr8	97506696	97507459	1.06E-09	4.29E-01	2.07E-04	1.00E+00

39595779.1

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr17	21415130	21415142	1.30E-06	7.23E-01	8.87E-02	8.91E-01
chr13	109792984	109792996	5.24E-07	3.45E-01	2.80E-02	9.25E-01
chr3	50312991	50313021	1.89E-10	3.16E-01	7.01E-03	9.78E-01
chr22	47070183	47070285	2.86E-09	3.75E-01	2.54E-02	9.37E-01
chr14	70654420	70654435	2.18E-08	6.67E-01	1.50E-03	9.98E-01
chr13	25745075	25745153	6.52E-07	3.13E-01	6.84E-04	9.98E-01
chr19	37498702	37498719	6.82E-07	3.95E-01	2.68E-02	9.36E-01
chr9	134148582	134148640	3.05E-19	6.79E-01	9.44E-03	9.86E-01
chr4	1210659	1210715	7.42E-16	4.74E-01	1.36E-02	9.72E-01
chr4	1210785	1210901	1.09E-05	4.87E-01	2.80E-02	9.46E-01
chr9	1045544	1045552	2.57E-04	3.75E-01	1.25E-02	9.68E-01
chr12	115109640	115109662	1.71E-13	5.53E-01	2.10E-02	9.64E-01
chr13	28498832	28498881	5.02E-08	4.38E-01	7.14E-03	9.84E-01
chr6	85483796	85483914	4.12E-05	3.46E-01	2.47E-02	9.33E-01
chr12	115109783	115109814	2.69E-06	5.32E-01	4.73E-02	9.18E-01
chr19	39826266	39826272	6.52E-20	4.91E-01	6.18E-03	9.88E-01
chr7	156798241	156798342	9.72E-04	3.33E-01	2.90E-04	9.99E-01
chr1	91183790	91184000	4.36E-10	4.55E-01	4.86E-03	9.89E-01
chr7	2672631	2672761	1.11E-04	4.05E-01	2.50E-03	9.94E-01
chr7	156795190	156795205	1.01E-05	3.93E-01	4.41E-02	8.99E-01
chr8	28637012	28637063	4.74E-05	6.42E-01	1.63E-01	7.97E-01
chr13	25320305	25320319	4.63E-09	3.42E-01	2.22E-03	9.94E-01
chr8	86350666	86350689	1.38E-04	3.85E-01	2.72E-02	9.34E-01
chr12	51718198	51718251	1.92E-06	6.90E-01	1.05E-01	8.68E-01
chr6	117086335	117086374	1.36E-09	5.62E-01	1.21E-02	9.79E-01
chr15	96904744	96904779	2.79E-10	3.39E-01	2.01E-02	9.44E-01
chr13	112758478	112758496	2.98E-05	5.00E-01	5.11E-03	9.90E-01
chr17	17295442	17295619	7.27E-12	4.01E-01	1.59E-02	9.62E-01
chr16	28996532	28996602	9.67E-07	5.09E-01	3.23E-02	9.40E-01
chr13	112758807	112758835	4.04E-08	3.10E-01	5.53E-03	9.82E-01
chr14	99706032	99706158	1.18E-04	3.69E-01	2.88E-02	9.28E-01
chr14	36987163	36987184	9.34E-11	3.42E-01	6.76E-03	9.81E-01
chr2	43398069	43398085	2.49E-11	4.51E-01	4.46E-02	9.10E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	90504290	90504359	1.81E-06	7.29E-01	1.66E-01	8.14E-01
chr13	28503165	28503235	9.87E-04	3.76E-01	3.78E-02	9.09E-01
chr3	154797891	154798152	1.40E-06	3.28E-01	1.28E-04	1.00E+00
chr6	28367069	28367128	9.50E-13	6.30E-01	9.85E-03	9.85E-01
chr1	40974762	40974913	1.66E-04	3.33E-01	1.10E-03	9.97E-01
chr22	50319158	50319194	3.73E-04	6.52E-01	1.15E-01	8.50E-01
chr12	115103722	115103767	7.11E-09	3.86E-01	1.94E-02	9.52E-01
chr5	122435250	122435295	9.37E-06	4.75E-01	3.19E-03	9.93E-01
chr10	134599210	134599245	4.40E-05	3.37E-01	2.06E-03	9.94E-01
chr7	117854567	117854818	1.23E-05	5.35E-01	2.16E-02	9.61E-01
chr15	96890743	96890812	7.68E-10	3.52E-01	1.92E-02	9.48E-01
chr11	45114836	45114845	8.70E-07	5.36E-01	6.09E-02	8.98E-01
chr1	151811102	151811203	1.24E-08	3.25E-01	7.38E-04	9.98E-01
chr7	633049	633101	6.28E-10	4.28E-01	2.77E-02	9.39E-01
chr4	148402315	148402431	3.63E-05	4.18E-01	5.27E-03	9.88E-01
chr3	50313747	50313762	1.37E-12	4.38E-01	9.92E-03	9.78E-01
chr6	42072513	42072527	5.78E-08	3.21E-01	1.58E-03	9.95E-01
chr1	32237851	32237877	5.70E-10	6.67E-01	9.60E-04	9.99E-01
chr11	60718854	60718885	8.52E-08	3.54E-01	3.38E-03	9.91E-01
chr4	5894624	5894634	1.68E-04	4.60E-01	5.81E-02	8.88E-01
chr2	105478832	105478873	2.70E-05	3.94E-01	6.33E-03	9.84E-01
chr15	96890847	96890880	5.34E-04	3.33E-01	1.03E-02	9.70E-01
chr9	101822013	101822096	1.36E-04	4.60E-01	6.09E-02	8.83E-01
chr11	65661484	65661564	1.68E-08	3.69E-01	1.70E-02	9.56E-01
chr5	177989059	177989077	2.47E-11	5.23E-01	3.12E-02	9.44E-01
chr1	40781681	40781778	3.33E-04	3.33E-01	3.30E-03	9.90E-01
chr13	28498910	28498956	3.03E-09	3.33E-01	2.40E-03	9.93E-01
chr14	52294765	52294902	1.93E-08	5.15E-01	3.38E-02	9.38E-01
chr3	47029500	47029628	1.24E-04	3.57E-01	2.57E-02	9.33E-01
chr7	151328862	151328895	6.68E-05	3.76E-01	0.00E+00	1.00E+00
chr19	55553070	55553087	6.55E-16	6.19E-01	3.45E-02	9.47E-01
chr17	40822427	40822434	5.35E-05	6.03E-01	1.03E-01	8.54E-01
chr19	3178517	3178545	2.45E-17	4.00E-01	1.59E-02	9.62E-01

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr11	76299657	76299682	3.10E-06	7.94E-01	5.62E-02	9.34E-01
chr13	29106898	29106944	4.99E-06	3.15E-01	1.04E-02	9.68E-01
chr20	36533265	3653330	3.02E-11	3.49E-01	1.03E-02	9.71E-01
chr2	102187418	102187570	1.31E-04	6.28E-01	1.22E-01	8.37E-01
chr11	32355271	32355288	4.37E-08	3.56E-01	1.54E-03	9.96E-01
chr1	2706134	2706161	1.10E-07	3.08E-01	3.31E-03	9.89E-01
chr16	58497351	58497369	1.26E-04	4.61E-01	5.66E-02	8.91E-01
chr17	56409588	56409742	1.53E-04	4.96E-01	5.48E-02	9.00E-01
chr6	50787949	50787984	4.81E-09	4.53E-01	3.52E-02	9.28E-01
chr20	61992292	61992315	5.29E-09	3.92E-01	1.86E-03	9.95E-01
chr16	68390251	68390361	7.18E-17	4.71E-01	2.22E-02	9.55E-01
chr22	19710841	19710855	3.13E-06	3.94E-01	4.28E-03	9.89E-01
chr13	109793122	109793152	6.05E-07	3.08E-01	2.56E-03	9.92E-01
chr1	4715727	4715734	6.53E-06	6.02E-01	3.67E-02	9.43E-01
chr16	88496790	88496819	2.91E-04	3.40E-01	3.41E-02	9.09E-01
chr8	105478845	105478857	7.40E-10	3.48E-01	1.02E-02	9.72E-01
chr11	67204881	67205042	6.69E-04	3.75E-01	5.62E-02	8.70E-01
chr10	12223668	12223715	4.13E-08	4.51E-01	7.50E-02	8.57E-01
chr5	169724354	169724470	2.48E-04	4.26E-01	3.10E-02	9.32E-01

Table 4b:
Differentially methylated MHB regions between lung cancer tissues (LCT) and normal plasma.

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr3	187676563	187676644	1.78E-18	3.70E-01	6.60E-04	9.98E-01
chr7	391332	391394	2.80E-05	4.15E-01	2.80E-02	9.37E-01
chr8	11565988	11565997	9.49E-06	3.59E-01	3.91E-02	9.02E-01
chr7	104581511	104581826	4.32E-09	6.73E-01	6.48E-02	9.12E-01
chr7	5371156	5371187	7.04E-10	4.23E-01	1.16E-02	9.73E-01
chr21	38080525	38080550	9.91E-07	3.67E-01	7.82E-03	9.79E-01
chr8	17509194	17509223	2.08E-12	6.94E-01	3.56E-02	9.51E-01
chr5	140207807	140207820	1.13E-10	5.27E-01	5.73E-02	9.02E-01
chr20	62369207	62369221	5.29E-06	4.11E-01	0.00E+00	1.00E+00
chr17	79406164	79406253	6.94E-09	4.67E-01	8.66E-02	8.44E-01

-100-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr16	29164276	29164303	5.35E-05	7.19E-01	1.71E-01	8.08E-01
chr2	177017367	177017378	1.11E-05	5.37E-01	4.78E-03	9.91E-01
chr21	47971700	47971722	6.33E-05	6.49E-01	9.40E-02	8.73E-01
chr1	218098391	218098398	7.78E-07	3.22E-01	1.47E-02	9.56E-01
chr1	19251999	19252025	1.40E-06	7.94E-01	2.08E-01	7.92E-01
chr8	96085373	96085420	6.00E-07	7.50E-01	7.16E-02	9.13E-01
chr19	18113497	18113578	1.06E-05	6.39E-01	4.65E-02	9.32E-01
chr3	65342396	65342412	4.86E-09	3.34E-01	3.36E-02	9.09E-01
chr14	105640730	105640770	2.39E-05	5.71E-01	6.76E-02	8.94E-01
chr15	96865945	96865955	3.15E-06	3.85E-01	7.33E-03	9.81E-01
chr12	125039261	125039301	4.07E-08	4.82E-01	2.69E-02	9.47E-01
chr2	112124096	112124151	7.06E-08	5.36E-01	3.50E-02	9.39E-01
chr2	232186706	232186761	4.03E-05	6.99E-01	1.17E-01	8.57E-01
chr4	2814622	2814718	5.99E-04	3.94E-01	4.03E-02	9.07E-01
chr2	171570261	171570269	8.68E-09	3.71E-01	1.07E-02	9.72E-01
chr3	138663980	138664002	4.72E-08	3.89E-01	7.93E-03	9.80E-01
chr22	36773395	36773570	2.26E-08	4.53E-01	3.57E-02	9.27E-01
chr17	40822613	40822653	5.50E-05	4.87E-01	1.05E-01	8.23E-01
chr19	12306143	12306166	1.34E-04	3.85E-01	4.09E-03	9.89E-01
chr22	20783998	20784048	2.77E-09	3.56E-01	5.28E-02	8.71E-01
chr1	154375706	154375773	5.12E-13	6.97E-01	1.67E-02	9.77E-01
chr17	76127253	76127408	1.33E-22	5.19E-01	1.32E-02	9.75E-01
chr1	45009127	45009259	4.98E-04	5.01E-01	1.03E-01	8.30E-01
chr18	72916692	72916705	3.71E-04	5.30E-01	9.87E-02	8.43E-01
chr2	109834010	109834148	1.32E-05	5.00E-01	4.07E-02	9.25E-01
chr4	8082086	8082177	1.19E-05	5.85E-01	1.48E-01	7.98E-01
chr6	41515921	41515938	3.84E-08	3.06E-01	2.81E-03	9.91E-01
chr19	10624949	10624963	3.00E-07	3.35E-01	8.91E-03	9.74E-01
chr14	94254076	94254092	1.90E-08	3.55E-01	2.20E-02	9.42E-01
chr16	31488890	31488912	1.14E-06	5.07E-01	4.44E-02	9.19E-01
chr17	74141537	74141566	1.61E-14	3.84E-01	1.43E-02	9.64E-01
chr21	44819099	44819123	2.08E-11	5.22E-01	4.98E-02	9.13E-01
chr15	27136915	27136924	3.07E-04	3.89E-01	3.66E-02	9.14E-01

39595779.1

-101-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	10393363	10393497	3.04E-04	7.22E-01	9.25E-02	8.86E-01
chr2	73147431	73147444	3.25E-06	4.33E-01	1.05E-02	9.76E-01
chr19	48946520	48946657	1.80E-05	3.24E-01	1.39E-03	9.96E-01
chr5	50694887	50694985	1.26E-04	3.83E-01	8.31E-03	9.79E-01
chr8	38757858	38757899	2.69E-14	3.74E-01	1.47E-02	9.62E-01
chr9	126785062	126785097	2.43E-05	5.00E-01	6.59E-03	9.87E-01
chr8	142276522	142276533	5.10E-10	4.30E-01	1.94E-02	9.57E-01
chr16	88758463	88758479	1.78E-11	8.05E-01	1.47E-02	9.82E-01
chr10	114136302	114136365	8.24E-11	5.72E-01	8.45E-03	9.85E-01
chr11	67804601	67804688	2.54E-05	3.32E-01	5.21E-03	9.85E-01
chr2	160761118	160761165	2.90E-08	5.97E-01	0.00E+00	1.00E+00
chr9	127265806	127265816	4.03E-06	4.42E-01	1.61E-02	9.65E-01
chr1	151870591	151870681	6.08E-04	8.44E-01	2.31E-01	7.85E-01
chr21	37757965	37758006	5.82E-04	3.33E-01	1.27E-03	9.96E-01
chr19	3179842	3179904	9.81E-04	3.67E-01	4.09E-02	9.00E-01
chr5	110406515	110406530	8.65E-10	3.81E-01	3.67E-02	9.12E-01
chr19	58629492	58629522	1.39E-04	3.66E-01	4.52E-02	8.90E-01
chr21	46334450	46334460	3.21E-05	5.04E-01	1.06E-01	8.26E-01
chr4	184828032	184828058	6.44E-15	3.83E-01	6.36E-03	9.84E-01
chr12	64784380	64784575	1.47E-10	5.89E-01	6.57E-04	9.99E-01
chr2	70057016	70057043	3.20E-06	5.84E-01	6.67E-02	8.97E-01
chr17	40822041	40822061	4.17E-05	6.03E-01	1.43E-01	8.08E-01
chr2	25425395	25425417	3.58E-05	7.18E-01	1.72E-01	8.07E-01
chr16	68028052	68028131	2.34E-12	5.09E-01	3.45E-02	9.37E-01
chr2	177016618	177016669	3.62E-04	3.58E-01	1.44E-02	9.61E-01
chr16	88769941	88770131	2.70E-07	3.11E-01	6.08E-03	9.81E-01
chr3	48520626	48520665	3.50E-14	6.36E-01	2.73E-02	9.59E-01
chr3	170137672	170137687	6.55E-07	5.19E-01	5.62E-03	9.89E-01
chr14	73712647	73712670	4.07E-07	3.61E-01	5.16E-03	9.86E-01
chr17	59482675	59482690	1.11E-04	4.21E-01	2.63E-02	9.41E-01
chr11	69229767	69229875	8.65E-07	4.01E-01	5.93E-02	8.71E-01
chr7	4248707	4248866	2.55E-07	5.66E-01	4.92E-02	9.20E-01
chr2	66666590	66666648	1.13E-10	3.44E-01	7.08E-03	9.80E-01

39595779.1

-102-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr8	145900816	145901119	3.14E-20	4.38E-01	1.22E-02	9.73E-01
chr5	139227605	139227613	5.76E-04	3.33E-01	1.59E-03	9.95E-01
chr17	72462820	72462912	5.48E-08	3.35E-01	1.52E-02	9.57E-01
chr2	235055127	235055264	6.84E-04	6.67E-01	6.33E-02	9.13E-01
chr5	2740503	2740527	6.31E-07	3.19E-01	2.80E-03	9.91E-01
chr17	1588162	1588384	8.00E-04	3.33E-01	8.14E-05	1.00E+00
chr5	3600942	3600993	2.60E-11	4.95E-01	2.38E-03	9.95E-01
chr2	129492325	129492360	7.82E-10	4.17E-01	2.83E-02	9.37E-01
chr3	44063876	44063942	7.92E-05	3.69E-01	2.83E-02	9.29E-01
chr11	69920045	69920065	2.23E-11	6.87E-01	4.04E-02	9.44E-01
chr19	55766296	55766320	3.46E-06	4.63E-01	3.13E-02	9.37E-01
chr2	239755874	239755896	5.19E-07	4.11E-01	2.82E-02	9.36E-01
chr12	62584918	62584946	1.95E-04	3.57E-01	5.97E-03	9.84E-01
chr6	157557373	157557415	9.37E-07	3.55E-01	1.69E-03	9.95E-01
chr19	58868331	58868347	8.87E-06	4.91E-01	1.11E-01	8.16E-01
chr12	115103838	115103858	4.15E-08	5.28E-01	1.08E-02	9.80E-01
chr7	150329304	150329321	3.02E-06	3.87E-01	9.41E-03	9.76E-01
chr6	29521498	29521580	2.88E-08	3.76E-01	2.16E-03	9.94E-01
chr3	184504534	184504593	1.14E-06	4.37E-01	3.32E-02	9.29E-01
chr6	45631538	45631561	1.78E-09	4.37E-01	1.35E-03	9.97E-01
chr6	137814560	137814567	2.67E-08	3.32E-01	2.44E-02	9.31E-01
chr2	220313636	220313657	5.91E-09	4.57E-01	5.31E-03	9.89E-01
chr12	111099505	111099564	8.70E-09	5.20E-01	1.34E-02	9.75E-01
chr10	130832367	130832381	4.88E-05	4.15E-01	6.27E-02	8.69E-01
chr7	155556566	155556600	8.23E-10	6.17E-01	5.61E-02	9.17E-01
chr20	62473619	62473818	3.02E-06	7.50E-01	6.71E-02	9.18E-01
chr6	27835219	27835269	1.85E-06	3.44E-01	6.72E-03	9.81E-01
chr6	158411773	158411853	2.67E-08	8.08E-01	1.75E-01	8.22E-01
chr8	142413012	142413121	1.55E-08	5.00E-01	1.42E-02	9.72E-01
chr5	1878568	1878601	5.87E-05	4.16E-01	3.34E-02	9.26E-01
chr12	113675732	113675977	4.20E-04	5.67E-01	7.33E-02	8.85E-01
chr16	68390618	68390663	1.58E-04	4.09E-01	5.70E-02	8.78E-01
chr20	22562975	22562988	3.53E-05	6.56E-01	5.58E-02	9.22E-01

39595779.1

-103-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr17	44847184	44847220	1.20E-04	5.04E-01	1.07E-01	8.26E-01
chr16	50715337	50715362	2.27E-05	5.83E-01	5.81E-02	9.09E-01
chr2	19556014	19556079	3.21E-06	4.09E-01	6.03E-02	8.72E-01
chr4	107774	107793	3.69E-04	3.42E-01	4.92E-03	9.86E-01
chr17	79315001	79315013	1.43E-08	5.28E-01	2.49E-02	9.55E-01
chr13	19918950	19918984	2.81E-07	3.96E-01	2.96E-02	9.30E-01
chr15	57072647	57072802	1.03E-05	8.20E-01	1.26E-01	8.67E-01
chr6	45631260	45631269	6.27E-09	3.54E-01	3.37E-03	9.91E-01
chr11	1779926	1779973	3.24E-06	4.44E-01	4.30E-02	9.12E-01
chr17	44846986	44847044	7.65E-06	3.85E-01	5.77E-02	8.70E-01
chr16	50308506	50308653	1.06E-05	6.67E-01	4.90E-02	9.32E-01
chr2	220313356	220313384	8.12E-12	4.40E-01	1.08E-02	9.76E-01
chr3	196367690	196367895	4.42E-16	3.92E-01	5.65E-03	9.86E-01
chr5	171605385	171605414	4.30E-05	4.69E-01	3.74E-02	9.26E-01
chr9	37903889	37903895	2.04E-10	4.45E-01	3.02E-02	9.36E-01
chr5	172659760	172659766	1.65E-09	5.74E-01	8.65E-03	9.85E-01
chr19	3369814	3369834	1.95E-05	4.64E-01	7.57E-02	8.60E-01
chr3	128211534	128211604	4.81E-06	4.33E-01	3.14E-02	9.32E-01
chr2	133403686	133403797	1.24E-06	5.08E-01	6.77E-02	8.82E-01
chr5	3594724	3594745	6.75E-06	5.00E-01	2.09E-02	9.60E-01
chr8	61764645	61764655	6.49E-04	5.04E-01	1.40E-01	7.83E-01
chr8	145106961	145106989	5.57E-05	3.89E-01	2.20E-02	9.46E-01
chr17	79450423	79450451	1.04E-04	4.18E-01	8.46E-02	8.31E-01
chr2	45231787	45231831	1.06E-05	4.20E-01	1.79E-02	9.59E-01
chr14	52534718	52534943	6.22E-04	3.33E-01	2.25E-03	9.93E-01
chr12	5997166	5997213	1.87E-05	4.83E-01	4.88E-02	9.08E-01
chr15	96886880	96886915	7.19E-04	5.00E-01	5.04E-02	9.08E-01
chr22	50050779	50050844	1.32E-06	4.67E-01	7.48E-02	8.62E-01
chr17	36105480	36105490	1.07E-06	3.29E-01	7.20E-03	9.79E-01
chr22	46367343	46367373	1.39E-08	4.42E-01	8.10E-03	9.82E-01
chr22	19753445	19753480	6.98E-07	3.17E-01	9.14E-03	9.72E-01
chr9	95821885	95821901	4.75E-17	4.25E-01	1.20E-02	9.73E-01
chr10	77167577	77167584	2.17E-07	3.47E-01	4.30E-04	9.99E-01

39595779.1

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr16	28505658	28505862	7.74E-15	3.98E-01	2.18E-02	9.48E-01
chr7	121950428	121950564	1.32E-06	3.47E-01	1.15E-02	9.68E-01
chr6	107956258	107956267	4.29E-11	3.91E-01	2.30E-02	9.44E-01
chr8	1058991	10589016	3.90E-13	3.27E-01	8.12E-03	9.76E-01
chr9	110350145	110350159	5.33E-05	4.62E-01	9.62E-02	8.28E-01
chr9	134465995	134466015	4.92E-05	5.38E-01	5.18E-02	9.12E-01
chr8	145901158	145901388	1.52E-11	6.37E-01	3.42E-02	9.49E-01
chr15	74832704	74833132	8.57E-04	3.33E-01	1.42E-03	9.96E-01
chr6	30431723	30431770	3.39E-11	7.14E-01	9.37E-03	9.87E-01
chr5	3606576	3606611	1.00E-04	4.21E-01	3.77E-02	9.18E-01
chr20	56247297	56247315	1.93E-08	5.30E-01	4.68E-02	9.19E-01
chr5	76373633	76373641	9.52E-06	3.73E-01	5.07E-02	8.80E-01
chr17	68164684	68164749	2.93E-04	5.00E-01	2.80E-02	9.47E-01
chr13	79181662	79181784	1.22E-05	6.39E-01	4.53E-02	9.34E-01
chr6	1384169	1384189	1.20E-06	4.57E-01	2.20E-02	9.54E-01
chr22	23558969	23559146	3.12E-04	3.87E-01	3.67E-02	9.13E-01
chr5	3600677	3600690	2.47E-12	3.26E-01	1.04E-02	9.69E-01
chr17	73749653	73749685	1.44E-09	3.09E-01	4.38E-03	9.86E-01
chr4	24801685	24801698	1.75E-13	3.88E-01	6.79E-03	9.83E-01
chr10	77167421	77167439	5.92E-06	3.92E-01	4.57E-02	8.96E-01
chr4	111543271	111543299	4.88E-05	4.03E-01	4.87E-03	9.88E-01
chr17	73483828	73483897	2.89E-10	3.47E-01	1.46E-02	9.60E-01
chr12	124941370	124941476	1.77E-04	3.93E-01	3.26E-02	9.23E-01
chr22	24181054	24181071	7.17E-07	4.84E-01	2.29E-02	9.55E-01
chr19	23185803	23185818	7.63E-04	3.47E-01	1.35E-02	9.62E-01
chr2	73114219	73114237	2.09E-06	4.40E-01	3.28E-02	9.31E-01
chr16	68318872	68318933	5.46E-18	3.84E-01	1.58E-02	9.60E-01
chr20	31215990	31216058	4.38E-10	7.35E-01	9.67E-02	8.84E-01
chr20	21486625	21486638	6.78E-04	3.37E-01	1.34E-02	9.62E-01
chr6	117869097	117869149	9.80E-07	6.53E-01	8.17E-02	8.89E-01
chr5	177210046	177210086	2.14E-08	6.67E-01	1.74E-02	9.74E-01
chr16	56669343	56669402	4.45E-11	3.35E-01	1.28E-02	9.63E-01
chr9	124982160	124982166	2.60E-11	3.65E-01	1.32E-02	9.65E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr1	9910785	9910800	2.92E-05	6.74E-01	1.20E-01	8.49E-01
chr8	22547919	22547980	1.22E-05	5.76E-01	1.03E-01	8.48E-01
chr21	38077472	38077482	5.22E-11	3.18E-01	2.52E-03	9.92E-01
chr4	2801307	2801417	2.43E-05	6.33E-01	5.75E-02	9.17E-01
chr2	10471620	10471684	1.69E-14	3.91E-01	2.08E-02	9.50E-01
chr1	179713235	179713255	1.10E-08	3.63E-01	2.14E-02	9.44E-01
chr7	45068249	45068308	1.41E-10	3.33E-01	2.03E-02	9.43E-01
chr2	177016575	177016595	7.60E-06	3.48E-01	1.31E-02	9.64E-01
chr1	109633869	109633984	4.97E-04	3.33E-01	4.50E-03	9.87E-01
chr12	66276079	66276191	3.58E-06	3.37E-01	7.66E-03	9.78E-01
chr1	246745796	246745918	5.00E-05	7.93E-01	2.06E-01	7.94E-01
chr21	35321105	35321239	8.41E-04	3.33E-01	7.56E-03	9.78E-01
chr17	36666183	36666252	1.94E-13	3.79E-01	1.76E-02	9.56E-01
chr1	161039226	161039378	3.02E-04	4.05E-01	5.05E-02	8.89E-01
chr11	114010500	114010527	2.54E-07	5.11E-01	1.20E-01	8.10E-01
chr19	55765372	55765427	4.24E-08	3.53E-01	5.98E-03	9.83E-01
chr3	5061424	5061583	3.40E-04	5.22E-01	8.07E-02	8.66E-01
chr7	102066465	102066505	1.89E-13	4.22E-01	8.68E-03	9.80E-01
chr17	27912311	27912415	6.66E-10	4.04E-01	5.00E-02	8.90E-01
chr17	80358785	80358913	3.51E-04	5.81E-01	9.93E-02	8.54E-01
chr1	169637173	169637186	3.30E-04	8.33E-01	1.18E-01	8.76E-01
chr12	22094366	22094373	7.59E-09	3.72E-01	2.13E-02	9.46E-01
chr6	10391058	10391093	1.93E-04	4.86E-01	3.30E-02	9.37E-01
chr19	4108798	4108818	2.16E-04	5.96E-01	9.90E-02	8.58E-01
chr1	78512134	78512213	6.75E-04	3.47E-01	1.55E-02	9.57E-01
chr2	105459552	105459571	1.32E-05	3.20E-01	6.70E-03	9.79E-01
chr13	28498669	28498711	1.21E-06	3.35E-01	4.67E-03	9.86E-01
chr17	76886714	76886754	1.27E-04	5.30E-01	1.20E-01	8.16E-01
chr7	155595704	155595716	1.75E-05	4.20E-01	7.48E-02	8.49E-01
chr4	186560264	186560311	4.04E-06	3.95E-01	3.89E-02	9.10E-01
chr19	49338972	49339092	2.41E-05	5.02E-01	1.74E-03	9.97E-01
chr5	54527322	54527355	4.00E-06	3.75E-01	3.64E-03	9.90E-01
chr17	59482595	59482615	5.63E-09	3.23E-01	2.20E-02	9.36E-01

-106-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr8	141849327	141849357	6.86E-08	7.37E-01	6.36E-02	9.21E-01
chr17	1960955	1960972	7.77E-05	4.36E-01	8.82E-02	8.32E-01
chr6	42072312	42072322	1.24E-06	3.51E-01	2.83E-03	9.92E-01
chr3	170137395	170137413	1.94E-08	4.11E-01	1.00E-02	9.76E-01
chr17	7287595	7287642	2.52E-04	4.10E-01	2.09E-02	9.51E-01
chr14	103568958	103569000	1.27E-05	3.86E-01	1.09E-02	9.73E-01
chr2	71116046	71116096	1.97E-09	3.11E-01	2.26E-04	9.99E-01
chr8	142046260	142046307	6.96E-23	4.62E-01	1.10E-02	9.77E-01
chr12	51717855	51717865	2.91E-13	4.38E-01	7.87E-03	9.82E-01
chr2	177030210	177030276	3.93E-09	3.23E-01	2.26E-03	9.93E-01
chr9	138999207	138999217	5.27E-08	3.93E-01	2.42E-02	9.42E-01
chr6	27114547	27114655	4.10E-06	4.50E-01	7.12E-02	8.63E-01
chr5	176858541	176858704	6.13E-05	3.82E-01	5.84E-02	8.67E-01
chr9	92099055	92099114	4.20E-12	3.95E-01	8.57E-03	9.79E-01
chr2	71116232	71116257	9.65E-10	3.18E-01	2.33E-03	9.93E-01
chr6	100895678	100895697	3.33E-11	3.61E-01	9.92E-03	9.73E-01
chr5	132158710	132158741	1.04E-07	4.24E-01	5.06E-02	8.93E-01
chr17	33815330	33815432	1.51E-06	4.03E-01	3.10E-03	9.92E-01
chr6	26123664	26123979	3.92E-06	4.44E-01	1.92E-04	1.00E+00
chr7	24324675	24324683	5.83E-06	4.17E-01	2.16E-02	9.51E-01
chr1	221050180	221050273	2.00E-07	3.21E-01	1.73E-02	9.49E-01
chr5	2748507	2748516	1.13E-09	3.44E-01	2.01E-02	9.45E-01
chr7	3134644	3134670	1.87E-05	5.31E-01	6.62E-02	8.89E-01
chr16	55362634	55362703	3.10E-14	3.68E-01	1.34E-02	9.65E-01
chr22	47070535	47070641	8.91E-19	6.30E-01	1.10E-02	9.83E-01
chr12	96564033	96564157	7.95E-06	6.65E-01	1.53E-01	8.13E-01
chr19	3672066	3672080	1.72E-14	4.49E-01	2.99E-02	9.38E-01
chr12	111845175	111845270	9.62E-10	4.32E-01	4.87E-02	8.99E-01
chr7	100465077	100465161	6.37E-05	3.08E-01	7.24E-03	9.77E-01
chr3	15314976	15315000	6.16E-04	8.61E-01	2.36E-01	7.85E-01
chr3	186490345	186490366	2.87E-05	4.72E-01	1.83E-02	9.63E-01
chr9	127258187	127258274	1.55E-12	4.75E-01	2.86E-02	9.43E-01
chr7	156797306	156797346	2.76E-05	4.79E-01	1.19E-02	9.76E-01

39595779.1

-107-

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr19	14089841	14089870	5.49E-04	3.85E-01	6.06E-02	8.64E-01
chr10	134046101	134046188	5.32E-10	3.87E-01	2.63E-02	9.36E-01
chr4	3304775	3304980	1.63E-05	5.79E-01	3.69E-02	9.40E-01
chr2	129037557	129037569	1.39E-05	6.98E-01	1.42E-01	8.31E-01
chr7	1959874	1959887	3.73E-14	3.58E-01	1.27E-02	9.66E-01
chr6	26216260	26216306	9.79E-04	3.33E-01	5.18E-04	9.98E-01
chr22	50986956	50986986	2.70E-08	3.27E-01	1.02E-02	9.70E-01
chr1	112058209	112058223	1.41E-07	4.88E-01	5.20E-02	9.04E-01
chr22	19510978	19511029	5.11E-09	3.18E-01	4.50E-04	9.99E-01
chr11	639650	639673	8.45E-07	3.32E-01	1.25E-03	9.96E-01
chr2	71134585	71134608	1.74E-09	3.81E-01	5.34E-03	9.86E-01
chr6	43612899	43612908	1.19E-07	3.75E-01	3.73E-03	9.90E-01
chr18	13641703	13641723	5.87E-05	4.62E-01	8.62E-02	8.43E-01
chr16	85296061	85296086	5.83E-07	4.38E-01	2.30E-02	9.50E-01
chr16	50308403	50308417	5.12E-08	4.96E-01	3.88E-02	9.27E-01
chr19	2252604	2252633	1.77E-08	3.14E-01	1.01E-02	9.69E-01
chr2	176948107	176948120	1.78E-09	4.05E-01	9.63E-03	9.77E-01
chr2	66673634	66673719	9.71E-06	4.40E-01	2.95E-02	9.37E-01
chr19	16197083	16197146	1.06E-14	4.42E-01	2.90E-02	9.38E-01
chr7	139553288	139553487	9.78E-06	4.81E-01	8.12E-02	8.56E-01
chr3	128210683	128210713	6.79E-09	3.10E-01	5.29E-03	9.83E-01
chr2	162280538	162280561	9.52E-06	3.82E-01	2.47E-02	9.39E-01
chr15	70767238	70767293	3.17E-05	3.38E-01	1.22E-02	9.65E-01
chr1	27160026	27160126	1.21E-06	4.16E-01	4.34E-02	9.05E-01
chr3	184405603	184405680	5.42E-12	3.81E-01	1.36E-02	9.65E-01
chr19	58220236	58220403	6.79E-09	4.90E-01	2.92E-03	9.94E-01
chr9	99639697	99639735	1.16E-06	3.84E-01	3.91E-03	9.90E-01
chr20	57224759	57224776	2.24E-08	5.00E-01	2.33E-03	9.95E-01
chr8	99961491	99961523	7.30E-05	3.99E-01	3.46E-02	9.20E-01
chr5	92908162	92908213	3.04E-10	4.59E-01	4.24E-03	9.91E-01
chr5	157346588	157346611	1.33E-04	5.51E-01	9.27E-02	8.56E-01
chr19	13210490	13210503	1.73E-17	4.84E-01	8.76E-03	9.82E-01
chr5	140871224	140871238	1.16E-05	6.25E-01	8.89E-02	8.75E-01

39595779.1

-108-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	20024224	20024319	1.55E-27	5.62E-01	1.49E-02	9.74E-01
chr4	186560344	186560355	6.39E-05	3.58E-01	3.22E-02	9.17E-01
chr9	971674	971703	3.57E-08	4.07E-01	1.71E-02	9.60E-01
chr11	44340718	44340785	6.21E-05	3.13E-01	1.19E-02	9.63E-01
chr19	35505322	35505399	1.68E-06	8.67E-01	1.14E-01	8.84E-01
chr8	145104394	145104455	5.65E-09	6.61E-01	2.24E-03	9.97E-01
chr16	31488816	31488849	1.16E-04	4.59E-01	5.63E-02	8.91E-01
chr12	131464602	131464618	7.53E-07	5.53E-01	2.99E-02	9.49E-01
chr12	47701048	47701077	3.17E-08	4.59E-01	5.69E-02	8.90E-01
chr12	115138274	115138280	7.75E-06	3.78E-01	4.41E-02	8.95E-01
chr1	112058303	112058337	1.91E-06	3.19E-01	1.43E-02	9.57E-01
chr17	7287359	7287475	4.51E-06	3.11E-01	7.76E-03	9.76E-01
chr7	156799150	156799168	2.16E-07	3.62E-01	1.38E-02	9.63E-01
chr19	13215304	13215327	1.47E-10	3.85E-01	3.02E-02	9.27E-01
chr3	196367554	196367677	6.53E-13	5.40E-01	1.73E-02	9.69E-01
chr11	67186574	67186616	5.15E-05	4.45E-01	4.76E-02	9.03E-01
chr17	35299478	35299512	6.24E-04	3.96E-01	4.70E-02	8.94E-01
chr9	87905314	87905325	2.86E-13	5.79E-01	4.55E-02	9.27E-01
chr13	109793272	109793287	2.83E-11	3.16E-01	1.25E-02	9.62E-01
chr17	47816851	47816941	4.67E-11	5.31E-01	3.06E-02	9.45E-01
chr1	170630065	170630080	1.01E-06	4.08E-01	3.23E-02	9.26E-01
chr5	180230752	180230844	1.58E-04	3.25E-01	1.79E-02	9.48E-01
chr8	38627595	38627699	6.68E-04	3.33E-01	7.09E-03	9.79E-01
chr21	38597974	38598066	1.66E-07	6.26E-01	4.48E-02	9.33E-01
chr20	24759511	24759539	8.74E-12	3.26E-01	2.22E-02	9.36E-01
chr1	47899663	47899729	1.11E-09	4.20E-01	4.39E-02	9.05E-01
chr2	73147905	73147942	7.49E-09	4.59E-01	1.10E-03	9.98E-01
chr6	152631073	152631120	4.98E-04	4.33E-01	1.08E-01	8.00E-01
chr9	137663233	137663244	2.27E-07	8.13E-01	7.45E-02	9.16E-01
chr8	898406	898416	2.44E-05	5.51E-01	5.73E-02	9.06E-01
chr20	22557745	22557769	3.71E-06	4.34E-01	6.63E-03	9.85E-01
chr19	3180080	3180104	6.72E-13	3.96E-01	1.17E-02	9.71E-01
chr19	3369701	3369714	2.94E-04	3.86E-01	6.95E-02	8.47E-01

39595779.1

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr2	10471490	10471518	1.87E-05	4.11E-01	3.32E-02	9.25E-01
chr3	67706421	67706705	3.76E-21	4.93E-01	1.30E-02	9.74E-01
chr11	20618888	20618918	7.45E-04	3.75E-01	1.64E-02	9.58E-01
chr2	63282949	63282956	6.75E-09	3.31E-01	3.88E-03	9.88E-01
chr21	38069675	38069705	3.68E-06	3.55E-01	2.42E-03	9.93E-01
chr2	119602826	119602837	1.92E-06	4.21E-01	1.97E-02	9.55E-01
chr8	22398142	22398178	4.91E-06	5.45E-01	2.24E-02	9.61E-01
chr19	46319060	46319085	1.23E-06	6.90E-01	8.18E-02	8.94E-01
chr7	105319545	105319563	2.55E-05	5.84E-01	5.15E-02	9.19E-01
chr14	70756090	70756120	4.49E-11	3.87E-01	2.14E-02	9.48E-01
chr2	242140582	242140589	1.62E-04	5.30E-01	1.12E-01	8.25E-01
chr15	65689324	65689361	2.30E-04	3.75E-01	3.57E-02	9.13E-01
chr19	1132141	1132179	5.08E-05	3.93E-01	2.87E-02	9.32E-01
chr8	101661869	101661916	5.00E-04	3.33E-01	4.99E-03	9.85E-01
chr8	134510944	134511107	1.72E-08	4.17E-01	9.85E-03	9.77E-01
chr9	136710340	136710362	7.64E-17	4.27E-01	1.48E-02	9.66E-01
chr20	62199474	62199549	5.19E-05	4.27E-01	4.94E-02	8.96E-01
chr15	78632860	78633002	4.41E-04	3.33E-01	7.41E-03	9.78E-01
chr12	121530481	121530644	9.49E-08	6.00E-01	9.17E-02	8.67E-01
chr10	13701335	13701356	5.10E-10	4.18E-01	4.02E-02	9.12E-01
chr8	61764675	61764770	9.93E-05	7.12E-01	1.62E-01	8.15E-01
chr3	195488911	195488924	5.58E-05	3.17E-01	1.23E-02	9.63E-01
chr1	200842879	200842912	6.99E-05	5.09E-01	3.64E-02	9.33E-01
chr3	24536601	24536664	3.17E-05	3.25E-01	7.85E-03	9.76E-01
chr10	21789267	21789303	1.68E-06	4.11E-01	4.30E-03	9.90E-01
chr8	66864145	66864243	1.03E-16	5.13E-01	1.40E-02	9.73E-01
chr7	45067930	45067978	2.72E-16	8.50E-01	2.75E-02	9.69E-01
chr5	37834828	37834844	4.04E-08	3.28E-01	1.46E-02	9.57E-01
chr6	166582771	166582783	4.20E-04	3.22E-01	1.44E-02	9.57E-01
chr15	76635175	76635200	2.75E-09	3.15E-01	2.22E-03	9.93E-01
chr19	57182975	57183010	4.12E-07	4.54E-01	7.38E-03	9.84E-01
chr8	84051808	84052120	1.59E-04	3.55E-01	4.31E-02	8.92E-01
chr4	85418694	85418729	1.11E-08	3.45E-01	5.22E-03	9.85E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr12	125039318	125039343	9.92E-04	3.91E-01	5.70E-02	8.73E-01
chr2	66673956	66673991	2.04E-04	5.39E-01	9.85E-02	8.45E-01
chr6	151004402	151004610	9.19E-13	4.51E-01	3.29E-02	9.32E-01
chr1	6054866	6054949	1.28E-05	4.51E-01	4.24E-02	9.14E-01
chr17	46691788	46691802	2.77E-08	5.83E-01	2.92E-04	9.99E-01
chr11	65408444	65408551	4.29E-11	5.41E-01	2.10E-02	9.63E-01
chr12	115103886	115103900	1.27E-23	3.69E-01	1.83E-03	9.95E-01
chr11	130271850	130271875	2.00E-04	5.73E-01	8.38E-02	8.72E-01
chr8	22409072	22409160	2.06E-04	5.56E-01	4.93E-02	9.18E-01
chr5	72715504	72715573	1.58E-08	4.21E-01	1.48E-02	9.66E-01
chr18	46307700	46307753	5.34E-06	3.83E-01	6.43E-02	8.56E-01
chr17	40441433	40441470	3.71E-13	4.62E-01	9.83E-03	9.79E-01
chr13	26625651	26625662	7.65E-04	3.33E-01	6.37E-03	9.81E-01
chr1	157164679	157164693	2.77E-06	4.02E-01	3.69E-02	9.16E-01
chr11	15134983	15135110	8.60E-04	4.93E-01	6.85E-02	8.78E-01
chr20	44037518	44037624	6.73E-16	6.65E-01	5.82E-02	9.20E-01
chr12	49627611	49627624	3.73E-04	7.93E-01	2.23E-01	7.80E-01
chr4	8263785	8263806	6.39E-04	4.65E-01	1.04E-01	8.17E-01
chr9	126776168	126776183	2.84E-08	4.29E-01	7.07E-03	9.84E-01
chr16	85644561	85644674	8.30E-04	3.33E-01	8.82E-04	9.97E-01
chr10	22542557	22542582	3.45E-06	5.81E-01	4.71E-02	9.25E-01
chr6	106434088	106434131	1.94E-13	3.33E-01	9.45E-04	9.97E-01
chr13	112707953	112708006	6.09E-06	3.18E-01	7.61E-03	9.77E-01
chr12	49689395	49689400	5.07E-13	4.45E-01	2.60E-02	9.45E-01
chr19	18760938	18760989	5.44E-05	4.69E-01	1.06E-01	8.15E-01
chr12	123518638	123518773	3.09E-08	5.40E-01	2.29E-02	9.59E-01
chr17	14203010	14203043	4.56E-04	3.33E-01	6.05E-03	9.82E-01
chr16	86321477	86321575	1.67E-06	4.58E-01	3.66E-02	9.26E-01
chr21	38081090	38081187	3.33E-15	7.28E-01	8.76E-04	9.99E-01
chr17	1960261	1960274	3.34E-09	5.56E-01	1.35E-02	9.76E-01
chr6	20024500	20024569	4.22E-05	6.27E-01	8.49E-02	8.81E-01
chr19	38705029	38705203	2.18E-04	4.26E-01	4.58E-02	9.03E-01
chr14	73712457	73712489	5.67E-09	4.72E-01	1.85E-02	9.62E-01

-111-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr3	186685648	186685718	5.86E-06	5.25E-01	6.56E-02	8.89E-01
chr9	98111682	98111739	4.95E-07	3.32E-01	6.94E-04	9.98E-01
chr6	43192429	43192519	5.70E-07	7.64E-01	3.74E-02	9.53E-01
chr2	45029881	45029891	2.09E-07	3.64E-01	1.09E-02	9.71E-01
chr7	105319436	105319486	3.03E-08	3.21E-01	1.55E-02	9.54E-01
chr17	56565349	56565374	2.46E-04	4.73E-01	7.89E-02	8.57E-01
chr13	21050140	21050156	1.78E-08	4.15E-01	1.94E-02	9.55E-01
chr10	11704590	11704643	2.22E-06	4.40E-01	6.40E-02	8.73E-01
chr8	143813256	143813294	5.40E-05	3.74E-01	1.92E-02	9.51E-01
chr12	52627050	52627068	1.16E-05	5.33E-01	3.87E-02	9.32E-01
chr6	27655609	27656143	1.40E-04	4.64E-01	3.74E-02	9.25E-01
chr7	129781353	129781457	5.91E-11	3.51E-01	5.57E-03	9.84E-01
chr16	85344383	85344563	8.47E-06	3.90E-01	5.91E-02	8.69E-01
chr1	36788051	36788066	3.37E-04	6.67E-01	1.78E-01	7.89E-01
chr10	126300538	126300694	4.40E-09	6.00E-01	4.38E-02	9.32E-01
chr2	242785210	242785261	2.72E-07	5.37E-01	1.92E-03	9.96E-01
chr7	42267805	42267822	7.40E-07	3.55E-01	2.20E-03	9.94E-01
chr5	122422290	122422309	6.93E-09	4.33E-01	3.23E-02	9.31E-01
chr12	52115509	52115548	7.49E-09	3.44E-01	1.67E-02	9.54E-01
chr10	131757046	131757057	6.65E-06	3.37E-01	8.49E-03	9.75E-01
chr8	57026239	57026410	1.93E-08	3.28E-01	4.56E-03	9.86E-01
chr22	20078423	20078511	6.16E-06	6.86E-01	1.44E-01	8.27E-01
chr2	220313255	220313272	4.03E-08	3.68E-01	3.21E-02	9.20E-01
chr11	67205095	67205113	3.79E-07	7.63E-01	8.17E-02	9.03E-01
chr14	37136001	37136019	9.85E-06	4.30E-01	2.22E-02	9.51E-01
chr16	29675880	29675899	4.12E-06	4.79E-01	4.11E-02	9.21E-01
chr9	14346053	14346085	4.08E-17	5.81E-01	2.77E-02	9.55E-01
chr16	66621433	66621441	2.35E-09	6.75E-01	5.34E-02	9.27E-01
chr3	128210809	128210826	6.33E-04	4.86E-01	4.03E-02	9.23E-01
chr8	55371721	55371728	4.25E-07	3.40E-01	1.02E-02	9.71E-01
chr10	22624363	22624510	1.05E-07	3.33E-01	2.06E-04	9.99E-01
chr22	20785437	20785462	1.26E-08	3.66E-01	1.84E-04	9.99E-01
chr5	42952051	42952113	2.23E-05	4.92E-01	3.72E-02	9.30E-01

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr9	1045691	1045747	2.51E-04	4.08E-01	3.90E-02	9.13E-01
chr5	177379756	177379834	2.91E-04	1.00E+00	2.08E-01	8.28E-01
chr1	225655594	225655619	4.66E-05	4.43E-01	2.76E-02	9.41E-01
chr3	184504253	184504484	6.14E-10	5.79E-01	3.15E-02	9.48E-01
chr22	19711226	19711257	3.96E-10	7.14E-01	4.03E-03	9.94E-01
chr12	110236002	110236058	8.63E-04	4.97E-01	3.84E-02	9.28E-01
chr10	106401268	106401303	4.47E-04	3.77E-01	3.02E-02	9.26E-01
chr18	19756867	19756897	2.80E-07	3.80E-01	3.11E-02	9.24E-01
chr22	30662059	30662077	6.90E-07	4.11E-01	3.86E-02	9.14E-01
chr10	105344492	105344526	2.12E-07	3.88E-01	5.35E-03	9.86E-01
chr20	46044275	46044343	2.36E-05	4.75E-01	4.31E-02	9.17E-01
chr16	88767062	88767096	1.34E-09	4.25E-01	9.07E-04	9.98E-01
chr14	60952393	60952405	6.09E-06	4.61E-01	7.13E-03	9.85E-01
chr2	8715562	8715637	5.45E-06	5.34E-01	4.17E-02	9.28E-01
chr7	2757237	2757317	1.64E-04	3.73E-01	6.19E-02	8.58E-01
chr8	145900685	145900715	8.83E-08	5.91E-01	6.76E-02	8.97E-01
chr3	128765137	128765158	1.05E-04	4.28E-01	5.16E-02	8.92E-01
chr6	137244466	137244484	7.76E-08	3.36E-01	1.89E-03	9.94E-01
chr6	170581242	170581310	2.46E-10	5.39E-01	2.58E-02	9.54E-01
chr4	103790673	103790911	2.67E-07	6.67E-01	7.22E-05	1.00E+00
chr4	37585982	37585993	9.95E-14	5.31E-01	6.66E-02	8.88E-01
chr12	2030316	2030331	1.86E-14	4.37E-01	1.73E-02	9.62E-01
chr22	39147922	39148082	4.71E-06	3.96E-01	1.43E-02	9.65E-01
chr19	18770323	18770328	3.19E-04	4.00E-01	5.61E-02	8.77E-01
chr16	56672432	56672464	2.32E-04	4.80E-01	2.86E-02	9.44E-01
chr5	37834847	37834913	7.29E-05	4.37E-01	2.96E-02	9.37E-01
chr6	27832225	27832236	2.55E-04	3.67E-01	2.26E-02	9.42E-01
chr5	72677925	72677938	2.97E-09	3.99E-01	2.66E-02	9.38E-01
chr5	134386305	134386316	3.59E-06	3.38E-01	2.43E-03	9.93E-01
chr3	138664166	138664281	5.46E-05	5.81E-01	2.87E-02	9.53E-01
chr6	1393814	1393833	7.38E-10	4.27E-01	1.51E-02	9.66E-01
chr5	92908827	92908853	4.44E-05	6.67E-01	3.43E-02	9.51E-01
chr19	4374698	4374799	6.90E-07	5.17E-01	9.54E-02	8.44E-01

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr22	37881408	37881495	4.03E-13	4.79E-01	1.83E-02	9.63E-01
chr12	7284161	7284199	2.20E-04	5.37E-01	3.59E-02	9.37E-01
chr20	56194053	56194195	7.05E-07	3.39E-01	2.48E-02	9.32E-01
chr11	47416527	47416633	2.70E-07	5.00E-01	4.61E-02	9.16E-01
chr16	11327169	11327181	1.05E-04	5.04E-01	5.02E-02	9.09E-01
chr2	162280416	162280438	9.67E-04	3.59E-01	2.67E-02	9.31E-01
chr17	36718317	36718336	5.27E-07	3.49E-01	1.43E-02	9.61E-01
chr14	102555012	102555047	2.08E-17	6.89E-01	1.60E-02	9.77E-01
chr19	16482409	16482500	6.18E-07	3.27E-01	1.56E-02	9.54E-01
chr5	72678105	72678117	1.78E-08	4.18E-01	2.73E-02	9.39E-01
chr2	118593925	118593956	1.39E-04	3.60E-01	1.81E-02	9.52E-01
chr19	36247848	36247862	2.78E-12	4.57E-01	1.49E-02	9.68E-01
chr21	45577764	45577826	8.01E-07	6.26E-01	3.40E-02	9.49E-01
chr18	77199431	77199443	2.12E-09	4.78E-01	4.53E-02	9.13E-01
chr7	632944	632992	5.74E-04	5.05E-01	9.13E-02	8.47E-01
chr6	170581003	170581096	4.18E-11	3.19E-01	1.64E-02	9.51E-01
chr1	1564676	1564690	6.52E-05	6.73E-01	1.11E-01	8.59E-01
chr19	3179929	3179985	9.94E-18	4.80E-01	2.15E-02	9.57E-01
chr8	687582	687636	7.09E-06	4.90E-01	1.22E-03	9.98E-01
chr3	47122364	47122371	1.86E-04	1.00E+00	3.51E-01	7.40E-01
chr9	96715493	96715509	1.50E-08	3.83E-01	5.39E-03	9.86E-01
chr1	221050478	221050491	9.83E-04	4.24E-01	6.61E-02	8.65E-01
chr2	219738428	219738541	2.68E-10	3.04E-01	5.43E-04	9.98E-01
chr2	73928331	73928354	2.38E-06	5.81E-01	5.85E-02	9.09E-01
chr16	85095505	85095535	2.38E-04	8.33E-01	1.39E-01	8.57E-01
chr2	47206751	47206855	2.92E-08	7.65E-01	9.86E-02	8.86E-01
chr3	186490394	186490603	2.13E-05	4.22E-01	1.49E-02	9.66E-01
chr20	21492917	21492931	1.71E-05	3.40E-01	2.56E-02	9.30E-01
chr5	45696095	45696108	2.15E-09	3.50E-01	8.03E-03	9.78E-01
chr21	44819374	44819432	5.77E-05	5.13E-01	4.71E-02	9.16E-01
chr1	200842812	200842842	4.08E-09	4.93E-01	3.95E-02	9.26E-01
chr2	220313472	220313492	1.71E-14	3.99E-01	2.99E-03	9.93E-01
chr9	134128181	134128299	2.18E-09	6.25E-01	1.31E-02	9.79E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr9	96713741	96713787	2.78E-05	3.96E-01	1.28E-02	9.69E-01
chr17	75231575	75231615	3.46E-06	3.95E-01	1.28E-02	9.69E-01
chr17	73839499	73839523	3.40E-10	6.60E-01	8.65E-03	9.87E-01
chr10	118900000	118900044	2.78E-07	3.85E-01	8.75E-03	9.78E-01
chr17	3814755	3814785	8.76E-07	3.79E-01	2.69E-02	9.34E-01
chr18	11948030	11948038	8.74E-12	5.24E-01	3.45E-02	9.38E-01
chr20	62200085	62200109	1.57E-05	6.19E-01	1.07E-01	8.52E-01
chr4	102711896	102711974	3.50E-06	4.00E-01	3.42E-02	9.21E-01
chr19	38886138	38886154	9.69E-11	6.95E-01	5.02E-02	9.33E-01
chr4	81189711	81189729	5.42E-10	3.67E-01	1.25E-02	9.67E-01
chr6	1620159	1620204	1.80E-11	3.09E-01	3.46E-03	9.89E-01
chr14	74981388	74981406	4.39E-07	4.44E-01	1.46E-02	9.68E-01
chr9	5841793	5841956	8.50E-06	5.20E-01	6.28E-02	8.92E-01
chr7	27197885	27197933	7.44E-06	5.00E-01	3.43E-02	9.36E-01
chr6	166970726	166970758	3.17E-04	3.63E-01	6.28E-02	8.53E-01
chr6	163670668	163670704	1.04E-08	4.62E-01	3.55E-02	9.29E-01
chr13	79182034	79182044	1.84E-08	3.68E-01	1.59E-02	9.59E-01
chr6	20024352	20024422	2.85E-11	4.82E-01	2.60E-02	9.49E-01
chr15	99993103	99993143	4.02E-05	4.33E-01	8.52E-02	8.36E-01
chr3	148336056	148336126	4.03E-04	4.44E-01	4.12E-02	9.15E-01
chr3	38081466	38081537	6.00E-07	7.88E-01	3.53E-02	9.57E-01
chr7	156796586	156796676	3.68E-06	4.53E-01	1.55E-03	9.97E-01
chr16	56313147	56313258	8.53E-05	4.20E-01	8.26E-02	8.36E-01
chr1	155265254	155265284	8.52E-05	3.47E-01	3.08E-02	9.19E-01
chr9	69197259	69197281	1.39E-04	4.56E-01	1.14E-01	8.00E-01
chr7	3134380	3134520	4.96E-05	6.92E-01	1.25E-01	8.47E-01
chr16	86613029	86613062	1.75E-05	3.55E-01	8.16E-03	9.78E-01
chr22	19711923	19711938	3.15E-05	4.44E-01	2.46E-02	9.48E-01
chr16	54324380	54324392	7.32E-06	3.70E-01	9.13E-03	9.76E-01
chr20	43945282	43945292	2.51E-06	3.42E-01	3.14E-02	9.16E-01
chr11	123355219	123355301	3.66E-07	4.78E-01	6.22E-02	8.85E-01
chr7	149487306	149487319	1.16E-04	6.15E-01	4.45E-02	9.32E-01
chr7	27209526	27209582	2.14E-08	3.02E-01	6.59E-04	9.98E-01

-115-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	19692010	19692112	1.87E-05	3.58E-01	2.68E-02	9.30E-01
chr2	43020374	43020394	1.67E-07	4.38E-01	6.83E-03	9.85E-01
chr3	157824208	157824217	3.09E-06	3.25E-01	2.11E-03	9.94E-01
chr19	23254232	23254336	2.70E-07	4.44E-01	3.97E-02	9.18E-01
chr2	69064579	69064611	1.60E-04	4.61E-01	4.45E-02	9.12E-01
chr20	62199644	62199673	8.34E-19	5.12E-01	1.60E-02	9.70E-01
chr6	20320098	20320141	2.19E-05	5.24E-01	9.09E-02	8.52E-01
chr8	121021099	121021235	6.37E-10	6.44E-01	2.89E-02	9.57E-01
chr8	57069774	57069786	8.00E-10	3.24E-01	5.67E-03	9.83E-01
chr21	46902922	46902964	5.56E-04	6.01E-01	1.36E-01	8.16E-01
chr17	36666562	36666605	5.63E-09	4.16E-01	5.53E-03	9.87E-01
chr19	3688173	3688181	3.58E-06	5.02E-01	7.75E-02	8.66E-01
chr7	139876374	139876418	2.37E-04	5.00E-01	4.02E-02	9.26E-01
chr9	87905154	87905186	4.51E-04	5.16E-01	4.44E-02	9.21E-01
chr8	10588937	10588948	9.19E-13	6.06E-01	1.99E-03	9.97E-01
chr6	100901494	100901618	1.48E-05	4.66E-01	7.50E-02	8.61E-01
chr11	1332088	1332114	4.20E-04	6.27E-01	1.76E-01	7.81E-01
chr10	6094775	6094827	6.74E-10	4.20E-01	4.55E-02	9.02E-01
chr22	50319330	50319373	5.31E-11	5.65E-01	3.88E-02	9.36E-01
chr8	8654900	8654936	1.41E-07	4.79E-01	4.60E-02	9.12E-01
chr3	194118678	194118738	1.32E-05	3.26E-01	2.16E-03	9.93E-01
chr22	40796311	40796475	4.96E-13	6.13E-01	1.07E-02	9.83E-01
chr1	22104959	221050007	8.05E-04	3.52E-01	4.38E-02	8.89E-01
chr11	63637374	63637410	8.65E-10	3.30E-01	1.54E-02	9.55E-01
chr9	36458494	36458734	2.10E-11	4.02E-01	2.69E-02	9.37E-01
chr8	142276277	142276297	1.64E-09	5.06E-01	4.18E-02	9.24E-01
chr7	1552786	1552837	2.77E-06	5.42E-01	4.44E-02	9.24E-01
chr3	133661417	133661556	1.52E-13	3.89E-01	2.48E-02	9.40E-01
chr15	64244503	64244684	6.72E-13	6.85E-01	7.20E-02	9.05E-01
chr4	184828385	184828403	6.43E-06	3.74E-01	5.40E-02	8.74E-01
chr8	98290100	98290333	4.63E-04	3.59E-01	4.52E-03	9.88E-01
chr14	38080503	38080552	1.33E-09	3.71E-01	2.16E-02	9.45E-01
chr13	31480961	31480979	5.67E-04	3.61E-01	5.41E-02	8.69E-01

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr11	67171584	67171666	3.87E-08	3.65E-01	1.59E-02	9.58E-01
chr12	22094594	22094657	3.92E-08	3.77E-01	1.52E-02	9.61E-01
chr11	65408370	65408410	2.33E-11	7.04E-01	5.45E-02	9.28E-01
chr18	13641584	13641601	2.56E-04	3.68E-01	6.46E-02	8.51E-01
chr19	4916881	4917023	1.31E-09	4.65E-01	2.71E-02	9.45E-01
chr1	210391015	210391076	2.04E-04	8.33E-01	2.28E-01	7.85E-01
chr14	38080557	38080579	1.56E-04	4.78E-01	8.62E-02	8.47E-01
chr17	46619163	46619183	1.23E-21	4.90E-01	1.36E-02	9.73E-01
chr12	52652454	52652489	9.69E-11	3.65E-01	4.24E-03	9.89E-01
chr2	171786280	171786350	1.69E-06	5.00E-01	1.06E-02	9.79E-01
chr2	66667886	66667913	2.02E-09	3.92E-01	2.72E-02	9.35E-01
chr3	197840391	197840554	1.17E-07	4.83E-01	1.61E-02	9.68E-01
chr17	62774587	62774620	1.06E-05	3.95E-01	6.01E-02	8.68E-01
chr11	2847667	2847695	6.03E-04	4.29E-01	3.50E-02	9.25E-01
chr2	19561301	19561322	3.77E-05	4.72E-01	3.90E-02	9.24E-01
chr5	132161708	132161721	2.11E-18	3.78E-01	2.62E-02	9.35E-01
chr3	127469356	127469468	6.90E-04	4.17E-01	7.37E-02	8.50E-01
chr4	24801800	24801810	4.87E-05	3.89E-01	3.08E-02	9.27E-01
chr14	73181172	73181203	1.91E-09	5.13E-01	3.07E-02	9.44E-01
chr6	20024006	20024141	1.45E-09	6.50E-01	5.05E-02	9.28E-01
chr3	170137456	170137470	1.86E-09	3.67E-01	5.69E-03	9.85E-01
chr9	96715356	96715374	7.45E-06	3.39E-01	2.65E-02	9.27E-01
chr9	139090266	139090459	2.60E-04	3.38E-01	1.05E-03	9.97E-01
chr22	44391923	44391940	1.10E-04	4.22E-01	7.15E-02	8.55E-01
chr1	145562732	145562809	8.46E-09	3.95E-01	1.05E-02	9.74E-01
chr19	6753492	6753523	1.12E-08	4.54E-01	3.22E-02	9.34E-01
chr4	8202378	8202403	3.87E-04	4.70E-01	8.39E-02	8.49E-01
chr5	115298830	115298848	2.94E-12	5.62E-01	1.42E-03	9.97E-01
chr2	105459597	105459623	1.08E-04	3.24E-01	6.49E-03	9.80E-01
chr9	139716146	139716160	9.16E-08	6.29E-01	1.08E-01	8.54E-01
chr6	42335118	42335231	4.31E-06	3.39E-01	3.88E-02	8.97E-01
chr3	122296676	122296709	8.84E-06	3.69E-01	1.35E-02	9.65E-01
chr2	131797613	131797636	8.87E-10	3.60E-01	2.87E-02	9.26E-01

-117-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr10	94822392	94822422	6.90E-10	3.57E-01	1.20E-02	9.68E-01
chr13	50701694	50701724	2.33E-11	3.64E-01	6.83E-03	9.82E-01
chr13	100630745	100630775	4.05E-05	3.67E-01	7.75E-03	9.79E-01
chr10	11212536	11212603	5.61E-07	3.63E-01	2.76E-02	9.29E-01
chr9	129377680	129377714	2.85E-10	3.07E-01	5.94E-04	9.98E-01
chr15	74658177	74658190	1.28E-07	4.13E-01	1.06E-03	9.97E-01
chr2	25500296	25500316	4.07E-04	4.30E-01	5.83E-02	8.81E-01
chr10	102986952	102986973	2.98E-04	3.75E-01	2.61E-02	9.35E-01
chr21	44819226	44819239	2.77E-05	6.30E-01	8.77E-02	8.78E-01
chr12	115109842	115109855	4.44E-09	3.63E-01	1.72E-02	9.55E-01
chr22	39148239	39148256	1.93E-05	5.06E-01	3.62E-02	9.33E-01
chr20	43726616	43726631	3.72E-04	4.38E-01	4.98E-02	8.98E-01
chr19	50058236	50058275	2.29E-13	5.75E-01	4.48E-02	9.28E-01
chr15	96903029	96903092	2.06E-05	6.64E-01	5.67E-02	9.21E-01
chr4	94755870	94755918	2.12E-08	3.41E-01	2.31E-02	9.37E-01
chr6	37551119	37551243	1.05E-06	6.51E-01	8.25E-02	8.88E-01
chr7	36679543	36679660	2.83E-04	4.77E-01	9.14E-02	8.39E-01
chr11	123066909	123066942	4.87E-04	3.11E-01	4.72E-03	9.85E-01
chr17	16570355	16570361	1.81E-05	6.09E-01	1.51E-01	8.01E-01
chr6	166976776	166976796	1.34E-05	5.02E-01	8.63E-02	8.53E-01
chr4	16084937	16084953	5.58E-06	3.80E-01	4.31E-02	8.98E-01
chr3	151936218	151936283	2.48E-05	6.28E-01	1.05E-01	8.57E-01
chr12	114878759	114878791	1.42E-08	3.23E-01	1.02E-02	9.69E-01
chr2	220299635	220299643	2.68E-08	3.31E-01	3.76E-03	9.89E-01
chr21	38076854	38076871	1.16E-06	3.28E-01	1.72E-02	9.50E-01
chr2	157186445	157186485	1.12E-05	5.66E-01	1.01E-01	8.48E-01
chr11	65661646	65661754	4.06E-04	4.26E-01	5.66E-02	8.83E-01
chr9	134128101	134128135	4.95E-09	3.78E-01	5.26E-03	9.86E-01
chr3	129694365	129694378	1.06E-07	3.41E-01	1.99E-03	9.94E-01
chr2	233925180	233925317	2.36E-13	3.65E-01	1.38E-02	9.63E-01
chr17	17415395	17415403	1.49E-04	6.37E-01	1.14E-01	8.48E-01
chr9	1045788	1045818	8.55E-08	3.44E-01	1.52E-02	9.58E-01
chr4	16085139	16085180	4.96E-05	3.04E-01	1.03E-03	9.97E-01

39595779.1

-118-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr6	157970252	157970307	4.22E-09	5.22E-01	1.07E-02	9.80E-01
chr21	38081477	38081492	6.68E-09	3.17E-01	8.98E-03	9.72E-01
chr19	16482519	16482633	8.91E-19	3.02E-01	5.29E-04	9.98E-01
chr18	18823158	18823243	1.48E-04	3.70E-01	3.26E-03	9.91E-01
chr7	158362709	158362731	2.24E-08	3.85E-01	1.28E-02	9.68E-01
chr1	178455974	178456004	2.91E-04	3.38E-01	2.46E-02	9.32E-01
chr2	27301795	27301826	6.09E-06	6.96E-01	1.10E-01	8.64E-01
chr13	79170266	79170271	5.04E-07	3.10E-01	2.99E-03	9.90E-01
chr10	121277212	121277385	3.51E-04	4.03E-01	4.23E-02	9.05E-01
chr13	79170302	79170323	8.63E-06	5.10E-01	1.16E-02	9.78E-01
chr16	68676460	68676466	1.56E-05	3.66E-01	2.69E-02	9.32E-01
chr17	1960297	1960323	6.24E-04	5.38E-01	6.90E-02	8.86E-01
chr2	25384521	25384555	3.92E-08	5.00E-01	7.83E-02	8.65E-01
chr19	3179353	3179432	1.23E-21	7.94E-01	2.26E-02	9.72E-01
chr2	43397997	43398041	3.50E-04	3.91E-01	6.99E-02	8.49E-01
chr11	134202273	134202302	2.01E-04	3.33E-01	1.32E-03	9.96E-01
chr21	44819309	44819337	2.25E-09	4.32E-01	5.37E-02	8.89E-01
chr19	47152996	47153006	5.84E-05	3.87E-01	3.17E-02	9.24E-01
chr3	128211282	128211310	2.87E-08	3.73E-01	1.07E-02	9.72E-01
chr11	63974228	63974540	5.83E-09	3.10E-01	9.13E-03	9.71E-01
chr5	115152519	115152539	4.81E-06	5.17E-01	6.13E-04	9.99E-01
chr7	56242832	56242853	7.21E-04	3.94E-01	5.09E-02	8.86E-01
chr5	172370231	172370335	1.56E-05	5.50E-01	9.01E-02	8.59E-01
chr4	102711857	102711879	4.93E-07	3.23E-01	1.62E-02	9.52E-01
chr1	203032620	203032687	1.94E-05	6.52E-01	6.82E-02	9.05E-01
chr5	176023813	176023852	3.71E-07	6.67E-01	2.31E-02	9.67E-01
chr4	158143143	158143233	4.15E-04	3.33E-01	1.51E-03	9.95E-01
chr2	119611407	119611447	1.42E-06	3.13E-01	8.30E-03	9.74E-01
chr3	192125879	192125900	2.69E-04	4.18E-01	2.92E-02	9.35E-01
chr18	9404557	9404633	1.79E-09	3.57E-01	3.34E-02	9.15E-01
chr1	47899900	47899973	3.19E-07	4.06E-01	5.61E-02	8.79E-01
chr14	105527019	105527137	8.63E-10	3.28E-01	1.18E-02	9.65E-01
chr3	50312991	50313021	2.27E-21	5.56E-01	7.01E-03	9.88E-01

39559779.1

-119-

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr22	47070183	47070285	9.36E-16	5.14E-01	2.54E-02	9.53E-01
chr14	27065959	27066053	6.38E-04	4.31E-01	5.19E-02	8.92E-01
chr7	27197712	27197726	3.37E-05	4.24E-01	1.53E-02	9.65E-01
chr16	2317491	2317523	6.22E-14	5.26E-01	3.55E-02	9.37E-01
chr7	557581	557595	1.30E-10	7.50E-01	1.08E-02	9.86E-01
chr9	134148582	134148640	3.73E-14	5.69E-01	9.44E-03	9.84E-01
chr14	38066924	38066944	9.06E-09	3.65E-01	1.99E-02	9.48E-01
chr4	1210659	1210715	5.63E-18	4.84E-01	1.36E-02	9.73E-01
chr4	1210785	1210901	5.66E-05	5.03E-01	2.80E-02	9.47E-01
chr3	181438113	181438133	2.31E-10	9.35E-01	7.17E-02	9.29E-01
chr2	63284055	63284066	2.24E-05	4.13E-01	8.01E-03	9.81E-01
chr10	13390041	13390139	8.83E-04	5.51E-01	7.01E-02	8.87E-01
chr18	11947874	11947897	9.55E-07	4.40E-01	2.68E-02	9.43E-01
chr9	6645173	6645184	6.65E-04	4.67E-01	6.59E-02	8.76E-01
chr8	142414178	142414227	2.36E-05	6.39E-01	4.76E-02	9.31E-01
chr13	28498832	28498881	5.00E-08	3.99E-01	7.14E-03	9.82E-01
chr11	63774985	63775047	3.55E-05	3.14E-01	5.80E-03	9.82E-01
chr7	87228807	87228922	2.54E-07	5.27E-01	1.88E-02	9.66E-01
chr11	128609074	128609245	9.64E-05	3.72E-01	3.64E-02	9.11E-01
chr22	47070690	47070785	3.33E-04	6.89E-01	9.78E-02	8.76E-01
chr13	96294381	96294499	2.52E-09	3.29E-01	5.24E-03	9.84E-01
chr5	179507584	179507627	1.78E-04	5.75E-01	8.02E-02	8.78E-01
chr8	28637012	28637063	5.18E-07	7.88E-01	1.63E-01	8.28E-01
chr10	79380018	79380032	1.83E-07	1.00E+00	1.47E-01	8.72E-01
chr12	51718198	51718251	2.70E-07	7.57E-01	1.05E-01	8.78E-01
chr7	155259609	155259665	2.94E-08	3.81E-01	1.55E-03	9.96E-01
chr12	52652502	52652523	5.67E-12	4.31E-01	7.01E-03	9.84E-01
chr12	125053451	125053464	4.82E-05	3.43E-01	5.20E-03	9.85E-01
chr2	97428602	97428718	5.20E-07	4.57E-01	4.71E-02	9.07E-01
chr17	17295442	17295619	1.03E-16	5.55E-01	1.59E-02	9.72E-01
chr16	28996532	28996602	2.14E-04	3.75E-01	3.23E-02	9.21E-01
chr16	2317602	2317631	4.83E-06	3.77E-01	4.62E-02	8.91E-01
chr7	95402297	95402335	1.36E-06	4.19E-01	7.12E-04	9.98E-01

39595779.1

-120-

Chrom	Start	End	q-value	Cancer-exp	Normal-exp	Weight
chr17	15244821	15244853	1.52E-04	3.20E-01	1.14E-02	9.66E-01
chr3	192126077	192126114	2.45E-08	3.16E-01	4.81E-03	9.85E-01
chr2	43398069	43398085	1.25E-09	3.89E-01	4.46E-02	8.97E-01
chr12	124877816	124877872	7.30E-08	1.00E+00	1.23E-01	8.90E-01
chr1	208132616	208132726	3.42E-08	3.04E-01	1.18E-03	9.96E-01
chr17	61512658	61512683	1.25E-10	3.96E-01	1.75E-02	9.58E-01
chr19	42071605	42071726	7.47E-04	3.70E-01	3.07E-02	9.23E-01
chr10	101292552	101292646	8.64E-05	3.61E-01	1.65E-03	9.95E-01
chr5	137795191	137795306	9.26E-06	5.07E-01	4.49E-02	9.19E-01
chr22	50319158	50319194	5.63E-04	6.23E-01	1.15E-01	8.44E-01
chr12	115103722	115103767	3.73E-17	4.94E-01	1.94E-02	9.62E-01
chr12	123131513	123131635	2.56E-04	5.00E-01	3.05E-02	9.42E-01
chr13	114875997	114876008	1.30E-05	3.91E-01	4.49E-02	8.97E-01
chr17	79134237	79134248	2.42E-06	4.99E-01	3.71E-02	9.31E-01
chr3	129295780	129295896	6.74E-09	6.11E-01	4.42E-02	9.33E-01
chr10	102881084	102881114	2.40E-08	3.06E-01	5.53E-03	9.82E-01
chr15	96890743	96890812	7.95E-14	6.41E-01	1.92E-02	9.71E-01
chr2	219847214	219847286	3.27E-06	5.00E-01	0.00E+00	1.00E+00
chr10	6213930	6214079	5.67E-05	5.83E-01	5.20E-02	9.18E-01
chr11	45114836	45114845	7.51E-11	7.75E-01	6.09E-02	9.27E-01
chr4	81952496	81952527	1.08E-05	4.33E-01	5.67E-02	8.84E-01
chr4	159091800	159091863	1.51E-04	4.09E-01	6.40E-02	8.65E-01
chr19	16197202	16197336	2.39E-04	4.71E-01	6.04E-02	8.86E-01
chr10	122708551	122708626	4.00E-06	3.59E-01	1.37E-03	9.96E-01
chr7	633049	633101	5.38E-07	3.63E-01	2.77E-02	9.29E-01
chr14	105531892	105532022	4.98E-04	5.08E-01	4.31E-02	9.22E-01
chr3	50313747	50313762	1.35E-22	5.09E-01	9.92E-03	9.81E-01
chr6	42072513	42072527	4.31E-04	3.39E-01	1.58E-03	9.95E-01
chr11	88242385	88242459	4.29E-08	4.12E-01	6.93E-03	9.83E-01
chr1	32237851	32237877	3.37E-04	3.33E-01	9.60E-04	9.97E-01
chr15	41795065	41795089	1.38E-06	3.63E-01	2.37E-02	9.39E-01
chr18	46298176	46298192	3.78E-05	8.79E-01	2.60E-01	7.72E-01
chr5	1495316	1495431	2.09E-12	6.25E-01	4.74E-02	9.29E-01

39595779.1

-121-

Chrom	Start	End	q-value	Cancer.exp	Normal.exp	Weight
chr3	138679355	138679366	7.38E-07	3.57E-01	5.73E-03	9.84E-01
chr5	134880200	134880327	1.16E-06	4.87E-01	2.63E-03	9.95E-01
chr2	66666698	66666809	2.72E-06	6.67E-01	1.20E-02	9.82E-01
chr6	42072328	42072373	3.87E-10	5.14E-01	7.48E-03	9.86E-01
chr9	101822013	101822096	2.70E-09	7.30E-01	6.09E-02	9.23E-01
chr6	19691985	19692003	3.43E-10	3.18E-01	5.76E-03	9.82E-01
chr5	177989059	177989077	4.41E-10	5.37E-01	3.12E-02	9.45E-01
chr2	177014492	177014518	6.48E-08	3.57E-01	2.05E-02	9.46E-01
chr19	55553070	55553087	3.65E-14	4.87E-01	3.45E-02	9.34E-01
chr17	40822427	40822434	2.19E-05	6.55E-01	1.03E-01	8.64E-01
chr19	3178517	3178545	1.09E-16	3.70E-01	1.59E-02	9.59E-01
chr10	124910708	124910730	9.40E-06	3.87E-01	3.09E-03	9.92E-01
chr13	29106180	29106363	9.39E-07	6.50E-01	3.71E-02	9.46E-01
chr9	126135891	126135950	3.10E-05	4.01E-01	1.60E-03	9.96E-01
chr17	56409588	56409742	2.12E-07	6.50E-01	5.48E-02	9.22E-01
chr6	50787949	50787984	3.26E-05	4.40E-01	3.52E-02	9.26E-01
chr16	2286867	2286879	3.18E-06	3.23E-01	2.11E-02	9.39E-01
chr7	129425931	129426013	2.12E-07	4.15E-01	3.15E-02	9.30E-01
chr16	68390251	68390361	3.65E-19	4.71E-01	2.22E-02	9.55E-01
chr22	19710841	19710855	2.14E-05	4.05E-01	4.28E-03	9.90E-01
chr1	1148346	1148354	1.67E-04	4.74E-01	5.68E-02	8.93E-01
chr16	88496790	88496819	1.25E-05	5.00E-01	3.41E-02	9.36E-01
chr8	134094859	134094975	2.26E-05	5.05E-01	1.37E-01	7.86E-01
chr5	138850540	138850606	5.90E-20	4.59E-01	1.74E-02	9.63E-01
chr12	6882964	6882984	2.49E-07	4.76E-01	6.43E-02	8.81E-01
chr7	155246915	155246940	5.33E-06	4.07E-01	2.30E-02	9.47E-01
chr10	12223668	12223715	2.97E-07	4.35E-01	7.50E-02	8.53E-01
chr22	20783780	20783802	2.62E-09	3.73E-01	8.99E-03	9.76E-01

Table 5a:
The sets of cancer specific markers derived from MARS based features selection on training data sets

Chrom	Start	End
chr10	103454457	103454477
chr10	104352036	104352196
chr10	10847503	10847691
chr10	27702716	27702729
chr10	30110010	30110051
chr11	106698602	106698622
chr11	121395074	121395118
chr11	2020027	2020065
chr11	33929206	33929320
chr11	36145336	36145373
chr11	47416330	47416377
chr1	156087833	156087873
chr1	16159932	16159950
chr11	62139527	62139609
chr11	62370212	62370238
chr1	201509291	201509346
chr1	207921424	207921528
chr12	110033418	110033548
chr12	25537415	25537441
chr1	245970850	245970891
chr12	52438106	52438159
chr12	52889045	52889173
chr12	54332947	54332963
chr12	54473544	54473561
chr12	6901902	6901997
chr12	96564033	96564157
chr13	19918783	19918794
chr13	20370973	20371021
chr13	34253491	34253585
chr14	64130579	64130727
chr14	91818657	91818672
chr14	98190804	98190884
chr15	58734402	58734645
chr15	86298535	86298604
chr16	69961433	69961449
chr17	17295442	17295619
chr17	21415130	21415142
chr17	26173343	26173438
chr17	37211873	37211969
chr17	64948116	64948158
chr17	74526150	74526161
chr17	75143193	75143219
chr17	76875977	76875998
chr18	72916692	72916705
chr19	1169117	1169138

Chrom	Start	End
chr19	16178526	16178600
chr1	92946665	92946768
chr19	3403636	3403711
chr19	3403810	3403840
chr19	47496561	47496586
chr19	50931391	50931435
chr19	52880780	52880980
chr19	58220626	58220669
chr20	30816216	30816330
chr20	34391305	34391406
chr20	41153703	41153752
chr2	101926462	101926479
chr2	111141088	111141148
chr2	11884471	11884591
chr21	35197108	35197284
chr2	197040996	197041007
chr2	202945219	202945279
chr2	232348683	232348704
chr2	232745254	232745369
chr22	37942946	37942985
chr22	38104696	38104738
chr22	39148367	39148393
chr2	242195200	242195217
chr22	42815670	42815694
chr22	45692585	45692725
chr22	46022276	46022300
chr2	25138739	25138897
chr2	36717718	36717897
chr2	46219462	46219480
chr2	48776047	48776154
chr2	66673054	66673077
chr2	85637199	85637376
chr2	99155846	99155916
chr3	126645737	126645771
chr3	141098240	141098373
chr3	18284952	18285049
chr3	184056458	184056470
chr3	52250610	52250789
chr4	1195669	1195719
chr4	177420126	177420335
chr4	186808246	186808296
chr4	6697966	6698085
chr4	8727069	8727091
chr5	138730608	138730648
chr5	139725539	139725550
chr5	156570589	156570747
chr5	176831296	176831309
chr6	107012070	107012253
chr6	108984666	108984700

Chrom	Start	End
chr6	133523044	133523327
chr6	149805994	149806024
chr6	158411773	158411853
chr6	159128335	159128361
chr6	163767808	163767820
chr6	170585982	170586062
chr6	170589857	170589872
chr6	21246436	21246559
chr6	37533092	37533142
chr6	43192429	43192519
chr6	43650729	43650959
chr6	43894465	43894595
chr6	43939175	43939212
chr6	71874735	71874758
chr7	100549051	100549245
chr7	100875668	100875696
chr7	101884716	101884822
chr7	140096553	140096723
chr7	157071940	157071954
chr7	17139387	17139535
chr7	33176934	33177005
chr7	36320659	36320695
chr7	39393611	39393625
chr7	41428132	41428227
chr7	56183788	56183856
chr7	601463	601503
chr7	614586	614611
chr7	64019659	64019780
chr7	75018128	75018144
chr8	125766626	125766669
chr8	131076680	131076724
chr9	116327621	116327718
chr9	130517759	130517991
chr9	130955063	130955167
chr9	95964149	95964200
chr10	81004561	81004609
chr1	162442502	162442619
chr1	212415595	212415691
chr1	235062685	235062859
chr15	77286562	77286576
chr16	17499260	17499310
chr16	78636379	78636435
chr17	16873955	16874009
chr17	41669725	41669744
chr17	45786237	45786298
chr1	8400056	8400143
chr19	49839124	49839203
chr19	7723032	7723142
chr20	57425892	57425937

-125-

Chrom	Start	End
chr20	62200085	62200109
chr2	174890243	174890269
chr2	5836838	5836852
chr3	196347229	196347361
chr4	68411077	68411095
chr6	155434075	155434119
chr6	2999401	2999433
chr7	50850278	50850660
chr14	28485394	28485506
chr19	55766296	55766320
chr22	46480824	46480971

Table 5b:

The sets of tissue specific markers derived from MARS based features selection on training data sets

5

Chrom	Start	End
chr10	12526101	12526198
chr10	134062557	134062567
chr10	2992831	2992986
chr10	71250913	71251031
chr10	894133	894176
chr11	10529517	10529548
chr11	10580059	10580207
chr11	117454728	117454907
chr11	125840213	125840218
chr11	126226966	126226994
chr11	20626148	20626166
chr1	120906033	120906056
chr11	2828703	2828805
chr1	154978354	154978371
chr11	61467518	61467601
chr1	164545706	164545781
chr11	66871960	66871973
chr1	168687549	168687721
chr1	172291575	172291727
chr11	85423851	85424008
chr11	91958190	91958214
chr1	19249130	19249147
chr11	94600637	94600757
chr1	20821827	20821848
chr12	121164886	121164929
chr12	123352751	123352815
chr12	124905745	124905758
chr1	219786354	219786374
chr12	214138	214168
chr1	24468442	24468457
chr1	25257913	25257952

-126-

Chrom	Start	End
chr12	58131153	58131181
chr12	63157011	63157117
chr12	6336078	6336124
chr12	63828393	63828590
chr12	65813665	65813787
chr12	69440660	69440851
chr12	92793926	92794091
chr13	24575928	24576139
chr13	75981727	75981971
chr1	37941197	37941228
chr1	4045074	4045329
chr1	41898133	41898168
chr14	24078834	24078895
chr14	35319476	35319527
chr14	55170355	55170479
chr14	62403215	62403374
chr14	64130579	64130727
chr14	65708245	65708394
chr14	73279457	73279546
chr14	75683126	75683335
chr14	77089391	77089627
chr14	78108418	78108650
chr15	23892535	23892608
chr15	31850650	31850711
chr15	61497667	61497716
chr15	66998995	66999028
chr15	67179342	67179441
chr15	67356599	67356824
chr15	68136062	68136196
chr15	70488068	70488255
chr15	76268973	76269064
chr15	90727362	90727449
chr16	11749588	11749704
chr16	12510137	12510244
chr16	1544949	1545020
chr16	16244310	16244423
chr16	67687178	67687277
chr16	69516112	69516139
chr16	71560275	71560421
chr16	75279686	75279773
chr16	84967640	84967652
chr16	86965384	86965400
chr16	88096279	88096362
chr17	16924593	16924617
chr17	17473311	17473523
chr17	25821621	25821649
chr17	29469666	29469741
chr17	43483074	43483225
chr17	43506900	43506997

Chrom	Start	End
chr17	71538984	71539080
chr17	79456619	79456772
chr17	950579	950597
chr18	13447880	13447894
chr18	14458513	14458531
chr18	45058265	45058273
chr18	5489107	5489117
chr18	74499570	74499587
chr18	74665755	74665892
chr19	14673070	14673115
chr19	18131156	18131182
chr19	2291605	2291625
chr19	33879823	33879853
chr19	35505322	35505399
chr19	35630532	35630541
chr19	53088342	53088389
chr1	9563349	9563412
chr19	58629652	58629666
chr20	17731262	17731276
chr20	32336264	32336332
chr20	3732686	3732700
chr20	3732852	3732883
chr20	47421687	47421875
chr20	50352688	50352798
chr20	62200085	62200109
chr20	62959294	62959302
chr20	8880040	8880231
chr2	106048152	106048370
chr2	109788225	109788402
chr2	110798982	110799019
chr2	14375315	14375355
chr21	44781718	44781727
chr2	175739843	175739898
chr22	20237340	20237395
chr2	231277281	231277433
chr22	32475876	32475936
chr22	34095314	34095472
chr2	236611862	236611872
chr22	38484978	38484992
chr22	45887405	45887583
chr22	46476111	46476156
chr22	46786620	46786697
chr22	50343067	50343097
chr22	50985445	50985457
chr2	39472294	39472340
chr2	86248350	86248380
chr2	97305261	97305286
chr3	106937697	106937906
chr3	113160421	113160434

Chrom	Start	End
chr3	122459452	122459655
chr3	128712749	128712769
chr3	130098135	130098295
chr3	170965322	170965449
chr3	195356987	195357011
chr3	29320866	29320988
chr3	38012448	38012522
chr3	51965036	51965246
chr4	1244005	1244024
chr4	129556483	129556502
chr4	151504097	151504165
chr4	169561292	169561324
chr4	185722597	185722644
chr4	2257619	2257679
chr4	26755387	26755468
chr4	965678	965736
chr5	10624799	10624898
chr5	164830490	164830637
chr5	170041185	170041290
chr5	5140660	5140708
chr5	95198558	95198771
chr6	115682789	115682877
chr6	135642658	135642719
chr6	139858289	139858355
chr6	150623168	150623365
chr6	155434075	155434119
chr6	169351059	169351075
chr6	170574715	170574763
chr6	21246436	21246559
chr6	24981775	24981874
chr6	3136197	3136313
chr6	42213769	42213962
chr6	47427927	47427990
chr7	100875668	100875696
chr7	102963572	102963623
chr7	151290079	151290342
chr7	1577891	1577924
chr7	47419767	47419848
chr8	126645983	126646069
chr8	129040440	129040473
chr8	129103440	129103463
chr8	141523764	141523791
chr8	142189542	142189581
chr8	143408017	143408047
chr8	144361314	144361344
chr8	144361364	144361393
chr8	145654572	145654583
chr8	27530948	27530986
chr8	41979371	41979431

Chrom	Start	End
chr8	53070170	53070416
chr8	8230665	8230840
chr9	116356861	116357046
chr9	130517759	130517991
chr9	130860588	130860681
chr9	130877668	130877724
chr9	134455414	134455516
chr9	136921088	136921169
chr9	138904160	138904187
chr9	140188472	140188494
chr9	6788700	6788803
chr9	86064195	86064250
chr10	10840659	10840699
chr10	135174002	135174025
chr10	6254206	6254244
chr10	75647560	75647648
chr10	77109265	77109419
chr10	80919230	80919249
chr1	151870591	151870681
chr1	156261302	156261327
chr11	64740009	64740016
chr1	205449277	205449449
chr12	123402308	123402427
chr1	2164204	2164384
chr14	106774181	106774432
chr14	65681455	65681535
chr1	47039258	47039357
chr14	91818657	91818672
chr15	79082104	79082154
chr1	6515867	6515873
chr16	78636379	78636435
chr17	17415395	17415403
chr17	19429425	19429525
chr17	8371613	8371640
chr1	8433676	8433790
chr19	1169117	1169138
chr19	13823667	13823731
chr19	22610862	22610883
chr20	50073805	50073870
chr20	57417351	57417392
chr20	60760760	60760972
chr21	30612916	30613022
chr21	34399242	34399284
chr21	45705693	45705715
chr21	47581404	47581439
chr2	158694393	158694462
chr2	161348608	161348781
chr2	219124246	219124271
chr2	219147468	219147621

Chrom	Start	End
chr3	123118882	123119159
chr3	13828799	13828855
chr3	48540230	48540315
chr4	154460014	154460117
chr4	19929150	19929272
chr4	54889382	54889454
chr5	112538998	112539022
chr5	1725410	1725419
chr5	17366324	17366430
chr5	66564412	66564447
chr7	158890101	158890118
chr7	79875750	79875982
chr8	22447199	22447304
chr8	81279608	81279654
chr9	102213820	102213912
chr9	95964655	95964697
chr10	104840711	104840824
chr10	3917268	3917367
chr11	120418670	120418806
chr11	57199801	57200049
chr12	126978451	126978675
chr1	28431421	28431512
chr13	107869633	107869814
chr14	34504893	34504986
chr16	23303806	23303916
chr17	38708456	38708524
chr19	1221753	1221771
chr19	48983562	48983583
chr19	57412139	57412201
chr22	20237221	20237247
chr22	29866375	29866425
chr2	237581359	237581546
chr22	38148519	38148762
chr3	58103818	58103844
chr3	65940157	65940167
chr3	88312669	88312786
chr4	113634382	113634491
chr4	187071070	187071120
chr4	2904348	2904411
chr4	68411077	68411095
chr5	149994824	149994881
chr6	57607668	57607879
chr8	10559717	10559846
chr9	132658755	132658787
chr1	2077786	2077884
chr13	21900491	21900523
chr15	79342608	79342818
chr17	46697530	46697565
chr18	60050875	60051139

-131-

Chrom	Start	End
chr2	187482101	187482155
chr3	194014878	194014894
chr9	129511230	129511375
chr9	78871348	78871392
chr10	121578130	121578169
chr10	88729962	88730008
chr15	23115150	23115163
chr22	37942946	37942985
chr3	126645737	126645771
chr5	170041887	170042039
chr6	31869020	31869031
chr11	128558360	128558418
chr14	93565785	93565854
chr10	30327101	30327236
chr10	6183283	6183437
chr10	30110010	30110051
chr17	45786237	45786298
chr3	13063400	13063494
chr2	105990523	105990613

[0088] The disclosures of all references listed herein, including the following references, are incorporated by reference herein in their entireties:

REFERENCES:

1. EP Patent No. 2,258,871
- 5 2. EP Patent No. 2,643,479
3. US Patent No. 8,026,067
4. US Patent No. 8,110,361
5. US Patent No. 8,288,100
6. US Patent No. 8,293,470
- 10 7. US Patent No. 8,586,313
8. US Patent No. 9,096,906
9. WO 2015/116837
10. Accession GSE79279, Gene Expression Omnibus (GEO)
11. The ENCODE Project Consortium. An integrated encyclopedia of DNA
- 15 elements in the human genome. *Nature* 489:57-74 (2012).
12. Bernstein et al. The NIH Roadmap Epigenomics Mapping Consortium. *Nat Biotechnol* 28:1045-1048 (2010).
13. Blattler et al. Global loss of DNA methylation uncovers intronic enhancers in genes showing expression changes. *Genome Biol* 15:469 (2014).
- 20 14. Chen et al. Loss of 5-hydroxymethylcytosine is linked to gene body hypermethylation in kidney cancer. *Cell Res* 26:103-118 (2016).
15. Diep et al. Library-free methylation sequencing with bisulfite padlock probes. *Nat Methods*. 9(3):270-272 (2012).
16. Dixon et al. Topological domains in mammalian genomes identified by analysis
- 25 of chromatin interactions. *Nature* 485:376-380 (2012).
17. Friedman Multivariate Adaptive Regression Splines. *Ann. Of Stat.* 19(1): 1-67 (1991).
18. Gong et al. DeconRNASeq: a statistical framework for deconvolution of heterogeneous tissue samples based on mRNA-Seq data. *Bioinformatics* 29:1083-1085 (2013).
- 30 19. Guelen et al. Domain organization of human chromosomes revealed by mapping of nuclear lamina interactions. *Nature* 453:948-951 (2008).

20. Guo et al. Single-cell methylome landscapes of mouse embryonic stem cells and early embryos analyzed using reduced representation bisulfite sequencing. *Genome Res* 23:2126-2135 (2013).
21. Hansen et al. Increased methylation variation in epigenetic domains across
5 cancer types. *Nat Genet* 43:768-775 (2011).
22. Heyn et al. Distinct DNA methylomes of newborns and centenarians. *Proc Natl Acad Sci USA* 109:10522-10527 (2012).
23. Heyn et al. Epigenomic analysis detects aberrant super-enhancer DNA
methylation in human cancer. *Genome Biol* 17:11 (2016).
- 10 24. Houseman et al. DNA methylation arrays as surrogate measures of cell mixture
distribution. *BMC Bioinformatics* 13:86 (2012).
25. Houseman et al. Reference-free deconvolution of DNA methylation data and
mediation by cell composition effects. *BMC Bioinformatics* 17:259 (2016).
26. Irizarry et al. The human colon cancer methylome shows similar hypo- and
15 hypermethylation at conserved tissue-specific CpG island shores. *Nat Genet* 41:178-186
(2009).
27. Johansson et al. Targeted resequencing of candidate genes using selector
probes. *Nucleic Acids Res.* 39(2):e8 (2011).
28. Johnson et al. Adjusting batch effects in microarray expression data using
20 empirical Bayes methods. *Biostatistics* 8:118-127 (2007).
29. Jones, P.A. & Martienssen, R. A blueprint for a Human Epigenome Project: the
AACR Human Epigenome Workshop. *Cancer Res* 65:11241-11246 (2005).
30. Komori et al. Application of microdroplet PCR for large-scale targeted bisulfite
sequencing. *Genome Res.* 21(10):1738-1745 (2011).
- 25 31. Landan et al. Epigenetic polymorphism and the stochastic formation of
differentially methylated regions in normal and cancerous tissues. *Nat Genet.* 44(11):1207-
1214 (2012).
32. Landau et al. Locally disordered methylation forms the basis of intratumor
methylome variation in chronic lymphocytic leukemia. *Cancer Cell* 26:813-825 (2014).
- 30 33. Lehmann-Werman et al. Identification of tissue-specific cell death using
methylation patterns of circulating DNA. *Proc Natl Acad Sci USA* 113, E1826-1834 (2016).

34. Leung et al. Integrative analysis of haplotype-resolved epigenomes across human tissues. *Nature* 518:350-354 (2015).
35. Lohk et al. DNA methylome profiling of human tissues identifies global and tissue-specific methylation patterns. *Genome Biology* 15(4):r54 (2014).
- 5 36. Meissner et al. Reduced representation bisulfite sequencing for comparative high-resolution DNA methylation analysis. *Nucleic Acids Res.* 33(18):5868-5877 (2005).
37. Mitsui et al. The homeoprotein Nanog is required for maintenance of pluripotency in mouse epiblast and ES cells. *Cell* 113:631-642 (2003).
38. Papageorgiou et al. Fetal-specific DNA methylation ratio permits noninvasive prenatal diagnosis of trisomy 21. *Nat Med.* 17(4):510-513 (2011).
- 10 39. Pujadas et al. Regulated noise in the epigenetic landscape of development and disease. *Cell* 148:1123-1131 (2012).
40. Saito et al. Linkage disequilibrium analysis of allelic heterogeneity in DNA methylation. *Epigenetics* 10:1093-1098 (2015).
- 15 41. Schultz et al. Human body epigenome maps reveal noncanonical DNA methylation variation. *Nature* 523:212-216 (2015).
42. Shao et al. Deciphering the heterogeneity in DNA methylation patterns during stem cell differentiation and reprogramming. *BMC Genomics* 15:978 (2014).
43. Shoemaker et al. Allele-specific methylation is prevalent and is contributed by CpG-SNPs in the human genome. *Genome Res* 20:883-889 (2010).
- 20 44. Shu et al. Induction of pluripotency in mouse somatic cells with lineage specifiers. *Cell* 153:963-975 (2013).
45. Slatkin Linkage disequilibrium -- understanding the evolutionary past and mapping the medical future. *Nat Rev Genet* 9:477-485 (2008).
- 25 46. Snyder et al. Cell-free DNA Comprises an In Vivo Nucleosome Footprint that Informs Its Tissues-Of-Origin. *Cell* 164:57-68 (2016).
47. Sun et al. Plasma DNA tissue mapping by genome-wide methylation sequencing for noninvasive prenatal, cancer, and transplantation assessments. *Proc Natl Acad Sci USA* 112, E5503-5512 (2015).
- 30 48. Takai et al. Comprehensive analysis of CpG islands in human chromosomes 21 and 22. *Proc Natl Acad Sci USA* 99:3740-3745 (2002).

-135-

49. Wen, B., Wu, H., Shinkai, Y., Izarray, R.A. & Feinberg, A.P. Large histone H3 lysine 9 dimethylated chromatin blocks distinguish differentiated from embryonic stem cells. *Nat Genet* 41:246-250 (2009).
50. Wigler et al. The somatic replication of DNA methylation. *Cell* 24:33-40
5 (1981).
51. Williams et al. TET1 and hydroxymethylcytosine in transcription and DNA methylation fidelity. *Nature* 473:343-348 (2011).
52. Xie, H. et al. Genome-wide quantitative assessment of variation in DNA methylation patterns. *Nucleic Acids Res.* 39(10):4099-4108 (2011).
- 10 53. Xie, W. et al. Epigenomic analysis of multilineage differentiation of human embryonic stem cells. *Cell* 153:1134-1148 (2013).
54. Ziller et al. Charting a dynamic DNA methylation landscape of the human genome. *Nature* 500:477-481 (2013).