



162162

ABSTRACT



G00002

5 **METHODS AND APPARATUS FOR PERFORMING ARRAY** **MICROCRYSTALLIZATIONS**

An apparatus for setting up crystallisation experiments comprises a mother
liquor drop station that removes mother liquor from a plurality of wells of a
10 multi well plate and delivers submicroliter volumes of mother liquor to a
plurality of different drop regions at the same time. The apparatus also
comprises a molecule drop station that delivers submicroliter volumes of a
solution containing a molecule to be crystallized to the drop regions.

Fig 6A

15