

INDIRECT HOMOGENEOUS MOBILITY SHIFT ASSAYS FOR THE DETECTION OF BIOLOGICS IN PATIENT SAMPLES

CROSS-REFERENCE TO RELATED APPLICATIONS

[0001] This application claims priority to U.S. Provisional Application No. 62/088,465, filed December 5, 2014, U.S. Provisional Application No. 62/113,317, filed February 6, 2015, and U.S. Provisional Application No. 62/158,791, filed May 8, 2015, the disclosures of which are hereby incorporated by reference in their entirety for all purposes.

BACKGROUND OF THE INVENTION

[0002] Inflammatory bowel disease (IBD), which occurs world-wide and afflicts millions of people, is the collective term used to describe three gastrointestinal disorders of unknown etiology: Crohn's disease (CD), ulcerative colitis (UC), and indeterminate colitis (IC). IBD, together with irritable bowel syndrome (IBS), will affect one-half of all Americans during their lifetime, at a cost of greater than \$2.6 billion dollars for IBD and greater than \$8 billion dollars for IBS. A primary determinant of these high medical costs is the difficulty of diagnosing digestive diseases and how these diseases will progress. The cost of IBD and IBS is compounded by lost productivity, with people suffering from these disorders missing at least 8 more days of work annually than the national average.

[0003] Despite the successes of anti-TNF α therapies in the treatment of IBD, a subpopulation of patients are refractory to treatment, highlighting an unmet medical need for new therapies. Vedolizumab is a gut-specific, $\alpha 4\beta 7$ integrin-neutralizing monoclonal antibody, which does not affect peripheral blood cell counts and appears to lack systemic effects. Vedolizumab is a new anti-inflammatory treatment option for the management of therapy-refractory patients. In addition, ustekinumab is a IL12p40 monoclonal antibody, which is another novel IBD therapeutic. However, the availability of diagnostic tests to accurately measure the levels of biologics such as vedolizumab and ustekinumab is necessary for the effective use of these novel therapeutics in IBD patients.

[0004] As such, there is a need in the art for assays to detect the presence or level of biologics such as vedolizumab and ustekinumab in a patient sample to monitor drug therapy and to guide treatment decisions. Such assays are particularly useful for the therapeutic

management of diseases such as ulcerative colitis and Crohn's disease using an individualized approach to monitor drug efficacy and optimize therapy accordingly, and can include assessing disease course and clinical parameters such as pharmacodynamics, disease activity indices, disease burden, and inflammatory biomarkers. The present invention satisfies this need and provides related advantages as well.

BRIEF SUMMARY OF THE INVENTION

[0005] The present invention provides novel indirect homogeneous mobility shift assays for detecting and measuring the presence or level of a biologic in a sample. The assays of the present invention are particularly advantageous for detecting and measuring the presence or level of biologics that target complex antigens including cell surface proteins, transmembrane proteins, heavily glycosylated proteins, multimeric proteins, and the like. As such, the present invention provides information for guiding treatment decisions for those subjects receiving therapy with a biologic agent and improves the accuracy of optimizing therapy, reducing toxicity, and/or monitoring the efficacy of therapeutic treatment to biologic therapy.

The present invention also provides isolated soluble $\alpha 4\beta 7$ integrin heterodimers and isolated soluble IL-12p40 monomers that are suitable for use in the assays described herein.

[0006] In certain aspects, the present invention provides a method for determining the presence or level of a biologic in a sample, the method comprising:

(a) contacting the sample with an unlabeled soluble antigen that binds to the biologic to form an unlabeled complex between the antigen and the biologic in the sample;

(b) contacting the sample from step (a) with a labeled form of the biologic to form a labeled complex between the antigen and the labeled biologic;

(c) subjecting the unlabeled and labeled complexes to size exclusion chromatography to separate the unlabeled and labeled complexes from free labeled biologic and to detect an amount of the free labeled biologic; and

(d) comparing the amount of the free labeled biologic detected in step (c) to a standard curve of known amounts of the biologic, thereby determining the presence or level of the biologic in the sample.

[0007] In some embodiments, the biologic includes antibodies, antibody fragments, proteins, polypeptides, peptides, fusion proteins, multivalent binding proteins, antibody-drug conjugates, vaccines, nucleic acids, sugars, recombinant forms thereof, engineered forms thereof, and combinations thereof.

[0008] In certain embodiments, the antigen is a soluble version (*e.g.*, a soluble fragment, variant, or monomeric form) of a membrane-bound protein, a glycosylated protein, a multimeric protein, an insoluble protein, a protein that is difficult to express or purify, and/or a large protein. In certain instances, the antigen is a soluble extracellular domain of a membrane-bound protein (*e.g.*, a soluble cytokine receptor extracellular domain). In certain other instances, the antigen is a soluble homodimer or heterodimer comprising the extracellular domains of two membrane-bound proteins (*e.g.*, a soluble integrin heterodimer). In yet other instances, the antigen is a soluble protein that does not multimerize and remains in monomeric form once isolated and/or purified (*e.g.*, a soluble cytokine variant with one or more cysteine residues mutated to minimize or eliminate the formation of multimers).

[0009] In other embodiments, the sample is a whole blood, serum, or plasma sample, *e.g.*, from a subject receiving biologic therapy. In preferred embodiments, the sample is serum. In particular embodiments, the subject has a disease or disorder such as, *e.g.*, an autoimmune disease (*e.g.*, rheumatoid arthritis), an inflammatory disease (*e.g.*, inflammatory bowel disease (IBD) such as Crohn's disease (CD) or ulcerative colitis (UC)), or cancer.

[0010] In particular embodiments, the standard curve is generated by incubating the antigen and the labeled biologic with a (*e.g.*, two-fold) serial dilution of known amounts of the biologic. In certain embodiments, the area under the curve (AUC) of the free labeled biologic is plotted against (*e.g.*, the logarithm of) known amounts of the biologic obtained from the standard curve, and the level of the biologic in the sample is calculated by interpolation, *e.g.*, based upon the size of the peak area of the free labeled biologic. In other embodiments, free label added to a stock solution of labeled biologic is used as a labeled biologic loading control. The ratio of the free labeled biologic to free label is plotted against known amounts of biologic.

[0011] In one particular embodiment, the presence and/or level of an anti- $\alpha 4\beta 7$ integrin drug (*e.g.*, vedolizumab) is determined with an indirect homogeneous mobility shift assay using size exclusion chromatography as described herein.

[0012] In another particular embodiment, the presence and/or level of an anti-IL12p40 drug (*e.g.*, ustekinumab) is determined with an indirect homogeneous mobility shift assay using size exclusion chromatography as described herein.

[0013] In other embodiments, the presence and/or level of anti-drug antibodies (ADA) (*e.g.*, autoantibodies including HACA, HAHA, *etc.*) that are generated against anti- $\alpha 4\beta 7$

integrin drugs and anti-IL12p40 drugs as well as other biologics is determined with a homogeneous mobility shift assay as described in, *e.g.*, U.S. Patent Nos. 8,574,855 and 8,865,417, and U.S. Patent Publication Nos. 2014/0051184 and 2014/0141983, the disclosures of which are hereby incorporated by reference in their entirety for all purposes.

5 **[0014]** In other aspects, the present invention provides an isolated soluble $\alpha 4$ integrin polypeptide comprising an amino acid sequence having at least 80% identity to SEQ ID NO:1 or SEQ ID NO:3. In yet other aspects, the present invention provides an isolated soluble $\beta 7$ integrin polypeptide comprising an amino acid sequence having at least 80% identity to SEQ ID NO: 2 or SEQ ID NO:4.

10 **[0015]** In particular embodiments, the present invention provides an isolated soluble $\alpha 4\beta 7$ integrin heterodimer comprising:

 (a) an $\alpha 4$ integrin polypeptide having an amino acid sequence that has at least 80% identity to SEQ ID NO:1, wherein the $\alpha 4$ integrin polypeptide is linked to a first member of a binding pair (*e.g.*, SEQ ID NO:3), and

15 (b) a $\beta 7$ integrin polypeptide having an amino acid sequence that has at least 80% identity to SEQ ID NO:2, wherein the $\beta 7$ integrin polypeptide is linked to a second member of the binding pair (*e.g.*, SEQ ID NO:4).

[0016] In certain other aspects, the present invention provides an isolated soluble IL-12p40 polypeptide comprising an amino acid sequence having at least 80% identity to SEQ ID
20 NOS:6, 7, 11, 12, or 13.

[0017] In particular embodiments, the unlabeled soluble antigen used in the indirect homogeneous mobility shift assays of the present invention comprise the isolated soluble $\alpha 4\beta 7$ integrin heterodimers or isolated soluble IL-12p40 polypeptides described herein.

[0018] In further aspects, the present invention provides expression vectors encoding the
25 soluble polypeptides described herein, host cells comprising the expression vectors, and methods for generating the soluble polypeptides described herein.

[0019] Other objects, features, and advantages of the present invention will be apparent to one of skill in the art from the following detailed description and figures.

BRIEF DESCRIPTION OF THE DRAWINGS

30 **[0020]** **FIG. 1** shows the principle of the indirect homogeneous mobility shift assays (HMSA) of the present invention. Using vedolizumab (“VLM” or “Vedo”) as a non-limiting

example, in the first step, serum is added to a 96 well plate along with integrin $\alpha 4\beta 7$ and diluent buffer. In the second step, labeled VLM (*e.g.*, Vedo Alexa 488) is added. The samples are then injected sequentially on an HPLC size exclusion column.

[0021] FIG. 2 shows the chromatograms of Vedo Alexa 488 (left) and Vedo Alexa 488 plus integrin $\alpha 4\beta 7$ antigen (right) with retention times (x-axis) and light units (y-axis) of one embodiment of the present invention. Both are in 4% normal human serum. Note that the antigen binds up the majority of the labeled VLM.

[0022] FIG. 3 shows the chromatograms from standard curve overlays with the retention time (x-axis) and light units (y-axis) of the various components of a VLM assay in one embodiment of the present invention. “Vedo Alexa 488/ $\alpha 4\beta 7$ ” = Alexa Fluor[®] 488-labeled VLM bound to a soluble $\alpha 4\beta 7$ heterodimer. “Alexa 488” = Blocked (*e.g.*, inactivated) Alexa Fluor[®] 488 loading control.” Note that the free Vedo Alexa 488 peak area gets larger when there is more therapeutic VLM present.

[0023] FIG. 4 shows a VLM standard curve in one embodiment of the present invention. The standard curve was generated using a serial dilution of VLM with a concentration range of between 0.15625 $\mu\text{g/ml}$ and 80 $\mu\text{g/ml}$.

[0024] FIG. 5 shows a VLM drug assay validation in one embodiment of the present invention. The limit of blank (LOB), limit of detection (LOD), lower limit of quantitation (LLOQ), and upper limit of quantitation (ULOQ) were calculated using a standard curve generated from a serial dilution of VLM with a concentration range of between 0.15625 $\mu\text{g/ml}$ and 80 $\mu\text{g/ml}$.

[0025] FIG. 6 shows the intra-assay precision and accuracy of the VLM drug assay in one embodiment of the present invention.

[0026] FIG. 7 shows the retention time (x-axis) and light units (y-axis) of an anti-vedolizumab autoantibody (ATV) assay in one embodiment of the present invention.

[0027] FIG. 8 shows the validation of the ATV assay in one embodiment of the present invention.

[0028] FIG. 9 shows the validation of the ATV assay in one embodiment of the present invention.

[0029] FIG. 10 shows the interassay precision of the ATV assay in one embodiment of the present invention.

[0030] FIG. 11 shows a ustekinumab (UTK) drug assay validation in one embodiment of the present invention. The limit of blank (LOB), limit of detection (LOD), lower limit of quantitation (LLOQ), and upper limit of quantitation (ULOQ) were calculated using a standard curve generated from a serial dilution of UTK with a concentration range of between 0.078 µg/ml and 40 µg/ml.

[0031] FIG. 12 shows the validation of the assay for autoantibodies to ustekinumab (ATU) in one embodiment of the present invention.

[0032] FIG. 13 shows a schematic diagram of the exemplary embodiments of the soluble $\alpha 4$ integrin antigen and soluble $\beta 7$ integrin antigen of the present invention. Left: Full-length proteins. Right: Truncated $\alpha 4 \beta 7$ integrin heterodimer with a cysteine bridge of the ACID-BASE peptide.

[0033] FIG. 14 shows the linearity of drug dilution in human serum for (A) VLM and (B) autoantibodies to vedolizumab (ATV).

[0034] FIG. 15 shows an analysis of common interfering agents in serum: (A) hemolyzed serum interference; (B) lipemic serum interference; and (C) RF serum interference.

[0035] FIG. 16 shows the principle of the ustekinumab (UTK) homogeneous mobility shift assay (HMSA).

[0036] FIG. 17 shows the linearity of dilution in normal human serum for (A) UTK and (B) autoantibodies to ustekinumab (ATU).

[0037] FIG. 18 shows an analysis of common interfering agents in serum: (A) hemolyzed serum interference; (B) lipemic serum interference; and (C) RF serum interference.

[0038] FIG. 19 shows integrin $\alpha 4 \beta 7$ substrate interference in one embodiment of the VLM assay of the present invention.

[0039] FIG. 20 shows standard curves using fixed amounts of labeled vedolizumab ("Vedo488") and integrin $\alpha 4 \beta 7$ antigen (top left), as well as 2-fold increased (top right) or 4-fold increased (bottom left) amounts of both reagents. Note that the top end of the curve saturates at higher VLM concentrations. Similarly, the bottom end flattens out at slightly higher levels of VLM.

DETAILED DESCRIPTION OF THE INVENTION

I. Introduction

[0040] The present invention is based in part on the discovery that an indirect homogeneous mobility shift assay (HMSA) using size exclusion chromatography is particularly advantageous for measuring the presence or level of biologics that target antigens having one or more of the following characteristics: cell surface or membrane-bound, (heavily) glycosylated, multimeric (*e.g.*, forms heterodimers, homodimers, *etc.*), insoluble, difficult to express, difficult to purify, large in size, and combinations thereof. In certain aspects, the use of a soluble form (*e.g.*, a soluble fragment, variant, or monomer) of the antigen overcomes the difficulties and limitations associated with antigens having one or more of the above characteristics and enables the precise and accurate measurement of any biologic of interest in a sample from a patient receiving therapy with that biologic.

[0041] The principle behind the indirect assays of the present invention is that the amount of (unlabeled) biologic in a sample (*e.g.*, serum) obtained from a patient receiving biologic therapy determines how much unlabeled antigen remains free to bind to a labeled form of the biologic. By tracking changes in the area of the free (unbound) labeled biologic, the presence or level of (unlabeled) biologic in the patient sample can be determined. More particularly, the relative amount (*e.g.*, ratio) of labeled and unlabeled biologic determines how much antigen is bound to each and determines the amount (*e.g.*, peak area) of free labeled biologic following size exclusion chromatography. The amount (*e.g.*, peak area) of the free labeled biologic can then be compared to a standard curve of known amounts of the biologic to provide an accurate measurement of biologic levels in the patient sample with high sensitivity and dynamic range. In certain embodiments, the size of the peak area of free labeled biologic following size exclusion chromatography is calculated and compared to the standard curve to interpolate the concentration of biologic in a patient sample.

[0042] The importance of measuring serum concentrations of biologics is illustrated by the fact that the FDA requires pharmacokinetic and tolerability (*e.g.*, immune response) studies to be performed during clinical trials. The present invention also finds utility in monitoring patients receiving these drugs to make sure they are getting the right dose, that the drug isn't being cleared from the body too quickly, and that they are not developing an immune response against the drug. Furthermore, the present invention is useful in guiding the switch between different drugs due to failure with the initial drug.

II. Definitions

[0043] As used herein, the following terms have the meanings ascribed to them unless specified otherwise.

[0044] The terms “a,” “an,” or “the” as used herein not only include aspects with one member, but also include aspects with more than one member. For instance, the singular forms “a,” “an,” and “the” include plural referents unless the context clearly dictates otherwise.

[0045] The terms “competition,” “competition-based,” and “indirect” are used interchangeably herein to refer to an assay of the present invention for determining the presence or level of (unlabeled) biologic in a sample that relies on detecting the amount of free (unbound) labeled biologic remaining in the sample after unlabeled antigen and labeled biologic are added (sequentially) to the sample.

[0046] The terms “VLM,” “VDZ,” and “Vedo” are used interchangeably herein to refer to vedolizumab.

[0047] The terms “biologic” or “biologic agent” or “biological drug” as used herein encompass products and substances produced from or extracted from a biological system (*e.g.*, a living organism). Non-limiting examples of biologics include antibodies, antibody fragments, proteins, polypeptides, peptides, fusion proteins (*e.g.*, Ig fusion proteins or Fc fusion proteins), multivalent binding proteins (*e.g.*, DVD Ig), antibody-drug conjugates, vaccines, nucleic acids, sugars, recombinant forms thereof, engineered forms thereof, and combinations thereof.

[0048] The term “antibody” includes large (150 kDa) “Y-shaped” molecules that consist of two identical light chains (~220 amino acids) and two identical heavy chains (~440 amino acids) that are held together by a combination of covalent (disulfide) and non-covalent interactions. Each light and heavy chain consists of repeating segments of constant or variable regions that contain one intrachain disulfide bond. The variable regions are located at the N-termini of the light and heavy chains, while the constant domains are located at the C-termini of the light and heavy chains. The N-termini of the light and heavy chains come together to form the antigen-binding site. The light chain is comprised of one variable domain and one constant domain and the heavy chain is comprised of one variable domain and three constant domains. Located at the ends of the “Y” are two identical (bivalent)

antigen-binding sites. The distance between the two antigen binding sites varies due to the flexible hinge region, and as a result, the antigen binding efficiency can be greatly increased. The formation of the antigen-binding region is caused by the pairing of the variable region from the heavy chain (V_H) with the variable region of the light chain (V_L). Variations in amino acid sequences of the variable regions are responsible for the vast diversity of antigen-binding sites, and the greatest variability occurs throughout three hypervariable regions, termed complementary determining regions (CDRs). The tail region of the antibody, known as the F_C region, is comprised of two constant domains (C_{H2} , and C_{H3}) from each of the heavy chains. The F_C region is responsible for recruiting effector functions through binding of F_C receptors on neutrophils and macrophages.

[0049] The term “antigen” includes any molecule, agent, or substance that (*e.g.*, specifically) binds to or interacts with a biologic. As one non-limiting example, the antigen comprises a soluble fragment, variant, or monomer of a membrane-bound protein, a glycosylated protein, a multimeric protein, an insoluble protein, a protein that is difficult to express or purify, and/or a large protein. As another non-limiting example, the antigen comprises a soluble fragment of a cell surface molecule such an integrin receptor (*e.g.*, $\alpha 4\beta 7$ integrin), wherein the soluble fragment contains one or more extracellular domains of the corresponding full-length molecule (*e.g.*, a soluble $\alpha 4\beta 7$ antigen heterodimer comprising extracellular domain sequences from the corresponding full-length $\alpha 4$ and $\beta 7$ proteins). As yet another non-limiting example, the antigen comprises a cytokine such as $TNF\alpha$ or a subunit thereof such as IL-12p40 that does not form homodimers or heterodimers.

[0050] The term “size exclusion chromatography” or “SEC” includes a chromatographic method in which molecules in solution are separated based on their size and/or hydrodynamic volume. It is applied to large molecules or macromolecular complexes such as proteins and their conjugates. Typically, when an aqueous solution is used to transport the sample through the column, the technique is known as gel filtration chromatography.

[0051] The term “complex” includes an antigen bound (*e.g.*, by non-covalent means) to a biologic (*e.g.*, an unlabeled or labeled biologic), and a biologic (*e.g.*, a labeled biologic) bound (*e.g.*, by non-covalent means) to an autoantibody against the biologic.

[0052] As used herein, an entity that is modified by the term “labeled” includes any antigen, molecule, protein, enzyme, antibody, antibody fragment, cytokine, or related species that is conjugated with another molecule or chemical entity that is empirically detectable.

Chemical species suitable as labels include, but are not limited to, fluorescent dyes, *e.g.* Alexa Fluor[®] dyes such as Alexa Fluor[®] 488, quantum dots, optical dyes, luminescent dyes, and radionuclides, *e.g.*, ¹²⁵I.

[0053] The phrase “fluorescence label detection” includes a means for detecting a fluorescent label. Means for detection include, but are not limited to, a spectrometer, a fluorimeter, a photometer, and a detection device commonly incorporated with a chromatography instrument such as, but not limited to, size exclusion-high performance liquid chromatography, such as, but not limited to, an Agilent-1200 HPLC System.

[0054] The term “subject,” “patient,” or “individual” typically includes humans, but also includes other animals such as, *e.g.*, other primates, rodents, canines, felines, equines, ovines, porcines, and the like.

[0055] The term “sample” includes any biological specimen obtained from an individual. Samples include, without limitation, whole blood, plasma, serum, red blood cells, white blood cells (*e.g.*, peripheral blood mononuclear cells (PBMC), polymorphonuclear (PMN) cells), ductal lavage fluid, nipple aspirate, lymph (*e.g.*, disseminated tumor cells of the lymph node), bone marrow aspirate, saliva, urine, stool (*i.e.*, feces), sputum, bronchial lavage fluid, tears, fine needle aspirate (*e.g.*, harvested by random periareolar fine needle aspiration), any other bodily fluid, a tissue sample such as a biopsy of a site of inflammation (*e.g.*, needle biopsy), cellular extracts thereof, and an immunoglobulin enriched fraction derived from one or more of these bodily fluids or tissues. In some embodiments, the sample is whole blood, a fractional component thereof such as plasma, serum, or a cell pellet, or an immunoglobulin enriched fraction thereof. One skilled in the art will appreciate that samples such as serum samples can be diluted prior to the analysis. In certain embodiments, the sample is obtained by isolating PBMCs and/or PMN cells using any technique known in the art. In certain other embodiments, the sample is a tissue biopsy such as, *e.g.*, from a site of inflammation such as a portion of the gastrointestinal tract or synovial tissue.

[0056] The term “isolated,” when applied to a nucleic acid or polypeptide, denotes that the nucleic acid or polypeptide is essentially free of other cellular components with which it is associated in the natural state. It is preferably in a homogeneous state although it can be in either a dry or aqueous solution. Purity and homogeneity are typically determined using analytical chemistry techniques such as polyacrylamide gel electrophoresis or high performance liquid chromatography. A polypeptide that is the predominant species present

in a preparation is substantially purified. In particular, an isolated gene is separated from open reading frames that flank the gene and encode a protein other than the gene of interest. The term “purified” denotes that a nucleic acid or polypeptide gives rise to essentially one band in an electrophoretic gel. Particularly, it means that the nucleic acid or polypeptide is at least about 80% pure, at least about 85% pure, at least about 90% pure, at least about 95% pure, or at least about 99% pure.

[0057] The term “soluble,” in the context of a polypeptide, refers to polypeptide that can be prepared in a soluble and functional form using a host cell or a cell-free protein synthesis system. For instance, a soluble polypeptide does not form insoluble aggregates comprising misfolded and/or functionally inactive polypeptides.

[0058] The term “nucleic acid” or “polynucleotide” includes deoxyribonucleotides or ribonucleotides and polymers thereof in either single- or double-stranded form. Unless specifically limited, the term encompasses nucleic acids containing known analogues of natural nucleotides that have similar binding properties as the reference nucleic acid and are metabolized in a manner similar to naturally occurring nucleotides. Unless otherwise indicated, a particular nucleic acid sequence also implicitly encompasses conservatively modified variants thereof (*e.g.*, degenerate codon substitutions), alleles, orthologs, SNPs, and complementary sequences as well as the sequence explicitly indicated. Specifically, degenerate codon substitutions may be achieved by generating sequences in which the third position of one or more selected (or all) codons is substituted with mixed-base and/or deoxyinosine residues (Batzer *et al.*, *Nucleic Acid Res.* 19:5081 (1991); Ohtsuka *et al.*, *J. Biol. Chem.* 260:2605-2608 (1985); and Cassol *et al.* (1992); Rossolini *et al.*, *Mol. Cell. Probes* 8:91-98 (1994)). The term nucleic acid is used interchangeably with gene, cDNA, and mRNA encoded by a gene.

[0059] The terms “polypeptide,” “peptide,” and “protein” are used interchangeably herein to include a polymer of amino acid residues. The terms apply to amino acid polymers in which one or more amino acid residue is an artificial chemical mimetic of a corresponding naturally occurring amino acid, as well as to naturally occurring amino acid polymers and non-naturally occurring amino acid polymers. As used herein, the terms encompass amino acid chains of any length, including full-length proteins, truncated forms, or fragments thereof, wherein the amino acid residues are linked by covalent peptide bonds.

[0060] The term “amino acid” refers to naturally-occurring and unnatural amino acids, as well as amino acid analogs and amino acid mimetics that function in a manner similar to the naturally occurring amino acids.

[0061] Amino acids may be referred to herein by either their name, their commonly known three letter symbols or by the one-letter symbols recommended by the IUPAC-IUB Biochemical Nomenclature Commission. Additionally, nucleotides may be referred to by their commonly accepted single-letter codes.

[0062] The steps of the methods of the present invention do not necessarily have to be performed in the particular order in which they are presented. A person of ordinary skill in the art would understand that other orderings of the steps of the methods of the invention are encompassed within the scope of the present invention.

[0063] Brackets, “[]” indicate that the species within the brackets are referred to by their concentration.

III. Description of the Embodiments

[0064] The present invention provides novel indirect homogeneous mobility shift assays for detecting and measuring the presence or level of a biologic in a sample. The assays of the present invention are particularly advantageous for detecting the presence or level of biologics that target complex or large antigens including cell surface proteins, transmembrane proteins, heavily glycosylated proteins, and multimeric proteins, as well as antigens that cannot be purified, impure antigens, and partially or substantially purified antigens. The present invention also provides isolated soluble $\alpha 4 \beta 7$ integrin heterodimers and isolated soluble IL-12p40 monomers that are suitable for use in the assays described herein.

[0065] In one aspect, the present invention provides a method for determining the presence or level of a biologic in a sample, the method comprising:

(a) contacting the sample with an unlabeled soluble antigen that binds to the biologic to form an unlabeled complex (*e.g.*, a plurality of unlabeled complexes) between the antigen and the biologic in the sample;

(b) contacting the sample from step (a) with a labeled form of the biologic (“labeled biologic”) to form a labeled complex (*e.g.*, a plurality of labeled complexes) between the antigen and the labeled biologic;

(c) subjecting the (*e.g.*, plurality of) unlabeled and labeled complexes to size exclusion chromatography to separate the (*e.g.*, plurality of) unlabeled and labeled complexes from free labeled biologic and to detect an amount of the free labeled biologic; and

(d) comparing the amount of the free labeled biologic detected in step (c) to a standard curve of known amounts of the biologic, thereby determining the presence or level of the biologic in the sample.

[0066] In some embodiments, the biologic includes antibodies, antibody fragments, proteins, polypeptides, peptides, fusion proteins, multivalent binding proteins, antibody-drug conjugates, vaccines, nucleic acids, sugars, recombinant forms thereof, engineered forms thereof, and combinations thereof. In particular embodiments, the biologic comprises an antibody (*e.g.*, a monoclonal antibody) or a fragment thereof (*e.g.*, an antigen-binding fragment of a monoclonal antibody) or a conjugate thereof (*e.g.*, an antibody-drug conjugate). Non-limiting examples of antibody-based biologics are shown in Table 1.

[0067] In particular embodiments, the method of the present invention detects the presence of and/or measures the level of unbound (free) biologic in a sample, *e.g.*, the population of biologic in a sample that is not bound to its (endogenous) target antigen or a fragment thereof.

[0068] In certain embodiments, the antigen is a soluble version (*e.g.*, a soluble fragment, variant, or monomeric form) of a membrane-bound protein, a (heavily) glycosylated protein, a multimeric protein, an insoluble protein, a protein that is difficult to express or purify, and/or a large protein. In certain instances, the antigen is a soluble extracellular domain of a membrane-bound protein (*e.g.*, a soluble cytokine receptor extracellular domain). In certain other instances, the antigen is a soluble homodimer or heterodimer comprising the extracellular domains of two membrane-bound proteins (*e.g.*, a soluble integrin heterodimer). In yet other instances, the antigen is a soluble protein that does not multimerize and remains in monomeric form once isolated and/or purified (*e.g.*, a soluble cytokine variant with one or more cysteine residues mutated to minimize or eliminate the formation of multimers).

[0069] In some embodiments, the antigen is a soluble fragment (*e.g.*, extracellular domain) of a cell surface molecule such as, *e.g.*, a cell adhesion molecule (CAM). Non-limiting examples of CAMs include immunoglobulin superfamily (IgSF) CAMs, integrins, cadherins, and selectins.

[0070] IgSF CAMs are any of a variety of polypeptides or proteins located on the surface of a cell that have one or more immunoglobulin-like fold domains, and which function in intercellular adhesion and/or signal transduction. In many cases, IgSF CAMs are transmembrane proteins. Non-limiting examples of IgSF CAMs include mucosal addressin cell adhesion molecule 1(MADCAM1), neural cell adhesion molecules (NCAMs; *e.g.*, NCAM-120, NCAM-125, NCAM-140, NCAM-145, NCAM-180, NCAM-185, *etc.*), intercellular adhesion molecules (ICAMs, *e.g.*, ICAM-1, ICAM-2, ICAM-3, ICAM-4, and ICAM-5), vascular cell adhesion molecule-1 (VCAM-1), platelet-endothelial cell adhesion molecule-1 (PECAM-1), L1 cell adhesion molecule (L1CAM), cell adhesion molecule with homology to L1CAM (close homolog of L1) (CHL1), sialic acid binding Ig-like lectins (SIGLECs; *e.g.*, SIGLEC-1, SIGLEC-2, SIGLEC-3, SIGLEC-4, *etc.*), nectins (*e.g.*, Nectin-1, Nectin-2, Nectin-3, *etc.*), and nectin-like molecules (*e.g.*, Necl-1, Necl-2, Necl-3, Necl-4, and Necl-5).

[0071] Integrins are transmembrane $\alpha\beta$ heterodimers and at least 18 α and eight β subunits are known in humans, generating 24 heterodimers. The α and β subunits have distinct domain structures, with extracellular domains from each subunit contributing to the ligand-binding site of the heterodimer. Non-limiting examples of integrins include $\alpha_1\beta_1$, $\alpha_2\beta_1$, $\alpha_3\beta_1$, $\alpha_4\beta_1$, $\alpha_5\beta_1$, $\alpha_6\beta_1$, $\alpha_7\beta_1$, $\alpha_8\beta_1$, $\alpha_9\beta_1$, $\alpha_{10}\beta_1$, $\alpha_{11}\beta_1$, $\alpha_v\beta_1$, $\alpha_v\beta_3$, $\alpha_v\beta_5$, $\alpha_v\beta_6$, $\alpha_v\beta_8$, $\alpha_{IIB}\beta_3$, $\alpha_4\beta_7$, $\alpha_E\beta_7$, $\alpha_6\beta_4$, $\alpha_L\beta_2$, $\alpha_M\beta_2$, $\alpha_X\beta_2$, and $\alpha_D\beta_2$.

[0072] In particular embodiments, the antigen is an $\alpha_4\beta_7$ integrin and the biologic is an anti- $\alpha_4\beta_7$ integrin drug such as vedolizumab (VLM). In certain instances, the soluble fragment of the $\alpha_4\beta_7$ integrin that binds to the anti- $\alpha_4\beta_7$ integrin drug comprises an α_4 fragment comprising an amino acid sequence having at least 80% identity to SEQ ID NO:1 or SEQ ID NO:3 and/or a β_7 fragment comprising an amino acid sequence having at least 80% identity to SEQ ID NO: 2 or SEQ ID NO:4.

[0073] In other embodiments, the antigen is an $\alpha_4\beta_1$ integrin and the biologic is an anti- $\alpha_4\beta_1$ integrin drug such as natalizumab. In certain instances, the soluble fragment of the $\alpha_4\beta_1$ integrin that binds to the anti- $\alpha_4\beta_1$ integrin drug comprises a heterodimer of the extracellular domains of the α_4 and β_1 subunits.

[0074] Cadherins are calcium-dependent transmembrane proteins that play important roles in cell adhesion, forming adherens junctions to bind cells within tissues together. Non-

limiting examples of cadherins include E-cadherin, N-cadherin, N-cadherin 2, and P-cadherin.

[0075] Selectins are heterophilic CAMs that bind fucosylated carbohydrates, *e.g.*, mucins. The three family members are E-selectin (endothelial), L-selectin (leukocyte), and P-selectin (platelet).

[0076] In other embodiments, the antigen is a soluble fragment (*e.g.*, extracellular domain) of a cell surface molecule such as, *e.g.*, a cytokine receptor.

[0077] Non-limiting examples of cytokine receptors include type I cytokine receptors, type II cytokine receptors, members of the immunoglobulin (Ig) superfamily, tumor necrosis factor receptors, chemokine receptors, and TGF β receptors. Examples of type I cytokine receptors include, but are not limited to, interleukin receptors (*e.g.*, IL-2 receptor, IL-3 receptor, IL-4 receptor, IL-5 receptor, IL-6 receptor, IL-7 receptor, IL-9 receptor, IL-11 receptor, IL-12 receptor, IL-13 receptor, IL-15 receptor, IL-21 receptor, IL-23 receptor, IL-27 receptor, *etc.*), colony stimulating factor receptors (*e.g.*, erythropoietin receptor, GM-CSF receptor, G-CSF receptor, *etc.*), hormone receptors or neuropeptide receptors (*e.g.*, growth hormone receptor, prolactin receptor, *etc.*), and other cytokine receptors such as oncostatin M receptor and leukemia inhibitory factor receptor. Examples of type II cytokine receptors include, but are not limited to, interferon receptors (*e.g.*, interferon-alpha/beta receptor, interferon-gamma receptor, *etc.*) and interleukin receptors (*e.g.*, IL-10 receptor, IL-20 receptor, IL-22 receptor, IL-28 receptor, *etc.*). Examples of immunoglobulin (Ig) superfamily receptors include, but are not limited to, IL-1 receptor, CSF1, c-kit receptor, and IL-18 receptor. Examples of tumor necrosis factor receptors include, but are not limited to, TNF receptor (CD120), lymphotoxin β receptor, CD134, CD40, FAS, TNFRSF6B, CD27, CD30, CD137, TNFRSF10A, TNFRSF10B, TNFRSF10C, TNFRSF10D, RANK, osteoprotegerin, TNFRSF12A, TNFRSF13B, TNFRSF13C, TNFRSF14, nerve growth factor receptor, TNFRSF17, TNFRSF18, TNFRSF19, TNFRSF21, TNFRSF25, and ectodysplasin A2 receptor. Examples of chemokine receptors include, but are not limited to, CXC chemokine receptors, CC chemokine receptors, C chemokine receptors, and CX3C chemokine receptors. Examples of TGF β receptors include, but are not limited to, TGF β receptor 1, TGF β receptor 2, and TGF β receptor 3.

[0078] In certain embodiments, the antigen is an IL-6 receptor and the biologic is an anti-IL-6 receptor drug such as tocilizumab. In certain instances, the soluble fragment of the IL-6

receptor that binds to the anti-IL-6 receptor drug comprises an extracellular domain of the IL-6 receptor.

[0079] In yet other embodiments, the antigen is a soluble fragment (*e.g.*, extracellular domain) of a cluster of differentiation (CD) molecule. Non-limiting examples of CD molecules include CD3, CD4, CD8, CD11a, CD11b, CD14, CD15, CD16, CD19, CD20, CD22, CD24, CD25, CD30, CD31, CD34, CD38, CD45, CD56, CD61, CD91, CD114, CD117, CD182, and the like. In certain instances, the biologic that binds to a soluble fragment of a CD molecule is a member selected from the group consisting of visilizumab, priliximab, rituximab, ofatumumab, obinutuzumab, ibritumomab tiuxetan, tositumomab, ocrelizumab, veltuzumab, daclizumab, and combinations thereof.

[0080] In some embodiments, the antigen is a cytokine or a monomer thereof (*e.g.*, a soluble cytokine variant with one or more cysteine residues mutated to minimize or eliminate the formation of multimers).

[0081] Non-limiting examples of cytokines include TNF α , TNF-related weak inducer of apoptosis (TWEAK), osteoprotegerin (OPG), IFN- α , IFN- β , IFN- γ , interleukins (*e.g.*, IL-1 α , IL-1 β , IL-1 receptor antagonist (IL-1ra), IL-2, IL-4, IL-5, IL-6, soluble IL-6 receptor (sIL-6R), IL-7, IL-8, IL-9, IL-10, IL-12, IL-13, IL-15, IL-17, IL-23, and IL-27), adipocytokines (*e.g.*, leptin, adiponectin, resistin, active or total plasminogen activator inhibitor-1 (PAI-1), visfatin, and retinol binding protein 4 (RBP4)), and the like.

[0082] In particular embodiments, the cytokine is a p40 subunit of IL-12 or IL-23 and the biologic is an anti-IL-12p40 drug such as ustekinumab (UTK). In certain instances, the cytokine is a p40 variant which comprises one or more cysteine residues mutated to minimize or eliminate the formation of multimers. In some instances, the p40 variant comprises an amino acid sequence having at least 80% identity to SEQ ID NOS:6, 7, 11, 12 or 13.

[0083] In other embodiments, the cytokine is TNF α and the biologic is an anti-TNF α drug. Non-limiting examples of anti-TNF α drugs include REMICADE[®] (infliximab), HUMIRA[®] (adalimumab), ENBREL[®] (etanercept), CIMZIA[®] (certolizumab pegol), SIMPONI[®] (golimumab), and combinations thereof.

[0084] The soluble antigens described herein can be produced by any method known to one of ordinary skill in the art, such as but not limited to, synthetic methods, such as solid phase and liquid phase synthesis, or recombinant biology methods.

[0085] In some embodiments, the sample is a whole blood, serum, or plasma sample, *e.g.*, obtained from a subject receiving biologic therapy. In preferred embodiments, the sample is serum. In particular embodiments, the subject has a disease or disorder such as, *e.g.*, an autoimmune disease (*e.g.*, rheumatoid arthritis), an inflammatory disease (*e.g.*, inflammatory bowel disease (IBD) such as Crohn's disease (CD) or ulcerative colitis (UC)), or cancer.

[0086] In particular embodiments, the standard curve is generated by incubating the antigen and the labeled biologic with a (*e.g.*, two-fold) serial dilution of known amounts of the biologic. In certain embodiments, the area under the curve (AUC) of the free labeled biologic is plotted against (*e.g.*, the logarithm of) known amounts of the biologic obtained from the standard curve, and the level of the biologic in the sample is calculated by interpolation, *e.g.*, based upon the size of the peak area of the free labeled biologic. In other embodiments, a ratio of the free labeled biologic to a loading control (*e.g.*, free label) is determined and used to normalize the level of the biologic in the sample from the standard curve.

[0087] In certain embodiments, the size exclusion chromatography (SEC) is size exclusion-high performance liquid chromatography (SE-HPLC). In particular embodiments, the (*e.g.*, plurality of) unlabeled and labeled complexes are eluted first through a stationary phase, followed by the free labeled biologic. The underlying principle of SEC is that molecules or complexes of different sizes will elute (filter) through a stationary phase at different rates. This results in the separation of a solution of molecules or complexes based on size. Provided that all the molecules or complexes are loaded simultaneously or near simultaneously, molecules or complexes of the same size elute together. Each size exclusion column has a range of molecular weights that can be separated. The exclusion limit defines the molecular weight at the upper end of this range and is where molecules or complexes are too large to be trapped in the stationary phase. The permeation limit defines the molecular weight at the lower end of the range of separation and is where molecules or complexes of a small enough size can penetrate into the pores of the stationary phase completely and all molecules or complexes below this molecular mass are so small that they elute as a single band.

[0088] In some instances, the eluent is collected in constant volumes, or fractions. The more similar the molecules or complexes are in size, the more likely they will be in the same fraction and not detected separately. Preferably, the collected fractions are examined by

spectroscopic techniques to determine the concentration of the molecules or complexes eluted. Typically, the spectroscopy detection techniques useful in the present invention include, but are not limited to, fluorometry, refractive index (RI), and ultraviolet (UV). In certain instances, the elution volume decreases roughly linearly with the logarithm of the molecular hydrodynamic volume (*i.e.*, heavier molecules or complexes come off first).

[0089] A biologic (*e.g.*, therapeutic antibody) can be labeled with any of a variety of detectable group(s). In certain embodiments, a biologic is labeled with a fluorophore or a fluorescent dye. In other embodiments, a biologic is labeled with a luminescent tag, a metal, a radionuclide, and the like. Specific immunological binding of an antigen to a labeled biologic or the amount of free labeled biologic can be detected directly or indirectly. A signal from the direct or indirect label can be analyzed, *e.g.*, using a spectrophotometer to detect color from a chromogenic substrate, a radiation counter to detect radiation such as a gamma counter for detection of ^{125}I , or a fluorometer to detect fluorescence in the presence of light of a certain wavelength.

[0090] Non-limiting examples of fluorophores or fluorescent dyes include those listed in the Molecular Probes Catalogue, which is herein incorporated by reference (*see*, R. Haugland, *The Handbook-A Guide to Fluorescent Probes and Labeling Technologies*, 10th Edition, Molecular probes, Inc. (2005)). Such exemplary fluorophores or fluorescent dyes include, but are not limited to, Alexa Fluor[®] dyes such as Alexa Fluor[®] 350, Alexa Fluor[®] 405, Alexa Fluor[®] 430, Alexa Fluor[®] 488, Alexa Fluor[®] 514, Alexa Fluor[®] 532, Alexa Fluor[®] 546, Alexa Fluor[®] 555, Alexa Fluor[®] 568, Alexa Fluor[®] 594, Alexa Fluor[®] 610, Alexa Fluor[®] 633, Alexa Fluor[®] 635, Alexa Fluor[®] 647, Alexa Fluor[®] 660, Alexa Fluor[®] 680, Alexa Fluor[®] 700, Alexa Fluor[®] 750, and/or Alexa Fluor[®] 790, as well as other fluorophores including, but not limited to, Dansyl Chloride (DNS-Cl), 5-(iodoacetamido)fluorescein (5-IAF), fluorescein 5-isothiocyanate (FITC), tetramethylrhodamine 5- (and 6-)isothiocyanate (TRITC), 6-acryloyl-2-dimethylaminonaphthalene (acrylodan), 7-nitrobenzo-2-oxa-1,3,-diazol-4-yl chloride (NBD-Cl), ethidium bromide, Lucifer Yellow, 5-carboxyrhodamine 6G hydrochloride, Lissamine rhodamine B sulfonyl chloride, Texas Red[™] sulfonyl chloride, BODIPY[™], naphthalamine sulfonic acids (*e.g.*, 1-anilinonaphthalene-8-sulfonic acid (ANS), 6-(p-toluidinyl)naphthalene-2-sulfonic acid (TNS), and the like), Anthroyl fatty acid, DPH, Parinaric acid, TMA-DPH, Fluorenyl fatty acid, fluorescein-phosphatidylethanolamine, Texas Red-phosphatidylethanolamine, Pyrenyl-phosphatidylcholine, Fluorenyl-phosphatidylcholine, Merocyanine 540, 1-(3-sulfonatopropyl)-4-[β -[2[(di-n-butylamino)-6

naphthyl]vinyl]pyridinium betaine (Naphtyl Styryl), 3,3'dipropylthiadicarbocyanine (diS-C₃-(5)), 4-(p-dipentyl aminostyryl)-l-methylpyridinium (di-5-ASP), Cy-3 Iodo Acetamide, Cy-5-N-Hydroxysuccinimide, Cy-7-Isothiocyanate, rhodamine 800, IR-125, Thiazole Orange, Azure B, Nile Blue, Al Phthalocyanine, Oxaxine 1, 4', 6-diamidino-2-phenylindole (DAPI), Hoechst 33342, TOTO, Acridine Orange, Ethidium Homodimer, N(ethoxycarbonylmethyl)-6-methoxyquinolinium (MQAE), Fura-2, Calcium Green, Carboxy SNARF-6, BAPTA, coumarin, phytofluors, Coronene, metal-ligand complexes, IRDye[®] 700DX, IRDye[®] 700, IRDye[®] 800RS, IRDye[®] 800CW, IRDye[®] 800, Cy5, Cy5.5, Cy7, DY 676, DY680, DY682, DY780, and mixtures thereof. Additional suitable fluorophores include enzyme-cofactors; lanthanide, green fluorescent protein, yellow fluorescent protein, red fluorescent protein, or mutants and derivatives thereof.

[0091] Typically, the fluorescent group is a fluorophore selected from the category of dyes comprising polymethines, pthalocyanines, cyanines, xanthenes, fluorenes, rhodamines, coumarins, fluoresceins and BODIPY[™].

[0092] In certain embodiments, the fluorescent group is a near-infrared (NIR) fluorophore that emits in the range of between about 650 to about 900 nm. Use of near infrared fluorescence technology is advantageous in biological assays as it substantially eliminates or reduces background from auto fluorescence of biosubstrates. Another benefit to the near-IR fluorescent technology is that the scattered light from the excitation source is greatly reduced since the scattering intensity is proportional to the inverse fourth power of the wavelength. Low background fluorescence and low scattering result in a high signal to noise ratio, which is essential for highly sensitive detection. Furthermore, the optically transparent window in the near-IR region (650 nm to 900 nm) in biological tissue makes NIR fluorescence a valuable technology for *in vivo* imaging and subcellular detection applications that require the transmission of light through biological components. Within aspects of this embodiment, the fluorescent group is preferably selected from the group consisting of IRDye[®] 700DX, IRDye[®] 700, IRDye[®] 800RS, IRDye[®] 800CW, IRDye[®] 800, Alexa Fluor[®] 660, Alexa Fluor[®] 680, Alexa Fluor[®] 700, Alexa Fluor[®] 750, Alexa Fluor[®] 790, Cy5, Cy5.5, Cy7, DY 676, DY680, DY682, and DY780. In certain embodiments, the near infrared group is IRDye[®] 800CW, IRDye[®] 800, IRDye[®] 700DX, IRDye[®] 700, or Dynamic DY676.

[0093] Fluorescent labeling can be accomplished using a chemically reactive derivative of a fluorophore. Common reactive groups include amine reactive isothiocyanate derivatives

such as FITC and TRITC (derivatives of fluorescein and rhodamine), amine reactive succinimidyl esters such as NHS-fluorescein, and sulfhydryl reactive maleimide activated fluors such as fluorescein-5-maleimide, many of which are commercially available. Reaction of any of these reactive dyes with a biologic results in a stable covalent bond formed between a fluorophore and a biologic.

[0094] In certain instances, following a fluorescent labeling reaction, it is often necessary to remove any non-reacted fluorophore from the labeled target molecule. This is often accomplished by size exclusion chromatography, taking advantage of the size difference between fluorophore and labeled protein.

[0095] Reactive fluorescent dyes are available from many sources. They can be obtained with different reactive groups for attachment to various functional groups within the target molecule. They are also available in labeling kits that contain all the components to carry out a labeling reaction. In certain instances, Alexa Fluor[®] 488 NHS ester is used from Life Technologies (Cat. No. A-10235).

IV. Indirect Homogeneous Mobility Shift Assays

[0096] The present invention provides novel indirect assays for detecting and measuring the presence or level of a biologic (“drug”) in a sample using size exclusion chromatography. The assays of the present invention are particularly advantageous for detecting the presence or level of drugs that target complex or large antigens including cell surface proteins, transmembrane proteins, heavily glycosylated proteins, and multimeric proteins, as well as antigens that cannot be purified, impure antigens, and partially or substantially purified antigens. The antigens are not labeled and thus patient drug/antigen complexes do not appear in the chromatogram. The principle behind the indirect assays is that the amount of patient drug determines how much unlabeled antigen remains free to bind to a labeled version of the drug. By tracking changes in the area of the free (unbound) labeled drug, one can determine how much patient drug is present.

[0097] In certain aspects, the first step of the indirect assays described herein comprises incubating a sample (*e.g.*, serum) containing therapeutic drug (*e.g.*, vedolizumab (VDZ)) with a fixed amount of antigen to the drug (*e.g.*, soluble $\alpha 4\beta 7$). In the second step, a fixed amount of labeled drug (*e.g.*, VDZ coupled with Alexa Fluor[®] 488) is added. The amount of therapeutic drug in the sample determines how much antigen remains free and available to bind the labeled drug. This, in turn, determines how much labeled drug is free. Since the

peak area of the free labeled drug is proportional to the amount of therapeutic drug in the sample, one can quantify the amount of therapeutic drug by interpolation against a standard curve containing known amounts of drug.

[0098] The following description of the principles of the indirect assays of the invention uses vedolizumab (VDZ) as the therapeutic drug for illustrative purposes only (*see*, Figure 1), and is not intended to limit the scope of the assay methodology for detecting or measuring the presence or level of other biologics in patient samples:

1. To each patient sample (*e.g.*, serum), fixed amounts of antigen (*e.g.*, soluble $\alpha 4\beta 7$) and labeled VDZ are added. The amount of antigen and labeled VDZ can be added to each sample in a controlled ratio. For example, adding an amount of antigen which would bind up about 75-80% of the labeled VDZ provides optimal sensitivity without limiting the dynamic range of the assay. The ratio of antigen to labeled VDZ was determined by titrating the antigen with a fixed amount of labeled VDZ so that when the antigen is added to the labeled VDZ, the peak of free labeled VDZ is reduced by about 75-80% (*see*, Figure 2).

2. Quantification can be performed by tracking the increase of the labeled VDZ peak area ($R_t = 7.5-8.5$ min). This area is proportional to the amount of therapeutic drug present. Tris-blocked Alexa488 can be added to labeled VDZ stock solutions as a loading control. Raw chromatograms can be collected in Agilent ChemStation and then exported to the program "R" for automated analysis. The standard curve can be generated by plotting the labeled VDZ peak area as a function of the log of known VDZ sample concentrations. A 10-point standard curve can be used and fitted with a 5-parameter logistic (5-PL) model to account for asymmetry. Unknowns can be determined from the standard curve by interpolation.

[0099] In certain embodiments, the ratio of antigen to labeled drug that is added to a sample is an amount of each reagent that provides the best compromise between the low-end sensitivity needed as well as a dynamic range that enables the measurement of drug in patient samples without requiring dilutions. As a non-limiting example, the ratio of antigen to labeled drug that is added to a sample is an amount of antigen that binds up about 75% to about 80% (*e.g.*, about 75%, 76%, 77%, 78%, 79%, or 80%) of the labeled drug. In some

instances, the ratio of antigen to labeled drug that is added to a sample is an amount of antigen that binds up at least about 50%, 55%, 60%, 65%, 70%, 75%, 80%, 85%, 90%, or 95% of the labeled drug. In other instances, the ratio of antigen to labeled drug that is added to a sample is an amount of antigen that binds up about 50% to about 90%, about 60% to about 90%, about 70% to about 90%, about 80% to about 90%, about 50% to about 80%, about 60% to about 80%, about 70% to about 80%, about 50% to about 70%, about 60% to about 70%, or about 50% to about 60% of the labeled drug. The ratio of antigen to labeled drug can be determined by titrating the antigen with a fixed amount of labeled drug so that when the antigen is added to the labeled drug, the peak of free labeled drug is reduced by a desired percent (*e.g.*, about 75-80%).

[0100] In certain other embodiments, the dynamic range of the indirect assays described herein can be improved by proportionately increasing the amount of both antigen and labeled drug that is added to a sample. In some instances, the amount of both antigen and labeled drug can be about 1.5-fold, 1.6-fold, 1.7-fold, 1.8-fold, 1.9-fold, 2-fold, 2.5-fold, 3-fold, 3.5-fold, 4-fold, 4.5-fold, 5-fold, 5.5-fold, 6-fold, 6.5-fold, 7-fold, 7.5-fold, 8-fold, 8.5-fold, 9-fold, 9.5-fold, or 10-fold more than a reference amount of antigen and labeled drug. As a non-limiting example, the reference amount of labeled VDZ can be about 75 ng and the increased amount of labeled VDZ can be about 120 ng (*i.e.*, 1.6-fold more than the reference amount).

[0101] In some embodiments, the lower limit of quantitation (LLOQ) of the indirect assays described herein is about 0.125 µg/mL, 0.25 µg/mL, 0.375 µg/mL, 0.5 µg/mL, 0.625 µg/mL, 0.75 µg/mL, 0.875 µg/mL, 1 µg/mL, 1.25 µg/mL, 1.5 µg/mL, 1.75 µg/mL, 2 µg/mL, 3 µg/mL, 4 µg/mL, or 5 µg/mL. In other embodiments, the upper limit of quantitation (ULOQ) of the indirect assays described herein is about 8 µg/mL, 9 µg/mL, 10 µg/mL, 11 µg/mL, 12 µg/mL, 13 µg/mL, 14 µg/mL, 15 µg/mL, 16 µg/mL, 17 µg/mL, 18 µg/mL, 19 µg/mL, 20 µg/mL, 21 µg/mL, 22 µg/mL, 23 µg/mL, 24 µg/mL, 25 µg/mL, 26 µg/mL, 27 µg/mL, 28 µg/mL, 29 µg/mL, 30 µg/mL, 35 µg/mL, 40 µg/mL, 45 µg/mL, or 50 µg/mL. In particular embodiments, the LLOQ is about 1 µg/mL and the ULOQ is about 25 µg/mL.

V. Soluble $\alpha 4\beta 7$ Integrin Polypeptide Antigens

[0102] In one aspect, the present invention provides an isolated soluble $\alpha 4$ integrin polypeptide comprising an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%,

94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:1. In some embodiments, the isolated soluble $\alpha 4$ integrin polypeptide comprises β -propeller repeats (*i.e.*, repeats 1-7) and a thigh domain of the human $\alpha 4$ integrin extracellular domain (*see*, Figure 13; “ $\alpha 4\Delta 620$ ”), or a fragment thereof. In other embodiments, the isolated soluble $\alpha 4$ integrin polypeptide comprises β -propeller repeats, a thigh domain, and one or both Calf domains (*i.e.*, Calf-1 and/or Calf-2) of the human $\alpha 4$ integrin extracellular domain, or a fragment thereof. In yet other embodiments, the isolated soluble $\alpha 4$ integrin polypeptide is a truncated receptor comprising the entire human $\alpha 4$ integrin extracellular domain. The isolated soluble $\alpha 4$ integrin polypeptide described herein includes a ligand binding domain or a portion thereof.

[0103] In some embodiments, the isolated soluble $\alpha 4$ integrin polypeptide also includes an ACID peptide. Such a peptide may have an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:8. An ACID peptide can form an α -helical coiled-coil conformation with a BASE peptide. In some embodiments, the ACID peptide includes a cysteine residue that can form a disulfide bridge with a cysteine residue on the BASE peptide. Acid coiled-coil region peptides (ACID peptides) and basic coiled-coil region peptides (BASE peptides) are described in, *e.g.*, O’Shea et al., Curr Biol, 1993, 3:658-667, Jun et al., Proc Natl Acad Sci U.S.A., 2001, 98(12):6830-6835, Takagi et al., Nat Struct. Biol., 2001, 8:412-416, Nishida et al., Immunity, 2006, 25:583-594, and Dong et al., Biochemistry, 2012, 51(44):8814-8828.

[0104] In some embodiments, the isolated soluble $\alpha 4$ integrin polypeptide includes a linker, such as one or more amino acid residues, *e.g.*, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10 or more amino acid residues. The linker can be located between the end of extracellular domain and the ACID peptide.

[0105] In some embodiments, the isolated soluble $\alpha 4$ integrin polypeptide comprises an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:3.

[0106] The soluble $\alpha 4$ integrin polypeptide can also include an affinity tag or epitope tag, such as a histidine tag, avidin tag, V5 tag, FLAG tag, HA tag, Myc tag, cleavable tag, and the like. In some instances, the soluble $\alpha 4$ integrin polypeptide can include a fluorescent tag,

such as GFP, DsRed, CFP, YFP, RFP, and the like, or other detectable tag, such as horseradish peroxidase, chloramphenicol acetyltransferase, beta-galactosidase, luciferase, and the like.

[0107] In another aspect, the present invention provides an isolated soluble $\beta 7$ integrin polypeptide comprising an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:2. In some embodiments, the isolated soluble $\beta 7$ integrin polypeptide comprises a PSI domain, an I-like domain, hybrid domains, and an I-EGF 1 domain of the human $\beta 7$ integrin extracellular domain (*see*, Figure 13; " $\beta 7\Delta 527$ "), or a fragment thereof. In other embodiments, the isolated soluble $\beta 7$ integrin polypeptide comprises one or more of the PSI domain, I-like domain, and one or both hybrid domains of the human $\beta 7$ integrin extracellular domain, or a fragment thereof. In yet other embodiments, the isolated soluble $\beta 7$ integrin polypeptide comprises a PSI domain, an I-like domain, hybrid domains, I-EGF domains (*i.e.*, domains 1-4), and optionally a β -tail of the human $\beta 7$ integrin extracellular domain, or a fragment thereof. In further embodiments, the isolated soluble $\beta 7$ integrin polypeptide is a truncated receptor comprising the entire human $\beta 7$ integrin extracellular domain. The soluble $\beta 7$ integrin polypeptide described herein includes a ligand binding domain or a portion thereof.

[0108] In some embodiments, the isolated soluble $\beta 7$ integrin polypeptide also includes a BASE peptide. Such a peptide may have an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:9. A BASE peptide can form an α -helical coiled-coil conformation with an ACID peptide.

[0109] In some embodiments, the isolated soluble $\beta 7$ integrin polypeptide includes a protease cleavage site. In some instances, the cleavage site is a tobacco etch virus (TEV) protease cleavage site. The TEV site can include the amino acid sequence EXXYXQ/S, wherein X is any amino acid residue (SEQ ID NO:10). The TEV site may be located upstream the BASE peptide.

[0110] In some embodiments, the isolated soluble $\beta 7$ integrin polypeptide includes a linker, such as one or more amino acid residues, *e.g.*, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10 or more amino acid

residue. The linker can be located between the end of the EGF-I domain of the $\beta 7$ integrin extracellular domain and the cleavage site.

[0111] In some embodiments, the isolated soluble $\beta 7$ integrin polypeptide further comprises an affinity tag or epitope tag. Useful affinity or epitope tags include, but are not limited to, a histidine tag, avidin tag, V5 tag, FLAG tag, HA tag, Myc tag, cleavable tag, and the like. In some instances, the soluble $\beta 7$ integrin polypeptide can include a fluorescent tag, such as GFP, DsRed, CFP, YFP, RFP, and the like, or other detectable tag, such as horseradish peroxidase, chloramphenicol acetyltransferase, beta-galactosidase, luciferase, and the like.

[0112] In some embodiments, the isolated soluble $\beta 7$ integrin polypeptide comprises an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:4.

[0113] The isolated soluble $\beta 7$ integrin polypeptide can associate with a soluble $\alpha 4$ integrin polypeptide described herein to form a complex, such as a covalently linked heterodimer. In some instances, the $\alpha 4\beta 7$ integrin complex can bind to $\alpha 4\beta 7$ ligands, such as but not limited to VCAM-1 and MAdCAM-1, and to antibodies directed against $\alpha 4\beta 7$ integrin, such as, but not limited to, vedolizumab, natalizumab, and etrolizumab.

[0114] In some aspects, the present invention provides an isolated soluble $\alpha 4\beta 7$ integrin polypeptide comprising a soluble $\alpha 4$ integrin polypeptide having an amino acid sequence that has at least 80% identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO: 1, wherein the $\alpha 4$ integrin polypeptide is linked to a first member of a binding pair, and a soluble $\beta 7$ integrin polypeptide having an amino acid sequence that has at least 80% identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO: 2, wherein the $\beta 7$ integrin polypeptide is linked to a second member of the binding pair.

[0115] In some embodiments, the binding pair can be any peptides, molecules, motifs, and compounds that can allow the $\alpha 4$ subunit and $\beta 7$ subunit to form a heterodimer, such as a covalently linked heterodimer. The $\alpha 4\beta 7$ heterodimer is capable of binding to $\alpha 4\beta 7$ ligands, such as, but not limited to VCAM-1 and MAdCAM-1, and to antibodies directed against $\alpha 4\beta 7$ integrin, such as, but not limited to, vedolizumab, natalizumab, and etrolizumab. The

soluble $\alpha 4$ integrin polypeptide and the soluble $\beta 7$ integrin polypeptide can heterodimerize via a cysteine bridge or a derivative thereof (*see*, Figure 13). In some embodiments, the binding pair is selected from the group consisting of coiled-coil peptides, leucine zipper peptides, dock-and-lock peptides, avidin-biotin, and derivatives thereof. In some instances, the coiled-coil peptides are ACID-BASE peptides.

[0116] The polypeptides described herein can be produced by any method known to one of ordinary skill in the art, such as but not limited to, synthetic methods, such as solid phase and liquid phase synthesis, or recombinant biology methods, such as those described herein.

[0117] In other aspects, the present invention provides an expression vector encoding a soluble $\alpha 4\beta 7$ integrin polypeptide comprising a first polynucleotide sequence comprising a nucleic acid sequence encoding a soluble $\alpha 4$ integrin polypeptide having an amino acid sequence that has at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:1 and a nucleic acid sequence encoding a first member of a binding pair, and a second polynucleotide sequence comprising a nucleic acid sequence encoding a soluble $\beta 7$ integrin polypeptide having an amino acid sequence that has at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:2 and a nucleic acid sequence encoding a second member of the binding pair. In some embodiments, the first member of the binding pair is an ACID peptide. In some embodiments, the first member of the binding pair is a BASE peptide. In some instances, the second polynucleotide sequence further comprises a nucleic acid sequence encoding an affinity tag, such as a histidine tag. In some instances, the second polynucleotide sequence further comprises a nucleic acid sequence encoding a protease cleavage site, such as a TEV site.

[0118] In some embodiments, the first polynucleotide sequence comprises a nucleic acid sequence encoding a polypeptide having an amino acid sequence that has at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:3. In some embodiments, the second polynucleotide sequence comprises a nucleic acid sequence encoding a polypeptide having an amino acid sequence that has at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%,

90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:4.

[0119] In some embodiments, the expression vector is capable of directing expression of the polynucleotide sequences preferentially in a particular cell type. The expression vector can be a plasmid, phage, phagemid, cosmid, bacteriophage, baculovirus vector, lentiviral vector, retroviral vector, adenoviral vector, yeast plasmid, and the like. The expression vector can also comprise a promoter. Useful promoters include constitutive promoters and inducible promoters. The first polynucleotide sequence and/or the second polynucleotide sequence of the expression vector may be operably linked to a promoter. The promoter can be selected depending on the host cell containing the expression vector or used to generate or produce the soluble $\alpha 4\beta 7$ integrin polypeptide encoded by the expression vector described herein. The expression vector may include regulatory elements, a selectable marker cassette, antibiotic resistance cassette, or any other component that facilitates the expression of the polypeptide by a host cell.

[0120] In some embodiments, the first and second polynucleotide sequences are found in a single expression vector. Such polynucleotide sequence can be located in a bicistronic expression vector such that an IRES sequence is located between the first and second polynucleotides sequences in the vector. A single promoter can drive the expression of both polynucleotides sequences. In some embodiments, the first polynucleotide sequence is operably linked to the promoter and is located immediately upstream from a nucleic acid sequence encoding a ribosomal skip, such as a viral 2A peptide, which is immediately upstream of the second polynucleotide sequence. In other embodiments, the second polynucleotide sequence is operably linked to the promoter and is located immediately upstream from a nucleic acid sequence encoding a ribosomal skip which is immediately upstream of the first polynucleotide sequence. The soluble $\alpha 4$ integrin polypeptide and soluble $\beta 7$ integrin polypeptide can be generated from one expression vector.

[0121] Methods for constructing an expression vector are known to those of ordinary skill in the art. Detailed descriptions of protocols and methods are described in, *e.g.*, Green, M.R., and Sambrook, J., eds., *Molecular Cloning: A Laboratory Manual*, 4th ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y. (2012), Ausubel, F. M., *et al.*, *Current Protocols in Molecular Biology* (Supplement 99), John Wiley & Sons, New York (2012),; Berger and Kimmel, *Guide to Molecular Cloning Techniques*, *Methods in Enzymology*:

Volume 152, Academic Press, Inc., San Diego, CA, (1987); and PCR Protocols: A Guide to Methods and Applications, Academic Press, San Diego, CA, (1990).

[0122] In other aspects, the present invention provides a host cell comprising any one of the expression vectors encoding a soluble $\alpha 4\beta 7$ integrin polypeptide described herein. The host cell can be a stable cell line, such as, but not limited to, HEK293 cells, CHO cells, COS cells, Jurkat cells, NIH3T3 cells, and derivatives thereof. The host cell can be a bacterial cell, yeast cell, fungal cell, algal cell, plant cell, insect cell, animal cell, mammalian cell, non-human cell, or human cell. Suitable host cells are described in Goeddel, Gene Expression Technology: Methods in Enzymology, 185, Academic Press, San Diego, CA, (1990).

[0123] The expression vector can be introduced into the host cell by methods including, but not limited to, transformation, transfection, lipofection, nucleofection, microinjection, electroporation, biolistics, virosomes, liposomes, immunoliposomes, polycation or lipid:polynucleotide conjugates, naked DNA, artificial virions, and agent-enhanced uptake of DNA.

[0124] In yet another aspect, the present invention provides a method of generating a soluble $\alpha 4\beta 7$ integrin polypeptide encoded by any one of the expression vectors described herein. The method comprises (a) introducing the expression vector encoding the soluble $\alpha 4\beta 7$ integrin polypeptide into a host cell, (b) culturing the resulting host cell under conditions to produce the soluble $\alpha 4\beta 7$ integrin polypeptide, and (c) isolating the soluble $\alpha 4\beta 7$ integrin polypeptide.

[0125] The cells containing the expression vector can be cultured under conditions that allow, promote or induce the production of the soluble $\alpha 4\beta 7$ integrin polypeptide.

[0126] The soluble $\alpha 4$ integrin polypeptide, soluble $\beta 7$ integrin polypeptide and soluble $\alpha 4\beta 7$ integrin polypeptide can be purified from, for example, a cell culture supernatant or soluble fraction of a cell extract, according to standard methods known in the art including, but not limited to, chromatography (*e.g.*, ion exchange, affinity, hydrophobic, chromatofocusing, and size exclusion), electrophoretic procedures (*e.g.*, preparative isoelectric focusing), differential solubility (*e.g.*, ammonium sulfate precipitation), SDS-PAGE, or extraction (*see, e.g.*, Protein Purification, J. C. Janson and Lars Ryden, eds., VCH Publishers, New York, (1989)) to obtain substantially pure polypeptides. Methods for protein purification, chromatography, electrophoresis, centrifugation, and crystallization are

described in, *e.g.*, Coligan *et al.*, Current Protocols in Protein Science, Vol. 1, John Wiley and Sons, Inc., New York, (2000).

[0127] The soluble recombinant $\alpha 4\beta 7$ integrin polypeptide can complex with its cognate ligand, such as a ligand that specifically binds to wild-type, full-length $\alpha 4\beta 7$ integrin. The soluble recombinant $\alpha 4\beta 7$ integrin polypeptide can be an antigen for an anti- $\alpha 4\beta 7$ integrin antibody.

VI. Soluble IL-12p40 Polypeptide Antigens

[0128] In one aspect, the present invention provides an isolated soluble IL-12p40 polypeptide comprising an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NOS:6, 11, 12, or 13. In some embodiments, the polypeptide further comprises an affinity tag. In other embodiments, the isolated soluble IL-12p40 polypeptide comprising an amino acid sequence having at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:7. In some embodiments, the polypeptide further comprises an affinity tag. In particular embodiments, the soluble IL-12p40 polypeptide is a monomer, and cannot dimerize or form multimers.

[0129] In another aspect, the present invention provides an expression vector encoding a soluble IL-12p40 polypeptide comprising a polynucleotide sequence comprising a nucleic acid sequence encoding an IL-12p40 polypeptide having an amino acid sequence that has at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NOS:6, 11, 12, or 13. The polynucleotide sequence can further comprise a nucleic acid sequence encoding an affinity tag. Such an affinity tag can be a histidine tag, such as hexahistidine. Other non-limiting examples of affinity tags include an avidin tag, V5 tag, FLAG tag, HA tag, Myc tag, cleavable tag, and the like. In some embodiments, the polynucleotide comprising a nucleic acid sequence encoding an IL-12p40 polypeptide having an amino acid sequence that has at least 80% sequence identity, *e.g.*, at least 80%, 81%, 82%, 83%, 84%, 85%, 86%, 87%, 88%, 89%, 90%, 91%, 92%, 93%, 94%, 95%, 96%, 97%, 98%, 99% or 100% sequence identity, to SEQ ID NO:7.

[0130] In some embodiments, the expression vector is capable of directing expression of the polynucleotide sequences preferentially in a particular cell type. The expression vector can be a plasmid, phage, phagemid, cosmid, bacteriophage, baculovirus vector, lentiviral vector, retroviral vector, adenoviral vector, yeast plasmid, and the like. The expression vector can also comprise a promoter. Useful promoters include constitutive promoters and inducible promoters. The polynucleotide sequence of the expression vector may be operably linked to a promoter. The promoter can be selected depending on the host cell selected to generate or produce the soluble IL-12p40 polypeptide encoded by the expression vector described herein. The expression vector may include regulatory elements, a selectable marker cassette, antibiotic resistance cassette, or any other component that facilitates the expression of the polypeptide.

[0131] Methods for constructing an expression vector are known to those of ordinary skill in the art. Detailed descriptions of protocols and methods are described in, *e.g.*, Green, M.R., and Sambrook, J., eds., *Molecular Cloning: A Laboratory Manual*, 4th ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y. (2012); Ausubel, F. M., *et al.*, *Current Protocols in Molecular Biology* (Supplement 99), John Wiley & Sons, New York (2012); Berger and Kimmel, *Guide to Molecular Cloning Techniques*, *Methods in Enzymology*, Volume 152, Academic Press, Inc., San Diego, CA (1987); and *PCR Protocols: A Guide to Methods and Applications*, Academic Press, San Diego, CA, (1990).

[0132] In other aspects, the present invention provides a host cell comprising any one of the expression vectors encoding a soluble IL-12p40 polypeptide described herein. The host cell can be a stable cell line, such as, but not limited to, HEK293 cells, CHO cells, COS cells, Jurkat cells, NIH3T3 cells, and derivatives thereof. The host cell can be a bacterial cell, yeast cell, fungal cell, algal cell, plant cell, insect cell, animal cell, mammalian cell, non-human cell, or human cell. Suitable host cells are described in Goeddel, *Gene Expression Technology: Methods in Enzymology* 185, Academic Press, San Diego, CA, 1990.

[0133] The expression vector can be introduced into the host cell by methods including, but not limited to, transformation, transfection, lipofection, nucleofection, microinjection, electroporation, biolistics, virosomes, liposomes, immunoliposomes, polycation or lipid:polynucleotide conjugates, naked DNA, artificial virions, and agent-enhanced uptake of DNA.

[0134] In yet another aspect, the present invention provides method of generating a soluble IL-12p40 polypeptide encoded by any one of the expression vectors described herein. The method comprises (a) introducing the expression vector encoding the soluble IL-12p40 polypeptide into a host cell, (b) culturing the resulting host cell under conditions to produce the soluble IL-12p40 polypeptide, and (c) isolating the soluble IL-12p40 polypeptide.

[0135] The cells containing the expression vector can be cultured under conditions that allow, promote or induce the production of the soluble IL-12p40 polypeptide.

[0136] The soluble IL-12p40 polypeptide can be purified from, for example, a cell culture supernatant or soluble fraction of a cell extract, according to standard methods known in the art including, but not limited to, chromatography (*e.g.*, ion exchange, affinity, hydrophobic, chromatofocusing, and size exclusion), electrophoretic procedures (*e.g.*, preparative isoelectric focusing), differential solubility (*e.g.*, ammonium sulfate precipitation), SDS-PAGE, or extraction (*see, e.g.*, Protein Purification, J. C. Janson and Lars Ryden, eds., VCH Publishers, New York, (1989)) to obtain substantially pure polypeptides. Methods for protein purification, chromatography, electrophoresis, centrifugation, and crystallization are described in, *e.g.*, Coligan *et al.*, Current Protocols in Protein Science, Vol. 1, John Wiley and Sons, Inc., New York, (2000).

[0137] Unlike wild-type IL-12p40, the soluble recombinant IL-12p40 polypeptide described herein cannot form a dimer, trimer, or oligomer. In particular embodiments, the presence of either two cysteine to alanine substitutions (SEQ ID NO:6) or two cysteine to serine substitutions (SEQ ID NO:11) in the wild-type IL-12p40 polypeptide sequence prevents the IL-12p40 antigen from oligomerizing. In other embodiments, the soluble IL-12p40 polypeptide monomer has the wild-type IL-12p40 polypeptide sequence with one cysteine to alanine substitution and one cysteine to alanine substitution (SEQ ID NOS:12 or 13). An anti-IL-12p40 antibody, such as ustekinumab, can specifically bind to the soluble recombinant IL-12p40 polypeptide.

VII. Biologic Therapy

[0138] The indirect homogeneous mobility shift assays of the present invention are suitable for detecting and/or measuring the presence or level of a biologic in a sample from a subject (*e.g.*, a subject receiving biologic therapy). Non-limiting examples of biologics include antibodies, antibody fragments, proteins, polypeptides, peptides, fusion proteins (*e.g.*, Ig fusion proteins or Fc fusion proteins), multivalent binding proteins (*e.g.*, DVD Ig), antibody-

drug conjugates, vaccines, nucleic acids, sugars, recombinant forms thereof, engineered forms thereof, and combinations thereof.

[0139] Examples of antibody-based biologics include, but are not limited to, diagnostic or therapeutic monoclonal antibodies and antigen-binding fragments or conjugates thereof. In certain embodiments, the antibody comprises an anti-integrin drug such as an anti- $\alpha 4\beta 7$ integrin drug (*e.g.*, vedolizumab (ENTYVIO[™]), etrolizumab) and/or an anti- $\alpha 4\beta 1$ integrin drug (*e.g.*, natalizumab (TYSABRI[®])). In other embodiments, the antibody comprises an anti-cytokine drug such as an anti-IL12p40 drug (*e.g.*, ustekinumab (STELARA[®])). In yet other embodiments, the antibody comprises an anti-cytokine receptor drug such as an anti-IL-6 receptor drug (*e.g.*, tocilizumab (ACTEMRA[®])). In further embodiments, the antibody comprises an anti-CD receptor drug such as an anti-CD3 receptor drug (*e.g.*, visilizumab), an anti-CD4 receptor drug (*e.g.*, priliximab), an anti-CD20 receptor drug (*e.g.*, rituximab (RITUXAN[®]), ofatumumab (ARZERRA[®]), obinutuzumab (GAZYVA[®]), ibritumomab tiuxetan (ZEVALIN[®]), tositumomab (BEXXAR[®]), ocrelizumab, veltuzumab, an anti-CD25 receptor drug (*e.g.*, daclizumab (ZENAPAX[®])), or combinations thereof. In other embodiments, the antibody comprises an anti-TNF α drug such as infliximab (REMICADE[®]), adalimumab (HUMIRA[®]), etanercept (ENBREL[®]), golimumab (SIMPONI[®]), certolizumab pegol (CIMZIA[®]), or combinations thereof. Additional examples of antibody-based biologics include antibody-drug conjugates such as brentuximab vedotin (ADCETRIS[®]).

[0140] Table 1 provides an exemplary and non-exhaustive list of diagnostic and therapeutic monoclonal antibodies which have either been approved or are currently in development. An extensive list of biologic medicines including monoclonal antibody-based therapeutics and diagnostics in clinical development and approved products is provided in the 2006 PhRMA Report entitled “418 Biotechnology Medicines in Testing Promise to Bolster the Arsenal Against Disease” and the 2013 PhRMA Report entitled “Medicines in Development – Biologics,” the disclosures of which are hereby incorporated by reference in their entirety for all purposes.

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
Digestive Disorders		
ABT 874	Abbott Laboratories	Crohn's disease

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
AMG 139/MEDI-2070	Amgen	Crohn's disease
AMG 181/MEDI-7183	Amgen	Crohn's disease, ulcerative colitis
anrukinzumab (IMA-638)	Pfizer	ulcerative colitis
anti-IP10	Bristol-Myers Squibb	Crohn's disease, ulcerative colitis
clazakizumab (anti-IL6)	Bristol-Myers Squibb, Alder Biopharmaceuticals	Crohn's disease
etrolizumab (rhuMAb- β 7)	Genentech	ulcerative colitis
GSK1070806 (IL-18 mAb)	GlaxoSmithKline	inflammatory bowel disease
Humira [®] (adalimumab)	AbbVie	Crohn's disease
MDX-1100	Millennium Pharmaceuticals	ulcerative colitis
Nuvion [®] (visilizumab)	PDL BioPharma	I.V. steroid-refractory ulcerative colitis and Crohn's disease
PF-00547659	Pfizer	Crohn's disease
PF-04236921	Pfizer	Crohn's disease
QAX576	Novartis Pharmaceuticals	Crohn's disease
Remicade [®] (infliximab)	Janssen Biotech	Crohn's disease
SAR252067 (anti-LIGHT mAb)	Sanofi US	Crohn's disease, ulcerative colitis
SAR339658 (VLA2 antagonist)	Sanofi US	inflammatory bowel disease
Simponi [®] (golimumab)	Janssen Biotech	ulcerative colitis
Stelara [®] (ustekinumab)	Janssen Biotech	Crohn's disease
tralokinumab	AstraZeneca, MedImmune	ulcerative colitis
Tysabri [®] (natalizumab)	Biogen Idec	Crohn's disease
vedolizumab (MLN0002)	Takeda Pharmaceuticals	Crohn's disease, ulcerative colitis
Autoimmune disorders		
ABT-122	AbbVie	rheumatoid arthritis
Actemra [®] (tocilizumab)	Genentech, Roche	early rheumatoid arthritis, systemic sclerosis
AGS-009	Argos Therapeutics	systemic lupus erythematosus (SLE)
alemtuzumab	Genzyme	relapsing-remitting multiple sclerosis
AME 527	Applied Molecular	rheumatoid arthritis
AMG 108	Amgen	rheumatoid arthritis
AMG 557/MEDI-5872	Amgen, AstraZeneca, MedImmune	SLE
AMG 714	Amgen	rheumatoid arthritis
AMG 729	Amgen	autoimmune diseases
AMG 811	Amgen	discoid lupus erythematosus, SLE
ART 874	Abbott Laboratories	multiple sclerosis
anti-CD16 mAb	MacroGenics	immune thrombocytopenic
anti-IL17 mAb (RG7624)	Genentech	autoimmune disorders
anti-LINGO (BIIB033)	Biogen Idec	Multiple sclerosis
Benlysta [®] (belimumab)	GlaxoSmithKline	rheumatoid arthritis, SLE, systemic sclerosis, scleroderma
BI-695500 (rituximab biosimilar)	Boehringer-Ingelheim Pharmaceuticals	rheumatoid arthritis
BI-695501 (adalimumab biosimilar)	Boehringer-Ingelheim Pharmaceuticals	rheumatoid arthritis
BT-061	AbbVie, Biotest	rheumatoid arthritis
Cimzia [®] (certolizumab pegol)	UCB	ankylosing spondylitis, juvenile rheumatoid arthritis
clazakizumab (anti-IL6)	Bristol-Myers Squibb, Alder	rheumatoid arthritis

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
	Biopharmaceuticals	
CNTO-136 (sirukumab)	Janssen Biotech	rheumatoid arthritis
CNTO-1959	Janssen Biotech	rheumatoid arthritis
daclizumab (anti-CD25 mAb)	AbbVie, Biogen Idec	multiple sclerosis
epratuzumab	Immunomedics, UCB	SLE
ETI-201	Elusys Therapeutics	SLE
GSK1223249 (NOGO-A mAb)	GlaxoSmithKline	multiple sclerosis
Humira [®] (adalimumab)	AbbVie	rheumatoid arthritis, ankylosing spondylitis, juvenile rheumatoid arthritis, psoriasis
HuZAF [®] (fontolizumab)	PDL BioPharma, Biogen Idec	rheumatoid arthritis
Ilaris [®] (canakinumab)	Novartis Pharmaceuticals	systemic juvenile idiopathic arthritis
IMMU-106 (hCD20)	Immunomedics	autoimmune diseases
mavrilimumab	AstraZeneca, MedImmune	rheumatoid arthritis
MEDI-545 (MDX-1103)	Medarex, MedImmune	lupus
MEDI-546 (anti-IFN-alphaR mAb)	AstraZeneca, MedImmune	SLE
MEDI-551 (anti-CD19 mAb)	AstraZeneca, MedImmune	scleroderma
MEDI-570 (anti-ICOS mAb)	AstraZeneca, MedImmune	SLE
MLN 1202	Millennium Pharmaceuticals	multiple sclerosis
NN8209 (anti-C5aR-151 mAb)	Novo Nordisk	rheumatoid arthritis
NN8210 (anti-C5aR-215 mAb)	Novo Nordisk	rheumatoid arthritis
NN8226 (anti-IL-20 mAb)	Novo Nordisk	rheumatoid arthritis
NN8765 (anti-NKG2 mAb)	Novo Nordisk	rheumatoid arthritis
NN8828 (anti-IL-21 mAb)	Novo Nordisk	rheumatoid arthritis
ocrelizumab (anti-CD20 mAb)	Biogen Idec, Genentech, Roche	multiple sclerosis, rheumatoid arthritis
ofatumumab	GlaxoSmithKline	multiple sclerosis, rheumatoid arthritis
OKT3-gamma-1	Johnson & Johnson	psoriatic arthritis
olokizumab	UCB	rheumatoid arthritis
otelixizumab (anti-CD3 mAb)	GlaxoSmithKline	rheumatoid arthritis
ozoralizumab (ATN-103)	Ablynx	rheumatoid arthritis
pateclizumab (anti-LT alpha mAb)	Genentech	rheumatoid arthritis
PD-360324	Pfizer	cutaneous lupus erythematosus
PF-04236921	Pfizer	SLE, rheumatoid arthritis
PF-05280586 (rituximab biosimilar)	Pfizer	rheumatoid arthritis
Prolia [®] (denosumab)	Amgen	rheumatoid arthritis
Remicade [®] (infliximab)	Janssen Biotech	rheumatoid arthritis
Rituxan [®] (rituximab)	Genentech, Biogen Idec	rheumatoid arthritis, lupus, primary progressive multiple sclerosis, SLE, relapsing-remitting multiple sclerosis
rontalizumab (RG7415)	Genentech	SLE
SAN-300 (anti-VLA-1 antibody)	Santarus	rheumatoid arthritis
SAR113244 (anti-CXCR5 mAb)	Sanofi US	SLE
sarilumab (SAR153191)	Regeneron Pharmaceuticals, Sanofi US	rheumatoid arthritis
secukinumab (AIN457)	Novartis Pharmaceuticals	ankylosing spondylitis, rheumatoid arthritis, multiple sclerosis
sifalimumab (anti-IFN-alpha mAb)	AstraZeneca, MedImmune	SLE

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
Simponi [®] (golimumab)	Janssen Biotech	rheumatoid arthritis, juvenile rheumatoid arthritis, sarcoidosis, ankylosing spondylitis, psoriatic arthritis
siplizumab (MEDI-507)	MedImmune	psoriasis
Soliris [®] (eculizumab)	Alexion Pharmaceuticals	severe or refractory myasthenia gravis
Stelara [®] (ustekinumab)	Janssen Biotech	rheumatoid arthritis, sarcoidosis, plaque psoriasis, multiple sclerosis
tabalumab (BAFF inhibitor)	Eli Lilly	SLE
TRX 1 (anti-CD4)	TolerRx	cutaneous lupus erythematosus
TRX 4	TolerRx	psoriasis
Tysabri [®] (natalizumab)	Biogen Idec	Multiple sclerosis
veltuzumab (IMMU-106)	Immunomedics, Takeda Pharmaceuticals USA	immune thrombocytopenic purpura, rheumatoid arthritis
VX15	Teva Pharmaceuticals, Vaccinex	multiple sclerosis
XmAb [®] 5871 (anti-CD19 mAb)	Xencor	autoimmune disorders
Musculoskeletal Disorders/Arthritis		
ABT-981	AbbVie	osteoarthritis
AMG 167	Amgen	metabolic bone diseases
AMG 745	Amgen	muscular atrophy
AMG 827 (brodalumab)	AstraZeneca, Amgen	psoriatic arthritis
blosozumab (LY2541546)	Eli Lilly	osteoporosis
BYM338	Novartis Pharmaceuticals	sporadic inclusion body myositis, muscular atrophy
Cimzia [®] (certolizumab pegol)	UCB	psoriatic arthritis
clazakizumab (anti-IL6)	Bristol-Myers Squibb, Alder Biopharmaceuticals	psoriatic arthritis
gevokizumab (IL-1B inhibitor mAb)	XOMA	osteoarthritis of the hand
Humira [®] (adalimumab)	AbbVie	spondylarthritis
Ilaris [®] (canakinumab)	Novartis Pharmaceuticals	gouty arthritis
ixekizumab (IL-17 antibody)	Eli Lilly	psoriatic arthritis
LY2495655 (anti-myostatin mAb)	Eli Lilly	disuse muscular atrophy
MCS110	Novartis Pharmaceuticals	synovitis
Prolia [®] (denosumab)	Amgen	male osteoporosis, postmenopausal osteoporosis
romosozumab (AMG 785)	Amgen	postmenopausal osteoporosis
SAR391786 (REGN1033)	Regeneron Pharmaceuticals, Sanofi US	treatment of muscle atrophy post-orthopedic surgery
secukinumab (AIN 457)	Novartis Pharmaceuticals	psoriatic arthritis, polymyalgia rheumatica
Stelara [®] (ustekinumab)	Janssen Biotech	psoriatic arthritis
tanezumab	Pfizer	osteoarthritis
Cancer and Related Conditions		
1311-huA33	Life Science Pharmaceuticals	colorectal cancer
1D09C3	GPC Biotech	relapsed/refractory B-cell lymphomas
8H9 mAb	United Therapeutics	metastatic brain cancer
212-Pb-TCMC-trastuzumab	AREVA Med	HER2-positive cancer metastasized to the abdominal region

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
AbGn-7	AbGenomics International	solid tumors
ABT-806	AbbVie	solid tumors
Actimab-A (M195 mAb)	Actinium Pharmaceuticals	acute myeloid leukemia (AML)
Adcetris [®] (brentuximab vedotin)	Seattle Genetics	cutaneous T-cell lymphoma, front-line Hodgkin lymphoma, post-transplant Hodgkin lymphoma relapse prevention, non-Hodgkin lymphoma, non-lymphoma malignancies, CD30-positive hematologic malignancies
AGS PSCA mAb	Agensys, Merck	prostate cancer
ALT-836	Altor BioScience, Genentech	cancer
AME-133v	MENTRIK Biotech	non-Hodgkin lymphoma
AMG 102	Amgen	cancer
AMG 479	Amgen	cancer
AMG 623	Amgen	B-cell chronic lymphocytic leukemia (CLL)
AMG 655	Amgen	cancer
AMG 706	Amgen	imatinib resistant GIST, advanced thyroid cancer
AMG 780	Amgen	solid tumors
AMG 820	Amgen	solid tumors
AMG 888 (U3-1287)	Amgen	non-small-cell lung cancer (NSCLC)
antibody-drug conjugate (RG7600)	Genentech	ovarian cancer, pancreatic cancer
anti-CD22 ADC (RG7593)	Genentech, Seattle Genetics	diffuse large B-cell lymphoma, non-Hodgkin lymphoma, chronic lymphocytic leukemia (CLL)
anti-CD23 MAb	Biogen Idec	CLL
anti-CD45 mAb	Actinium Pharmaceuticals	AML
anti-CD79b ADC (RG7596)	Genentech, Seattle Genetics	diffuse large B-cell lymphoma,, non-Hodgkin lymphoma, CLL
anti-CD80 MAb	Biogen Idec	non-Hodgkin B-cell lymphoma
anti-CXCR4	Bristol-Myers Squibb	hematological malignancies
anti-EGFL7 mAb (RG7414)	Genentech	metastatic colorectal cancer, NSCLC, solid tumors
anti-FGFR3 mAb (RG7444)	Genentech	solid tumors
anti-HER3/EGFR DAF mAb (RG7597)	Genentech	colorectal cancer, head and neck cancer
anti-idiotypic cancer vaccine	Viventia Biotech	malignant melanoma
anti-lymphotoxin beta receptor mAb	Biogen Idec	solid tumors
anti-PD-L1	Bristol-Myers Squibb	cancer
anti-PD-L1 mAb (RG7446)	Genentech	melanoma, solid tumors
anti-PEM MAb	Somanta Pharmaceuticals	cancer
anti-STEAP1 ADC (RG7450)	Genentech, Seattle Genetics	prostate cancer
anti-Tac(Fv)-E38 immunotoxin	National Cancer Institute	leukemia, lymphoma
APN301 (hu14.18-IL2)	Apeiron Biologics	malignant melanoma, neuroblastoma in children
Archexin [™] (RX-0201)	Rexahn Pharmaceuticals	pancreatic cancer
Arzerra [®] (ofatumumab)	GlaxoSmithKline	CLL, diffuse large B-cell lymphoma, follicular lymphoma

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
ASG-5ME	Agensys, Seattle Genetics	pancreatic cancer, castration-resistant prostate cancer
ASG-22ME	Agensys, Seattle Genetics	solid tumors
AV-203	AVEO Oncology	solid tumors
Avastin [®] (bevacizumab)	Genentech, Roche	ovarian cancer, HER-2 negative-breast cancer, HER-2 positive breast cancer, high-risk carcinoid tumors, glioblastoma multiforme, metastatic ovarian cancer, NSCLC, metastatic colorectal cancer
AVE 9633 maytansin-loaded anti-CD33 mAb	Sanofi Aventis	AML
bavituximab	Peregrine Pharmaceuticals	NSCLC, pancreatic cancer, breast cancer, liver cancer, prostate cancer, rectal adenocarcinoma
BAX-69	Baxter International	solid tumors
BAY 79-4620	Bayer HealthCare Pharmaceuticals	solid tumors
BAY 94-9343	Bayer HealthCare Pharmaceuticals	solid tumors
BAY 20-10112	Amgen, Bayer HealthCare Pharmaceuticals	solid tumors
Bexxar [®] (tositumomab)	GlaxoSmithKline	non-Hodgkin lymphoma
BHQ880	Novartis Pharmaceuticals	multiple myeloma
BI-505	BioInvent International	multiple myeloma
BI-836845	Boehringer Ingelheim Pharmaceuticals	solid tumors
bivatuzumab	Boehringer Ingelheim Pharmaceuticals	cancer
blinatumomab	Amgen	leukemia and lymphoma
BrevaRex [™]	ViRexx	breast cancer, multiple myeloma
BT-062 (indatuximab ravtansine)	Biotest	multiple myeloma
BYM338	Novartis Pharmaceuticals	cancer-related cachexia
Campath [®] (alemtuzumab)	Berlex Laboratories, Genzyme	B-cell chronic lymphocytic leukemia, lymphoma
catumaxomab	Fresenius Biotech	malignant ascites, ovarian cancer
CAT 3888	Cambridge Antibody Technology	hairy cell leukemia
CDX-011 (glembatumumab vedotin)	Celldex Therapeutics	breast cancer, malignant melanoma
CDX-1127	Celldex Therapeutics	hematological malignancies, solid tumors
CEP-37250/KNK-2804	Teva North America, Kyowa Hakko Kirin Pharma	adenocarcinoma
ch14.18 mAb	United Therapeutics	neuroblastoma
chimeric mAb	National Cancer Institute	neuroblastoma
Cixutumumab (LY3012217)	Eli Lilly, Imclone Systems	NSCLC
CNTO-328 (siltuximab)	Janssen Biotech	giant lymph node hyperplasia, multiple myeloma, myelodysplastic syndromes, prostate cancer, renal cancer,
Cotara [™] mAb TNT	Peregrine Pharmaceuticals	recurrent glioblastoma
CP-751871 (figitumumab)	Pfizer	adrenocortical carcinoma, non-small cell lung cancer
CS-1008 (tigatuzumab)	Daiichi Sankyo	pancreatic cancer, colorectal

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
		cancer, non-small cell lung cancer, ovarian cancer
CSF-1R mAb (IMC-CS4)	Eli Lilly	solid tumors
CT-011 (pidilizumab)	CureTech	AML, colorectal cancer, diffuse large B-cell lymphoma, follicular lymphoma, malignant melanoma
dalotuzumab (MK-0646)	Merck	breast cancer, neuroendocrine tumors, NSCLC
daratumumab	Janssen Biotech	multiple myeloma
DEDN6526A	Genentech	malignant melanoma
demcizumab (OMP-21M18)	GlaxoSmithKline, OncoMed Pharmaceuticals	solid tumors
DFRF4539A	Genentech	multiple myeloma
DI17E6 (anti-integrin mAb)	EMD Serono	colorectal cancer, prostate cancer
DKN-01	Dekkun	solid tumors
ecromeximab (KW-2871)	Life Science Pharmaceuticals	metastatic melanoma
elotuzumab	Bristol-Myers Squibb, AbbVie	multiple myeloma
EMD 273063	EMD Lexigen	solid tumors, malignant melanoma, neuroblastoma, SCLC
enavatuzumab	AbbVie	solid tumors
ensituximab (NPC-1C)	Neogenix Oncology	colorectal cancer, pancreatic cancer
epratuzumab Y-90/veltuzumab combination	Immunomedics	non-Hodgkin lymphoma
Erbix [®] (cetuximab)	Bristol-Myers Squibb, Eli Lilly, ImClone Systems	esophageal cancer, colorectal cancer, squamous cell cancer of the head and neck
farletuzumab (MORAb-003)	Eisai	platinum-sensitive ovarian cancer, NSCLC
FG-3019	FibroGen	pancreatic cancer
ficlatuzumab	AVEO Oncology	glioblastoma, lymphoma, multiple myeloma, solid tumors
flanvotumab (TYRP1 protein)	Eli Lilly	malignant melanoma
Fzd7 (vantictumab)	Bayer HealthCare Pharmaceuticals, OncoMed Pharmaceuticals	solid tumors
ganitumab	Amgen	pancreatic cancer, breast cancer, colorectal cancer, sarcoma
GC-33/RG7686	Chugai Pharma USA, Roche	liver cancer
GMK	Progenies Pharmaceuticals	prevention of recurrence following surgery to remove primary melanoma in high-risk patients
GS-6624 (simtuzumab)	Gilead Sciences	colorectal cancer, pancreatic cancer
HCD122 (anti-CD40 mAb)	Novartis Pharmaceuticals, XOMA	lymphoma
Herceptin [®] (trastuzumab)	Genentech	HER2-overexpressing early stage or metastatic breast cancer
HGS-ETR1 (mapatumumab)	GlaxoSmithKline	liver cancer, multiple myeloma, NSCLC, hematologic and solid tumors
HGS-TR2J	Human Genome Sciences	advanced solid tumors
HuC242-DM4	ImmunoGen	colorectal, gastrointestinal, NSCLC, pancreatic cancers
HuL2G7	Galaxy Biotech	solid tumors
HuMax-CD4 (zanolimimumab)	Genmab, Serono	cutaneous T-cell lymphoma, non-

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
		cutaneous T-cell lymphoma
HuMax CD20 (ofatumumab)	Gemnab	CLL, non-Hodgkin lymphoma
HuMax-EGFr	Genmab	head and neck cancer
HuM195/rGel	Targa Therapeutics	AML, CML, myelodysplastic syndromes
huN901-DM1 (lorvotuzumab mertansine)	ImmunoGen	SCLC, multiple myeloma, solid tumors
icrucumab (LY3012212)	Eli Lilly, ImClone Systems	bladder cancer, breast cancer, colorectal cancer
IMC-TR1 (LY3022859)	Eli Lilly, ImClone Systems	solid tumors
IMGN529	ImmunoGen	non-Hodgkin lymphoma
IMGN853	ImmunoGen	solid tumors
IMMU-102 (epratuzumab Y-90)	Immunomedics	non-Hodgkin lymphoma
inotuzumab ozogamicin (CMC-544)	Pfizer, UCB	aggressive non-Hodgkin lymphoma, ALL
Iomab-B (anti-CD45 mAb)	Actinium Pharmaceuticals	AML
J 591 Lu-177	BZL Biologics	prostate cancer
KB004	KaloBios Pharmaceuticals	hematological malignancies
LFA102	Novartis Pharmaceuticals, XOMA	breast cancer, prostate cancer
lirilumab (anti-KIR)	Bristol-Myers Squibb	cancer
LY2495655 (anti-myostatin mAb)	Eli Lilly	cancer cachexia
LY2875358 (c-met mAb)	Eli Lilly	cancer
M195-bismuth 213 conjugate	Actinium Pharmaceuticals	AML
M200 (volociximab)	PDL BioPharma, Biogen Idec	advanced solid tumors
MAb HeFi-1	National Cancer Institute	lymphoma, non-Hodgkin lymphoma
MABp1	XBioTech	cancer-related cachexia, advanced cancer, leukemia
MDX-060 (iratumumab)	Medarex	Hodgkin's disease, anaplastic large-cell-lymphoma
MDX-070	Medarex	prostate cancer
MDX-214	Medarex	EGFR-expressing cancers
MEDI-522	MedImmune	T-cell lymphoma, melanoma, prostate cancer, solid tumors
MEDI-551 (anti-CD19 mAb)	AstraZeneca, MedImmune	hematological malignancies
MEDI-573 (anti-IGF mAb)	AstraZeneca, MedImmune	solid tumors
MEDI-575 (anti-PDGFR α mAb)	AstraZeneca, MedImmune	glioblastoma, NSCLC
MEDI-0639 (anti-DLL-4 mAb)	AstraZeneca, MedImmune	solid tumors
MEDI-3617 (anti-ANG-2 mAb)	AstraZeneca, MedImmune	solid tumors
MEDI-4736 (anti-CD274 mAb)	AstraZeneca, MedImmune	cancer
MEDI-6469 (anti-OX40 mAb)	AgonOx, AstraZeneca, MedImmune	solid tumors
MGA271 (anti-B7-H3)	MacroGenics	solid tumors
MGAH22 (anti-HER2)	MacroGenics	solid tumors
milatuzumab	Immunomedics	CLL
milatuzumab-DOX	Immunomedics	multiple myeloma
MINT1526A	Genentech	solid tumors
MK-3475	Merck	malignant melanoma, NSCLC
MLN0264 (GCC antibody drug conjugate)	Millennium Pharmaceuticals	gastrointestinal cancer
mogamulizumab	Kyowa Hakko Kirin Pharma	cutaneous T-cell lymphoma, adult T-cell lymphoma, T-cell leukemia

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
MORAb 003	Eisai	ovarian cancer
MORAb-004	Eisai	colorectal cancer, melanoma, sarcoma
MORAb-009 (amatuximab)	Eisai	mesothelioma
moxetumomab pasudotox	AstraZeneca, MedImmune	hematological malignancies
Mylotarg [™] (gemtuzumab ozogamicin)	Wyeth	acute myeloid leukemia
necitumumab	Bristol-Myers Squibb, Eli Lilly, ImClone Systems	NSCLC
neuradiab	Bradmer Pharmaceuticals	glioblastoma
nimotuzumab	InnoMab PTE	glioma, squamous cell carcinomas of the head and neck, recurrent or refractory high grade malignant glioma, anaplastic astrocytomas, glioblastomas and diffuse intrinsic pontine glioma
nivolumab (anti-PD1)	Bristol-Myers Squibb	melanoma, NSCLC, renal cell carcinoma, solid tumors
obinutuzumab (GA101)	Biogen Idec, Genentech	CLL, diffuse large B-cell lymphoma, non-Hodgkin lymphoma
olaratumab (LY3012207)	Eli Lilly, ImClone Systems	glioblastoma
OMP-52M51 (anti-Notch 1)	GlaxoSmithKline, OncoMed Pharmaceuticals	hematological malignancies
OMP-59R5 (anti-Notch 2/3)	GlaxoSmithKline, OncoMed Pharmaceuticals	pancreatic cancer
onartuzumab (anti-c-met- mAb)	Genentech	metastatic NSCLC
oregovomab	Quest Pharmatech	ovarian cancer
PAM 4	Merck	pancreatic cancer
panitumumab (rHuMab EGFr)	Abgenix	colorectal cancer
Perjeta [™] (pertuzumab)	Genentech	early HER2-positive breast cancer, HER2-positive metastatic breast cancer, HER2-positive gastric cancer, ovarian cancer
PF-03446962	Pfizer	solid tumors
PF-04605412	Pfizer	solid tumors
PF-05082566	Pfizer	cancer, lymphoma
PF-05280014 (trastuzumab biosimilar)	Pfizer	metastatic breast cancer
PSMA-ADC	Progenics Pharmaceuticals	prostate cancer
R1550 RadioTheraCIM	Roche, YM BioSciences	metastatic breast cancer, glioma
ramucirumab (LY3009806)	Eli Lilly, ImClone Systems	breast cancer, colorectal cancer, gastric cancer
RAV 12	Raven Biotechnologies	cancer
Redectane [®] (girentuximab I-124)	Wilex AG	diagnosis of kidney cancer
REGN1400	Regeneron Pharmaceuticals	cancer
Rencarex [®] G250	Wilex AG	renal cancer
RG7116	Roche	solid tumors
RG7155	Roche	solid tumors
RG7160 (humAb EGFR)	Roche	colorectal cancer
RG7212	Roche	solid tumors
RG7356 (anti-CD44 mAb)	Roche	AML, solid tumors

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
RG7458 (antibody-drug conjugate)	Genentech, Seattle Genetics	ovarian cancer
RIGScan [™]	Navidea Biopharmaceuticals	diagnosis of colorectal cancer
rilotumumab	Amgen	colorectal cancer, gastric cancer, prostate cancer, SCLC
Rituxan [®] (rituximab)	Genentech	diffuse large B-cell lymphoma, B-cell non-Hodgkin lymphoma, indolent non-Hodgkin lymphoma induction therapy, relapsed or refractory CLL
RON8 mAb	Eli Lilly, ImClone Systems	cancer
SAR3419 (maytansin-loaded anti-CD19 mAb)	Sanofi US	ALL, non-Hodgkin lymphoma
SAR153192 (REGN 421) (anti-DLL4 mAb)	Regeneron Pharmaceuticals, Sanofi US	cancer
SAR256212 (MM-121) (anti-ErbB3 mAb)	Merrimack Pharmaceuticals, Sanofi US	breast cancer, solid tumors
SAR307746 (REGN910) (anti-angiopoietin-2 mAb)	Regeneron Pharmaceuticals, Sanofi US	solid tumors
SAR566658 (maytansin-loaded anti-DS6)	Sanofi US	DS6-positive solid tumors
SAR650984 (anti-CD38 naked mAb)	Sanofi US	hematological malignancies
SGN30	Seattle Genetics	cutaneous anaplastic large-cell lymphoma, systemic anaplastic large-cell lymphoma, Hodgkin's disease
SGN-33 (lintuzumab)	Seattle Genetics	AML, myelodysplastic syndromes CLL multiple myeloma, non Hodgkin lymphoma
SGN-40	Seattle Genetics	AML, myelodysplastic syndromes CLL multiple myeloma, non Hodgkin lymphoma
SGN-75 (vorsetuzumab mafodotin)	Seattle Genetics	non-Hodgkin lymphoma, renal cancer
sibroturtumab	Life Science Pharmaceuticals	colorectal, head and neck, lung cancers
Sym004	EMD Serono, Symphogen	head and neck cancer, solid tumors
tabalumab (BAFF inhibitor)	Eli Lilly	multiple myeloma
Tarvacin [™] (bavituximab)	Peregrine Pharmaceuticals	solid tumors
TF2	Immunomedics	diagnosis of colorectal cancer
TG-1101 (ublituximab)	TG Therapeutics	CLL, non-Hodgkin lymphoma
tigatuzumab	Daiichi Sankyo	breast cancer, liver cancer, ovarian cancer, pancreatic cancer
TNX-650	Tanox	refractory Hodgkin lymphoma
trastuzumab emtansine (T-DM1)	Genentech, Roche	HER2-positive metastatic breast cancer, early HER2-positive breast cancer, advanced HER2-positive gastric cancer
TRC105	TRACON Pharmaceuticals	bladder cancer, liver cancer, ovarian cancer, prostate cancer, solid tumors
tremelimumab (anti-CTLA4 mAb)	AstraZeneca, MedImmune	solid tumors, metastatic melanoma, prostate cancer

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
tumor immunotherapy mAb	Genentech	solid tumors
U3-1565	Daiichi Sankyo	solid tumors
urelumab (anti-CD137)	Bristol-Myers Squibb	cancer
VAY736	MorphoSys, Novartis Pharmaceuticals	CLL
VB4-845	Viventia Biotechnologies	bladder cancer
Vectibix [®] (panitumumab)	Amgen	colorectal cancer
VEGFR3 mAb (IMC-3C5)	Eli Lilly, ImClone Systems	cancer
veltuzumab (IMMU-106)	Immunomedics	CLL, non-Hodgkin lymphoma
VGX-100	Circadian Technologies	solid tumors
volociximab	AbbVie	NSCLC
VX15	Vaccinex	solid tumors
Xgeva [®] (denosumab)	Amgen	delay or prevention of bone, metastases in prostate cancer or breast cancer, giant cell tumor of the bone
XmAb [®] (high ADCC mAb)	Boehringer Ingelheim Pharmaceuticals, Xencor	cancer
XmAb [®] 2513 (anti-CD30 mAb)	Xencor	Hodgkin disease, T-cell lymphoma
XmAb [®] 5574 (anti-CD19 mAb)	MorphoSys, Xencor	CLL
Y-90 hPAM 4 (IMMU-107)	Immunomedics	pancreatic cancer
Yervoy [™] (ipilimumab)	Bristol-Myers Squibb	adjuvant melanoma, NSCLC, prostate cancer, SCLC, gastric cancer, ovarian cancer, leukemia, lymphoma, renal cell cancer
zanolimumab	Emergent BioSolutions	peripheral T-cell lymphoma
Zevalin [®] (ibritumomab tiuxetan)	Spectrum Pharmaceuticals	diffuse large B-cell lymphoma, non-Hodgkin lymphoma
Infectious diseases		
ABthrax [™] (raxibacumab)	Human Genome Sciences	anthrax
Anthim [™] (ETI-204)	Elusys Therapeutics	anthrax
anthrax immune globulin	Cangene	anthrax
anti-HIV-1 mAb	Polymun Scientific	HIV infection
anti-hsp90 mAb	NeuTec Pharma	candidiasis
anti-PD-L1	Bristol-Myers Squibb	hepatitis B
anti-staph mAb	MedImmune	prevention of staphylococcal infections
Aurexis (tefibazumab)	Inhibitex	prevention and treatment of <i>S. aureus</i> bacteremia
bavituximab	Peregrine Pharmaceuticals	hepatitis C
CCR5 MAb	Human Genome Sciences	HIV infection
Cytolin [®] (anti-CD8 mAb)	CytoDyn	HIV infection
FGI-101-1A6	Functional Genetics	influenza
foravirumab	Crucell	post-exposure prevention of rabies
ibalizumab (TMB-355)	TaiMed Biologics USA	HIV-1 infection
KB001-A (antibody-fragment product)	KaloBios Pharmaceuticals, Sanofi Pasteur	<i>Pseudomonas</i> infections in cystic fibrosis patients, prevention of ventilator-associated pneumonia
KD-247	Kaketsuken	HIV-1 infection
MBL-HCV1	MassBiologics	Hepatitis C
MDX-066 (CDA-1)	Medarex	<i>C. difficile</i> disease

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
MDX-1303	Medarex, PharmAthene	anthrax
MEDI-557 (RSV mAb-extended half-life)	AstraZeneca, MedImmune	prevention of respiratory syncytial virus (RSV) infections
MK-3415A (actoxumab/bezlotoxumab)	Merck	<i>Clostridium difficile</i> infections
NM01	SRD Pharmaceuticals	HIV infection
Numax TM (motavizumab)	MedImmune	RSV
PRO 140	CytoDyn	HIV-1 infection
SAR279356 (anti-PNAG mAb)	Sanofi US	prevention of bacterial infections
Soliris [®] (eculizumab)	Alexion Pharmaceuticals	Shiga toxin <i>E. coli</i> -related hemolytic uremic syndrome (STEC-HUS)
streptococcal B vaccine conjugate	Novartis Vaccines	prevention of streptococcal B infections
Synagis [®] (palivizumab)	MedImmune	prevention of RSV infections
Tarvacin TM	Peregrine Pharmaceuticals	hepatitis C
TCN-032 (IgG mAb)	Theraclone Sciences	influenza A virus infections
TCN-202	Theraclone Sciences	CMV infections
Thraxiva TM (fully human anthrax mAb)	Emergent BioSolutions	post-exposure treatment of anthrax
TNX-355	Tanox	HIV infection
UB-421	United Biomedical	HIV-1 infection
XOMA 3AB	XOMA, National Institute of Allergy and Infectious Diseases	botulism
XTL 6865	XTL Biopharmaceuticals	hepatitis C
Blood disorders		
afelimomab	Abbot Laboratories	sepsis, septic shock
Benlysta [®] (belimumab)	GlaxoSmithKline	vasculitis
BI-655075	Boehringer Ingelheim Pharmaceuticals	blood coagulation disorders
eculizumab	Alexion Pharmaceuticals	paroxysmal nocturnal hemoglobinuria
ferroportin mAb	Eli Lilly	anemia
hepcidin mAb	Eli Lilly	anemia
ReoPro [®] (abciximab)	Eli Lilly	adjunct to percutaneous coronary intervention for the prevention of cardiac ischemic complications
SelGI	Selexys Pharmaceuticals	sickle cell anemia
urtioxazumab	Teijin Pharma	hemolytic uremic
Cardiovascular disease		
AMG 145	Amgen	Hypercholesterolemia, hyperlipoproteinemia type IIa
anti-fibrin mAb (3B6/22 Tc-99m)	Agenix	diagnosis of deep vein thrombosis, diagnosis of pulmonary embolism
anti-oxLDL (BI-204/RG7418)	BioInvent International, Genentech	atherosclerosis
anti-PCSK9 mAb (RG7652)	Genentech	cardiovascular disease
GSK249320	GlaxoSmithKline	stroke
IL-1 β antibody	Eli Lilly	cardiovascular disease
Ilaris [®] (canakinumab)	Novartis Pharmaceuticals	secondary prevention of cardiovascular events
inclacumab (RG1512)	Roche	peripheral vascular disease

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
MABp1	XBioTech	vascular restinosis
MLN 1202	Millennium Pharmaceuticals	atherosclerosis
pexelizumab	Alexion Pharmaceuticals, Procter & Gamble Pharmaceuticals	acute myocardial infarction, cardiopulmonary bypass
PF-04950615 (RN316)	Pfizer	hypercholesterolemia
SAR236553/REGN727 (anti-PCSK-9 mAb)	Regeneron Pharmaceuticals, Sanofi US	hypercholesterolemia
Diabetes and Related Conditions		
anti-CD3 mAb	MacroGenics	type 1 diabetes
gevokizumab (IL-1B inhibitor mAb)	XOMA	type 1 diabetes, type 2 diabetes
GSK1070806 (IL-18 mAb)	GlaxoSmithKline	type 2 diabetes
Ilaris [®] (canakinumab)	Novartis Pharmaceuticals	type 1 diabetes, type 2 diabetes
MABp1	XBioTech	type 2 diabetes
OKT3-gamma-1	Johnson & Johnson	type 1 diabetes
teplizumab	MacroGenics	type 1 diabetes
TRX 4 (anti-CD3)	TolerRx	type 1 diabetes
Eye Conditions		
anti-factor D (RG7417)	Genentech	geographic atrophy associated with age-related macular degeneration
anti-LINGO (BIIB033)	Biogen Idec	optic neuritis
gevokizumab (IL-1B inhibitor mAb)	XOMA	intermediate or posterior uveitis
GSK933776A (anti-B amyloid mAb)	GlaxoSmithKline	age-related macular degeneration
Humira [®] (adalimumab)	AbbVie	uveitis
iSONEP [™] (sonenpizumab)	Lpath	wet age-related macular degeneration
Lucentis [®] (ranibizumab)	Genentech, Roche	age-related macular degeneration
PF-04382923 (RN6G)	Pfizer	age-related macular degeneration
secukinumab (AIN457)	Alcon Labs, Novartis Pharmaceuticals	uveitis
Soliris [®] (eculizumab)	Alexion Pharmaceuticals	severe or refractory neuromyelitis optica
Genetic Disorders		
KRN-23	Kyowa Hakko Kirin Pharma	X-linked dominant hypophosphatemic rickets
Soliris [®] (eculizumab)	Alexion Pharmaceuticals	paroxysmal nocturnal hemoglobinuria (PNH)
Neurological Disorders		
AAB-002	Janssen Alzheimer Immunotherapy, Pfizer	Alzheimer's disease
AAB-003/PF-05236812	Janssen Alzheimer Immunotherapy, Pfizer	Alzheimer's disease
ABT-110	AbbVie	chronic pain
ALD403	Alder Biopharmaceuticals	prevention of migraine
ATI355 (anti-Nogo-A mAb)	Novartis Pharmaceuticals	spinal cord injury
BAN2401 (amyloid beta-protein inhibitor)	BioArctic Neuroscience, Eisai	Alzheimer's disease

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
bapineuzumab	Janssen Alzheimer Immunotherapy, Pfizer	Alzheimer's disease
crenezumab (anti-Aβeta)	Genentech	Alzheimer's disease
fulranumab	Janssen Research & Development	cancer pain
GSK1223249 (NOGO-A mAb)	GlaxoSmithKline	amyotrophic lateral sclerosis (ALS)
GSK933776A (anti-B amyloid mAb)	GlaxoSmithKline	Alzheimer's disease
LY2951742 (CGRP peptide)	Arteaus Therapeutics, Eli Lilly	migraine prevention
MEDI-5117 (anti-IL-6 mAb)	AstraZeneca, MedImmune	osteoarthritis pain
RG1450 (gantenerumab)	Roche	prodromal Alzheimer's disease
RN-307 (anti-CGRP mAb)	Labrys Biologics	migraine
RN624	Rinat Neuroscience	osteoarthritis pain
RN1219	Rinat Neuroscience	Alzheimer's disease
SAR228810 (anti-protofibrillar AB mAb)	Sanofi US	Alzheimer's disease
solanezumab (LY2062430)	Eli Lilly	Alzheimer's disease
tanezumab	Pfizer	chronic pain
Respiratory Disorders		
ABN 912	Novartis Pharmaceuticals	asthma, chronic obstructive pulmonary disorders (COPD)
ABX-IL8	Amgen	COPD
ALT-836	Altor BioScience, Genentech	acute lung injury, adult respiratory distress syndrome
AMG 157/MEDI-9929	Amgen, AstraZeneca	asthma
AMG 317	Amgen	asthma
AMG 761	Amgen	asthma
AMG 827 (brodalumab)	Amgen, AstraZeneca	asthma
benralizumab	AstraZeneca, MedImmune	asthma, chronic obstructive pulmonary disease (COPD)
Bosatria™ (mepolizumab)	GlaxoSmithKline	asthma
carlumab	Janssen Biotech	pulmonary fibrosis
CNTO-3157	Janssen Biotech	asthma
CNTO-5825	Janssen Biotech	allergic asthma
daclizumab (anti-CD25 MAb)	Protein Design Labs, Roche	asthma
FG-3019	FibroGen	idiopathic pulmonary fibrosis
GS-6624 (simtuzumab)	Gilead Sciences	idiopathic pulmonary fibrosis
KB003	KaloBios Pharmaceuticals	severe asthma
MEDI-528 (anti-TL-9 mAb)	MedImmune	asthma
MEDI-4212 (anti-IgE mAb)	AstraZeneca, MedImmune	asthma
MEDI-7814 (anti-C5/C5a mAb)	AstraZeneca, MedImmune	COPD
MEDI-8968 (anti-IL-1R mAb)	AstraZeneca, MedImmune	COPD
mepolizumab (anti-TL5 mAb)	GlaxoSmithKline	asthma and nasal polyposis
QAX576	Novartis Pharmaceuticals	asthma, idiopathic pulmonary fibrosis
QBX258	Novartis Pharmaceuticals	asthma
QGE031	Novartis Pharmaceuticals	allergic asthma
quilizumab (anti-M1 prime mAb)	Genentech	allergic asthma, allergic rhinitis
reslizumab	Cephalon	asthma, eosinophilic esophagitis
RG3637 (lebrikizumab)	Genentech, Roche	severe asthma
SAR156597 (bispecific interleukin-4/interleukin-13 mAb)	Sanofi US	idiopathic pulmonary fibrosis

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
SAR231893 (anti-IL4 mAb)	Regeneron Pharmaceuticals, Sanofi US	asthma
STX-100	Biogen Idec	idiopathic pulmonary fibrosis
TNX-832	Tanox	respiratory diseases
tralokinumab	AstraZeneca, MedImmune	asthma
Xolair® (omalizumab)	Genentech, Novartis Pharmaceuticals	pediatric asthma
Skin Diseases		
AbGn-168H	AbGenomics International	plaque psoriasis
AMG 827 (brodalumab)	Amgen, AstraZeneca	psoriasis
BT-061	AbbVie, Biotest	plaque psoriasis
CNTO-1959	Janssen Biotech	plaque psoriasis
gevokizumab (IL-1B inhibitor mAb)	XOMA	acne vulgaris
Humira® (adalimumab)	AbbVie	hidradenitis suppurativa
ixekizumab (IL-17 antibody)	Eli Lilly	psoriasis
MABp1	XBiotech	acne, psoriasis
MK-3222 (tildrakizumab)	Merck	plaque psoriasis
QGE031	Novartis Pharmaceuticals	atopic dermatitis
REGN846	Regeneron Pharmaceuticals	atopic dermatitis
SAR231893 (anti-IL4 mAb)	Regeneron Pharmaceuticals, Sanofi US	atopic dermatitis
secukinumab (AIN457)	Novartis Pharmaceuticals	plaque psoriasis
Xolair® (omalizumab)	Genentech, Novartis Pharmaceuticals	chronic idiopathic urticaria
Transplantation		
ASKP-1240	Astellas Pharma US, Kyowa Hakko Kirin Pharma	prevention of organ transplant rejection
Benlysta® (belimumab)	GlaxoSmithKline	immunosuppression
ORTHOCLONE OKT® 3 (muromomab-CD3)	Ortho Biotech	acute kidney transplant rejection, reversal of heart and liver transplant rejection
OKT3-gamma-1	Protein Design Labs, Johnson & Johnson	renal transplant rejection
Simulect® (basiliximab)	Novartis Pharmaceuticals	prevention of renal transplant rejection
Soliris® (eculizumab)	Alexion Pharmaceuticals	presensitized kidney transplant (acute humoral rejection)
TOL101	Tolera Therapeutics	prevention of transplant rejection
Zenapax® (daclizumab)	Roche	prophylaxis of acute kidney transplant rejection
Other		
anti-IL31	Bristol-Myers Squibb	immunology
anti-TWEAK (BIIB 023)	Biogen Idec	lupus nephritis
BAX-69	Baxter International	lupus nephritis
CNTO-5	Janssen Biotech, MorphoSys	inflammation
CNTO-136 (sirukumab)	Janssen Biotech	lupus nephritis
CR 0002	CuraGen	kidney inflammation
FB 301	Cytovance Biologics, Fountain Biopharma	hypersensitivity (IgE-mediated allergic diseases)

TABLE 1		
<u>Monoclonal Antibodies (mAb)</u>		
Product Name	Company	Indication
FG-3019	FibroGen	liver fibrosis due to chronic hepatitis B infection
fresolimumab (TGF β antagonist)	Genzyme	fibrosis
GS-6624 (simtuzumab)	Gilead Sciences	liver fibrosis, myelofibrosis
GSK1070806 (anti-interleukin 18 mAb)	GlaxoSmithKline	metabolic disorders
Humira [®] (adalimumab)	AbbVie	interstitial cystitis
LY2382770 (TGF- β antibody)	Eli Lilly	diabetic nephropathy
mAb	Genentech	metabolic disorders
mepolizumab (anti-IL5 mAb)	GlaxoSmithKline	hypereosinophilic syndrome, eosinophilic esophagitis
Meth-mAb	InterveXion Therapeutics	methamphetamine abuse
Stelara [®] (ustekinumab)	Janssen Biotech	primary biliary cirrhosis
VAY736	MorphoSys, Novartis Pharmaceuticals	inflammation
Xolair [®] (omalizumab)	Genentech, Tanox	peanut allergy

[0141] Non-limiting examples of protein-based or polypeptide-based biologics include cytokines (*e.g.*, interleukins), chemokines, growth factors, blood-production stimulating proteins (*e.g.*, erythropoietin), hormones (*e.g.*, Elonva[®] (follicle stimulating hormone), growth hormone), enzymes (*e.g.*, Pulmozyme[®] (dornase alfa)), clotting factors, insulin, albumin, fragments thereof, conservatively modified variants thereof, analogs thereof, and combinations thereof.

[0142] Examples of cytokines include, but are not limited to, TNF α , TNF-related weak inducer of apoptosis (TWEAK), osteoprotegerin (OPG), IFN- α , IFN- β , IFN- γ , interleukins (*e.g.*, IL-1 α , IL-1 β , IL-1 receptor antagonist (IL-1ra), IL-2, IL-4, IL-5, IL-6, soluble IL-6 receptor (sIL-6R), IL-7, IL-8, IL-9, IL-10, IL-12, IL-13, IL-15, IL-17, IL-23, and IL-27), adipocytokines (*e.g.*, leptin, adiponectin, resistin, active or total plasminogen activator inhibitor-1 (PAI-1), visfatin, and retinol binding protein 4 (RBP4)), and combinations thereof. In particular embodiments, the interleukin comprises IL-2 such as Proleukin[®] (aldesleukin; recombinant IL-2).

[0143] Examples of chemokines include, but are not limited to, CXCL1/GRO1/GRO α , CXCL2/GRO2, CXCL3/GRO3, CXCL4/PF-4, CXCL5/ENA-78, CXCL6/GCP-2, CXCL7/NAP-2, CXCL9/MIG, CXCL10/IP-10, CXCL11/I-TAC, CXCL12/SDF-1, CXCL13/BCA-1, CXCL14/BRAK, CXCL15, CXCL16, CXCL17/DMC, CCL1, CCL2/MCP-1, CCL3/MIP-1 α , CCL4/MIP-1 β , CCL5/RANTES, CCL6/C10, CCL7/MCP-3,

CCL8/MCP-2, CCL9/CCL10, CCL11/Eotaxin, CCL12/MCP-5, CCL13/MCP-4, CCL14/HCC-1, CCL15/MIP-5, CCL16/LEC, CCL17/TARC, CCL18/MIP-4, CCL19/MIP-3 β , CCL20/MIP-3 α , CCL21/SLC, CCL22/MDC, CCL23/MPIF1, CCL24/Eotaxin-2, CCL25/TECK, CCL26/Eotaxin-3, CCL27/CTACK, CCL28/MEC, CL1, CL2, CX₃CL1, and

5 combinations thereof.

[0144] Non-limiting examples of growth factors include epidermal growth factor (EGF), heparin-binding epidermal growth factor (HB-EGF), vascular endothelial growth factor (VEGF), pigment epithelium-derived factor (PEDF; also known as SERPINF1), amphiregulin (AREG; also known as schwannoma-derived growth factor (SDGF)), basic
 10 fibroblast growth factor (bFGF), hepatocyte growth factor (HGF), transforming growth factor- α (TGF- α), transforming growth factor- β (TGF- β 1, TGF- β 2, TGF- β 3, *etc.*), endothelin-1 (ET-1), keratinocyte growth factor (KGF; also known as FGF7), bone morphogenetic proteins (*e.g.*, BMP1-BMP15), platelet-derived growth factor (PDGF), nerve growth factor (NGF), β -nerve growth factor (β -NGF), neurotrophic factors (*e.g.*, brain-derived neurotrophic
 15 factor (BDNF), neurotrophin 3 (NT3), neurotrophin 4 (NT4), *etc.*), growth differentiation factor-9 (GDF-9), granulocyte-colony stimulating factor (G-CSF), granulocyte-macrophage colony stimulating factor (GM-CSF), myostatin (GDF-8), erythropoietin (EPO), thrombopoietin (TPO), and combinations thereof.

[0145] Examples of receptor construct-based or fusion protein-based biologics include, but
 20 are not limited to, naturally-occurring receptors linked to an immunoglobulin frame (*e.g.*, Orencia[®] (abatacept; immunoglobulin CTLA-4 fusion protein), Amevive[®] (alefacept; IgG1 fusion protein), ENBREL[®] (etanercept; recombinant human TNF-receptor fusion protein), engineered proteins combining two different polypeptide species (*e.g.*, Ontak[®] (denileukin diftitox; engineered protein comprising interleukin-2 and diphtheria toxin), and combinations
 25 thereof.

[0146] The present invention can therefore be used in methods for detecting and measuring the presence or level of a biologic in a sample from a subject receiving biologic therapy for one or more of the diseases or disorders referred to herein and Table 1, including one or more of the following:

30 [0147] Inflammatory diseases, such as inflammatory bowel disease (IBD) (*e.g.*, Crohn's disease (CD) and ulcerative colitis (UC)), uveitis, sarcoidosis, Wegener's granulomatosis, and other diseases with inflammation as a central feature;

[0148] Autoimmune diseases, such as rheumatoid arthritis (RA), multiple sclerosis (MS), systemic lupus erythematosus (SLE), ankylosing spondylitis (Bechterew's disease), lupus, psoriatic arthritis, juvenile idiopathic arthritis, psoriasis, erythematosus, and celiac disease;

[0149] Cancer, such as digestive and gastrointestinal cancers (*e.g.*, colorectal cancer, small intestine (small bowel) cancer; gastrointestinal stromal tumors, gastrointestinal carcinoid tumors, colon cancer, rectal cancer, anal cancer, bile duct cancer, gastric (stomach) cancer; esophageal cancer; appendix cancer; and the like); gallbladder cancer; liver cancer; pancreatic cancer; breast cancer; lung cancer (*e.g.*, non-small cell lung cancer); prostate cancer; ovarian cancer; renal cancer (*e.g.*, renal cell carcinoma); cancer of the central nervous system; skin cancer; choriocarcinomas; head and neck cancers; hematological malignancies (*e.g.*, leukemia, lymphoma such as B-cell non-Hodgkin's lymphoma); osteogenic sarcomas (*e.g.*, Ewing sarcoma); soft tissue sarcomas (*e.g.*, Dermatofibrosarcoma Protuberans (DFSP), rhabdomyosarcoma); other soft tissue malignancies, and papillary thyroid carcinomas;

[0150] Infectious diseases, such as *C. difficile* disease, respiratory syncytial virus (RSV), HIV, anthrax, candidiasis, staphylococcal infections, and hepatitis C;

[0151] Blood disorders, such as sepsis, septic shock, paroxysmal nocturnal hemoglobinuria, and hemolytic uremic syndrome;

[0152] Cardiovascular disease, such as atherosclerosis, acute myocardial infarction, cardiopulmonary bypass, and angina;

[0153] Metabolic disorders, such as diabetes, *e.g.*, type 1 diabetes mellitus and type 2 diabetes;

[0154] Genetic disorders, such as paroxysmal nocturnal hemoglobinuria (PNH);

[0155] Neurological disorders, such as osteoarthritis pain and Alzheimer's disease;

[0156] Respiratory disorders, such as asthma, chronic obstructive pulmonary disorders (COPD), nasal polyposis, and pediatric asthma;

[0157] Skin diseases, such as psoriasis, including chronic moderate to severe plaque psoriasis;

[0158] Transplant rejection, such as acute kidney transplant rejection, reversal of heart and liver transplant rejection, prevention of renal transplant rejection, prophylaxis of acute kidney transplant rejection, and renal transplant rejection; and/or

[0159] Other disorders, such as kidney inflammation, postmenopausal osteoporosis (bone disorders), hypereosinophilic syndrome, eosinophilic esophagitis, and peanut allergy.

[0160] In particular embodiments, the subject has an inflammatory disease (*e.g.*, inflammatory bowel disease (IBD) such as Crohn's disease (CD) or ulcerative colitis (UC)) or an autoimmune disease (*e.g.*, rheumatoid arthritis).

VIII. Examples

[0161] The following examples are offered for illustrative purposes, and are not intended to limit the invention in any manner. Those of skill in the art will readily recognize a variety of non-critical parameters which can be changed or modified to yield essentially the same results.

Example 1

[0162] This example illustrates the validation of a homogeneous mobility shift assay (HMSA) for the measurement of vedolizumab (VLM) and anti-VLM antibodies (ATV) in inflammatory bowel disease (IBD) patient serum.

[0163] Vedolizumab (VLM), an $\alpha 4\beta 7$ integrin antagonist, is a therapeutic monoclonal antibody recently approved for use in moderate to severe ulcerative colitis or Crohn's disease patients that have failed to demonstrate adequate response to conventional therapies or TNF alpha antagonists. $\alpha 4\beta 7$ integrin is a gut-specific heterodimeric glycoprotein that is important for leukocyte homing to sites of inflammation within the intestine, via its interaction with mucosal addressin cell adhesion molecule-1 (MadCam-1) expressed on the intestinal vascular endothelium. Availability of diagnostic tests to accurately measure VLM drug levels and anti-VLM antibodies is necessary for the effective use of this novel therapeutic in IBD patients. In succession with our anti-TNF alpha therapeutic drug monitoring assays, we have now developed and validated assays to measure VLM and anti-VLM levels in patient serum.

Methods

[0164] Soluble $\alpha 4\beta 7$ protein heterodimer was expressed and purified from mammalian cells (*e.g.*, CHO cells). In particular, the integrin $\alpha 4$ subunit was truncated after the thigh domain ($\alpha 4\Delta 620$) and integrin $\beta 7$ after the I-EGF 1 domain ($\beta 7\Delta 527$) (*see*, Figure 13). Acidic/Basic α -helical coiled coil peptides containing one disulfide-bridge forming Cys residue were attached to the C-terminus of the subunits to stabilize the heterodimer (Takagi *et al.*, *Embo J.*, 18:4607-4615 (2003); O'Shea *et al.*, *Curr. Biol.*, 3:658-667 (1993)). A

hexahistidine tag was fused to the C-terminus of the basic peptide to facilitate purification. This recombinant heterodimer was used in a “competition” based homogeneous mobility shift assay (HMSA) format to measure VLM levels in sera of patients on VLM therapy. Patient serum was combined with $\alpha 4\beta 7$ to allow therapeutic VLM to bind with $\alpha 4\beta 7$.
5 Subsequently, fluorescently labeled VLM competed with unlabeled VLM in patient sera for binding to its target, $\alpha 4\beta 7$ followed by separation on HPLC size exclusion chromatography. The amount of “free” VLM-alexa fluor determined VLM levels in patient sera. For the anti-VLM assay, standard curves were created by incubating normal human serum containing known amounts of rabbit anti-VLM antibodies with fluorescently labeled VLM; bound and
10 free VLM were then separated by SEC-HPLC. Method validation was determined according to industry recommendations.

Results

[0165] Sensitivity for VLM was 0.35 $\mu\text{g/mL}$. The lower and upper limits of quantitation (LLOQ and ULOQ) were determined to be 0.625 $\mu\text{g/mL}$ and 14 $\mu\text{g/mL}$, respectively. For
15 anti-VLM, the LLOQ and ULOQ were 3.13 U/ml and 150 U/ml, respectively. The standard curves generated for each assay showed high reproducibility and sensitivity. Inter- and intra-assay precision showed less than 10% CV and accuracy was within 20%. There was no significant interference from lipemic, hemolyzed, or rheumatoid factor (Rf) serum. *See*, Figures 4-10.

Conclusions

[0166] A sensitive and specific assay has been developed and validated to measure VLM and anti-VLM levels in patients undergoing treatment for IBD. The assay format required a unique approach owing to the complexity of a large heterodimeric, heavily glycosylated membrane bound antigen as drug target. Both the drug and anti-drug assays demonstrate
25 high accuracy and precision with high sensitivity with a high tolerance to known interfering agents. The assays are useful for clinical monitoring and drug optimization in individual patients.

Example 2

[0167] This example illustrates the validation of a homogeneous mobility shift assay (HMSA) and for the measurement of ustekinumab (UTK) and antibodies-to-ustekinumab (ATU) in inflammatory bowel disease (IBD) patient serum.
30

[0168] Ustekinumab is a therapeutic monoclonal antibody which has potential utility in the treatment of IBD in patients that have failed to respond to conventional therapies or TNF alpha antagonists. Ustekinumab is a monoclonal antibody with specificity for interleukin 12 and interleukin 23 via their common p40 subunit and blocks inflammation through these pathways. Here we describe the analytical validation of a HPLC-based high mobility shift assay to measure UTK as well as antibodies-to-ustekinumab.

Methods

[0169] Recombinant, soluble IL12p40 subunit (*e.g.*, SEQ ID NO:7) was expressed and purified from mammalian cells. The assay takes advantage of fluorescently labeled UTK competing with unlabeled UTK in patient sera for binding to IL12p40. After incubation of patient sera with recombinant IL12p40, UTK-alexa fluor 488 was added, before running samples on an HPLC size exclusion column which separates “free” UTK-alexa fluor 488 from UTK-alexa fluor 488 which is bound to IL12p40. The amount of “free” UTK-alexa fluor 488 is a measure of the amount of therapeutic UTK in patient sera. The area under the curve (AUC) of the “free” UTK-alexa fluor 488 is plotted against the log of the UTK concentration in known standard samples and UTK concentration in patient sera are calculated by interpolation. For the ATU assay, standard curves were created by incubating normal human serum containing known amounts of rabbit ATU with fluorescently labeled UTK; bound and free GLM were then separated by SEC-HPLC. Method validation was determined according to industry recommendations.

Results

[0170] Sensitivity for UTK was 0.15 µg/ml. The upper limit of quantitation (ULOQ) was determined to be 8 µg/mL and the lower limit of quantitation (LLOQ) was determined to be 0.625 µg/ml. For ATU, the LLOQ and ULOQ were 3.13 U/ml and 150 U/ml respectively. The standard curves generated for each assay showed high reproducibility and sensitivity. Inter- and intra-assay precision showed less than 10% CV and accuracy was within 20%. There was no significant interference from lipemic, hemolyzed, or rheumatoid factor (Rf) serum. *See*, Figures 11-12.

Conclusions

[0171] This study describes the validation of a novel “competition” based assay to measure ustekinumab levels in patients undergoing treatment for IBD. The assay format required a unique approach owing to complexity of working with this antigen, due to its tendency to

form homodimers and heterodimers in serum making usage of a conventional HMSA assay difficult.

[0172] The competition-based HMSA assay showed high accuracy and precision with high sensitivity. In addition, the assay has high tolerance to known interfering agents. These assays are useful for clinical monitoring and drug optimization in individual patients.

Example 3

[0173] This example illustrates an exemplary vedolizumab (VLM) competition-based assay methodology of the present invention. One of ordinary skill in the art will appreciate that the assay methodology described in this example is applicable to determining the presence or level of ustekinumab (UTK) as well as other biologics in a sample in situations where the complexity of working with the antigen that binds to the biologic (*e.g.*, the antigen is a membrane-bound protein, a glycosylated protein, a multimeric protein, an insoluble protein, a protein that is difficult to express or purify, and/or a large protein) necessitate the use of a soluble form (*e.g.*, a soluble fragment, variant, or monomer) of the antigen.

[0174] Normal human serum (NHS) samples spiked with known amounts of vedolizumab (VLM) are serially diluted. Two-fold serial dilutions starting from 80 µg/ml VLM are diluted in NHS to make 10-point curve (*i.e.*, 80 µg/ml to 0.15625 µg/ml). NHS spiked with 12, 4, and 1 µg/ml VLM are used as positive controls. Standard serum samples, positive controls, and patient samples are added to a 96 well plate. Patient samples are added undiluted or diluted 4x or 8x in NHS to increase the assay dynamic range. Soluble α4β7 antigen and assay diluent are added. The plate is placed on a shaker and allowed to incubate at room temperature for 1 hour. After 1 hour, labeled VLM (*e.g.*, VLM-Alexa Fluor[®] 488) is added. The plate is again placed on a shaker for a 1 hour incubation. Samples are filtered using a 0.2 µm filter plate. Samples are loaded onto an HPLC autosampler and run sequentially through a size exclusion chromatography column (*e.g.*, a Phenomenex BioSep-SEC-s3000 column) which separates free labeled VLM (*e.g.*, VLM-Alexa Fluor[®] 488) from labeled VLM bound to the soluble α4β7 antigen.

[0175] Software written in R-programming language is used to identify the peak representing free labeled VLM (*e.g.*, VLM-Alexa Fluor[®] 488) and to determine the area of the peak. The area of this peak gets larger when there is VLM in the patient's sera. By comparing the size of this peak to the standard curve, one can interpolate the patient's VLM levels.

[0176] Prism (*e.g.*, GraphPad Prism 6) is used to generate a standard curve by plotting the area of free labeled VLM (*e.g.*, VLM-Alexa Fluor[®] 488) as a function of serum VLM levels. By comparing the size of this peak to the standard curve, one can interpolate the patient's VLM drug concentration.

- 5 [0177] The assay described in this example is premised on the competition between the VLM in a sample from a patient receiving VLM therapy and the labeled VLM added to the sample reaction for binding to the soluble $\alpha 4\beta 7$ antigen. The relative ratios of labeled and unlabeled VLM determines how much $\alpha 4\beta 7$ antigen is bound to each and determines the free labeled VLM (*e.g.*, VLM-Alexa Fluor[®] 488) peak area. The more drug present in the patient
10 sample, the more the labeled VLM remains free as opposed to bound to the $\alpha 4\beta 7$ antigen.

Example 4

[0178] This example illustrates the validation of a homogeneous mobility shift assay (HMSA) for the measurement of vedolizumab (VLM) and anti-VLM antibodies in inflammatory bowel disease (IBD) patient serum.

15 **Background and Aims**

- [0179] Vedolizumab, an $\alpha 4\beta 7$ integrin antagonist, is a therapeutic monoclonal antibody recently approved for use in moderate to severe ulcerative colitis and Crohn's disease patients that have failed to demonstrate adequate response to conventional therapies or TNF α antagonists. Availability of diagnostic tests to accurately measure serum VLM and anti-
20 VLM (ATV) levels is necessary for the effective use of this novel therapeutic in IBD patients. Here we describe the analytical validation of the HMSA developed to measure VDM and ATV levels in patient serum as well as its clinical utility.

Methods

- [0180] Soluble $\alpha 4\beta 7$ heterodimer (*e.g.*, $\alpha 4\Delta 620/\beta 7\Delta 527$ heterodimer; *see*, Figure 13) was
25 expressed and purified in mammalian cells. This recombinant heterodimer was used in a "competition" based HMSA format to measure serum VLM levels in patients on VLM therapy. Patient serum was combined with $\alpha 4\beta 7$ to allow therapeutic VLM to bind with $\alpha 4\beta 7$. Subsequently, fluorescently labeled VLM competed with unlabeled VLM in patient sera for binding to its target, $\alpha 4\beta 7$, followed by separation on HPLC size exclusion
30 chromatography (*see*, Figure 1). For the anti-VLM assay, standard curves were created by incubating normal human serum containing known amounts of rabbit anti-VLM antibodies

with fluorescently labeled VLM; bound and free VLM were then separated by SEC-HPLC. Validation was performed according to industry recommendations (Shankar, G., et al. 2008).

Results

[0181] VLM and ATV assays show high intra-assay precision and accuracy. Intra-assay precision is less than 10% and intra-assay accuracy is less than 15% error. Run to run, instrument to instrument, and analyst to analyst variability are less than 15% CV and less than 20% error in almost all cases. (Table 2). Sensitivity for VLM was 0.348 µg/mL with a dynamic range of 0.625-14 µg/mL (Table 3). The limit of detection for the ATV assay is <1.56 U/mL. A precise value could not be determined as it is too low to interpolate our curve. The dynamic range for ATV was 3.13-150 U/mL in undiluted serum (Table 3). Normal human serum spiked with VLM or ATV showed good linearity and recovery across serial dilutions (Figure 14). ATV shows high drug tolerance at levels of 20 µg/mL of VLM (Table 4). Common interfering agents in patient serum were tested and did not interfere at levels seen in IBD patients (Figure 15).

Table 2. Validation data for accuracy and precision of the assays.

Accuracy and Precision of the VLM-HMSA

	Intra-Assay Precision						Inter-Assay Precision					
	(n=10)			Run to Run (n=10)			Analyst to Analyst (n=5)			Instrument to Instrument (n=10)		
	Hgh	Med	Low	High	Med	Low	High	Med	Low	High	Med	Low
Expected (µ/mL)	12	4	1	12	4	1	12	4	1	12	4	1
Measured (Mean, U/mL)	10.34	4.36	0.89	11.22	4.16	1.13	10.51	4.19	1.05	12.30	4.18	1.29
SD	0.56	0.21	0.10	1.11	0.17	0.18	0.47	0.33	0.16	0.86	0.19	0.08
CV%	5.45	4.82	11.64	9.87	4.07	14.53	4.46	7.97	14.91	6.97	4.44	6.07
Accuracy (%Error)	13.83	9.07	11.00	6.48	3.91	12.72	12.44	4.76	5.30	2.46	4.44	29.04

Accuracy and Precision of the ATV-HMSA

	Intra-Assay Precision						Inter-Assay Precision					
	(n=10)			Run to Run (n=10)			Analyst to Analyst (n=10)			Instrument to Instrument (n=10)		
	Hgh	Med	Low	High	Med	Low	High	Med	Low	High	Med	Low
Expected (U/mL)	50	25	12.5	50	25	12.5	50	25	12.5	50	25	12.5
Measured (Mean, U/mL)	47.31	23.32	11.72	51.53	25.01	13.54	51.52	25.56	13.51	47.31	23.34	11.75
SD	2.82	1.51	0.79	5.00	3.39	1.76	5.00	2.43	1.86	2.73	1.46	0.78
CV%	5.96	6.46	6.75	9.71	13.56	13.00	9.71	9.52	13.78	5.77	6.25	6.62
Accuracy (%Error)	5.37	6.71	6.27	7.93	0.04	8.32	3.04	2.25	8.06	5.38	6.65	5.96

Table 3. Limits of quantitation for each assay.

Assay	Limit of Blank	Limit of Detection	Lower Limit of Quantitation	Upper Limit of Quantitation
VLM-HMSA	n=30	n=30	n=36	n=36
	0.054µg/mL	0.348µg/mL	0.625µg/mL	14µg/mL
ATV-HMSA	n=25	n=25	n=36	n=36
	<1.56U/mL	<1.56U/mL	3.13U/mL	150U/mL

Table 4. Drug interference in the anti-VLM assay.

Vedo [µg/mL]	Total ATV [U/mL]
20	51.6
0	49.2
20	26.4
0	24.8
20	13
0	14

5 Conclusions

[0182] A sensitive and specific assay has been developed and validated to measure VLM and anti-VLM levels in patients undergoing treatment for IBD. The assay format required a unique approach owing to the complexity of a large heterodimeric, heavily glycosylated membrane protein as drug target. Both the drug and anti-drug assays demonstrate high accuracy and precision with tolerance to known interfering agents. The development of VLM and anti-VLM assays is useful for clinical monitoring and drug optimization in individual patients.

Example 5

[0183] This example illustrates the validation of a homogeneous mobility shift assay (HMSA) for the measurement of ustekinumab (UTK) and antibodies-to-ustekinumab (ATU) in inflammatory bowel disease (IBD) patient serum.

Background and Aims

[0184] Ustekinumab is a therapeutic monoclonal antibody which has potential utility in the treatment IBD patients that have failed to respond to conventional therapies or TNFα antagonists. Ustekinumab is specific for IL-12 and IL-23 via their common p40 subunit and blocks inflammation through these pathways. Availability of diagnostic tests to accurately measure serum UST and anti-UST (ATU) levels is necessary for the effective use of this

novel therapeutic in IBD patients. Here we describe the analytical validation of a “competition” based HMSA developed to measure UTK levels as well as conventional HMSA to measure ATU levels in patient serum.

Methods

5 [0185] Recombinant, soluble IL12p40 (*e.g.*, SEQ ID NO:7) was expressed and purified from mammalian cells. Recombinant IL12p40 was used in a “competition” based HMSA format to measure serum UTK levels in patients on UTK therapy. Patient serum was combined with rIL12p40 to allow therapeutic UTK to bind with rIL12p40. Subsequently, fluorescently labeled UTK competed with unlabeled UTK in patient sera for binding to its
10 target, rIL12p40, followed by separation on HPLC size exclusion chromatography (*see*, Figure 16). For the anti-UTK assay, standard curves were created by incubating normal human serum containing known amounts of rabbit anti-UTK antibodies with fluorescently labeled UTK; bound and free UTK were then separated by SEC-HPLC.

Results

15 [0186] UTK and ATU assays show high intra-assay precision and accuracy. Intra-assay precision is less than 10% and intra-assay accuracy is less than 15% error. Run to run, instrument to instrument, and analyst to analyst variability are less than 15% CV and less than 20% error. (Table 5). Sensitivity for UTK is 0.224 µg/mL with a dynamic range of 0.625-10 µg/mL (Table 6). The limit of detection for the ATU assay is <1.56 U/mL. A
20 precise value could not be determined as it is too low to interpolate from the standard curve. The dynamic range for ATU is 3.13-150 U/mL in undiluted serum (Table 6). Normal human serum spiked with UTK or ATU showed good linearity and recovery across serial dilutions within the assay’s dynamic range (Figure 17). ATU shows high drug tolerance. Levels of 20 µg/mL of UTK do not significantly interfere with ATU detection (Table 7). Common
25 interfering agents in patient serum were tested and did not interfere at levels seen in IBD patients. Only highly hemolyzed serum may potentially interfere but not at levels normally seen in patient sera (Figure 18).

Table 5. Validation data for precision and accuracy of the assays.

Accuracy and Precision of the UTK-HMSA

	Intra-Assay Precision						Inter-Assay Precision					
	(n=5)			Run to Run (n=5)			Analyst to Analyst (n=5)			Instrument to Instrument (n=5)		
	Hgh	Med	Low	High	Med	Low	High	Med	Low	High	Med	Low
Expected (µg/mL)	10	5	2.5	10	5	2.5	10	5	2.5	10	5	2.5
Measured (Mean, µg/mL)	8.89	4.52	2.40	10.59	5.03	2.59	10.09	4.93	2.66	10.29	5.26	2.89
SD	0.62	0.08	0.05	0.97	0.29	0.13	0.52	0.27	0.18	0.49	0.29	0.09
CV%	6.96	1.76	2.16	9.18	5.67	5.06	5.18	5.39	6.95	4.79	5.44	3.15
Accuracy (%Error)	11.14	9.55	4.03	5.88	0.69	3.71	0.94	1.42	6.23	2.87	5.28	15.43

Accuracy and Precision of the ATU-HMSA

	Intra-Assay Precision						Inter-Assay Precision					
	(n=5)			Run to Run (n=5)			Analyst to Analyst (n=5)			Instrument to Instrument (n=5)		
	Hgh	Med	Low	High	Med	Low	High	Med	Low	High	Med	Low
Expected (U/mL)	50	20	10	50	20	10	50	20	10	50	20	10
Measured (Mean, U/mL)	49.61	21.08	9.47	50.63	21.51	9.20	51.24	21.48	9.19	49.61	20.39	8.78
SD	0.33	0.83	0.78	1.16	0.73	0.60	1.80	0.73	0.61	1.85	0.69	0.20
CV%	0.67	3.95	8.27	2.30	3.39	6.52	3.51	3.39	6.60	3.74	3.37	2.30
Accuracy (%Error)	0.79	5.42	5.30	1.25	7.54	7.96	2.47	7.41	8.11	0.79	1.97	12.23

Table 6. Limits of quantitation for each assay.

Limits of Quantitation

Assay	Limit of Blank	Limit of Detection	Lower Limit of Quantitation	Upper Limit of Quantitation
UTK-HMSA	n=30	n=30	n=36	n=36
	0.077µg/mL	0.224µg/mL	0.625µg/mL	10µg/mL
ATU-HMSA	n=30	n=30	n=30	n=30
	Too low to interpolate	< 1.56U/mL	3.13U/mL	150U/mL

5

Table 7. Drug interference in the anti-UTK assay.

UTK [µg/mL]	Total ATU [U/mL]
20	58.8
0	51.7
20	16.7
0	22.1
20	5.6
0	9.2

Conclusions

[0187] HMSA for UTK and ATU showed high accuracy and precision across a wide dynamic range. The HMSA platform allowed detection of UTK and ATU even in the presence of interfering agents which are known to limit the utility of ELISA/ECLIA methods.

- 5 The development of a new “competition” based HMSA platform allows for the measurement of therapeutic drug levels when a conventional HMSA is not feasible.

Example 6

[0188] This example illustrates experiments performed to improve the dynamic range of the homogeneous mobility shift assay (HMSA) for the measurement of vedolizumab (VDZ).

- 10 [0189] In an effort to improve the assay’s dynamic range, the assay was modified by increasing the amount of VDZ-Alexa488 used in the assay 1.6-fold. Increasing the amount of labeled VDZ and proportionately increasing the amount of $\alpha 4\beta 7$ antigen (*e.g.*, $\alpha 4\Delta 620/\beta 7\Delta 527$ heterodimer; *see*, Figure 13) increased the dynamic range of the assay. In addition, changing the amount of $\alpha 4\beta 7$ relative to labeled VDZ affected the lower limit of
- 15 quantification. Excess antigen relative to labeled VDZ made the assay relatively insensitive to drug in patient sera and raised the lower limit of quantification. Conversely, too little antigen relative to labeled VDZ had the effect of lowering the upper limit of quantification. Increasing the amount of labeled VDZ from 75 ng/well to 120 ng/well (1.6-fold) and titrating the antigen such that the presence of antigen binds up 75-80% of the labeled VDZ provided
- 20 the best compromise between the low-end sensitivity needed as well as an improved dynamic range that would enable the measurement of drug in patient sera without requiring dilutions in most cases.

- [0190] In addition, the data is plotted such that the area of the VDZ-Alexa488 peak (without dividing by the Blocked-Alexa 488 control peak area) is plotted against the log of
- 25 the VDZ concentration.

1. Assay Limits

Limit of Blank

- [0191] The Limit of Blank (LOB) was determined from 30 replicates of the standard curve blank. The standard curve blank (negative control) in all assays consisted of 4% normal
- 30 human serum + 40ng VDZ-Alexa488 + 165ng $\alpha 4\beta 7$ per 100 μ L injection. The average (mean) + 1.645SD of the VDZ-Alexa488 peak area was calculated and then used for calculation of the LOD.

Table 8. Assay Limit of Blank (LOB).

N	30
Mean (VDZ-Alexa488 peak area)	1.91
SD (VDZ-Alexa488 peak area)	0.09
The average (mean) + 1.645SD	2.06
Interpolated LOB	0.156/ml*

*The average (mean) + 1.645SD was just below the lowest point of the curve and therefore the lowest point on the curve was taken for the LOB.

Limit of Detection

[0192] The Limit of Detection (LOD) was determined by utilizing the measured LOB and replicates of serum containing VDZ at a low concentration which is approaching the LOB. Standard 10 was chosen because it is the lowest point on the curve and nearest the LOB. The LOD was calculated using the equation: $LOD = LOB + 1.645(SD_{\text{low concentration sample}})$ (Armbruster et al., 2008). The value was then interpolated from the averaged standard curve of the experiments used in the calculation to yield the concentration in $\mu\text{g/mL}$.

Table 9. Assay Limit of Detection.

N	30
Mean (VDZ-Alexa488 peak area)	1.91
SD (VDZ-Alexa488 peak area)	0.061
1.645*SD	0.100
Interpolated (1.645*SD)	Too low to interpolate
LOD=The average (mean) + 3*SD	0.457 $\mu\text{g/mL}$

Limit of Quantitation

[0193] The Lower Limit of Quantitation (LLOQ) was determined by analyzing interpolated concentrations of 30 replicates of a low concentration VDZ positive sample. In this case, standard 8 (effective serum concentration of 0.625 $\mu\text{g/mL}$) was chosen. The upper limit of quantitation (ULOQ) was determined by analyzing 30 replicates of a high concentration VDZ positive sample (effective serum concentration equal to 14 $\mu\text{g/mL}$). LLOQ was defined as the concentration that results in a $CV \leq 20\%$ with $\text{Error} \leq 25\%$ and thus measures the assay's precision and accuracy at a low analyte concentration. The ULOQ was also qualified by $CV \leq 20\%$ with $\text{Error} \leq 25\%$.

Table 10. LLOQ and ULOQ.

	LLOQ	ULOQ
N	30	30
Expected (µg/mL)	1	25
Mean (µg/mL)	1.05	20.95
SD (µg/mL)	0.12	1.44
CV (%)	11.39	6.86
Error (%)	5.13	-16.19

These criteria resulted in the following values for the Assay Limits:

5 Table 11. Assays Limits of Quantitation.

LOD	0.457 µg/mL
LLOQ	1 µg/mL
ULOQ	25 µg/mL

2. Interference

Antibody to Vedolizumab (ATV) Interference

- 10 [0194] In this experiment, 2.5, 5, 10 and 20 µg/mL VDZ was added to either normal human serum or various concentrations of rabbit ATV positive serum. These samples were then analyzed using the assay and the % recovery calculated.

15 Table 12. ATV Interference.

VDZ spike (µg/mL)	Measured VDZ (µg/mL) Control (No ATV)	% Recovery					
		0 U/mL	3.13 U/mL	6.25 U/mL	12.5 U/mL	25 U/mL	50 U/mL
2.5	2.23	100	90	82	62	0	0
5	4.47	100	93	88	77	60	22
10	9.15	100	91	85	77	71	55
20	17.53	100	87	81	73	69	59

[0195] The assay tolerance to ATV is up to 6.25 U/ml in this experiment. The ATV interference is expected as neutralizing ATV will compete with $\alpha 4\beta 7$ for binding to patient VDZ. Only VDZ bound to non-neutralizing ATV is expected to be detected.

20 Integrin $\alpha 4\beta 7$ Substrate Interference

[0196] To test the interference of $\alpha 4\beta 7$, titrations were performed in the range of 1 to 1000 ng/mL. VDZ concentrations were plotted on the intercept of the X-axis correspond to the zero concentration of each interfering agent (Figure 19).

Table 13. $\alpha 4\beta 7$ Interference (100 ng/mL).

VDZ Positive Control	High Control (12 $\mu\text{g/mL}$)	Medium Control (4 $\mu\text{g/mL}$)	Low Control (1 $\mu\text{g/mL}$)
Mean ($\mu\text{g/mL}$)	12.57	3.77	1.22
SD ($\mu\text{g/mL}$)	0.23	0.01	0.04
CV (%)	1.8	0.3	3.5
Recovery (%)	109.6	108.1	106.7

[0197] There was no significant interference from levels of $\alpha 4\beta 7$ far exceeding levels one would expect to see in sera from patients being treated with VDZ.

3. Antigen Stability

[0198] The stability of $\alpha 4\beta 7$ antigen was assessed through storage at 4 degrees (4°C) for 1 week. After 1 week, the $\alpha 4\beta 7$ was removed from 4°C storage and analyzed against a fresh aliquot stored at -70°C. The data from the 4°C storage samples was then compared to that obtained before storage and the percent error calculated. The aliquots were deemed stable at 4°C if the error was within 25%.

Table 14. Accelerated Stability.

Time Point		Week 0			Week 1		
VDZ Pos. Ctrl		High	Med	Low	High	Med	Low
4°C	Mean ($\mu\text{g/mL}$)	10.00	3.15	0.88	11.06	4.42	1.12
	SD ($\mu\text{g/mL}$)	0.21	0.45	0.06	0.81	0.21	0.06
	CV (%)	2.11	14.24	6.65	7.32	4.81	5.21
	Error (%)	-16.64	-21.35	-12.34	-7.81	10.50	12.05

[0199] Accelerated stability of $\alpha 4\beta 7$. Antigen stored at -70°C was used in the assay to test VDZ positive controls. Controls came out within specifications (Error $\leq 25\%$ and CV $\leq 25\%$). Antigen stored at 4°C was used in the assay to test VDZ positive controls. These controls also came out within specifications. In addition, there was no obvious loss of antigen potency as a result of being stored at 4°C versus -70°C.

4. Additional Dynamic Range Experiments

[0200] Three standard curves were generated using a 1X amount of labeled VDZ (e.g., 75 ng/well of VDZ-Alexa488), as well as 2X and 4X amounts. Figure 20 shows that proportionately increasing the concentrations of both VDZ-Alexa488 and integrin $\alpha 4\beta 7$ used in the assay increases the assay's dynamic range. It also increases anti-drug antibody tolerance. A maximum LLOQ of 1 $\mu\text{g/mL}$ was the criteria for this assay. In particular,

increasing the amount of labeled VDZ and $\alpha 4\beta 7$ by 1.6-fold enabled a maximum LLOQ of 1 $\mu\text{g/mL}$ while providing an improved top-end range of 25 $\mu\text{g/mL}$ (*see*, Table 11).

[0201] Although the foregoing invention has been described in some detail by way of illustration and example for purposes of clarity of understanding, one of skill in the art will appreciate that certain changes and modifications may be practiced within the scope of the appended claims. In addition, each reference provided herein is incorporated by reference in its entirety to the same extent as if each reference was individually incorporated by reference.

INFORMAL SEQUENCE LISTING

SEQ ID NO:1

Human alpha 4 integrin fragment

MAWEARREPGPRRAAVRETVMLLLCLGVPTGRPYNVDTESALLYQGPHTLFGYSVVLHSHGANRWLL
 5 VGAPTANWLANASVINPGAIRYCRIGKNPGQTCEQLQLGSPNGEPCGKTCLEERDNQWLGVTLRSRQPG
 ENGSIVTCGHRWKNIFYIKNENKLPTGGCYGVPPDLRTELKRIAPCYQDYVKKFGENFASCQAGISS
 FYTKDLIVMGAPGSSYWTGSLFVYNITTNKYKAFLDKQNQVKFGSYLGYSVGAGHFRSQHTTEVVGGA
 PQHEQIGKAYIFSIDEKELNILHEMKGKKLGSYFGASVCAVDLNADGFSDLLVGAPMQSTIREEGRVF
 10 VYINSGSGAVMNAMETNLVGSCKYARFGESIVNLGIDNDGFEDVAIGAPQEDDLQGAIIYINGRAD
 GISSTFSQRIEGLQISKSLSMFGQSISGQIDADNNGYVDVAVGAFRSDSAVLLRTRPVVIVDASLSHP
 ESVNRTKFDCEVNGWPSVCIDLTLCFSYKGKEVPGYIVLFYNMSLDVNRKAESPFRFYSSNGTSDVI
 TGSIQVSSREANCRTHQAFMRKDVRDILTPIQIEAAYHLGPHVISKRSTEEFPPLQPILOQKKEKDIM
 KKTINFAR

SEQ ID NO:2

Human beta 7 integrin fragment

MVALPMVLVLLLVLSSRGESELDKIPSTGDATEWRNPHLSMLGSCQPAPSCQKCILSHPSCAWCKQLN
 FTASGEAEARRCARREELLARGCPLEELPEPRGQQEVLQDQPLSQGARGEGATQLAPQVRVTLRPGE
 20 PQQQLQVRFLRAEGYPVDLYYLMDSLSSMKDDLERVRQLGHALLVRLQEVTHSVRIGFGSFVDKTVLPF
 VSTVPSKLRHPCPTRLERCQSPFSFHHVLSLTGDAQAFEREVGRQSVSGNLDSPEGGFDAIQAALCQ
 EQIGWRNVSRLLVFTSDDTFHTAGDGKLGFI FMPSDGHCHLDSNGLYSRSTEFDPYPSVGQVAQALSAA
 NIQPIFAVTSAAALPVYQELSKLIPKSAVGELSEDSSNVVQLIMDAYNSLSSTVTLEHSSLPPGVHISY
 ESQCEGPEKREGKAEDRGQCNHVRINQTVTFWVSLQATHCLPEPHLLRLRALGFSEELIVELHTLCLDC
 25 NCSDTQPQAPHCSDGQGHLCGVCSCAPGRLGRLCECSVAELSSPDLESGC

SEQ ID NO:3

Human alpha 4 integrin fragment with acidic peptide

MAWEARREPGPRRAAVRETVMLLLCLGVPTGRPYNVDTESALLYQGPHTLFGYSVVLHSHGANRWLL
 VGAPTANWLANASVINPGAIRYCRIGKNPGQTCEQLQLGSPNGEPCGKTCLEERDNQWLGVTLRSRQPG
 30 ENGSIVTCGHRWKNIFYIKNENKLPTGGCYGVPPDLRTELKRIAPCYQDYVKKFGENFASCQAGISS
 FYTKDLIVMGAPGSSYWTGSLFVYNITTNKYKAFLDKQNQVKFGSYLGYSVGAGHFRSQHTTEVVGGA
 PQHEQIGKAYIFSIDEKELNILHEMKGKKLGSYFGASVCAVDLNADGFSDLLVGAPMQSTIREEGRVF
 VYINSGSGAVMNAMETNLVGSCKYARFGESIVNLGIDNDGFEDVAIGAPQEDDLQGAIIYINGRAD
 GISSTFSQRIEGLQISKSLSMFGQSISGQIDADNNGYVDVAVGAFRSDSAVLLRTRPVVIVDASLSHP
 35 ESVNRTKFDCEVNGWPSVCIDLTLCFSYKGKEVPGYIVLFYNMSLDVNRKAESPFRFYSSNGTSDVI
 TGSIQVSSREANCRTHQAFMRKDVRDILTPIQIEAAYHLGPHVISKRSTEEFPPLQPILOQKKEKDIM
 KKTINFARTGGLAQCEKELQALEKENAQLEWELQALEKELAQ

SEQ ID NO:4

Human beta 7 integrin fragment with the BASE-p1 sequence of peptide
 Velcro containing a Cys at the "d" positions of the heptad repeat
 peptide with TEV cleavage site and His6 tag

MVALPMVLVLLLVLSSRGESELDKIPSTGDATEWRNPHLSMLGSCQPAPSCQKCILSHPSCAWCKQLN
 FTASGEAEARRCARREELLARGCPLEELPEPRGQQEVLQDQPLSQGARGEGATQLAPQVRVTLRPGE
 45 PQQQLQVRFLRAEGYPVDLYYLMDSLSSMKDDLERVRQLGHALLVRLQEVTHSVRIGFGSFVDKTVLPF
 VSTVPSKLRHPCPTRLERCQSPFSFHHVLSLTGDAQAFEREVGRQSVSGNLDSPEGGFDAIQAALCQ
 EQIGWRNVSRLLVFTSDDTFHTAGDGKLGFI FMPSDGHCHLDSNGLYSRSTEFDPYPSVGQVAQALSAA
 NIQPIFAVTSAAALPVYQELSKLIPKSAVGELSEDSSNVVQLIMDAYNSLSSTVTLEHSSLPPGVHISY
 ESQCEGPEKREGKAEDRGQCNHVRINQTVTFWVSLQATHCLPEPHLLRLRALGFSEELIVELHTLCLDC
 50 NCSDTQPQAPHCSDGQGHLCGVCSCAPGRLGRLCECSVAELSSPDLESGCGGLENGYFQGGKNAQCK
 KKLQALKKKNAQLKWKQLQALKKKLAQGGHHHHHH

SEQ ID NO:5

Human IL-12p40 wild-type

MCHQQLVISWFSVLFLASPLVAIWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEV
 LGSGKTLTIQVKEFGDAGQYTCHKGGSEVLSHSLLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNY
 SGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVGRDNKEYEYSVEQCQEDSACPAAEE
 SLPIEVMVDAVHKLKYENYTSSFFIRDIIKPDPPKNLQKPLKNSRQVEVSWEYPDTWSTPHSYFSLT
 5 FCVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDRYYSSSWSEWASVPC

SEQ ID NO:6

Human IL-12p40 variant (C199A, C274A)

MCHQQLVISWFSVLFLASPLVAIWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEV
 10 LGSGKTLTIQVKEFGDAGQYTCHKGGSEVLSHSLLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNY
 SGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVGRDNKEYEYSVEQCQEDSACPAAEE
 SLPIEVMVDAVHKLKYENYTSSFFIRDIIKPDPPKNLQKPLKNSRQVEVSWEYPDTWSTPHSYFSLT
 FAVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDRYYSSSWSEWASVPC

15 SEQ ID NO:7

Human IL-12p40 variant (C199A, C274A) with hexahistidine tag

MCHQQLVISWFSVLFLASPLVAIWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEV
 LGSGKTLTIQVKEFGDAGQYTCHKGGSEVLSHSLLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNY
 20 SGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVGRDNKEYEYSVEQCQEDSACPAAEE
 SLPIEVMVDAVHKLKYENYTSSFFIRDIIKPDPPKNLQKPLKNSRQVEVSWEYPDTWSTPHSYFSLT
 FAVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDRYYSSSWSEWASVPCHHHHHH

SEQ ID NO:8

ACID peptide with cysteine residue (3)

25 AQCEKELQALEKENAQLEWELQALEKELAQ

SEQ ID NO:9

BASE peptide with cysteine residue (16), TEV cleavage site (3-9),
 and hexahistidine tag (46-51)

30 GGLENGYFQGGKNAQCKKKLQALKKKNAQLKWKLQALKKKLAQGGHHHHHH

SEQ ID NO:10

TEV cleavage site

EXXYXQ/S

35 X is any amino acid residue

SEQ ID NO:11

Human IL-12p40 variant (C199S, C274S)

MCHQQLVISWFSVLFLASPLVAIWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEV
 40 LGSGKTLTIQVKEFGDAGQYTCHKGGSEVLSHSLLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNY
 SGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVGRDNKEYEYSVEQCQEDSASPAAEE
 SLPIEVMVDAVHKLKYENYTSSFFIRDIIKPDPPKNLQKPLKNSRQVEVSWEYPDTWSTPHSYFSLT
 FSVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDRYYSSSWSEWASVPC

45 SEQ ID NO:12

Human IL-12p40 variant (C199A, C274S)

MCHQQLVISWFSVLFLASPLVAIWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEV
 LGSGKTLTIQVKEFGDAGQYTCHKGGSEVLSHSLLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNY
 50 SGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVGRDNKEYEYSVEQCQEDSACPAAEE
 SLPIEVMVDAVHKLKYENYTSSFFIRDIIKPDPPKNLQKPLKNSRQVEVSWEYPDTWSTPHSYFSLT
 FSVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDRYYSSSWSEWASVPC

SEQ ID NO:13

55 Human IL-12p40 variant (C199S, C274A)

MCHQQLVISWFSVLFLASPLVAIWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLDQSSEV
LGSGKTLTIQVKEFGDAGQYTCHKGGEVLSHSLLLHKKEDGIWSTDILKDQKEPKNKTFLRCEAKNY
SGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVVRGDNKEYEYSVEQCEDSASPAAEE
5 SLPIEVMVDAVHKLKYENYTSSFFIRDIIKPDPPKNLQLKPLKNSRQVEVSWEYPDTWSTPHSYFSLT
FAVQVQGKSKREKKDRVFTDKTSATVICRKNASISVRAQDRYYSSSWSEWASVPC

WHAT IS CLAIMED IS:

1 1. A method for determining the presence or level of a biologic in a
2 sample, the method comprising:

3 (a) contacting the sample with an unlabeled soluble antigen that binds to
4 the biologic to form an unlabeled complex between the antigen and the biologic in the
5 sample;

6 (b) contacting the sample from step (a) with a labeled form of the biologic
7 to form a labeled complex between the antigen and the labeled biologic;

8 (c) subjecting the unlabeled and labeled complexes to size exclusion
9 chromatography to separate the unlabeled and labeled complexes from free labeled biologic
10 and to detect an amount of the free labeled biologic; and

11 (d) comparing the amount of the free labeled biologic detected in step (c)
12 to a standard curve of known amounts of the biologic, thereby determining the presence or
13 level of the biologic in the sample.

1 2. The method of claim 1, wherein the antigen comprises a soluble
2 fragment, variant, or monomer of a membrane-bound protein, a glycosylated protein, a
3 multimeric protein, an insoluble protein, a protein that is difficult to express or purify, and/or
4 a large protein.

1 3. The method of claim 1 or 2, wherein the antigen is a soluble fragment
2 of a cell surface molecule.

1 4. The method of claim 3, wherein the cell surface molecule is a cell
2 adhesion molecule (CAM).

1 5. The method of claim 4, wherein the cell adhesion molecule (CAM) is a
2 member selected from the group consisting of an Ig superfamily CAM, an integrin, a
3 cadherin, and a selectin.

1 6. The method of claim 5, wherein the integrin is an $\alpha 4\beta 7$ integrin.

1 7. The method of claim 6, wherein the soluble fragment comprises an $\alpha 4$
2 fragment comprising an amino acid sequence having at least 80% identity to SEQ ID NO:1 or

3 SEQ ID NO:3 and/or a $\beta 7$ fragment comprising an amino acid sequence having at least 80%
4 identity to SEQ ID NO: 2 or SEQ ID NO:4.

1 8. The method of claim 6 or 7, wherein the biologic is vedolizumab.

1 9. The method of claim 1 or 2, wherein the antigen is a cytokine or a
2 monomer thereof.

1 10. The method of claim 9, wherein the cytokine is a p40 subunit of IL-12
2 or IL-23.

1 11. The method of claim 10, wherein the p40 subunit comprises an amino
2 acid sequence having at least 80% identity to SEQ ID NOS:6, 7, 11, 12, or 13.

1 12. The method of claim 10 or 11, wherein the biologic is ustekinumab.

1 13. The method of any one of claims 1 to 12, wherein the standard curve is
2 generated by incubating the antigen and the labeled biologic with a serial dilution of known
3 amounts of the biologic.

1 14. The method of any one of claims 1 to 13, wherein the area under the
2 curve (AUC) of the free labeled biologic is plotted against the log of known amounts of the
3 biologic and the level of the biologic in the sample is calculated by interpolation.

1 15. The method of any one of claims 1 to 14, wherein the sample is serum.

1 16. The method of any one of claims 1 to 15, wherein the labeled biologic
2 is a fluorophore-labeled biologic.

1 17. The method of any one of claims 1 to 16, wherein the sample is
2 obtained from a subject receiving therapy with the biologic.

1 18. An isolated soluble $\alpha 4$ integrin polypeptide comprising an amino acid
2 sequence having at least 80% identity to SEQ ID NO:1.

1 19. The isolated soluble $\alpha 4$ integrin polypeptide of claim 18, further
2 comprising an ACID peptide.

1 20. The isolated soluble $\alpha 4$ integrin polypeptide of claim 18 or 19, further
2 comprising a linker.

1 21. The isolated soluble $\alpha 4$ integrin polypeptide of any one of claims 18 to
2 20, wherein the amino acid sequence has at least 80% identity to SEQ ID NO:3.

1 22. The isolated soluble $\alpha 4$ integrin polypeptide of any one of claims 18 to
2 21, further comprising an affinity tag.

1 23. An isolated soluble $\beta 7$ integrin polypeptide comprising an amino acid
2 sequence having at least 80% identity to SEQ ID NO:2.

1 24. The isolated soluble $\beta 7$ integrin polypeptide of claim 23, further
2 comprising a BASE peptide.

1 25. The isolated soluble $\beta 7$ integrin polypeptide of claim 23 or 24, further
2 comprising a protease cleavage site.

1 26. The isolated soluble $\beta 7$ integrin polypeptide of any one of claims 23 to
2 25, further comprising an affinity tag.

1 27. The isolated soluble $\beta 7$ integrin polypeptide of any one of claims 23 to
2 26, further comprising a linker.

1 28. The isolated soluble $\beta 7$ integrin polypeptide of any one of claims 23 to
2 27, wherein the amino acid sequence has at least 80% identity to SEQ ID NO:4.

1 29. An isolated soluble $\alpha 4\beta 7$ integrin polypeptide comprising:

2 (a) an $\alpha 4$ integrin polypeptide having an amino acid sequence that has at
3 least 80% identity to SEQ ID NO:1, wherein the $\alpha 4$ integrin polypeptide is linked to a first
4 member of a binding pair, and

5 (b) a $\beta 7$ integrin polypeptide having an amino acid sequence that has at
6 least 80% identity to SEQ ID NO:2, wherein the $\beta 7$ integrin polypeptide is linked to a second
7 member of the binding pair.

1 30. The isolated soluble $\alpha 4\beta 7$ integrin polypeptide of claim 29, wherein
2 the binding pair is selected from the group consisting of coiled-coil peptides, leucine zipper
3 peptides, dock-and-lock peptides, avidin-biotin, and derivatives thereof.

1 31. The isolated soluble $\alpha 4\beta 7$ integrin polypeptide of claim 30, wherein
2 the coiled-coil peptides are ACID-BASE peptides.

1 32. The isolated soluble $\alpha 4\beta 7$ integrin polypeptide of any one of claims 29
2 to 31, wherein the $\alpha 4$ integrin polypeptide and the $\beta 7$ integrin polypeptide heterodimerize via
3 a cysteine bridge or a derivative thereof.

1 33. An isolated soluble IL-12p40 polypeptide comprising an amino acid
2 sequence having at least 80% identity to SEQ ID NOS:6, 11, 12 or 13.

1 34. The isolated soluble IL-12p40 polypeptide of claim 33, wherein the
2 polypeptide further comprises an affinity tag.

1 35. The isolated soluble IL-12p40 polypeptide of claim 33 or 34, wherein
2 the amino acid sequence has at least 80% identity to SEQ ID NO:7.