




- (51) **International Patent Classification:**
G02C 7/02 (2006.01) *G02C 7/04* (2006.01)
- (21) **International Application Number:**
PCT/IB2013/000607
- (22) **International Filing Date:**
22 January 2013 (22.01.2013)
- (25) **Filing Language:** English
- (26) **Publication Language:** English
- (30) **Priority Data:**
61/594,799 3 February 2012 (03.02.2012) US
- (71) **Applicant:** COOPERVISION INTERNATIONAL HOLDING COMPANY, LP [GB/BB]; Suite #2, Edgill House, Wildey Business Park, St. Michael (BB).
- (72) **Inventors:** BACK, Arthur; C/o CooperVision, Inc., 5870 Stoneridge Drive, Suite 1, Pleasanton, CA 94588 (US).
TUAN, Kuang-mon, Ashley; C/o CooperVision, Inc., 5870 Stoneridge Drive, Suite 1, Pleasanton, CA 94588 (US).
- (74) **Agents:** CRITTEN, Matthew, Peter et al.; Abel & Imray, 20 Red Lion Street, London WC1R 4PQ (GB).
- (81) **Designated States** (*unless otherwise indicated, for every kind of national protection available*): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IS, JP, KE, KG, KM, KN, KP, KR, KZ, LA, LC, LK, LR, LS, LT, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.
- (84) **Designated States** (*unless otherwise indicated, for every kind of regional protection available*): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

[Continued on next page]

(54) **Title:** MULTIFOCAL CONTACT LENSES AND RELATED METHODS AND USES TO IMPROVE VISION OF PRESBY-
OPIC SUBJECTS



(57) **Abstract:** Multifocal contact lenses and methods and uses are described. The multifocal contact lenses include an optic zone. The optic zone has an aspheric power profile that provides a near vision refractive power and a distance vision refractive power, and provides an Add power that corresponds to the difference between the near vision refractive power and the distance vision refractive power. The multifocal contact lenses can improve binocular vision of presbyopic subjects by being prescribed such that the non-dominant eye contact lens is over-corrected for distance vision, and both multifocal contact lenses are under-corrected for the Add power requirement of the subject. Batches and sets of multifocal contact lenses are also described.



Published:

— with international search report (Art. 21(3))

— before the expiration of the time limit for amending the claims and to be republished in the event of receipt of amendments (Rule 48.2(h))