

HEARING AID SYSTEM WITH LOST PARTNER FUNCTIONALITY

1. FIELD OF THE INVENTION

The present invention relates to hearing aids. The invention, more particularly relates to hearing aids with a wireless link to transmit and receive data. The invention still
5 more particularly relates to binaural systems, i.e. systems including a pair of hearing aids for aiding both ears. The invention further relates to a method of monitoring the operation of a hearing aid.

A hearing aid is here understood to be an electronic device, for helping the hard of hearing, including an input transducer for receiving an acoustic input signal and
10 convert said input signal into an electrical signal, a signal processing unit for amplifying the signal according to the user's needs and an output transducer for converting the electrical signal into an acoustic output signal. The hearing aid can be any type of hearing aid e.g. In-The-Ear (ITE), Completely-In-the-Canal (CIC), Behind-The-Ear (BTE) or hybrids thereof.

2. THE PRIOR ART

DE 102007039455 A1 shows a hearing aid system, where each hearing aid in a pair includes a magnetic field sensor and a wireless connection between the two hearing aids. If the magnetic field sensor in one hearing aid detects an increased magnetic field without the contralateral hearing aid detecting the same magnetic field, this is
20 interpreted as being due to the magnetic field from a telephone and the hearing aids adjust accordingly. If both hearing aids detect a magnetic field of the same strength, this is interpreted as being due to both hearing aids residing in the hearing aid battery charger, and both hearing aids will automatically deactivate.

Similarly DE 102007046437 A1 describes a hearing aid system with two hearing aids
25 being wirelessly connected to each other. If the two hearing aids are placed in a storage box in close proximity of each other, the field strength of the wireless connection will increase due to the decrease in distance, and the both hearing aids will power down or go in to sleep mode.

US 20050238190 A1 describes a hearing aid with a proximity sensor. Said proximity
30 sensor detects whether a hand is in proximity of the hearing aid, and based on the

signal from said proximity sensor, the hearing aid may change program or turn the volume up or down.

WO 2009006889 A1 describes a hearing aid with means for identifying the transducers in said hearing aid by measuring the impedance of said transducers.

EP 1276349 B1 describes a way for the hearing aid to test parts of the electronics or the full function of the hearing aid, to establish which parts need service.

SUMMARY OF THE INVENTION

Modern hearing aids are so light and comfortable to wear, that the user might not notice wearing them. Consequently, if a hearing aid falls off its position at or in the ear, the user might not notice losing one of the hearing aids. It is a feature of the invention to better monitor the placement of the hearing aids.

Embodiments of the present invention provide a failure detection and notification system for detecting defective or missing parts of a hearing aid system and especially for informing the user of said hearing aid system of such errors.

According to a first aspect, there is provided a hearing aid system including a first hearing aid and an external device, each device having link means for establishing a short-range wireless link between the first hearing aid and the external device, where both the first hearing aid and the external device have means to monitor whether the link to the external device is active, characterised in that both the first hearing aid and the external device have means to issue an alert based on the detection of no traffic at the link within a defined period of time, following a period of time with monitored active traffic on the link to the external device.

With wireless inter-ear connection i.e. connection between the users left and right hearing aid, the wireless connection can be used to control whether the contralateral hearing aid is located at or in the opposite ear and if not, give an alarm to the user that he might have lost his left or right hearing aid, or it is not operating.

The wireless inter-ear connection is normally designed with a short range, e.g. no more than 0,5 m, sufficient for communication while the hearing aids are located at the respective ears. If a hearing aid is moved out of range, the connection is lost.

By periodically communicating with the contralateral hearing aid through the wireless link while constantly monitoring the traffic received, it can be verified whether an active connection to the opposite hearing aid exists. If the connection is lost for a longer period of time, it is likely that one of the following situations has occurred.

Either the battery of the opposite hearing aid is flat and there is no transmission from the contralateral hearing aid, or the opposite hearing aid is located too far away for wireless link to be able to keep an active wireless connection with the contralateral hearing aid. When the system has detected a loss of wireless connection to the second hearing aid over a predetermined period of time, for example 10 seconds, both hearing aids, if still powered, will alert the user to check the status of the opposite hearing aid. If one hearing aid is not powered, the other one will alert the user. The right hearing aid will tell the user to check the status of the left hearing aid and vice versa. As one of the hearing aids is either turned off because of a flat battery or is not located within reach of the wireless link and thereby not located at the correct position at or in the ear of the user, the user will only hear the message in one ear and hereby be alerted that he or she has lost the hearing aid at the opposite ear or the battery of the opposite hearing aid is flat.

In order to ensure that the user will not receive an alarm for no connection between the two hearing aids before the wireless link has been able to initially establish this connection after turning on the hearing aids, the system must preferably log an active connection for some period of time before initiating the monitoring of whether the wireless link is active.

According to a second aspect, there is provided a method for notifying a user of a defect or error in a hearing aid system, comprising establishing a link for short range wireless connection between a first hearing aid and an external device, characterised by comprising monitoring whether the link to the external device is active by starting a first counter once traffic has commenced; if the connection is not active, and the first counter has reached a preset value, initiating a second counter; if the connection is reestablished, resetting the second counter ; if the second counter reaches a predetermined count, issuing an alert in at least the first hearing aid.

BRIEF DESCRIPTION OF THE DRAWINGS

The invention will now be described in greater detail based on preferred embodiments and the drawings, wherein

Fig. 1 shows a binaural hearing aid system connected with a wireless link, according to an embodiment of the invention;

Fig. 2 shows a block diagram of a detection unit according to a first embodiment of the invention; and

Fig. 3 shows a block diagram of a detection unit according to a second embodiment of the invention.

DETAILED DESCRIPTION

After power-up of either one of the hearing aids, it will start looking for its 'partner'
5 on the wireless link, attempting to establish a connection to the opposite hearing aid. Once the connection has been established between the two hearing aids, a first counter will start in order to log the time the connection has been active. If the connection is broken, after it has been confirmed that there has been logged a stable connection for a predetermined length of time, a second counter will initiate. If the wireless
10 connection is reestablished before the second counter reaches a predetermined count, the first counter will reset and start over, and monitor whether there is a stable connection. If the wireless connection is not reestablished before the second counter reaches its predetermined count, this will trigger an alarm. The alarm may be given in several different ways. It could be a spoken message given by the hearing aid to the
15 user, it could be turning on or flashing an LED on the hearing aid, or it could be an alarm sent to a remote control, or another external unit convenient for alerting the user, a relative or an attendant.

Fig. 1 shows two hearing aids 1, 11, each with a microphone 2, 12, a signal processor 3, 13, a receiver 4, 14, a diode 6, 16, an antenna 5, 15 and a wireless link 20
20 connecting the two hearing aids to form a binaural hearing aid system. The control unit 30 may be connected to each of the hearing aids through wireless links 21. In one embodiment the control unit acts as a relay station of the connection between the hearing aids through wireless links 21.

Fig. 2 shows a system adapted to detect whether there is an active connection to the
25 opposite hearing aid, and, if there is no longer an active connection, give an alarm to the user of the hearing aid. A system like this is sometimes referred to as a watchdog. The wireless link activity detector 100 detects whether there is an active wireless connection. Terminal 101 is connected to the wireless link. As long as the wireless link receives data from the opposite hearing aid, the output 104 from the flip-flop or
30 similar device 102 will confirm that the opposite hearing aid is present. Terminal 103 is a boot input to the flip-flop 102, to make sure that the initial state of flip-flop 102 is correct after power-up. Block 110 is an alarm logic. The output 180 from alarm logic

110 will trigger an alarm to the user if the wireless connection has been inactive for a preset amount of time. Box 111 is a dead time counter. Once the wireless link is no longer detecting an active connection to the opposite hearing aid, the dead time counter will start counting the clock cycles received through terminal 112. The circuit
5 120 is a reset circuit to reset the dead time counter 111, the moment the link has been found active again. The comparator 113 will compare the number of clock cycles with a preset value. When the preset value of comparator 113 is reached, it will send a signal to decision block 115, which will change the state of the output from “opposite hearing aid present” to “opposite hearing aid not present”.

10 Fig. 3 shows the wireless link activity detector 100 and the alarm logic 110 like in figure 2. Further figure 3 shows an initialization circuit 200 that detects whether there has been an active wireless connection between the two hearing aids prior to the loss of connection. Circuit 200 comprises a live time counter 201 and comparator circuit 202, with input 101 from the wireless link and the clock terminal 112, or a decimated
15 version of the clock terminal 112 to the counter 201. Similarly the live time counter may be reset through the boot terminal 103. This will be done for instance at start-up to ensure a correct initial state. Block 202 is a comparator, which compares the live time with a preset length of time, to delay the start of the link activity detector 100.

To ensure that an alarm is not triggered due to small interruptions in the connection,
20 the system will include a counter and comparator to determine when the wireless link has been ‘silent’ for too long. The dead time accepted before alarm is set as convenient. Setting it too short may incur a risk of alarm at spurious interruptions of the link, setting it too long implies that a desired alarm may be delayed too much. A suitable length of time may be 10 to 20 seconds. One way to determine the time is to
25 count the clock cycles and compare it with a number of clock cycles preset in the comparator.

If the connection between the two hearing aids is lost and is not reestablished within a period specified, the output of the detection system might enable a light emitting diode (LED) to start flashing on both hearing aids or play an audible message or send
30 an alarm to the remote control, or other external unit. This will help the user locate the hearing aid if lost. If the connection is lost because of a flat battery on one hearing aid, the LED will only flash on the hearing aid still powered by the battery. This will be very helpful in order to determine which hearing aid has a flat battery, especially for

people who need help from relatives or care takers to change their battery, as the hearing aid with a flat battery will not be flashing the LED. While the spoken message may stop after three messages the LED will keep flashing until the connection is reestablished or the hearing aids are disabled or the battery is flat. Disabling the LED
5 and/or the audio alarm might be done by a constant press on the program switch for 5 seconds or similar input to the hearing aid or the remote control.

As some users may occasionally only wear a hearing aid in one ear, or may not want the alert feature, it is important to be able to turn the system on and off. This is preferably done with the appropriate fitting software at the initial fitting by the
10 audiologist. The fitting software includes the option to enable or disable this feature of the hearing aid. Further, it should be possible to set the number of times a message should be repeated and the interval by which the message should be repeated. The frequency of repeating the message may not be constant. One desirable way to repeat the message could be to repeat the message at a 2 minute interval for the first three
15 times. Repeating the message from 4 to 10 times at an interval of 10 minutes, after that the message will be repeated at an interval of half an hour. Even though the audiologist has adjusted these settings according to the users wish, the user might still want to be able to control these settings him or herself. This can be done by giving the hearing aid user a 'light' version of the fitting software, where he or she is able to
20 adjust these settings at home at their own PC.

A problem sometimes arising with hearing aids during use is that a transducer stops working, partly or totally. This may occur because the transducer has been exposed to water or moist. While the transducers may be replaced, once the failure has been detected, it is not always easy for the user to identify the problem as it may arise
25 slowly from just a slight degree of malfunctioning, barely audible to the user, to a severe degree of malfunctioning. If this happens over several weeks or months, the user might just think that it is due to a natural increase in his or her hearing loss. However it may be possible to measure this error in the transducer electrically. WO 2009006889 A1 teaches a way to identify whether it is the correct receiver attached to
30 the hearing aid in a receiver in the ear (RITE) type of hearing aid. Similarly, the impedance of a broken transducer may change and indicate that it is no longer working properly. If the receiver of the first hearing aid is broken or malfunctioning, a message can be sent via the wireless connection to the second hearing aid, indicating a

problem with the receiver and informing the user to seek professional assistance to solve the problem. Further EP 1276349 B1 teaches a way for the hearing aid itself to test parts of the hearing aid or the full function of the hearing aid. Such a test may be triggered in one hearing aid and if it is not possible the alarm to user of a defect by the
5 respective hearing aid, the alarm may be given in the opposite hearing aid.

AMENDED CLAIMS:

1. A hearing aid system including a first hearing aid and an external device, each device having link means for establishing a short-range wireless link between the first hearing aid and the external device, where both the first hearing aid and the external device have means to monitor whether the link to the external device is active, C H A R A C T E R I S E D I N T H A T both the first hearing aid and the external device have means to issue an alert based on the detection of no traffic at the link within a defined period of time, following a period of time with monitored active traffic on the link to the external device.
2. The system according to claim 1, wherein the means to monitor whether the link to the external device is active comprises, a first timer to count the clock cycles elapsed since the onset of traffic, the first timer having a first comparator adapted to delay the start of active link detector until a first predetermined time interval is reached, a watchdog timer adapted to count the clock cycles elapsed since the tripping by the first comparator, the watchdog timer having a second comparator to determine if a second predetermined interval is reached, and alarm logic for issuing an alert on the conditions of reaching the predetermined time interval by the second comparator.
3. The system according to claim 1, wherein the external device has means for conducting a self-test and means for sending a message on the link to the first hearing aid about the result of the self-test.
4. The system according to claims 1 or 2, wherein the alert means is adapted to issue an audible message or tone.
5. The system according to claims 1 or 2, wherein the alert means comprises a visual indicator.
6. The system according to claims 1 to 5, wherein the external device is a second hearing aid.
7. The system according to claims 1 to 5, wherein the external device is a remote control.
8. A method for notifying a user of a defect or error in a hearing aid system, comprising establishing a link for short range wireless connection between a first hearing aid

and an external device, C H A R A C T E R I S E D by comprising monitoring whether the link to the external device is active by starting a first counter once traffic has commenced;

if the connection is not active, and the first counter has reached a preset value, initiating a second counter;

if the connection is reestablished, resetting the second counter;

if the second counter reaches a predetermined count, issuing an alert in at least the first hearing aid.

9. The method according to claim 8, wherein the external device is a second hearing aid.

10. The method according to claim 8, wherein the alert is a spoken message, a visual indicator, or a message sent to an external unit.

11. The method according to claim 8, wherein the external device is a remote control.