

Fig 1

```

Met Ala Gln Val Leu Arg Gly Thr Val Thr Asp Phe Pro Gly Phe Asp
1           5           10           15

Glu Arg Ala Asp Ala Glu Thr Leu Arg Lys Ala Met Lys Gly Leu Gly
          20           25           30

Thr Asp Glu Glu Ser Ile Leu Thr Leu Leu Thr Ser Arg Ser Asn Ala
          35           40           45

Gln Arg Gln Glu Ile Ser Ala Ala Phe Lys Thr Leu Phe Gly Arg Asp
          50           55           60

Leu Leu Asp Asp Leu Lys Ser Glu Leu Thr Gly Lys Phe Glu Lys Leu
65           70           75           80

Ile Val Ala Leu Met Lys Pro Ser Arg Leu Tyr Asp Ala Tyr Glu Leu
          85           90           95

Lys His Ala Leu Lys Gly Ala Gly Thr Asn Glu Lys Val Leu Thr Glu
          100          105          110

Ile Ile Ala Ser Arg Thr Pro Glu Glu Leu Arg Ala Ile Lys Gln Val
          115          120          125

Tyr Glu Glu Glu Tyr Gly Ser Ser Leu Glu Asp Asp Val Val Gly Asp
          130          135          140

Thr Ser Gly Tyr Tyr Gln Arg Met Leu Val Val Leu Leu Gln Ala Asn
145          150          155          160

Arg Asp Pro Asp Ala Gly Ile Asp Glu Ala Gln Val Glu Gln Asp Ala
          165          170          175

Gln Ala Leu Phe Gln Ala Gly Glu Leu Lys Trp Gly Thr Asp Glu Glu
          180          185          190

Lys Phe Ile Thr Ile Phe Gly Thr Arg Ser Val Ser His Leu Arg Lys
          195          200          205

Val Phe Asp Lys Tyr Met Thr Ile Ser Gly Phe Gln Ile Glu Glu Thr
          210          215          220

```

Ile Asp Arg Glu Thr Ser Gly Asn Leu Glu Gln Leu Leu Leu Ala Val
225 230 235 240

Val Lys Ser Ile Arg Ser Ile Pro Ala Tyr Leu Ala Glu Thr Leu Tyr
245 250 255

Tyr Ala Met Lys Gly Ala Gly Thr Asp Asp His Thr Leu Ile Arg Val
260 265 270

Met Val Ser Arg Ser Glu Ile Asp Leu Phe Asn Ile Arg Lys Glu Phe
275 280 285

Arg Lys Asn Phe Ala Thr Ser Leu Tyr Ser Met Ile Lys Gly Asp Thr
290 295 300

Ser Gly Asp Tyr Lys Lys Ala Leu Leu Leu Leu Cys Gly Glu Asp Asp
305 310 315 320

Fig 2




Fig 3

Process Step	Parameter	3 L process	100 L process
AX ▼	Resin:	Q Sepharose XL	Q Sepharose XL
	Column:	XK 50	BPG 300
	Column diameter (cm):	5	30
	Bed height (cm):	13.5	13.5
	Column volume (L):	0.265	9.5
	Flow rate:	200 cm/h	170 cm/h
	Buffer temperature:	Ambient temperature	Ambient temperature
EQ ▼	Pump:	with WATSON-MARLOW Tubing pump 323S	BioProcess system
	Buffer A:	20 mM Tris pH 7.4, 0.1 % Tween80, 25 mM NaCl	20 mM Tris pH 7.4, 0.1 % Tween80, 25 mM NaCl
Load ▼	Vol:	2 CV	3 CV
			Online dilution (1+1) with 1% (v/v) Tween 80, 4 mM EDTA, pH 8.0
Wash ▼	Buffer:	20 mM Tris pH 7.4, 0.1 % Tween80, 25 mM NaCl	20 mM Tris pH 7.4, 0.1 % Tween80, 25 mM NaCl
	Vol:	10 CV	10 CV
Step Elution ▼	100 % Buffer B:	20 mM Tris pH 7.4, 0.1 % Tween80, 300 mM NaCl	20 mM Tris pH 7.4, 0.1 % Tween80, 300 mM NaCl
	Vol:	9 CV	9 CV
Pooling ▼	Criteria:	Main peak Pool (0.1 AU-0.2 AU)	Main peak Pool (target 0.2 AU-0.4 AU ; 5 mm light path)
			(within 3 weeks)
CIP ▼	Step 1:	2 M NaCl	2 M NaCl
	Vol:	3 CV	3 CV
	Step 2:	1 M NaOH	1 M NaOH
Storage ▼	Vol:	3 CV	3 CV
	Incubation:	15 h	15 h
Storage ▼	Storage:	20 mM NaOH	10 mM NaOH

Fig 4

Process Step	Parameter	3 L Process	100 L Process
<div>AF</div> <div>▼</div>	Resin:	Heparin HyperD	Heparin HyperD
	Column:	XK 50	BPG300
	Column diameter (cm):	5	30
	Bed height (cm):	25	25
	Column volume (L):	0.491	17.7
	Flow rate:	60-100 cm/h	60-100 cm/h
	Buffer temperature:	Ambient temperature	Ambient temperature
<div>EQ</div> <div>▼</div>	Buffer A:	20 mM Tris pH 7.4; 0.1 % Tween80; 2 mM CaCl ₂ , 25 mM NaCl	20 mM Tris pH 7.4; 0.1 % Tween80; 2 mM CaCl ₂ , 25 mM NaCl
	Vol:	2 CV	3 CV
<div>Load</div> <div>▼</div>	Flow rate:	100 cm/h	60-100 cm/h; online 1:8 dilution of AX eluate with buffer A
<div>Wash 1</div> <div>▼</div>	Buffer:	20 mM Tris pH 7.4; 0.1 % Tween80; 2 mM CaCl ₂ , 25 mM NaCl	20 mM Tris pH 7.4; 0.1 % Tween80; 2 mM CaCl ₂ , 25 mM NaCl
	Vol:	15 CV	15 CV
	Flow rate:	100 cm/h	60-100 cm/h
<div>Wash 2</div> <div>▼</div>	Buffer:	20 mM Tris pH 7.4; 0.1 % Tween80, 25 mM NaCl	20 mM Tris pH 7.4; 0.1 % Tween80, 25 mM NaCl
	Vol:	2 CV	2 CV
	Flow rate:	100 cm/h	60-100 cm/h
<div>Step Elution</div> <div>▼</div>	100 % Buffer B:	20 mM Tris pH 7.4; 0.1 % Tween80; 10 mM EDTA; 25 mM NaCl	20 mM Tris pH 7.4; 0.1 % Tween80; 10 mM EDTA; 25 mM NaCl
	Vol:	2-3 CV	2-3 CV
	Flow rate:	40-60 cm/h	40-60 cm/h
<div>Pooling</div> <div>▼</div>	Criteria:	Complete peak Pool (0.1 AU-0.1 AU)	Complete peak Pool (0.2 AU-0.2 AU ; 5 mm light path)
			(within 3 weeks)
<div>CIP</div> <div>▼</div>	Step 1:	2 M NaCl	2 M NaCl
	Vol:	3 CV	3 CV
	Step 2:	0.1 M NaOH	0.1 M NaOH
<div>Storage</div> <div>▼</div>	Vol:	3 CV	3 CV
	Incubation:	15 h	≥15 h (contact time, no pumping)
	Storage:	1 M NaCl, 25% EtOH	1 M NaCl, 25% EtOH

Fig 5

Process Step	Parameter	3 L process	100 L process
AX ▼	Resin:	Source15 Q	Source15 Q
	Column:	-	Fineline200
	Column diameter (cm):	-	20
	Bed height (cm):	15	15
	Column volume (L):	0.4	4.7
	Flow rate:	100 cm/h	100 cm/h
	Buffer temperature:	Ambient temperature	Ambient temperature
EQ ▼	Buffer A:	20 mM Bis-Tris pH 7; 25 mM NaCl	20 mM Bis-Tris pH 7; 25 mM NaCl
	Vol:	2 CV	3 CV
Load ▼	Flow rate:	100 cm/h	100 cm/h AF pool splitting (1/4-1/5)
Wash ▼	Buffer:	20 mM Bis-Tris pH 7; 25 mM NaCl	20 mM Bis-Tris pH 7; 25 mM NaCl
	Vol:	3 CV	3 CV
	Flow rate:	100 cm/h	100 cm/h
Gradient Elution ▼	Buffer B:	20 mM Bis-Tris pH 7; 180 mM NaCl 0-100% B in 33 CV	20 mM Bis-Tris pH 7; 180 mM NaCl 0-100% B in 33 CV
	Gradient:		
Pooling ▼	Criteria:	major peak starting at 0.05 AU to the valley between peak 1 and peak 2	major peak starting at 0.1 AU to the valley between peak 1 and peak 2 (5 mm light path)
CIP ▼	Step 1:	2 M NaCl	(within 3 weeks) 2 M NaCl
	Vol:	3 CV	3 CV
	Step 2:	1 M NaOH	1 M NaOH
	Vol:	3 CV	3 CV
Storage ▼	Incubation:	15 h	≥15 h (contact time, no pumping)
	Storage:	25 mM NaOH	10 mM NaOH

Fig 6

Process Step	Parameter	3 L process	100 L process
Concentration	cassette: area: Target concentration: Processing time: Parameter: Concentration factor	Hydrosart 10 K 0.1 m ² 5 g/L ~ 2h TMP 0.9-1.1 Entrance pressure 1 bar Retentate back pressure 0.8 bar ~ 6	Hydrosart 10 K 2.4 m ² 5 g/L - TMP 0.9-1.1 Entrance pressure 1 bar Retentate back pressure 0.8 bar 6-8
▼			
Buffer change	► Dilution buffer: Diafiltration volumes: Concentration factor Target concentration: Dilution	20 mM Bis-Tris, 150 mM NaCl, 1mM CaCl ₂ pH 7 10 ~ 6 12 g/L Addition of cassette wash	20 mM Bis-Tris, 150 mM NaCl, 1mM CaCl ₂ pH 7 10 ~ 6 12 g/L Addition of cassette wash
▼			
Addition of Tween80	► Buffer: Target concentration Tween80: Target concentration bulk:	20 mM Bis-Tris, 150 mM NaCl, 1mM CaCl ₂ , 10 % Tween80 pH 7 0.05 % 10 g/L	20 mM Bis-Tris, 150 mM NaCl, 1mM CaCl ₂ , 10 % Tween80 pH 7 0.05 % 10 g/L
▼			
Sterile filtration	► Device: Cut off: Target concentration bulk:	Sartopore2 0.45-0.2 µm 10 g/L	Sartopore2 0.45-0.2 µm 10 g/L
▼			
Storage of bulk	► Temperature:	2-8°C	2-8°C




FIG. 7A




FIG. 7B