


FIG. 1







Fig. 2

Fr V Susp + 20 mM Caprylate				
Protein	Temp	Symbol	LRV	
			pH 4.5	pH 5.1
25 AU	5C	○	5.7	2.8
	12C	□	5.7	5.7
	20C	△	5.7	5.7
65 AU	5C	●	2.6	0.5
	12C	■	4.1	0.6
	20C	▲	5.7	2.6


Albumin Susp + 20 mM Caprylate				
Protein	Temp	Symbol	LRV	
			pH 4.5	pH 5.1
25 AU	5C	○	5.7	3.7
	12C	□	5.7	5.7
	20C	△	5.7	5.7
65 AU	5C	●	0.3	1.6
	12C	■	0.5	2.6
	20C	▲	1.2	5.7

Fig. 3

A



B



Fig, 4





FIG. 5

A



B





FIG. 6

A



B




FIG. 7

