



G00002

14

ABSTRACT

PROCESS AND DEVICE FOR TRANSFER OF STERILE PRODUCTS
BETWEEN A CONTAINER AND AN ISOLATOR

Sterile products contained in a container (10) are transferred into a sterile isolator (12) by connecting the container onto a lock (28) fitted on the outside of the isolator (12) and equipped with pulsed light sterilization means (28). The connection between the lock (28) and the isolator (12) is made by a sealed double door transfer device. The throwaway container is equipped with a plug type closer (16), cutout from the inside of the isolator (12) after sterilizing the volume (40) inside the lock (28) and opening the double door (24, 38).

Fig. 1



162162