



(51) International Patent Classification:

A01N 37/16 (2006.01) A01N 25/02 (2006.01)
A01N 37/02 (2006.01) A01N 25/00 (2006.01)
A01N 25/30 (2006.01) A01N 59/00 (2006.01)

(21) International Application Number:

PCT/AU2020/050101

(22) International Filing Date:

07 February 2020 (07.02.2020)

(25) Filing Language:

English

(26) Publication Language:

English

(30) Priority Data:

2019900479 15 February 2019 (15.02.2019) AU

(71) Applicant: WHITELEY CORPORATION PTY LTD

[AU/AU]; Suite 501, Level 5, 12 Mount Street, North Sydney, New South Wales 2060 (AU).

(72) Inventors: GLASBEY, Trevor; 26 Pershing Place, Taniba

Bay, New South Wales 2319 (AU). ROBERTS, Nicholas Alan; 30 Endeavour Close, Woodrising, New South Wales

2284 (AU). WHITELEY, Gregory Stuart; 155 Queenscliff Street, Queenscliff, New South Wales 2096 (AU).

(81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ, CA, CH, CL, CN, CO, CR, CU, CZ, DE, DJ, DK, DM, DO, DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN, HR, HU, ID, IL, IN, IR, IS, JO, JP, KE, KG, KH, KN, KP, KR, KW, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME, MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ, OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA, SC, SD, SE, SG, SK, SL, ST, SV, SY, TH, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, WS, ZA, ZM, ZW.

(84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).

(54) Title: IMPROVED ENDOSCOPE DISINFECTANT

Comparison of Inventive formulations with prior art examples




FIGURE 3

(57) Abstract: A working disinfectant solution comprising peracetic acid and at least one surfactant, prepared by the dilution of either a single-part or two-part concentrate for use in the disinfection or sterilisation of reusable medical devices such as endoscopes. The working solution exhibits rapid wetting characterised by a dynamic surface tension of less than 50mN/m at 250ms and 46mN/m at 500ms.

[Continued on next page]

Published:

— with international search report (Art. 21(3))