



Figure 1: Overlay of H1H17203P, H1H17139P, and H1H17161P RP-UPLC Chromatograms



RP-UPLC chromatograms of 50 mg/mL H1H17203P (dashes — —), 50 mg/mL H1H17139P (solid —), and 50 mg/mL H1H17161P (circle dots).

Figure 2: Overlay of H1H17203P, H1H17139P and H1H17161P SE-UPLC Chromatograms



SE-UPLC chromatograms of H1H17203P (dashes—), H1H17139P (solid —), H1H17161P (circle dots.....).