






Fig. 2. Concentration test
Spores of *Bacillus subtilis* ATCC 6633 (from H₂O) in ethanol




log **Fig. 3. Concentration test**
Spores of *Clostridium sporogenes* ATCC 19404 (from H₂O) in ethanol



100


Fig. 4. Concentration test

Spores of *Bacillus subtilis* ATCC 6633 (from EtOH) in isopropanol




log

Fig. 5. Concentration test

Spores of *Bacillus subtilis* ATCC 6633 (from H₂O) in isopropanol


log

Fig. 6. Concentration test
Spores of *Clostridium sporogenes* ATCC 19404 (from H₂O) in isopropanol




log

Fig. 7. Concentration test
Spores of *Bacillus subtilis* ATCC 6633 (from EtOH) in n-propanol



log

Fig. 8. Concentration test
Spores of *Bacillus subtilis* ATCC 6633 (from H₂O) in n-propanol



log

Fig. 9. Concentration test
Spores of *Clostridium sporogenes* ATCC 19404 (from H₂O) in n-propanol




Fig. 10. Spore count in composition X




Fig. 11. Spore count in composition Y

