

















[Fig. 1]


[Fig. 2]

Buffer Composition	Without SA Solution	With SA Solution (2.5 mM KCl)	With SA Solution (20 mM KCl)
50 mM Citric acid (pH 2.5)			
2.5 mM Citric acid (pH 2.9)			
0.7 mM Citric acid (pH 3.3)			
0.7 mM Citric acid (pH 4.5)	Precipitated	Precipitated	Precipitated
0.7 mM Citric acid (pH 5.5)			


[Fig. 3]




[Fig. 4]




[Fig. 5]




[Fig. 6]






[Fig. 7]


[Fig. 8]




[Fig. 9]



[Fig. 10]



[Fig. 11]



[Fig. 12]

