

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property
Organization
International Bureau

(43) International Publication Date
22 May 2020 (22.05.2020)



(10) International Publication Number
WO 2020/101680 A1

(51) International Patent Classification:

A61B 8/06 (2006.01) *A61B 17/11* (2006.01)
A61B 8/08 (2006.01) *A61B 8/12* (2006.01)
A61B 8/00 (2006.01) *A61B 5/00* (2006.01)

(21) International Application Number:

PCT/US2018/061191

(22) International Filing Date:

15 November 2018 (15.11.2018)

(25) Filing Language:

English

(26) Publication Language:

English

(71) Applicants: **BAXTER INTERNATIONAL INC.**
[US/US]; One Baxter Parkway, Deerfield, Illinois 60015
(US). **BAXTER HEALTHCARE SA** [CH/CH]; Thur-
gauerstrasse 130, CH-8152 Glattpark (Opfikon) (CH).

(72) Inventors: **WITKOWSKI, Shannon**; 2035 Berkeley Av-
enue, St. Paul, Minnesota 55105 (US). **SCHEIDNES,**
Michael; 8193 Shadyview Lane N, Maple Grove, Minneso-
ta 55311 (US). **KWON, Sung**; 2410 University Ave. W
#512, St. Paul, Minnesota 55114 (US). **STUDER, James**;
W7296 770th Ave., Beldenville, Wisconsin 54003 (US).

(74) Agent: **WEBER, Brett J.**; K&L Gates LLP, c/o Foreign
Patents, P.O. Box 1135, Chicago, Illinois 60690-1135 (US).

(81) Designated States (*unless otherwise indicated, for every
kind of national protection available*): AE, AG, AL, AM,
AO, AT, AU, AZ, BA, BB, BG, BH, BN, BR, BW, BY, BZ,
CA, CH, CL, CN, CO, CR, CU, CZ, DE, DJ, DK, DM, DO,
DZ, EC, EE, EG, ES, FI, GB, GD, GE, GH, GM, GT, HN,
HR, HU, ID, IL, IN, IR, IS, JO, JP, KE, KG, KH, KN, KP,
KR, KW, KZ, LA, LC, LK, LR, LS, LU, LY, MA, MD, ME,
MG, MK, MN, MW, MX, MY, MZ, NA, NG, NI, NO, NZ,
OM, PA, PE, PG, PH, PL, PT, QA, RO, RS, RU, RW, SA,

(54) Title: VASCULAR MONITORING SYSTEM




FIG. 1B

(57) Abstract: A Doppler blood flow monitoring device (150) includes a signal generation module, a signal reception module, a signal filtration module, a signal conversion module, at least one speaker, and a user interface (304). The signal generation module is configured to send a signal to a probe (106) positioned in a probe receptacle on a vascular coupler positioned about a patient's vessel. The signal reception module is configured to receive a return signal from the probe. The signal filtration module is configured to filter the return signal. The signal conversion module is configured to convert the filtered signal into an audible indication and a visual indication corresponding to a characteristic of blood flow in the patient's vessel. The at least one speaker (214) is configured to emit the first audible indication. Additionally, the user interface is configured to display the visual indication.

[Continued on next page]

SC, SD, SE, SG, SK, SL, SM, ST, SV, SY, TH, TJ, TM, TN,
TR, TT, TZ, UA, UG, US, UZ, VC, VN, ZA, ZM, ZW.

(84) Designated States (*unless otherwise indicated, for every kind of regional protection available*): ARIPO (BW, GH, GM, KE, LR, LS, MW, MZ, NA, RW, SD, SL, ST, SZ, TZ, UG, ZM, ZW), Eurasian (AM, AZ, BY, KG, KZ, RU, TJ, TM), European (AL, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR, GB, GR, HR, HU, IE, IS, IT, LT, LU, LV, MC, MK, MT, NL, NO, PL, PT, RO, RS, SE, SI, SK, SM, TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, KM, ML, MR, NE, SN, TD, TG).

Published:

— with international search report (Art. 21(3))