

SEQUENCE LISTING

<110> National University Corporation Okayama University

<120> A method for differentiating organ composed mature cells from immature cells, reconstituting 3-dimensional organ structure, and medical reconstituting and regenerating organs for therapy

<130> OGTP1002F

<150> JP 2009-188578

<151> 2009-8-17

<160> 22

<170> PatentIn version 3.1

<210> 1

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized forward primer for detecting auaporin-1 gene

<400> 1

cctccaggca cagtcttctc 20

<210> 2

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized reverse primer for detecting auaporin-1 gene

<400> 2

cagtggcctc ctgactcttc 20

<210> 3

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized forward primer for detecting Na-Ca transporter gene

<400> 3

tggaagctgg tctgtctcct 20

<210> 4

<211> 25

<212> DNA

<213> Artificial

<220>

<223> synthesized reverse primer for detecting Na-Ca transporter gene

<400> 4

tgatgacatc aagaaggtgg tgaag 25

<210> 5

<211> 17

<212> DNA

<213> Artificial

<220>

<223> synthesized forward primer for detecting nephrin gene

<400> 5

gctccacacca tccgtgc 17

<210> 6

<211> 22

<212> DNA

<213> Artificial

<220>

<223> synthesized reverse primer for detecting nephrin gene

<400> 6

gactatgtcc acacaacccc ca 22

<210> 7

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized forward primer for detecting pax-2 gene

<400> 7

agggcacatctg cgataatgac 20

<210> 8

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized reverse primer for detecting pax-2 gene

<400> 8

ctcgcgtttc ctcttctcac 20

<210> 9

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized forward primer for detecting gdnf gene

<400> 9

cccgaagatt atcctgacca 20

<210> 10

<211> 20

<212> DNA

<213> Artificial

<220>

<223> synthesized reverse primer for detecting gdnf gene

<400> 10

tagcccaaac ccaagtcagt 20

<210> 11
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Synthesized forward primer for detecting ret gene

<400> 11
gcgtcagggga gatggtaaag 20

<210> 12
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Synthesized reverse primer for detecting ret gene

<400> 12
catcagggaa acagttgcag 20

<210> 13
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Synthesized forward primer for detecting ncam gene

<400> 13
acgtccgggt catagtcctg 20

<210> 14
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Synthesized reverse primer for detecting ncam gene

<400> 14
ctatgggttc cccatccttt 20

<210> 15
<211> 20
<212> DNA
<213> Artificial

<220>
<223> Synthesized forward primer for detecting wt-1 gene

<400> 15
acccaggctg caataagaga 20

<210> 16
<211> 20
<212> DNA
<213> Artificial

<220>

<223> synthesized reverse primer for detecting wt-1 gene

<400> 16
gctgaagggc ttttcacttg 20

<210> 17
<211> 20
<212> DNA
<213> Artificial

<220>
<223> synthesized forward primer for detecting AQP-2 gene

<400> 17
gtggctgccc agctgctggg 20

<210> 18
<211> 20
<212> DNA
<213> Artificial

<220>
<223> synthesized reverse primer for detecting AQP-2 gene

<400> 18
agctccaccg actgccgccg 20

<210> 19
<211> 20
<212> DNA
<213> Artificial

<220>
<223> synthesized forward primer for detecting bf-2 gene

<400> 19
aggagacaga catcgacgtg 20

<210> 20
<211> 19
<212> DNA
<213> Artificial

<220>
<223> synthesized reverse primer for detecting bf-2 gene

<400> 20
tgacgaagca gtcgttgag 19

<210> 21
<211> 25
<212> DNA
<213> Artificial

<220>
<223> synthesized forward primer for detecting gapdh gene

<400> 21
tgatgacatc aagaaggtgg tgaag 25

<210> 22
<211> 23

<212> DNA
<213> Artificial

<220>
<223> synthesized reverse primer for detecting gapdh gene

<400> 22
tccttggagg ccatgtaggc cat 23