




Fig. 1

Addition of lig sequences to the ends of single stranded cDNA templates.





Fig. 2

Overview of Level-Up normalization process.



Front and Back K-linker structures as shown annealed to single strand cDNA template.

Figure 3




Figure 3A




Fig 3B




Fig. 4

Length distribution from gel electrophoresis of input cDNA.




Fig. 5

Length distribution from gel electrophoresis of normalized cDNA.




Fig. 6

Saturation curves for input cDNA and Level-Up normalized cDNA using Nanopore cDNA sequencing.